

M E T E R S

Product Series	SKU	Model	Product Description	public price
ZenTest Smart Testers	AI3711	PH60-Z	ZenTest Smart Pocket pH Tester Kit	\$ 181.29
ZenTest Smart Testers	AI3715	PH60F-Z	ZenTest Smart Pocket Flat pH Tester Kit	\$ 207.54
ZenTest Smart Testers	AI3713	PH60S-Z	ZenTest Smart Pocket Spear pH Tester Kit	\$ 343.82
ZenTest Smart Testers	AI3717	ORP60-Z	ZenTest Smart Pocket ORP TesterKit	\$ 222.55
ZenTest Smart Testers	AI3718	EC60-Z	ZenTest Smart Pocket Conductivity/TDS/Salinity/Resistivity Tester Kit	\$ 181.29
ZenTest Smart Testers	AI3719	PC60-Z	ZenTest Smart Pocket pH/Conductivity/TDS/Salinity/Resistivity Tester Kit	\$ 243.80
ZenTest Smart Testers	AI3803	PCO60-Z	ZenTest Smart Pocket pH/Conductivity/TDS/Salinity/Resistivity/ORP Tester Kit	\$ 325.07
Zentest Smart LabSen Testers	AI3737	PH60Z-WW	PH60Z-WW Smart pH Meter Kit for Wastewater and Solids-Containing Solutions	\$ 266.83
Zentest Smart LabSen Testers	AI3739	PH60Z-PW	PH60Z-PW Smart pH Meter Kit for Pure Water	\$ 404.98
Zentest Smart LabSen Testers	AI3741	PH60Z-MS	PH60Z-MS Smart pH Meter Kit for Small-Volume Samples	\$ 377.89
Zentest Smart LabSen Testers	AI3743	PH60Z-SA	PH60Z-SA Smart pH Meter Kit for Strong Alkaline and High Salinity Solutions	\$ 350.80
Zentest Smart LabSen Testers	AI3745	PH60Z-HF	PH60Z-HF Smart pH Meter Kit for HF-Containing and Strong Acidic Solutions	\$ 350.80
Zentest Smart LabSen Testers	AI3747	PH60Z-HT	PH60Z-HT Smart pH Meter Kit for High Temperature and Caustic Solutions	\$ 377.89
Zentest Smart LabSen Testers	AI3749	PH60Z-VS	PH60Z-VS Smart pH Meter Kit for Viscous Samples	\$ 364.35
Zentest Smart LabSen Testers	AI3751	PH60Z-MT	PH60Z-MT Smart pH Meter Kit for Meats	\$ 494.48
60 Premium Series Testers	AI311	PH60	Premium Series Pocket pH Tester Kit	\$ 131.28
61 Premium Series Testers	AI312	PH60F	Premium Series Pocket Flat pH Tester Kit	\$ 161.28
62 Premium Series Testers	AI313	PH60S	Premium Series Pocket Spear pH Tester Kit	\$ 310.06
63 Premium Series Testers	AI315	ORP60	Premium Series Pocket ORP Tester Kit	\$ 162.53
64 Premium Series Testers	AI314	EC60	Premium Series Pocket Conductivity/TDS/Salinity Tester Kit	\$ 133.78
65 Premium Series Testers	AI316	PC60	Premium Series Pocket pH/Conductivity/TDS/Salinity Tester Kit	\$ 218.79
66 Premium Series Testers	AI3801	PCO60	Premium Series Pocket pH/Conductivity/TDS/Salinity/ORP Tester Kit	\$ 287.56
60 Premium Series LabSen Testers	AI3721	PH60-WW	PH60-WW pH Meter Kit for Wastewater and Solids-Containing Solutions	\$ 212.65
61 Premium Series LabSen Testers	AI3723	PH60-PW	PH60-PW pH Meter Kit for Pure Water	\$ 350.80
62 Premium Series LabSen Testers	AI3725	PH60-MS	PH60-MS pH Meter Kit for Small-Volume Samples	\$ 323.71



63 Premium Series LabSen Testers	AI3727	PH60-SA	PH60-SA pH Meter Kit for Strong Alkaline and High Salinity Solutions	\$ 296.62
64 Premium Series LabSen Testers	AI3729	PH60-HF	PH60-HF pH Meter Kit for HF-Containing and Strong Acidic Solutions	\$ 296.62
65 Premium Series LabSen Testers	AI3731	PH60-HT	PH60-HT pH Meter Kit for High Temperature and Caustic Solutions	\$ 323.71
66 Premium Series LabSen Testers	AI3733	PH60-VS	PH60-VS pH Meter Kit for Viscous Samples	\$ 310.17
67 Premium Series LabSen Testers	AI3735	PH60-MT	PH60-MT pH Meter Kit for Meats	\$ 404.98
20 Value Series Testers	AI209	PH20	Value Series Pocket pH Tester Kit	\$ 82.52
21 Value Series Testers	AI212	EC20	Value Series Pocket Conductivity Tester Kit	\$ 87.52
22 Value Series Testers	AI213	TDS20	Value Series Pocket TDS Tester Kit	\$ 87.52
23 Value Series Testers	AI214	Salt20	Value Series Pocket Salinity Tester Kit	\$ 87.52
GroStar Series Pen Testers for Growers	AI101G	GS1	GroStar pH Pen Tester Kit	\$ 145.03
GroStar Series Pen Testers for Growers	AI102G	GS2	GroStar Soil pH Pen Tester Kit	\$ 321.32
GroStar Series Pen Testers for Growers	AI103G	GS3	GroStar EC/500ppm/700ppm Pen Tester Kit	\$ 142.53
GroStar Series Pen Testers for Growers	AI104G	GS4	GroStar pH/EC/500ppm/700ppm Combo Pen Tester Kit	\$ 223.80
GroStar Series Pen Testers for Growers	AI201G	GS1-P	GroStar pH/ORP Pen Tester Kit with ORP Probe	\$ 220.05
GroStar Series Pen Testers for Growers	AI204G	GS4-P	GroStar pH/EC/500ppm/700ppm/ORP Combo Pen Tester Kit with ORP probe	\$ 298.81
Eco Series Compact Pens	AI231	5021	Compact Conductivity Pen	\$ 62.51
Eco Series Compact Pens	AI232	5031	Compact TDS Pen	\$ 62.51
Eco Series Compact Pens	AI239	5041	Compact ORP Pen	\$ 87.52
Eco Series Compact Pens	AI236	5051	Compact Salinity Pen	\$ 62.51
Eco Series Compact Pens	AI235	5052	Compact Saltwater Salinity Pen	\$ 106.27
Eco Series Compact Pens	AI237	5053	Compact NaCl Salinity Pen	\$ 106.27
pH Pen Testers	AI221	SX610	pH Pen Tester Kit (± 0.1 pH)	\$ 93.77
pH Pen Testers	AI1223	SX620	pH Pen Tester Kit (± 0.01 pH)	\$ 118.77

850/8500 Series Portable pH/Conductivity Meters	AI5510	PH850	Value Series Portable pH Meter Kit	\$ 412.58
850/8500 Series Portable pH/Conductivity Meters	AI5512	EC850	Value Series Portable Conductivity/TDS Meter Kit	\$ 450.09
850/8500 Series Portable pH/Conductivity Meters	AI5514	PC850	Value Series Portable pH/Conductivity/TDS Meter Kit	\$ 612.63
850/8500 Series Portable pH/Conductivity Meters	AI5520	PH8500	Premium Series Portable pH Meter Kit	\$ 500.10
850/8500 Series Portable pH/Conductivity Meters	AI5522	EC8500	Premium Series Portable Conductivity/TDS/Salinity Meter Kit	\$ 575.12
850/8500 Series Portable pH/Conductivity Meters	AI5524	PC8500	Premium Series Portable pH Conductivity/TDS/Salinity Meter Kit	\$ 737.65
PH850 Portable pH Meters for Specific Applications	AI5543	PH850-SS	Value Series Portable Spear pH Meter Kit for Food	\$ 575.12
PH850 Portable pH Meters for Specific Applications	AI5545	PH850-MT	Value Series Portable Blade Spear pH Meter Kit for Meat	\$ 600.12
PH850 Portable pH Meters for Specific Applications	AI5535	PH850-DP	Value Series Portable pH Meter Kit for Liquid Food (milk, yogurt, sauce...)	\$ 493.85
PH850 Portable pH Meters for Specific Applications	AI5541	PH850-PW	Value Series Portable pH Meter Kit for Purified Water	\$ 537.61
PH850 Portable pH Meters for Specific Applications	AI5531	PH850-WW	Value Series Portable pH Meter Kit for Wastewater	\$ 431.34
PH850 Portable pH Meters for Specific Applications	AI5537	PH850-SL	Value Series Portable pH Meter Kit for Soil	\$ 493.85
PH850 Portable pH Meters for Specific Applications	AI5539	PH850-MS	Value Series Portable pH Meter Kit for Micro Samples	\$ 493.85
PH850 Portable pH Meters for Specific Applications	AI5533	PH850-BR	Value Series Portable pH Meter Kit for Beverage Making	\$ 450.09
PH850 Portable pH Meters for Specific Applications	AI5529	PH850-HF	Value Series Portable pH Meter Kit for Strong Acid Solutions	\$ 493.85
PH850 Portable pH Meters for Specific Applications	AI5527	PH850-SA	Value Series Portable pH Meter Kit for Strong Alkali and High Salinity Solutions	\$ 493.85

PH850 Portable pH Meters for Specific Applications	AI5564	PH850-HT	Value Series PH850-HT Portable pH Meter for High Temperature Liquid and Caustic Solutions	\$ 493.85
PH850 Portable pH Meters for Specific Applications	AI5560	PH850-FT	Value Series PH850-FT Portable pH Meter Kit with Flat pH Electrode	\$ 493.85
PH850 Portable pH Meters for Specific Applications	AI5590	PH850-A	Value Series PH850 Portable pH Meter with Temperature Probe Only	\$ 311.31
PH8500 Portable pH Meters for Specific Applications with Data Logger	AI5557	PH8500-SS	Premium Series Portable Spear pH Meter for Food with GLP Data Logger	\$ 712.65
PH8500 Portable pH Meters for Specific Applications with Data Logger	AI5558	PH8500-MT	Premium Series Portable Blade Spear pH Meter for Meat with GLP Data Logger	\$ 818.92
PH8500 Portable pH Meters for Specific Applications with Data Logger	AI5552	PH8500-DP	Premium Series Portable pH Meter for Liquid Food with GLP Data Logger	\$ 612.63
PH8500 Portable pH Meters for Specific Applications with Data Logger	AI5556	PH8500-PW	Premium Series Portable pH Meter for Purified Water with GLP Data Logger	\$ 662.64
PH8500 Portable pH Meters for Specific Applications with Data Logger	AI5559	PH8500-WW	Premium Series Portable pH Meter for Wastewater with GLP Data Logger	\$ 537.61
PH8500 Portable pH Meters for Specific Applications with Data Logger	AI5553	PH8500-SL	Premium Series Portable pH Meter for Soil with GLP Data Logger	\$ 612.63
PH8500 Portable pH Meters for Specific Applications with Data Logger	AI5554	PH8500-MS	Premium Series Portable pH Meter for Micro Samples with GLP Data Logger	\$ 612.63
PH8500 Portable pH Meters for Specific Applications with Data Logger	AI5551	PH8500-BR	Premium Series Portable pH Meter for Beverage Making with GLP Data Logger	\$ 562.62
PH8500 Portable pH Meters for Specific Applications with Data Logger	AI5525	PH8500-HF	Premium Series Portable pH Meter for Strong Acid with GLP Data Logger	\$ 612.63



PH8500 Portable pH Meters for Specific Applications with Data Logger	AI5523	PH8500-SA	Premium Series Portable pH Meter for Alkali and High Salinity Solutions with GLP Data Logger	\$ 612.63
PH8500 Portable pH Meters for Specific Applications with Data Logger	AI5566	PH8500-HT	Premium Series PH8500-HT Portable pH Meter for High Temperature Liquid and Caustic Solutions with GLP Data Logger	\$ 687.64
PH8500 Portable pH Meters for Specific Applications with Data Logger	AI5562	PH8500-FT	Premium Series PH8500-FT Portable pH Meter Kit with Flat pH Electrode with GLP Data Logger	\$ 612.63
PH8500 Portable pH Meters for Specific Applications with Data Logger	AI5592	PH8500-A	Premium Series PH8500 Portable pH Meter with Temperature Probe Only with GLP Data Logger	\$ 436.34
910/950 Value Series Benchtop Meters	AI5610	PH910	Value Series Benchtop pH Meter Kit	\$ 487.60
910/950 Value Series Benchtop Meters	AI5612	EC910	Value Series Benchtop Cond./TDS Meter Kit	\$ 537.61
910/950 Value Series Benchtop Meters	AI5614	PC910	Value Series Benchtop pH/Cond./TDS Meter Kit	\$ 787.66
910/950 Value Series Benchtop Meters	AI5620	PH950	Value Series Benchtop pH Meter Kit with TestBench	\$ 612.63
910/950 Value Series Benchtop Meters	AI5622	EC950	Value Series Benchtop Cond./TDS/Salinity Meter Kit with TestBench	\$ 687.64
910/950 Value Series Benchtop Meters	AI5624	PC950	Value Series Benchtop pH/Cond./TDS/Salinity Meter Kit with TestBench	\$ 950.19
9500 Premium Series GLP Benchtop Meters	AI5630	PH9500	Premium Series Research-grade GLP Benchtop pH Meter Kit with TestBench (printer sold separately)	\$ 812.67
9501 Premium Series GLP Benchtop Meters	AI5632	EC9500	Premium Series Research-grade GLP Benchtop Cond./TDS/Salinity/Resistivity Meter Kit with TestBench (printer sold separately)	\$ 912.69
9502 Premium Series GLP Benchtop Meters	AI5634	PC9500	Premium Series Research-grade GLP Benchtop pH/Cond./TDS/Salinity/Resistivity Meter Kit with TestBench (printer sold separately)	\$ 1,375.28

700 Series Benchtop Meters	AI501	PH700	Benchtop pH/mV Meter Kit	\$ 367.58
701 Series Benchtop Meters	AI502	EC700	Benchtop Conductivity/TDS Meter Kit	\$ 412.58
800 Series Benchtop Meters	AI521	PH800	Benchtop pH Meter Kit	\$ 537.61
800 Series Benchtop Meters	AI522	EC800	Benchtop Conductivity/TDS/Salinity/Resistivity Meter Kit	\$ 600.12
800 Series Benchtop Meters	AI523	PC800	Benchtop pH/Conductivity/TDS/Salinity/Resistivity Meter Kit	\$ 862.68
820 Series Precision Benchtop Meters	AI531	PH820	High Precision Benchtop pH Meter Kit	\$ 737.65
821 Series Precision Benchtop Meters	AI532	EC820	High Precision Benchtop Conductivity/TDS/Salinity/Resistivity Meter Kit	\$ 812.67
822 Series Precision Benchtop Meters	AI533	PC820	High Precision Benchtop pH/Conductivity/TDS/Salinity/Resistivity Meter Kit	\$ 1,112.73
400 Series Portable Meters	AI411	PH400	Portable pH Meter Kit	\$ 375.08
400 Series Portable Meters	AI412	EC400	Portable Conductivity/TDS Meter Kit	\$ 412.58
400 Series Portable Meters	AI413	PC400	Portable pH/Conductivity/TDS Meter Kit	\$ 612.63
400S Series Portable Meters	AI421	PH400S	Portable pH Meter Kit with Datalogger and USB Data Output	\$ 481.35
400S Series Portable Meters	AI422	EC400S	Portable Conductivity/TDS/Salinity/Resistivity Meter Kit with Datalogger and USB Data Output	\$ 541.36
400S Series Portable Meters	AI423	PC400S	Portable pH/Conductivity/TDS/Salinity/Resistivity Meter Kit with Datalogger and USB Data Output	\$ 821.42
Optical DO Meters	AI480	DO850	Portable Optical Dissolved Oxygen Meter	\$ 1,351.74
Optical DO Meters	AI485	DO8500	Portable Optical Dissolved Oxygen Meter with Auto. Salinity Compensation and datalogger	\$ 1,622.62
Turbidity Meters	AI481	TN400	Portable Turbidity Meter (infrared light source, Compliant with ISO7027)	\$ 912.69
Turbidity Meters	AI483	TN500	Portable White Light Turbidity Meter with GLP Data Logger, Compliant with EPA 180.1	\$ 1,475.30
Turbidity Meters	AI486	TN480	Portable Infrared Turbidity Meter with GLP Data Logger, Compliant with ISO7027	\$ 1,275.26
Turbidity Meters	AI482	TN420	Portable White Light Turbidity Meter, Compliant with EPA 180.1	\$ 975.20
Water Hardness Meters	AI4731	YD300	Portable Water Hardness Meter Kit	\$ 850.17
Fluoride Meters	AI4725	WS100	Portable Fluoride/pH Meter Kit	\$ 837.67
Fluoride Meters	AI4721	WS200	Portable Fluoride/pH/Conductivity/TDS/Salinity/Resistivity Meter Kit	\$ 987.70
SX700 Series Portable Meters	AI469	SX751	Handheld pH/ORP/DO/Conductivity/TDS/Salinity/Resistivity Multiparameter Meter Kit	\$ 1,076.78
SX700 Series Portable Meters	AI468	SX736	Handheld pH/DO/Conductivity/TDS/Salinity/Resistivity Multiparameter Meter Kit	\$ 1,015.83

SX700 Series Portable Meters	AI467	SX731	Handheld pH/ORP/Conductivity/TDS/Salinity/Resistivity Multiparameter Meter Kit	\$ 652.84
SX700 Series Portable Meters	AI465	SX725	Handheld pH/DO Meter Kit	\$ 784.22
SX700 Series Portable Meters	AI463	SX721	Handheld pH/ORP Meter Kit	\$ 399.56
SX700 Series Portable Meters	AI461	SX716	Handheld DO Meter Kit	\$ 625.75
SX700 Series Portable Meters	AI477	SX723	Handheld pH/Conductivity/TDS/Salinity/Resistivity Multiparameter Meter Kit	\$ 474.06
SX700 Series Portable Meters	AI475	SX713	Handheld Conductivity/TDS/Salinity/Resistivity Meter Kit	\$ 311.52
SX700 Series Portable Meters	AI473	SX712	Handheld ORP Meter Kit	\$ 338.61
SX700 Series Portable Meters	AI460	SX711	Handheld pH Meter Kit	\$ 311.52
SX700 Series Portable Meters	AI471	SX713-02	Handheld Pure Water Conductivity/TDS/Salinity/Resistivity Meter Kit	\$ 528.23

Electrodes & Accessories

Product Series	SKU	Model	Product Description	public price
PC-body pH Electrodes	AI1301	201-C	pH Electrode (gel, BNC connector)	\$ 40.01
PC-body pH Electrodes	AI3115	201DJ-C	pH Electrode (double-junction, BNC connector)	\$ 53.76
PC-body pH Electrodes	AI1302	201T-F	pH Electrode (ATC, BNC/RCA connector)	\$ 87.52
PC-body pH Electrodes	AI1312	201T-S	pH Electrode (ATC, 8-pin connector)	\$ 87.52
PC-body pH Electrodes	AI1319	201DJ-F	pH Electrode (double-junction, ATC, BNC/RCA connector)	\$ 106.27
ORP Electrodes	AI1315	301Pt-S	ORP Electrode (platinum pin, 8-pin connector)	\$ 103.77
ORP Electrodes	AI1303	301Pt-C	ORP Electrode (platinum ring, BNC connector)	\$ 106.27
ORP Electrodes	AI1321	301DJ-C	Double-junction ORP Electrode (platinum ring, BNC connector)	\$ 125.03
ORP Electrodes	AI1323	301DJ-CG	Double-junction ORP Electrode (Gold pin, BNC connector)	\$ 150.03
Conductivity Electrodes	AI1304	2301-C	Conductivity Electrode (k=1, BNC connector)	\$ 72.51
Conductivity Electrodes	AI1308	2401-C	Conductivity Electrode (k=1, BNC connector)	\$ 111.27
Conductivity Electrodes	AI1305	2301T-F	Conductivity Electrode (k=1, ATC, BNC/RCA connector)	\$ 87.52
Conductivity Electrodes	AI1313	2301T-S	Conductivity Electrode (k=1, ATC, 8-pin connector)	\$ 87.52
Conductivity Electrodes	AI1309	2401T-F	Conductivity Electrode (k=1,ATC/RCA connector)	\$ 117.52
Conductivity Electrodes	AI1306	2310-C	Conductivity Electrode (k=10, BNC connector)	\$ 111.27
Conductivity Electrodes	AI1307	2310T-F	Conductivity Electrode (k=10,ATC, BNC/RCA connector)	\$ 152.53
Conductivity Electrodes	AI1310	DJS-0.1-C	Conductivity Electrode (k=0.1, with flow cell, BNC connector)	\$ 225.05
Conductivity Electrodes	AI1311	DJS-0.1-F	Conductivity Electrode (k=0.1,ATC, with flow cell, BNC/RCA connector)	\$ 275.06
Temperature Electrode	AI1350	MP500	Temperature Electrode (30K thermister)	\$ 71.26
LabSen Combination pH Electrodes (BNC connector)	AI3101	LabSen 211	Glass-body pH Electrode for General Water Solutions (TRIS compatible)	\$ 143.78
LabSen Combination pH Electrodes (BNC connector)	AI3111	LabSen 231	Glass-body pH Electrode for Wastewater, Suspension, Emulsion	\$ 136.28

LabSen Combination pH Electrodes (BNC connector)	AI3102	LabSen 221	Glass-body pH Electrode for Viscous and Low Ion Concentration Samples	\$ 237.55
LabSen Combination pH Electrodes (BNC connector)	AI3131	LabSen 331	POM-body pH Electrode for Wastewater Treatment	\$ 100.02
LabSen Combination pH Electrodes (BNC connector)	AI3107	LabSen 371	PVC-body Flat pH Electrode	\$ 188.79
LabSen Combination pH Electrodes (BNC connector)	AI3103	LabSen 241-6	Glass-body pH Electrode for Small Samples (>0.2mL)	\$ 235.05
LabSen Combination pH Electrodes (BNC connector)	AI3104	LabSen 241-3	Glass-body pH Electrode for Micro Samples (>30 microliters)	\$ 276.31
LabSen Combination pH Electrodes (BNC connector)	AI3151	LabSen 241-3SP	Glass-body Micro pH Electrode for Micro-Volume Solutions Containing Proteins	\$ 311.31
LabSen Combination pH Electrodes (BNC connector)	AI3161	LabSen 241-180	Glass-body pH Electrode for NMR Tubes (3mm diameter, 180mm length)	\$ 368.83
LabSen Combination pH Electrodes (BNC connector)	AI3105	LabSen 251	Glass-body Spear pH Electrode for Semi-Solid Samples	\$ 201.29
pH Electrodes	AI3106	LabSen 551	PVC-body Spear pH Electrode	\$ 212.54
LabSen Combination pH Electrodes (BNC connector)	AI3157	LabSen 751	Stainless Steel Spear pH Electrode for Solid Food Samples	\$ 268.81
LabSen Combination pH Electrodes (BNC connector)	AI3159	LabSen 761	Stainless Steel Blade Spear pH Electrode for Meat	\$ 436.34

LabSen Combination pH Electrodes (BNC connector)	AI3108	LabSen 801	Glass-body pH Electrode for Purified Water& Low Ion-Concentration Solutions	\$ 257.55
LabSen Combination pH Electrodes (BNC connector)	AI3141	LabSen 821	Glass-body pH Electrode for Dairy Products and Liquid Food	\$ 161.28
LabSen Combination pH Electrodes (BNC connector)	AI3109	LabSen 831	Glass-body pH Electrode for Strong Acid Samples and HF (Hydrofluoric) Acid	\$ 186.29
LabSen Combination pH Electrodes (BNC connector)	AI3145	LabSen 841	Glass-body pH Electrode for Strong Alkali and High Salinity Solutions	\$ 186.29
LabSen Combination pH Electrodes (BNC connector)	AI3110	LabSen 851-1	Glass-body pH Electrode for Viscous Solutions	\$ 210.04
LabSen Combination pH Electrodes (BNC connector)	AI3155	LabSen 851-S	Glass-body pH Electrode for Highly Viscous Solutions	\$ 353.82
LabSen Combination pH Electrodes (BNC connector)	AI3149	LabSen 881	Glass-body pH Electrode for Low Temperature Solutions	\$ 218.79
LabSen Combination pH Electrodes (BNC connector)	AI3167	LabSen 861	Glass-body pH Electrode for High-temp. and Caustic Solutions	\$ 236.30
LabSen Combination pH Electrodes (BNC connector)	AI3168	LabSen 851-H	Pre-Pressurized Glass-body pH Electrode for Highly Viscous Samples with Strong Alkali or High Temperature	\$ 301.31
LabSen Combination pH Electrodes (BNC connector)	AI3147	LabSen 871	Glass-body Double-junction pH Electrode for non-aquaous/Organic Solutions	\$ 361.32
LabSen pH/Temp. Electrodes (BNC/RCA connector)	AI3121	LabSen 213	Glass-body pH/Temp. Electrode for General Water Solutions (TRIS compatible)	\$ 166.28

LabSen pH/Temp. Electrodes (BNC/RCA connector)	AI3123	LabSen 223	Glass-body pH/Temp. Electrode for Viscous and Low Ion Concentration Samples	\$ 265.05
LabSen pH/Temp. Electrodes (BNC/RCA connector)	AI3133	LabSen 333	POM-body pH/Temp. Electrode for Wastewater Treatment	\$ 146.28
LabSen pH/Temp. Electrodes (BNC/RCA connector)	AI3137	LabSen 803	Glass-body pH/Temp. Electrode for Purified Water& Low Ion-Concentration Solutions	\$ 281.31
LabSen pH/Temp. Electrodes (BNC/RCA connector)	AI3201	LabSen 753	Titanium Spear pH/Temp. Electrode for Solid Food Samples	\$ 312.56
LabSen pH/Temp. Electrodes (BNC/RCA connector)	AI3202	LabSen 763	Stainless Steel Blade Spear pH/Temp. Electrode for Meat	\$ 476.35
LabSen pH/Temp. Electrodes (BNC/RCA connector)	AI3113	LabSen 553	PVC-body Spear pH/Temp. Electrode	\$ 231.30
LabSen pH/Temp. Electrodes (BNC/RCA connector)	AI3153	LabSen 823	Glass pH/Temp. Electrode for Dairy Products and Liquid Food	\$ 170.03
LabSen pH/Temp. Electrodes (BNC/RCA connector)	AI3117	LabSen 243-6	Glass pH/Temp. Electrode for Small Samples (>0.2mL)	\$ 248.80
LabSen pH/Temp. Electrodes (BNC/RCA connector)	AI3179	LabSen 843	Glass-body pH/Temp. Electrode for Strong Alkali and High Salinity Solutions	\$ 211.29
LabSen pH/Temp. Electrodes (BNC/RCA connector)	AI3177	LabSen 833	Glass-body pH/Temp. Electrode for Strong Acid Samples and HF (Hydrofluoric) Acid	\$ 223.80
LabSen pH/Temp. Electrodes (BNC/RCA connector)	AI3165	LabSen 853-S	Glass-body pH/Temp. Electrode for Highly Viscous Solutions	\$ 372.58
LabSen pH/Temp. Electrodes (BNC/RCA connector)	AI3176	LabSen 863	Glass-body pH Electrode for High-temp. and Caustic Solutions	\$ 273.81
LabSen pH/Temp. Electrodes (BNC/RCA connector)	AI3172	LabSen 373	PVC-body Flat pH Electrode	\$ 225.05
LabSen pH/Temp. Electrodes (BNC/RCA connector)	AI3219	LabSen 873	Glass-body Double-junction pH/Temp. Electrode for non-aquaeous/Organic Solutions	\$ 411.33
LabSen pH/Temp. Electrodes (8-pin connector)	AI3231	LabSen 335	POM-body pH/Temp. Electrode for Wastewater Treatment	\$ 141.49
LabSen pH/Temp. Electrodes (8-pin connector)	AI3233	LabSen 805	Glass-body pH/Temp. Electrode for Purified Water& Low Ion-Concentration Solutions	\$ 290.27

LabSen pH/Temp. Electrodes (8-pin connector)	AI3235	LabSen 246-5	Glass pH/Temp. Electrode for Small Samples (>0.2mL)	\$ 261.10
LabSen pH/Temp. Electrodes (8-pin connector)	AI3237	LabSen 845	Glass-body pH/Temp. Electrode for Strong Alkali and High Salinity Solutions	\$ 231.92
LabSen pH/Temp. Electrodes (8-pin connector)	AI3239	LabSen 835	Glass-body pH/Temp. Electrode for Strong Acid Samples and HF (Hydrofluoric) Acid	\$ 231.92
LabSen pH/Temp. Electrodes (8-pin connector)	AI3241	LabSen 865	Glass-body pH Electrode for High-temp. and Caustic Solutions	\$ 261.10
LabSen pH/Temp. Electrodes (8-pin connector)	AI3243	LabSen 855	Pre-Pressurized Glass-body pH/Temp. Electrode for Highly Viscous Samples with Strong Alkali or High Temperature	\$ 246.51
LabSen pH/Temp. Electrodes (8-pin connector)	AI3245	LabSen 765	Titanium Blade Spear pH/Temp. Electrode for Meat	\$ 348.61
Water Hardness Electrode	AI4739	601-S	Water Hardness Electrode (built-in temperature sensor)	\$ 362.57
Fluoride Electrode	AI4729	F501-S	3-in-1 Fluoride Electrode (built-in temperature sensor)	\$ 170.03
60 Series Pocket Testers Replacement Probes	AI1201	PH60-E	Replacement pH/temp. Electrode	\$ 53.76
61 Series Pocket Testers Replacement Probes	AI1203	PH60F-E	Replacement Flat pH/temp. Electrode	\$ 78.77
62 Series Pocket Testers Replacement Probes	AI1205	PH60S-E	Replacement Spear pH/temp. Electrode	\$ 218.79
63 Series Pocket Testers Replacement Probes	AI1207	ORP60-E	Replacement ORP/temp. Electrode for ORP60-Z	\$ 75.02
64 Series Pocket Testers Replacement Probes	AI1209	EC60-E	Replacement Conductivity/temp. Electrode	\$ 57.51
65 Series Pocket Testers Replacement Probes	AI1211	PC60-E	Replacement pH/Conductivity/temp. Electrode	\$ 100.02
60-Z Smart Testers Replacement pH Probes	AI3711-E	PH60-DE	Replacement pH/temp. Electrode (double-junction, gel inner solution)	\$ 68.76
60-Z Smart Testers Replacement pH Probes	AI3715-E	PH60F-DE	Replacement Flat pH/temp. Electrode (double-junction, gel inner solution)	\$ 90.02
60-Z Smart Testers Replacement pH Probes	AI3713-E	PH60S-DE	Replacement Spear pH/temp. Electrode (double-junction, gel inner solution)	\$ 222.55

60-Z Smart Testers Replacement pH Probes	AI3717-E	ORP60-DE	Replacement ORP/temp. Electrode for ORP60-Z (double-junction)	\$ 86.27
60-Z Smart Testers Replacement pH Probes	AI3718-E	EC60-DE	Replacement Conductivity/temp. Electrode	\$ 57.51
60-Z Smart Testers Replacement pH Probes	AI3719-E	PC60-DE	Replacement pH/Conductivity/temp. Electrode (double-junction, gel inner solution)	\$ 118.77
60-Z Smart Testers Replacement pH Probes	AI3717-DA	ORP60-DA	Replacement ORP Electrode for PC60-Z, PH60-Z, PH60S-Z, PH60F-Z (double-junction)	\$ 86.27
pH Pen Tester Replacement Probes	AI1213	SX615	Replacement pH Electrode for SX610 pH Pen Tester	\$ 47.51
pH Pen Tester Replacement Probes	AI1214	SX625	Replacement pH Electrode for SX620 pH Pen Tester	\$ 60.01
Replacement ORP Probe for 5041	AI238	301-G	Replacement ORP Electrode for 5041 Compact ORP Pen Tester	\$ 37.51
GroStar Pens Replacement Probes	AI801G	GS1-E	Replacement pH Probe for GS1 pH Pen	\$ 63.76
GroStar Pens Replacement Probes	AI802G	GS2-E	Replacement Soil pH Probe for GS2 Soil pH Pen	\$ 216.71
GroStar Pens Replacement Probes	AI803G	GS3-E	Replacement EC Probe for GS3 EC/ppm Pen	\$ 67.51
GroStar Pens Replacement Probes	AI804G	GS4-E	Replacement pH/EC Combo Probe for GS4 pH/EC/ppm Combo Pen	\$ 100.02
GroStar Pens Replacement Probes	AI805G	GS5-E	Replacement ORP Electrode for PC60-Z, PH60-Z, PH60S-Z, PH60F-Z (double-junction)	\$ 86.27
Accessories for 9500 Series GLP Benchtop Meters	AI5671	TH192G	Pin-Type Mini Printer Set for 9500 Series Benchtop Meter, including printer×1, power adapter×1, data cable×1, ink ribbon×2, printing paper×2	\$ 312.56
Accessories for 9500 Series GLP Benchtop Meters	AI5672	TH5740	Printing Paper (14 rolls per box) for TH192G Printer	\$ 25.01
Accessories for 9500 Series GLP Benchtop Meters	AI5673	ORO-09	Ink ribbon (black, 5 pieces per box) for TH192G Printer	\$ 25.01
Accessories for 9500 Series GLP Benchtop Meters	AI5674	K1000	Mini USB Keyboard	\$ 37.51

Accessories for Optical DO Meters	AI4101	DO803	Optical Dissolved Oxygen Probe (3-meter cable)	\$ 897.06
Accessories for Optical DO Meters	AI4102	DO810	Optical Dissolved Oxygen Probe (10-meter cable)	\$ 1,037.09
Accessories for Optical DO Meters	AI4104	2301T-F-3M	Salinity Probe (3-meter cable)	\$ 102.52
Accessories for Optical DO Meters	AI4105	2301T-F-10M	Salinity Probe (10-meter cable)	\$ 117.52
Accessories for Optical DO Meters	AI4103	DO8032	Luminescent Sensor Cap	\$ 218.79
Accessories for Turbidity Meters	AI4201	TN400-S1	AMCO Turbidity Standard Calibration Solution Kit (0,20,100,800NTU) for TN400 and TN500i	\$ 189.62
Accessories for Turbidity Meters	AI4202	TN400-S2	Sample Vials (pack of 3)	\$ 33.55
Accessories for Turbidity Meters	AI4203	TN400-S3	Silicone Oil for wiping the glass vials	\$ 17.50
Accessories for Turbidity Meters	AI4205	T500-S1	AMCO Turbidity Standard Calibration Solution Kit (20/100/400/800NTU) for TN500	\$ 189.62
Accessories for Turbidity Meters	AI4204	T200-1	T200-1 Turbidity Standard Solution Set (20/100/400/800 NTU) for TN400 and TN480	\$ 189.62
Accessories for Turbidity Meters	AI4206	TN500-S2	AMCO 0.0NTU Standard Calibration Solution 100ml	\$ 109.40
Accessories for Turbidity Meters	AI4207	TN500-S3	Sample Vials (pack of 6)	\$ 67.10
Accessories for Turbidity Meters	AI4208	TN500-S4	Rechargeable Lithium Batteray	\$ 36.47
Accessories for Turbidity Meters	AI4209	TN500-S5	Replaceable Tungsten Filament Lamp for TN500	\$ 160.45
Accessories for Turbidity Meters	AI4211	TN500-S6	USB charger with cable	\$ 29.17
Accessories for SX700 DO Meters	AI4110	DO500	DO500 Polarographic DO Electrode with 2-meter cable	\$ 325.07

Accessories for SX700 DO Meters	AI4110-M	DO500-M	DO500 Polarographic DO Electrode with 5-meter cable	\$ 362.57
Accessories for SX700 DO Meters	AI4110-L	DO500-L	DO500 Polarographic DO Electrode with 10-meter cable	\$ 425.09
Accessories for SX700 DO Meters	AI4113	DO503	Polarographic DO Electrode Replacement Membrane Caps	\$ 60.01
Accessories for SX700 DO Meters	AI4112	DO502	DO Electrode Inner Solution Refill (30ml)	\$ 25.01
Magnetic Stirrers	AI2801	801	Powerful Magnetic Stirrer	\$ 107.52
Magnetic Stirrers	AI2901	901	Intelligent Magnetic Stirrer	\$ 175.04
Magnetic Stirrers	AI2802	802	Heating Magnetic Stirrer	\$ 185.04
Others	AI2803	601	Magnetic Stirrer with Electrode Holder	\$ 120.02
Others	AI2810	602	Flexible Electrode Holder	\$ 72.51
Others	AI2910	N/A	CalPod Solution Organizer	\$ 18.75
Others	AI7103	N/A	8-pin to BNC/RCA Converter	\$ 18.75
Others	AI7109	N/A	100-240V Switching Power Adapter	\$ 25.01
Others	AI7110	N/A	100-240V Switching Power Adapter with Universal Plugs	\$ 37.51
Others	AI7112	N/A	8-pin Connection Adapter for 60 Series Testers	\$ 21.88

Solutions

Product Series	SKU	Product Description	public price
Solutions	AI1101	pH 7.00 Calibration Buffer Solution 8 oz.	\$ 12.50
Solutions	AI1102	pH 4.00 Calibration Buffer Solution 8 oz.	\$ 12.50
Solutions	AI1103	pH 10.01 Calibration Buffer Solution 8 oz.	\$ 12.50
Solutions	AI1101-L	pH 7.00 Calibration Buffer Solution 16 oz.	\$ 18.75
Solutions	AI1101-L	pH 4.00 Calibration Buffer Solution 16 oz.	\$ 18.75
Solutions	AI1101-L	pH 10.01 Calibration Buffer Solution 16 oz.	\$ 18.75
Solutions	AI1130-S	1413 μ S/cm Conductivity Calibration Solution 4 oz.	\$ 12.50
Solutions	AI1132-S	84 μ S/cm Conductivity Calibration Solution 4 oz.	\$ 12.50
Solutions	AI1131-S	12.88 mS/cm Conductivity Calibration Solution 4 oz.	\$ 18.75
Solutions	AI1134-S	111.8 mS/cm Conductivity Calibration Solution 4 oz.	\$ 25.01
Solutions	AI1130	1413 μ S/cm Conductivity Calibration Solution 8 oz.	\$ 15.00
Solutions	AI1132	84 μ S/cm Conductivity Calibration Solution 8 oz.	\$ 15.00
Solutions	AI1131	12.88 mS/cm Conductivity Calibration Solution 8 oz.	\$ 22.50
Solutions	AI1134	111.8 mS/cm Conductivity Calibration Solution 8 oz.	\$ 30.01
Solutions	AI1130-L	1413 μ S/cm Conductivity Calibration Solution 16 oz.	\$ 22.50
Solutions	AI1132-L	84 μ S/cm Conductivity Calibration Solution 16 oz.	\$ 22.50
Solutions	AI1131-L	12.88 mS/cm Conductivity Calibration Solution 16 oz.	\$ 33.76
Solutions	AI1134-L	111.8 mS/cm Conductivity Calibration Solution 16 oz.	\$ 45.01
Solutions	AI1170	pH Electrode Maintenance Probe Care Kit, including pH 4/7/10 buffers (4oz), 3M KCL Soaking Solution (4oz), Probe Cleaning Solution (8oz), a CalPod Solution Organizer	\$ 62.51
Solutions	AI1114	pH Calibration Solution Kit 7.00, 4.00 (16oz), plus Apera CalPod Solution Organizer	\$ 37.51
Solutions	AI1116	pH Calibration Solution Kit 7.00, 4.00, 10.01 (16oz), plus CalPod Solution Organizer	\$ 50.01
Solutions	AI1112	pH/Cond. Calibration Solution Kit pH 7.00, 4.00, 10.01 (8oz), 1413 μ S, 12.88mS, 84 μ S (8oz)+3M KCL Soaking solution (4oz)+CalPod Solution Organizer	\$ 100.02
Solutions	AI1108	pH Calibration Solution Kit 7.00, 4.00 (8oz)+3M KCL Soaking solution (4oz)	\$ 31.26
Solutions	AI1109	pH Calibration Solution Kit 7.00, 4.00, 10.01 (8oz)+3M KCL Soaking solution (4oz)	\$ 37.51
Solutions	AI1104	pH Calibration Solution Kit (7.00, 4.00, 10.01) 16oz/bottle+3M KCL storage solution 4oz+CalPod Solution Organizer	\$ 68.76
Solutions	AI1105	pH/Cond. Calibration solution kit pH 7.00, 4.00, 1413 μ S, 12.88 mS (8 oz)+3M KCL soaking solution (4oz)+CalPod Solution Organizer	\$ 75.02

Solutions		AI1106	3M KCl Soaking Solution 8 oz.	\$ 18.75
Solutions		AI1107	3M KCl Soaking Solution 4 oz.	\$ 12.50
Solutions		AI1172	1M LiCl Electrolyte Refill Solution 8 oz.	\$ 37.51
Solutions		AI1166	Electrode Cleaning Solution 8oz.	\$ 25.01
Solutions		AI1141	ORP Standard Solution (222mv) 4 oz.	\$ 18.75
Solutions		AI1141-L	ORP Standard Solution (650mv) 8 oz.	\$ 31.26
Solutions		AI1190	Protelyte Soaking Solution 4oz	\$ 25.01
Solutions		AI1121	pH 1.68 Calibration Buffer Solution 8oz	\$ 18.75
Solutions		AI1121-L	pH 1.68 Calibration Buffer Solution 16oz	\$ 31.26
Solutions		AI1123	pH 12.45 Calibration Buffer Solution 8oz	\$ 18.75
Solutions		AI1123-L	pH 12.45 Calibration Buffer Solution 16oz	\$ 31.26