



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

Item#	Description	P.LISTA
00024Q	9 v. Battery Alkaline	\$ 44.33
001215	JUMPER, 2 POS. LO PRO	\$ 14.65
002819	BATT, 3V COIN CELL	\$ 20.62
00709318	MANUAL, WDM PIPESONDE PROBE,CHLORINE	\$ 281.37
007144	SEAL, MCDONALD, PS, ENGRAVED	\$ 406.64
007147	SEAL, FORD, PS, ENGRAVED	\$ 402.78
007160	MANUAL, WDM PIPESONDE PROBE	\$ 267.88
007179	SEAL, TEFLON, MUELLER, ENGRAVED	\$ 402.78
007183	HST REFRIGED SAMPLER FOR EM	\$ 35,653.20
007192	KIT, SEAL, FORD CORP.STOP, PIPESONDE	\$ 402.78
007193	KIT, SEAL, MCDONALD CORP., PIPESONDE	\$ 406.64
007194	KIT, SEAL, MUELLER CORP., PIPESONDE	\$ 406.64
010293	SCR. 6-32 X 3/4 SKT. CAP	\$ 5.59
010721	STANDOFF 6/32 X 1 1/8 MF	\$ 41.43
011458	TUBE, ACID, 1500	\$ 229.34
011583	MANIFOLD BOTTOM SPARGER TOC ANALYZER	\$ 3,488.23
011622	WEIGHT, REAGENT, PVC, 1/8 NPT	\$ 206.21
012039	CABLE DB9F/DB9F 6 FT. PIN/PIN	\$ 204.28
013090-01	COVER, UV LAMP FINISHED	\$ 1,031.05
013509	COMPRESSOR, POC FURNACE 2001	\$ 94.63
013827	BARB ADAPTOR 1/4-28 x 0.164 (16)	\$ 68.42
013829	BARB ADAPTOR 1/4-28 x 0.102 (13)	\$ 68.42
013831	BARB, REDUCER 1/8 x 1/4 (16 x 24)	\$ 38.93
013832	FITTING, 1/4 NPT MALE x 1/4 BARB, PP	\$ 30.45
013849	FERRULE, 1/8 SF, ETFE YEL P359	\$ 71.88
013861	FITTING,1/4X1/8T,PP	\$ 75.93
013956	BARB,ADAPTOR 1/4-28 / 1/8"	\$ 191.37
014016	UNION, 1/8MPT X 1/8T (JACO)	\$ 53.38
014029	UNION, 1/2T X 1/8T (PFA)	\$ 1,108.14
014038	FTG, BARB, KYNAR, 1/8 X 3/16 (16X15)	\$ 60.71
014084	ELBOW, 1/8T X 1/8MP (JACO)	\$ 119.00
014657	BARB, 1/4-28 X 3/32 (7014)	\$ 75.93
014659	PLUG, 1/4-28	\$ 91.35
014699	REGULATOR PRESSURE, 0-30 PSI	\$ 7,227.00
014704	FLOW CONTROLLER, BK/GN	\$ 4,933.63
015804	FUSE, 1A, 5x20, ANTISURGE	\$ 102.14
015805	FUSE, T, 2A, 250V 5x20, ANTISURGE	\$ 16.19
016163	SWITCH,MICRO PCB (MINITURE)	\$ 123.73
020008	GUARD,FAN,92mm,HT	\$ 87.69
020011	FAN, 2.36", 24V	\$ 676.45
021049	MEMBRANE, HYDROPHOBIC, .63 X .87	\$ 217.77
021135	COMB AMMONIA ION SELECTIVE ELECTRODE	\$ 3,006.43
021136	COMB NITRATE ION SELECTIVE ELECTRODE	\$ 3,102.79
021140	COMB COPPER ION SELECTIVE ELECTRODE	\$ 2,601.72
024045	MOTOR, W/MTG. BOSS 1RPM/115V	\$ 2,139.19



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

024505	PUMP HEAD, 7013-20	\$	3,179.88
024506	PUMP HEAD, 7014-20	\$	3,179.88
024507	PUMP HEAD, 7015-20	\$	3,179.88
024508	PUMP HEAD, 7016-20	\$	3,179.88
026009	TUBING, 1/8" OD X .06" ID PFA PER FT	\$	58.97
026013	TUBING TYGON 7015 (PER FT)	\$	53.58
026043	aa TUBING NORPRENE 7013 P/FT, 50 FT	\$	33.73
026044	aa TUBING NORPRENE 7014 P/FT, 50 FT	\$	39.70
026045	aa TUBING NORPRENE 7015 P/FT, 50 FT	\$	50.88
026046	aa TUBING NORPRENE7016 P/FT 50FT MIN	\$	33.73
026059	TUBE PTFE .062 X .03 (PER INCH)	\$	14.16
026067	TUBING, PTFE 1/8" OD X .03"ID PER FT	\$	89.04
026071	aa TUBING NORPRENE 7024 P/FT, 50 FT	\$	59.36
08310=A=0000	SENSOR, 2-ELECTRODE COND. 8310	\$	1,667.03
08311=A=0000	SENSOR, 2-ELECTRODE COND. 8311	\$	2,004.29
08312=A=0000	SENSOR, 2-ELECTRODE COND. 8312	\$	2,042.83
08313=A=0001	PVC FLOWCHAMBER,8313.1,W/3X3/4 FNPT	\$	1,031.05
08315=A=0000	SENSOR, 2-ELECTRODE COND. 8315	\$	4,490.38
08315=A=0002	2 electrodes conductivity sensor k=0,01 (8315.2 for Yokogawa flow chamber)	\$	2,929.34
08315=A=1111	2 electrodes conductivity sensor k=0,01 3/4"G thread (8315.1) with test certificate	\$	4,509.65
08316=A=0000	SENSOR, 2-ELECTRODE COND., 8316	\$	5,550.34
08317=A=0000	2 electrodes conductivity sensor k=1 3/4 NPT Thread with test certificate	\$	6,012.86
08318=A=0001	FLOW-CHAMBER, SST 8318.1 WITH	\$	1,734.48
08319=A=0000	female connector 6+T with connexion drawing	\$	211.99
08319=A=0005	CABLE, 5 M + IP 65 CONNECTOR	\$	805.57
08319=A=0010	CABLE, 10 M + IP 65 CONNECTOR	\$	1,566.81
08319=A=0020	CABLE, 20M + IP 65 CONNECTOR	\$	1,597.65
08319=A=1115	Sensor cable with IP65 plug, 5mFor 83XX sensors and gateways	\$	576.23
08346=A=8302	MOTOR, SYNCHRONOUS GEAR COMPLE	\$	9,713.09
08346=C=0418	MOTOR, FOR PROBE 8346	\$	3,892.94
08346=C=8302	Gearbox for probe 8346	\$	4,143.48
08350=A=0003	8350 pH sensor (10meter cable) with special glass for HF applications (8350.3)	\$	1,293.15
08350=A=1005	TUBE, (PP) WITH ADJUSTABLE FL	\$	822.91
08350=A=1010	TUBE, (PP) WITH ADJUSTABLE FL	\$	824.84
08350=A=1015	tube (PP) with adjustable flange PN10 ND32 (PE 42) L= 1500 MM	\$	859.53
08350=A=1105	tube (PP) 500 mm, with 2 fixing clips for 8350/8351	\$	969.38
08350=A=1110	tube (PP) 1000 mm with 2 fixing clips for 8350/8351	\$	1,121.63
08350=A=1115	tube (PP) 1500 mm with 2 fixing tubes dor 8350/8351	\$	1,056.11
08350=A=6010	TUBE WITH FIXED FLANGE, 1000MM LONG	\$	1,682.45
08350=A=6015	TUBE WITH FIXED FLANGE, 1500MM LONG	\$	1,248.83
08350=A=7000	KIT FOR CHEMICAL CLEANING PP	\$	1,075.02



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

08350=A=8000	ENCLOSURE W/ BUILT-IN PREAMP	\$	2,081.38
08350=A=9500	PVC DN40 tee for 8350/8351 probes with thread union	\$	2,717.35
08350=C=0004	tee SENSOR, PH 8350.4 W/10M. HI CBL.	\$	3,372.60
08350=C=0005	tee 8350 FLAT PH SENSOR W/CPVC BODY	\$	1,643.90
08351=C=0000	SENSOR, 8351 REDOX W/ 10M HIGH	\$	2,466.82
08361=C=0000	PROBE, PH FOR 8361	\$	3,411.14
08361=C=1000	ELECTRODE, PT 100 FOR 8361	\$	1,085.01
08362=A=0000	tee 8362 HI PURITY WATER PH SENSOR ASSY	\$	5,319.07
08362=A=0010	tee ULTRAPURE WATER PH SYSTEM 8362 W/10M CABLE	\$	6,899.38
08362=A=1001	TEMPERATURE SENSOR (8362)	\$	1,081.16
08362=A=1111	8362 HI PURITY WATER ORP SENSOR ASSY	\$	5,319.07
08362=A=2000	ELECTRODE, PH (8362)	\$	2,447.54
08362=A=2100	ELECTRODE, ORP (8362)	\$	2,370.46
08362=A=3001	CABLE FOR TEMP SENSOR, 3M (8362)	\$	422.06
08362=A=3002	CABLE FOR TEMP SENSOR, 10M (8362)	\$	423.98
08362=A=3003	CABLE 20m temperature sensor	\$	551.18
08362=A=8300	KIT OF 5X GASKETS DN40 FOR 8362 PH	\$	514.56
08362=C=4000	CHAMBER, CLEAR FLOW THROUGH (8362)	\$	1,734.48
08365=A=0010	Probe	\$	2,235.55
08394=A=0510	Kit for 8394 with EPDM gasket, clamp and 316L SS welding ferrule (H=13mm) for 2"(38mm) clamp installation	\$	1,243.04
08394=A=1500	COND PROBE, 2 ELECTRODE (9125)	\$	3,468.96
08394=A=2000	SENSOR, 2-ELECTRODE COND. 8314	\$	6,070.68
08394=A=2011	SONDE 8314+CERTIFAT RUGOSITE/MATIERE	\$	8,614.58
08398=A=3000	8398 DN50 MODEL	\$	4,933.63
08398=A=5000	8398 1"G thread inductive cond.probe	\$	5,531.06
08398=A=5004	8398 replacement kit for 8335,4 with 8398.5, adapter and nut	\$	5,319.07
08398=A=5005	8398 replacement kit for 8335,5 with 8398.5, adapter and nut	\$	5,781.60
08398=A=6000	Kit (316LSS tee, gasket, nut) for 8398.3 in DN50 pipes	\$	2,409.00
08398=A=8300	Kit (316LSS chamber, gasket, nut)-8398.3 in derivation	\$	5,261.26
08398=C=6300	HANGER, PROBE, SST	\$	352.68
08544=A=0002	AIR/WATER CLEANING SYSTEM 110 VAC	\$	982.87
08544=A=0004	NETTOYAGE CHIMIQUE 220 VAC	\$	6,494.66
08544=A=0005	8544.5 incl.air-powered hose pump, electrovalve 110 VAC	\$	11,833.01
08544=A=0006	NETTOYAGE CHIMIQUE 24 VAC	\$	6,494.66
085G0021	OXY 1100 25mm Membrane Crtg Single	\$	928.91
085G0028	OXY 1100 Sensor 50 um KBr	\$	1,724.84
085G3309	OXY Membrane Cartridge Holder	\$	144.35
085G4085	Mounting Bracket For 1.5" PVC Pipe	\$	1,073.45
08810=A=3001	EEPROM 8810 ISE VERSION	\$	1,339.40
08810=A=4000	BOARD, MEASUREMENT	\$	1,838.55
08810=A=4600	8810 V2 USER CONNECTION BOX IN METAL	\$	2,100.65
08810=A=5300	AUTO SAMPLE DILUTION, INTERNAL	\$	10,657.42
08810=A=6250	8810 reagent/calibration level detector	\$	1,160.17
08810=A=7000	CARTE AFFI EQUI 8810 SANS 4/20	\$	8,730.22



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

08810=A=7050	Conditionning peristaltic pump 50 Hz	\$	12,950.78
08810=A=7060	POMPE PERISTAL 60 Hz-20 Tr/mn-1K7	\$	7,419.72
08810=A=7160	24 VAC 60 Hz micro pump only for 8810	\$	9,809.45
08810=A=7200	Auto.calibr.pump for 8810 premounted w/o canister/level detector	\$	10,445.42
08810=A=8300	Assembly ready-mounting: Electrolyte bridge + Reference-Elect.8429 + electrolyte reservoir-8810	\$	2,967.89
08810=A=8308	Electrolyte bridge ready-mounting on 8810	\$	1,106.21
08810=C=5062	PISTON DE DILUTEUR VERSION 3 8810	\$	653.32
08810=C=5063	PISTON (HYDRAULIC SIDE) FOR 8810 DIL	\$	2,062.10
08810=C=5066	SEAL FOR 8810 DIL. 10X14,3X1,4 EPDM	\$	368.10
08810=C=5067	SHEET GASKET 5X9,3X1,4 FOR 8810 DIL.	\$	368.10
08811=A=0000	ELECTRONIC UNIT COMPLETE	\$	23,550.38
08811=A=1000	CARTE MICRO EQUIPEE SEQUENCEUR 8811	\$	8,498.95
08811=A=2323	3--ch_8811 discontinuous 110 V with inches connectors (Sample = Drain = 3/8-1/2" tube, Instrument air = 1/6-1/4" tube)	\$	61,111.51
08811=A=2326	6--ch_8811 discontinuous 110 V with inches connectors (Sample = Drain = 3/8-1/2" tube, Instrument air = 1/6-1/4" tube)	\$	89,672.62
08811=A=6200	EPROM STANDARD SEQUENCEUR 8811	\$	749.68
08878=A=1500	vv ProbeIMMERSION, 8878.4, 1.5 METER	\$	911.57
08878=C=3600	MEMBRANE PROTECTION CAP	\$	256.32
08891=A=6000	PUMP TUBE SET (ISMAPRENE TUBES	\$	2,524.63
09073=C=0340	POROUS CARTRIDGE	\$	632.12
09073=C=0350	DIPA BOTTLE CAP ADAPTER	\$	487.58
09078=A=2000	CELLULE OXISTAT 9078 COMPLETE	\$	1,626.56
09078=C=1010	PROBE, BODY (9078)	\$	2,254.82
09078=C=1020	WASHER, SEALING FOR FILL PLUG (9078)	\$	35.46
09078=C=1030	PLUG, ELECTROLYTE FILL (9078)	\$	42.01
09120=A=1500	MEASURING BOARD, DUAL CONDUCTIVITY	\$	2,717.35
09123=A=7010	Resin kit (includes 2 filters, 2L sachet_resin, funnel, instructions)	\$	1,135.12
09123=C=0200	CATIONIC RESIN EMPTY CARTRIDGE	\$	1,468.53
09123=C=0500	Flat gasket, Viton	\$	135.87
09123=C=0600	Resin cartridge filter	\$	358.46
09125=A=2485	RS485 JBUS/MODBUS W/BOARD	\$	585.87
09125=A=4000	RELAY BOARD, 9125, 9135	\$	1,778.81
09126=A=1000	Yearly recalibration_the standard 9126 PURECAL in our factory. For optional COFRAC certificate see P/N TE9126=A=000	\$	2,926.94
09135=A=1501	pH-measurement-board compatible with CPU-after-Feb2004	\$	2,177.74
09180=A=0501	Equipped shield plate with earth terminal and fixing system for amperometric transmitters (9181, 9182, 9183, 9184 all versions, 9185, 9186 and 9187)	\$	670.67
09180=A=1501	AMPEROMETRIC BRD.,REV.E OR LATER	\$	3,256.97
09180=A=8010	CLARK CELL 10m CONNEXION CABLE	\$	867.24
09181=A=3600	ELECTROLYTE GEL FOR OXISTAT D.O. MTR	\$	254.75
09181=A=3600-LM	ELECTROLYTE GEL FOR OXISTAT D.O. MTR	\$	254.75



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

09181=A=4001	SYSTEME DE NETTOYAGE 9181 240 VAC	\$	3,353.33
09182=A=1000	OXYGEN ELECTRODE W/NTC SENSOR FOR 9500 AMPERO ANALYZER	\$	5,762.33
09182=A=1200	Caps-calibration (in air)-9182 sensor (and 9078 - 8878)	\$	389.29
09183=C=4400	Adaptor F 1/4" NPT vs. M 1/2"	\$	543.47
09184=A=1001	ELECTRODE W/O CABLE (9184, 9187)	\$	6,533.21
09184=A=3500	BOX OF 4 MEMBRANES 9184 CHLOROMAT	\$	1,270.02
09184=A=3600	ELECTROLYTE 100 ML FOR 9184 CHLOROMT	\$	433.62
09184=A=3600-LM	ELECTROLYTE 100 ML FOR 9184 CHLOROMT	\$	433.62
09184=A=4020	1/4 IN. TUBING ADAPTER	\$	167.09
09184=A=4300	ELECTRODE CABLE, DIGITAL AMP	\$	1,023.34
09184=A=4400	CABLE EQUIPE AS7 AMPERO	\$	1,052.25
09184=C=1030	ELECTROLYTE FILLING PLUG (9184)	\$	53.00
09184=C=4010	FLOW REGULATOR	\$	333.41
09184=C=4100	REPLACEMENT PROBE BODY	\$	389.29
09185=A=1000	SENSOR, FOR 9185 OZONE ANALYZER	\$	6,089.95
09185=A=3500	MEMBRANES FOR POLYMETRON O2 /O3	\$	1,360.60
09185=A=3600	ELECTROLYTE 100 ML FOR 9185 OZONMAT	\$	319.79
09185=A=3600-LM	ELECTROLYTE 100 ML FOR 9185 OZONMAT	\$	319.79
09185=A=6211	KIT, REPLACEMENT CPU W/ 9185-03 FIRMWARE FOR 91XX METERS	\$	4,683.10
09186=A=0100	MEASURING CELL (9186)	\$	7,747.34
09186=A=0200	CONDITIONING BOTTLE W/TAP&CART 9186	\$	2,524.63
09186=A=0300	WORKING ELECTRODE (9186 HYDRASTAT)	\$	4,721.64
09186=A=0400	FLOW CONTROLLER (9186)	\$	13,683.12
09186=A=0500	CONNECTION CABLE FOR 9186 HYDRASTAT	\$	1,449.25
09186=A=8000	SPARE PARTS KIT, 2 YR. 9186	\$	8,306.23
09186=C=0110	VENTURI BLK & BODY W/O INJ NOZL,9186	\$	4,297.66
09186=C=0120	SS COUNTER ELECTRODE (9186)	\$	1,085.01
09186=C=0360	ORION, MONOETHYLAMINE BOTTLE CAP ADPT	\$	892.29
09187=A=3500	REPLACEMENT MEMBRANES FOR 9187 ELEC	\$	1,352.89
09187=A=3600	ELECTROLYTE 100 ML FOR 9187 CHLOR	\$	310.76
09200=A=1040	BOARD, CPU W/DISPLAY	\$	4,933.63
09200=A=1045	CPU BOARD W/DISPLAY, SODIUM	\$	5,627.42
09200=A=1110	CPU board with display for Monec E9200, specifically loaded with 9210-silica phase1	\$	6,803.02
09200=A=2000	BOARD, POWER SUPPLY, MONEC E9200	\$	3,353.33
09200=A=5012	INTERFACE BOARD, SODIUM	\$	2,177.74
09200=A=5036	BOARD, USER CONNECTIONS INTERFACE, 9240 SODIUM (MULTI) ANALYZER (1-4ch)	\$	3,044.98
09200=A=5500	921X Internal bus communication module	\$	2,620.99
09200=A=5510	INTERNAL BUS COMM MODULE, SODIUM	\$	1,131.27
09200=A=5511	MODULE, INTERNAL BUS COMMUNICATION, 9240 (CONDUCTIVITY)	\$	1,152.47
09210=A=0640	CAPSULE-FILTERS (4X)	\$	210.06



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

09210=A=0650	Weights (4x), 9210 Phase 2 reagent pump	\$	312.21
09210=A=8540	PROGRAMMING INTERFACE FOR 9100/9200	\$	6,340.49
09210=C=0330	Spacer ring-921X mirror	\$	123.73
09210=C=0340	ee Mirror-921X-NLZs (SILKOSTAT, PHOSPHAMAT,...)	\$	609.00
09210=C=7010	9210 SILKOSTAT DRY CHEMICALS FOR 50 DAYS	\$	1,559.10
09210=C=7010-LM	9210 SILKOSTAT DRY CHEMICALS FOR 50 DAYS	\$	1,559.10
09211=A=3122	9211 PHOSPHAMAT 0. .5ppm PO4, ½blue method, with wall-mount enclosure, mit Profibus DP Option, 2 -ch_, including solid reagents kit used to start up the-NLZ.	\$	118,368.62
09211=C=7001	Set_dry chemicals-9211 PHOSPHAMAT 0. .50 ppm PO4 (except H2SO4), 45 days worth_operation	\$	1,649.56
09240=A=0260	Level detector, floating design, 924x	\$	867.24
09240=A=0320	CABLE AS7, LENGTH, 1M, SODIUM	\$	260.17
09240=A=0510	TAP FOR DIPA BOTTLE W/QUICK CONNECTOR	\$	160.15
09240=A=0700	ELECTROLYTE RESERVOIR, SODIUM	\$	637.90
09240=A=1500	POTENTIOMETRIC MODULE, SODIUM	\$	1,559.10
09240=A=7170	ANTI-BACKFLUSH VALVE	\$	362.31
09240=A=8000	1 YEAR SPARE PART KIT SODIUM	\$	2,794.44
09240=A=8001	2 YRS SPARE KIT FOR 924X SOLENOID VALVES REPL.	\$	2,582.45
09240=A=8010	SERVICE KIT SODIUM	\$	4,432.56
09240=A=8020	REAGENT BOTTLES, SODIUM	\$	177.42
09240=A=8102	Second (2nd) Extra-Channel complete for 9240 analyser -- To proceed this P/N , please provide the analyzer-S/N and the CPU-S/N of the instrument it will be install onto --	\$	7,169.18
09240=A=8103	Third (3rd) Extra--ch_ complete-9240-NLZ -- To proceed this P/N , please provide the analyzer-S/N and the CPU-S/N_the instrument it will be install onto --	\$	7,227.00
09240=A=8104	Fourth (4th) Extra--ch_ complete-9240-NLZ -- To proceed this P/N , please provide the analyzer-S/N and the CPU-S/N_the instrument it will be install onto --	\$	7,227.00
09240=A=8315	KIT, INSTALL FOR 9240 SODIUM ANALZR	\$	5,627.42
09240=A=9110	Complete sampling block (capable 1-4 -ch_3/2way valve) equipped-1-ch_ - Mounted on 9240-1ch	\$	7,573.90
09240=A=9120	Complete sampling block (capable 1-4 channels_3/2way valve) equipped for 2 channels - Mounted on 9240-2ch	\$	10,618.87
09240=A=9140	Complete sampling block (capable 1-4 -ch_3/2way valve) equipped-4-ch_ - Mounted on 9240-4ch	\$	11,042.86
09240=A=9171	IP65 CONNECTOR FOR 9240 ELECTROVALVE	\$	343.04
09240=A=9172	CONNECTOR, IP 65, ELECTROVALVE W/MOLDED CABLE	\$	211.99
09240=A=9201	OverFlowVessel (complete) for 924x (with inlet and outlet quick fittings - Service only)	\$	3,623.14
09240=A=9300	MEASURING CELL FOR SODIUM ANALYZER	\$	1,478.16
09240=A=9510	DIPA BOTTLE W/POROUS CARTRIDE, TAP	\$	1,649.56



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

09240=A=9665	Locking -shaft (hook model)-924x-NLZs (all version SN above or equals to 86117). Includes all necessary items-mounting (on frame-side). For locking-hook (on panel side), see 09240=A=9660.	\$ 734.26
09240=A=9840	BRUSH, CLEANING, STATIC HEAT EXCHANGER SYSTEM	\$ 275.59
09240=C=0201	DOOR, 9240/9245 SODIUM ANALYZER	\$ 1,821.20
09240=C=0230	LOCKING SCREW FOR TUBE 3.2MM-1/8"	\$ 10.60
09240=C=0256	CONDUCTIVITY SENSOR, SODIUM	\$ 1,456.96
09240=C=0258	-Conduc. sensor in 9240-cell_1035mm (for tightening, see seal P/N 429=120=030 ; -mounting , see threaded sleeves P/N 359015,20035)	\$ 1,428.06
09240=C=0300	SODIUM MEASUREMENT CELL	\$ 260.17
09240=C=0310	REFERENCE ELECTRODE, SODIUM	\$ 1,092.72
09240=C=0320	SODIUM ION SELECTIVE ELECTRODE	\$ 1,522.49
09240=C=0340	Pump connection board for 924X	\$ 468.31
09245=A=8300	US (FITTINGS & OM) KIT, PK/2	\$ 362.31
09245=A=8310	CALIBRATION, AUTOMATIC, OPTIONAL	\$ 14,222.74
09245=A=8315	KIT, INSTALLATION FOR 9245 SODIUM	\$ 5,627.42
09245=A=9100	COMPLETE SAMPLING BLOCK FOR 1 CHANNEL	\$ 3,777.31
09523=A=7000	POLYMETRON 9X23 RESIN CARTRIDGE	\$ 2,736.62
09523=A=7010	POLYMETRON 9X23 RESIN KIT	\$ 1,167.88
09523=A=9000	POLYMETRON 9X23 Hydraulic kit	\$ 8,845.85
09523=C=0100	POLYMETRON 9X23 MEASUREMENT CELL	\$ 2,659.54
09523=C=0200	POLYMETRON 9X23 EMPTY RESIN CARTRIDGE	\$ 1,478.16
09523=C=0300	POLYMETRON 9X23 RESIN CARTRIDGE CAP	\$ 354.60
09523=C=0400	POLYMETRON 9X23 RESIN TUBE	\$ 256.32
1000A3025	ASSY-CABLE,WITH BNC CONN.	\$ 310.28
1000A3316-102	ASSY-FLOW CHAMBER,316 SS	\$ 1,566.81
1000A3335-001	ASSY, LCP SPRAY WASHER, 1-1/2"	\$ 703.43
1000A3335-004	ASSY, PHD SPRAY WASHER,1"	\$ 489.51
1000A3398-001	AIR COMPR SUB ASY 115V	\$ 3,989.30
1000A3398-002	AIR COMPR SUB ASY 230V	\$ 6,128.50
1000A3430-001	KIT,RS485-TO-RS232 CONVER	\$ 2,004.29
1000A3434-001	KIT,HVIEW SOFTWARE	\$ 6,571.75
1000A3455-001	ASY-PANEL WALL MOUNT	\$ 229.34
1000A3457-001	KIT-UNIVERSAL ADAPTER	\$ 398.93
1000A4A1048	CABLE, .1"CL,5.4"LG,15CND	\$ 277.52
1000A4A3146	ASSY-PANEL, BACKLIGHTING	\$ 512.64
1000A7F1037	TOOL, DEPTH GAGE	\$ 61.28
1000C4F3070	COVER FOR X/P ENCLOS.,LEV	\$ 294.86
1000C4F3217-101	FAB-BRACKET,MOUNTING,	\$ 25.63
1000D4F1067	FAB-ENCL,REAR,670, (LTB)	\$ 462.53
1000D4F3206	si FAB,SUB-PANEL,NGA	\$ 229.34
1000F1058-001	MODEL 33 MOLDED ENCLOSURE	\$ 416.28
1000F1236-111	si FAB-TUBE FITTNG,KYNAR,1/2	\$ 189.25
1000F1236-122	si FAB-TUBE FITTNG,KYNAR,3/4	\$ 292.93



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

1000F3293-007	si FAB-SEALING HUB,PVDF	\$	711.86
1000F3321-001	FAB-ADAPTER pHOX DO SENSO	\$	275.59
1000F3374-002	1' SENSOR GUARD, PEEK	\$	194.65
1000F3374-003	1in. SENSOR GUARD, RYTON	\$	235.12
1000F3382-001	FAB-FITTING,ADPTR,INGOLD,K	\$	229.34
1000F3389-001	si FAB-SWIVEL FRONT	\$	56.27
1000F3399-001	Enclosure, Airblaster	\$	3,989.30
1000F3407-00882	si PIPE,1-1/2",PVC80,SLOT,	\$	296.79
1000F3407-00888	si PIPE,1-1/2",PVC80.88L,NPT	\$	260.17
1000F3407-01202	PIPE,1-1/2",PVC80,SLOT,	\$	634.05
1000F3508-001	COLLAR FOR 1/4" SHAFT	\$	235.12
1000F6045-04802	PIPE 1"PVC.SCH.80,48"THD.	\$	48.76
1000F6047-04802	PIPE,1 1/2, PVC80, 48.0L	\$	92.31
1000F6243-01852	si PIPE,1/2,316SS 80,18.5LG,	\$	553.11
1000F6343-01852	PIPE, 1/2, CPVC80, 18.5L	\$	271.74
1000F6343-04802	si PIPE, 1/2, CPVC80, 48.0L	\$	213.92
1000G1064	ASSY,DOOR,MODEL 670	\$	300.64
1000G1103	ASSY-CASE & DOOR;MOD 670	\$	1,131.27
1000G1107-101	ASSY-ENCL,DOOR,670 LOGO	\$	601.29
1000G1110	ASSY-KIT,ENCL MTG;670 TYP	\$	210.06
1000G1113	ASSY-BARRIER & BRACKET	\$	5,897.23
1000G1212	BOM-PROGRM EPROM FOR 95P(\$	1,248.83
1000G3073-030	BOM-CABLE,LOW TEMP,30'	\$	703.43
1000G3074-060	CABLE,Hi Temp Triax 60ft	\$	3,353.33
1000G3075-101	X-P ENCL KIT,ALUMINUM,BAS	\$	528.05
1000G3088-001	Sun Shield with Hardware	\$	1,285.44
1000G3296	KIT-FERRITE TORROID CORE	\$	876.88
1000G3316-101	ASSY-FLOW CHAMBER,KYNAR	\$	616.70
1000G3501-002	HARDWARE, BULK CHEMICAL	\$	48.76
10049	ACETIC ACID, GLACIAL ACS 500ML	\$	169.11
100735	rr Arsenic Dual Range Test Kit	\$	321.84
100735RGT	rr REAGENT SET, DUAL RANGE ARSENIC REPL.	\$	296.79
100744	ee CLEAN & CAL pH MOUNT, ACCESSORY, US	\$	2,871.53
100745	ee CLEAN & CAL pH MOUNT, ACCESSORY, EU/UK	\$	2,871.53
100747	ee CLEAN & CAL pH MOUNT, ACCESSORY, APAC/LATAM	\$	3,314.78
100748	CLEAN & CAL pH MOUNT, ACCESSORY, BATTERY OPTION	\$	2,871.53
100749	BATTERY, POWER BANK, USB, 5000mAH	\$	89.04
100763	PANEL, MOUNTING, ACCESSORY, CLEAN & CAL	\$	612.85
100866	FLUORESCENCE SAMPLE VIALS PK/6	\$	26.21
100930	BRACKET, CDx440sc	\$	101.37
100934-05	CABLE ASSY, RS232, 5 METERS, FOR H2S	\$	252.46
100934-10	CABLE ASSY, RS232, 10 METERS, FOR H2S	\$	308.35
100934-30	CABLE ASSY, RS232, 30 METERS, FOR H2S	\$	456.75
100935-05	CABLE ASSY, POWER, 5 METERS, FOR H2S	\$	252.46
100935-10	CABLE ASSY, POWER, 10 METERS, FOR H2S	\$	308.35



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

100935-30	CABLE ASSY, POWER, 30 METERS, FOR H2S	\$	456.75
100936	MODIFIED, CONN, RCPT 5 CKT CIRCULAR, GATEWAY	\$	144.15
101335	CABLE PULL-THROUGH TOOL FOR CABLE GUARD	\$	61.48
102133	TITRAVER TITRANT, 0.100M 29ML	\$	39.19
102133-LM	TITRAVER TITRANT, 0.100M 29ML	\$	39.19
102149	EDTA STD SOLN, 0.20N 500ML	\$	101.35
102149-LM	EDTA STD SOLN, 0.20N 500ML	\$	101.35
102233	MAGNESIUM SULFATE STD SOLN 29ML	\$	66.13
1030414	O-RING NIT .070W X 1.114ID	\$	15.61
1032511	SCREW 4-40X .375PAN SST PHIL	\$	9.83
103769	FERROUS IRON RGT PP 25ML PK/100	\$	93.04
103769-LM	FERROUS IRON RGT PP 25ML PK/100	\$	93.04
104128	KTO: MOLYBDATE RGT HR PWD PW PK/1000	\$	906.99
104128-LM	KTO: MOLYBDATE RGT HR PWD PW PK/1000	\$	906.99
104199	MOLYBDATE RGT HR PP 25ML PK/100	\$	125.03
104199-LM	MOLYBDATE RGT HR PP 25ML PK/100	\$	125.03
104228	KTO: ACID REAGENT PP 25ML PK/1000	\$	1,167.16
104228-LM	KTO: ACID REAGENT PP 25ML PK/1000	\$	1,167.16
104299	ACID REAGENT PP 25ML PK/100	\$	158.80
104299-LM	ACID REAGENT PP 25ML PK/100	\$	158.80
1043811	SCREW 10-32X2.5 PAN SST PHIL	\$	9.06
1043911	SCREW 10-32X2.25 PAN SST PHIL	\$	8.48
104399	CHLORIDE 2 IND PWD PLWS PK/100	\$	71.19
104399-LM	CHLORIDE 2 IND PWD PLWS PK/100	\$	71.19
104517	SODIUM HYDROXIDE 1.000N, 3.60L	\$	234.88
104532	SODIUM HYDROXIDE,1.00N 100ML MDB	\$	57.82
104542	SODIUM HYDROXIDE, 1.0N 100ML	\$	50.23
104542-LM	SODIUM HYDROXIDE, 1.0N 100ML	\$	50.23
104553	SODIUM HYDROXIDE, 1.00N 1000ML	\$	89.78
104553-LM	SODIUM HYDROXIDE, 1.00N 1000ML	\$	89.78
104633	FILTRATION AID SOLN 29ML DB	\$	56.01
104732	SODIUM ARSENITE SOLN, 5G/L 100ML	\$	50.23
104732-LM	SODIUM ARSENITE SOLN, 5G/L 100ML	\$	50.23
10501H	AMMONIUM CHLORIDE AR 454G	\$	137.48
105599	SULFAMIC ACID PWD PLWS PK/100	\$	62.32
105599-LM	SULFAMIC ACID PWD PLWS PK/100	\$	62.32
105766	CHLORIDE 2 IND PP 100ML PK/50	\$	77.33
105766-LM	CHLORIDE 2 IND PP 100ML PK/50	\$	77.33
105932	DETERGENTS REAGENT 100 TESTS MDB	\$	83.47
105932-LM	DETERGENTS REAGENT 100 TESTS MDB	\$	83.47
10604	AMMONIUM HYDROXIDE(28% AS NH3) 2.5L	\$	184.29
1061717	MOUNT, ADHESIVE BACKED	\$	16.57
106199	PHOSPHATE 2 RGT PWD PLWS PK/100	\$	86.72
10649	AMMONIUM HYDROXIDE(28% AS NH3) 500ML	\$	117.97
10682401	CON CRIMP Pin .040"dia.	\$	43.17



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

107166	MANGANOUS SULFATE PWD PLWS PK/50	\$	49.31
107166-LM	MANGANOUS SULFATE PWD PLWS PK/50	\$	49.31
107266	ALKALINE IODIDE-AZIDE PP PK/50	\$	50.03
107266-LM	ALKALINE IODIDE-AZIDE PP PK/50	\$	50.03
107399	SULFAMIC ACID PP 300ML PK/100	\$	107.85
1075411	SCREW 1/4 -20X .5 PAN SST PHIL	\$	5.59
1076416	MAGNET, STIR BAR 1.500L X .313D	\$	53.19
107728	KTO: POTASSIUM IODIDE PWD PLW PK/1000	\$	775.10
107760	KI REFILL VIAL, 50 grams	\$	193.32
107799	POTASSIUM IODIDE PP PK/100	\$	80.03
107799-LM	POTASSIUM IODIDE PP PK/100	\$	80.03
108041	BEAKER, PP LOW FORM 50ML	\$	16.96
108042	BEAKER, PP LOW FORM 100ML	\$	19.46
108044	BEAKER, PP LOW FORM 150ML	\$	23.51
108046	BEAKER, PP LOW FORM 250ML	\$	32.38
108048	BEAKER, PP LOW FORM 400ML	\$	47.41
108052	BEAKER, PP LOW FORM 600ML	\$	57.24
108053	BEAKER, PP LOW FORM 1000ML	\$	60.71
108071	BEAKER, PP LOW FORM 50ML PK/12	\$	123.15
108072	BEAKER, PP LOW FORM 100ML PK/12	\$	157.64
108074	BEAKER, PP LOW FORM 150ML PK/12	\$	198.50
108076	BEAKER, PP LOW FORM 250ML PK/6	\$	139.14
108078	BEAKER, PP LOW FORM 400ML PK/6	\$	198.50
108082	BEAKER, PP LOW FORM 600ML PK/4	\$	172.10
108083	BEAKER, PP LOW FORM 1000ML PK/3	\$	162.27
108138	CYLINDER, GRADUATED PP 10ML	\$	60.13
108140	CYLINDER, GRADUATED PP 25ML	\$	65.52
108141	CYLINDER, GRADUATED PP 50ML	\$	72.85
108142	CYLINDER, GRADUATED PP 100ML	\$	89.04
108146	CYLINDER, GRADUATED PP 250ML	\$	112.93
108149	CYLINDER, GRADUATED PP 500ML	\$	146.47
108153	CYLINDER, GRADUATED PP 1000ML	\$	181.93
108246	FLASK, ERLLENMEYER PP 250ML	\$	53.77
108273	FLASK, ERLLENMEYER PP 125ML PK/12	\$	422.06
108276	FLASK, ERLLENMEYER PP 250ML PK/6	\$	271.74
108367	FUNNEL, ANALYTICAL PP 65MM	\$	25.63
108368	FUNNEL, ANALYTICAL 75MM POLPRPLN	\$	20.62
108370	FUNNEL, ANALYTICAL 100MM PP	\$	35.08
1085524	oo FUSE SLO-BLOW 1 AMP 3AG	\$	21.58
1097752	SEAL-CUP	\$	34.11
1101526	OUTLET 3 WIRE CONVENIENCE	\$	51.26
11032H	AMMONIUM MOLYBDATE RGT 100ML MDB	\$	50.23
110649	SILICA STD SOLN, 1MG/L 500ML	\$	99.73
111729	SILICA STD SOLN, 50MG/L 200ML	\$	77.33
1124224	N NAPHTHYL DIHYDROCHLORIDE 25G	\$	440.85



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

1124224-LM	N NAPHTHYL DIHYDROCHLORIDE 25G	\$	440.85
1125614	FERROUS AMMONIUM SULFATE, 113G	\$	113.46
1125614-LM	FERROUS AMMONIUM SULFATE, 113G	\$	113.46
11401	BARIUM CHLORIDE, 454G	\$	220.42
1142922	HOLDER, COLORIMETER CELL	\$	13.49
1143322	BOTTLE HOLDER-SPARE	\$	48.76
1145732	WINDOW PHOTOCELL	\$	16.96
1147911	SCREW 6-32X .438FLAT SST PHIL	\$	1.93
1148811	SCREW 4-40X .188PAN SST PHIL	\$	1.93
11523	MP CASE GASKET	\$	116.79
1157051	TRANSFORMER	\$	1,719.06
1160124	FUSE, 3AG, 1/8 AMP SLO BLO	\$	73.81
1160324	FUSE, 1/ 16A 250V SLOW-BLOW	\$	42.78
1170300	MICROSCOPE SLIDES PK/72	\$	95.78
1181724	FUSE, 15 A 250V SLOW-BLOW	\$	29.10
1182726	PLUG, 15 AMP, 125V AC OR DC	\$	41.43
12001	CALCIUM CARBONATE PRI STD 454G	\$	408.33
12014	CALCIUM CARBONATE, PRI STD 113G	\$	147.79
12014701	si UNION BULKHEAD 1/4-28 TO 10/32	\$	112.93
1202556	HEAT LIMITER-ASSY	\$	177.88
1206000	HUB 1/2 RIGID CONDUIT -SCRUTITE	\$	72.85
1206499	BARIVER 4 POWDER PILLOWS PK/100	\$	139.48
1206499-LM	BARIVER 4 POWDER PILLOWS PK/100	\$	139.48
1206514	SULFAVER 4, SULFATE RGT PWD 113G	\$	563.71
1206528	KTO: SULFAVER 4 PWD PLWS PK/1000	\$	968.42
1206528-LM	KTO: SULFAVER 4 PWD PLWS PK/1000	\$	968.42
1206599	SULFAVER 4 POWDER PILLOWS PK/100	\$	118.88
1206599-LM	SULFAVER 4 POWDER PILLOWS PK/100	\$	118.88
1206628	KTO: CHROMAVER 3 PWD PLW 25ML PK/1000	\$	1,253.88
1206628-LM	KTO: CHROMAVER 3 PWD PLW 25ML PK/1000	\$	1,253.88
1206699	CHROMAVER 3 PP 25ML PK/100	\$	164.05
1206699-LM	CHROMAVER 3 PP 25ML PK/100	\$	164.05
12153	CALCIUM STD SLN, 1000MG/L 1000ML	\$	122.50
1216607	CABLE, 3 COND 20 AWG STRD 2SHLD	\$	19.46
1218629	COD STD SOLN, 300MG/L 200ML	\$	80.03
1218629-LM	COD STD SOLN, 300MG/L 200ML	\$	80.03
1218649	COD STD SOLN, 300MG/L 500ML	\$	146.35
1218649-LM	COD STD SOLN, 300MG/L 500ML	\$	146.35
1218900	PIPETTE FILLER, 3-VALVE	\$	104.07
1222149	BUFFER SOLN, PH10.01 500ML	\$	42.09
1222156	BUFFER SOLN, PH10.01 4L	\$	172.53
1222220	BUFFER SOLN, PH7.00 25ML	\$	57.45
1222249	BUFFER SOLN, PH7.00 500ML	\$	42.09
1222256	BUFFER SOLN, PH7.00 4L	\$	172.17
1222349	BUFFER SOLN, PH4.01 500ML	\$	42.09



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

1222356	BUFFER SOLN, PH4.01 4L	\$	172.17
1225600	SPATULA, MICRO STAINLESS STEEL	\$	23.51
1225700	SCOOP, DOUBLE ENDED	\$	17.54
1226300	SYRINGE, 1CC LUER SLIP TIP PKGD	\$	8.67
1227642	ISOPROPYL ALCOHOL, TECH 100ML	\$	81.47
1228949	POTASSIUM IODIDE SOLN, 10% 500ML	\$	122.30
1228949-LM	POTASSIUM IODIDE SOLN, 10% 500ML	\$	122.30
1232120	CALCIUM SULFATE POWDER, 50 TESTS	\$	69.02
125=000=003	ELECTRODE, AMMONIA W/ELECT. & MEMB.	\$	9,327.65
125=000=004	ELECTRODE, COPPER	\$	12,141.36
125=000=007	Low level chloride ion-SELECTive electrode with AS9 connector	\$	11,505.38
125=000=203	ELECTROLYTE, INTRNL, AMM (8810)	\$	447.11
125=000=204	MEMBRANES, AMM, PK/20 (8810)	\$	678.37
125=010=003	STICK SILVER ELECTRODE	\$	2,447.54
125=010=004	MEASURING ELECTRODE	\$	2,832.98
1255149	SILVER NITRATE STD 0.141N 500ML	\$	147.79
1255149-LM	SILVER NITRATE STD 0.141N 500ML	\$	147.79
1261699	DITHIVER METALS RGT PP PK/100	\$	354.12
1261699-LM	DITHIVER METALS RGT PP PK/100	\$	354.12
1261742	LEAD STD SOLN, 100MG/L 100ML	\$	95.40
1261742-LM	LEAD STD SOLN, 100MG/L 100ML	\$	95.40
1268001	SCALE, INTERFERENCE WEDGE X	\$	33.34
126936	SILICONE OIL 15ML SCDB	\$	42.28
127032	SULFURIC ACID, 1.000N 100ML MDB	\$	48.23
127053	SULFURIC ACID, 1.0N 1000ML	\$	89.78
1271028	KTO: CHROMAVER 3 PP 5 OR 10ML PK/1000	\$	939.51
1271099	ee CHROMAVER 3 PP 5 OR 10ML PK/100	\$	117.97
1271099-LM	ee CHROMAVER 3 PP 5 OR 10ML PK/100	\$	117.97
1279142	MANGANESE STD, 1000MG/L 100ML	\$	62.32
1279142-LM	MANGANESE STD, 1000MG/L 100ML	\$	62.32
1279249	NITROGEN, NITRATE 1000MG/L 500ML	\$	99.73
1279249-LM	NITROGEN, NITRATE 1000MG/L 500ML	\$	99.73
1279400	TERM, PIN MALE, .062 DIA.	\$	12.72
1279642	LEAD STD SOLN, 1000MG/L 100ML	\$	95.93
1279642-LM	LEAD STD SOLN, 1000MG/L 100ML	\$	95.93
12842	COPPER STD SOLN, 100MG/L 100ML	\$	73.16
12842-LM	COPPER STD SOLN, 100MG/L 100ML	\$	73.16
12932	COPPER STD SOLN,10MG/L 100ML MDB	\$	91.23
1307611	md BOX, RSC 23X16-5/8X11 No. 23B	\$	16.57
1320100	TBG PVC 1" OD X 1/8" W	\$	91.35
1321A2F9994-101	vv FAB-SHAFT,PADDLE	\$	204.28
1321B2F9986-104	vv ADAPTOR,2" NPT FLOW-THRU,	\$	840.26
1321B2F9988-107	ADAPTOR,2" NPT INSERTION,	\$	928.91
13224	ERIOCHROME BLACK T, PRAC 25G	\$	119.78
13224-LM	ERIOCHROME BLACK T, PRAC 25G	\$	119.78



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

1324200	O-RING NIT .070W X .801ID	\$	5.59
13406	HYDROCHLORIC ACID, ACS 2.5L	\$	211.39
13449	HYDROCHLORIC ACID, ACS 500ML	\$	120.33
13449-LM	HYDROCHLORIC ACID, ACS 500ML	\$	120.33
135129	D-Fructose, 99%, 250 g	\$	98.09
1352900	300 ML MAGNETIC FILTER FUNNEL 47MM BASE	\$	969.38
1353000	FILTER, MEMB PLAIN 47MM PK/100	\$	449.04
1353000-LM	FILTER, MEMB PLAIN 47MM PK/100	\$	449.04
1353001	FILTER, MEMB GRID 47MM PK/200 GELMAN	\$	553.11
1353002	FILTER, MEMB GRID 47MM PK/1000	\$	2,023.56
1353702	SAMPLE CELL, 1 INCH SQ PK/2	\$	368.10
1361700	WINDOW, PHOTOCELL	\$	108.89
1377400	BRDG RECT 10.0 A. 200V. PIV	\$	105.23
1381401	FAUCET ADAPTER ASSY, W/HOSE FTG	\$	208.14
13949	IRON STD SOLN, 1MG/L 500ML	\$	90.87
1400153	POTASSIUM IODIDE-IODATE 1000ML	\$	152.13
1400153-LM	POTASSIUM IODIDE-IODATE 1000ML	\$	152.13
1402442	CADMIUM STD SOLN, 100MG/L 100ML	\$	83.47
1402442-LM	CADMIUM STD SOLN, 100MG/L 100ML	\$	83.47
1403268	ZINCOVER 5 ZINC REAGENT PP PK/25	\$	60.35
1403268-LM	ZINCOVER 5 ZINC REAGENT PP PK/25	\$	60.35
1403332	CYCLOHEXANONE 100ML MDB	\$	59.98
1403332-LM	CYCLOHEXANONE 100ML MDB	\$	59.98
1403353	CYCLOHEXANONE, ACS 1000ML	\$	196.94
1403353-LM	CYCLOHEXANONE, ACS 1000ML	\$	196.94
1403428	KTO: NITRAVER 5 PP 25ML PK/1000	\$	1,291.83
1403428-LM	KTO: NITRAVER 5 PP 25ML PK/1000	\$	1,291.83
140349	SILICA STD SOLN, 10MG/L 500ML	\$	109.30
1403499	NITRAVER 5 PWD PLWS 25ML PK/100	\$	173.62
1403499-LM	NITRAVER 5 PWD PLWS 25ML PK/100	\$	173.62
140349-LM	SILICA STD SOLN, 10MG/L 500ML	\$	109.30
1403546	NITRAVER 5 PWD PLWS 5ML PK/50	\$	56.18
1403599	NITRAVER 5 PWD PLWS 5ML PK/100	\$	107.68
1403599-LM	NITRAVER 5 PWD PLWS 5ML PK/100	\$	107.68
1403969	CYANIVER 3 RGT PP 25ML PK/100	\$	129.36
1403969-LM	CYANIVER 3 RGT PP 25ML PK/100	\$	129.36
1404099	CYANIVER 4 RGT PP 25ML PK/100	\$	129.36
1404099-LM	CYANIVER 4 RGT PP 25ML PK/100	\$	129.36
1404169	CYANIVER 5 RGT PP 25ML PK/100	\$	129.36
1404169-LM	CYANIVER 5 RGT PP 25ML PK/100	\$	129.36
14049	IRON STD SOLN, 10MG/L 500ML	\$	109.30
14049-LM	IRON STD SOLN, 10MG/L 500ML	\$	109.30
1405566	BUFFER PWD PLWS, PH6.00 PK/50	\$	75.33
1405940	BURET, STRAIGHT ACRYLIC 25ML	\$	331.48
1405941	BURET, STRAIGHT ACRYLIC 50ML	\$	373.88



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

1406041	FLASK, VOLUMETRIC PP 50ML	\$	95.01
1406042	FLASK, VOLUMETRIC PP 100ML	\$	101.76
1406049	FLASK, VOLUMETRIC PP 500ML	\$	151.09
1406428	DPD TOT CHLORINE PP 25ML PK/1000	\$	867.24
1406428-LM	DPD TOT CHLORINE PP 25ML PK/1000	\$	867.24
1406499	DPD TOT CHLORINE PP 25ML PK/100	\$	100.62
1406499-LM	DPD TOT CHLORINE PP 25ML PK/100	\$	100.62
1406528	KTO: NITRIVER 3, NITRITE RGT PK/1000	\$	1,145.48
1406528-LM	KTO: NITRIVER 3, NITRITE RGT PK/1000	\$	1,145.48
1406599	NITRIVER 3 PWD PLWS 25ML PK/100	\$	158.80
1406599-LM	NITRIVER 3 PWD PLWS 25ML PK/100	\$	158.80
1407028	DPD FREE CHLORINE 25ML PK/1000	\$	867.24
1407028-LM	DPD FREE CHLORINE 25ML PK/1000	\$	867.24
1407053	DPD FREE CHLORINE RGT PWD, 1 KG	\$	4,300.07
1407099	DPD FREE CHLORINE PP 25ML PK/100	\$	100.62
1407099-LM	DPD FREE CHLORINE PP 25ML PK/100	\$	100.62
1407628	DPD TOT CHLORINE PP 5ML PK/1000	\$	652.24
1407628-LM	DPD TOT CHLORINE PP 5ML PK/1000	\$	652.24
1407699	DPD TOT CHLORINE PP 5ML PK/100	\$	74.97
1407699-LM	DPD TOT CHLORINE PP 5ML PK/100	\$	74.97
1407728	DPD FREE CHLORINE PP 5ML PK/1000	\$	652.24
1407728-LM	DPD FREE CHLORINE PP 5ML PK/1000	\$	652.24
1407799	DPD FREE CHLORINE PP 5ML PK/100..	\$	74.97
1407799-LM	DPD FREE CHLORINE PP 5ML PK/100..	\$	74.97
1407828	KTO: NITRIVER 3 PWD PLWS 5ML PK/1000	\$	740.77
1407899	NITRIVER 3 PWD PLWS 5ML PK/100	\$	87.09
1407899-LM	NITRIVER 3 PWD PLWS 5ML PK/100	\$	87.09
1407995	BUFFER PWD PLWS, PH8.00 PK/15	\$	45.70
1408026	CDTA, MAGNESIUM SALT 100G	\$	581.77
1408026-LM	CDTA, MAGNESIUM SALT 100G	\$	581.77
1408099	CDTA MAGNESIUM SALT PP PK/100	\$	139.29
1408099-LM	CDTA MAGNESIUM SALT PP PK/100	\$	139.29
1408100	TEST KIT, NI-12 NITRATE-NITRITE	\$	572.38
1408100-LM	TEST KIT, NI-12 NITRATE-NITRITE	\$	572.38
1409417	OXALIC ACID 2 REAGENT 3.78L	\$	579.97
1409895	BUFFER PWD PLWS, PH6.86 PK/15	\$	45.70
1410415	EC MEDIUM MPN TUBES PK/15	\$	108.12
1410766	BUFFER, PH9.00 PP 50ML PK/50	\$	70.63
1411100	TEST KIT, CN-67 CHLORINE DPD PH	\$	418.20
1411100-LM	TEST KIT, CN-67 CHLORINE DPD PH	\$	418.20
1411928	KTO: NITRAVER 6 RGT PP 30ML PK/1000	\$	1,463.47
1411999	NITRAVER 6 RGT PP 30ML PK/100	\$	209.58
1411999-LM	NITRAVER 6 RGT PP 30ML PK/100	\$	209.58
1412099	NITRAVER 6 PWD PLWS 5ML PK/100	\$	137.48
1412099-LM	NITRAVER 6 PWD PLWS 5ML PK/100	\$	137.48
141410	STD COLOR SOLN 500 UNITS, 10ML	\$	215.00



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

141410-LM	STD COLOR SOLN 500 UNITS, 10ML	\$	215.00
141453	STD COLOR SOLN 500 UNITS, 1000ML	\$	635.98
141453-LM	STD COLOR SOLN 500 UNITS, 1000ML	\$	635.98
1414669	MOLYVER 1 REAGENT PP 25ML PK/100	\$	154.10
1414669-LM	MOLYVER 1 REAGENT PP 25ML PK/100	\$	154.10
1414869	MOLYVER 2 REAGENT PP 25ML PK/100	\$	154.10
1414869-LM	MOLYVER 2 REAGENT PP 25ML PK/100	\$	154.10
1415142	CHROMIUM, TRIVALENT 50MG/L 100ML	\$	89.06
1415142-LM	CHROMIUM, TRIVALENT 50MG/L 100ML	\$	89.06
1416066	NUTRIENT BFR SOLN PLW .5ML PK/50	\$	39.75
1416100	TEST KIT, NI-14 NITRATE LOW	\$	339.19
1416100-LM	TEST KIT, NI-14 NITRATE LOW	\$	339.19
1416369	LITHIUM HYDROXIDE PWD PLW PK/100	\$	117.97
1416369-LM	LITHIUM HYDROXIDE PWD PLW PK/100	\$	117.97
1416442	POTASSIUM PERMANGANATE SOL 100ML	\$	160.08
1417099	BOROVER 3 RGT PWD PLWS PK/100	\$	285.47
1417099-LM	BOROVER 3 RGT PWD PLWS PK/100	\$	285.47
1417442	ALUMINUM STD SOLN, 100MG/L 100ML	\$	89.06
1417442-LM	ALUMINUM STD SOLN, 100MG/L 100ML	\$	89.06
1417542	IRON STD SOLN, 100.0MG/L 100ML	\$	56.01
1417542-LM	IRON STD SOLN, 100.0MG/L 100ML	\$	56.01
1417642	NICKEL STD SOLN, 1000MG/L 100ML	\$	79.30
1417642-LM	NICKEL STD SOLN, 1000MG/L 100ML	\$	79.30
1417742	ZINC STD SOLN, 1000MG/L 100ML	\$	95.40
1417742-LM	ZINC STD SOLN, 1000MG/L 100ML	\$	95.40
1417869	MOLYVER 3 REAGENT PP 25ML PK/100	\$	154.10
1417869-LM	MOLYVER 3 REAGENT PP 25ML PK/100	\$	154.10
141832	HYDROCHLORIC ACID,2.5N 100ML MDB	\$	57.82
141832-LM	HYDROCHLORIC ACID,2.5N 100ML MDB	\$	57.82
1418642	MOLYBDENUM STD, 1000MG/L 100ML	\$	90.87
1418742	MOLYBDENUM STD SOLN,10MG/L 100ML	\$	82.92
1418828	KTO: CUVER 1, COPPER REAGENT PK/1000	\$	1,378.55
1418899	CUVER 1 COPPER PP 25ML PK/100	\$	184.29
1418899-LM	CUVER 1 COPPER PP 25ML PK/100	\$	184.29
1419301	TEST KIT, MO-2 MOLYBDATE, 0-40 mg/L	\$	790.15
1419301-LM	TEST KIT, MO-2 MOLYBDATE, 0-40 mg/L	\$	790.15
1419700	DROPPER, TUBE 0.5 & 1.0ML	\$	11.95
1419705	DROPPER, GLASS .5&1ML PK/5	\$	55.70
1420299	BUFFER PWD PLW, CITRATE PK/100	\$	263.79
1420299-LM	BUFFER PWD PLW, CITRATE PK/100	\$	263.79
1420416	PHOSPHATE STD SOLN, 10MG/L 946ML	\$	106.04
1420542	ACETATE STD SOLN, 1000MG/L 100ML	\$	122.86
1420703	aa MOLYBDOVANADATE REAGENT 2.9L	\$	419.17
1420703-LM	aa MOLYBDOVANADATE REAGENT 2.9L	\$	419.17
1420704	KTO: MOLYBDOVANADATE RGT, PK/4X946ML	\$	402.91
1420716	MOLYBDOVANADATE RGT-P04, 946ML	\$	115.44



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

1421199	CUVER 1 RGT PWD PLWS 5ML PK/100	\$	76.77
1421199-LM	CUVER 1 RGT PWD PLWS 5ML PK/100	\$	76.77
1421300	TEST KIT, CU-5 COPPER, CUVER I	\$	310.28
1421300-LM	TEST KIT, CU-5 COPPER, CUVER I	\$	310.28
1422249	COBALT CHLORIDE SOLN, 500ML	\$	196.94
1423188	INSTRUCTIONS, MICRO. PIPET	\$	25.05
1424210	PHOSPHATE 500MG/L 10ML AMP PK/16	\$	198.74
1424232	PHOSPHATE STD, 500MG/L 100ML MDB	\$	172.17
1424342	PHOSPHATE STD SOLN, 15MG/L 100ML	\$	110.02
1424610	ZINC, 25MG/L STD 10ML AMP PK/16	\$	216.81
1424610-LM	ZINC, 25MG/L STD 10ML AMP PK/16	\$	216.81
1425010	CHLORIDE, 12500 10ML AMP PK/16	\$	198.74
1425249	SULFATE STD SOLN, 2500MG/L 500ML	\$	182.48
1425249-LM	SULFATE STD SOLN, 2500MG/L 500ML	\$	182.48
1425310	IRON, 25MG/L STD 10ML AMP PK/16	\$	222.23
1425310-LM	IRON, 25MG/L STD 10ML AMP PK/16	\$	222.23
1425610	CHROMIUM 6, 12.5MG/L 10ML PK/16	\$	218.62
1426010	NITROGEN 500MG/L 10ML AMP PK/16	\$	204.16
1426810	CHLORINE, 50-75MG/L PK/16 10ML	\$	187.90
1426820	CHLORINE, 50-75MG/L 2ML PK/20	\$	173.62
1426901	SILICA GEL W/INDICATOR 454G	\$	191.52
1426901-LM	SILICA GEL W/INDICATOR 454G	\$	191.52
1427010	VOLATILE ACIDS, 10ML AMP PK/16	\$	374.00
1427110	DETERGENTS 60MG/L 10ML AMP PK/16	\$	202.36
1427810	ALKALINITY 0.500N 10ML AMP PK/16	\$	205.97
1428200	TWEEZERS, BALANCE WEIGHT	\$	31.80
1428300	PUMP, VACUUM HAND OPERATED	\$	447.11
1428400	DESSICATOR PLATE, 230MM	\$	312.21
1428500	DESSICATOR, 250MM	\$	350.75
1429099	ALUVER 3 POWDER PILLOWS PK/100	\$	285.47
1429099-LM	ALUVER 3 POWDER PILLOWS PK/100	\$	285.47
1429449	BLEACHING 3 RGT PWD PLWS PK/100	\$	83.47
1429449-LM	BLEACHING 3 RGT PWD PLWS PK/100	\$	83.47
1429525	SODIUM PYROPHOSPHATE, 10-HYD 50G	\$	91.23
1429600	ff FURNACE, MUFFLE 1093C 120V 50/60HZ	\$	4,875.82
1429900	DEMINERALIZER RESIN BOTTLE, 177ML	\$	73.81
1429900-LM	DEMINERALIZER RESIN BOTTLE, 177ML	\$	73.81
1430598	Buffered Dilution Water, (MgCl2/KH2PO4) 99mL,PK/25	\$	252.46
143201	nn TEST KIT, DE-2 DETERGENT SUBASSEM	\$	1,786.51
143201-LM	nn TEST KIT, DE-2 DETERGENT SUBASSEM	\$	1,786.51
143203	KTO: TEST KIT, DE-2 DETERGENTS	\$	1,189.08
143203-LM	KTO: TEST KIT, DE-2 DETERGENTS	\$	1,189.08
1432198	POTASSIUM 1 RGT PWD PLWS PK/25	\$	52.76
1432198-LM	POTASSIUM 1 RGT PWD PLWS PK/25	\$	52.76
1432298	POTASSIUM 2 RGT SOLN PLWS, PK/25	\$	46.06
1432298-LM	POTASSIUM 2 RGT SOLN PLWS, PK/25	\$	46.06



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

1432328-LM	KTO: POTASSIUM 3 RGT PWD PLWS PK/1000	\$	2,023.56
1432399	POTASSIUM 3 RGT PWD PLWS PK/100	\$	258.37
1432399-LM	POTASSIUM 3 RGT PWD PLWS PK/100	\$	258.37
143601	TEST KIT, CA-23 CARBON DIOXIDE	\$	174.22
143601-LM	TEST KIT, CA-23 CARBON DIOXIDE	\$	174.22
1436120	TUBE, CULT W/CAP 16X150MM PK/10	\$	123.15
1436369	DECHLORINATING REAGENT PP PK/100	\$	61.86
1436369-LM	DECHLORINATING REAGENT PP PK/100	\$	61.86
1436401	TITRAVER EDTA 0.0800M, CARTRIDGE	\$	70.10
1436401-LM	TITRAVER EDTA 0.0800M, CARTRIDGE	\$	70.10
1436703	aa PHOSPHATE STD SOLN, 30MG/L 2.9L	\$	126.64
1436716	PHOSPHATE STD SOLN, 30MG/L 946ML	\$	135.51
1436832	PHOSPHATE STD, 100MG/L 100ML MDB	\$	57.45
1437297	BAG, WHIRL-PAK* 5-1/2X9 PK/10	\$	26.98
1437701	SODIUM HYDROXIDE 0.1600N, CRTG	\$	76.05
1437701-LM	SODIUM HYDROXIDE 0.1600N, CRTG	\$	76.05
1437801	SODIUM HYDROXIDE 0.3636N, CRTG	\$	74.08
1437801-LM	SODIUM HYDROXIDE 0.3636N, CRTG	\$	74.08
1437901	SODIUM HYDROXIDE 1.600N, CRTG	\$	74.08
1437901-LM	SODIUM HYDROXIDE 1.600N, CRTG	\$	74.08
1438001	SODIUM HYDROXIDE 3.636N, CRTG	\$	74.08
1438001-LM	SODIUM HYDROXIDE 3.636N, CRTG	\$	74.08
143801	TEST KIT, CA-10WR DO PH DISC M	\$	599.36
1438101	SODIUM HYDROXIDE 8.00N, CRTG	\$	76.77
1438801	SULFURIC ACID 0.1600N, CARTRIDGE	\$	67.93
1438801-LM	SULFURIC ACID 0.1600N, CARTRIDGE	\$	67.93
1438901	SULFURIC ACID 1.600N, CARTRIDGE	\$	69.02
1438901-LM	SULFURIC ACID 1.600N, CARTRIDGE	\$	69.02
1439001	HYDROCHLORIC ACID 8.00N, CRTG	\$	77.33
1439101	SULFURIC ACID 8.00N, CARTRIDGE	\$	74.08
1439101-LM	SULFURIC ACID 8.00N, CARTRIDGE	\$	74.08
1439301	rr MERCURIC NITRATE 0.2256N, CRTG	\$	81.83
1439301-LM	rr MERCURIC NITRATE 0.2256N, CRTG	\$	81.83
1439601	SILVER NITRATE STD 0.2256N, CRTG	\$	76.77
1439601-LM	SILVER NITRATE STD 0.2256N, CRTG	\$	76.77
1439701	SILVER NITRATE STD 1.128N, CRTG	\$	74.08
1439701-LM	SILVER NITRATE STD 1.128N, CRTG	\$	74.08
1439901	TITRAVER EDTA 0.800M, CARTRIDGE	\$	70.10
1439901-LM	TITRAVER EDTA 0.800M, CARTRIDGE	\$	70.10
144001	TEST KIT, 8-P CHLORIDE LOW	\$	210.06
144001-LM	TEST KIT, 8-P CHLORIDE LOW	\$	210.06
1440026	SODIUM CHLORIDE, 1000 S/cm, 50ML	\$	35.22
1440042	SODIUM CHLORIDE, 491 MG/L 100ML	\$	58.18
1440049	SODIUM CHLORIDE, 1000uS/cm, 500ML	\$	146.88
1440053	SODIUM CHLORIDE, 491MG/L 1L	\$	137.84
1440101	SODIUM THIOSULFATE 2.00N, CRTG	\$	74.08



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

1440201	HEXAVER CDTA 0.0800M, CARTRIDGE	\$	74.97
1440301	HEXAVER CDTA 0.800M, CARTRIDGE	\$	67.02
1440301-LM	HEXAVER CDTA 0.800M, CARTRIDGE	\$	67.02
14411	HYDROGEN PEROXIDE 30%, ACS 473ML	\$	213.20
1442326	GALLIC ACID SOLUTION 50ML SCDB	\$	61.60
1442917	ACETONE, ACS 4L	\$	355.93
1442949	ACETONE, ACS 500ML	\$	135.51
1445817	CHLOROFORM, ACS 4L	\$	773.29
1445817-LM	CHLOROFORM, ACS 4L	\$	773.29
1445849	CHLOROFORM, ACS 500ML	\$	193.32
1445849-LM	CHLOROFORM, ACS 500ML	\$	193.32
1445853	CHLOROFORM, ACS 1L	\$	325.22
1445853-LM	CHLOROFORM, ACS 1L	\$	325.22
1445917	ISOPROPYL ALCOHOL, ACS 4L	\$	328.83
1445949	ISOPROPYL ALCOHOL, ACS 500ML	\$	205.97
1446417	METHYL ALCOHOL, ACS 4L	\$	294.50
1446449	METHYL ALCOHOL, ACS 500ML	\$	118.88
1446449-LM	METHYL ALCOHOL, ACS 500ML	\$	118.88
1446949	PYRIDINE, ACS 500ML	\$	316.18
1446949-LM	PYRIDINE, ACS 500ML	\$	316.18
1447017	TOLUENE, ACS 4L	\$	364.96
1447624	SILVER DIETHYLDITHIOCARBAMATE, 25G	\$	484.21
1447624-LM	SILVER DIETHYLDITHIOCARBAMATE, 25G	\$	484.21
1447849	HEXANES, ACS 500ML	\$	143.09
1448000	STOPPER, HOLLOW 11X17X26MM PK/6	\$	43.94
1448001	STOPPER, HOLLOW 14X20X26MM PK/6	\$	53.00
1449501	CARTRIDGE EMPTY, DGTL TITRATOR	\$	36.81
1451520	PIPET, VOLUMETRIC CLASS A 20ML	\$	90.39
1451534	PIPET, VOLUMETRIC CLASS A 0.5ML	\$	43.94
1451535	PIPET, VOLUMETRIC CLASS A 1ML	\$	43.94
1451536	PIPET, VOLUMETRIC CLASS A 2ML	\$	43.94
1451537	PIPET, VOLUMETRIC CLASS A 5ML	\$	43.94
1451538	PIPET, VOLUMETRIC CLASS A 10ML	\$	44.90
1451540	PIPET, VOLUMETRIC CLASS A 25ML	\$	82.48
1451541	PIPET, VOLUMETRIC CLASS A 50ML	\$	114.48
1451542	PIPET, VOLUMETRIC CLASS A 100ML	\$	138.18
1451550	PIPET, VOL CLASS A 20ML PK/6	\$	379.66
1451571	PIPET, VOL CLASS A 50ML PK/6	\$	533.83
145201	TEST KIT, HA-71A HARDNESS TOTAL	\$	252.46
145201-LM	TEST KIT, HA-71A HARDNESS TOTAL	\$	252.46
145300	TEST KIT, 5-B HARDNESS 0-30GPG	\$	75.93
145300-LM	TEST KIT, 5-B HARDNESS 0-30GPG	\$	75.93
1453300	ff EFFERVESCENT TABLET, QTY. 36	\$	44.45
1453300S	ff EFFERVESCENT TABLET, QTY 36, MIN LABELED	\$	45.89
1453700	FORCEPS, POINTED TIPS	\$	41.82
145400	TEST KIT, 5-EP HARDNESS 0-30GPG	\$	123.15



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

145400-LM	TEST KIT, 5-EP HARDNESS 0-30GPG	\$	123.15
145401	TEST KIT, 5-EP MG/L HARDNESS	\$	123.15
145401-LM	TEST KIT, 5-EP MG/L HARDNESS	\$	123.15
1454200	TEST KIT, CN-70 CHLORINE F&T	\$	323.77
1454200-LM	TEST KIT, CN-70 CHLORINE F&T	\$	323.77
1454201	TEST KIT, CN-70F CHLORINE FREE	\$	323.77
1454201-LM	TEST KIT, CN-70F CHLORINE FREE	\$	323.77
1454202	TEST KIT, CN-70T CHLORINE TOTAL	\$	337.26
1454599	ACID RGT PP (HR) 5ML PK/100	\$	92.67
1454599-LM	ACID RGT PP (HR) 5ML PK/100	\$	92.67
1454699	MOLYBDATE RGT PP 5ML PK/100	\$	79.30
1454699-LM	MOLYBDATE RGT PP 5ML PK/100	\$	79.30
1454828	CITRIC ACID PWD PW 25ML PK/1000	\$	834.72
1454828-LM	CITRIC ACID PWD PW 25ML PK/1000	\$	834.72
1454899	CITRIC ACID PP 25ML PK/100	\$	103.51
1454899-LM	CITRIC ACID PP 25ML PK/100	\$	103.51
1454999	CITRIC ACID PP 5ML PK/100	\$	58.71
1454999-LM	CITRIC ACID PP 5ML PK/100	\$	58.71
1455099	BROMPHENOL BLUE PP 50ML PK/100	\$	87.09
1455099-LM	BROMPHENOL BLUE PP 50ML PK/100	\$	87.09
1455232	BROMPHENOL BLUE IND, 100ML MDB	\$	56.90
1455232-LM	BROMPHENOL BLUE IND, 100ML MDB	\$	56.90
1455353	BOTTLE, MR AMB HDPE 28-410 1L	\$	4.05
1455400	TEST KIT, SI-5 SILICA HIGH 0-30	\$	491.44
1455400-LM	TEST KIT, SI-5 SILICA HIGH 0-30	\$	491.44
1455523	ff AMMONIA NITROGEN 1 SOLN, 25ML	\$	49.31
1456842	POTASSIUM IODIDE SOLN, 20% w/v 100ML	\$	136.95
1456842-LM	POTASSIUM IODIDE SOLN, 20% w/v 100ML	\$	136.95
1456942	STANNOUS CHLORIDE SOLN, 100ML	\$	148.32
1456942-LM	STANNOUS CHLORIDE SOLN, 100ML	\$	148.32
145700	TEST KIT, HA-4P HARDNESS CALCIUM	\$	339.19
145700-LM	TEST KIT, HA-4P HARDNESS CALCIUM	\$	339.19
145701	TEST KIT, HA-4P MG/L HARDNESS	\$	339.19
145701-LM	TEST KIT, HA-4P MG/L HARDNESS	\$	339.19
1457142	ARSENIC STD SOLN, 1000MG/L 100ML	\$	65.93
1457142-LM	ARSENIC STD SOLN, 1000MG/L 100ML	\$	65.93
1457439	FLASK, VOLUMETRIC CLASS A, 10ML	\$	108.50
1457440	FLASK, VOLUMETRIC CLASS A 25ML	\$	114.48
1457441	FLASK, VOLUMETRIC CLASS A 50ML	\$	122.18
1457442	FLASK, VOLUMETRIC CLASS A 100ML	\$	126.81
1457445	FLASK, VOLUMETRIC CLASS A 200ML	\$	145.31
1457446	FLASK, VOLUMETRIC CLASS A 250ML	\$	157.64
1457449	FLASK, VOLUMETRIC CLASS A 500ML	\$	192.14
1457453	FLASK, VOLUMETRIC CLASS A 1000ML	\$	217.77
1457799	ASCORBIC ACID PWD PLWS PK/100	\$	80.22
1458042	LEAD ACETATE SOLN, 10% 100ML	\$	87.09



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

1458042-LM	LEAD ACETATE SOLN, 10% 100ML	\$	87.09
1461142	BARIUM STD SOLN, 1000MG/L 100ML	\$	83.64
1461342	ff SILVER STD SOLN, 1000MG/L 100ML	\$	174.34
1462300	M-ENDO BROTH MF, DHYD 500G	\$	435.55
1462434	M-FC BROTH BASE, DHYD 500G	\$	491.44
146300	TEST KIT, IR-20 IRON MANGANESE	\$	689.94
146300-LM	TEST KIT, IR-20 IRON MANGANESE	\$	689.94
146400	TEST KIT, IR-18 IRON 0-4 mg/L	\$	260.17
146400-LM	TEST KIT, IR-18 IRON 0-4 mg/L	\$	260.17
146401	TEST KIT, IR-18B IRON 0-7 mg/L	\$	254.39
146401-LM	TEST KIT, IR-18B IRON 0-7 mg/L	\$	254.39
146430	TEST KIT, IR-10 IRON 0-10 MG/L	\$	235.12
146500	TEST KIT, IR-18A IRON 0-1 MGL	\$	291.01
146500-LM	TEST KIT, IR-18A IRON 0-1 MGL	\$	291.01
1465100	SAFETY BULB	\$	128.54
1466442	CHROMIUM STD, 1000MG/L 100ML	\$	83.64
1466442-LM	CHROMIUM STD, 1000MG/L 100ML	\$	83.64
146700	TEST KIT, MN-5 MANGANESE 0-3	\$	562.74
146700-LM	TEST KIT, MN-5 MANGANESE 0-3	\$	562.74
146803	TEST KIT, NI-11 NITRATE 0-40	\$	339.19
146803-LM	TEST KIT, NI-11 NITRATE 0-40	\$	339.19
146900	TEST KIT, OX-2P DISSOLVED OXYGEN	\$	294.86
146900-LM	TEST KIT, OX-2P DISSOLVED OXYGEN	\$	294.86
1469600	COUNTER, HAND TALLY	\$	71.31
147006	TEST KIT, 17-F PH 5.6-8.4	\$	335.33
147006-LM	TEST KIT, 17-F PH 5.6-8.4	\$	335.33
147008	TEST KIT, 17-H PH 6.6-8.4	\$	335.33
147009	TEST KIT, 17-J PH 7.4-9.6	\$	335.33
147011	TEST KIT, 17-N PH 4.0-10.0	\$	335.33
147011-LM	TEST KIT, 17-N PH 4.0-10.0	\$	335.33
1470332	ZINC ACETATE SOLN, 2N 100ML MDB	\$	129.89
1470953	BORATE BUFFER SOLN, 1000ML	\$	141.65
1471799	PETRI DISH W/PAD 9 x 50 mm, pk 100, (PALL)	\$	208.14
1472456	BOTTLE W/CAP, WM PP 250ML PK/12	\$	145.31
1472459	BOTTLE W/CAP, WM PP 500ML PK/12	\$	198.50
1472463	BOTTLE W/CAP, WM PP 1000ML PK/6	\$	161.11
1473632	AMMONIUM HYDROXIDE,10% 100ML MDB	\$	76.77
1473632-LM	AMMONIUM HYDROXIDE,10% 100ML MDB	\$	76.77
1473649	AMMONIUM HYDROXIDE SLN 10% 500ML	\$	189.71
1473885	WEIGHING PAPER, 100X100MM PK/500	\$	85.37
1474949	SODIUM STD SOLN, 1000MG/L 500ML	\$	101.35
147500	TEST KIT, PO-14 PHOSPHATE	\$	262.10
1476253	MAGNESIUM CHLORIDE SOLN, 51% W/V	\$	195.13
1476353	SODIUM HYDROXIDE, 0.250N 1000ML	\$	84.72
1476353-LM	SODIUM HYDROXIDE, 0.250N 1000ML	\$	84.72



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

1476453	SODIUM HYDROXIDE, 0.050N 1000ML	\$	125.39
1476932	PHOSPHORIC ACID, 10% 100ML MDB	\$	80.03
1477149	ff BOTTLE, REAGENT GL STPR 500ML	\$	126.81
1479010	POTASSIUM 250MG/L 10ML AMP PK/16	\$	209.58
1479110	NITROGEN, 50MG/L 10ML AMP PK/16	\$	204.16
1479120	NITROGEN, 50MG/L 2ML AMP PK/20	\$	175.98
1479442	MAGNESIUM STD, 1000MG/L 100ML	\$	97.20
1479442-LM	MAGNESIUM STD, 1000MG/L 100ML	\$	97.20
148002	TEST KIT, SU-5 SULFITE HIGH LOW	\$	383.51
148002-LM	TEST KIT, SU-5 SULFITE HIGH LOW	\$	383.51
1480407	ff STOPPER, NPRN no 7 12MM HOLE PK/6	\$	47.41
1480474	STOPPER W/HOLE NEOPRENE NO 1	\$	12.33
1480477	ff STOPPER, NEOPRENE no 8 7.9MM HOLE	\$	18.31
1480801	STOPPER, NEOPRENE SOLID #1 PK/12	\$	36.81
1480802	STOPPER, NEOPRENE SOLID #2 PK/12	\$	37.77
1480806	ff STOPPER, NEOPRENE SOLID no 6 PK/6	\$	46.25
1481000	TUBING CONN, QDISC 9.5-6.4 PK/12	\$	129.12
1481253	HYDROCHLORIC ACID, 0.1N 1000ML	\$	67.21
1481799	BORIC ACID PWD PLWS PK/100	\$	178.87
1481799-LM	BORIC ACID PWD PLWS PK/100	\$	178.87
1481866	COPPER SULFATE PWD PLWS, PK/50	\$	69.02
1481866-LM	COPPER SULFATE PWD PLWS, PK/50	\$	69.02
14832	METHYL ORANGE, 0.50G/L 100ML MDB	\$	40.83
1483531	BOILING CHIPS, CARBON 227G	\$	59.79
148410	TEST KIT, AD-16E WATER CONDTNG	\$	2,447.54
1484201	SODIUM HYDROXIDE, 0.9274N CRTG	\$	81.83
14849	METHYL ORANGE IND SOLN, 500ML	\$	65.57
1485200	DISH, PETRI DSPL F/ 47MM PK/500	\$	607.07
1485299	DISH, PETRI DSPL FOR 47MM PK/100	\$	170.36
1486166	NUTRIENT BFR SOLN PLWS 3ML PK/50	\$	128.45
1486166-LM	NUTRIENT BFR SOLN PLWS 3ML PK/50	\$	128.45
1486266	NUTRIENT BFR SOLN PLWS 6ML PK/50	\$	162.97
1486266-LM	NUTRIENT BFR SOLN PLWS 6ML PK/50	\$	162.97
1486398	NUTRIENT BFR SOLN PWS 19ML PK/25	\$	215.00
1486398-LM	NUTRIENT BFR SOLN PWS 19ML PK/25	\$	215.00
1486510	BOD STD, DILUTION 10ML AMP PK/16	\$	149.77
1486610	BOD STD, MANOMETRIC 10ML PK/16	\$	148.32
1486817	BOTTLE, CARBOY W/SPIGOT 4L	\$	323.77
1486858	BOTTLE, CARBOY W/SPIGOT 10L	\$	470.24
1486860	BOTTLE, CARBOY W/SPIGOT 20L	\$	572.38
1488600	RACK, TEST TUBE POLYETHYLENE	\$	37.77
1490817	SODIUM CITRATE, SILICA ANAL 3.78L	\$	191.52
1490849	SODIUM CITRATE FOR SILICA 500ML	\$	78.58
1490932	ACETATE BUFFER SOLN 100TESTS MDB	\$	66.13
1490932-LM	ACETATE BUFFER SOLN 100TESTS MDB	\$	66.13
1490953	ACETATE BUFFER SOLN, pH 4.0, 1 L	\$	229.46



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

1491800	PAD, ABSORBENT W/DISP PK/1000	\$	770.88
1495901	EDTA, 0.714M CRTG	\$	75.33
1495901-LM	EDTA, 0.714M CRTG	\$	75.33
1496001	EDTA, 0.1428M CRTG	\$	75.33
1496001-LM	EDTA, 0.1428M CRTG	\$	75.33
1496101	IODATE-IODIDE 0.3998N, CARTRIDGE	\$	67.93
1496101-LM	IODATE-IODIDE 0.3998N, CARTRIDGE	\$	67.93
1497054	VIAL, SPECIMEN W/CAP 24ML PK/10	\$	63.79
14G1005-203	BOM-ELE MDL,LEVEL,CAP SW,	\$	3,006.43
151065,08699	PTFE tubing 0.8x1.6 mm (per meter)	\$	158.22
1513700	DISC, WHITE OPAL GLASS	\$	112.16
151399,90002	TUBING, TYGON DIAM. = 1.6/3.2	\$	55.60
151399,90006	TUBING, TYGON DIAM. = 2.4/4.0	\$	126.42
151400,22387	TUYAU POLYETHYLENE 6X8 MM	\$	58.01
151400,22390	TUBING, POLYETHYLEN 12MM X 10M	\$	86.34
1523200	STIR BAR RETRIEVER, MAGNETIC PE	\$	55.70
15249	NITRIC ACID, ACS 500ML	\$	180.12
15249-LM	NITRIC ACID, ACS 500ML	\$	180.12
1531400	CAP. ELECTLT 1000 MFD 50 VDC	\$	20.24
15349	NITROGEN, AMMONIA 10MG/L 500ML	\$	81.47
1577800	O-RING NIT .103W X .490ID	\$	7.71
1590800	IC AMP,OP LOW POWER DUAL 358	\$	42.78
15934	BRILLIANT GRN BILE, DHYD 500G	\$	393.15
1601700	BRACKET, RELAY HOLDER	\$	52.61
1607400	JACK, POWER, PANEL MOUNTING	\$	38.16
1608000	IC REFERENCE VOLT 2.5V	\$	77.67
1608400	PLUG, PHONE 1/4 IN, 3 CONDUCTOR	\$	75.93
1608600	CLAMP, CABLE	\$	41.43
16232	PHENOLPHTHALEIN, 5G/L 100ML MDB	\$	48.42
16232-LM	PHENOLPHTHALEIN, 5G/L 100ML MDB	\$	48.42
16236	PHENOLPHTHALEIN, 5G/L 15ML SCDB	\$	44.45
1624100	FTG,ELBOW STREET 1/8 TO 1/8 NPT	\$	150.51
16253	PHENOLPHTHALEIN IND SOLN, 1000ML	\$	120.50
16253-LM	PHENOLPHTHALEIN IND SOLN, 1000ML	\$	120.50
1634700	O-RING NIT .070W X .551ID	\$	13.49
1636100	MEMBRANE, DO PROBE STD PK/15	\$	1,297.01
1636400	DIAPHRAGM KIT, DO PROBE	\$	635.98
1646800	CLAMP, WORM DRIVE HOSE, SS	\$	47.02
1650600	KNOB, ROUND, .250 SHAFT DIA	\$	25.63
1659400	WASHER SEALING STAINLESS STEEL	\$	9.25
16649H	POTASSIUM CHROMATE SOLN 500ML	\$	103.51
16649H-LM	POTASSIUM CHROMATE SOLN 500ML	\$	103.51
16701H	POTASSIUM IODIDE AR, 454G	\$	420.97
16726H	POTASSIUM IODIDE, ACS 100G	\$	173.98
16801H	POTASSIUM PERMANGANATE, 454G	\$	218.62
1685900	THERMISTOR 30K HONEYWELL	\$	141.65



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

1687800	FTG, BARB .062TBG X10-32TPI PLYP	\$	19.85
1690001	ee DIGITAL TITRATOR	\$	707.28
1690002	TITRATOR ASSY, DIGITAL	\$	724.63
1690008	MANUAL DGTL TITRATOR 16900-01	\$	70.34
1691600	O-RING BUTYL .103W X .206ID	\$	9.83
1695300	COUNTER, VEEDOMETER, DIGITAL TITRATOR	\$	158.03
1698500	LED,RED WITH 7.5 IN LEADS	\$	47.02
17001H	POTASSIUM PHOSPHATE MONO 454G	\$	179.40
170234	Di-Sodium-Hydrogen Phosphate, 500 g	\$	162.97
1704200	BOTTLE, SQUARE - W/ 25 ML MARK	\$	14.07
1704500	DROPPER, 2.50 ML	\$	39.31
1707400	BATTERY, ALKALINE C 1.5V	\$	36.04
17110	PHOSPHATE, 50MG/L 10ML AMP PK/16	\$	178.31
17120H	PHOSPHATE, 50MG/L 2ML AMP PK/20	\$	160.08
17149	PHOSPHATE STD SOLN, 50MG/L 500ML	\$	96.48
1716900	DIODE ZN 1N5347B 10V 5W 10%	\$	21.78
1720500	DELIVERY TUBE, 180 DEG HOOK PK/5	\$	37.58
172533	ROCHELLE SALT SOLUTION 29ML DB	\$	49.31
172549	ROCHELLE SALT SOLN, 500ML	\$	256.56
173006	TUBE, VIEWING 5-10ML-LP PK/6	\$	82.87
173106	STOPPER FOR 18MM TUBE PK/6	\$	31.22
173125	STOPPER FOR 18MM TUBE PK/25	\$	97.32
17314H	SILVER NITRATE, 113G	\$	1,244.85
17314H-LM	SILVER NITRATE, 113G	\$	1,244.85
173200	COLOR COMPARATOR BOX, HACH LOGO	\$	44.33
173250	COLOR COMPARATOR BOX, ACCUVAC	\$	109.46
1733	SCREW, THUMB, 1/4-20 X .5 BLACK CAP	\$	63.40
1734400	TUBE, PP .085 OD .025 ID 4.76	\$	7.52
1734433	TUBE, DELIVERY 4-3/4" STR PK/100	\$	44.90
173606	VIEWING TUBE W/5&20ML MARK PK/6	\$	119.49
1748000	O-RING 1.051 ID X .070W NITRILE PTFE	\$	11.18
1749400	SWITCH, SP3POS ROTARY SELECTOR,	\$	65.72
1757700	SAMPLE COOLER, STAINLESS STEEL	\$	4,606.01
177001	ROD, GLASS STIRRING PK/3	\$	48.95
178600	BULB, RUBBER	\$	12.72
1787500	FLOWMETER, AIR 0.2 TO 2.0 SCFH	\$	141.07
1787700	STOPPER ASSY	\$	71.69
179003	HYDRAVER 2, HYDRAZINE RGT 2.9L	\$	330.64
179032	HYDRAVER 2 REAGENT 100 TESTS MDB	\$	67.93
179032-LM	HYDRAVER 2 REAGENT 100 TESTS MDB	\$	67.93
179049	HYDRAVER 2 HYDRAZINE RGT, 500ML	\$	170.73
179049-LM	HYDRAVER 2 HYDRAZINE RGT, 500ML	\$	170.73
179053	HYDRAVER 2 HYDRAZINE RGT, 1000ML	\$	215.00
179053-LM	HYDRAVER 2 HYDRAZINE RGT, 1000ML	\$	215.00
179073	KTO: HYDRAVER 2 HYDRAZINE RGT PK/3X1L	\$	570.93
1790900	SCREW 10-32X1.875 PAN SST PHIL	\$	14.26
1800600	SWITCH, SNAP MOUNT ROCKER DPDT	\$	70.34

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

1801000	CORD, 18/3 SVT 7.5' 10A-125V	\$	60.90
180202	TEST KIT, AL-36B ECOLOGY COMBINE	\$	886.51
1810600	TUBING, LATEX 3/16ID 3/32W BLACK	\$	74.39
181133NC	FERROUS AMMONIUM SULFATE, 0.250N 29ML DB	\$	87.45
181233	FERROIN INDICATOR 29ML DB	\$	53.65
181233-LM	FERROIN INDICATOR 29ML DB	\$	53.65
181369	NITRIVER 2, NITRITE RGT PK/100	\$	107.68
181369-LM	NITRIVER 2, NITRITE RGT PK/100	\$	107.68
18149	SODIUM CARBONATE, 0.02N 500ML	\$	107.49
1816300	GREASE, SILICONE LUBRICANT	\$	93.66
181632	SULFIDE 1 REAGENT, 100ML MDB	\$	69.19
181632-LM	SULFIDE 1 REAGENT, 100ML MDB	\$	69.19
181649	SULFIDE 1 REAGENT, 500ML	\$	211.39
181649-LM	SULFIDE 1 REAGENT, 500ML	\$	211.39
181732	SULFIDE 2 REAGENT, 100ML MDB	\$	69.19
181732-LM	SULFIDE 2 REAGENT, 100ML MDB	\$	69.19
181749	SULFIDE 2 REAGENT, 500ML	\$	211.39
181749-LM	SULFIDE 2 REAGENT, 500ML	\$	211.39
18201H	SODIUM CHLORIDE, ACS 454G	\$	87.09
18201H-LM	SODIUM CHLORIDE, ACS 454G	\$	87.09
183400	TEST KIT, CH-8 CHROMIUM HEX	\$	260.17
18349	CHLORIDE STD SOLN 1000MG/L 500ML	\$	107.32
183699	PHENOL 2 RGT PWD PLWS PK/100	\$	139.65
183699-LM	PHENOL 2 RGT PWD PLWS PK/100	\$	139.65
183700	TEST KIT, HA-62 HDNS IRON 0-4	\$	614.78
183701	TEST KIT, HA-62A HDNS IRON 0-7	\$	614.78
183702	TEST KIT, HA62B HDNS IRON 0-7	\$	714.99
1838001	HEAT SHIELD, GLASS CYLINDRICAL	\$	316.06
1842100	GOGGLES, SAFETY CLEAR	\$	28.14
1843200	CABLE, 3 COND 18 AWG STRD SVT	\$	8.09
1847300	SCREW 1/4 -20X .75 HEX SST CAP	\$	5.78
184900	TEST KIT, HY-2 HYDRAZINE 0-1	\$	423.98
184900-LM	TEST KIT, HY-2 HYDRAZINE 0-1	\$	423.98
1850600	STRAINER ASSY, SAMPLE LINE	\$	653.32
185102	TEST KIT, CY-3 CYANURIC ACID	\$	354.60
185189	INSTRUCTIONS, CY-3 TEST KIT	\$	14.65
1856501	HEATER ELEMENT ASSY 120V	\$	562.74
1856502	HEATER ELEMENT ASSY 240V	\$	562.74
1857100	O-RING FLCAR .070W X .926ID	\$	21.97
1862400	BOTTLE, RND W/20&25ML MARK & CAP	\$	31.80
1863800	HEAT LIMITER 216 DEG C	\$	262.10
1864100	RACK, COD TEST TUBE COOLING	\$	310.28
186441	COPPER TUBE FOR SIPHON	\$	90.39
1870040	INSERT, MOLDED SATURATION EXTRACT	\$	80.75
18734	SODIUM HYDROXIDE, ACS 500G	\$	125.56
18765	ff SODIUM HYDROXIDE, ACS 2.5KG	\$	310.76
187701	THERMOMETER, POCKET -10 TO 110C	\$	79.02



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

187866	HEXAMETHYLENETETRAMINE PP PK/50	\$	99.90
187866-LM	HEXAMETHYLENETETRAMINE PP PK/50	\$	99.90
1882800	FRO POT TRM BOURNS 3299X-1-104	\$	33.34
188703	TEST KIT, STPL-WRT W/PH & TEMP	\$	2,813.71
188703-LM	TEST KIT, STPL-WRT W/PH & TEMP	\$	2,813.71
1889500	IC REFERENCE VOLT 1.23V	\$	57.24
189149	NITROGEN, AMMONIA 1MG/L 500ML	\$	81.47
1892300	SWITCH, SNAP MT RKR SPST ON/OFF	\$	46.25
1893300	TOOL, TRIMMER POT ADJUSTMENT	\$	124.50
189457	FILTER PAPER, FOLD 12.5CM PK/100	\$	198.50
189458	FILTER PAPER, FOLDED 15CM PK/100	\$	265.95
189459	FILTER PAPER, FOLD 18.5CM PK/100	\$	227.41
189460	FILTER PAPER, FOLDED 24CM PK/100	\$	265.95
189500	THERMOMETER, POCKET -30 TO 120F	\$	79.02
1895000	si LAMP ASSY, 1720 D/E	\$	233.19
189640	CYLINDER, GRAD MIX 25ML LOW FORM	\$	113.51
189641	CYLINDER, GRAD MIXING 50 ML	\$	141.26
189736	PHENOLPHTHALEIN SOLN 15ML SCDB	\$	39.19
189753	PHENOLPHTHALEIN SOLN, 1000ML	\$	102.62
1902000	JACK, BUSHING MOUNT, PHONE, MINI	\$	16.00
1902200	WASHER, FLAT, BAKELITE, .281 ID	\$	19.66
1906000	NIPPLE, OVERFLOW TUBE	\$	53.38
1906100	STAND ASSY, FLOW-THRU	\$	314.13
1906101	BASE, STAND	\$	66.10
1906102	UPRIGHT, STAND	\$	84.41
1906200	TUBE ASSY, OVERFLOW	\$	472.16
1906500	FTG, BARB .125 TBG X10-32 PDLYP	\$	16.19
1907100	FTG, PLUG HEX 10-32 PORT POLYPRO	\$	7.52
190901	STOPPER, BOD BOTTLE GLASS PENNYHEAD	\$	29.68
190902	BOTTLE, BOD 60ML W/30ML LINE	\$	97.32
1913400	SWITCH, DIP, 8 STATION SPST	\$	68.03
1914200	IC AMP, TRANSCONDUCTANCE 13600	\$	25.82
191442	BORON STD SOLN, 1000MG/L 100ML	\$	145.79
191442-LM	BORON STD SOLN, 1000MG/L 100ML	\$	145.79
19153	SODIUM HYDROXIDE, 0.100N 1000ML	\$	92.67
19153-LM	SODIUM HYDROXIDE, 0.100N 1000ML	\$	92.67
192049	irr MERCURIC NITRATE STD .141N 500ML	\$	155.02
192200	HYDROMETER, 1.000-1.225	\$	119.49
1922800	BUMPER, BUTTON, .375 DIA	\$	10.02
1924700	JACK, PHONE, PC BD MOUNT, 3 COND	\$	49.91
19253	SODIUM HYDROXIDE, 0.0227N 1000ML	\$	90.14
19253-LM	SODIUM HYDROXIDE, 0.0227N 1000ML	\$	90.14
192L0125	Model 531 Loop Powered Display	\$	1,468.53
192LX020	OXY 1100 Moisture Barrier Paste	\$	149.55
1930008	MANUAL, AMPEROMETRIC TITRATOR	\$	177.30
1930011	DIGITAL TITRATOR FORWARD TITRATN	\$	738.12
193317	AMMONIUM MOLYBDATE FOR ANLZR, 3.78L	\$	243.91

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

193332	AMMONIUM MOLYBDATE RGT 100ML MDB	\$	61.60
193332-LM	AMMONIUM MOLYBDATE RGT 100ML MDB	\$	61.60
193417	AMINO ACID REAGENT, 3.78L	\$	231.26
193417-LM	AMINO ACID REAGENT, 3.78L	\$	231.26
193432	AMINO ACID RGT, 100ML MDB	\$	58.18
193432-LM	AMINO ACID RGT, 100ML MDB	\$	58.18
193449	AMINO ACID REAGENT, 500ML	\$	106.04
193449-LM	AMINO ACID REAGENT, 500ML	\$	106.04
19353	SODIUM HYDROXIDE, 0.02N 1000ML	\$	83.11
193701	TEST KIT, TA-3 TANNIN LIGNIN	\$	585.87
193701-LM	TEST KIT, TA-3 TANNIN LIGNIN	\$	585.87
1938004	BATTERY, ALKALINE AA 1.5V PK/4	\$	56.27
1940088	MANUAL, TITRSTIR	\$	23.13
1940100	HOLDER ASSY, TITRATOR/ELECTRODE	\$	1,127.41
1940400	si BLOCK, CLAMP TITRSTIR	\$	510.71
1940700	SCREW, SHOULDER TITRSTIR	\$	166.70
1940900	KNOB, ROSETTE 1/4-20 X 1.25	\$	39.70
1942700	ROD ASSY	\$	547.32
19449	SILICA STD SOLN, 1000MG/L 500ML	\$	125.39
19449-LM	SILICA STD SOLN, 1000MG/L 500ML	\$	125.39
194749	NITROGEN, NITRATE 100MG/L 500ML	\$	77.50
19501H	SODIUM SULFITE, ANHYD ACS 454G	\$	114.36
1954900	PHOTOCELL, SILICON	\$	535.76
195555	RACK FOR ODI TEST	\$	360.39
1956400	WASHER #6 EXT LOCK SST	\$	7.71
1961800	MIRROR, SPHERICAL	\$	271.74
19634	LACTOSE BROTH, DEHYDRATED 500G	\$	466.38
1969900	DIODE MOTOROLA MPTE-5	\$	56.27
1970001	TENSETTE PIPET 0.1-1.0 ML	\$	1,073.45
1970002	TENSETTE PIPET 0.1-1.0ML REPL	\$	753.54
1970008	MANUAL, TENSETTE PIPET	\$	37.58
1970010	ee TENSETTE PIPET 1.0-10.0 ML	\$	1,073.45
1970088	MANUAL, TENSETTE PIPET	\$	42.40
1970200	ee TENSETTE PIPET FINAL ASSY 1-10ML	\$	1,073.45
1970700	TURRET, METERING .1 - 1 ML	\$	395.08
1970800	SCREW, TURRET ADJUSTMENT	\$	211.99
1970900	PISTON ASSY, TENSETTE 0.1-1.0ML	\$	240.90
1971300	BALL .125 DIA 440 STAINLESS	\$	13.49
1971400	SPRING,CPRSN .088 DIA X.500 MW	\$	20.43
1971500	NUT, RETAINER 1.0-10.0 ML	\$	29.87
1971600	PIN,DOWEL 0.0627 DIA. 0.750 LONG	\$	25.05
1971700	RING,RTNG EXTERNAL .156 SHAFT ST	\$	11.37
1971800	SPRING CPRSN .352 ODX 3.821 SST	\$	56.85
1971900	SEAL QUAD RING .250 NOM ID	\$	13.30
1972000	SCREW 6-32X .125SET BRS TIP	\$	21.78
1972400	WASHER, .258-.265 I.D. FLAT SST	\$	11.37
1972800	PISTON 1.0-10.0 ML	\$	82.29

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

1972900	COUPLING 1.0-10.0 ML	\$	53.00
1973000	TIP, EXTENSION 1.0-10.0 ML	\$	82.29
19734	LAURYL TRYPTOSE BROTH, DHYD 500G	\$	294.86
1983500	IC REGULATOR VOLT ADJ POS 317T	\$	34.11
1989800	TUBING, GLASS, PRECISION BORE	\$	200.43
1989900	SEAL QUAD RING .562 NOM ID	\$	26.60
1990000	O-RING NITRILE .097W X .755ID	\$	21.97
1990100	RING,RTNG WIRE .250 SHAFT CS	\$	12.14
1990600	HOUSING ASSY	\$	140.30
1990700	RETAINER, O-RING	\$	21.97
1990800	COUPLING, TENSETTE PIPET	\$	68.61
1990900	EXTENSION, HOUSING	\$	46.83
1993500	SAMPLE CELL, MATCHED PR DR/3000	\$	414.35
199503	aa MOLYBDATE 3 REAGENT, 2.9L	\$	388.45
199503-LM	aa MOLYBDATE 3 REAGENT, 2.9L	\$	388.45
199504	KTO: MOLYBDATE 3 REAGENT 1L PK/4	\$	457.11
199526	MOLYBDATE 3 RGT SOLN, 50ML SCDB	\$	65.57
199526-LM	MOLYBDATE 3 RGT SOLN, 50ML SCDB	\$	65.57
199532	MOLYBDATE 3 RGT SOLN, 100ML	\$	73.88
199532-LM	MOLYBDATE 3 RGT SOLN, 100ML	\$	73.88
199549	MOLYBDATE 3 RGT SOLN, 500ML	\$	102.62
199549-LM	MOLYBDATE 3 RGT SOLN, 500ML	\$	102.62
199553	MOLYBDATE 3 RGT SOLN, 1000ML	\$	126.28
199553-LM	MOLYBDATE 3 RGT SOLN, 1000ML	\$	126.28
199849	PAO STD SOLN, 0.000564N 500ML	\$	99.90
199901	PAO STD SOLN 0.00564N, CRTG	\$	70.63
199917	PAO STD SOLN, 0.00564N 3.78L	\$	303.53
199942	PAO STD SOLN, 0.00564N 100ML	\$	51.67
199953	PAO STD SOLN, 0.00564N 1000ML	\$	160.08
1999800	FILTER ASSY, INTERFERENCE 455NM	\$	1,628.48
1999900	FILTER ASSY, INTERFERENCE 860NM	\$	795.93
1D1044-104	PLASTIC METER FACE,4+1/2"	\$	194.65
1D1062	METER,2.5",0-1MADC,.2"	\$	1,032.98
1D1069	BLANK METER SCALE 3+1/2	\$	194.65
1D1072	MTR,ANALG,4-20ma,2%,3+1/2	\$	1,287.37
1E0565	ELECTRODE, PH,COMB. 30" W	\$	528.05
1E1095	ELECTRODE,PH COMB.PG13.5,	\$	1,135.12
1E1097	ELECTRODE,J GLASS W/TC, P	\$	1,667.03
1F1036	FUSE, 1/4 AMP,250V,SLOW	\$	31.41
1H1200	KNOB WITH STUD,3/8-16 X 2	\$	86.92
1H1201	KNOB WITH STUD 1/4-20 X 1	\$	40.66
1H1204	GROMMET 7/16"ID X 15/16"	\$	10.65
1P1109-327	HEX BUSH, CPVC, 80, S/T, 1-1/4x1/2	\$	65.14
1P1109-329	HEX BUSH, CPVC, 80, S/T, 1-1/4x1	\$	65.14
1P2101-107	PIPE/FT, PVC, 80,1+1/2	\$	3.66
1P2102-107	UNION, PVC, 80,THD,1+1/	\$	102.91
1P2104-107	TEE, PVC, 80,THD,1+1/	\$	240.90

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

1P2104-108	TEE, PVC, 80,THD, 2	\$	258.24
1P2109-238	HEX BUSH,PVC,80,SOC,2X1.5	\$	65.14
1P2116-107	END CAP, PVC, 80,THD,1+1/	\$	70.34
1P2120-105	COUPLING,PVC, 80,THD, 1	\$	33.34
1P2120-107	COUPLING,PVC,80,THD,1+1/2	\$	102.91
1P2120-132	Reducer Coupling, 1+1/2" x 1" PVC	\$	204.28
1W0905	CBL, MONITOR EXTEN, HART	\$	283.30
1W0979	CABLE 4 COND.	\$	8.29
1W0980	CABLE, 2 COND. + SHIELD	\$	10.41
1W1000	BUSS WIRE, 20 GAUGE	\$	11.37
1W1055	CABLE, 5 COND. WO/BLU WIR	\$	118.80
1W1098	CABLE,5 COND.+SHLD 1RR.24	\$	11.37
1W1100	si !! CABLE,5C.24GAPOLY150CSHLD	\$	15.03
1W1121	CABLE,PAIRED SHEILDED,	\$	10.21
1W3000	CABLE 3 PR XLPE OVERALL S	\$	11.37
1W5002	CABLE,10 COND,PVC W/2SHLD	\$	19.66
200096	KIT,MM TO IN 8810 W/O SAMPLE VALVE	\$	335.33
200097	KIT,MM TO IN, 8810 STD SAMPLE VALVE	\$	393.15
200-1002	SCREW, 4-40 X 1/2" PHMS,SS	\$	8.48
200117	mc KIT, 7013 PUMP	\$	3,391.87
200118	mc KIT, 7014 PUMP	\$	3,391.87
200119	mc KIT, 7015 PUMP	\$	3,391.87
200120	mc KIT, 7016 PUMP	\$	3,391.87
200138	KIT, 7024 PUMP, BARB REDUCED	\$	3,391.87
200140	KIT, 7024 PUMP, BARB REDUCED	\$	3,391.87
200142	KIT, 7015 PUMP, BARB REDUCED	\$	3,391.87
201002	TEST KIT, CYN-3 CYANIDE 0-0.3	\$	724.63
201002-LM	TEST KIT, CYN-3 CYANIDE 0-0.3	\$	724.63
201104	KTO: COPPER 1 REAGENT, 946ML PK/4	\$	355.93
201204	KTO: COPPER 2 REAGENT, 946ML PK/4	\$	359.54
202-0005	SET SCREW, 10-32 X 1/4"	\$	9.25
202-0012	*SET SCREW, 4-40 X 1/8" CUP PT,	\$	12.14
2022R0	SNSR-ORP,RYTON,5W,10'	\$	3,777.31
2022R060	ORP Sensor, Ryton, 60 Foot Cable	\$	5,338.34
202300	TEST KIT, HA-77 HARDNESS & IRON	\$	398.93
202300-LM	TEST KIT, HA-77 HARDNESS & IRON	\$	398.93
20232H	SULFURIC ACID, 0.100N 100ML MDB	\$	70.63
20253	SULFURIC ACID, 0.100N 1000ML	\$	67.93
20253-LM	SULFURIC ACID, 0.100N 1000ML	\$	67.93
2028G2000-A1A4X	SNSR-ORP,LCP,5W,40',NPT,T	\$	96.94
2028R0	Convertible ORP Sensor	\$	3,565.32
2028R010020000N	SNSR-ORP,5W,LCP,10',TNLM	\$	3,603.86
2028R010020000X	SNSR-ORP,5W,LCP,10',TNLM	\$	96.94
2028R012	SENSOR-ORP,5WIRE,LCP,12'	\$	4,586.74
2028R015	SENSOR-ORP,5WIRE,LCP,15'	\$	4,586.74
2028R020	ORP SNSR,LCP,CONVRTBLE,	\$	4,625.28
2028R025	SNSR-ORP,LCP,5W,25',CONV.	\$	3,777.31

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

2028R030	ORP SNSR,LCP,CONVRTBLE,	\$	4,740.91
2028R030000010N	ORP SNSR,LCP,CONVRTBLE,	\$	4,933.63
2028R035	ORP SNSR,LCP,CONV,35'CBL,	\$	4,818.00
2028R050	SENSOR-ORP,5WIRE,LCP,50'	\$	4,066.39
2028R075	SNSR-ORP,5W,LCP,75',STAND	\$	4,355.47
2028R1	SNSR-ORP,5W,LCP,GOLD,10'	\$	3,892.94
2028R120	SNSR-ORP,5W,LCP,GOLD,20'	\$	4,952.90
20342	SULFURIC ACID STD, 0.020N 100ML	\$	41.19
20353	SULFURIC ACID STD, 0.020N 1000ML	\$	70.63
20353-LM	SULFURIC ACID STD, 0.020N 1000ML	\$	70.63
2035400	KTO: 5500SC PHOSPHATE LR RGT SET	\$	1,468.89
2035400-LM	KTO: 5500SC PHOSPHATE LR RGT SET	\$	1,468.89
20356	SULFURIC ACID STD, 0.020N 4L	\$	156.27
2035600	KTO: 9610sc SILICA REAGENT SET	\$	1,468.89
2035601	KTO: 9610SC REAGENT SET SILICA, EU	\$	1,468.89
20356-LM	SULFURIC ACID STD, 0.020N 4L	\$	156.27
2035702	9610sc REAGENT 1, SILICA, 2,1-L INT	\$	623.33
2035702-LM	9610sc REAGENT 1, SILICA, 2,1-L INT	\$	623.33
2035802	9610sc REAGENT 2 SILICA 2L	\$	290.89
2035802-LM	9610sc REAGENT 2 SILICA 2L	\$	290.89
2035902	9610sc STANDARD 1 SILICA 2L	\$	180.31
2036002	KTO: 9610sc REAGENT 3 SILICA PKG	\$	527.57
2036002-LM	KTO: 9610sc REAGENT 3 SILICA PKG	\$	527.57
2036100	KTO: 9611SC PHOSPHATE LR RGT SET	\$	1,468.89
2036101	KTO: 9611SC PHOSPHATE LR RGT SET, EU	\$	1,468.89
2036202	9611sc REAGENT 1, LR PHOSPHATE, 2, 1-L INT	\$	374.00
2036302	9611SC REAGENT 2 LR PHOSPHATE 2L	\$	213.20
2036402	9611SC STANDARD 2 LR PHOSPHATE 2L	\$	292.69
2036502	KTO: 9611SC REAGENT 3 LR PHOSPHATE PKG	\$	574.55
2036602	9611SC STANDARD 1 LR PHOSPHATE 2L	\$	292.69
2036700	KTO: 9611SC PHOSPHATE HR RGT SET	\$	1,468.89
2036701	KTO: 9611SC PHOSPHATE HR RGT SET, EU	\$	1,468.89
2036802	9611sc REAGENT 1, HR PHOSPHATE, 2, 1-L INT	\$	771.48
2036902	9611SC REAGENT 2 HR PHOSPHATE 2L	\$	213.20
2037002	9611sc REAGENT 3, HR PHOSPHATE, 2, 1-L INT	\$	420.97
2037102	9611SC STANDARD 1 HR PHOSPHATE 2L	\$	204.16
2037502	9610sc REAGENT 4 SILICA 2L	\$	152.49
2037600	KIT;5 EMPTY BOTTLES,5500/9611 HLR P, 9610 SILICA	\$	383.51
2037601	KIT; 4 EMPTY BOTTLES,5500/9611 HR P, 5500 SILICA	\$	312.21
203832	SULFURIC ACID, 19.2N 100ML MDB	\$	76.77
203832-LM	SULFURIC ACID, 19.2N 100ML MDB	\$	76.77
203849	SULFURIC ACID STD, 19.2N 500ML	\$	105.88
203849-LM	SULFURIC ACID STD, 19.2N 500ML	\$	105.88
2038700	9610SC REAGENT 3 SILICA 1.8L	\$	205.97
2038700-LM	9610SC REAGENT 3 SILICA 1.8L	\$	205.97
2038800	9610SC REAGENT 3 SILICA POWDER	\$	334.25
2038800-LM	9610SC REAGENT 3 SILICA POWDER	\$	334.25

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

2038900	9611SC REAGENT 3 LR PHOSPHATE POWDER	\$	252.95
2039000	9611SC REAGENT 3 LR PHOSPHATE 1.8L	\$	419.17
203953	ETHYLENE GLYCOL 1000ML	\$	139.48
203953-LM	ETHYLENE GLYCOL 1000ML	\$	139.48
204053	SODIUM HYDROXIDE STD 4.5N 1000ML	\$	89.06
204053-LM	SODIUM HYDROXIDE STD 4.5N 1000ML	\$	89.06
204059	TOT INORG NITRGN BASE CONC 50ML	\$	69.02
204059-LM	TOT INORG NITRGN BASE CONC 50ML	\$	69.02
204253	FERRIC CHLORIDE-SULFURIC ACID	\$	92.14
204253-LM	FERRIC CHLORIDE-SULFURIC ACID	\$	92.14
20426	TITRAVER HARDNESS RGT, ACS 100G	\$	121.05
20426-LM	TITRAVER HARDNESS RGT, ACS 100G	\$	121.05
20434	TITRAVER HARDNESS RGT, ACS 500G	\$	276.43
20434-LM	TITRAVER HARDNESS RGT, ACS 500G	\$	276.43
204399	CHROMIUM 1 RGT PWD PLWS PK/100	\$	114.36
204399-LM	CHROMIUM 1 RGT PWD PLWS PK/100	\$	114.36
204499	CHROMIUM 2 RGT PWD PLWS PK/100	\$	114.36
204499-LM	CHROMIUM 2 RGT PWD PLWS PK/100	\$	114.36
204649	NITROGEN, NITRATE 1.0MG/L 500ML	\$	81.47
2048300	SILICA STD. SOLN, 5.0MG/L SIO2 2L	\$	-
204839-6	mc Accessory Kit,WATER,BASIC	\$	1,219.92
2048500	PHOSPHATE STD. SOLN,45.0 MG/L PO4 2L	\$	-
205-0008	5/16-18 HEX NUT, SS	\$	7.52
2052626	CHARCOAL ACTIVATED POWDER, 100G	\$	111.83
20529	EDTA STD SOLN, 0.0200N 200ML	\$	44.45
2052R1	ORP Sensor, Ryton Body, Au Electrode	\$	4,586.74
2054849	FERROUS AMM SULFATE 0.125N 500ML	\$	100.62
2054959	STIR BAR, MICRO 3MM DIA X 10MM	\$	44.90
2054960	STIR BAR, MICRO 2MM DIA X 7MM	\$	47.41
2055133	FERROIN IND SOLN 29ML DB	\$	52.76
20553	EDTA STD SOLN, 0.0200N 1000ML	\$	67.21
20556	EDTA STD SOLN, 0.0200N 4L	\$	146.35
2055734	BOILING CHIPS, SILICON CARBIDE	\$	193.32
2057799	SULFATE EXTRACTANT 1 PP PK/100	\$	87.09
2057799-LM	SULFATE EXTRACTANT 1 PP PK/100	\$	87.09
2058016	HARD STD SOLN 0.5 MG/L 946ML	\$	252.95
2058116	CALCIUM CHLORIDE 2MG/L 946ML	\$	179.95
2058R0	KTO: SNSR-ORP,5W,LCP	\$	3,989.30
2058R1	FINISH-SNSR,ORP,LCP,GOLD	\$	4,394.02
205932	FORMALDEHYDE, ACS 100ML MDB	\$	52.76
205932-LM	FORMALDEHYDE, ACS 100ML MDB	\$	52.76
2059703	aa PHOSPHATE STD SOLN, 3MG/L 2.9L	\$	170.73
2059716	PHOSPHATE STD SOLN, 3MG/L 946ML	\$	144.18
2059996	POTASSIUM IODIDE PWD PLWS PK/50	\$	76.77
2059996-LM	POTASSIUM IODIDE PWD PLWS PK/50	\$	76.77
2062501	MAGNESIUM CHLORIDE 0.0800M, CRTG	\$	74.08
2063300	TEST KIT, SU-DT SULFITE DIG TITR	\$	1,029.12

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

2063500	rr TEST KIT, CD-DT CHLORIDE	\$	1,065.74
2063500-LM	rr TEST KIT, CD-DT CHLORIDE	\$	1,065.74
2063600	TEST KIT, HA-DT HARDNESS	\$	1,002.14
2063600-LM	TEST KIT, HA-DT HARDNESS	\$	1,002.14
2063700	TEST KIT, AL-DT ALK W/DIG TITR	\$	1,015.63
2063700-LM	TEST KIT, AL-DT ALK W/DIG TITR	\$	1,015.63
2063900	TEST KIT, HAC-DT CAL/MAGNESIUM	\$	1,135.12
2063900-LM	TEST KIT, HAC-DT CAL/MAGNESIUM	\$	1,135.12
2064000	TEST KIT, AC-DT ACIDITY DIG TITR	\$	1,032.98
2064100	TEST KIT, CA-DT CO2 DIG TITRATOR	\$	988.65
2064342	CYLINDER, GRDTD LOW FORM 100ML	\$	145.70
2064400	TEST KIT, OZ2 OZONE	\$	352.68
2064400-LM	TEST KIT, OZ2 OZONE	\$	352.68
2065801	BLADE, REPLACEMENT F/CLIPPER	\$	35.65
206640	PIPET, SEROLOGICAL 25ML	\$	82.48
206670	PIPET, SEROLOGICAL 25ML PK/12	\$	873.02
206700	IMHOFF CONE, PC W/SCREW CLOSURE	\$	148.78
2067299	BROMTHYMOL BLUE IND PP PK/100	\$	90.14
206800	STAND, IMHOFF CONE FOR TWO CONES	\$	597.43
206801	STAND, IMHOFF CONE, HOLDS 4	\$	283.30
207-0004	1/4" FLAT WASHER,SS	\$	10.02
2070042	CHLORINE DIOXIDE RGT 1 100 TESTS	\$	56.18
2070142	CHLORINE DIOXIDE RGT 2 100 TESTS	\$	56.18
2070242	CHLORINE DIOXIDE RGT 3, 100ML	\$	56.18
2073549	ff AMMONIA RGT 1 PWD PLWS PK/100	\$	66.13
2073900	BURET KIT,AUTO SELF-ZEROING	\$	106.77
2074037	TUBING, CLEAR 6.4MMID 1.6W 50FT	\$	472.16
2074038	TUBING, CLEAR 3.2MMID 0.8W 50FT	\$	126.42
2074039	TUBING, CLEAR 9.5MMID 1.6W 50FT	\$	707.28
2074040	TUBING, CLEAR 12.7MMID 1.6W 50FT	\$	838.33
2074041	TUBING, CLEAR 7.9MMID 1.6W 50FT	\$	574.31
2074042	TUBING, CLEAR 3.2MMID 1.6W 50FT	\$	291.01
2074043	TUBING, CLEAR 4.8MMID 1.6W 50FT	\$	439.40
2074146	TUBING, VAC 9.5MMID 6.4W 50FT	\$	2,678.81
2074200	DEMINERALIZER CRTG, ORG REMOVAL	\$	368.10
2074846	PHENOL RED IND PWD PLWS PK/50	\$	45.70
2075325	ff BAG, 100 ML STERILE W/THIOSULFATE, PK/25	\$	43.55
2075333	BAG, STERILE W/THIOSULFATE PK/100	\$	123.15
2075333-LM	BAG, STERILE W/THIOSULFATE PK/100	\$	123.15
2075408	CARBOY, LOWBOY HDPE W/SPIGOT 8L	\$	466.38
2075460	CARBOY, LOWBOY HDPE W/SPIGOT 15L	\$	529.98
2075533	BAG, STERILE W/THIO, STANDUP, PK100	\$	132.78
2076004	KTO: MOLYBDOVANADATE RGT 1000ML PK/4	\$	419.17
2076017	MOLYBDOVANADATE RGT, 3.78L	\$	475.18
2076026	nn MOLYBDOVANADATE RGT, 25ML	\$	41.00
2076026-LM	nn MOLYBDOVANADATE RGT, 25ML	\$	41.00
2076032	MOLYBDOVANADATE RGT 100TESTS MDB	\$	56.90

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

2076032-LM	MOLYBDOVANADATE RGT 100TESTS MDB	\$	56.90
2076049	MOLYBDOVANADATE RGT, 500ML	\$	99.73
2076049-LM	MOLYBDOVANADATE RGT, 500ML	\$	99.73
2076053	MOLYBDOVANADATE RGT, 1000ML	\$	134.42
2076053-LM	MOLYBDOVANADATE RGT, 1000ML	\$	134.42
2076266	SULFAMIC ACID PP PK/50	\$	118.69
2076266-LM	SULFAMIC ACID PP PK/50	\$	118.69
2079100	DEMINERALIZER CRTG, ULTRAPURE HN	\$	470.24
2079200	DEMINERALIZER CRTG, HIGH CAP HN	\$	470.24
2079400	WALL BRACKET, DEMINERALIZER CRTG	\$	175.95
2080838-4	CABLE ASSY, 2 INCH	\$	86.34
2081335-1	WATER WEIR	\$	1,258.46
2081424-3	DWS	\$	8,286.96
2081599	BUFFER, CITRATE PP F/IRON PK/100	\$	160.08
2081599-LM	BUFFER, CITRATE PP F/IRON PK/100	\$	160.08
2081627-1	FLOW-CELL, COATING, QUART	\$	711.14
2081653-01	ASSY,P.S.,UNIV IN,NEMA US/JAPAN	\$	8,614.58
2081653-02	NEMA POWER SUPPLY ASY, UNIVERSAL ROW	\$	4,182.02
2081669	SULFOSALICYLIC ACID PP PK/100	\$	147.79
2081699GP	Brush W/SS HANDLE (LONG)	\$	319.92
2081701	TITRAVER, 0.0716M CRTG	\$	73.52
2081701-LM	TITRAVER, 0.0716M CRTG	\$	73.52
2081801	TITRAVER, 0.716M CRTG	\$	77.50
2081801-LM	TITRAVER, 0.716M CRTG	\$	77.50
2082127	ACCESSORY KIT, MOD. PCX,PCM	\$	262.10
2082411	si CABLE ASSY,5 PIN/5 PIN, W1	\$	13.49
2082412	CABLE ASSY, 7 PIN/7 PIN,W2	\$	89.04
2082414	CABLE ASSY, W5	\$	562.74
2082800	rr UV LAMP KIT, 115V	\$	1,985.02
2082802	rr UV LAMP KIT, 230V 50HZ	\$	2,139.19
2082803	rr UV LAMP KIT, W/ 230V EURO PLUG	\$	2,197.01
2082904	KTO: CHROMATE IND SOLN, 500ML PK/4	\$	350.51
2082949	CHROMATE INDICATOR SOLN, 500ML	\$	129.36
2083517-03	Cable, FEMALE SMB 16"	\$	427.84
2083838-01	Cable Assembly, FEMALE,Q-DISC,2FT	\$	418.20
2084766	POTASSIUM PERSULFATE PP PK/50	\$	63.60
2084766-LM	POTASSIUM PERSULFATE PP PK/50	\$	63.60
2084769	POTASSIUM PERSULFATE PP PK/100	\$	114.36
2084769-LM	POTASSIUM PERSULFATE PP PK/100	\$	114.36
2084900	SAMPLE CELL, LAB TURB PK/6	\$	188.29
2084953	Assembly, Pump 12VDC W/CONN	\$	4,683.10
2085859	BOTTLE W/CAP, AMB PE 500ML PK/6	\$	123.15
208601	TEST KIT, CD-51 CHLORIDE HIGH	\$	443.26
208601-LM	TEST KIT, CD-51 CHLORIDE HIGH	\$	443.26
2086047	BOTTLE, BLACK BOD	\$	84.60
2086256	ff BOTTLE, DROP DISP PE 250ML PK/6	\$	163.62
2086273	ff BOTTLE, DROP DISP PE 125ML PK/12	\$	231.26

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

2086280	BOTTLE, DROP DISP PE 60ML PK/12	\$	198.50
2086773	BOTTLE W/CAP, RECT 125ML PK/12	\$	88.46
2086776	BOTTLE W/CAP, RECT 250ML PK/12	\$	126.81
2086779	BOTTLE W/CAP, RECT 500ML PK/12	\$	198.50
2086783	BOTTLE W/CAP, RECT 1000ML PK/6	\$	229.34
2086784	BOTTLE W/CAP, RECT 2000ML PK/4	\$	252.46
2086871	BOTTLE W/CAP, SQ PE 60ML PK/12	\$	91.35
2086874	BOTTLE W/CAP, SQ PE 175ML PK/12	\$	139.14
2086879	BOTTLE W/CAP, SQ PE 500ML PK/12	\$	262.10
2086883	BOTTLE W/CAP, SQ PE 1000ML PK/6	\$	198.50
2086946	ff BOTTLE, WASH SAFETY 250ML	\$	56.85
2086949	ff BOTTLE, WASH SAFETY 500ML	\$	63.79
2087050	BOTTLE W/CAP, WM PE 30ML PK/12	\$	82.48
2087076	BOTTLE W/CAP, WM PE 250ML PK/12	\$	149.17
2087079	BOTTLE W/CAP, WM PE 500ML PK/12	\$	206.21
2087083	BOTTLE W/CAP, WM PE 1000ML PK/6	\$	185.97
2087117	BOTTLE, WM LARGE HDPE 4L	\$	95.01
2087154	BOTTLE, WM LARGE HDPE 2L	\$	63.79
2087300	BUCKET, W/GRADUATIONS 9.5L	\$	237.05
2087640	BURET, AUTO SELF-ZEROING 25ML	\$	383.51
2087742	BURNER, ALCOHOL 100ML	\$	141.26
2087801	CLAMP, BURET DOUBLE	\$	61.48
2087901	CLAMP, TUBING LARGE PK/12	\$	33.15
2088000	ALCONOX DETERGENT PWD #1104 1.8KG	\$	236.68
2088117	DETERGENT, LIQUI-NOX* SOLN 3.78L	\$	164.94
2088153	DETERGENT, LIQUI-NOX* SOLN .946L	\$	53.29
2088200	DETERGENT, ALCOTABS* PK/100	\$	79.30
2088300	CLEANER, ALNOCHROMIX* PK/10	\$	314.37
2088638	CYLINDER, GRAD MIX 10ML HDPE STPR	\$	113.51
2088640	CYLINDER, GRAD MIX 25ML TALL FORM	\$	119.49
2088641	CYLINDER, GRAD MIX 50ML TALL FORM	\$	146.47
2088642	CYLINDER, GRAD MIX 100ML TALL FORM	\$	165.55
2088701	DESICCANT, DRIERITE* W/IND 454G	\$	67.02
2088701-LM	DESICCANT, DRIERITE* W/IND 454G	\$	67.02
2088800	DESICCATOR, VACUUM	\$	529.98
2089000	RACK, DRYING GLASSWARE	\$	379.66
2089010	RACK, DRYING GLASSWARE, LARGE	\$	537.69
2089641	FLASK, ERLLENMEYER SCR CAP 50ML	\$	84.80
2089643	FLASK, ERLLENMEYER SCR CAP 125ML	\$	82.48
2089649	FLASK, ERLLENMEYER SCR CAP 500ML	\$	95.01
2089653	FLASK, ERLLENMEYER SCR CAP 1000ML	\$	104.07
2089673	FLASK, ERLLEN SCR CAP 125ML PK/6	\$	481.80
2089679	FLASK, ERLLEN SCR CAP 500ML PK/6	\$	514.56
2089743	FLASK, ERLLENMEYER PC 125ML	\$	55.70
2089746	FLASK, ERLLENMEYER PC 250ML	\$	79.02
2089843	FLASK, ERLLENMEYER PMP 125ML	\$	60.13
2089846	FLASK, ERLLENMEYER PMP 250ML	\$	66.49

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

2089873	FLASK, ERLN PMP 125ML PK/6	\$	304.50
2089876	FLASK, ERLN PMP 250ML PK/4	\$	254.39
2089879	FLASK, ERLN PMP 500ML PK/4	\$	383.51
2089953	FLASK, FILTERING PP 1000ML	\$	123.15
2091800	NOTEBOOK, FIELD	\$	356.53
2091900	NOTEBOOK, LABORATORY	\$	250.54
2092000	MARKER, LABORATORY BROAD LINE	\$	18.50
2092625	PIPET, SER STERILE 1ML PK/50	\$	120.45
2092626	PIPET, SER STERILE 2ML PK/50	\$	101.76
2092627	PIPET, SER STERILE 5ML PK/50	\$	129.12
2092628	PIPET, SER STERILE 10ML PK/50	\$	148.97
2092655	PIPET, SER STERILE 1ML PK/1000	\$	2,062.10
2092656	PIPET, SER STERILE 2ML PK/1000	\$	1,817.35
2092667	PIPET, SER STERILE 5ML PK/200	\$	487.58
2092668	PIPET, SER STERILE 10ML PK/200	\$	506.85
2092826	PIPET, SER STERILE 2ML PK/50	\$	87.30
2092828	PIPET, SER STERILE 10ML PK/50	\$	109.27
2092835	PIPET, SER STERILE 1ML PK/50	\$	71.50
2092855	PIPET, SER STERILE 1ML PK/1000	\$	1,385.66
2092868	PIPET, SER STERILE 10ML PK/500	\$	1,131.27
2093127	PIPET, SER STERILE 5X0.1ML PK/25	\$	116.60
2093128	PIPET, SER STERILE 10X.1ML PK/25	\$	137.22
2093135	PIPET, SER STERILE 1X.01ML PK/50	\$	88.65
2093136	PIPET, SER STERILE 2X.01ML PK/35	\$	66.49
2093158	PIPET, SER STERILE 10ML PK/400	\$	1,370.24
2093166	PIPET, SER STERILE 2ML PK/175	\$	429.77
2093438	PIPET, MOHR MEAS 10X0.1ML ORN	\$	43.94
2094100	PUMP, VACUUM HAND-OPR W/GAUGE	\$	537.69
2094200	BOD BOTTLE RACK 12 BOTTLE	\$	184.05
2094300	BOTTLE CARRIER, PINT	\$	346.90
2094700	RACK, PIPET	\$	119.49
2095000	SAMPLE CELL, MATCHED PR DR/2000	\$	472.16
2095100	SAMPLE CELL, 1CM MATCHED PAIR	\$	890.37
2095348	STIR BAR, OCTAGONAL 12.7X7.9MM	\$	18.31
2095349	STIR BAR, OCTAGONAL 15.9X7.9MM	\$	18.50
2095350	STIR BAR, OCTAGONAL 22.2X7.9MM	\$	22.55
2095351	STIR BAR, OCTAGONAL 25.4X7.9MM	\$	19.46
2095352	STIR BAR, OCTAGONAL 28.6X7.9MM	\$	20.62
2095355	STIR BAR, OCTAGONAL 50.8X7.9MM	\$	25.63
2095661	STOPCOCK PLUG, PTFE 2MM	\$	202.36
2096208	TUBE, CULT W/CAP 13X100MM PK/10	\$	124.50
2096230	TUBE, CULT W/CAP 20X150MM PK/10	\$	122.18
2096250	TUBE, CULT W/CAP 25X150MM PK/10	\$	141.65
2096900	WIPER, DISPOSABLE 11X22CM PK/280	\$	19.66
2097000	WIPER, DISPOSABLE 11X22CM PK/280	\$	20.62
209738	PIPET, DSPL WIDE-TIP 10-11ML	\$	10.60
2097810	WICK, FOR ALCOHOL BURNER PK/10	\$	61.86



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

209798	PIPET, STERILE 10-11ML PK/25	\$	87.30
2098061	CARBOY, RECT W/SPIGOT 9L	\$	470.24
2098080	CARBOY, RECT W/SPIGOT 20L	\$	597.43
2098640	CRUCIBLE, GOOCH 25 ML	\$	61.86
2098641	CRUCIBLE, GOOCH 40 ML, EA	\$	64.37
2099542	FLASK, VOLUMETRIC PMP 100ML	\$	106.38
2100300	SAMPLE CELL, 2100A PK/6	\$	208.14
2100803	aa SILICA STD SOLN 0.50MG/L 2.9L	\$	110.02
2100817	SILICA STD SOLN 0.50 MG/L 3.78L	\$	172.53
2100Q01	2100Q PORTABLE TURBIDIMETER	\$	4,875.82
2100Q01-LM	2100Q PORTABLE TURBIDIMETER	\$	4,875.82
2100Q01NB	2100Q PORTABLE TURBIDIMETER - NO BATTERY	\$	5,376.89
2100Q01USB	KTO: 2100Q w/USB + POWER MODULE	\$	5,723.78
2100Q01USB-LM	KTO: 2100Q w/USB + POWER MODULE	\$	5,723.78
2100QIS01	2100Q IS PORTABLE TURBIDIMETER	\$	4,875.82
2100QIS01-LM	2100Q IS PORTABLE TURBIDIMETER	\$	4,875.82
2100QIS01USB	KTO: 2100QIS W/USB + POWER MODULE	\$	5,723.78
2100QIS01USB-LM	KTO: 2100QIS W/USB + POWER MODULE	\$	5,723.78
2101100	WEIGHT SET, 100G	\$	196.57
2101415	LAURYL TRYPTOSE MPN TUBES PK/15	\$	116.21
2103200	BRUSH, CAMEL'S HAIR	\$	10.60
210542	SODIUM CHLORIDE, 1000MG/L 100ML	\$	51.84
2105503	KTO: DPD FREE CL PP PK/100X3	\$	263.79
2105503-LM	KTO: DPD FREE CL PP PK/100X3	\$	263.79
2105528	DPD FREE CHLORINE RGT 10 ML, PK/1000	\$	776.90
2105528-LM	DPD FREE CHLORINE RGT 10 ML, PK/1000	\$	776.90
210553	SODIUM CHLORIDE STD, 1000MG/L 1L	\$	138.76
2105560	DPD FREE REFILL VIAL	\$	202.36
2105560-LM	DPD FREE REFILL VIAL	\$	202.36
2105569	DPD FREE CHLORINE RGT PP, 10 ML, PK/100	\$	90.87
2105569-LM	DPD FREE CHLORINE RGT PP, 10 ML, PK/100	\$	90.87
2105603	KTO: DPD TOTAL CL PP PK/100X3	\$	263.79
2105603-LM	KTO: DPD TOTAL CL PP PK/100X3	\$	263.79
2105628	DPD TOTAL CHLORINE RGT 10 ML PK/1000	\$	776.90
2105628-LM	DPD TOTAL CHLORINE RGT 10 ML PK/1000	\$	776.90
2105660	DPD TOTAL REFILL VIAL	\$	202.36
2105660-LM	DPD TOTAL REFILL VIAL	\$	202.36
2105669	DPD TOT CHLORINE RGT PP 10 ML PK/100	\$	90.87
2105669-LM	DPD TOT CHLORINE RGT PP 10 ML PK/100	\$	90.87
2105728	KTO: FERROVER IRON RGT PP PK/1000	\$	715.47
2105728-LM	KTO: FERROVER IRON RGT PP PK/1000	\$	715.47
2105769	FERROVER, IRON REAGENT PK/100	\$	84.72
2105769-LM	FERROVER, IRON REAGENT PK/100	\$	84.72
2105828	KTO: CUVER 1, COPPER REAGENT PK/1000	\$	1,266.53
2105828-LM	KTO: CUVER 1, COPPER REAGENT PK/1000	\$	1,266.53
2105869	CUVER 1, COPPER REAGENT PK/100	\$	163.69
2105869-LM	CUVER 1, COPPER REAGENT PK/100	\$	163.69

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

2106028	KTO: PHOSVER 3 PHOSPHATE RGT PK/1000	\$	952.16
2106028-LM	KTO: PHOSVER 3 PHOSPHATE RGT PK/1000	\$	952.16
2106046	PHOSVER 3, PHOSPHATE PP PK/50	\$	61.60
2106046-LM	PHOSVER 3, PHOSPHATE PP PK/50	\$	61.60
2106069	ee PHOSVER 3 PWD PLWS 10ML PK/100	\$	118.69
2106069-LM	ee PHOSVER 3 PWD PLWS 10ML PK/100	\$	118.69
2106103	KTO: NITRAVER 5, NITRATE RGT PK/300	\$	419.17
2106103-LM	KTO: NITRAVER 5, NITRATE RGT PK/300	\$	419.17
2106169	NITRAVER 5 PWD PLWS 10ML PK/100	\$	150.13
2106169-LM	NITRAVER 5 PWD PLWS 10ML PK/100	\$	150.13
2106228	KTO: CITRIC ACID, FOR SILICA PK/1000	\$	690.18
2106228-LM	KTO: CITRIC ACID, FOR SILICA PK/1000	\$	690.18
2106269	CITRIC ACID, FOR SILICA PK/100	\$	80.22
2106269-LM	CITRIC ACID, FOR SILICA PK/100	\$	80.22
210637	PIPET, PP MOHR MEASURING 5ML	\$	208.14
210638	PIPET, PP MOHR MEASURING 10ML	\$	215.85
2106469	ff AMMONIA 1 RGT PP 10ML PK/100	\$	98.46
2106569	rr ff AMMONIA 2 RGT PP 10ML PK/100	\$	98.46
2106669	ZINCOVER 5, ZINC REAGENT PK/100	\$	176.51
2106669-LM	ZINCOVER 5, ZINC REAGENT PK/100	\$	176.51
210667	PIPET, PP GRAD MEAS 5ML PK/12	\$	2,409.00
210668	PIPET, PP GRAD MEAS 10ML PK/12	\$	2,466.82
2106769	SULFAVER 4, SULFATE RGT PK/100	\$	112.02
2106769-LM	SULFAVER 4, SULFATE RGT PK/100	\$	112.02
2106869	CYANIVER 3 CYANIDE RGT PP PK/100	\$	107.68
2106869-LM	CYANIVER 3 CYANIDE RGT PP PK/100	\$	107.68
2106969	CYANIVER 4 CYANIDE RGT PP PK/100	\$	107.68
2106969-LM	CYANIVER 4 CYANIDE RGT PP PK/100	\$	107.68
2107069	CYANIVER 5 CYANIDE RGT PP PK/100	\$	107.68
2107069-LM	CYANIVER 5 CYANIDE RGT PP PK/100	\$	107.68
2107169	NITRIVER 3 PWD PLW 10ML PK/100	\$	132.42
2107169-LM	NITRIVER 3 PWD PLW 10ML PK/100	\$	132.42
2107249	NITRAVER 6 RGT PP 10ML PK/100	\$	162.61
2107249-LM	NITRAVER 6 RGT PP 10ML PK/100	\$	162.61
2107328	KTO: MOLYBDATE RGT, SILICA HR PK/1000	\$	858.21
2107328-LM	KTO: MOLYBDATE RGT, SILICA HR PK/1000	\$	858.21
2107369	MOLYBDATE RGT, SILICA HR PK/100	\$	107.32
2107369-LM	MOLYBDATE RGT, SILICA HR PK/100	\$	107.32
2107428	KTO: ACID RGT, SILICA HR PK/1000	\$	888.92
2107428-LM	KTO: ACID RGT, SILICA HR PK/1000	\$	888.92
2107469	ACID RGT, SILICA HIGH RNG PK/100	\$	105.88
2107469-LM	ACID RGT, SILICA HIGH RNG PK/100	\$	105.88
2107569	NITRIVER 2, NITRITE RGT PK/100	\$	143.09
2107569-LM	NITRIVER 2, NITRITE RGT PK/100	\$	143.09
2107669	BUFFER, CITRATE TYPE PK/100	\$	141.10
2107669-LM	BUFFER, CITRATE TYPE PK/100	\$	141.10
2107769	SODIUM PERIODATE PP 10ML PK/100	\$	92.14



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

2107769-LM	SODIUM PERIODATE PP 10ML PK/100	\$	92.14
2109210	PHOSPHORUS, 25MG/L 10ML PK/16	\$	202.36
2110042	POLYVINYL ALCOHOL 20G/L 100TESTS	\$	57.98
2110042-LM	POLYVINYL ALCOHOL 20G/L 100TESTS	\$	57.98
2111132	Glycine, ACS, 100 g	\$	111.83
2112100	INOCULATING LOOP AND HANDLE	\$	46.25
211220	PHENOL SOLN 30 G/L 29ML	\$	48.23
211220-LM	PHENOL SOLN 30 G/L 29ML	\$	48.23
21126	PHENOL RED IND SOLN, 50ML SCDB	\$	46.06
2112820	MANGANESE, 25MG/L 2ML AMP PK/20	\$	180.12
21132	PHENOL RED IND SOLN, 100ML MDB	\$	53.29
2113300	KTO: TEST KIT, PN-10 115V	\$	2,755.90
2113300-LM	KTO: TEST KIT, PN-10 115V	\$	2,755.90
2113301	TEST KIT, PHOSPHONATE W/O UV	\$	572.38
2113301-LM	TEST KIT, PHOSPHONATE W/O UV	\$	572.38
2113302	KTO: TEST KIT, PN-10 PHOSPHONATE 230V	\$	2,755.90
2113400	GOGGLES, SPECTACLES UV CONTROL	\$	95.01
21136	PHENOL RED IND SOLN, 15ML SCDB	\$	46.06
2114500	CLAMP, 2-PRONG EXTENSION 38MM	\$	114.48
21149	PHENOL RED IND SOLN, 500ML	\$	112.91
2115700	SETTLOMETER KIT	\$	572.38
2115839	ff TUBE, CENTRIFUGE PC 15ML PK/10	\$	155.91
2117940	CYLINDER, GRADUATED PC 25ML	\$	79.02
2117941	CYLINDER, GRADUATED PC 50ML	\$	88.65
2117942	CYLINDER, GRADUATED PC 100ML	\$	98.29
2117946	CYLINDER, GRADUATED PC 250ML	\$	119.10
2117949	CYLINDER, GRADUATED PC 500ML	\$	151.09
2117953	CYLINDER, GRADUATED PC 1000ML	\$	183.85
211801	STOPPER, RUBBER SOLID #1 PK/12	\$	24.67
211802	ff STOPPER, RUBBER SOLID no 2 PK/12	\$	24.67
211806	STOPPER, RUBBER SOLID #6 PK/12	\$	44.90
2118669	COPPER RGT, FREE COPPER, 25ml pk/100	\$	79.30
2118869	HYDROSULFITE RGT PWD PLWS PK/100	\$	70.99
2118869-LM	HYDROSULFITE RGT PWD PLWS PK/100	\$	70.99
211907	STOPPER, RUBBER 1 HOLE #7 PK/6	\$	73.81
211908	STOPPER, RUBBER 1 HOLE #8 PK/6	\$	161.69
2119432	rr NESSLER RGT EX ALK, 100TESTS MDB	\$	89.78
2119432-LM	rr NESSLER RGT EX ALK, 100TESTS MDB	\$	89.78
2119449	rr NESSLER RGT EX ALK, 500ML	\$	173.98
2119449-LM	rr NESSLER RGT EX ALK, 500ML	\$	173.98
2119649	HYDROGEN PEROXIDE 50%, 490ML	\$	131.34
212132	SULFURIC ACID, 0.500N 100ML MDB	\$	66.85
212132-LM	SULFURIC ACID, 0.500N 100ML MDB	\$	66.85
2122326	ALKALINE CYANIDE RGT, 50ML SCDB	\$	64.68
2122326-LM	ALKALINE CYANIDE RGT, 50ML SCDB	\$	64.68
2122332	ALKALINE CYANIDE RGT 100ML MDB	\$	89.78
2122332-LM	ALKALINE CYANIDE RGT 100ML MDB	\$	89.78



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

2122349	ALKALINE CYANIDE REAGENT, 500ML	\$	238.49
2122349-LM	ALKALINE CYANIDE REAGENT, 500ML	\$	238.49
2122426	PAN IND SOLN, 0.1% 50ML SCDB	\$	61.60
2122426-LM	PAN IND SOLN, 0.1% 50ML SCDB	\$	61.60
2122432	PAN IND SOLN, 0.1% 100ML MDB	\$	91.23
2122432-LM	PAN IND SOLN, 0.1% 100ML MDB	\$	91.23
2122449	PAN INDICATOR SOLN, 0.1% 500ML	\$	243.91
2122449-LM	PAN INDICATOR SOLN, 0.1% 500ML	\$	243.91
2122531	SILICA STD SOLN, 25MG/L 236ML	\$	164.94
2122800	SAMPLE CELL, 25X54MM 10ML W/CAP PKGD	\$	28.14
212368	NICKEL 1 REAGENT PWD PLWS PK/25	\$	78.58
212468	NICKEL 2 REAGENT PWD PLWS PK/25	\$	114.55
212468-LM	NICKEL 2 REAGENT PWD PLWS PK/25	\$	114.55
2124720	DROPPER, LDPE 0.5&1.0ML PK/20	\$	23.51
212528	KTO: PHOSVER 3 PWD PLWS 25ML PK/1000	\$	1,089.47
212528-LM	KTO: PHOSVER 3 PWD PLWS 25ML PK/1000	\$	1,089.47
2125704	KTO: WATER, STANDARD ZERO, 946ML PK/4	\$	195.13
2125716	WATER, STANDARD ZERO 946ML	\$	52.76
2125815	rr COD DIGESTION VIAL, LR HW PK/150	\$	912.07
2125815-LM	rr COD DIGESTION VIAL, LR HW PK/150	\$	912.07
2125825	rr COD DIGESTION VIAL, LR HW PK/25	\$	166.80
2125825-LM	rr COD DIGESTION VIAL, LR HW PK/25	\$	166.80
2125915	rr COD DIGESTION VIAL, HR HW PK/150	\$	912.07
2125915-LM	rr COD DIGESTION VIAL, HR HW PK/150	\$	912.07
2125925	rr COD DIGESTION VIAL, HR HW PK/25	\$	166.80
2125925-LM	rr COD DIGESTION VIAL, HR HW PK/25	\$	166.80
212599	ee PHOSVER 3 PWD PLW 25ML PK/100	\$	137.67
212599-LM	ee PHOSVER 3 PWD PLW 25ML PK/100	\$	137.67
212699	ACID RGT PWD PLWS PK/100	\$	138.76
212699-LM	ACID RGT PWD PLWS PK/100	\$	138.76
2127716	CALCIUM CHLORIDE STD SOLN 946ML	\$	420.97
2128410	NITROGEN, 150MG/L 10ML AMP PK/16	\$	204.16
2128800	VIEWING TUBE INSERT, OPTICAL	\$	96.75
2129000	TEST KIT, CN-80 CHLORINE	\$	526.13
2129000-LM	TEST KIT, CN-80 CHLORINE	\$	526.13
21301H	UNIVER 3 POWDER, 454G	\$	184.29
21301H-LM	UNIVER 3 POWDER, 454G	\$	184.29
213100	ASPIRATOR, WITH CHECK VALVE	\$	168.82
213102	ASPIRATOR, W/O CHECK VALVE	\$	190.02
21320H	UNIVER 3 POWDER, 28.3G	\$	43.72
21320H-LM	UNIVER 3 POWDER, 28.3G	\$	43.72
2135100	REAGENT SET, PHOSVER 3 HHC AND ASCORBIC ACID	\$	109.13
2140800	STIR BAR, CUVETTE 9X8MM	\$	65.52
2141100	FORCEPS, MEMBRANE FILTER	\$	126.62
2141299	TRIAZOLE REAGENT PWD PLWS PK/100	\$	186.10
2141299-LM	TRIAZOLE REAGENT PWD PLWS PK/100	\$	186.10
2141342	BENZOTRIAZOLE, 500MG/L 100ML	\$	167.47

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

2141342-LM	BENZOTRIAZOLE, 500MG/L 100ML	\$	167.47
2142001	PAO STD SOLN, 0.0451N CARTRIDGE	\$	56.90
2142210	BOX, TAB-LOCK CORR 11.5X3.5X3.5	\$	4.24
2143166	BUFFER DILUTION WATER PLWS PK/25	\$	119.10
2150166	PHTHALATE-PHOSPHATE RGT PP PK/50	\$	99.54
2150166-LM	PHTHALATE-PHOSPHATE RGT PP PK/50	\$	99.54
2150232	PAN INDICATOR 0.3% 100ML MDB	\$	136.95
2150232-LM	PAN INDICATOR 0.3% 100ML MDB	\$	136.95
2150342	COBALT STD SOLN, 1000MG/L 100ML	\$	87.45
2150342-LM	COBALT STD SOLN, 1000MG/L 100ML	\$	87.45
2155332	BUFFER SOLN, PHOSPHATE 100ML MDB	\$	61.60
2155353	BUFFER SOLN, PHOSPHATE PH7.0 1L	\$	154.10
2158706	CAP, WHT PP 28-400 FOAM LNR PK/6	\$	16.96
2158806	CAP, WHT PP 28-410 FOAM LNR PK/6	\$	17.34
2162315	LAURYL TRYPTOSE 1X BROTH, PK/15	\$	113.13
2162700	BURNER, BUNSEN PROPANE	\$	105.61
2162800	TUBE, FERMENTATION 0.7ML PK/72	\$	66.10
2162921	ROSOLIC ACID, PWD 1G	\$	114.09
2163600	STIRRER, MAGNETIC BATT POWERED	\$	903.86
2164000	DISH, ALUMINUM 57 X16 MM PK/100	\$	151.09
2164100	FILTER FUNNEL, ACRYLIC PLATE 7CM	\$	932.76
2164200	FILTER PLATE, PTFE 7CM FOR 21641	\$	402.78
2164300	FILTER, GLASS FBR 70MM 2.7UM PK/100	\$	292.93
2164400	ff HAND PROTECTOR, HOT-HAND	\$	82.48
2165600	Algal Toxin Analyzer	\$	18,424.03
2165700	Algal Toxin Test Kit	\$	1,476.24
2166506	CAP, WHT PP 22-400 FOAM LNR PK/6	\$	16.00
2166706	CAP, WHT PP 24-400 FOAM LNR PK/6	\$	16.00
2167500	KTO: TEST KIT, TZ-1 115V	\$	2,717.35
2167502	rr KTO: TEST KIT, TZ-1 230/50 UV LAMP	\$	2,755.90
2167570	TEST KIT, TZ-1 W/O UV LAMP KIT	\$	516.49
2167969	DEHA REAGENT 1 PWD PLWS PK/100	\$	150.13
2168042	DEHA REAGENT 2 SOLN, 100ML	\$	69.19
2168042-LM	DEHA REAGENT 2 SOLN, 100ML	\$	69.19
2168049	DEHA REAGENT 2 SOLN, 500ML	\$	209.58
2168049-LM	DEHA REAGENT 2 SOLN, 500ML	\$	209.58
2169349	Octyl Alcohol, 99%, 500 mL	\$	108.93
2171820	FERRIC CHLORIDE RGT 5ML PK/20	\$	411.94
2171820-LM	FERRIC CHLORIDE RGT 5ML PK/20	\$	411.94
217240	CYLINDER, GRADUATED PMP 25ML	\$	68.99
217241	CYLINDER, GRADUATED PMP 50ML	\$	73.81
217242	CYLINDER, GRADUATED PMP 100ML	\$	88.65
217246	CYLINDER, GRADUATED PMP 250ML	\$	113.51
217249	CYLINDER, GRADUATED PMP 500ML	\$	167.67
217253	CYLINDER, GRADUATED PMP 1000ML	\$	190.21
217500	DROPPER, .583 ML	\$	18.31
2175749	SULFATE STD SOLN, 1000MG/L 500ML	\$	99.73

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

2177200	SECCHI DISC, 9 INCHES W/ ROPE	\$	254.39
2178800	GLOVES, HOT MILL 10-INCH LARGE	\$	175.95
2178900	DISH, PETRI 100X15MM DSPL PK/600	\$	728.48
2178901	DISH, PETRI 100X15MM DSPL PK/24	\$	36.81
217900	COVER FOR COOKIT	\$	43.94
2179000	WEIGHING BOAT, 8.9CM SQ PK/500	\$	406.64
2179020	WEIGHING BOAT, 8.9CM SQ PK/20	\$	52.81
2179100	WEIGHING BOAT, 14CM SQ PK/500	\$	784.37
218049	SODIUM HYDROXIDE SOLN 50% 500ML	\$	209.58
218049-LM	SODIUM HYDROXIDE SOLN 50% 500ML	\$	209.58
2182000	TEST KIT, NI-15 NITRITE LR	\$	339.19
2182000-LM	TEST KIT, NI-15 NITRITE LR	\$	339.19
2182115	LAURYL TRYPTOSE/MUG SOLN PK/15	\$	116.60
2182369	FREE COPPER RGT PWD PLWS PK/100	\$	84.72
2182369-LM	FREE COPPER RGT PWD PLWS PK/100	\$	84.72
2182999	FORMALDEHYDE RGT 1 PP PK/100	\$	118.52
2182999-LM	FORMALDEHYDE RGT 1 PP PK/100	\$	118.52
2183032	SULFURIC ACID, 1.90N 100ML MDB	\$	52.20
2183032-LM	SULFURIC ACID, 1.90N 100ML MDB	\$	52.20
2183100	TEST KIT, FORMALDEHYDE 0-0.5-10%	\$	323.77
2183100-LM	TEST KIT, FORMALDEHYDE 0-0.5-10%	\$	323.77
2183152	INSERT, MOLDED FORMALDEHYDE HR	\$	75.93
2183188	INSTRUCTIONS, FORMALDEHYDE HR	\$	14.65
2183668	GLYCOL TEST RGT 1 PWD PLWS PK/25	\$	31.08
2183668-LM	GLYCOL TEST RGT 1 PWD PLWS PK/25	\$	29.08
2183768	GLYCOL TEST RGT 2 PWD PLWS PK/25	\$	31.61
2183768-LM	GLYCOL TEST RGT 2 PWD PLWS PK/25	\$	31.61
2183825	GLYCOL TEST RGT 3 VIALS, PK/25	\$	256.56
2183825-LM	GLYCOL TEST RGT 3 VIALS, PK/25	\$	240.30
2183936	ETHYLENE GLYCOL, 0.05% 15ML SCDB	\$	42.28
2184300	LAMP, UV LONGWAVE 6W 115V 60HZ	\$	720.77
2184302	rr LAMP, UV LONGWAVE 230V 50/60HZ	\$	678.37
2184311	BULB, REPLACEMENT F/ULTRAVIOLET	\$	281.37
2184600	DEMINERALIZER RESIN ASSY, 473ML	\$	146.15
2184600-LM	DEMINERALIZER RESIN ASSY, 473ML	\$	146.15
2185336	THYMOLPHTHALEIN, 1G/L 15ML SCDB	\$	66.13
2185336-LM	THYMOLPHTHALEIN, 1G/L 15ML SCDB	\$	66.13
2185628	PIPET TIP, FOR 19700-01 PK/1000	\$	551.18
2185696	PIPET TIP, FOR 19700-01 PK/50	\$	48.76
2186400	TEST KIT, EG-1 ETHYLENE GLYCOL	\$	450.96
2186400-LM	TEST KIT, EG-1 ETHYLENE GLYCOL	\$	423.98
2186489	INSTRUCTIONS, EG-1 TEST KIT HACH	\$	36.81
218710	CALCIUM 10000MG/L 10ML AMP PK/16	\$	180.12
218710-LM	CALCIUM 10000MG/L 10ML AMP PK/16	\$	180.12
2187201	CAP ASSY, FLIP-SPOUT W/SCREENS	\$	19.27
2187399	COPPER MASKING RGT PP PK/100	\$	136.95
2187469	PORPHYRIN 1 PWD PLWS PK/100	\$	128.28

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

2187569	PORPHYRIN 2 PWD PLWS PK/100	\$	78.58
2187569-LM	PORPHYRIN 2 PWD PLWS PK/100	\$	78.58
2188226	CUVER 2, COPPER REAGENT PWD 100G	\$	359.54
2188299	CUVER 2 PWD PLWS 25ML PK/100	\$	204.16
2188299-LM	CUVER 2 PWD PLWS 25ML PK/100	\$	204.16
2190143	BOTTLE, WASH 125ML	\$	14.07
2190304	CALVER B, CONC & WTR 946ML PK/4	\$	1,555.61
2190304-LM	CALVER B, CONC & WTR 946ML PK/4	\$	1,555.61
219300	VIAL 2,5,10,15,20,25 ML	\$	16.96
2193432	METHYL PURPLE IND SOLN,100ML MDB	\$	60.35
2193449	METHYL PURPLE IND SOLN, 500ML	\$	118.52
2194100	TEST KIT, CU-6 TOT COPPER DISC	\$	300.64
2194100-LM	TEST KIT, CU-6 TOT COPPER DISC	\$	300.64
2194300	SAMPLER, WATER CORE 15 FOOT	\$	437.47
2194301	CASE, WATER CORE SAMPLER	\$	163.62
2194804	KTO: SULFURIC ACID, 2.50N 946ML PK/4	\$	216.81
2195900	MILLIPORE SDI TESTER W/FILTER PK	\$	4,220.57
2195900-LM	MILLIPORE SDI TESTER W/FILTER PK	\$	4,220.57
2196800	AMPULE BREAKER KIT	\$	84.80
2197153	POTASSIUM DICHROMATE, 1.000N 1L	\$	134.23
2197846	DPD FREE CHLORINE PP 5ML PK/50	\$	40.28
2197846-LM	DPD FREE CHLORINE PP 5ML PK/50	\$	40.28
2198246	DPD TOTAL CHLORINE PP 5ML PK/50	\$	40.28
2198246-LM	DPD TOTAL CHLORINE PP 5ML PK/50	\$	40.28
2198720	FILTER, GLASS FBR 20MM 1.5UM PK/100	\$	89.04
2198725	FILTER, GLASS FBR 25MM 1.5UM PK/100	\$	92.89
2198755	FILTER, GLASS FBR 90MM 1.2UM PK/100	\$	240.90
2199432	ROCHELLE SALT-PVA SOLN 100ML MDB	\$	44.45
2199449	ROCHELLE SALT-PVA REAGENT 500ML	\$	95.76
2199725	PIPET TIP, FOR 19700-10 PK/250	\$	170.75
2199796	PIPET TIP, FOR 19700-10 PK/50	\$	41.82
2202400	SYRINGE, 10CC LUER-LOK* TIP	\$	10.41
220399	SULFITE 1 RGT PWD PLWS PK/100	\$	113.63
220399-LM	SULFITE 1 RGT PWD PLWS PK/100	\$	113.63
220600	STOVE, FOLDING W/FUEL TABLETS	\$	43.94
220600-LM	STOVE, FOLDING W/FUEL TABLETS	\$	43.94
220700	HEATAB* DRY FUEL TABLET, PK/21	\$	35.22
220700-LM	HEATAB* DRY FUEL TABLET, PK/21	\$	35.22
220928	KTO: PHOSVER 3 PWD PLWS 5ML PK/1000	\$	785.94
220928-LM	KTO: PHOSVER 3 PWD PLWS 5ML PK/1000	\$	785.94
2209525	FILTER, MEMB 0.45UM 25MM PK/25	\$	108.12
220999	ee PHOSVER 3 PWD PLWS 5ML PK/100	\$	95.93
220999-LM	ee PHOSVER 3 PWD PLWS 5ML PK/100	\$	95.93
221120	BROMINE WATER, 30G/L 29ML	\$	46.42
221120-LM	BROMINE WATER, 30G/L 29ML	\$	46.42
221149	BROMINE WATER, 30G/L 500ML	\$	182.48
2212129	rr MERCURIC THIOCYANATE SOLN, 200ML	\$	164.41

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

2212129-LM	rr MERCURIC THIOCYANATE SOLN, 200ML	\$	164.41
2212242	FERRIC ION SOLUTION, 100ML	\$	73.88
2212242-LM	FERRIC ION SOLUTION, 100ML	\$	73.88
221500	COLIFORM TUBE RACK FOR 25 TUBES	\$	352.68
2215500	FILTER, GLASS FBR 24MM PK/100	\$	124.30
2217515	LAURYL TRYPTOSE/MUG PK/15	\$	116.60
221800	TUBE, GLASS F/DRAW OFF PIPETTE	\$	125.27
221969	NITRIVER 2 PP 25ML PK/100	\$	162.97
221969-LM	NITRIVER 2 PP 25ML PK/100	\$	162.97
2220038	ff DISPENSER W/FLASK, POUR-OUT 10ML	\$	620.56
222301	TEST KIT, AC-6 ACIDITY	\$	198.50
222301-LM	TEST KIT, AC-6 ACIDITY	\$	198.50
222338	INSTRUCTIONS, AC-6 TEST KIT	\$	16.38
222403	TEST KIT, AL-94 MG/L BOILER WTR	\$	946.26
2224553	HARDNESS BUFFER SOLN, 1,2&5MG/L 1L	\$	179.40
2224553-LM	HARDNESS BUFFER SOLN, 1,2&5MG/L 1L	\$	179.40
2224653	HARDNESS BUFFER SOLN, 10MG/L 1L	\$	186.10
2224653-LM	HARDNESS BUFFER SOLN, 10MG/L 1L	\$	186.10
2224753	HARDNESS BUFFER SOLN, 20MG/L 1L	\$	186.10
2225053	HARDNESS IND RGT SOLN, 2MG/L 1L	\$	122.30
2225053-LM	HARDNESS IND RGT SOLN, 2MG/L 1L	\$	122.30
2225153	HARDNESS IND RGT, 5-100MG/L 1L	\$	122.30
2225800	SYRINGE, 30CC LUER-LOK* TIP	\$	18.50
2226964	BUFFER, PH4.01 PWD PLWS R PK/250	\$	292.69
2226964-LM	ff BUFFER, PH4.01 PWD PLWS R PK/250	\$	292.69
2226966	BUFFER, PH4.01 PWD PLWS R PK/50	\$	70.63
2226995	BUFFER, PH4.01 PWD PLWS R PK/15	\$	37.75
2226995-LM	BUFFER, PH4.01 PWD PLWS R PK/15	\$	37.75
222702	TEST KIT, CH-14 CHROMATE	\$	260.17
2227064	BUFFER, PH7.00 PWD PLWS Y PK/250	\$	292.69
2227064-LM	BUFFER, PH7.00 PWD PLWS Y PK/250	\$	292.69
2227066	BUFFER, PH7.00 PWD PLWS Y PK/50	\$	70.63
2227066-LM	BUFFER, PH7.00 PWD PLWS Y PK/50	\$	70.63
2227095	BUFFER, PH7.00 PWD PLWS Y PK/15	\$	37.75
2227095-LM	BUFFER, PH7.00 PWD PLWS Y PK/15	\$	34.47
2227164	BUFFER, PH10.01 PWD PLW BLUE PK/250	\$	292.69
2227164-LM	BUFFER, PH10.01 PWD PLW BLUE PK/250	\$	292.69
2227166	BUFFER, PH10.01 PWD PLWS BLUE PK/50	\$	70.63
2227166-LM	BUFFER, PH10.01 PWD PLWS BLUE PK/50	\$	70.63
2227195	BUFFER, PH10.01 PWD PLWS BLUE PK/15	\$	37.75
2227195-LM	BUFFER, PH10.01 PWD PLWS BLUE PK/15	\$	37.75
222800	TEST KIT, CH-12 CHROMIUM TOTAL	\$	728.48
222800-LM	TEST KIT, CH-12 CHROMIUM TOTAL	\$	728.48
2228725	SYRINGE, 2ML ALL PP/PE PK/25	\$	53.77
2229653	BUFFER-MASKING SOLUTION, 1000ML	\$	136.95
2229653-LM	BUFFER-MASKING SOLUTION, 1000ML	\$	136.95
2229726	EGTA SOLUTION, 50ML SCDB	\$	70.99

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

2229726-LM	EGTA SOLUTION, 50ML SCDB	\$	70.99
223002	TEST KIT, CN-39WR WTR TRT COMB	\$	847.97
223002-LM	TEST KIT, CN-39WR WTR TRT COMB	\$	847.97
2230356	CALMAGITE SOLUTION, 0.003% 4L	\$	240.30
2230356-LM	CALMAGITE SOLUTION, 0.003% 4L	\$	240.30
223101	TEST KIT, CN-66 CL2 FREE & TOTAL	\$	235.12
223101-LM	TEST KIT, CN-66 CL2 FREE & TOTAL	\$	235.12
223102	TEST KIT, CN-66F FREE CHLORINE	\$	235.12
223102-LM	TEST KIT, CN-66F FREE CHLORINE	\$	235.12
223103	TEST KIT, CN-66T TOTAL CHLORINE	\$	235.12
223103-LM	TEST KIT, CN-66T TOTAL CHLORINE	\$	235.12
2232512	PIPET, TRANSFER STERILE PK/15	\$	21.78
2233100	BAG, WHIRL-PAK* STER 3X7 PK/500	\$	238.97
2233125	BAG, WHIRL-PAK* 3X7 207ML PK/25	\$	23.70
2233199	BAG, WHIRL-PAK* 3X7 207ML PK/100	\$	90.39
2233200	CRUCIBLE HOLDER, WALTER	\$	55.70
223400	TEST KIT, CO-1 COLOR 0-500	\$	323.77
223604	KTO: MOLYBDATE REAGENT, 1000ML PK/4	\$	366.77
223617	MOLYBDATE RGT F CR 2 PHOS, 3.78L	\$	377.61
223632	MOLYBDATE REAGENT, 100ML MDB	\$	59.98
223632-LM	MOLYBDATE REAGENT, 100ML MDB	\$	59.98
223653	MOLYBDATE REAGENT, 1000ML	\$	109.66
223653-LM	MOLYBDATE REAGENT, 1000ML	\$	109.66
2236945	ff WATER IN OIL RGT, TNT 4.5ML PK50	\$	581.77
2237400	SYRINGE, 20ML ECCENTRIC LUER TIP	\$	21.58
223801	TEST KIT, HS-WR HYDROGEN SULFIDE	\$	784.37
223801-LM	TEST KIT, HS-WR HYDROGEN SULFIDE	\$	784.37
2238300	DESICCATOR, VACUUM W/140MM PLATE	\$	227.41
223840	INSERT, MOLDED HS-6 TEST KIT	\$	45.10
2238600	SUPPORT, SEPARATORY FUNNEL PP	\$	114.48
2239514	EYE WASH BOTTLE, 120ML	\$	43.94
2239800	HANDBOOK, NIOSH/OSHA CHEM HAZ	\$	196.57
224000	TEST KIT, NI-6 NITRITE HR	\$	466.38
224000-LM	TEST KIT, NI-6 NITRITE HR	\$	466.38
2240349	CALCIUM STD SOLN, 1000MG/L 500ML	\$	105.51
224100	rr TEST KIT, NI-8 AMMONIA NITROGEN	\$	339.19
224100-LM	rr TEST KIT, NI-8 AMMONIA NITROGEN	\$	339.19
2241106	CAP, WHITE PP 15-415 PK/6	\$	13.88
2241125	CAP, WHITE PP 15-415 PK/25	\$	31.80
2241646	LAB COAT, WHITE POPLIN SIZE 46	\$	237.05
2241732	ALKALI SOLUTION, 100ML MDB	\$	66.13
2241732-LM	ALKALI SOLUTION, 100ML MDB	\$	66.13
2241749	ALKALI SOLN FOR CA/MG, 500ML	\$	170.00
2241749-LM	ALKALI SOLN FOR CA/MG, 500ML	\$	170.00
2241832	CALCIUM/MAGNESIUM IND, 100ML MDB	\$	74.08
2241832-LM	CALCIUM/MAGNESIUM IND, 100ML MDB	\$	74.08
2241849	CALCIUM/MAGNESIUM IND, 500ML	\$	189.71



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

2241849-LM	CALCIUM/MAGNESIUM IND, 500ML	\$	189.71
2241926	EDTA SOLN, 1M 50ML SCDB	\$	49.67
2241926-LM	EDTA SOLN, 1M 50ML SCDB	\$	49.67
2242000	RGT SET, ALUMINUM ALUMINON	\$	440.85
2242000-LM	RGT SET, ALUMINUM ALUMINON	\$	440.85
2242200	KTO: REAGENT SET, CADMIUM 100 TESTS	\$	2,168.10
2242300	REAGENT SET, CHLORINE DIOXD	\$	254.75
2242500	REAGENT SET, TOTAL CHROMIUM	\$	502.28
2242500-LM	REAGENT SET, TOTAL CHROMIUM	\$	502.28
2242600	REAGENT SET, COBALT OR NICKEL	\$	847.37
2242600-LM	REAGENT SET, COBALT OR NICKEL	\$	847.37
2242700	REAGENT SET, COPPER 25ML	\$	518.54
2242700-LM	REAGENT SET, COPPER 25ML	\$	518.54
2242800	REAGENT SET, CYANIDE 25ML	\$	381.22
2242800-LM	REAGENT SET, CYANIDE 25ML	\$	381.22
2243100	KTO: REAGENT SET, LEAD 90 TESTS	\$	2,402.98
2243101	KTO: RGT SET, LEAD DITHIZONE METHOD	\$	3,197.95
2243200	REAGENT SET, MANGANESE HR 25ML	\$	243.91
2243300	REAGENT SET, MANGANESE LR 25ML	\$	500.47
2243300-LM	KTO: REAGENT SET, MANGANESE LR 25ML	\$	500.47
2243400	REAGENT SET, MOLYBDENUM/DATE	\$	444.46
2243400-LM	REAGENT SET, MOLYBDENUM/DATE	\$	444.46
2243500	KTO: REAGENT SET, FOR NICKEL	\$	930.48
2243500-LM	KTO: REAGENT SET, FOR NICKEL	\$	930.48
2243700	REAGENT SET, NITROGEN AMNA	\$	1,114.76
2243700-LM	REAGENT SET, NITROGEN AMNA	\$	1,114.76
2243900	KTO: REAGENT SET, FOR PHENOLS	\$	2,529.45
2243900-LM	KTO: REAGENT SET, FOR PHENOLS	\$	2,529.45
2243901	KTO: RGT SET, PHENOLS 0-2MG/L	\$	3,378.62
2243901-LM	KTO: RGT SET, PHENOLS 0-2MG/L	\$	3,378.62
2244000	REAGENT SET, PHOSPHONATES	\$	377.61
2244000-LM	REAGENT SET, PHOSPHONATES	\$	377.61
2244100	REAGENT SET, PHOSPHORUS REAC	\$	114.36
2244100-LM	REAGENT SET, PHOSPHORUS REAC	\$	114.36
2244300	REAGENT SET 25 ML, SILICA HR/UHR	\$	364.96
2244300-LM	REAGENT SET 25 ML, SILICA HR/UHR	\$	364.96
2244500	REAGENT SET, FOR SULFIDE	\$	193.32
2244500-LM	REAGENT SET, FOR SULFIDE	\$	193.32
2244501	REAGENT SET, HYDROGEN SULFIDE	\$	198.74
2244600	REAGENT SET, TANNIN & LIGNIN	\$	317.99
2244600-LM	REAGENT SET, TANNIN & LIGNIN	\$	317.99
2244700	KTO: REAGENT SET, VOLATILE ACIDS	\$	464.33
2244700-LM	KTO: REAGENT SET, VOLATILE ACIDS	\$	464.33
2244800	RGT SET, ZINC TEST KIT	\$	285.47
2244800-LM	RGT SET, ZINC TEST KIT	\$	285.47
2245100	REAGENT SET, NITROGEN NITRATE LR	\$	354.12
224800	TEST KIT, PO-19 PHOSPHATE HIGH	\$	389.29



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

224800-LM	TEST KIT, PO-19 PHOSPHATE HIGH	\$	389.29
224801	TEST KIT, PO-19A PO4 HIGH TURBID	\$	529.98
224801-LM	TEST KIT, PO-19A PO4 HIGH TURBID	\$	529.98
224902	TEST KIT, PO-23 PHOSPHATE	\$	537.69
224902-LM	TEST KIT, PO-23 PHOSPHATE	\$	537.69
225001	TEST KIT, PO-24 PHOSPHATE TOTAL	\$	757.39
225001-LM	TEST KIT, PO-24 PHOSPHATE TOTAL	\$	757.39
2250800	CAP, WHITE PP 15-415WT COCT	\$	1.73
225100	TEST KIT, SF-1 SULFATE	\$	350.75
225100-LM	TEST KIT, SF-1 SULFATE	\$	350.75
225188	INSTRUCTIONS, SF-1 TEST KIT	\$	16.38
2251926	TKN INDICATOR SOLN, 50ML SCDB	\$	47.51
2252102	COLONY COUNTER, DARKFIELD 220V	\$	6,571.75
2253534	TRYPTIC SOY BROTH, DHYD 500G	\$	163.43
2253828	KTO: AMINO ACID F PP 25ML PK/1000	\$	691.99
2253828-LM	KTO: AMINO ACID F PP 25ML PK/1000	\$	691.99
2253869	AMINO ACID F PP 25ML PK/100	\$	86.72
2253869-LM	AMINO ACID F PP 25ML PK/100	\$	86.72
2253929	COD STD SOLN, 1000MG/L 200ML	\$	70.10
2253929-LM	COD STD SOLN, 1000MG/L 200ML	\$	70.10
225401	TEST KIT, CN-65 CHLORINE	\$	335.33
225401-LM	TEST KIT, CN-65 CHLORINE	\$	335.33
2254069	AMINO ACID F, PP FOR 10ML PK/100	\$	74.97
2254069-LM	AMINO ACID F, PP FOR 10ML PK/100	\$	74.97
2254117	KTO: AMINO ACID F REAGENT PACKAGE	\$	460.72
2254117-LM	KTO: AMINO ACID F REAGENT PACKAGE	\$	460.72
2254217	CITRIC ACID RGT SOLN, 3.78L MODIFIED	\$	179.23
2254217-LM	CITRIC ACID RGT SOLN, 3.78L MODIFIED	\$	179.23
2254232	CITRIC ACID F RGT SOLN 100ML MDB	\$	69.19
2254232-LM	CITRIC ACID F RGT SOLN 100ML MDB	\$	69.19
2254249	CITRIC ACID RGT SOLN, 500ML	\$	109.30
2254249-LM	CITRIC ACID RGT SOLN, 500ML	\$	109.30
2254400	FILTER HOLDER, WITH 250ML RCVR	\$	529.98
2255000	TEST KIT, SI-7 SILICA LR 0-1MG/L	\$	491.44
2255000-LM	TEST KIT, SI-7 SILICA LR 0-1MG/L	\$	491.44
2255453	HARDNESS INDICATOR RGT, 1MG/L	\$	122.30
2255653	CALIB STD SOLN, DR CELL 1000ML	\$	170.73
2255734	LAURYL TRYPTOSE-MUG BROTH, 500G	\$	1,416.49
2256900	REAGENT SET, SILICA 100 TESTS	\$	229.46
2256900-LM	REAGENT SET, SILICA 100 TESTS	\$	229.46
2257169	MBTH PWD PLWS PK/100	\$	292.69
2257169-LM	MBTH PWD PLWS PK/100	\$	292.69
2257249	DEVELOPING SOLN, LR FORMALDEHYDE	\$	95.76
2257249-LM	DEVELOPING SOLN, LR FORMALDEHYDE	\$	95.76
2257310	FORMALDEHYDE 4000MG/L 10ML AMP PK/16	\$	215.00
2257310-LM	FORMALDEHYDE 4000MG/L 10ML AMP PK/16	\$	215.00
2257700	REAGENT SET, FORMALDEHYDE LR	\$	661.27

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

2257700-LM	REAGENT SET, FORMALDEHYDE LR	\$	661.27
2258700	SYRINGE, 60CC LUER-LOK* TIP	\$	42.40
2259901	PAO STD SOLN, 0.00451N CRTG	\$	56.90
226=004=013	TEFLON COATED MAG. STIRRER 4.5X15MM	\$	48.76
2260549	PHOSPHORUS 1 RGT SOLN, 500ML	\$	113.46
2260549-LM	PHOSPHORUS 1 RGT SOLN, 500ML	\$	113.46
2260649	PHOSPHORUS 2 RGT SOLN, 500ML	\$	218.62
2260649-LM	PHOSPHORUS 2 RGT SOLN, 500ML	\$	218.62
226142	CARBON DIOXIDE STD 100MG/L 100ML	\$	172.17
2261432	BURET, AUTO 10ML TFE STPCK ONLY	\$	317.99
2261438	BURET, AUTO 10ML TFE STPCK COMPL	\$	358.46
2261440	BURET, AUTO 25ML TFE STPCK COMPL	\$	393.15
2261442	BURET, AUTO 25ML TFE STPCK ONLY	\$	341.11
2261451	STOPCOCK PLUG, PTFE 2MM REPL	\$	81.91
2263356	FILTER PAPER, FOLDED 11CM PK/500	\$	516.49
2263411	CHLORINE (TOTAL) IND SOLN 473ML	\$	74.97
2263411-LM	CHLORINE (TOTAL) IND SOLN 473ML	\$	74.97
2263511	BUFFER SOLN, TOT CHLORINE 473ML	\$	74.97
2263511-LM	BUFFER SOLN, TOT CHLORINE 473ML	\$	74.97
2264467	FUNNEL, POWDER PP 65MM	\$	16.96
2264472	FUNNEL, POWDER PP 150MM	\$	106.38
2264649	FLASK, DISTILLATION 500ML	\$	790.15
2264806	CONNECTOR, 20-400 TO 20-400 PK/6	\$	572.38
2264900	J-TUBE, 20-400 SCREW THREADS	\$	487.58
2265000	CONDENSER, 235MM 20-400 SCR THD	\$	930.84
2265100	DRIP TUBE, 20-400 SCREW THREAD	\$	167.09
2265200	THISTLE TUBE, 20-400 SCR THREAD	\$	368.10
2265300	DISTILLATION APPARATUS SET	\$	2,524.63
2265308	MANUAL, DISTILLATION PROCEDURES	\$	85.37
2265400	ARSENIC DISTILLATION APPAR SET	\$	1,144.76
2265500	GAS SCRUBBER, 20-400 SCR THREADS	\$	331.48
2265600	GAS BUBBLER, MDM DIFFUSION TUBE	\$	653.32
2265700	CYLINDER, SINGLE LINE 24-400	\$	292.93
2265800	CYANIDE DISTILLATION APPAR SET	\$	1,655.46
2266349	MEHLICH 2 SOIL EXTRACTANT, 500ML	\$	81.83
2266349-LM	MEHLICH 2 SOIL EXTRACTANT, 500ML	\$	81.83
2267301	SODIUM THIOSULFATE, 0.113N CRTG	\$	70.99
2267410	SULFITE EQUIV, 10ML AMP PK/16	\$	204.16
2267501	SODIUM THIOSULFATE, 0.2000N CRTG	\$	73.16
2267601	SODIUM THIOSULFATE, 0.2068N CRTG	\$	74.08
2268732	ALKALINE EDTA SOLN, 100ML MDB	\$	70.63
2268732-LM	ALKALINE EDTA SOLN, 100ML MDB	\$	70.63
2270701	aa CERIC STD SOLN, 0.5N CARTRIDGE	\$	75.33
2270800	HANDBOOK, STD METHODS WATER &	\$	1,268.10
2270803	ff AWWA WATER DISTRIBUTION OPERATOR TRAINING HANDBOOK	\$	360.39
2270804	AWWA WATER TREATMENT OPERATOR HANDBOOK	\$	354.60

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

2270806	AWWA WASTEWATER OPERATOR FIELD GUIDE	\$	360.39
2270900	TEST KIT, UNIV DIGITAL TITRATOR	\$	1,243.04
2270900-LM	TEST KIT, UNIV DIGITAL TITRATOR	\$	1,243.04
227142	IRON STD SOLN, 1000MG/L 100ML	\$	75.33
227142-LM	IRON STD SOLN, 1000MG/L 100ML	\$	75.33
2271900	REAGENT SET, ALKALINITY F/DT KIT	\$	249.33
2271900-LM	REAGENT SET, ALKALINITY F/DT KIT	\$	249.33
2272000	REAGENT SET, TOT HDNS F/DT KIT	\$	236.68
2272000-LM	REAGENT SET, TOT HDNS F/DT KIT	\$	236.68
2272100	REAGENT SET, CAL MAG TOT HDNS	\$	402.91
2272100-LM	REAGENT SET, CAL MAG TOT HDNS	\$	402.91
2272200	REAGENT SET, DO FOR DT KIT	\$	317.99
2272200-LM	REAGENT SET, DO FOR DT KIT	\$	317.99
2272300	REAGENT SET, SULFITE FOR DT KIT	\$	202.36
2272300-LM	REAGENT SET, SULFITE FOR DT KIT	\$	202.36
2272500	KTO: RGT SET, CHLORINE F/DT KIT	\$	271.01
2272600	rr REAGENT SET, CHLORIDE F/DT KIT	\$	310.76
2272600-LM	rr REAGENT SET, CHLORIDE F/DT KIT	\$	310.76
2272700	REAGENT SET, CARBON DIOXIDE	\$	200.55
2272700-LM	REAGENT SET, CARBON DIOXIDE	\$	200.55
2272800	REAGENT SET, ACIDITY FOR DT KIT	\$	285.47
2272800-LM	REAGENT SET, ACIDITY FOR DT KIT	\$	285.47
2273702	HEATER/STIRRER, 230V 50/60HZ	\$	5,627.42
2273800	SUPPORT ROD, ALUMINUM 1/2X24 IN	\$	88.65
2273900	CAP, BLK PHEN 13-425 PTFE/RUBBER	\$	10.21
2274000	CONNECTOR, THERM TO 13-425	\$	140.69
2274106	CONNECTOR, 24-400 TO 24-400 PK/6	\$	1,017.56
2274402	HEATER SUPPORT & APPAR, 230/50	\$	8,325.50
2275699	TPTZ IRON RGT PP 25ML PK/100	\$	120.33
2275699-LM	TPTZ IRON RGT PP 25ML PK/100	\$	120.33
2275806	TUBE, CULTURE DSPL 16X100MM PK/6	\$	25.05
2276145	BEAKER, BERZELIUS 200ML PK/1	\$	41.82
2277700	MEDIA SET, TOTAL BACTERIA COUNT, 15	\$	134.71
2277800	PRIMARY STD SET, FOR NITROGEN	\$	563.71
2277800-LM	PRIMARY STD SET, FOR NITROGEN	\$	563.71
2277924	AMMONIUM P-TOLUENESULFONATE, 25G	\$	176.70
2277924-LM	AMMONIUM P-TOLUENESULFONATE, 25G	\$	176.70
2278024	GLYCINE P-TOLUENESULFONATE, 25G	\$	187.90
2278124	NICOTINIC ACID P-TOLUENE-SULF 25G	\$	178.87
2279100	CAP, SAFETY VENT FOR ALNOCHROMIX*	\$	94.63
2280734	nn m-HPC AGAR, 500G	\$	553.11
2280934	PRESENCE-ABSENCE BROTH, 500G	\$	377.73
2280934-LM	PRESENCE-ABSENCE BROTH, 500G	\$	377.73
2281034	nn R2A AGAR, 500G	\$	553.11
2281234	nn M-E AGAR, 500G	\$	597.43
2282415	EC MEDIUM/MUG TUBES PK/15	\$	110.43
2283355	AMINO ACID F, PWD 410G	\$	303.53

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

2283426	BUFFER SOLN, RED PH 4.01, 50ML	\$	31.08
2283436	BUFFER SOLN, RED pH4.01, 15ML	\$	19.15
2283449	BUFFER SOLN, RED PH4.01 500ML	\$	45.17
2283453	ff BUFFER SOLN, RED PH4.01 1L	\$	87.81
2283456	BUFFER SOLN, RED PH4.01 4L	\$	178.31
2283461	BUFFER SOLN, RED PH 4.01, 20 LITER	\$	482.40
2283526	BUFFER SOLN, YELLOW pH 7.00 50ML	\$	31.08
2283536	BUFFER SOLN, YELLOW pH7.00, 15ML	\$	19.15
2283549	BUFFER SOLN, YELLOW PH7.00 500ML	\$	45.17
2283553	BUFFER SOLN, YELLOW PH7.00 1L	\$	87.81
2283556	BUFFER SOLN, YELLOW PH7.00 4L	\$	178.31
2283561	BUFFER SOLN, YELLOW PH 7.0, 20 LITER	\$	482.40
2283626	BUFFER SOLN, BLUE pH 10.01, 50ML	\$	31.08
2283636	ff BUFFER SOLN, BLUE pH10.01, 15ML	\$	19.15
2283649	BUFFER SOLN, BLUE PH10.01 500ML	\$	45.17
2283653	BUFFER SOLN, BLUE PH10.01 1L	\$	87.81
2283656	BUFFER SOLN, BLUE PH10.01 4L	\$	178.31
2283661	BUFFER SOLN, BLUE PH 10.01, 20 LITER	\$	482.40
2284799	METHYLTHYMOL BLUE IND PP PK/100	\$	73.52
2284799-LM	METHYLTHYMOL BLUE IND PP PK/100	\$	73.52
2285520	FRACTIONATING HEAD, W/O SIDEARM	\$	695.72
2285600	FUNNEL, CAPILLARY FLOW 3.25 IN	\$	402.78
2285901	DESICCANT, DRIERITE* 454G	\$	23.66
2287255	CALVER B CALCIUM INDICATOR	\$	451.69
2287255-LM	CALVER B CALCIUM INDICATOR	\$	451.69
2288000	REAGENT SET, CHLORIDE FOR	\$	292.69
2288000-LM	REAGENT SET, CHLORIDE FOR	\$	292.69
2289000	COUPLING, FEM 1/4 CPRSN NUT	\$	204.28
2291700	TEST KIT, HYDROGEN PEROXIDE DUAL	\$	314.13
2291700-LM	TEST KIT, HYDROGEN PEROXIDE DUAL	\$	314.13
2292301	FERROUS ETHYLENEDIAMONIUM SULF-	\$	59.79
2293566	SILVER 1 RGT PP PK/50	\$	218.62
2293566-LM	SILVER 1 RGT PP PK/50	\$	218.62
2293666	SILVER 2 RGT SOLN PLWS PK/50	\$	191.52
2293666-LM	SILVER 2 RGT SOLN PLWS PK/50	\$	191.52
2293766	SODIUM THIOSULFATE PWD PW, PK/50	\$	56.37
2294549	BARIUM CHLORIDE SOLN, 30% 500ML	\$	118.69
2296323	DIPHENYLAMINE SULPHONIC ACID IND	\$	61.60
2296600	REAGENT SET, SILVER CLRM METHOD	\$	422.78
2296600-LM	REAGENT SET, SILVER CLRM METHOD	\$	422.78
2297255	DPD COMPOUND, FREE & TOTAL 24G	\$	74.97
2297255-LM	DPD COMPOUND, FREE & TOTAL 24G	\$	74.97
2298504	KTO: ALKALI BUFFER FOR BRINE, PK/4X	\$	292.69
2298516	ALKALI BUFFER FOR BRINE, 946ML	\$	91.23
2298516-LM	ALKALI BUFFER FOR BRINE, 946ML	\$	91.23
2298604	HNB REAGENT PKG FOR LR CALCIUM/	\$	1,682.08
2298655	HNB PWD FOR LR CALCIUM/MAGNESIUM	\$	243.91

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

2299264	BUFFER, PH4.01 (RED) & PH7.00	\$	48.23
2299264-LM	BUFFER, PH4.01 (RED) & PH7.00	\$	48.23
2299300	TEST KIT, TPTZ IRON	\$	377.73
2299300-LM	TEST KIT, TPTZ IRON	\$	377.73
2299442	BEAKER, GRIFFIN DSPL PP 100ML	\$	14.07
2300051	ff DIGITAL REFRACTOMETER, BRIX 0 to 65%	\$	1,537.91
2300052	ff DIGITAL REFRACTOMETER BRIX 40-95 %	\$	1,537.91
2300054	ff DIGITAL REFRACTOMETER, SALINITY 0 to 280PPT	\$	1,863.60
2300055	ff DIGITAL REFRACTOMETER, AUTOMOTIVE	\$	1,863.60
230149	FERROZINE IRON RGT SOLN, 500ML	\$	256.56
230149-LM	FERROZINE IRON RGT SOLN, 500ML	\$	256.56
230166	FERROZINE RGT SOLN PLWS PK/50	\$	76.05
230166-LM	FERROZINE RGT SOLN PLWS PK/50	\$	76.05
23032	POTASSIUM HYDROXIDE 12N 100ML	\$	49.50
23032-LM	POTASSIUM HYDROXIDE 12N 100ML	\$	49.50
230424	FERROZINE POWDER 25G	\$	619.72
230424-LM	FERROZINE POWDER 25G	\$	619.72
230434	FERROZINE IRON REAGENT 500G	\$	6,197.15
23049H	POTASSIUM HYDROXIDE, 12N 500ML	\$	116.72
2305016	HNB DILUTION SOLVENT, 946ML	\$	94.48
2305042	HNB DILUTION SOLVENT 100ML	\$	81.83
2305042-LM	HNB DILUTION SOLVENT 100ML	\$	81.83
2305442	CALCIUM STD SOLN, 10.0MG/L 100ML	\$	85.81
2305726	AMMONIUM FLUORIDE 30% 50ML SCDB	\$	78.58
2305726-LM	AMMONIUM FLUORIDE 30% 50ML SCDB	\$	78.58
2305842	ALUMINUM STD SOLN, 10MG/L 100ML	\$	69.19
2305842-LM	ALUMINUM STD SOLN, 10MG/L 100ML	\$	69.19
2307242	ALKALI BUFFER SOLN 3.0N 100ML	\$	61.96
2307242-LM	ALKALI BUFFER SOLN 3.0N 100ML	\$	61.96
2307300	aa KTO: HNB REAGENT PKG	\$	99.90
2307300-LM	aa KTO: HNB REAGENT PKG	\$	99.90
2307355	HNB POWDER, ONE BOTTLE	\$	139.48
2307355-LM	HNB POWDER, ONE BOTTLE	\$	124.67
2307442	SODIUM CHLORIDE, 10246MG/L 100ML	\$	50.40
2307449	SODIUM CHLORIDE, 10246MG/L 500ML	\$	103.35
2307542	SODIUM CHLORIDE, 85.47MG/L 100ML	\$	50.23
2307553	SODIUM CHLORIDE, 85.47MG/L 1L	\$	138.76
2308100	TEST KIT, ASC-1 ASCORBIC ACID	\$	302.57
2308100-LM	TEST KIT, ASC-1 ASCORBIC ACID	\$	302.57
2308232	ASCORBIC ACID TITRANT 100ML MDB	\$	64.68
2308232-LM	ASCORBIC ACID TITRANT 100ML MDB	\$	64.68
2308249	ASCORBIC ACID TITRANT 500ML	\$	229.46
2308249-LM	ASCORBIC ACID TITRANT 500ML	\$	229.46
2308300	TUBE, CENTRIFUGE PS 15ML	\$	6.17
2309934	A-1 MEDIUM, BACTO 500G	\$	422.06
2310542	MASKING RGT, AL BRINE 100 ML	\$	119.78
2310542-LM	MASKING RGT, AL BRINE 100 ML	\$	119.78



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

2310629	ECR RGT SOLN, AL BRINE 200ML	\$	59.98
2310742	BUFFER RGT, AL BRINE 100 ML	\$	66.85
2310742-LM	BUFFER RGT, AL BRINE 100 ML	\$	66.85
2311149	HARDNESS BFR, AMINOETHANOL 500ML	\$	91.06
2312502	FLASK, DIGESDAHL W/STPR UNPACKED	\$	489.51
2312542	FLASK, DIGESDAHL W/STOPPER PKG	\$	522.27
2313018	MANUAL, DIGESDAHL APPARATUS COMPLIANT	\$	151.09
2313120	FRACTIONATING COLUMN KIT, W/O	\$	1,803.86
2314011	CHLORINE (FREE) IND, CL17 473ML	\$	74.97
2314011-LM	CHLORINE (FREE) IND, CL17 473ML	\$	74.97
2314111	CHLORINE (FREE) BFR, CL17 473ML	\$	69.02
2314111-LM	CHLORINE (FREE) BFR, CL17 473ML	\$	69.02
2314426	POTASSIUM HYDROXIDE STD, 1.0N	\$	41.19
2314426-LM	POTASSIUM HYDROXIDE STD, 1.0N	\$	41.19
2314500	TEST KIT, AL-TA ALKALINITY	\$	191.18
2314500-LM	TEST KIT, AL-TA ALKALINITY	\$	191.18
2315001	TEST KIT, POOL MASTER MGL	\$	728.48
2315001-LM	TEST KIT, POOL MASTER MGL	\$	728.48
2315040	INSERT, MOLDED-POOL MASTER	\$	63.98
2316949	ORP STD SOLN, ZOBELLS, 500ML	\$	186.10
2318069	CHLORIDE ADJ BFR PWD PLWS PK/100	\$	114.36
2318069-LM	CHLORIDE ADJ BFR PWD PLWS PK/100	\$	114.36
2318149	SODIUM STD SOLN, 100MG/L 500ML	\$	97.56
2318153	SODIUM STD SOLN, 100MG/L 1000ML	\$	112.02
2318406	BOTTLE, OBLONG HDPE 118ML PK/6	\$	19.85
2318506	DROPPER ASSY, CHILD RES PK/6 GLASS	\$	51.46
2319109	BUTTERFIELDS BUFFRD PHOSPHATE 90ML	\$	25.63
2319110	KTO: BUTTERFIELD'S BUFFERED PHOSPHATE	\$	462.53
2319125	KTO: BUTTERFIELD'S BUF PHOS 90ML PK/25	\$	582.01
2319145	BUTTERFIELDS BUFFERED PHOSPHATE	\$	51.46
2319800	rr KTO: REAGENT SET, CLRM CHLORIDE	\$	236.68
2319800-LM	rr KTO: REAGENT SET, CLRM CHLORIDE	\$	236.68
2319900	REAGENT SET, CA/MG HARDNESS	\$	249.33
2319900-LM	REAGENT SET, CA/MG HARDNESS	\$	249.33
2321142	PHOSPHOROUS STD, 1000MG/L 100ML	\$	70.63
2321353	HYDROCHLORIC ACID SOLN, 1.0M 1L	\$	69.02
2321353-LM	HYDROCHLORIC ACID SOLN, 1.0M 1L	\$	69.02
2321453	IODINE SOLN, 0.1M 1000ML	\$	93.95
2321753	POTASSIUM PERMANGANATE, 0.1N 1L	\$	56.90
2321953	SILVER NITRATE STD 0.100N 1000ML	\$	95.40
2321953-LM	SILVER NITRATE STD 0.100N 1000ML	\$	95.40
2323212	PRESENCE-ABSENCE TEST DSPL PK/12	\$	202.36
2323250	PRESENCE-ABSENCE TEST, PK/50	\$	686.32
2323250-LM	PRESENCE-ABSENCE TEST, PK/50	\$	686.32
2323810	FILTER, STERILIZING 0.2UM PK/10	\$	350.75
2324333	BOTTLE W/CAP, SQ PC 125ML PK/6	\$	175.95
2324373	nn BOTTLE W/CAP, SQ PC 125ML PK/48	\$	1,204.50

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

2324539	BOTTLE W/CAP, SQ PC 500ML PK/4	\$	175.38
2324579	BOTTLE W/CAP, SQ PC 500ML CASE/24	\$	969.38
2324663	ff BOTTLE W/CAP, SQ PC 1000ML PK/4	\$	235.12
23249	FLUORIDE STD SOLN, 100MG/L 500ML	\$	128.64
232706	BOTTLE, MXG 10,15,20,23ML MKS PK/6	\$	63.79
2327337	TUBING, C-FLEX* 1/4ID 1/16W 25FT	\$	283.30
2327367	TUBING, C-FLEX* 1/4ID 1/16W WHT	\$	12.33
2329032	FORMATE RGT SOLN, 100ML MDB	\$	141.65
2329232	BROMCRESOL GREEN-METHYL R, 100ML	\$	60.35
2329232-LM	BROMCRESOL GREEN-METHYL R, 100ML	\$	60.35
2329249	BROMCRESOL GREEN-METHYL R, 500ML	\$	99.54
2329249-LM	BROMCRESOL GREEN-METHYL R, 500ML	\$	99.54
2329332	WIDE RANGE-4 PH IND, 100ML MDB	\$	69.19
2329332-LM	WIDE RANGE-4 PH IND, 100ML MDB	\$	69.19
2329353	WIDE RANGE-4 PH IND SOLN, 1000ML	\$	261.98
2329353-LM	WIDE RANGE-4 PH IND SOLN, 1000ML	\$	261.98
2329549	BROMCRESOL GREEN IND SOLN, 500ML	\$	193.32
2329549-LM	BROMCRESOL GREEN IND SOLN, 500ML	\$	193.32
2329649	WIDE RANGE PH IND SOLN, 500ML	\$	139.48
2330242	PHENOLPHTHALEIN IND SOLN, 100ML	\$	45.70
2330242-LM	PHENOLPHTHALEIN IND SOLN, 100ML	\$	45.70
2330353	HYDROCHLORIC ACID, 0.02N 1000ML	\$	51.67
2332453	SODIUM HYDROXIDE STD SOLN, 6N 1L	\$	90.70
2333353	IODINE STD SOLN, 0.0282N 1000ML	\$	116.72
2333353-LM	IODINE STD SOLN, 0.0282N 1000ML	\$	116.72
2335549	rr ff BUFFER SOLN, PHOSPHATE 500ML	\$	50.03
2335649	DPD IND SOLN, 500ML	\$	122.86
2339053	BORIC ACID SOLN, 20G/L, 1000ML	\$	136.95
2339349	SULFURIC ACID 0.04N, 500ML	\$	56.90
234000	FILTER HOLDER, GLASS, 47MM	\$	1,840.48
2340249	STOCK NITRITE SOLN, 250 mg/L, 500ML	\$	119.78
2340249-LM	STOCK NITRITE SOLN, 250 mg/L, 500ML	\$	119.78
2340501	POTASSIUM SULFATE, 454G	\$	106.04
2340753	ff SODIUM HYDROXIDE THIOSULFATE RGT	\$	121.77
2343267	ff HOLMIUM TRICHLORIDE PP PK/20	\$	112.02
2343333	SYRINGE, 5CC LUER-LOK* PK/125	\$	175.95
234401	SULFAMIC ACID, 454G	\$	143.09
2344042	GLYCINE REAGENT, 100ML	\$	70.99
2345852	vv COD VIALS HG-F LR EEC PK/25	\$	191.52
2345952	vv MSDS: COD VIALS HG-F HR EEC PK/25	\$	184.29
2346800	REAGENT SET, SILICA LR DR/2000	\$	336.06
2346800-LM	REAGENT SET, SILICA LR DR/2000	\$	336.06
2347003	aa CITRIC ACID-SURFACTANT, 2.9L	\$	179.40
2347003-LM	aa CITRIC ACID-SURFACTANT, 2.9L	\$	179.40
2347900	STEAMBATH, CIRCULAR 5-RING 8 IN	\$	1,456.96
2348000	TIMER, 3-CHANNEL 1 SEC TO 99 HR	\$	364.24
2349332	SODIUM HYDROXIDE,.035N 100ML MDB	\$	66.85

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

2349332-LM	SODIUM HYDROXIDE, 035N 100ML MDB	\$	66.85
2349432	CHLORIDE TITRANT, HR 100ML MDB	\$	314.37
2349432-LM	CHLORIDE TITRANT, HR 100ML MDB	\$	314.37
2349500	TEST KIT, OS-1 OXYGEN SCAVENGER	\$	770.88
2349732	SULFURIC ACID, 0.035N 100ML MDB	\$	49.50
2349732-LM	SULFURIC ACID, 0.035N 100ML MDB	\$	49.50
2349832	SILVER NITRATE TITRANT 0.0493N 100ML	\$	97.56
2349832-LM	SILVER NITRATE TITRANT 0.0493N 100ML	\$	97.56
2349932	EDTA STD SOLN, 0.035N 100ML MDB	\$	57.98
2350206	CAP, WHT PP 20-400 FOAM LNR PK/6	\$	12.33
2350300	TEST KIT, BSC-1 BOILER FEED &	\$	447.11
2350340	INSERT, MOLDED BSC-1	\$	59.17
2350700	rr TEST KIT, PBC-DT PROFESSIONAL	\$	3,353.33
2350708	MANUAL, PROFESSIONAL BOILER	\$	135.67
2350800	TEST KIT, MN-PAN MANGANESE LR	\$	753.54
2350800-LM	TEST KIT, MN-PAN MANGANESE LR	\$	753.54
2351300	STANDARDIZATION KIT, GELEX* SSX	\$	1,007.93
2351400	REAGENT SET, DO F/BURET METHOD	\$	316.18
2351600	POTASSIUM IODIDE STARCH TEST	\$	31.22
2351749	POTASSIUM STD, 100MG/L 500ML	\$	156.27
2352449	MOLYBDENUM 1 LR PWD PLW PK/100	\$	103.51
2352512	MOLYBDENUM 2 RGT SOLN, 50ML MDB	\$	121.77
2352512-LM	MOLYBDENUM 2 RGT SOLN, 50ML MDB	\$	121.77
2352523	MOLYBDENUM 2 RGT SOLN 25ML SCDB	\$	118.69
2352532	MOLYBDENUM 2 RGT SOLN, 100ML MDB	\$	125.39
2352532-LM	MOLYBDENUM 2 RGT SOLN, 100ML MDB	\$	125.39
2352549	MOLYBDENUM 2 RGT SOLN, 500ML	\$	240.30
2352769	MOLYBDENUM 1 PP 50ML PK/100	\$	161.16
2353003	nn AMINO ACID F DILTN SOLVENT 2.7L	\$	178.68
2353003-LM	nn AMINO ACID F DILTN SOLVENT 2.7L	\$	178.68
2353011	AMINO ACID F DILTN SOLVENT 475ML	\$	87.45
2353011-LM	AMINO ACID F DILTN SOLVENT 475ML	\$	87.45
2353103	aa KTO: AMINO ACID F RGT PKG	\$	261.98
2353103-LM	aa KTO: AMINO ACID F RGT PKG	\$	261.98
2353255-LM	nn AMINO ACID F RGT PWD, 308G	\$	73.88
2354153	AMMONIA STD SOLN, 1000MG/L 1L	\$	79.30
2355237	DROPPER ASSY, CHILD RES PLASTIC TUBE	\$	10.60
2356901	BAG, PE ZIPPER 8X10X .006 IN	\$	5.40
2359299	MOLYBDENUM 1 RGT PP PK/100	\$	118.88
2359300	TEST KIT, MO-LR MOLYBDATE LR	\$	516.49
2365831	Tri-n-butyl phosphate, 250 mL	\$	70.10
236616	FERROZINE IRON RGT SOLN 946ML	\$	272.82
2367800	BATTERY, ZINC-AIR 1.4V PK/4	\$	59.94
2368531	PPB-1, ACID PRESERVATIVE 236ML	\$	58.71
2368531-LM	PPB-1, ACID PRESERVATIVE 236ML	\$	58.71
2368655	PPB-2, FIXER RGT SOLN 43ML	\$	99.54
2368655-LM	PPB-2, FIXER RGT SOLN 43ML	\$	99.54

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

2368749	PPB-3, ELUANT 500ML	\$	107.68
2368749-LM	PPB-3, ELUANT 500ML	\$	107.68
2368855	PPB-4, NEUTRALIZER 22ML	\$	114.55
2368855-LM	PPB-4, NEUTRALIZER 22ML	\$	114.55
2368964	PPB-5, IND REAGENT PWD PLW PK/20	\$	102.43
2368964-LM	PPB-5, IND REAGENT PWD PLW PK/20	\$	102.43
2369400	CLIPPER, (SHEARS) 7-1/4 INCH	\$	68.99
2370853	CHLORIDE STD SOLN, 0.00282N 1L	\$	101.35
2371006	CAP, PHENOLIC 33-430 PTFE/PE LNR PK/6	\$	33.15
2373220	M-FC AMPULES, GLASS 2ML PK/20	\$	106.38
2373250	M-FC, PLASTIC 2ML PK/50	\$	294.86
2373250-LM	M-FC, PLASTIC 2ML PK/50	\$	294.86
2373520	M-ENDO AMPULES, GLASS 2ML PK/20	\$	106.38
2373542	M-ENDO BROTH, 100ML	\$	152.63
2373550	M-ENDO PLASTIC AMPULES PK/50	\$	294.86
2373550-LM	M-ENDO PLASTIC AMPULES PK/50	\$	294.86
2373850	M-TGE AMPULES, PLASTIC 2ML PK/50	\$	294.86
2374755	PPB-6, DECOLORIZER 10ML SCDB	\$	81.83
2374755-LM	PPB-6, DECOLORIZER 10ML SCDB	\$	81.83
2374820	LEAD STD SOLN, 10MG/L 25ML	\$	45.89
2374820-LM	LEAD STD SOLN, 10MG/L 25ML	\$	45.89
2374900	Lab Accessory (Leadtrak PPB Extractor) No Chemistry	\$	13.30
2375000	REAGENT SET, LEADTRAK KIT	\$	814.84
2375000-LM	REAGENT SET, LEADTRAK KIT	\$	814.84
2375503	aa ANIONIC SURFACTANT SOLN 2.9L	\$	132.78
2375503-LM	aa ANIONIC SURFACTANT SOLN 2.9L	\$	132.78
2376400	SYRINGE PLUNGER, INTEGRATED SEAL	\$	14.07
2376526	POLYVINYL ALCOHOL DISPRSNG RGT 50ML	\$	45.70
2376549	POLYVINYL ALCOHOL DISPERSING	\$	116.52
2376626	MINERAL STABILIZER 50ML SCDB	\$	51.67
2376649	MINERAL STABILIZER 500ML	\$	238.49
237842	ZINC STD SOLN, 100MG/L 100ML	\$	61.60
237842-LM	ZINC STD SOLN, 100MG/L 100ML	\$	61.60
2380123	ECR MASKING RGT, 25ML SCDB	\$	50.40
2380123-LM	ECR MASKING RGT, 25ML SCDB	\$	50.40
2380268	ECR REAGENT PWD PLWS, PK/25	\$	69.02
2380268-LM	ECR REAGENT PWD PLWS, PK/25	\$	69.02
2382710	Acetylacetone 100 mL	\$	167.47
2386026	SODIUM SULFITE, MICRO GRADE 100G	\$	165.35
2386420	AMINO ACID F SOLN, 1.2ML AMPULES	\$	88.89
2386432	AMINO ACID F RGT SOLN, 100ML MDB	\$	79.50
2386432-LM	AMINO ACID F RGT SOLN, 100ML MDB	\$	79.50
2386442	AMINO ACID F RGT SOLN, 100ML	\$	79.50
2386442-LM	AMINO ACID F RGT SOLN, 100ML	\$	79.50
2387200	ff FILTER, MEMB 0.2UM 25MM PK/100	\$	489.51
2387300	FORCEPS, FILTER SST BENT TIP	\$	110.24
2387400	COVERSLIP, MICROSCOPE SLIDE 124/pk	\$	79.98



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

2393701	rr MERCURIC NITRATE CRTG 2.570N	\$	84.92
2395266	AMMONIA SALICYLATE RGT PP PK/50	\$	124.67
2395266-LM	AMMONIA SALICYLATE RGT PP PK/50	\$	124.67
2395366	AMMONIA SALICYLATE PP PK/50	\$	158.08
2395366-LM	AMMONIA SALICYLATE PP PK/50	\$	158.08
2395466	AMMONIA CYANURATE RGT PP PK/50	\$	77.50
2395466-LM	AMMONIA CYANURATE RGT PP PK/50	\$	77.50
2395566	AMMONIA CYANURATE RGT PP PK/50	\$	129.36
2395566-LM	AMMONIA CYANURATE RGT PP PK/50	\$	129.36
2401066	QAC REAGENT 1 PWD PLW PK/50	\$	78.58
2401268	QAC REAGENT 2, PWD PLWS PK/25	\$	39.75
2401268-LM	QAC REAGENT 2, PWD PLWS PK/25	\$	39.75
2401612	PRESENCE/ABSENCE W/MUG, PK/12	\$	244.75
2401650	PRESENCE-ABSENCE W/MUG, PK/50	\$	848.21
2401650-LM	PRESENCE-ABSENCE W/MUG, PK/50	\$	848.21
2401800	CAP, BLACK PP 22-400 PTFE/RUBBER	\$	7.13
2401812	CAP, BLK PP 22-400 PTFE/RBR PK/12	\$	65.52
2401906	SAMPLE CELL, 25X95MM 10-20-25ML PK/6	\$	183.47
2405200	ACCUVAC SNAPPER	\$	56.27
2406042	TTC SOLN 1%, STERILE 100ML GLASS	\$	158.03
2406549	NITROGEN, AMMONIA 100MG/L 500ML	\$	79.50
2406634	M-ENDO AGAR LES, DHYD 500G	\$	572.38
2406720	PLATE COUNT AGAR TUBES PK/20	\$	314.13
2407000	RACK, TEST TUBE, WHITE POLYPROPYLENE	\$	81.33
2407124	NALIDIXIC ACID, 25G	\$	420.13
2408449	SULFITE STD SOLN, 15MG/L 500ML	\$	176.15
2408532	SODIUM THIOSULFATE STD, 0.00246N	\$	56.01
2408549	SODIUM THIOSULFATE, .00246N 500ML	\$	99.54
2408732	SODIUM THIOSULFATE TITRANT, 100ML	\$	102.62
2408853	SODIUM THIOSULFATE, 0.00564N 1L	\$	120.33
2408932	SODIUM THIOSULFATE STD, 0.0109N	\$	56.90
2408949	SODIUM THIOSULFATE, .0109N 500ML	\$	77.50
2409032	SODIUM THIOSULFATE STD, 0.0113N	\$	57.82
2409101	SODIUM THIOSULFATE STD, 0.02256N	\$	70.99
2409232	SODIUM THIOSULFATE STD, 0.0246N	\$	56.01
2409249	SODIUM THIOSULFATE, .0246N 500ML	\$	70.99
2409301	SODIUM THIOSULFATE STD, 0.0250N	\$	73.16
2409317	SODIUM THIOSULFATE STD, 0.0250N	\$	256.56
2409353	SODIUM THIOSULFATE, .025N 1000ML	\$	91.06
2409453	SODIUM THIOSULFATE, .0375N 1000ML	\$	156.27
2410103	GLOVES, CHEMICAL RESIST 8-8-1/2	\$	25.63
2410104	GLOVES, CHEMICAL RESIST 9-9-1/2	\$	25.63
2410105	GLOVES, CHEMICAL RESIST SIZE 10	\$	25.63
2410106	GLOVES, CHEMICAL RESIST SIZE 11	\$	28.14
2410212	SAMPLE CELL, 1" SQ PK/12 PS W/COVER	\$	279.44
2410222	SAMPLE CELL, PS 1"SQ PK/2	\$	76.32
2412200	KIT, LONG PATH VIEWING ADAPTOR	\$	38.93

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

2415132	NITROGEN STOCK SOLN, 100ML MDB	\$	42.28
2415200	rr UV LAMP, PORTABLE	\$	223.56
2415342	QAC STANDARD, 100MG/L CTAB 100ML	\$	95.20
2415582	FUNNEL, POLYPROPYLENE 89MM	\$	14.26
2415600	FILTER, ACRODISC* 0.2UM 25MM	\$	21.78
2415815	rr COD DIGESTION VIAL, ULR PK/150	\$	912.07
2415815-LM	rr COD DIGESTION VIAL, ULR PK/150	\$	912.07
2415825	rr COD DIGESTION VIAL, ULR PK/25	\$	166.80
2415825-LM	rr COD DIGESTION VIAL, ULR PK/25	\$	166.80
2415915	rr COD DIGESTION VIAL, HR+ PK/150	\$	912.07
2415915-LM	rr COD DIGESTION VIAL, HR+ PK/150	\$	912.07
2415925	rr COD DIGESTION VIAL, HR+ PK/25	\$	166.80
2415925-LM	rr COD DIGESTION VIAL, HR+ PK/25	\$	166.80
241899	SULFIDE INHIBITOR PWD PW PK/100	\$	263.79
241899-LM	SULFIDE INHIBITOR PWD PW PK/100	\$	263.79
241906	CAP, BOD BOTTLE SNAP-OVER PK/6	\$	25.82
2425300	SAMPLE CELLS, 25ML MATCHED PAIR	\$	485.65
2426400GM	BACTO M-TGE BROTH, DYHD 500G	\$	784.37
2427601	nn SAMPLE CELL, 25X60MM 10ML W/CAP	\$	23.32
2427606	SAMPLE CELL, 25X60MM PK/6 10ML MARK	\$	135.67
2428350	M-GREEN Y&M, PLASTIC 2ML PK/50	\$	294.86
2428350-LM	M-GREEN Y&M, PLASTIC 2ML PK/50	\$	294.86
2428450	M-TGE - Tryptone Glucose Broth W/TTC indicator 2mL plastic (50 pkg)	\$	294.86
2428450-LM	M-TGE - Tryptone Glucose Broth W/TTC indicator 2mL plastic (50 pkg)	\$	294.86
2428520	M-FC/ROSOLIC ACID AMP PK/20 GLAS	\$	106.38
2428550	M-FC ROSOLIC ACID, PLASTIC 2ML PK/50	\$	294.86
2428550-LM	M-FC ROSOLIC ACID, PLASTIC 2ML PK/50	\$	294.86
2428700	TEST KIT, NI-SA AMMONIA	\$	352.68
2428700-LM	TEST KIT, NI-SA AMMONIA	\$	352.68
2428R0	SNSR-ORP,2W,LCP	\$	4,143.48
2428R025	ORP SNSR,ENCPS,CONV,LCP,	\$	4,394.02
2428R1	SNSR-ORP,2W,LCP,GOLD	\$	4,471.10
2428R105	SNSR-ORP,2W,LCP,GOLD,5FT.	\$	5,376.89
2428R125	SENSOR-ORP, 2 WIRE, LCP, GOLD ELTRD.	\$	4,740.91
2428U0	SNSR-ORP,2W,LCP,UL	\$	4,355.47
2429042	GLYCINE REAGENT, 10% SOLN 100ML	\$	70.99
2429104	PROCESS PH ELECTROLYTE PK/4	\$	359.54
2429149	PROCESS PH ELECTROLYTE 500ML	\$	113.63
2429300	RGT SET, ZINC 20ML	\$	220.42
2429300-LM	RGT SET, ZINC 20ML	\$	220.42
2429600	REAGENT SET, SILICA HIGH RANGE 10ML	\$	278.24
2429600-LM	REAGENT SET, SILICA HIGH RANGE 10ML	\$	278.24
2429700	REAGENT SET, PHOSPHONATE	\$	294.50
2429700-LM	REAGENT SET, PHOSPHONATE	\$	294.50
2429800	REAGENT SET, NITRATE NITROGEN LR	\$	285.47



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

2429800-LM	REAGENT SET, NITRATE NITROGEN LR	\$	285.47
2429900	rr ff REAGENT SET, NITROGEN, AMMONIA	\$	191.52
24-3000	SINGLE CAP. DEFOAM CONTRO	\$	6,070.68
2430000	REAGENT SET, MANGANESE	\$	225.84
2430000-LM	REAGENT SET, MANGANESE	\$	225.84
243001	rr TEST KIT, FF-2 FISHERIES/FISH	\$	2,389.73
243002	rr TEST KIT, FF-1A FISH FARMING	\$	1,329.77
243002-LM	rr TEST KIT, FF-1A FISH FARMING	\$	1,329.77
243003	rr TEST KIT, FF-3 FISH FARMING	\$	2,620.99
2430200	REAGENT SET, CYANIDE 10ML	\$	314.37
2430200-LM	REAGENT SET, CYANIDE 10ML	\$	314.37
243134	GLYCERIN, GLYCEROL, USP 500G	\$	170.56
2432309	BART TESTER, IRB PK/9	\$	481.80
2432327	KTO: BART TESTER, IRB PK/27	\$	1,403.00
2432409	BART TESTER, SRB PK/9	\$	481.80
2432427	KTO: BART TESTER, SRB PK/27	\$	1,403.00
2432509	BART TESTER, SLYM PK/9	\$	481.80
2432527	KTO: BART TESTER, SLYM PK/27	\$	1,403.00
2432609	BART TESTER, FLOR PK/9	\$	481.80
2432709	BART TESTER, ALG PK/9	\$	481.80
2432727	KTO: BART TESTER, ALG PK/27	\$	1,403.00
2434501	BISMUTH NITRATE, 0.0200M CRTG	\$	76.77
2434501-LM	BISMUTH NITRATE, 0.0200M CRTG	\$	76.77
2434702	SAMPLE CELL, PORTABLE TURB, PK/2	\$	98.87
2434706	SAMPLE CELL 2100 PORTABLE TURB, PK/6	\$	126.42
2434809	BART TESTER SET, IRB-SRB-SLYM	\$	481.80
2436466	NUTRIENT BFR SOLN PLWS 4ML PK/50	\$	128.64
2436466-LM	NUTRIENT BFR SOLN PLWS 4ML PK/50	\$	128.64
2437306	NUTRIENT AGAR W/MUG, 14ML PK/6	\$	158.22
2437401	REAGENT SET, ALKALINITY 0-400	\$	168.20
2437401-LM	REAGENT SET, ALKALINITY 0-400	\$	168.20
2437500	REAGENT SET, CARBON DIOXIDE	\$	92.14
2437600	REAGENT SET, CHLORIDE 0-400 MG/L	\$	164.05
2437600-LM	REAGENT SET, CHLORIDE 0-400 MG/L	\$	164.05
2437700	REAGENT SET, CHLORINE TOTAL	\$	186.10
2437700-LM	REAGENT SET, CHLORINE TOTAL	\$	186.10
2437800	REAGENT SET, HARDNESS 0-20MG/L	\$	149.77
2437900	REAGENT SET, HARDNESS 0-30 GPG	\$	97.20
2438600	REAGENT SET, CHLORINE TOTAL	\$	218.62
2438800	REAGENT SET, CHLORINE (F & T)	\$	144.90
2438800-LM	REAGENT SET, CHLORINE (F & T)	\$	144.90
2439200	REAGENT SET, COPPER 0-5 MG/L	\$	147.24
2439200-LM	REAGENT SET, COPPER 0-5 MG/L	\$	147.24
2439500	REAGENT SET, MOLYBDENUM	\$	229.46
2439802	REAGENT SET, DISSOLVED OXYGEN	\$	281.85
24399	HEXAVER POWDER PILLOWS PK/100	\$	71.55
24399-LM	HEXAVER POWDER PILLOWS PK/100	\$	71.55

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

2440000	REAGENT SET, SULFITE 0-200 MG/L	\$	209.58
2440000-LM	REAGENT SET, SULFITE 0-200 MG/L	\$	209.58
2440200	REAGENT SET, MOLYBDENUM	\$	281.85
2440400	REAGENT SET, CHLORINE (F&T)	\$	191.52
2440400-LM	REAGENT SET, CHLORINE (F&T)	\$	191.52
2443400	RGT SET, ALUMINUM 0-0.220 TEST KIT	\$	500.47
2443400-LM	RGT SET, ALUMINUM 0-0.220 TEST KIT	\$	500.47
2444100	REAGENT SET, CHLORIDE 10-100MG/L	\$	294.50
2444100-LM	REAGENT SET, CHLORIDE 10-100MG/L	\$	294.50
2444300	TEST KIT, AL-AP ALKALINITY P&M	\$	191.18
2444300-LM	TEST KIT, AL-AP ALKALINITY P&M	\$	191.18
2444301	TEST KIT, AL-AP MG/L ALKALINITY	\$	191.18
2444301-LM	TEST KIT, AL-AP MG/L ALKALINITY	\$	191.18
2444400	TEST KIT, CN-21P CHLORINE HIGH	\$	225.48
2444400-LM	TEST KIT, CN-21P CHLORINE HIGH	\$	225.48
2444800	KTO: RGT SET, CHLORINE 2000-70000	\$	388.45
2445300	KTO: RGT SET, CHLORINE (F & T)	\$	242.10
2446600	REAGENT SET, DEHA 0-450 UG/L	\$	271.01
2446600-LM	REAGENT SET, DEHA 0-450 UG/L	\$	271.01
2446732	Diacetyl, 99%, 100 g	\$	126.28
2446900	KTO: REAGENT SET, GLYCOL	\$	1,282.79
2447000	REAGENT SET, HARDNESS 0-5000MG/L	\$	187.90
2447000-LM	REAGENT SET, HARDNESS 0-5000MG/L	\$	187.90
2447200	REAGENT SET, HARDNESS (CALCIUM) 10-160MG/L	\$	189.71
2447200-LM	REAGENT SET, HARDNESS (CALCIUM) 10-160MG/L	\$	189.71
2447500	REAGENT SET, HARDNESS (CALCIUM) 100-4000	\$	189.71
2447500-LM	REAGENT SET, HARDNESS (CALCIUM) 100-4000	\$	189.71
2447600	REAGENT SET, HARDNESS 0-2500MG/L	\$	176.51
2447600-LM	REAGENT SET, HARDNESS 0-2500MG/L	\$	176.51
2447700	KTO: REAGENT SET, HARDNESS 100-25000	\$	209.58
2448000	REAGENT SET, HARDNESS 10-160	\$	180.12
2448000-LM	REAGENT SET, HARDNESS 10-160	\$	180.12
2448100	REAGENT SET, HARDNESS 100-4000	\$	180.12
2448100-LM	REAGENT SET, HARDNESS 100-4000	\$	180.12
2448200	REAGENT SET, HARDNESS 0-2500MG/L	\$	354.12
2448200-LM	REAGENT SET, HARDNESS 0-2500MG/L	\$	281.85
2448600	REAGENT SET, HARDNESS DT 10-160	\$	352.32
2448600-LM	REAGENT SET, HARDNESS DT 10-160	\$	352.32
2448700	REAGENT SET, HARDNESS DT 100-	\$	328.83
2448700-LM	REAGENT SET, HARDNESS DT 100-	\$	328.83
244903	aa SULFURIC ACID STD SOLN 2.9L	\$	242.10
244903-LM	aa SULFURIC ACID STD SOLN 2.9L	\$	242.10
2449100	REAGENT SET, HYDROGEN PEROXIDE	\$	283.30
2449200	KTO: RGT SET, IRON 10-100 MG/L	\$	464.33
2449300	KTO: RGT SET, IRON 100-1000 MG/L	\$	464.33
244932	SULFURIC ACID 5.25N 100ML MDB	\$	45.70
244932-LM	SULFURIC ACID 5.25N 100ML MDB	\$	45.70

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

2449400	REAGENT SET, MOLYBDENUM	\$	216.81
2449400-LM	REAGENT SET, MOLYBDENUM	\$	216.81
244953	SULFURIC ACID STD 5.25N, 1000ML	\$	76.05
244953-LM	SULFURIC ACID STD 5.25N, 1000ML	\$	76.05
245026	SODIUM HYDROXIDE 5.0N 50ML SCDB	\$	45.89
245026-LM	SODIUM HYDROXIDE 5.0N 50ML SCDB	\$	45.89
245032	SODIUM HYDROXIDE 5.0N 100ML MDB	\$	52.20
245032-LM	SODIUM HYDROXIDE 5.0N 100ML MDB	\$	52.20
245053	SODIUM HYDROXIDE 5.0N 1L	\$	90.14
245053-LM	SODIUM HYDROXIDE 5.0N 1L	\$	90.14
245199	POTASSIUM PERSULFATE PP PK/100	\$	111.66
245199-LM	POTASSIUM PERSULFATE PP PK/100	\$	111.66
2458200	rr REAGENT SET, NITROGEN-AMMONIA	\$	263.79
2458200-LM	REAGENT SET, NITROGEN-AMMONIA	\$	263.79
2458R0	SNSR-ORP,2W,LCP,UNION	\$	4,509.65
2459000	REAGENT SET, PHOSPHORUS (TOTAL)	\$	198.74
2459000-LM	REAGENT SET, PHOSPHORUS (TOTAL)	\$	198.74
2459100	REAGENT SET, POTASSIUM 0.25-2.50	\$	652.24
2459100-LM	REAGENT SET, POTASSIUM 0.25-2.50	\$	652.24
2459200	REAGENT SET, QAC 0-5 MG/L	\$	587.19
2459300	REAGENT SET, SILICA LR 10ML	\$	238.49
2459300-LM	REAGENT SET, SILICA LR 10ML	\$	238.49
2459500	REAGENT SET, CHLORINE (F&T)	\$	186.10
2459500-LM	REAGENT SET, CHLORINE (F&T)	\$	186.10
2459800	REAGENT SET, SULFITE 0-500+MG/L	\$	281.85
2459800-LM	REAGENT SET, SULFITE 0-500+MG/L	\$	281.85
246066	CYANURIC ACID 2 PWD PLWS PK/50	\$	87.09
246066-LM	CYANURIC ACID 2 PWD PLWS PK/50	\$	87.09
246102	nn FORMAZIN CALIB AMPULE 4000NTU	\$	135.51
24614	HYDROXYLAMINE HYDROCHLORIDE 113G	\$	155.55
246142	FORMAZIN TURB STD, 4000NTU 100ML	\$	111.46
246142-LM	FORMAZIN TURB STD, 4000NTU 100ML	\$	111.46
246149	FORMAZIN TURB STD, 4000NTU 500ML	\$	256.56
246149-LM	FORMAZIN TURB STD, 4000NTU 500ML	\$	256.56
2462300	REAGENT SET, SILICA HIGH RANGE 5ML	\$	339.67
2462300-LM	REAGENT SET, SILICA HIGH RANGE 5ML	\$	339.67
2463100	AUTOCLAVE TAPE, 1 INCH X 60 YD	\$	55.31
2463300	ff BAG, DECONTAMINATION PK/200	\$	489.51
2464101	TURBIDITY STD, GELEX 0-10NTU	\$	209.58
2464102	TURBIDITY STD, GELEX 0-100NTU	\$	209.58
2464103	TURBIDITY STD, GELEX 0-1000NTU	\$	209.58
2464105	GELEX STD SET, 2100 PORTABLE TURBS	\$	579.97
246800	FILTER HOLDER, SYRINGE LUER 25MM, PK 6	\$	466.38
246801	FILTER HOLDER, SYRINGE LUER 25MM, 1 EA	\$	90.58
2471100	TEST KIT, CN-DT CL2 W/DIG TITR	\$	1,050.32
2471100-LM	TEST KIT, CN-DT CL2 W/DIG TITR	\$	1,050.32
2471515	EC/MUG W/O DURHAM TUBES, PK/15	\$	110.43

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

247632	M-NITROPHENOL IND SOLN 100ML MDB	\$	81.47
247632-LM	M-NITROPHENOL IND SOLN 100ML MDB	\$	81.47
2478409	BART TESTER, POOL PK/9	\$	481.80
2481300	KTO: TEST KIT, STORM WATER	\$	1,877.09
2481569	PHENOL REAGENT PWD PLW PK/100	\$	137.48
2481569-LM	PHENOL REAGENT PWD PLW PK/100	\$	137.48
2483600	TEST KIT, PHENOL	\$	537.69
2483600-LM	TEST KIT, PHENOL	\$	537.69
2484600	AMPULE BREAKER, POURITE	\$	60.13
2486100	MANIFOLD, 3-PLACE ALUMINUM	\$	6,571.75
2488349	NITRATE IONIC STR ADJ SOLN 500ML	\$	117.44
2489443	FLASK, ERLLENMEYER WM 125ML	\$	35.08
2489446	FLASK, ERLLENMEYER WM 250ML	\$	35.65
2489449	FLASK, ERLLENMEYER WM 500ML	\$	37.77
2489454	FLASK, ERLLENMEYER WM 2000ML	\$	89.04
2489500	ULTRA-SONIC BATH, .75 GAL W/HEAT	\$	5,473.25
2489737	PIPET, SER DSPL WM 5ML PK/25	\$	79.02
2489738	PIPET, SER DSPL WM 10ML PK/20	\$	84.80
2490300	FILTER, MEMB BLACK 47MM PK/200	\$	616.70
2490409	BART TESTER, HAB PK/9	\$	481.80
2490427	KTO: BART TESTER, HAB PK/27	\$	1,403.00
2493023	BLANKING REAGENT, 29ML	\$	43.36
2493120	BUFFER SOLN, 1.5ML AMP PK/20	\$	83.11
2493220	DPD IND SOLN, 1.5ML AMP PK/20	\$	76.77
2494925	PRESENCE-ABSENCE BROTH AMPULES,PK/25	\$	273.66
2495050	BOTTLE W/CAP, 100ML STER PK/50	\$	458.67
2495300	rr KTO: RGT SET, COLORIMETRIC TKN	\$	635.98
2495402	SAMPLE CELL, 10ML MATCHED PK/2	\$	520.34
2495402S	SAMPLE CELL SUPER, 10ML MATCHED PK/2	\$	487.58
2495408	SAMPLE CELL, 10ML MATCHED PK/8	\$	1,576.45
2495525	PRESENCE-ABSENCE W/MUG AMPULES PK/25	\$	356.53
2495525-LM	PRESENCE-ABSENCE W/MUG AMPULES PK/25	\$	356.53
2495900APR	APPARATUS SET FOR TEST KIT, NPK-1 SOIL NPK	\$	1,965.74
2495989	MANUAL, TEST INTERPRETATION NPK1	\$	126.81
2496000RGT	ff REAGENT SET, SIW 1 SOIL/IRRIGATION	\$	4,952.90
2496088	Manual, SIW-1 TEST KIT	\$	133.36
2496400	STIRRER, TEST TUBE HAND HELD	\$	85.37
2496510	STIRRING RODS, HAND STIR PK/10	\$	47.41
2497900	RACK, TEST TUBE 13MM 6X15 ROWS	\$	89.04
2497901	RACK, TEST TUBE 16MM 5X12 ROWS	\$	89.04
2497902	RACK, TEST TUBE 20MM 4X10 ROWS	\$	89.04
2497903	RACK, TEST TUBE 25MM 3X8 ROWS	\$	89.04
2497904	RACK, TEST TUBE 30MM 3X7 ROWS	\$	89.04
2497912	RACK, TEST TUBE 20MM 4X5 ROWS	\$	89.04
2498032	SODIUM HYDROXIDE, .075N 100ML MDB	\$	75.33
2498132	EDTA STD SOLN, 0.0075N 100ML MDB	\$	67.02
2501025	DISSOLVED OXY LR ACCUVAC, PK/25	\$	111.66



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

2501025-LM	DISSOLVED OXY LR ACCUVAC, PK/25	\$	111.66
2502025	DPD FREE CHLORINE ACCUVAC, PK/25	\$	97.20
2502025-LM	DPD FREE CHLORINE ACCUVAC, PK/25	\$	97.20
2502749	25 PPM DIVO LE	\$	178.31
2502849	50 PPM DIVO LE	\$	178.31
2502949	100 PPM DIVO LE	\$	178.31
2503025	DPD TOT CHLORINE ACCUVAC, PK/25	\$	97.20
2503025-LM	DPD TOT CHLORINE ACCUVAC, PK/25	\$	97.20
2503400	CLAMP, SWIVEL 2 PRONG, LARGE	\$	208.14
25045000	DIGITAL COUNTER PEN	\$	105.03
2504810	CHLORINE, 1000-1500MG/L, PK/16	\$	182.48
2506025	SPADNS FLUORIDE ACCUVAC, PK/25	\$	111.66
2506025-LM	SPADNS FLUORIDE ACCUVAC, PK/25	\$	111.66
2506101	CLEANING BRUSH FOR SLUDGE JUDGE, 6FT	\$	133.94
2506102	SAMPLER, WATER CORE, 15', LARGE	\$	753.54
25061020	MODBUS RS485 CARD, AF7000	\$	1,325.91
2506103	CLEANING BRUSH, SLUDE JUDGE, LARGE	\$	133.94
25061030	AUTOFLUSH CARD, 110V S VALVE, AF7000	\$	2,967.89
25061040	AUTOFLUSH CARD, 220V S VALVE, AF7000	\$	2,967.89
25061050	TRANSMITTER DISPLAY PANEL, AF7000	\$	3,584.59
25062000	MOTOR GEARBOX ASSEMBLY 230V, AF7000	\$	7,824.43
25066000	OPTICAL SENSOR, REPLACEMENT, AF7000	\$	830.62
25067000	TRANSMITTER UNIT, AF7000	\$	10,773.05
25068000	KIT, FASTNERS, 79 PIECES, AF7000	\$	383.51
25069000	KIT, ELBOW CAP, AF7000	\$	830.62
2507025	FERROVER IRON ACCUVAC, PK/25	\$	113.63
2507025-LM	FERROVER IRON ACCUVAC, PK/25	\$	113.63
25071000	KIT, O RING, AF7000	\$	115.25
2507149	SODIUM NITRATE, 0.5M, 500 ML	\$	86.72
25071490	KIT, COLLAR CLAMP ASSEMBLY, AF7000	\$	265.95
2507200	KTO: PH BUFFER SOLUTION KIT 4L	\$	354.12
25072000	POWER CORD, RPLCMNT, AF7000, EUROPE	\$	190.21
2507300	KTO: PH BUFFER SOLUTION 20L KIT	\$	953.96
25073000	POWER CORD, REPLACEMENT, AF7000, USA	\$	143.19
25074010	KIT, PISTON HOUSING BRACKET, AF7000	\$	593.58
2507500	SPEC COLOR STD KIT, MONOCHLOR/F-NH3	\$	563.71
25075000	KIT, SENSOR SPRING & WASHER, AF7000	\$	106.57
2507600	GLASS FIBER PADS, 200/PK	\$	238.97
25076000	KIT, ACCESSPORT w/HEX BOLT, AF7000	\$	147.05
25078000	FUSE KIT, AF7000	\$	167.47
25079000	3 CONDUCTOR TERMINAL BLOCK, AF7000	\$	40.28
2508025	PHOSVER 3 ACCUVAC, PK/25	\$	112.91
25081000	GEARBOX COUPLER, AF7000	\$	279.44
25083000	BACK PLATE GASKET, AF7000	\$	537.69
25086000	OPTICAL SENSOR DISC, AF7000	\$	235.12
25088000	WATER CONNECTION W/AUTO FLUSH, 220V, AF7000	\$	6,167.04
25089000	SOLENOID FOR AUTO FLUSH CARD, 110V, AF7000	\$	2,813.71

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

25091000	GRIT FILTER, 500 MICRON, AF7000	\$	925.06
25092000	SOLENOID FOR AUTO FLUSH CARD, 220V, AF7000	\$	2,987.16
25093000	Main Board, AF7000	\$	4,683.10
2510253	N-HEXANE, ULTRA PURE, 1-L BOTTLE	\$	337.86
2510600	REAGENT SET GB VERIFICATION KIT	\$	1,739.90
2511025	oo NITRAVER 5 NITRATE ACCUVAC PK/25	\$	112.91
2511201	FILTER, GLASS FBR, 21MM PK/100	\$	103.68
2511202	FILTER, GLASS FBR, 24MM PK/100	\$	117.37
2511203	FILTER, GLASS FBR, 25MM PK/100	\$	117.37
2511204	FILTER, GLASS FBR, 35MM PK/100	\$	150.51
2511205	FILTER, GLASS FBR, 37MM PK/100	\$	150.51
2511206	FILTER, GLASS FBR, 42.5MM PK/100	\$	171.52
2511207	FILTER, GLASS FBR, 47MM PK/100	\$	208.14
2511208	FILTER, GLASS FBR, 55MM PK/100	\$	235.12
2511209	FILTER, GLASS FBR, 70MM PK/100	\$	260.17
2511210	FILTER, GLASS FBR, 90MM PK/100	\$	462.53
2511211	FILTER, GLASS FBR, 110MM PK/100	\$	526.13
2511212	FILTER, GLASS FBR, 125MM PK/100	\$	616.70
2511301	FILTER, GLASS FBR, VSS 21MM PK/100	\$	92.51
2511302	FILTER, GLASS FBR, VSS 24MM PK/100	\$	104.07
2511303	ff FILTER, GLASS FBR, VSS 25MM PK/100	\$	104.07
2511306	FILTER, GLASS FBR, VSS 42.5MM PK/100	\$	151.09
2511307	FILTER, GLASS FBR, VSS 47MM PK/100	\$	183.66
2511309	FILTER, GLASS FBR, VSS 70MM PK/100	\$	229.34
2511310	FILTER, GLASS FBR, VSS 90MM PK/100	\$	414.35
2511311	FILTER, GLASS FBR, VSS 110MM PK/100	\$	458.67
2511312	FILTER, GLASS FBR, VSS 125MM PK/100	\$	551.18
25118026	KTO: SAT. KCL SOLUTION, 59 ML, WITH DISPENSING KIT	\$	60.80
25118049	SAT. KCL SOLUTION, 500 ML	\$	126.28
25118059	SAT. KCL SOLUTION, 59 ML	\$	43.72
2511900	PIPET TIPS, 200-1000UL, 1000/PK	\$	155.91
2512100	WATER PURIF SYS, 100-230VAC, 50/60HZ	\$	17,421.89
2512200	rr WATER PURIF SYS, UV 100-230VAC, 50/60HZ	\$	20,197.06
2512201	PURIFICATION PACK FOR DQ3 WATER SYS	\$	2,293.37
2512202	0.22UM FILTER FOR DQ3 WATER SYSTEMS	\$	724.63
2512203	ff VENT FILTER, SIMPLICITY WATER SYSTEM	\$	235.12
2512204	WALL MOUNT BRACKET FOR DQ3 SYSTEM	\$	838.33
2512210	ff KTO: REAGENT SET, AMTAX II 0.2-12.0MG/L	\$	607.07
2512211	AMTAX II STD SOLN 10MG/L 237ML	\$	215.85
2512212	ff KTO: REAGENT SET, AMTAX II 0.2-30.0MG/L	\$	799.79
2512213	AMTAX II STD SOLN 25MG/L 237ML	\$	215.85
2512214	ff KTO: REAGENT SET, AMTAX II 2-120MG/L	\$	607.07
2512215	AMTAX II STD SOLN 100MG/L 237ML	\$	215.85
2512216	ff KTO: REAGENT SET, AMTAX II 20-1200MG/L	\$	607.07
2512217	AMTAX II STD SOLN 1000MG/L 237ML	\$	215.85
2512300	WATER SYS, 1-5L, 100-230VAC, 50/60HZ	\$	14,665.99
2512400	rr WATER SYS, 1-5L, UV 90.255VAC 50/60HZ	\$	16,940.09



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

2512401	PURIFICATION PACK 1, WATER PUR SYS	\$	1,032.98
2512402	PURIFICATION PACK 2, WATER PUR SYS	\$	1,032.98
2512403	FILTER FOR WATER SYSTEMS 1-5L	\$	516.49
251242	KTO: MEL POTABLE WATER LAB	\$	11,659.56
2512500	WATER SYS, 1-20L, 100-230 V/50-60 HZ	\$	15,591.05
2512600	rr WATER SYS, 1-20L UV 100-230V/50-60HZ	\$	19,117.82
2512602	ff PURIFICATION PACK 2, SYNERGY WATER SYS	\$	1,491.65
2512604	PRESSURE GAUGE/REGULATOR	\$	1,085.01
2512700	rr UV LAMP 185NM	\$	1,584.16
2512800	BIOPAK PURIFICATION CARTRIDGE	\$	1,032.98
2512900	30L PE RESERVOIR, DQ3 SYSTEM	\$	4,355.47
2512902	ff TANK WALL BRACKET, 30L RESERVOIR	\$	481.80
2513603	PVDF SYRINGE FILTER, .45UM, 33MM, 50/PK	\$	790.15
2513605	NYLON SYRINGE FILTER, .2UM 33MM, 50/PK	\$	873.02
2513608	PVDF SYRG FLTR, .22UM, 33MM, 50/PK	\$	788.22
2513801	FILTER, GLASS FIBER 0.7UM 47MM, PK/100	\$	238.97
2514101	FILTER, MEMB 0.45UM 25MM PK/100	\$	422.06
2514200	PREFILTERS, 5UM PP, 45MM, 48/PK	\$	618.63
2514300	SPE SOLVENT RECOVERY KIT	\$	5,164.90
2514420	EZ GGA, BOD, NO PIPETTE, 2ML PK/20	\$	149.77
2514654	0.5 MG/L NH4-N STD SOLN, 2000 ML	\$	314.37
2514754	2.5 MG/L NH4-N STD SOLN, 2000 ML	\$	334.25
2514836	ELECTROLYTE AMTAX, (LOW-RANGE), 11 ML	\$	196.94
2515025	DISSOLVED OXY HR ACCUVAC, PK/25	\$	111.66
2515025-LM	DISSOLVED OXY HR ACCUVAC, PK/25	\$	111.66
2515050	TEST KIT, ACCUVAC DO HIGH RANGE	\$	491.44
2515050-LM	TEST KIT, ACCUVAC DO HIGH RANGE	\$	491.44
2516025	OZONE ACCUVAC, 0-0.25MG/L PK/25	\$	116.72
2516025-LM	OZONE ACCUVAC, 0-0.25MG/L PK/25	\$	116.72
2517025	OZONE ACCUVAC, 0-0.75MG/L PK/25	\$	116.72
2517025-LM	OZONE ACCUVAC, 0-0.75MG/L PK/25	\$	116.72
2517150	MI BROTH, 50/PK	\$	516.49
2517150-LM	MI BROTH, 50/PK	\$	516.49
2517600	EMPTY COD VIALS WITH CAPS, PK/25	\$	117.37
251764	PARAFILM 4 IN X 125 FT ROLL	\$	106.77
2518025	OZONE ACCUVAC, 0-1.5MG/L PK/25	\$	116.72
2518025-LM	OZONE ACCUVAC, 0-1.5MG/L PK/25	\$	116.72
2518050	TEST KIT, ACCUVAC OZONE HR	\$	466.38
2518050-LM	TEST KIT, ACCUVAC OZONE HR	\$	466.38
25189000	SILICONE SPEC CHECK, HR CHLORINE, STD 3	\$	0.00
2519900	ff REPLACEMENT SENSOR	\$	354.60
25207059	1.0M KNO3 SOLUTION, 59 ML	\$	39.56
252074	GLASS WOOL FOR DETERGENTS 5 GRAM	\$	38.93
2521100	FUNNEL FILTER, .45UM, STERILE, BULK	\$	526.13
2521101	FUNNEL FILTER, .45UM, STERILE, INDV	\$	826.77
25221000	STABLCAL 600 NTU, 500 ML	\$	283.66
25224000	5500sc AMC Colorimeter Cleaning Kit	\$	49.53



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

25225010S	5500sc Ammonia Monochloramine Cleaning Soln, 15mL	\$	40.83
25227000	nn 10 NTU TU5200 StablCal (Production Use Only)	\$	832.91
25228000	nn 20 NTU TU5200 StablCal (Production Use Only)	\$	832.91
2522900	COTTON SWABS FOR CL-17 ANALYZER PK/1000	\$	158.22
25229000	nn 600 NTU TU5200 StablCal (Production Use Only)	\$	834.72
2522905	COTTON SWAB FOR CL-17, 5/PK	\$	8.09
25233000	KTO: 5500SC AMMONIA MONOCHLORAMINE REAGENT SET	\$	682.95
2523500	POWER INLET CORD	\$	22.55
25235000	5500SC REAGENT 2, AMMONIA/MONOCHLORAMINE	\$	87.45
2523600	STIRRER BASE	\$	262.10
25236000	5500SC REAGENT 3, AMMONIA/MONOCHLORAMINE	\$	124.47
2523700	WASHER GB/T96-1985	\$	1.54
25237000	5500SC STANDARD 1, AMMONIA MONOCHLORAMINE	\$	77.14
25238000	5500SC STANDARD 2, AMMONIA MONOCHLORAMINE	\$	96.48
2523900	RS232 CABLE	\$	21.20
25239000	5500SC ACIDIC SURFACTANT WASH	\$	92.14
2524100	Processor plate with hoses	\$	3,989.30
2524300	aa KTO: RGT SET COD BULK LR, 300 TESTS	\$	395.68
2524400	aa KTO: RGT SET COD BULK HR, 300 TESTS	\$	395.68
25248000	EO STERILIZED PETRI DISH W/ PAD	\$	400.86
25249000	SET OF 5X RUBBER SEALS, 5500SC/961X	\$	269.81
2525400	DIELECTRIC GREASE	\$	936.62
2527025	SPADNS2 (ARSENIC FREE) FLUORIDE RGT, ACCUVAC,25/PK	\$	128.64
2527025-LM	SPADNS2 (ARSENIC FREE) FLUORIDE RGT, ACCUVAC,25/PK	\$	128.64
2527100	KTO: CONSTRUCTION EFFLUENT KIT	\$	5,338.34
25277000	PROGNOSYS 5500sc AMC LICENSE KIT	\$	3,738.77
2527900	KTO: HIGH RANGE BORON STARTER SET	\$	1,337.48
2528300	MAGNETIC STIRR BAR 4MM X 12MM	\$	8.67
25286000	KTO: 5500SC SILICA REAGENT SET, R2 MODIFIED	\$	1,396.62
25286000INT	KTO: 5500SC SILICA REAGENT SET, R2 MODIFIED INT	\$	1,355.06
25286000-LM	KTO: 5500SC SILICA REAGENT SET, R2 MODIFIED	\$	1,396.62
25293000	KTO: ELECTRODELESS-PEEK, SWP W ADAPTOR	\$	3,758.04
25294000	KTO; ORP SENSOR-RYTON, SWP W ADAPTOR	\$	3,835.13
25295000	KTO: pH SENSOR-RYTON, SWP W ADAPTOR	\$	3,603.86
25296000	KTO: FP360sc 500 ppb, w ADAPTOR	\$	71,768.93
25297000	KTO: ASSY, PROBE, LDO MODEL 2, HACH, SWP W ADAPTOR	\$	7,381.18
25298000	KTO: T-LINE sc/IMMERS 4000NTU WIPER PVC, SWP w ADAPTOR	\$	14,030.02
25299000	KTO: TS-LINE sc/IMMERS 50g/l WIPER PVC, SWP w Adaptor	\$	17,209.90
253000	FILTER, GLASS FBR 47MM PK/100	\$	148.59
25300000	KTO: TSS sc, INSITU PROBE, W/O WIPER, SWP w Adaptor	\$	22,297.70
253051	FILTER, GLASS FBR 42.5MM PK/100	\$	157.45
253053	FILTER, GLASS FBR 70MM PK/100	\$	235.12
253055	FILTER, GLASS FBR 55MM PK/100	\$	235.12
25318000	CITRIC ACID RGT SOLN, 2L MODIFIED	\$	263.79
25318000-LM	CITRIC ACID RGT SOLN, 2L MODIFIED	\$	263.79
2532800	BOTTLE W/CAP, F/HYDROGEN SULFIDE	\$	16.38
253334	NITRIFICATION INHIBITOR 500G	\$	614.30

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

253334-LM	NITRIFICATION INHIBITOR 500G	\$	614.30
253335	NITRIFICATION INHIBITOR 35G	\$	78.58
253335-LM	NITRIFICATION INHIBITOR 35G	\$	78.58
2537600	1 1/2" SS WELD-IN SOCKET	\$	275.59
2537733	HYDROGEN SULFIDE PPR, CUT PK/100	\$	79.30
2537733-LM	HYDROGEN SULFIDE PPR, CUT PK/100	\$	79.30
2537800	TEST KIT, HS-C HYDROGEN SULFIDE	\$	192.14
2537801	TEST KIT, HS-C1 HYDROGEN SULFIDE, NO ALK	\$	143.19
2537801-LM	TEST KIT, HS-C1 HYDROGEN SULFIDE, NO ALK	\$	143.19
2537900	CHART, HYDROGEN SULFIDE HS-C	\$	14.65
254049	NITRIC ACID SOLN, 1:1 500ML	\$	109.30
254049-LM	NITRIC ACID SOLN, 1:1 500ML	\$	109.30
2541249	AMMONIA ELCTD CONDITIONING SOLN 500ML	\$	87.45
2541259	AMMONIA ELCTD CONDITIONING SOLN 50ML	\$	37.39
25432	BROMCRESOL PURPLE, 100ML MDB	\$	67.21
2543768	FERROMO IRON RGT 1 PWD PLW PK/25	\$	87.45
2543768-LM	FERROMO IRON RGT 1 PWD PLW PK/25	\$	87.45
2543866	FERROMO IRON RGT 2 PWD PLW PK/50	\$	84.92
2543866-LM	FERROMO IRON RGT 2 PWD PLW PK/50	\$	84.92
2544800	REAGENT SET, FERROMO	\$	493.24
2544800-LM	REAGENT SET, FERROMO	\$	493.24
2545049	SODIUM HYDROXIDE SOLN, 10N 500ML	\$	82.92
2546100	FILTER PAPER, PREWEIGHED PK/100	\$	695.72
2546110	FILTER, 42.5MM, PREWEIGHED PK/100	\$	714.99
2546120	FILTER, 55MM, PREWEIGHED PK/100	\$	695.72
2546135	FILTER, 70MM, PREWEIGHED PK/100	\$	784.37
2546140	FILTER, 90MM, PREWEIGHED PK/100	\$	884.58
2546150	ff FILTER, 110MM, PREWEIGHED PK/100	\$	961.67
2546200	FILTER, 47MM, PREWEIGHED PK/400	\$	1,480.09
2546240	FILTER, 90MM, PREWEIGHED PK/400	\$	2,351.18
2546250	FILTER, 110MM, PREWEIGHED PK/400	\$	2,601.72
2546259	ff FILTER PAPER, 18.5CM DIA PK/100	\$	379.66
2546902	POT CL REF ELCTRLT GEL CRTG PK/2	\$	51.84
2546902-LM	POT CL REF ELCTRLT GEL CRTG PK/2	\$	51.84
2550203	APRON, LONG SLEEVE SIZE LARGE	\$	439.40
2550501	GLOVES, NITRILE SMALL PK/100 DISP	\$	82.48
2550502	GLOVES, NITRILE MEDIUM PK/100 DISP	\$	82.48
2550503	GLOVES, NITRILE LARGE PK/100 DISP	\$	82.48
2550504	GLOVES, NITRILE X-LRG PK/100 DISP	\$	82.48
2550700	GOGGLES, SAFETY VENTED	\$	33.15
2550800	RESPIRATOR, DISPOSABLE PK/20	\$	157.64
2551000	FOAM PAD, PE FOR CASE 46601-00	\$	3.85
2551400	FILTER, GLASS FBR 47MM A/E PK/100	\$	190.21
2551456	FILTER, GLASS FBR 47MM A/E PK/100	\$	200.43
2551458	FILTER, GLASS FBR 70MM A/E PK/100	\$	304.50
2551459	FILTER, GLASS FBR 90MM A/E PK/100	\$	304.50
2551461	FILTER, GLASS FBR 110MM A/E PK/100	\$	489.51

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

2551462	FILTER, GLASS FBR 125MM A/E PK/100	\$	637.90
2551600	FILTER, GLASS FBR 21MM 1.1UM PK/100	\$	337.26
2551801	RACK, BOTTLE 20 COMPARTMENTS	\$	162.27
2551910	DISH, ALUMINUM 90MM PK/50	\$	341.11
2552114	SAMPLE CONTAINER, 45ML PK/600	\$	643.68
2552200	SAMPLE CONTAINER, 120ML HIGH PF	\$	8.67
2552212	SAMPLE CONTAINER,120ML HP PK/250	\$	360.39
2552411	SAMPLE CONTAINER, 240ML PK/140	\$	281.37
2552511	SAMPLE CONTAINER, 300ML PK/100	\$	260.17
25532	BROMTHYMOL BLUE, 100ML MDB	\$	56.90
25532-LM	BROMTHYMOL BLUE, 100ML MDB	\$	56.90
2553500	REAGENT SET, SILICA ULR BULK RGTS	\$	350.51
2553500-LM	REAGENT SET, SILICA ULR BULK RGTS	\$	350.51
2554300	SWABS, COTTON-TIP APPLCTR, STERILE, PK/100	\$	88.46
25549	BROMTHYMOL BLUE IND SOLN, 500ML	\$	112.91
25549-LM	BROMTHYMOL BLUE IND SOLN, 500ML	\$	112.91
255600	TEST KIT, IR-24 IRON FERROZINE	\$	398.93
255600-LM	TEST KIT, IR-24 IRON FERROZINE	\$	398.93
2556256	BOTTLE, RND JUG HDPE 38-400 4L	\$	15.61
2556900	aa REAGENT SET, CHLORINE FREE CL17	\$	195.13
2556900-LM	aa REAGENT SET, CHLORINE FREE CL17	\$	195.13
2557000	aa KTO: REAGENT SET, CHLORINE TOTAL CL17	\$	204.16
2557000G	vv REAGENT SET, CHLORINE TOTAL CL17	\$	187.90
2557000-LM	aa KTO REAGENT SET, CHLORINE TOTAL CL17	\$	204.16
25587000	9610sc REAGENT 2 SILICA 2L MODIFIED	\$	263.79
25588000	KTO: 9610sc SILICA REAGENT SET MODIFIED	\$	1,355.06
25588001	KTO: 9610sc SILICA REAGENT SET MODIFIED, EU	\$	1,355.06
2558996	PIPET TIP, STERILE 1-10ML PK/50	\$	283.30
2559701	TENSETTE NUMBER RING KIT 0-1	\$	78.44
2559710	TENSETTE NUMBER RING KIT 1-10	\$	78.44
2559800	TEST KIT, SURFACE WATER	\$	1,298.93
2559800-LM	TEST KIT, SURFACE WATER	\$	1,298.93
2559833	TEST KIT, EDUC SURFACE WATER	\$	1,304.71
2559833-LM	TEST KIT, EDUC SURFACE WATER	\$	1,304.71
2559836	WATER: A COMPREHENSIVE GUIDE FOR BREWERS	\$	104.26
256032	TANNIVER 3 REAGENT, 100ML MDB	\$	91.23
256032-LM	TANNIVER 3 REAGENT, 100ML MDB	\$	91.23
256049	TANNIVER 3 REAGENT, 500ML	\$	298.11
2560766	TOXTRAK REAGENT PWD PLW PK/50	\$	155.19
2560766-LM	TOXTRAK REAGENT PWD PLW PK/50	\$	155.19
2560836	TOXTRAK ACCELERATOR SOLN, 15ML	\$	86.72
2560836-LM	TOXTRAK ACCELERATOR SOLN, 15ML	\$	86.72
2560915	A-1 MEDIUM MPN TUBES, PK/15	\$	113.13
2560915-LM	A-1 MEDIUM MPN TUBES, PK/15	\$	113.13
256200	TEST KIT, IR-23 IRON 0-7 PH4-10	\$	699.57
2563000	REAGENT SET, ULTRA LOW CHLORINE	\$	177.42
2563137	DISPENSER, DIGITAL 1.0-5.0 ML	\$	1,580.30

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

2563140	DISPENSER, DIGITAL 1.0-10.0 ML	\$	1,803.86
25632200	KTO SL1000 PPA INSTRUMENT, NITRIFICATION KIT	\$	16,439.02
25638018	KTO: BOD LAB SYSTEM, HQ430D, LBOD101	\$	8,537.50
2563899	CAL STD, 25PPB, 40ML, PAH500	\$	476.98
2563901	CAL STD, 50PPB, 40ML, PAH500	\$	476.98
2563906	5500SC REAGENT 1, AMMONIA MONOCHLORAMINE, 1L	\$	204.16
2564325	TRYPTIC SOY BROTH, 100ML PK/25	\$	1,032.98
2565025	COD DIGEST VIAL, COD2 LR PK/25	\$	172.00
2565025-LM	COD DIGEST VIAL, COD2 LR PK/25	\$	172.00
2565115	COD DIGEST VIAL, COD2 HR PK/150	\$	939.51
2565115-LM	COD DIGEST VIAL, COD2 HR PK/150	\$	939.51
2565125	COD DIGEST VIAL, COD2 HR PK/25	\$	172.00
2565125-LM	COD DIGEST VIAL, COD2 HR PK/25	\$	172.00
2567620	ff FILTER & BARREL ASSEMBLY PK/20	\$	481.80
2567729	SOIL EXTRACTANT, 200ML	\$	118.52
2568500	COOLER, 9 QT.	\$	145.70
2569100	FILTER, GLASS FBR 110MM PK/100	\$	356.53
2569320	SPATULA, DISPOSABLE PK/20	\$	29.68
256949	PHOSPHATE STD SOLN, 1MG/L 500ML	\$	82.92
2569600	KTO: MEL P/A SAFE DRINKING WATER LAB	\$	10,734.50
2569600APR	CASE/APP MEL P/A SAFE DRINKING WATER LAB	\$	8,980.75
2569600RGT	REAGENT SET, MEL P/A SAFE DRINKING WATER LAB (2569600)	\$	449.04
2569688	MANUAL, MEL P/A	\$	312.21
2569693	PROCEDURES, MEL P/A, SPANISH	\$	64.37
2569694	PROCEDURES, MEL P/A, FRENCH	\$	42.21
2569695	PROCEDURES, MEL P/A, GERMAN	\$	42.21
2569700	KTO: MEL/MF LAB	\$	15,456.14
2569700A	MEL/MF LAB components A	\$	7,573.90
2569700A-LM	MEL/MF LAB components A	\$	7,573.90
2569700B	MEL/MF LAB components B	\$	5,184.17
2569700-LM	KTO: MEL/MF LAB	\$	15,456.14
2569800	Microbiological Lab/Most Probable Number	\$	8,441.14
2569800-LM	Microbiological Lab/Most Probable Number	\$	8,441.14
2569889	PACKAGING GUIDE, MEL/MPN, ENG	\$	11.18
2569893	PACKAGING GUIDE, MEL/MPN,SPANISH	\$	11.18
2569894	PACKAGING GUIDE, MEL/MPN, FRENCH	\$	5.40
2569895	PACKAGING GUIDE, MEL/MPN, GERMAN	\$	5.40
2569900	INCUBATOR, PORTABLE 12VDC	\$	4,586.74
2569922	oo PROCEDURES, MEL/MF & MEL/MPN	\$	71.69
2569988	MANUAL, MEL/P-A, MEL/MF & MEL/MPN	\$	283.30
2569993	PROCEDURES, MEL/MF&MPN, SPANISH	\$	69.57
2569994	PROCEDURES, MEL/MF&MPN, FRENCH	\$	42.21
2569995	PROCEDURES,MEL.MF&MPN,GERMAN	\$	42.21
2570305	FILTER ELEMENT, 0.75UM PK/5	\$	871.09
257201	COTTON BALLS, ABSORBENT PK/100	\$	33.15
25732	THYMOL BLUE IND SOLN, 100ML MDB	\$	70.63
25732-LM	THYMOL BLUE IND SOLN, 100ML MDB	\$	70.63



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

2575010	PUMP, SUBMERSIBLE SMALL	\$	826.77
257849	SULFATE STD SOLN, 50MG/L 500ML	\$	82.92
257924	N-SERVE*, TCMP 25G	\$	334.25
257924-LM	N-SERVE*, TCMP 25G	\$	334.25
2580100	KTO: MF APPARATUS REPLACEMENT SET	\$	3,603.86
2580200	KTO: MEL/MPN Consumables Set	\$	260.17
2580500	RACK, P/A BOTTLE	\$	70.34
2580501	RACK, MPN TUBES	\$	70.34
2580502	RACK, MF PLATES & GEN PURPOSE	\$	70.34
2580600	TEST KIT, CDS-DT CHLORIDE	\$	1,189.08
2580600-LM	TEST KIT, CDS-DT CHLORIDE	\$	1,189.08
2581400	REAGENT SET, SILICA ULR AMPULED RGT	\$	296.31
2581802	CAPS, FLIP SPOUT PK/2	\$	25.25
2584202	FORMAZIN STD, 7500NTU AMPULE PKG	\$	265.59
2584202-LM	FORMAZIN STD, 7500NTU AMPULE PKG	\$	265.59
2584335	FUNNEL, MICRO PP 35MM TOP ID	\$	12.72
2584600	LAMP, REPL FOR 24152-00 4W UVLW	\$	211.99
2585000	FUSE, 15A REPL F/ PORT INCUBATOR	\$	74.58
2585220	ff CAP, CENTRIFUGE TUBE 16MM PK/20	\$	106.77
2585400	MAGNIFIER, ILLUMINATED 3X	\$	67.84
2586201	MILLIPORE EZ FIT MANIFOLD	\$	5,376.89
2586300	FUNNEL, PUSH FIT STERILE PK/150	\$	1,108.14
2586400	KTO: MEL FIELD FILTRATION KIT, COLIFORM & E. COLI	\$	6,398.30
2586400RGT	RGT SET, MEL FIELD FILT KIT, COLIFORM & E. COLI (2586400)	\$	736.19
2587200	TEST KIT, GT-1 GLUTARALDEHYDE	\$	1,283.52
2587200-LM	TEST KIT, GT-1 GLUTARALDEHYDE	\$	1,283.52
2587300	RACK, TEST TUBE GRIP STYLE PP	\$	70.34
2589000	GELEX SEC TURB STD SET, 0-4000 NTU	\$	870.85
258901	FLUORIDE ADJ BUFFER 454G	\$	310.76
2589100	GELEX SEC TURB STD, STRAY LIGHT	\$	209.58
2589101	GELEX SEC TURB STD, 0-2NTU	\$	209.58
2589102	GELEX SEC TURB STD, 0-20NTU	\$	209.58
2589103	GELEX SEC TURB STD 0-200NTU	\$	209.58
2589104	GELEX SEC TURB STD, 200-4000NTU	\$	209.58
2589105	GELEX SEC TURB STD, 4000-10,000	\$	209.58
2589200	GELEX SEC TURB STD SET, 0-10,000 NTU	\$	1,038.88
2589549	CHLOROPHOSPHONAZO IND SOLN 500ML	\$	2,565.59
2589549-LM	CHLOROPHOSPHONAZO IND SOLN 500ML	\$	2,565.59
2589599	CHLOROPHOSPHONAZO SOLN PLWS PK/100	\$	776.90
2589599-LM	CHLOROPHOSPHONAZO SOLN PLWS PK/100	\$	776.90
2589723	CDTA SOLUTION, 30ML SCDB	\$	74.97
2589723-LM	CDTA SOLUTION, 30ML SCDB	\$	74.97
258999	FLUORIDE ADJ 25ML PWD PW PK/100	\$	193.32
258999-LM	FLUORIDE ADJ 25ML PWD PW PK/100	\$	193.32
2590401	GLOVES, LATEX SMALL PK/100	\$	151.09
2590402	ff GLOVES, LATEX MEDIUM PK/100	\$	151.09
2590403	GLOVES, LATEX LARGE PK/100	\$	151.09

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

2592920	SOIL SAMPLE CONTAINER PK/20	\$	838.33
259342	COPPER STD SOLN, 1000MG/L 100ML	\$	119.05
259342-LM	COPPER STD SOLN, 1000MG/L 100ML	\$	119.05
2594025	FILTER, VERSAPOR 3UM 25MM PK/25	\$	97.13
2594642	WATER, ULTRA PURE 100ML	\$	32.88
2594649	WATER, ULTRA PURE 500ML	\$	143.26
2595600	APPARATUS SET, ULTRA LR CHLORINE	\$	527.57
2595802	AMMONIUM CHLORIDE GEL CRTG, PK/2	\$	51.84
2595802-LM	AMMONIUM CHLORIDE GEL CRTG, PK/2	\$	51.84
259600	GLASS BEADS, FOR DO PK/100	\$	63.98
2597200	aa REAGENT SET, TOXTRAK KIT	\$	527.57
2597200-LM	aa REAGENT SET, TOXTRAK KIT	\$	527.57
2598500	CARTRIDGE, PRETREATMENT DI FEED	\$	350.75
2598501	CARTRIDGE, PRETREATMENT DSTL/RO	\$	335.33
2598502	CARTRIDGE, HIGH PURITY LOW TOC	\$	437.47
2598503	CARTRIDGE, ULTRAPURE	\$	350.75
2598505	FILTER, FINAL 0.2 MICRON	\$	291.01
2599738	ff SYRINGE FILTERS, 25MM 1.5U 934-AH	\$	2,293.37
2599803	aa MOLYBDATE RGT SOLN, 2.9L FOR LR	\$	231.26
2599803-LM	aa MOLYBDATE RGT SOLN, 2.9L FOR LR	\$	231.26
2599849	MOLYBDATE RGT SOLN, 500ML	\$	89.06
2599903	ASCORBIC ACID RGT DILUTION SOLN,	\$	261.98
2599949	ASCORBIC ACID DILUTN SOLN, 450ML	\$	89.78
25M1A1001-115	PH7.0 BUFFER, STD CELL 500ML	\$	200.55
25M1A1001-115AR	PH7.0 BUFFER, STD CELL 500ML	\$	215.00
25M1A1001-115CN	PH7.0 BUFFER, STD CELL 500ML	\$	215.00
25M1A1001-115LM	PH7.0 BUFFER, STD CELL 500ML	\$	200.55
25M1A1001-115TH	PH7.0 BUFFER, STD CELL 500ML	\$	215.00
25M1A1001-123	PH7.0 BUFFER, STD CELL 3.5L	\$	776.90
25M1A1014-123	PH4.0 BUFFER, 3.5L CALIB SOLN	\$	460.72
25M1A1016-123	PH7.0 BUFFER, 3.5L CALIB SOLN	\$	460.72
25M1A1017-123	PH10.0 BUFFER, 3.5L CALIB SOLN	\$	460.72
25M1A1025-115	EQUITRANSFERRANT PH7 BUFFER, 500ML	\$	229.46
25M1A1025-115BR	EQUITRANSFERRANT PH7 BUFFER, 500ML	\$	243.91
25M1A1025-115CN	EQUITRANSFERRANT PH7 BUFFER, 500ML	\$	243.91
25M1A1025-115KR	EQUITRANSFERRANT PH7 BUFFER, 500ML	\$	243.91
25M1A1025-115LM	EQUITRANSFERRANT PH7 BUFFER, 500ML	\$	229.46
25M1A1025-115VN	EQUITRANSFERRANT PH7 BUFFER, 500ML	\$	243.91
25M1A1025-123	EQUITRANSFERRANT PH7 BUFFER, 3.5L	\$	888.92
25M2A1001-115	200MV ORP SOLN, 500ML	\$	200.55
25M2A1001-115KR	200MV ORP SOLN, 500ML	\$	215.00
25M2A1001-119	200MV ORP SOLN, 1L	\$	229.46
25M2A1001-123	200MV ORP SOLN, 3.5L	\$	599.84
25M2A1002-115	600MV ORP SOLN, 500ML	\$	200.55
25M2A1002-115AR	600MV ORP SOLN, 500ML	\$	215.00
25M2A1002-115BR	600MV ORP SOLN, 500ML	\$	215.00
25M2A1002-115LM	600MV ORP SOLN, 500ML	\$	200.55

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

25M2A1002-115TW	600MV ORP SOLN, 500ML	\$	215.00
25M2A1002-119	600MV ORP SOLN, 1L	\$	229.46
25M2A1002-119CN	600MV ORP SOLN, 1L	\$	243.91
25M2A1002-119LM	600MV ORP SOLN, 1L	\$	229.46
25M2A1002-123	600MV ORP SOLN, 3.5L	\$	599.84
25M2A1002-123CN	600MV ORP SOLN, 3.5L	\$	655.85
25M2A1002-123LM	600MV ORP SOLN, 3.5L	\$	599.84
25M3A2000-100	COND SOLN, 100US/CM 1L	\$	419.17
25M3A2000-1000	COND SOLN, 1000US/CM 1L	\$	419.17
25M3A2000-200	COND SOLN, 200US/CM 1L	\$	419.17
25M3A2000-400	COND SOLN, 400US/CM 1L	\$	419.17
25M3A2000-500	COND SOLN, 500US/CM 1L	\$	419.17
25M3A2000-600	COND SOLN, 600US/CM 1L	\$	419.17
25M3A2050-1500	COND SOLN, 1500US/CM 1L	\$	458.91
25M3A2050-2000	COND SOLN, 2000US/CM 1L	\$	458.91
25M3A2100-100K	COND SOLN, 100,000US/CM 1L	\$	484.21
25M3A2100-10K	COND SOLN, 10,000US/CM 1L	\$	484.21
25M3A2100-150K	COND SOLN, 150,000US/CM 1L	\$	484.21
25M3A2100-2000	COND SOLN, 2000US/CM 1L	\$	484.21
25M3A2100-2500	COND SOLN, 2500US/CM 1L	\$	484.21
25M3A2100-3000	COND SOLN, 3000US/CM 1L	\$	484.21
25M3A2100-5000	COND SOLN, 5000US/CM 1L	\$	484.21
25M3A2100-50K	COND SOLN, 50,000US/CM 1L	\$	484.21
25M3A2200-200K	COND SOLN 200,000US/CM 1L	\$	608.87
25M3A2200-300K	COND SOLN 300,000US/CM 1L	\$	608.87
25M3A2300-350K	COND SOLN 350,000US/CM 1L	\$	691.99
25M3A2300-450K	COND SOLN 450,000US/CM 1L	\$	691.99
25M3A2300-500K	COND SOLN 500,000US/CM 1L	\$	691.99
25M7A1002-105	ELECTROLYTE FOR 5440D, 60ML	\$	285.47
25M8A1002-101	PROPRIETARY GEL, 2 GRAMS GLI	\$	200.55
25M8A1002-106	CHEM.PROPRIETARY GEL, 25	\$	444.46
25M8A1002-107	CHEM.PROPRIETARY GEL, 30	\$	603.45
25M9A1004-107	CLEANING SOLN FOR 5440D, 60ML	\$	285.47
2600-0252-001	nn ff FLEXY EWON UNIT, FOR REMOTE ACCESS	\$	6,725.93
2600103	aa WATER, ZERO STD 2.9L F/LR PHOS	\$	87.25
2600303	aa KTO: ASCORBIC ACID RGT PKG	\$	352.32
2600800	SPATULA, STERILE SAMPLER 6 INCH	\$	8.67
2600900	SCOOP, STERILE SAMPLER 60ML	\$	9.64
2601300	PH TEST STRIP, 0-14 PH PK/100	\$	78.77
2602853	STD COLOR SOLN 15 UNITS, 1000ML	\$	336.06
2602853-LM	STD COLOR SOLN 15 UNITS, 1000ML	\$	336.06
2603100	RGT SET, ULR HARDNESS TOTAL	\$	836.53
2603100-LM	RGT SET, ULR HARDNESS TOTAL	\$	836.53
2603101	REAGENT SET, ULR HARDNESS BULK	\$	2,547.52
2603101-LM	REAGENT SET, ULR HARDNESS BULK	\$	2,547.52
2603300	REAGENT SET, COPPER 10ML	\$	706.44
2603300-LM	REAGENT SET, COPPER 10ML	\$	706.44

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

2603449	COPPER MASKING RGT PP PK/100	\$	202.36
2603549	PORPHYRIN 1 RGT PP 10ML PK/100	\$	209.58
2603649	PORPHYRIN 2 RGT PP 10ML PK/100	\$	220.42
2603649-LM	PORPHYRIN 2 RGT PP 10ML PK/100	\$	220.42
2603700	REAGENT SET, ALUMINUM ECR 10ML	\$	496.86
2603700-LM	REAGENT SET, ALUMINUM ECR 10ML	\$	496.86
2603849	ECR REAGENT PWD PLW 20ML PK/100	\$	238.49
2603849-LM	ECR REAGENT PWD PLW 20ML PK/100	\$	238.49
2603999	HEXAMETHYLENETETRAMINE PK/100	\$	231.26
2603999-LM	HEXAMETHYLENETETRAMINE PK/100	\$	231.26
2604100	REAGENT SET, MOLYBDENUM 10ML	\$	354.12
2604100-LM	REAGENT SET, MOLYBDENUM 10ML	\$	354.12
2604299	MOLYVER 1 RGT PP 10ML PK/100	\$	122.50
2604299-LM	MOLYVER 1 RGT PP 10ML PK/100	\$	122.50
2604299NZ	MOLYVER 1 RGT PP, 10ML, PK/100, NEW ZEALAND	\$	126.28
2604399	MOLYVER 2 RGT PP 10ML PK/100	\$	122.50
2604399-LM	MOLYVER 2 RGT PP 10ML PK/100	\$	122.50
2604499	MOLYVER 3 RGT PP 10ML PK/100	\$	122.50
2604499-LM	MOLYVER 3 RGT PP 10ML PK/100	\$	122.50
2604545	RGT SET, TNT AMVER LR	\$	368.58
2604545-LM	RGT SET, TNT AMVER LR	\$	368.58
2604945	RGT SET, TNT TIN PRETREATMENT	\$	364.96
2604945-LM	RGT SET, TNT TIN PRETREATMENT	\$	364.96
2605150	TOT INORG NITROGEN REDUCT AMP PK/50	\$	229.46
2605150-LM	TOT INORG NITROGEN REDUCT AMP PK/50	\$	229.46
2605345	RGT SET, TNT NITRAVER X 50 TESTS	\$	283.66
2605345-LM	RGT SET, TNT NITRAVER X 50 TESTS	\$	283.66
2605546	NITRAVER X REAGENT B PP PK/50	\$	103.88
2605546-LM	NITRAVER X REAGENT B PP PK/50	\$	103.88
2605620	CHROMIUM 6, 5.0MG/L 2ML PK/20	\$	222.23
2605820	MANGANESE, 10MG/L 2ML AMP PK/20	\$	193.32
2605820-LM	MANGANESE, 10MG/L 2ML AMP PK/20	\$	193.32
2606945	RGT SET, TNT AMVER HR 50 TESTS	\$	337.86
2606945-LM	RGT SET, TNT AMVER HR 50 TESTS	\$	337.86
2608345	RGT SET, TNT NITRIVER 3 PK/50	\$	205.97
2608345-LM	RGT SET, TNT NITRIVER 3 PK/50	\$	205.97
2608420	M-COLIBLUE24 AMP, PK/20 GLASS	\$	144.11
2608420-LM	M-COLIBLUE24 AMP, PK/20 GLASS	\$	144.11
2608442	M-COLIBLUE24 BROTH, 100ML	\$	177.78
2608442-LM	M-COLIBLUE24 BROTH, 100ML	\$	177.78
2608450	M-COLIBLUE24 PLASTIC AMPULES PK/50	\$	338.22
2608450-LM	M-COLIBLUE24 PLASTIC AMPULES PK/50	\$	338.22
2608799	TPTZ IRON RGT PP 10ML PK/100	\$	122.50
2608799-LM	TPTZ IRON RGT PP 10ML PK/100	\$	122.50
2610696	PATHOSCREEN PWD PLW, 100ML PK/50	\$	204.28
2610696-LM	PATHOSCREEN PWD PLW, 100ML PK/50	\$	204.28
2610796	PATHOSCREEN PP 20ML PK/100	\$	144.35

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

2610796-LM	PATHOSCREEN PP 20ML PK/100	\$	144.35
2610810	PADDLE TSTR,TOT BAC/YEAST/MOLD PK/10	\$	174.60
2610910	PADDLE TSTR,TOT BACT/TOT COLIF PK/10	\$	174.60
2612520	ORP STANDARD, 20ML AMPULE PK/25	\$	522.15
2612520-LM	ORP STANDARD, 20ML AMPULE PK/25	\$	522.15
2612602	SAMPLE CELL, W/STOPPER MATCHED PAIR	\$	830.62
2612849	PITCHER, GRDTD TALL FORM 500ML	\$	33.92
2612853	PITCHER, GRDTD TALL FORM 1000ML	\$	35.65
2612854	PITCHER, GRDTD TALL FORM 2000ML	\$	37.77
2612953	CYLINDER, GRDTD W/HANDLE 1000ML	\$	87.30
2612954	CYLINDER, GRDTD W/HANDLE 2000ML	\$	169.40
2613049	BOTTLE, WASH SAFETY ASSORT 500ML	\$	119.49
2613212	MARKER, BLACK TIP F/COLONY CNTR	\$	72.46
2615199	PHTHALATE-PHOSPHATE PP, PK/100	\$	101.90
2615199-LM	PHTHALATE-PHOSPHATE PP, PK/100	\$	101.90
2615436	ff CDTA SOLUTION, 0.4M 15ML SCDB	\$	46.06
2615714	ff ELECTRODE HOLDER, F/FLEX LAB ARM	\$	132.98
2615900	STILL, WATER AUTO 2L/HR 120V	\$	16,997.90
2616100	BLENDER, 1.2 LITER 2 SPEED 120V	\$	2,158.46
2616300	WATER BATH, CIRCULATING 22L 110V	\$	11,505.38
2616302	WATER BATH, CIRCULATING 22L 220V	\$	11,042.86
2617300	KIT, SOLENOID VALVE MAINTENANCE	\$	867.24
2617501	POTASSIUM PERSULFATE, ACS 454G	\$	195.13
2617614	ff WEIGHT, BALANCE 100GM ASTM CL 1	\$	506.85
2617700	HAND SOAP, ANTIBACTERIAL 32OZ	\$	139.14
2618100	DESSICATOR CABINET, NON-VACUUM	\$	639.83
2618800	FILTER, MEMB 0.2UM 47MM PK/200	\$	728.48
2618901	FILTER, MEMB 0.45UM 47MM PK/100	\$	514.56
2619000	FILTER FUNNEL, MICRO .45UM PK/50	\$	969.38
2619200	INCUBATOR, CULT LOW PROFILE 110V	\$	10,869.41
2619202	INCUBATOR, CULT LOW PROFILE 230V	\$	10,869.41
2619309	BART TESTER, DN PK/9	\$	481.80
2619407	BART TESTER, N PK/7	\$	481.80
2619510	PADDLE TESTER, DISINFECTION	\$	174.60
2619793	MANUAL, BODTRAK SPANISH	\$	80.36
2620000	GELEX SEC TURB STD SET, 2100NIS	\$	870.85
2620532	SULFURIC ACID STD, 0.030N 100ML	\$	45.89
2620532-LM	SULFURIC ACID STD, 0.030N 100ML	\$	45.89
2620553	SULFURIC ACID STD, 0.030N 1000ML	\$	66.13
2621500	GABLE COVER, WATERBATH 26PC	\$	973.24
2622500	OXIDASE REAGENT, 0.5ML PK/50	\$	337.26
2623415	COD VIALS, MANGANESE III PK/150	\$	939.51
2623415-LM	COD VIALS, MANGANESE III PK/150	\$	939.51
2623425	COD VIALS, MANGANESE III PK/25	\$	172.00
2623425-LM	COD VIALS, MANGANESE III PK/25	\$	172.00
2624401	SAMPLE CELL, 100MM QUARTZ	\$	2,158.46
2624410	SAMPLE CELL, 10MM QUARTZ	\$	713.06

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

2624450	SAMPLE CELL, 50MM QUARTZ	\$	1,333.62
2626000	MANUAL, SENTRY SAMPLE MANIFOLD	\$	333.41
2628000	POWER SUPPLY, 12VDC 120VAC 60HZ	\$	487.58
2629201	SAMPLE CELL, 100MM OPTICAL GLASS	\$	969.38
2629250	SAMPLE CELL, 50MM OPTICAL GLASS	\$	529.98
2629500	SAMPLE CELL, 10MM 1.5ML PK/100	\$	248.61
2630020	CHLORINE, 25-30MG/L 2ML PK/20	\$	164.94
2630020-LM	ff CHLORINE, 25-30MG/L 2ML PK/20	\$	164.94
2630088	INSTRUCTIONS, CHLORINE 25-30MG/L	\$	8.31
2630200	STOPWATCH, WATER RESISTANT 24HR	\$	103.68
2630300	TIMER, COUNT UP/DOWN 99MIN LCD	\$	329.55
2630400	TIMER, 5 CHANNEL 1SEC TO 100HRS	\$	346.90
2630600	THERMOMETER, DIGITAL CERTIFIED	\$	967.45
2630601	PROBE, SS 3-PURPOSE OPTION F/2630600	\$	373.88
2630605	KTO: THERMOMETER WITH PROBE, DIGITAL CERTIFIED	\$	1,056.11
2631740	ff LAB STIRRER, PROGRAMMABLE 4-PADDLE JAR TESTER 2L	\$	12,642.43
2634001	SAMPLE VIAL, 120ML PK/200 THIO & TIE	\$	423.98
2634002	SAMPLE VIAL, 120ML PP PK/200 STERILE	\$	418.20
2634800	CONNECTOR, 1/16" NPT TO 1/8" OD	\$	102.91
2635300	SPEC COLOR STD, DPD CHLORINE - LR	\$	576.35
2635700	THERMOMETER -10 TO 225 C 405MM	\$	238.97
2635701	THERMOMETER -10 TO 260 C 405MM	\$	281.37
2635702	THERMOMETER -20 TO 110 C 305MM	\$	175.95
2635704	THERMOMETER 0 TO 300 F 305MM	\$	161.69
2636166	PIPET, VOL CLASS A CERT 2ML PK/6	\$	379.66
2636168	ff PIPET, VOL CLASS A CERT 10ML PK/6	\$	537.69
2636170	ff PIPET, VOL CLASS A CERT 25ML PK/6	\$	695.72
2636202	LOW TEMPERATURE BOD INCUBATOR, SHELDON MODEL SRI6-2	\$	19,965.79
2636240	CYLINDER, GRAD CLASS A 25ML	\$	389.29
2636241	CYLINDER, GRAD CLASS A 50ML	\$	346.90
2636242	CYLINDER, GRAD CLASS A 100ML	\$	362.31
2636246	CYLINDER, GRAD CLASS A 250ML	\$	418.20
2636249	CYLINDER, GRAD CLASS A 500ML	\$	462.53
2636253	CYLINDER, GRAD CLASS A 1000ML	\$	601.29
2636300	HRI3P-HACH B.O.D. 120V INCUBATOR 4.65CFT PELTIER	\$	14,608.18
2636302	HRI3P-2 240V HACH B.O.D. INCUBATOR 4.65CFT PELTIER	\$	14,608.18
2636340	CYLINDER, GRAD MIX CLASS A 25ML	\$	358.46
2636341	CYLINDER, GRAD MIX CLASS A 50ML	\$	373.88
2636342	CYLINDER, GRAD MIX CLASS A 100ML	\$	398.93
2636538	ff BURET, CERT PTFE PLUG 10ML CLASS A	\$	869.17
2636540	BURET, CERT PTFE PLUG 25ML CLASS A	\$	724.63
2636541	BURET, CERT PTFE PLUG 50ML CLASS A	\$	747.75
2636642	FLASK, VOL CERT CLASS A 100ML	\$	210.06
2636646	ff FLASK, VOL CERT CLASS A 250ML	\$	238.97
2636649	FLASK, VOL CERT CLASS A 500ML	\$	300.64
2638300	ABSORBENT PAD, PP 11X13X3/8 in	\$	9.06
2640800	FILTER, MEMB PLAIN 47MM PK/100	\$	528.05



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

2641549	WATER, ORGANIC FREE 500ML	\$	184.29
2651155	AMINO ACID F RGT PWD, 55G	\$	108.04
2651155-LM	AMINO ACID F RGT PWD, 55G	\$	108.04
2651255	ASCORBIC ACID PWD, 48G	\$	90.87
2651600	REAGENT SET, COBALT/NICKEL 10ML	\$	527.57
2651600-LM	REAGENT SET, COBALT/NICKEL 10ML	\$	527.57
2651700	REAGENT SET, LR MANGANESE 10ML	\$	196.94
2651700-LM	REAGENT SET, LR MANGANESE 10ML	\$	196.94
2653199	AMMONIA CYANURATE RGT PP PK/100	\$	165.50
2653199-LM	AMMONIA CYANURATE RGT PP PK/100	\$	165.50
2653299	AMMONIA SALICYLATE PP PK/100	\$	245.72
2653299-LM	AMMONIA SALICYLATE PP PK/100	\$	245.72
2655200	ADAPTER, 90 DEGREE W/HOSE ADPTR	\$	134.90
2655342	FLASK, ERLLENMEYER 100ML	\$	155.91
2655438	DISTILLING RECIEVER, 10ML	\$	335.33
2655500	MERCURY VAPOR ABSORBER	\$	450.96
2655510	MERCURY ABSORBER COLUMN	\$	1,092.72
2655800	MERCURY SCRUBBER	\$	47.41
2655900	STOPPER, PENNY HEAD	\$	153.02
2656200	CLAMP, EXTENSION THERMOMETER	\$	163.62
2656300	RING, SUPPORT PLAIN 4.6"OD / 4.0"ID	\$	65.52
2656600	FILTER UNIT, ANAL 0.45UM PK/12	\$	398.93
2656822	BOTTLE, B ROUND LDPE 59ML PK/6	\$	20.62
2657200	SOIL SCOOP SET, PLASTIC 1,2,5,10	\$	294.86
2657201	SOIL SCOOP, 1 GRAM, 0.85CC	\$	73.62
2657202	si SOIL SCOOP, 2 GRAM, 1.70CC	\$	73.62
2657205	SOIL SCOOP, 5 GRAM, 4.25CC	\$	73.62
2657210	SOIL SCOOP, 10 GRAM, 8.50CC	\$	73.62
2657364	PH STORAGE SOLN PP PK/20	\$	116.52
2657364-LM	PH STORAGE SOLN PP PK/20	\$	116.52
2657512	PHENOL RED SPEC GRADE IND 50ML	\$	48.23
2657512-LM	PHENOL RED SPEC GRADE IND 50ML	\$	48.23
2657888	RAPID LIQUID ANALYSIS MANUAL	\$	73.81
2658200	CABLE ASSEMBLY, BODTRAK PRINTER	\$	551.18
2658300	KTO: MERCURY RGT SET, COLD VAPOR EXT	\$	661.27
2658448	HGEX3 ALKALINE RGT PP PK/25	\$	71.19
2658548	HGEX4 INDICATOR PP PK/25	\$	56.37
2658548-LM	HGEX4 INDICATOR PP PK/25	\$	56.37
2658625	Ammonium Iron(III) Citrate (Green), 125 g	\$	350.51
2658636	HGEX5 HYDROXIDE SOLN 10ML SCDB	\$	64.68
2658636-LM	HGEX5 HYDROXIDE SOLN 10ML SCDB	\$	64.68
2658748	HGEX6 COMPLEXING RGT PP PK/25	\$	65.93
2658748-LM	HGEX6 COMPLEXING RGT PP PK/25	\$	65.93
2658825	HGEX RGT A 20 ML AMPULES PK/25	\$	278.24
2658825-LM	HGEX RGT A 20 ML AMPULES PK/25	\$	278.24
2658949	HGEX REAGENT B SOLN 500ML	\$	167.47
2658949-LM	HGEX REAGENT B SOLN 500ML	\$	167.47



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

2659059	HGEX RGT C SOLN 55ML	\$	57.98
2659059-LM	HGEX RGT C SOLN 55ML	\$	57.98
2659105	STABLCAL AMPULE KIT, 2100A	\$	968.42
2659105-LM	STABLCAL AMPULE KIT, 2100A	\$	968.42
2659110	STABLCAL CALIB SET, 2100A 100ML	\$	646.82
2659110-LM	STABLCAL CALIB SET, 2100A 100ML	\$	646.82
2659400	STABLCAL CALIB KIT, 2100P	\$	984.68
2659400-LM	STABLCAL CALIB KIT, 2100P	\$	984.68
2659405	STABLCAL AMPULE KIT, 2100P	\$	672.11
2659405-LM	STABLCAL AMPULE KIT, 2100P	\$	672.11
2659410	STABLCAL CALIB SET, 2100P	\$	576.35
2659410-LM	STABLCAL CALIB SET, 2100P	\$	576.35
2659500	STABLCAL CALIB KIT, 0 to 7500 NTU	\$	1,443.59
2659500-LM	STABLCAL CALIB KIT, 0 to 7500 NTU	\$	1,443.59
2659505	STABLCAL AMPULE KIT, 0-7500 NTU	\$	1,116.57
2659505-LM	STABLCAL AMPULE KIT, 0-7500 NTU	\$	1,116.57
2659510	STABLCAL CAL SET, 0 to 7500 NTU	\$	870.85
2659510-LM	STABLCAL CAL SET, 0 to 7500 NTU	\$	870.85
2659600	STABLCAL CALIB SET, 1720 SERIES	\$	2,023.56
2659600-LM	STABLCAL CALIB SET, 1720 SERIES	\$	2,023.56
2659701	nn STABLCAL SOLN, <0.1 NTU AMPULE	\$	139.29
2659742	STABLCAL SOLN, <0.1 NTU 100ML	\$	97.56
2659742-LM	STABLCAL SOLN, <0.1 NTU 100ML	\$	97.56
2659749	STABLCAL SOLN, <0.1 NTU 500ML	\$	169.11
2659749-LM	STABLCAL SOLN, <0.1 NTU 500ML	\$	169.11
2659753	STABLCAL SOLN, <0.1 NTU 1000ML	\$	225.84
2659753-LM	STABLCAL SOLN, <0.1 NTU 1000ML	\$	225.84
2659801	nn STABLCAL STD AMPULE, 1.0 NTU	\$	167.67
2659802	STABLCAL STD, 1.0 NTU AMPULE PKG	\$	173.98
2659842	STABLCAL STD, 1.0 NTU 100ML	\$	161.33
2659842-LM	STABLCAL STD, 1.0 NTU 100ML	\$	161.33
2659849	STABLCAL STD, 1.0 NTU 500ML	\$	285.47
2659849-LM	STABLCAL STD, 1.0 NTU 500ML	\$	285.47
2659853	STABLCAL STD, 1.0 NTU 1000ML	\$	411.94
2659853-LM	STABLCAL STD, 1.0 NTU 1000ML	\$	411.94
2659901	nn STABLCAL STD AMPULE, 10 NTU	\$	173.98
2659902	STABLCAL VERIF. AMPULE, 10 NTU	\$	97.93
2659942	STABLCAL STD, 10 NTU 100ML	\$	161.33
2659942-LM	STABLCAL STD, 10 NTU 100ML	\$	161.33
2659949	STABLCAL STD, 10 NTU 500ML	\$	285.47
2659949-LM	STABLCAL STD, 10 NTU 500ML	\$	285.47
2659953	STABLCAL STD, 10 NTU 1000ML	\$	478.79
2660142	STABLCAL STD, 20NTU 100ML	\$	169.11
2660142-LM	STABLCAL STD, 20NTU 100ML	\$	169.11
2660149	STABLCAL STD, 20 NTU 500ML	\$	285.47
2660149-LM	STABLCAL STD, 20 NTU 500ML	\$	285.47
2660153	STABLCAL STD, 20 NTU 1000ML	\$	444.46



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

2660153-LM	STABLCAL STD, 20 NTU 1000ML	\$	444.46
2660201	nn STABLCAL STD AMPULE, 100 NTU	\$	139.29
2660242	STABLCAL STD, 100 NTU 100ML	\$	169.11
2660242-LM	STABLCAL STD, 100 NTU 100ML	\$	169.11
2660249	STABLCAL STD, 100 NTU 500ML	\$	285.47
2660249-LM	STABLCAL STD, 100 NTU 500ML	\$	285.47
2660253	STABLCAL STD, 100 NTU 1000ML	\$	478.79
2660442	STABLCAL STD, 200NTU 100ML	\$	169.11
2660442-LM	STABLCAL STD, 200NTU 100ML	\$	169.11
2660449	STABLCAL STD, 200 NTU 500ML	\$	285.47
2660449-LM	STABLCAL STD, 200 NTU 500ML	\$	285.47
2660453	STABLCAL STD, 200 NTU 1000ML	\$	478.79
2660542	STABLCAL STD, 800 NTU 100ML	\$	169.11
2660542-LM	STABLCAL STD, 800 NTU 100ML	\$	169.11
2660549	STABLCAL STD, 800 NTU 500ML	\$	285.47
2660549-LM	STABLCAL STD, 800 NTU 500ML	\$	285.47
2660553	STABLCAL STD, 800 NTU 1000ML	\$	478.79
2660642	STABLCAL STD, 1000 NTU 100ML	\$	169.11
2660642-LM	STABLCAL STD, 1000 NTU 100ML	\$	169.11
2660649	STABLCAL STD, 1000 NTU 500ML	\$	285.47
2660649-LM	STABLCAL STD, 1000 NTU 500ML	\$	285.47
2660653	STABLCAL STD, 1000 NTU 1000ML	\$	478.79
2662100	STABLCAL CALIB KIT, 0 to 4000 NTU	\$	1,217.75
2662100-LM	STABLCAL CALIB KIT, 0 to 4000 NTU	\$	1,217.75
2662105	STABLCAL AMPULE KIT, 0-4000 NTU	\$	858.21
2662105-LM	STABLCAL AMPULE KIT, 0-4000 NTU	\$	858.21
2662110	STABLCAL CAL SET, 0 to 4000 NTU	\$	740.77
2662110-LM	STABLCAL CAL SET, 0 to 4000 NTU	\$	740.77
2662200	GAS WASHING BOTTLE, 1200ML	\$	1,834.69
2662300	STOPPER, F/GAS WASHING BOTTLE	\$	151.09
2665034	SILICA GEL, 100-200 MESH, 500GM	\$	1,476.11
2665902	SAMPLE CELL SET, 1 INCH SQ (2PK)	\$	682.23
2665904	SAMPLE CELL SET, 1 INCH SQ (4PK)	\$	1,354.82
2665908	SAMPLE CELL SET, 1 INCH SQ (8PK)	\$	2,524.63
2667200	TEST KIT, IR-18C FE2+ 0-10MGL	\$	238.97
2667200-LM	TEST KIT, IR-18C FE2+ 0-10MGL	\$	238.97
2668000	REAGENT SET, NITROGEN AMMONIA	\$	388.45
2668000-LM	REAGENT SET, NITROGEN AMMONIA	\$	388.45
2669748	IRON TEST RGT PK/25	\$	17.15
2669748-LM	IRON TEST RGT PK/25	\$	17.15
2670500	GROUNDWATER SAMPLING CAPSULE	\$	172.10
2670600	FIELD PACK, FOR PORT UV LAMP	\$	1,015.63
2670700	PWR SUPPLY F/UV LAMPS 115V/60HZ	\$	790.15
2670702	PWR SUPPLY F/UV LAMPS 220V/50HZ	\$	1,056.11
2670800	rr UV LAMP, SHORTWAVE, PENCIL TYPE	\$	506.85
2671000	rr UV LAMP, SHORTWAVE, PENCIL TYPE	\$	481.80
2671745	Total Nitrogen Test N Tube Hydroxide Reagent	\$	484.21

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

2671745-LM	TNT, TN HYDROXIDE RGT PK/50	\$	484.21
2671846	Total Nitrogen Persulfate Reagent Pillows pack 50	\$	53.29
2671846-LM	TN PERSULFATE PP PK/50	\$	53.29
2671946	Bisulfite Reagent A Powder Pillows PK 50	\$	53.29
2671946-LM	TN RGT A BISULFITE PP PK/50	\$	53.29
2672046	ee TN RGT B INDICATOR PP PK/50	\$	41.56
2672046-LM	ee TN RGT B INDICATOR PP PK/50	\$	41.56
2672145	ee nn Total Nitrogen Test N Tube Acid Reagent (Reagent C)	\$	263.79
2672145-LM	ee nn TNT, TN RGT C ACID SOLN 50 TESTS	\$	263.79
2672245	ee Total Nitrogen Reagent Set, Low Range Test N Tube	\$	493.24
2672245-LM	ee RGT SET, TNT LR TOTAL NITROGEN	\$	493.24
2672629	COD STANDARD SOLN, 800MG/L 200ML	\$	115.08
2672629-LM	COD STANDARD SOLN, 800MG/L 200ML	\$	115.08
2674400	APPARATUS SET, COLD VAPR MERCURY	\$	5,511.79
2676100	THERMOMETER, POCKET -30 TO 120F	\$	90.58
2676200	THERMOMETER, POCKET 0 TO 220F	\$	114.48
2676300	THERMOMETER, POCKET -5 TO 45C	\$	110.43
2676400	THERMOMETER, POCKET -10 TO 110C	\$	139.14
2676500	CENTRIFUGE, 8 X 15ML 115VAC	\$	4,471.10
2676600	FILTER, GLASS FBR 47MM PK/100	\$	599.36
2676700	FILTER, EXTRA THICK, 25MM PK/100	\$	167.09
2676800	ff PIPET FILLER, DROPPER SIZE PK/12	\$	66.49
2677925	ff ACCUVAC AMPULE, FOR BLANKS PK25	\$	88.89
2678300	FILTER HOLDER, MODIFIED TO 5MM	\$	177.30
2678500	KTO RGT SET, ULR SILICA RAPID LIQUID	\$	355.93
2678500-LM	KTO:RGT SET, ULR SILICA RAPID LIQUID	\$	355.93
2678600	RGT SET, PHOSPHORUS RAPID LIQUID	\$	252.95
2683318	MEMBRANE, VENTING PTFE PK/18	\$	1,029.12
2683533	BOTTLE, B ROUND LDPE 29ML PK/6	\$	19.46
2683706	PLUG, DROPPER 20MM PK/6	\$	14.65
2683806	CAP, NASAL 20MM WHITE PK/6	\$	12.33
2685300	INSTRUCTIONS, SAFETY & DISPOSAL	\$	12.72
2686100	ENVIROCHECK FILTER F/CRYPTOSPORIDIUM	\$	562.74
2686832	SODIUM THIOSULFATE 2.12N 100ML	\$	125.56
2686901	SODIUM THIOSULFATE STD, 2.26N	\$	70.99
2687000	KTO: RGT SET, HR HYPOCHLORITE (BLEACH)	\$	381.22
2687000-LM	KTO: RGT SET, HR HYPOCHLORITE (BLEACH)	\$	381.22
2687100	HR HYPOCHLORITE (BLEACH) DT	\$	1,021.42
2687100-LM	HR HYPOCHLORITE (BLEACH) DT	\$	1,021.42
2687200	TEST KIT, HYPOCHLORITE HR (BLEACH)	\$	402.78
2687200-LM	TEST KIT, HYPOCHLORITE HR (BLEACH)	\$	402.78
2687300	OIL CLOTH, WHITE 9" SQUARE	\$	70.34
2689400	RGT SET, CEL/820 WATER COND LAB	\$	1,223.17
2689700	APP SET, CEL BASIC DW LAB	\$	506.85
2689800	RGT SET, CEL BASIC DW LAB	\$	2,565.59
2689900	APP SET, CEL BASIC WASTE WTR	\$	813.28
2690000	RGT SET, CEL BASIC WASTE WTR	\$	1,208.72

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

2690100	APP SET, CEL ADV DW/GB VERIF	\$	1,643.90
2690200	REAGENT SET, CEL ADV DW	\$	3,505.10
2690300	APP SET, CEL ADV WASTE WTR	\$	1,605.36
2690400	RGT SET, CEL AD WASTE WTR	\$	3,252.15
2690400-LM	RGT SET, CEL AD WASTE WTR	\$	3,252.15
2690500	APP SET, CEL ADV PORT LAB	\$	1,574.52
2690600	RGT SET, CEL ADV PORT LAB	\$	4,083.26
2690600-LM	RGT SET, CEL ADV PORT LAB	\$	4,083.26
2690700	APP SET, CEL ENV WTR QUAL LAB	\$	1,923.35
2690800	RGT SET, CEL ENV WTR QUAL LAB	\$	3,288.29
2690800-LM	RGT SET, CEL ENV WTR QUAL LAB	\$	3,288.29
2691100	RGT SET, MEL POTABLE WTR	\$	2,177.74
2691200	APPARATUS SET, MEL/850	\$	1,707.50
2691400	CASE, ADV PORT LAB	\$	928.91
2691600	APPARATUS SET, AQUACULTURE LAB	\$	1,574.52
2691700	RGT SET, CEL AQUACULTURE LAB	\$	2,710.13
2694400	HUMIDITY/TEMP MONITOR	\$	283.30
2694500	BAROMETER, 4" DIAL READS IN HPA	\$	157.64
2695753	APA RGT 2 HARDNESS BFR SOLN, 1L	\$	89.06
2695753-LM	APA RGT 2 HARDNESS BFR SOLN, 1L	\$	89.06
2695853	APA RGT 1 HARDNESS IND SOLN, 1L	\$	245.72
2696253	APA STD 1, LR HDNS, 0 UG/L 1L	\$	109.30
2696353	APA STD 2, LR HDNS 5 MG/L 1L	\$	109.30
2696553	APA RGT 1 ALKLNITY ACID SLN 1L	\$	345.09
2696553-LM	APA RGT 1 ALKLNITY ACID SLN 1L	\$	345.09
2696653	APA RGT 2 ALKLNITY IND SOLN 1L	\$	193.32
2696753	APA STD 1 ALKALINITY ZERO STD 1L	\$	109.30
2697053	WASH SOLN, APA 6000 1L	\$	84.19
2697453	APA CLEANING SOLN, SURFACTANT	\$	123.03
2697942	STABLCAL CERT STD, 0.30 NTU 100ML	\$	161.33
2697942-LM	STABLCAL CERT STD, 0.30 NTU 100ML	\$	161.33
2697953	STABLCAL CERT STD, 0.30 NTU 1L	\$	482.40
2697953-LM	STABLCAL CERT STD, 0.30 NTU 1L	\$	482.40
2698042	STABLCAL CERT STD, 0.50 NTU 100ML	\$	161.33
2698042-LM	STABLCAL CERT STD, 0.50 NTU 100ML	\$	161.33
2698049	STABLCAL CERT STD, 0.50 NTU, 500ML	\$	204.16
2698049-LM	STABLCAL CERT STD, 0.50 NTU, 500ML	\$	204.16
2698053	nn STABLCAL CERT STD, 0.50 NTU 1L	\$	411.94
2698053-LM	nn STABLCAL CERT STD, 0.50 NTU 1L	\$	411.94
2700349	PAO TITRANT, 0.00564N 500ML	\$	77.50
2700653	EDTA TITRANT, 0.01M 1L	\$	75.33
2700742	BROMOCRESOL GREEN IND, 100ML	\$	87.09
2700832	SODIUM THIOSULFATE, 0.1N 100ML	\$	86.72
2701349	MUREXIDE INDICATOR PP, PK/100	\$	87.09
2701349-LM	MUREXIDE INDICATOR PP, PK/100	\$	87.09
27032	SILICA 1 RGT SOLN, 100ML MDB	\$	57.82
27032-LM	SILICA 1 RGT SOLN, 100ML MDB	\$	57.82

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

2703810	ff LAB STIRRER, PROGRAMMABLE 6-PADDLE JAR TESTER NB	\$	12,642.43
2704303	SOLVENT, FOR HI TEMP CITRIC ACID	\$	103.71
2704303-LM	SOLVENT, FOR HI TEMP CITRIC ACID	\$	103.71
2704420	SURFACTANT, F/HI TMP CITRIC ACID	\$	110.02
2704420-LM	SURFACTANT, F/HI TMP CITRIC ACID	\$	110.02
2704500	rr UV LAMP & POWER SUPPLY 115V	\$	1,194.86
2704502	rr UV LAMP & POWER SUPPLY 230V	\$	1,270.02
2704503	rr UV LAMP & P-SUPPLY, 230V EURO PLUG	\$	1,225.70
2704700	ELECTRODE WASHER, 500ML	\$	487.58
27053	SILICA 1 RGT SOLN, 1000ML	\$	154.47
2708000	SPEC COLOR STD, OZONE 0-0.75MG/L	\$	581.77
2709600	BODTRAK TUBING, EXTERNAL	\$	614.78
2709700	BODTRAK TUBING, INTERNAL	\$	335.33
2712000	TEST KIT, STREAM SURVEY	\$	1,574.52
2712500	SPEC COLOR STD KIT, FLUORIDE	\$	581.77
2714045	nn TNT, TN HYDROXIDE HR RGT PK/50	\$	263.79
2714045-LM	nn TNT, TN HYDROXIDE HR RGT PK/50	\$	263.79
2714100	ee RGT SET, TNT HR TOTAL NITROGEN	\$	518.54
2714100-LM	ee RGT SET, TNT HR TOTAL NITROGEN	\$	518.54
2714349	SALINITY STD SLN, 35.0 PPT 500ML	\$	141.46
2714600	ULR STABLCAL VERIFICATION KIT	\$	433.62
27169	SILICA 3 REAGENT PWD PLWS PK/100	\$	120.33
27169-LM	SILICA 3 REAGENT PWD PLWS PK/100	\$	120.33
2718	FILTER,AIR,10 MICRON	\$	126.81
2722000	CASE, SOFT 11.5H X 2.5D X 11.5W	\$	556.96
2723342	STABLCAL CERT STD, 0.10 NTU 100ML	\$	161.33
2723342-LM	STABLCAL CERT STD, 0.10 NTU 100ML	\$	161.33
2723353	STABLCAL CERT STD, 0.10 NTU 1L	\$	411.94
2723353-LM	STABLCAL CERT STD, 0.10 NTU 1L	\$	411.94
2723501	REPLACEMENT ELECTRODE, PHTESTR2	\$	283.30
2724106	R2A-I-MF AGAR 14ML PK/6	\$	123.15
27242	Deionized (Demineralized) Water 100 mL	\$	31.08
27248	WATER, DEIONIZED, 500ML	\$	40.28
27249	WATER, DEIONIZED 500ML (DEMINRLZD)	\$	39.19
27253	WATER, DEIONIZED 1L (DEMINERALIZED)	\$	31.44
27256	WATER, DEIONIZED 4L (DEMINERALIZED)	\$	95.20
2728800	BAG, FOIL 5 X 9.5 W/ZIPPER	\$	2.70
2742545	TNT, REACTIVE PHOSPHATE 50 TESTS	\$	166.94
2742545-LM	TNT, REACTIVE PHOSPHATE 50 TESTS	\$	166.94
2742645	TNT, TOTAL PHOSPHORUS 50 TESTS	\$	229.46
2742645-LM	TNT, TOTAL PHOSPHORUS 50 TESTS	\$	229.46
2742745	TNT, PHOSPHATE TOTAL & ACID HYDROLYZ	\$	245.72
2742745-LM	TNT, PHOSPHATE TOTAL & ACID HYDROLYZ	\$	245.72
2743042	SODIUM HYDROXIDE, 1.54N 100ML	\$	46.06
2743042-LM	SODIUM HYDROXIDE, 1.54N 100ML	\$	46.06
2743400	SAMPLE CELL, 1CM SQ PS PK/100	\$	208.14
2743811	FLUORIDE STD, 0.5MG/L 473ML	\$	84.19

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

2744850	TEST STRIP, TOT ALKALINITY BTL/50	\$	45.87
2744940	TEST STRIP, LR CHLORIDE BTL/40 (QT)	\$	183.85
2745050	TEST STRIP, CHLORINE BTL/50	\$	66.49
2745125	TEST STRIP, COPPER 0-3 PPM BTL/25	\$	88.27
2745250	TEST STRIP, TOTAL HARDNESS BTL/50	\$	45.87
2745325	TEST STRIP, IRON KIT W/ BTL/25 TESTS	\$	88.27
2745325-LM	TEST STRIP, IRON KIT W/ BTL/25 TESTS	\$	88.27
2745348	TEST STRIP IRON BTL/25	\$	64.95
2745425	TEST STRIP, NITRATE/NITRITE BTL/25	\$	80.36
2745650	TEST STRIP, PH LOW RANGE BTL/50	\$	45.87
2746300	ann STABLCAL STD, 40 NTU	\$	2,095.83
2746353	STABLCAL STD, 40NTU 1L	\$	444.46
2746353-LM	STABLCAL STD, 40NTU 1L	\$	444.46
2746356	STABLCAL STD, 40NTU 3.78L	\$	1,250.27
2749125	INOCULATING LOOP, PS 10UL PK/10	\$	17.15
2749220	DROPPER, POLY 3ML DISP PK/20	\$	18.50
2749320	DROPPER, POLY 5ML PK/20	\$	18.50
2750701	DROPPER, PE TRANSFER 15ML	\$	240.90
2751000	TEST STRIP, PH 7.5-14 PK/100	\$	100.09
2751200	THERMOMETER, POCKET DIAL 0-220F	\$	151.09
2751340	TEST STRIP, HR CHLORIDE BTL/40 (QT)	\$	183.85
2752300	FILTER PAPER, MIXED CELLULOSE	\$	375.80
2753200	STIR BAR, EGG SHAPE 3/4 X 3/8 INCH	\$	35.65
2753929	TRihalOMETHANES RGT 1, SCDB	\$	430.01
2753929-LM	TRihalOMETHANES RGT 1, SCDB	\$	430.01
2754048	TRihalOMETHANES RGT 2 SOLN 330ML	\$	475.18
2754048-LM	TRihalOMETHANES RGT 2 SOLN 330ML	\$	475.18
2754142	TRihalOMETHANES RGT 3 SOLN	\$	475.18
2754142-LM	TRihalOMETHANES RGT 3 SOLN	\$	475.18
27549	MANGANOUS SULFATE SOLN, 500ML	\$	64.68
2755250	TEST STRIP, 5-IN-1 WATER QUAL BTL/50	\$	72.46
2755325	TEST STRIP, AMMONIA SET W/ BTL/25	\$	88.27
2755953	APA RGT 1 COPPER BUFFER SOLN, 1L	\$	143.46
2756053	APA RGT 2 COPPER IND SOLN, 1L	\$	113.63
2756253	APA STD 1 COPPER, 0 MG/L 1L	\$	196.94
2756353	APA STD 2 (LR) COPPER 1MG/L 1L	\$	271.01
2756453	APA STD 2 (HR) COPPER 10MG/L 1L	\$	271.01
2756526	pH STORAGE SOLUTION, 50ML	\$	41.56
2756549	PH STORAGE SOLN, 500ML	\$	119.78
2756559	PH STORAGE SOLN, 50ML, 3.0 MOLAR	\$	59.98
2756699	TRihalOMETHANES RGT 4 PP PK/100	\$	888.92
2756707	THM STD, 10PPM AS CHLOROFORM PK/7	\$	334.25
2757150	TEST STRIP, PHOSPHATE BTL/50	\$	90.58
2758100	JAR, WM PP 500ML W/CAP	\$	61.86
2758101	JAR, WM PP 500ML W/CAP PKGD	\$	62.63
2758900	THERMOSTAT, HIGH LIMIT	\$	470.24
2759123	DO ELECTRODE FILL SOLN, 30 ML	\$	40.28



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

2759126	ELECT FILL SOLN, DO 50 ML SCDB	\$	52.76
2759126-LM	ELECT FILL SOLN, DO 50 ML SCDB	\$	52.76
2760345	TOC, LR TNT RGT SET	\$	1,524.90
2760345-LM	TOC, LR TNT RGT SET	\$	1,524.90
2760445	TOC, HR TNT RGT SET	\$	1,524.90
2760445-LM	TOC, HR TNT RGT SET	\$	1,524.90
2762133	GLYCINE RGT, 10% SOLN, 29ML SCDB	\$	73.88
2763900	DR/CHECK ABSORBANCE STANDARD KIT	\$	614.30
2764700	EVAPORATING DISH, 125MM X 65MM	\$	146.47
2767245	TNT KIT TOTAL PHOSPHATE HR	\$	339.67
2767245-LM	TNT KIT TOTAL PHOSPHATE HR	\$	339.67
2767345	TNT KIT, REACTIVE PHOSPHATE HR	\$	245.72
2767345-LM	TNT KIT, REACTIVE PHOSPHATE HR	\$	245.72
2768549	BUFFER SOLN 0.3,1,2,5MG/L HDNS 500ML	\$	163.69
2768549-LM	BUFFER SOLN 0.3,1,2,5MG/L HDNS 500ML	\$	163.69
2768649	BUFFER SOLN, 10MG/L HARDNESS	\$	163.69
2768649-LM	BUFFER SOLN, 10MG/L HARDNESS	\$	163.69
2768749	BUFFER SOLN, 20MG/L HARDNESS	\$	163.69
2768749-LM	BUFFER SOLN, 20MG/L HARDNESS	\$	163.69
2768849	BUFFER SOLN, 50MG/L HARDNESS 500ML	\$	163.69
2768849-LM	BUFFER SOLN, 50MG/L HARDNESS 500ML	\$	163.69
2768949	BUFFER SOLN, 100MG/L HARDNESS 500ML	\$	163.69
2769049	INDICATOR SOLN, 1MG/L HARDNESS	\$	136.95
2769049-LM	INDICATOR SOLN, 1MG/L HARDNESS	\$	136.95
2769149	INDICATOR SOLN, 2MG/L HARDNESS	\$	162.97
2769149-LM	INDICATOR SOLN, 2MG/L HARDNESS	\$	162.97
2769249	IND SOLN, 5-100 MG/L HARDNESS	\$	136.95
2769249-LM	IND SOLN, 5-100 MG/L HARDNESS	\$	136.95
2769820	SINGLET COMBO, PH 7.0 & 10.01, 10EA	\$	119.78
2769920	SINGLET COMBO, PH 4.01 & 7.0, 10EA	\$	119.78
276A3A1203	ASSY. MOUNTING BRACKET	\$	458.67
276A3M1210	vv PACKAGE OF 12 DISPOSABLE	\$	7.13
276C3A1200	aa SWIVEL, PIVOT, & CLAMP AS	\$	1,210.28
276G1200	ASSEMBLY-SWIVEL, & CLAMP	\$	1,616.92
276G1206	ASSEMBLY-SERVICE SUPPORT	\$	354.60
276G1226	KIT-MOUNTING HARDWARE	\$	225.48
276G1230	WASHER HEAD AIR BLAST CLE	\$	1,717.14
276G1233	ASSY-WASHER HEAD SUB-ASSY	\$	1,009.85
276M1210	PACKAGE OF 12 DISP.BAGS	\$	143.19
2770020	SINGLET, PH 4.01 PK/20	\$	119.78
2770120	SINGLET, PH 7.0 PK/20	\$	119.78
2770220	SINGLET, PH10.01 PK/20	\$	119.78
2770900	CL02 DPD PP/GLYCINE METHOD RGT SET	\$	155.02
2770900-LM	CL02 DPD PP/GLYCINE METHOD RGT SET	\$	155.02
2771320	SINGLET, COND 147US/CM, PK/20	\$	162.61
2771420	SINGLET, COND 1413US/CM, PK/20	\$	162.61
2771520	SINGLET, COND 12.88MS/CM, PK/20	\$	162.61

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

27749	ALKALINE IODIDE-AZIDE RGT 500ML	\$	95.40
2776000	CAP, IND SEAL WHITE PP 28-410	\$	1.16
2776353	RGT 1, AMMONIA/MONOCHLORAMINE 1L	\$	200.55
2776453	RGT 2, AMMONIA/MONOCHLORAMINE 1L	\$	80.22
2776553	RGT 3, AMMONIA/MONOCHLORAMINE 1L	\$	114.36
2776653	STD 1, AMMONIA/MONOCHLORAMINE 1L	\$	70.99
2776731	APA AMMONIA STD 2, 237 mL	\$	31.08
2776753	STD 2, AMMONIA/MONOCHLORAMINE 1L	\$	89.06
2779200	KTO: M-CB24 MF KIT, GLASS AMP 200 TESTS	\$	2,620.99
2779201	KTO: M-CB24 MF KIT, BROTH 200 TESTS	\$	1,852.04
2779202	KTO: M-CB24 MF KIT, BROTH 1000 TESTS	\$	7,997.88
2788400	nn STABLCAL CERT STD, 800MNTU BULK	\$	419.17
2788453	STABLCAL CERT STD 800 MNTU 1L	\$	482.40
2788453-LM	STABLCAL CERT STD 800 MNTU 1L	\$	482.40
2790166	TOC PERSULFATE PP, PK/50	\$	130.45
2790166-LM	TOC PERSULFATE PP, PK/50	\$	130.45
2790800	KTO: RGT SET, TRIHALOMETHANES PLUS	\$	2,204.24
2790800-LM	KTO: RGT SET, TRIHALOMETHANES PLUS	\$	2,204.24
2791505	TOC STD AMP, 1000MG/L C, 20ML PK/5	\$	157.19
2793053	vv RGT 2 BUFFER, APA HR HARDNESS 1L	\$	116.72
2793053-LM	vv RGT 2 BUFFER, APA HR HARDNESS 1L	\$	116.72
2793135	vv RGT 2 IND PWD, APA HR HDNS 0.5GM	\$	285.47
2793135-LM	vv RGT 2 IND PWD, APA HR HDNS 0.5GM	\$	285.47
2793253	STD 1, APA HR HARDNESS 0MG/L 1L	\$	116.72
2793353	STD 2, APA HR HDNS 1000MG/L 1L	\$	116.72
2793553	RGT 1, APA HR HARDNESS 1L	\$	131.34
2793553-LM	RGT 1, APA HR HARDNESS 1L	\$	131.34
2793600	KTO: RGT 2 KIT, APA HR HARDNESS	\$	314.37
2793600-LM	KTO: RGT 2 KIT, APA HR HARDNESS	\$	314.37
2793753	RGT 3, APA HR HARDNESS 1L	\$	99.54
2793828	SOFCHK HARDNESS TST STRIP pk/1000	\$	1,125.48
2793844	TEST STRIPS, SOFCHEK HDNS pk/250	\$	360.39
2793944	TST STRIPS, AC F/T CL1 pk/250	\$	695.72
2794005	VIAL, GLASS 40ML W/SEPTA CAP PK/5	\$	82.48
2794649	HARDNESS IND SOLN, 0.3MG/L, SP-510	\$	359.54
2794649-LM	HARDNESS IND SOLN, 0.3MG/L, SP-510	\$	359.54
2795000	ff PIPETTE TIPS, PK/1000 REPL F/2794900	\$	126.81
2795800	BEAKER, DISPOSABLE 50ML	\$	5.97
2796900	BALANCE, OHAUS MODEL YA501	\$	616.70
2797111	FLUORIDE STD, 5.0MG/L 473ML	\$	88.89
2797149	FLUORIDE STD SOLN, 5.0MG/L 500ML	\$	91.06
2797799	ARSENIC RGT #2 PP, PK/100	\$	112.02
2799400	SAFETYPETTE PIPETTE FILLER	\$	462.53
2799403	SILCONE PIPETTE HOLDER	\$	17.73
2799900	rr RGT SET, ARSENIC KIT REFILL	\$	443.26
2800000	rr TEST KIT, ARSENIC	\$	618.63
2800100	rr TEST STRIPS, ARSENIC BOTTLE/100	\$	208.14

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

2800200	REACTION VESSEL, ARSENIC KIT	\$	33.15
2800275	REACTION VESSEL, ARSENIC KIT, 75ML	\$	35.27
28014	MANVER 2 HARDNESS IND PWD, 113G	\$	136.95
28014-LM	MANVER 2 HARDNESS IND PWD, 113G	\$	136.95
28019BB	GLI pH/DO Sensor Extension Cable, 50ft, 1W1100, 6 Wire - USABB Only	\$	803.64
2802210	KTO: MONOCHLOR F REAGENT PK/1000	\$	1,703.77
2802246	MONOCHLORAMINE F RGT PP, PK/50	\$	107.49
2802246-LM	MONOCHLORAMINE F RGT PP, PK/50	\$	107.49
2802299	MONOCHLOR F REAGENT PK/100	\$	215.00
2802299-LM	MONOCHLOR F REAGENT PK/100	\$	215.00
2802300	DPD FREE BULK DISPENSER + POWDER	\$	272.82
2802300-LM	DPD FREE BULK DISPENSER + POWDER	\$	272.82
2802400	DPD TOTAL BULK DISPENSER + POWDER	\$	272.82
2802400-LM	DPD TOTAL BULK DISPENSER + POWDER	\$	272.82
2805145	ff RGT SET, HR MONOCHLORAMINE TNT	\$	395.68
2805145-LM	ff RGT SET, HR MONOCHLORAMINE TNT	\$	395.68
2805215	M-COLIBLUE24 AGAR PLATE, PK/15	\$	352.68
2805342	STABLCAL CERT STD 0.060NTU 100ML	\$	178.87
2805342-LM	STABLCAL CERT STD 0.060NTU 100ML	\$	178.87
2807710	FAN, CIRCULATING, 26162XX S/N >04 REPLACEMENT	\$	909.64
2808100	MICROPROC TEMP CONTROL 120V #1750573	\$	2,659.54
28101H	CALVER 2 CALCIUM INDICATOR, 454G	\$	169.11
28101H-LM	CALVER 2 CALCIUM INDICATOR, 454G	\$	169.11
28111H	CALVER 2, CALCIUM IND PWD 11.3KG	\$	2,818.53
2811200	ENVIROCHEK, HV FOR FINISHED WATER	\$	562.74
2811225	ENVIROCHEK HV FILTER, BULK PK/25	\$	12,835.15
2811415	M-HPC AGAR PLATES PK/15	\$	341.11
2811415-90	KTO: M-HPC AGAR PLATES PK/15x6	\$	1,738.33
28114H	CALVER 2 CALCIUM INDICATOR, 113G	\$	90.70
28114H-LM	CALVER 2 CALCIUM INDICATOR, 113G	\$	90.70
2811515	M-FC AGAR PLATES PK/15	\$	341.11
2811615	M-ENDO AGAR PLATES PK/15	\$	341.11
2811715	M-EI AGAR PLATES, PK/15	\$	341.11
2811815	M-TEC, MODIFIED AGAR PLATES, PK/15	\$	341.11
2812115	NUTRIENT AGAR W/MUG PLATES PK/15	\$	341.11
2812250	PSEUDO BROTH AMP, PLAS PK/50	\$	271.74
2812250-LM	PSEUDO BROTH AMP, PLAS PK/50	\$	271.74
2812350	R2A AMP, PLASTIC PK/50	\$	294.86
2812450	HPC BROTH AMP, PLAS PK/50	\$	314.13
2812650	TRYPTIC SOY BROTH AMP, PLASTIC PK/50	\$	294.86
2812811	RGT 1, TISAB SOLN 473ML	\$	122.30
2812811-LM	RGT 1, TISAB SOLN 473ML	\$	122.30
2813935	INDOXYL-BETA-D-GLUCOSIDE 2G	\$	1,817.35
2814215	R2A AGAR PLATES PK/15	\$	335.33
2815134	nn PROTEOSE PEPTONE #3 500G	\$	643.68
2815945	MR TOC RGT SET	\$	1,524.90

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

2815945-LM	MR TOC RGT SET	\$	1,524.90
2816900	KTO: CA610 FLUORIDE REAGENT SET	\$	278.24
2816900-LM	KTO: CA610 FLUORIDE REAGENT SET	\$	278.24
2818300	nn SYRINGE, 140CC W/CATHETER TIP	\$	48.76
282-101-0-1001N	RCDR,2-12PH,SNG SURF MNT,	\$	6,340.49
2822999	vv EZ ARSENIC RGT #1 PP PK/100 ROW	\$	66.85
2823099	vv EZ ARSENIC RGT #2 PP PK/100 ROW	\$	16.62
28232H	POTASSIUM HYDROXIDE 8N 100ML MDB	\$	43.00
28232H-LM	POTASSIUM HYDROXIDE 8N 100ML MDB	\$	43.00
2824349	AMMONIA ISA SOLN 500ML	\$	92.67
2824400	KTO: RGT SET, SERIES 5000 SILICA	\$	834.72
2824800	PUMP, VACUUM 1.2CFM 115V 60HZ	\$	2,409.00
2824801	PUMP, VACUUM 1.2CFM 230V 60HZ NAPLUG	\$	2,563.18
2824802	PUMP, VACUUM 1.2CFM 230V 50HZ EURO	\$	2,505.36
28249	POTASSIUM HYDROXIDE, 8N 500ML	\$	81.83
2824901	SERVICE KIT FOR 2824800	\$	240.90
28249-LM	POTASSIUM HYDROXIDE, 8N 500ML	\$	81.83
2825131	AMTAX STD SOLN 5MG/L 237ML	\$	178.84
2825253	PHOSPHAX RGT SOLN 1L	\$	191.52
2825254	REAGENT F/PHOSPHAX ANALYZER, 2000ML	\$	359.54
2825254-LM	REAGENT F/PHOSPHAX ANALYZER, 2000ML	\$	410.49
2825352	CLEANING SOLN F/PHOSPHAX, 1000ML	\$	180.68
2825352-LM	CLEANING SOL'N F/PHOSPHAX, 1000ML	\$	191.52
2825353	PHOSPHAX CLEANING SOLN 1L	\$	198.74
2825431	AMTAX EXPULSION SOLN 237ML	\$	210.06
2825553	AMTAX INDICATOR SOL. 0.2-12 & 0.2-30MG/L 1L	\$	173.83
2825653	AMTAX INDICATOR 5-120MG/L 1L	\$	184.05
2825753	AMTAX INDICATOR 50-1200MG/L 1L	\$	184.05
2825831	AMTAX STD SOLN 50MG/L 237ML	\$	213.92
2825854	50 MG/L NH4-N STD SOLN, 2000ML	\$	296.31
2825931	AMTAX STD SOLN 500MG/L 237ML	\$	213.92
2825954	500 MG/L NH4-N STD SOLN, 2000ML	\$	296.31
2826153	APA RGT 1 ALK ACID 0.080N 1L	\$	108.93
2826153-LM	APA RGT 1 ALK ACID 0.080N 1L	\$	108.93
2826253	APA STD 2 ALKALINITY 500MG/L 1L	\$	109.30
2829017	FLUORIDE ISA SOLN, 3.78L	\$	229.46
2829053	FLUORIDE ISA SOLN, 1L	\$	61.60
2830700	ff KTO: REAGENT SET, AMTAX 0.5-12.0MG/L	\$	562.74
2830800	ff KTO: REAGENT SET, AMTAX 5-120MG/L	\$	562.74
2830900	ff KTO: REAGENT SET, AMTAX 50-1200MG/L	\$	562.74
2831000	ff KTO: REAGENT SET, PHOSPHAX compact	\$	364.96
2831409	BART TESTER, APB PK/9	\$	481.80
2831500	FUNNEL ASSY, W/FILTER, DISH/PAD PK50	\$	1,125.48
2831600	MANIFOLD, FIELD VACUUM	\$	2,832.98
2833049	DW INORGANICS QC STD, 500ML	\$	146.35
2833149	WW INFLUENT INORG/MIX PARA STD 500ML	\$	294.50
2833249	WW EFFLUENT INORG/MIX PARA STD 500ML	\$	294.50

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

2833349	HR HARDNESS QC STD, 500ML	\$	176.15
2833349-LM	HR HARDNESS QC STD, 500ML	\$	176.15
2833449	LR HARDNESS QC STD, 500ML	\$	176.15
2833510	OXYGEN DEM MIXED QC STD 10ML PK/16	\$	204.16
2833649	METALS DW/HR QC STD, 500ML	\$	120.86
2833749	METALS DW/LR QC STD, 500ML	\$	120.86
2834100	KI BULK DISPENSER	\$	330.64
2834253	SODIUM STD, 100MG/L NA CL AS NA, 1L	\$	152.49
2834325	COD DIGEST VIAL, COD2 UHR PK/25	\$	172.00
2834325-LM	COD DIGEST VIAL, COD2 UHR PK/25	\$	172.00
2834453	DIISOPROPYLAMINE, 99% 1L	\$	285.47
2834453-LM	DIISOPROPYLAMINE, 99% 1L	\$	285.47
2835153	SODIUM STD, 10MG/L NA CL AS NA, 1L	\$	180.68
2835801	CHLORINE STD, 100-130MG/L 20ML AMP	\$	95.76
2835900	KTO: CL17 CALIB/VERF REFILL KIT	\$	373.88
2838300	SAMPLE CELLS, PK/6 (4X10ML & 2X25ML)	\$	152.83
2841700	3 M KCL SILVER SATD 28ML	\$	44.45
2841700-LM	3.0M KCL, SILVER SATURATED 28ML	\$	33.97
2841950	CAPSULE, GROUNDWATER 0.45UM PK/50	\$	5,781.60
2842100	FILTER, GL FBR XTHICK 1UM 47MM PK/100	\$	246.68
2842800	TEST KIT, CLD-2 CHLORINE DIOXIDE	\$	389.29
2842800-LM	TEST KIT, CLD-2 CHLORINE DIOXIDE	\$	389.29
2845900	REPLACEMENT SENSOR, OAKTON TDS & EC	\$	269.81
2850027	ff ANEMOMETER, MINI, ADVANCED	\$	1,032.98
2852800	ACCESSORY KIT, DO/CO	\$	52.61
2852900	PH ACCESSORY KIT	\$	147.24
28553	rr MERCURIC NITRATE, 0.0141N 1000ML	\$	112.02
2856300	KOHLRAUSCH 200ML CLASS A VOL FLASK	\$	291.01
28564026	KTO: 1.0M KNO3 SOLUTION, 59 ML	\$	66.85
2859100	PATHOSCREEN FIELD KIT	\$	194.65
2859100-LM	PATHOSCREEN FIELD KIT	\$	194.65
2859226	BLEACH SOLUTION, 5-6%	\$	31.61
2873200	SEQUENCER IV REAR HOUSING ENCLOSURE	\$	1,212.21
2873400	SEQUENCER IV FRONT HOUSING GASKET	\$	291.01
2873500	SEQUENCER IV FRONT HOUSING SCREWS	\$	58.59
2873800	SEQUENCER IV LCD	\$	3,931.49
2873900	SEQUENCER IV FUSE 2A FOR CONTACTS	\$	37.39
2874200	SEQUENCER IV RELAY BOARD	\$	8,614.58
2874400	SEQUENCER IV ANALOG INPUT TERMINAL	\$	225.48
2874500	SEQUENCER IV DISCRETE INPUT TERMINAL	\$	291.01
2874600	SEQUENCER IV VALVE OUTPUT TERMINAL	\$	420.13
2874700	SEQUENCER IV RS485 NETWORK TERMINAL	\$	137.99
2874800	SEQUENCER IV TRACK & HOLD BOARD	\$	14,126.38
2874900	SEQUENCER IV CPU BOARD	\$	15,860.86
2875100	SEQUENCER IV POWER CABLE	\$	89.04
2875200	SEQUENCER IV SIGNAL CABLE	\$	362.31
2875300	SEQUENCER IV MEMBRANE EXT. CABLE	\$	271.74

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

2875400	SEQUENCER IV ENCLOSURE BRACKET	\$	420.13
2875500	BAG, POLY ZIPPER 8 X 8 X .002	\$	0.58
2876453	APA ACIDIC SURFACTANT WASH, 1L	\$	89.25
2876453-LM	APA ACIDIC SURFACTANT WASH, 1L	\$	89.25
2877436	FA CHLORINATING SOLUTION 4ML SCDB	\$	164.41
2877436-LM	FA CHLORINATING SOLUTION 4ML SCDB	\$	164.41
2879200	FA REAGENT SET	\$	359.54
2879200-LM	FA REAGENT SET	\$	359.54
2881200	TOP STIRRER 4X4, CIMAREC DGT 115V	\$	1,391.44
2881300	TOP STIRRER 7X7, CIMAREC DIG 115V	\$	1,466.60
2881400	HOT PLATE 4X4, CIMAREC DIGTL, 115V	\$	1,129.34
2881500	HOT PLATE 7X7, CIMAREC DIGTL, 115V	\$	1,264.24
2881600	HOT PLATE/STIR 7X7, CIMAREC DIG 115V	\$	2,177.74
2881602	DIGITAL STIRRING/HOT PLATE, CIMAREC, 7X7, 230V	\$	2,698.08
2881900	KIT, CHAIN MOUNT, 316SS	\$	356.53
2882053	PH 7.0 BUFFER - LOW TEMP, 1L	\$	460.72
2882200	NITRIFICATION CONTROL KIT	\$	6,725.93
2882346	HARDNESS TREATMENT RGT PP, PK/50	\$	79.50
2883100	SCISSORS, 5" LENGTH, BLUNT TIP	\$	14.65
2883702	KTO: REPLACEMENT KTO FOR UNIVER 2	\$	182.48
2884000	RADALERT 100 RADIATION MONITOR	\$	2,389.73
2884200	INSPECTOR ALERT RADIATION MONITOR	\$	3,468.96
2890200	HACH AC HR FCHLORINE (100)	\$	90.58
2893300	SPEC COLOR STD, DPD CHLORINE - HR	\$	576.35
2893700	CAP, BOTTLE PP 38-439 PTFE LINER	\$	2.31
2894154	1 MG/L NH4-N STD SOLN, 2000ML	\$	314.37
2894154-LM	1 MG/L NH4-N STD SOLN, 2000ML	\$	314.37
2894246	CLEANING SOLN F/AMTAX ALNZR, 250ML	\$	359.54
2894246-LM	CLEANING SOL'N F/AMTAX ALNZR, 250ML	\$	359.54
2894354	10 MG/L NH4-N STD SOLN, 2000ML	\$	314.37
2894452	REAGENT F/AMTAX ANALYZER, 2500 ML	\$	359.54
2894452-LM	REAGENT F/AMTAX ANALYZER, 2500 ML	\$	359.54
2894536	ELECTROLYTE F/AMTAX ANALYZER, 11ML	\$	196.94
2894700	47mm, 0.45micron HYDROPHILIC POLY pk 50	\$	572.38
2895405	rr EZ COD RECYCLING KIT, 5-GAL	\$	1,000.22
2895405P	EZ COD RECYCLING, 5-GAL PICK-UP	\$	1,468.53
2895414	rr EZ COD RECYCLING KIT, 14-GAL	\$	1,884.80
2895420	EZ COD RECYCLING KIT, 20-GAL	\$	2,505.36
2895455	EZ COD RECYCLING KIT, 55 GAL	\$	6,379.03
2895805	REACTOR ADAPTER, 16MM TO 13MM, 5/PK	\$	61.28
2896742	STABLCAL REF SUSP I 3 NTU 100ML	\$	155.19
2896742-LM	STABLCAL REF SUSP I 3 NTU 100ML	\$	155.19
2896842	STABLCAL REF SUSP II 6 NTU 100ML	\$	374.00
2896842-LM	STABLCAL REF SUSP II 6 NTU 100ML	\$	374.00
2896942	STABLCAL REF SUSP III 18NTU 100ML	\$	374.00
2896942-LM	STABLCAL REF SUSP III 18NTU 100ML	\$	374.00
2897042	STABLCAL REF SUSP IV 30 NTU 100ML	\$	374.00

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

2897042-LM	STABLCAL REF SUSP IV 30 NTU 100ML	\$	374.00
2897100	StablCal Ref Suspension Amp Kit	\$	2,204.24
2897200	StablCal Ref Suspension Set	\$	1,373.13
2898700	BOD BOTTLE, PK/24 SERIALIZED (#1-24)	\$	1,032.98
2898800	BOD BOTTLE, PK/24 SERIALIZED, #25-48	\$	1,032.98
2898900	BOD BOTTLE, PK/24 SERIALIZED, #49-72	\$	1,032.98
2899300	TUBE, CULTURE, GLASS, PK/396, 16MM	\$	824.84
29028.0	Spacer, EPDM o-rings	\$	493.36
2907-0644-001	inn FLEXY EWON UNIT, FOR REMOTE ACCESS 4G	\$	5,974.32
29149	FLUORIDE STD SOLN, 1.0MG/L 500ML	\$	73.16
29153	FLUORIDE STD SOLN,1.0MG/L 1000ML	\$	102.43
2915900	FILTER MEMBRANE EZ-PAK .45 UM PK/600	\$	1,285.44
2916800	FILTER MEMBRAN EA-PAK .45 UM PK/600	\$	1,285.44
2916900	FILTER MEMBRANE EZ-PAK .22 UM PK/600	\$	1,528.27
2918553	25PPM TOTAL ALKALINITY STD 1L	\$	115.08
2918700	POLYSEED BOD INNOCULUM PK/50	\$	542.03
2918700-LM	POLYSEED BOD INNOCULUM PK/50	\$	542.03
2923200	RGT SET, DREL COMPL WATER QUAL	\$	3,179.88
2923200-LM	RGT SET, DREL COMPL WATER QUAL	\$	3,179.88
2923300	RGT SET, DREL INDUST WTR QUAL	\$	3,378.62
2923400	APPARATUS SET, DREL WATER QUAL	\$	2,409.00
2924050	PRY PLASTIC AMPULES PK/50	\$	262.10
29240500	MOTOR GEARBOX ASSEMBLY 110V, AF7000	\$	8,344.78
2925100	KTO: NUTRIENTS, ADV RGT SET TNT+	\$	856.40
2927202	BOTTLE, WASH, WIDE MOUTH 250ML PK/12	\$	178.65
2927203	BOTTLE, WASH, WIDE MOUTH 500ML PK/6	\$	108.50
2927204	BOTTLE, WASH, WIDE MOUTH 1L PK/3	\$	82.48
2927401	ff GLOVES, NITRILE, w/Neu-Thera SMALL PK/100	\$	84.80
2927402	ff GLOVES, NITRILE, w/Neu-Thera, MED PK/100	\$	84.80
2927404	ff GLOVES, NITRILE, w/Neu-Thera, X-LG PK/100	\$	84.80
2927701	GLASSES, SAFETY CLEAR UD LENS	\$	45.10
2927702	ff GLASSES, SAFETY CLEAR XTR LENS	\$	47.41
2927802	ff GLASSES, SAFETY SKYLINE CLEAR XTR	\$	43.55
2927901	GOGGLES, SPLASH CHEMICAL FUTURA	\$	50.69
2927902	GOGGLES, SPLASH CHEMICAL CLASSIC	\$	41.82
2928002	ff GOGGLES, SAFETY, FOGARD PLUS	\$	39.31
2928101	FACESHIELD	\$	91.54
2928202	BOTTLE, WASH, NALGENE BLUE PK/6	\$	162.66
2928203	ff BOTTLE, WASH, NALGENE GREEN PK/6	\$	162.66
2928204	BOTTLE, WASH, NALGENE RED PK/6	\$	162.66
2928301	CYLINDER, GRTD, CLS A DBL MTRC 50ML	\$	76.32
2928302	CYLINDER, GRTD, CLS A DBL MTRC 100ML	\$	90.58
2928303	CYLINDER, GRTD, CLS A DBL MTRC 250ML	\$	126.42
2928304	CYLINDER, GRTD, CLS A DBL MTRC 500ML	\$	162.66
2928305	CYLINDER, GRTD, CLS A DBL MTRC 1L	\$	252.46
2928401	CYLINDER, LOW-FORM DBL METRIC 50ML	\$	109.27
2928402	CYLINDER, LOW-FORM DBL METRIC 100ML	\$	126.42

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

2928501	ff MAGNETIC STIRRING BAR ARRESTER	\$	68.99
2928601	STIR BAR, MAGNETIC, "+" SHAPE 3/8"	\$	35.85
2928602	STIR BAR, MAGNETIC, "+" SHAPE 3/4"	\$	41.82
2928603	STIR BAR, MAGNETIC, "+" SHAPE 1"	\$	42.98
2928604	STIR BAR, MAGNETIC, "+" SHAPE 1 1/4"	\$	49.14
2928605	STIR BAR. MAGNETIC. "+" SHAPE 1 1/2"	\$	55.12
2929101	ff BOTTLE, WASH, RED 250 ML NARROW PK/6	\$	120.45
2929103	BOTTLE, WASH RED 500ML WIDE PK/6	\$	120.45
2929301	SAMPLER DIPPER 3 FT HANDLE	\$	164.78
2929302	SAMPLER DIPPER 6 FT HANDLE	\$	198.50
2929401	SAMPLER, EXTRA CUP, 500ML	\$	44.71
2929501	SAMPLER, WATER DIPPER 3'- 9' HANDLE	\$	533.83
2929601	BOTTLE BRUSH, SWIVEL HEAD	\$	191.18
2930003	BOTTLE, WRT-ON LBL 500ML WIDE PK/12	\$	169.98
2930004	BOTTLE, WRITE-ON LABEL 1L WIDE PK/6	\$	169.98
2930101	ff BOTTLE, WASH NARROW MTH 125ML PK/12	\$	178.65
2930102	ff BOTTLE, WASH NARROW MTH 250ML PK/12	\$	178.65
2930103	BOTTLE, WASH, NARROW MTH 500ML PK/12	\$	178.65
2930201	BOTTLE, WASH, AUTOCLAVE, 500ML, PK.6	\$	260.17
2930202	BOTTLE, WASH, AUTOCLAVE, 1-L, PK/4	\$	260.17
2932810	DISPOSABLE, 30X30 CM,196/BX	\$	20.43
2936100	FILTER, MEMBRANE GRID 47MM PK/150	\$	302.57
2936200	FILTER, STER 0.7UM 47MM PK/150	\$	599.36
2936300	PETRI DISH W/PAD 47MM PK/150	\$	275.59
2936500	FILTER, MEMB GRID 47MM PK/600 MILLIP	\$	1,165.96
2936600	FILTER, STER 0.7UM 47MM PK/600	\$	1,433.84
2936701	BALANCE, ANALYTICAL 80G X 0.1MG 100-240V 50/60HZ	\$	8,286.96
2936801	BALANCE, ANALYTICAL 120GX0.1MG 100-240V 50/60HZ	\$	9,134.93
2936901	BALANCE, ANALYTICAL 210GX0.1MG 100-240V 50/60HZ	\$	9,674.54
2937406	KTO: BOTTLE DROP 118ML/CAP, 6/PK	\$	29.39
2937506	KTO: BOTTLE DISP CAP PLUG 29ML, 6/PK	\$	41.72
2937606	KTO: BOTTLE DISP CAP PLUG 59ML, 6/PK	\$	42.88
2938200	ALCOHOL PREP PADS, STERILE, 200/PK	\$	48.76
2938300	CARRIER TRAY FOR BOTTLES & REAGENTS	\$	85.37
2938500	KTO: CUBITAINER, 20 L W/SPIGOT	\$	53.38
2939103	BOTTLE, WM, PP, AUTOCLV 125ML, PK/12	\$	91.35
2942100	TUBE, CULTURE 16X100MM NO PRTG T1	\$	6.75
2942700	MICROSCOPE, STEREO 10X W/8X-35X ZOOM	\$	5,685.24
2942800	MICROSCOPE, STEREO, 12.8X-56X ZOOM	\$	6,012.86
2943000	SAMPLER DIPPER 12 FT HANDLE	\$	533.83
2943100	aa BOTTLE, BOD, 300ML DISPOSABLE 100/BX	\$	422.06
2943101	BOTTLE, BOD, 300ML DISPOSABLE	\$	9.83
2943200	EPA 1664A SPE APPARATUS, STARTER KIT	\$	5,164.90
2943200NS	ff SPE STARTER KIT, NO STANDARD	\$	3,776.11
2943211	FLASK, TS 40/35, 1L	\$	784.37
2943212	FLASK, RND FLAT BOT, TS 24/25, 125ML	\$	181.93
2943213	FLASK, FLAT SID, TS 24/25, 125ML	\$	238.97

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

2943214	ELUTER, GLASS, W/#2HB	\$	973.24
2943215	CONDENSER, CROSSOVER W/#2HB	\$	1,505.14
2943216	HOLDER, STARTER, 47MM	\$	973.24
2943217	FUNNEL, 300ML, 47MM	\$	466.38
2943218	COUPLER, PTFE, 47MM	\$	156.10
2943219	SUPPORT, STAINLESS STEEL, 47MM	\$	165.35
2943220	CLAMP, ALUMINUM, 47MM	\$	466.38
2943221	TUBING, TYGON(R), 3/8" ID X 5/8" OD	\$	181.93
2943222	TUBE, COLLECTION, 30X200MM	\$	129.12
2943223	O-RINGS, COUPLER, 47MM	\$	61.86
2943224	CLAMP, PXME ACETAL, TS 24/25	\$	31.80
2943225	STOPCOCK, GLASS, 3 WAY, 4MM, ID	\$	616.70
2943226	PAN, ALUMINUM, 105MM, 100/PK	\$	208.14
2943227	PREFILTERS, 3XGMF, 42MM, 48/PK	\$	672.59
2943231	KIT STARTER W/O SR, EPA 1664A	\$	3,603.86
2943231NS	KIT STARTER W/O SR, EPA 1664A NO STANDARD	\$	2,563.18
2943232	SPE STARTER KIT, MANIFOLD, 3 PLACE	\$	5,164.90
2943233	FLASK, FLAT, TS 24/25, 12/PK	\$	2,601.72
2943235	HOLDER, MANIFOLD, 47MM	\$	973.24
2943237	PANS, ALUMINUM, 105MM, 24/PK	\$	61.86
2943238	STOPPER, SILICONE W/9.5MM HOLE, #8	\$	78.05
2943239	CARBOY, VACCUM, PP, 10L	\$	1,237.26
2943241	STOPPER, #13.5, 2HOLE/ADAPTER	\$	129.12
2943300	SPE FILTERS, 1664A, 47MM, 48/PK	\$	988.65
2943301	SPE FILTERS, 1664A, 47MM, 24/PK	\$	526.13
2943400	PIPETS, TRANSFER, PE, 15ML, 48/PK	\$	80.36
2943401	ff PIPETS, TRANSFER, PE, 15ML, 24/PK	\$	46.25
2943500	INLINE COLUMN W/SODIUM SULFATE 48/PK	\$	679.34
2943501	INLINE COLUMN W/SODIUM SULFATE 24/PK	\$	359.54
2943701	EPA 1664A STD 1X30, 4MG/ML	\$	57.98
2943702	EPA 1664A STD, 1X10ML, 4MG/ML, EA	\$	24.02
2943800	SPE CONSUMABLES KIT	\$	847.97
2943900	PLASTIC BOD BOTTLE STOPPERS, 25/PK	\$	373.88
2943901	PLASTIC BOD BOTTLE STOPPER	\$	20.43
2945029	PLAN OBJECTIVE 40X, MICROSCOPE LENS	\$	626.34
2945030	PLAN OBJECTIVE 100X, MICROSCOPE LENS	\$	869.17
2945035	LEICA, DM500, ACHROMATIC OUTFIT	\$	7,747.34
2945040	U.S. DETACHED POWER CORD ONLY	\$	69.38
2946001	FIXED VOLUME PIPETTE, 5UL	\$	93.66
2946002	FIXED VOLUME PIPETTE, 10UL	\$	93.66
2946003	FIXED VOLUME PIPETTE, 20UL	\$	93.66
2946005	FIXED VOLUME PIPETTE, 50UL	\$	93.66
2946006	FIXED VOLUME PIPETTE, 100UL	\$	93.66
2946007	FIXED VOLUME PIPETTE, 200UL	\$	103.30
2946008	FIXED VOLUME PIPETTE, 250UL	\$	103.30
2946009	FIXED VOLUME PIPETTE, 500UL	\$	103.30
2946010	FIXED VOLUME PIPETTE, 1000UL	\$	103.30

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

2946103	ff STOPPER, RUBBER 1HOLE 7/8", no 8	\$	11.95
2946401	MOISTURE ANALYZER, 90G 100-115VAC	\$	14,338.37
2946402	MOISTURE ANALYZER, 120G 100-115VAC	\$	17,421.89
2946701	ff ASSEMBLY, DISTILLATION KIT	\$	1,946.47
2946704	ASSEMBLY, SPE, STARTER, 47MM	\$	2,351.18
2946706	ASSEMBLY, FILTRATION, 47MM	\$	1,491.65
2946707	ff ASSEMBLY, MANIFOLD, 3 PLACE, SS	\$	5,164.90
2946900	2GB USB FLASH DRIVE	\$	78.05
2947000	PUMP, VACUUM 3.5CFM 115V 60HZ	\$	6,147.77
2947050	KTO: PLAN ACHROMATIC MICROSCOPE	\$	6,186.31
2947517	SPADNS2 (ARSENIC-FREE) FLUORIDE REAGENT, 3.78 L	\$	319.79
2947549	SPADNS2 (ARSENIC-FREE) FLUORIDE REAGENT, 500 ML	\$	91.23
2947549-LM	SPADNS2 (ARSENIC-FREE) FLUORIDE REAGENT, 500 ML	\$	91.23
2947553	SPADNS2 (ARSENIC-FREE) FLUORIDE REAGENT, 1 L	\$	130.98
2947600	KTO: PH BUFFER SOLUTION KIT 500mL	\$	96.84
2947801	SPE FILTERS, 47MM, 1-LAYER, 50/PK	\$	616.70
2947802	SPE FILTERS, 90MM, 1-LAYER, 50/PK	\$	1,636.19
2947808	SPE FILTERS, 47MM, 1-LAYER, 300/PK	\$	3,353.33
2947809	SPE FILTERS, 90MM, 1-LAYER, 300/PK	\$	7,843.70
2947820	FILTERS, SPE, 47MM, 50/PK	\$	323.77
2947822	FILTERS, SPE, 90MM, 50/PK	\$	971.31
2947901	ff PREFILTERS, 5UM PP, 47MM, 25/PK	\$	108.69
2947902	PREFILTERS, 5UM PP, 90MM, 25/PK	\$	217.77
2948000	PANS, ALUMINUM, 70MM, 100/PK	\$	116.60
2948600	SAMPLE PANS, ALUMINUM, 50/BX	\$	133.94
2948700	PIPET TIPS, 200 UL, 1000/PK	\$	155.91
2948710	PIPET TIP, FOR 2946006 PK/10	\$	15.61
295=100=203	FUSIBLE TEMPORISE 5 X 20 2 A	\$	41.63
2950830	aa COD SOLUTION 1, 300 TESTS	\$	91.06
2950915	aa COD SOLN 2, HR 0-1500, 150 TESTS	\$	148.15
2951015	aa COD SOLN 2, LR 0-150, 150 TESTS	\$	148.15
2952400	KTO: BOD TRAK II	\$	12,912.24
2952400-LM	KTO: BOD TRAK II	\$	12,912.24
2952401	nn vv BODTRAK II INSTRUMENT	\$	11,274.12
2952440	BODTRAK II, REAGENT SET W/BOTTLES	\$	328.83
2952440-LM	BODTRAK II, REAGENT SET W/BOTTLES	\$	328.83
2952450	aa BODTRAK II, INSTRUMENT, ASSEMBLY	\$	12,083.54
2952500	POWER SUPPLY EXT, 100-240VAC IN 24V OUT	\$	953.96
2955300	KTO: AMTAX MEMBRANE CAPS & ELECT (LR)	\$	794.97
2956200	CKT BD ASSY, PROCESSOR BOARD , BOD	\$	3,179.88
2956400	SENSOR COVER FOR BODTRAK 2	\$	292.93
2956600	LCD PCBA, 480X272, 4.3" TFT BODT 2	\$	2,062.10
2956700	SCREW, M3*14 HELIX, NYLOK CTD BODT 2	\$	1.54
2957000	FACEPLATE FOR BODTRAK II	\$	111.39
2957100	TUBING FOR BODTRAK 2	\$	110.43
2957200	CAP FOR BODTRAK 2	\$	16.57
2957300	STEM FOR BODTRAK 2	\$	7.71



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

2957400	NUMBERED PIPE SLEEVE, BODTRAK 2	\$	16.00
2957500	WASHER, PTFE, BODTRAK 2	\$	4.63
2957600	BOTTLE FENCE FOR BODTRAK II	\$	17.54
2957800	BATTERY COVER FOR BODTRAK II	\$	137.99
2957900	METAL HOUSING, BOD	\$	1,237.26
2958000	FOOT FOR BODTRAK II	\$	9.25
2958100	SCREW, PT5X10, FOR BODTRAK II	\$	3.28
2958200	2958200, SCREW, M4X20, 3 IN 1 TYPE, WITH NUT	\$	2.70
2958300	FAN FOR BODTRAK II	\$	152.44
2958400	LONG GROMMET, BODTRAK II	\$	6.17
2958500	ROUND GROMMET FOR BODTRAK II	\$	4.05
2958600	POWER SWITCH, BOD	\$	8.48
2958700	STIRRER FOR BODTRAK II	\$	1,499.36
2958800	BLOCKER, M3X2	\$	1.93
2958900	M3X8 SCREW, BODTRAK II	\$	1.73
2959000	RUBBER RING, BODTRAK II	\$	9.64
2959300	RS232 Cable, BODTrak II	\$	321.84
2959400	STIR BAR	\$	50.69
2959500	SEAL CUP, BODTRAK II	\$	11.76
2959700	BOX, BODTRAK II POWER ADAPTOR	\$	33.34
2959800	EPE CUSHION	\$	77.86
2959900	PLASTIC BAG FOR BODTRAK II	\$	5.40
2961701	STABLCAL VERIFICATION AMPULE, 10 NTU	\$	73.52
2961701-LM	STABLCAL VERIFICATION AMPULE, 10 NTU	\$	73.52
2961742	STABLCAL VERIF STD, 10 NTU, 100 ML	\$	167.67
2961742-LM	STABLCAL VERIF STD, 10 NTU, 100 ML	\$	167.67
2961749	STABLCAL VERIF STD, 10 NTU, 500 ML	\$	294.50
2961749-LM	STABLCAL VERIF STD, 10 NTU, 500 ML	\$	294.50
2961801	nn STABLCAL VERF AMPULE, 10 NTU	\$	120.86
2962266	RESPIROMETRIC BOD BUFFER SOLUTION PLW PK/50	\$	81.83
2962266-LM	ff RESPIROMETRIC BOD BUFFER SOLUTION PLW PK/50	\$	81.83
2964449	PROBE RINSE SOLN, NON-IONIC, 500ML	\$	88.70
2964926	FREECHLOR F REAGENT SOLUTION	\$	48.42
2965026	KTO: ELECTRODE FILLING SOLN, 2.44M KCl	\$	62.06
2965050	ELECTRODE FILLING SOLN, 2.44M KCl	\$	44.45
2965126	KTO: ELECTRODE FILLING SOLN, 0.02M NH4Cl	\$	77.23
2965249	ELECTRODE CLEANING SOLUTION	\$	121.77
2965249-LM	ELECTRODE CLEANING SOLUTION	\$	121.77
2965349	REDUCING SOLN, ORP ELECTRODES	\$	97.20
2965349-LM	REDUCING SOLN, ORP ELECTRODES	\$	97.20
2965503	DISPENSING TIP AND CAP, ASSY	\$	41.43
2965700	DISPENSING CAP KIT	\$	27.94
2966000	REAR PLATE, BOD	\$	75.93
2966100	BOTTOM PLATE, BOD	\$	105.23
2966500	THIMBLE, CELLULOSE, 25X80MM, 25/PK, 100% COTTON		
2966500	CELLULOSE MATERIAL	\$	443.26



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

	THIMBLE, CELLULOSE, 22X100MM, 25/PK, 100% COTTON		
2966600	CELLULOSE MATERIAL	\$	443.26
2967000	WTR CRTG FLTR - ORG REMOVAL	\$	360.39
2967100	WTR CRTG FLTR - UNIVERSAL USE	\$	466.38
2967710	aa FAST DIGEST COD SOLUTION 2 LR, 100 TESTS	\$	146.35
2968800	DO METER, PEN STYLE	\$	1,796.15
2968810	REPLACEMENT PROBE HEAD, DO PEN	\$	61.86
2968820	REPLACEMENT ELECTROLYTE, DO METER	\$	82.87
2971200	STABLCAL 500ML CALIBRATION KIT, 2100Q	\$	1,017.20
2971200-LM	STABLCAL 500ML CALIBRATION KIT, 2100Q	\$	1,017.20
2971205	STABLCAL AMPULE CALIBRATION KIT, 2100Q	\$	681.14
2971205-LM	STABLCAL AMPULE CALIBRATION KIT, 2100Q	\$	681.14
2971210	STABLCAL 100 ML CALIBRATION KIT, 2100Q	\$	587.19
2971210-LM	STABLCAL 100 ML CALIBRATION KIT, 2100Q	\$	587.19
2971304	BATTERY, NiMH AA, PK/4	\$	126.81
2971401	aa nn TURBIDITY STD, GELEX 1 NTU	\$	220.42
2971500	2100Q CASE ASSY	\$	314.13
2971826	SODIUM CHLORIDE, 100uS/cm, 50ML	\$	35.22
2971849	SODIUM CHLORIDE, 100uS/cm, 500ML	\$	146.88
2972226	SODIUM CHLORIDE, 10,000 S/cm, 50ML	\$	35.22
2972249	SODIUM CHLORIDE, 10,000 S/cm, 500ML	\$	146.88
2974226	0.001M KCI CONDUCTIVITY SOLN, 50ML	\$	35.22
2974249	0.001M KCI CONDUCTIVITY SOLN, 500ML	\$	141.10
2974326	0.01M KCI CONDUCTIVITY SOLN, 50ML	\$	35.22
2974349	0.01M KCI CONDUCTIVITY SOLN, 500ML	\$	141.10
2974426	0.1M KCI CONDUCTIVITY SOLN, 50ML	\$	35.22
2974449	0.1M KCI CONDUCTIVITY SOLN, 500ML	\$	141.10
2974549	442-30 TDS STANDARD SOLN, 500ML	\$	141.10
2974649	442-300 TDS STANDARD SOLN, 500ML	\$	141.10
2974749	442-1000 TDS STANDARD SOLN, 500ML	\$	141.10
2974826	442-3000 TDS STANDARD SOLN, 50ML	\$	35.22
2974849	442-3000 TDS STANDARD SOLN, 500ML	\$	141.10
2975149	ACID ELECTRODE CLEANING SOLN, 500ML	\$	114.55
2975149-LM	ACID ELECTRODE CLEANING SOLN, 500ML	\$	114.55
2980500	SPEC COLOR STD, DPD CHLORINE-MR	\$	655.85
2980699	AMMONIUM IONIC STR ADJ PP, PK/100	\$	196.94
2980699-LM	AMMONIUM IONIC STR ADJ PP, PK/100	\$	196.94
2982200-LM	ff KTO: CLF10sc, sc200 DUAL INPUT, COMBO pH	\$	19,349.09
2983400	ff KTO: SOLITAX sc PVC WIPER W/sc200	\$	19,773.07
2984799	NITRATE IONIC STRG ADJ PP, PK/100	\$	193.32
2984900	BACKPACK, SMALL, NYLON TEXTILE	\$	580.09
2985000	BACKPACK, LG, WHEELS, NYLON TEXTILE	\$	753.54
2985100	BRUSH, BOD BOTTLE	\$	19.08
2985900	aa KTO: HALOGEN LAMP W/FILTER PAD	\$	869.17
2986000	aa KTO: DEUTERIUM LAMP W/FILTER PAD	\$	5,627.42
2990700	KTO: SM BACKPACK, WITH CASES	\$	643.68
2990800	KTO: LG BACKPACK, WHEELS, WITH CASES	\$	803.64



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

2990900	KTO: DW STARTER KIT, SM BACKPACK	\$	13,124.23
2991300	CASE, BACKPACK, WATERPROOF	\$	157.07
2991400	CASE, COMPARTMENT, SM BACKPACK	\$	43.55
2991500	CASE, COMPARTMENT, LG BACKPACK	\$	53.58
2997100	REAGENT SET, CEL HYDRAULIC FRACTURING LAB	\$	3,017.27
2997200	APPARATUS SET, CEL/890 HYDRAULIC FRACTURING LAB	\$	3,122.06
2999200	THERMOMETER, 30 TO 50 C 250 MM	\$	314.13
2999831	1735 TRO IND SOLN 210ML	\$	77.50
2999931	1735 TRO BUFFER SOLN 237ML	\$	77.50
2A1014	ABRASIVE CLEANING PAPER	\$	76.32
2B1078	BEZEL, SPARE FOR 2B1076	\$	162.46
2D0960	MOTOR, 115VAC/50HZ/2RPM	\$	618.63
2D0961	MOTOR, 230VAC/60HZ/2RPM	\$	819.06
2D1068	TAKE UP ROLLER	\$	491.44
2D1069	ROLLER,SUPPLY,FOR 2"	\$	258.24
2D1088	GEAR TRAIN SPRING FOR	\$	42.98
2E0551	ELECTRODE, LIQUID LEVEL,	\$	316.06
2H1229	FASTENER,SCREW,4-40,SLF-C	\$	11.37
2K0086	CONNECTOR,CRIMP TERMINAL,	\$	6.36
2K0087	CONNECTOR,SOCKET TERMINAL	\$	6.36
2K0731	CONN, PLUG, 7 PIN	\$	443.26
2K0732	CONN, SOCKET, 7 PIN	\$	341.11
2K1120	CONTACT SOCKET	\$	11.37
2K1271	QUICK DISCONNECT,AMPHENOL	\$	153.02
2M1031-025-025	25' POLYETHY.TUBE,1/4"OD,	\$	304.50
2M1031-025-050	50' POLYETHY.TUBE,1/4"OD,	\$	304.50
2M1033-025-043	TUBING,1/4"ID X 7/16" OD	\$	40.66
2M1033-043	TUBING,PVC,1/4 ID.,7/16 O	\$	24.48
2M1033-043-020	PVC TUBING,20Ft.,1/4"ID,	\$	707.28
2M1033-043-025	PVC TUBING,25Ft.,1/4"ID,	\$	763.17
2M1034-025	TUBING,TEFLON,1/4"OD,1/8"	\$	21.39
2M1313	GASKET, ENCLOSURE, 8220	\$	83.83
2P2204-107	TEE,316 S.S.,150# THD. 1+1/2	\$	599.36
2P2204-108	TEE, SS.,316,THD, 2	\$	1,067.67
2P2215-105	BALL VLV,SS.,316,THD, 1	\$	1,285.44
2P2225-108	3"LG NPL,316SS,THD,2 DIA,	\$	196.57
2T1001	RELAY,REED SPST 5VDC .5A	\$	229.34
30001	ROVER, RUST REMOVER 454G	\$	57.45
30012	ROVER, RUST REMOVER 454G PK/12	\$	603.45
30040	ROVER, RUST REMOVER 18KG PAIL	\$	953.96
3004F0010-001	BRACKET,U53 SNSR SUPPORT	\$	204.28
3020800	FILTER, INTERFERENCE 880 NM	\$	1,077.30
3029100	CORD, POWER 12/3 SJT 15A-300V 6'	\$	148.97
3030700	FUSE, 1.6 A 250V 5MM X 20MM	\$	34.40
3030900	EOL FILTER, IR, LAB TURB 2100AN	\$	168.34
3031200	FILTER ASSY, EPA,LAB TURB 2100AN	\$	1,031.05
30314	FERROVER POWDER, 113G	\$	364.96

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

30353	FERROVER IRON RGT PWD, 1KG	\$	1,707.38
3039800	ASSY, FILTER HOLDER, EMPTY	\$	233.19
3046900	TOOL, COLORIMETER ALIGNMENT	\$	120.84
3047900	FLAP-DRAIN BLOCK, CL17	\$	45.67
3048400	TOOL, PLCC EXTRACTION	\$	501.07
3050001	SOFTENER, WATER PREMIUM	\$	928.91
3050051	SOFTENER, WATER PREMIUM EMPTY	\$	736.19
3050400	SCREEN, WATER SOFTENER	\$	53.58
3051100	STREET ELBOW, 1/8 NPT NYLON BLACK	\$	42.98
3060288	MANUAL, WASTEAWAY SAMPLE COND KIT	\$	124.30
3061100	TEE ASSY , W/CHK VAL, AIR INTRO	\$	88.65
3061600	TBG POLY .250 OD .040 W BLACK	\$	12.53
3061600A	TBG POLY .250 OD .040 W BLACK, NOT FDA CONTROLLED	\$	7.52
3061900	FILTER ELEMENT, 450 MESH, BROWN	\$	300.64
3062000	STRAINER, SEDIMENT 450 MESH	\$	377.73
3062900	VALVE, SOLENOID, 2WAY, AIR, SS	\$	516.49
3063900	PUMP, SUBMERSIBLE 1/6 HP MOTOR	\$	1,985.02
3064000	PUMP, AIR 0.76 CFM 115VAC 60 HZ	\$	3,931.49
3064100	PUMP, BELLOWS METERING 115V 60HZ	\$	10,773.05
3064200	PUMP BELLOWS ASSY, REPLACEMENT	\$	1,073.45
30749	NITROGEN, NITRATE 10MG/L 500ML	\$	65.57
30756	NITROGEN STD SOLN, 10MG/L 1GAL	\$	148.32
3100200	FTG, BARB 1/16-1/16 TUBE POLYPRO	\$	25.44
3100300	LATCH, PAWL CHROME PLATED	\$	58.78
3110400	VALVE, 3 WAY DISTRIBUTION	\$	1,857.82
3114600	FTG, BARB .062-.125 ID TUBE POLP	\$	25.82
3119800	TRANSMITTAL, CONDUCTIVITY INST.	\$	82.87
31249H	POTASSIUM BIODATE, 0.025N 500ML	\$	216.81
3125600	FILTER, INTERFERENCE 810NM	\$	1,106.21
3142100	IC COMPARATOR W/1.3 V REF 7665	\$	94.63
31425	POTASSIUM HYDROXIDE PELLETS, 25G	\$	53.29
31425-LM	POTASSIUM HYDROXIDE PELLETS, 25G	\$	53.29
31434	POTASSIUM HYDROXIDE, ACS 500G	\$	144.90
3146500	LENS, ASPHERIC CONDENSER 5X18 MM	\$	368.10
3148200	O-RING BUTYL .070W X .801ID	\$	20.24
3149200	IC COUPLER OPTICAL H11L2	\$	28.14
3149300	IC COUPLER OPTICAL 4N38	\$	20.24
3151600	O-RING, PORT	\$	32.96
315314	O-RING NIT .103W X .424ID	\$	12.33
31534	POTASSIUM ACID PHTHALATE, 500G	\$	296.31
31534-LM	POTASSIUM ACID PHTHALATE, 500G	\$	296.31
3153700	FASTENER, COMPRESSION R .187 ID	\$	4.43
3154100	TBG TYGON 5/32" OD X 1/32" W	\$	18.50
3155100	NIPPLE, 3/4X3/4NPT POLYETHYLENE	\$	39.70
31653	SILVER NITRATE STD .0141N 1000ML	\$	101.35
31653-LM	SILVER NITRATE STD .0141N 1000ML	\$	101.35
3176900	IC EPROM KIT T.P.A.	\$	651.39

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

3178800	PEN, DISPOSABLE CARTRIDGE, BLUE	\$	612.85
3180300	STRAW, PLASTIC ABLEX 105 RED	\$	16.96
3182700	MANUAL CHESSELL MOD NO 300 RCDR	\$	302.57
32215	BRILLIANT GREEN TUBES PK/15	\$	113.13
32332	SODIUM THIOSULFATE SOLN, 0.1N	\$	49.31
32353	SODIUM THIOSULFATE SOLN, 1000ML	\$	81.83
32600	CLAMP HOLDER	\$	96.75
32800	CLAMP, BURET DOUBLE	\$	213.92
32900	SUPPORT, BASE & ROD	\$	314.13
3305100	si SAMPLE CELL, PLASTIC SERIES 5000	\$	439.40
33109	Smart Chip for X1100 (include PCB and socket)	\$	1,161.14
33414	SILVER SULFATE 113G	\$	1,167.16
33414-LM	SILVER SULFATE 113G	\$	1,167.16
3422A1A	Conductivity Sensor,K0.05	\$	1,965.74
3422A1A34N	ASSY-SNSR,COND,K0.05,1/2"	\$	3,141.34
3422A1A50	ASSY-SNSR,COND,K.05,1/2"	\$	2,466.82
3422A1B	ASSY-SNSR,COND.,0.05K,1/2	\$	2,158.46
3422A1C	ASSY-SNSR,COND,0.05K,1/2"	\$	2,235.55
3422A2A	ASSY-SNSR,COND,K0.05, 1/2	\$	2,351.18
3422A2A50N	ASSY-SNSR,COND,K0.05, 1/2	\$	2,890.80
3422A2AX	CONDUCTIVITY SENSOR,K0.05, 1/2	\$	53.58
3422A2B	ASSY-SNSR,COND,0.05,1/2"S	\$	2,543.90
3422A2D	ASSY-SNSR,COND.,0.05K,1/2	\$	2,890.80
3422A3A	ASSY-SNSR,COND.K0,05,3/4"	\$	1,965.74
3422A3B	ASSY-SNSR,COND,0.05K,3/4"	\$	2,158.46
3422A3C	ASSY-SNSR,COND.0.05K,3/4"	\$	2,235.55
3422A4A	ASSY-SNSR,COND,K0.05,3/4"	\$	2,466.82
3422A4C	ASSY-SNSR,COND.,0.05K,3/4	\$	2,736.62
3422A5A	SNSR-COND,.05,SS-T,20'	\$	2,466.82
3422A5A30P	SNSR-COND,.05,SS-T,30'CBL	\$	3,815.86
3422A5C	ANSR-COND,.05,SS-T,AL HD	\$	2,235.55
3422B3A	Conductivity Sensor,K0.5	\$	1,965.74
3422B3A20G	Generic 0.5K Sensor with 20' Cable	\$	1,965.74
3422B3B	ASSY-SNSR,COND,K0.5,3/4"	\$	2,158.46
3422B3C	ASSY-SNSR,COND.,K0.5,3/4"	\$	2,235.55
3422B4A	ASSY-SNSR,COND,K0.5, 3/4"	\$	2,466.82
3422B4B	ASSY-SNSR,COND.,K0.5,3/4"	\$	2,659.54
3422B4C	ASSY-SNSR,COND,K0.5,3/4"	\$	2,736.62
3422B4D	ASSY-SNSR,COND,K0.5,3/4"	\$	3,006.43
3422B6A	SNSR-COND, .5,SS-T,20'	\$	2,466.82
3422B6B	SNSR-COND,.5,SS-T,PP HD	\$	2,659.54
3422C3A	Conductivity Sensor,K1.0	\$	1,965.74
3422C3A50	ASSY-SNSR,COND., K1.0, 3/	\$	2,466.82
3422C3AX	Conductivity Sensor,K1.0	\$	53.58
3422C3B	ASSY-SNSR,COND., K1.0, 3/	\$	2,158.46
3422C3C	ASSY-SNSR, COND., K1.0, 3	\$	2,235.55
3422C4A	ASSY-SNSR,COND.,K1.0, 3/4	\$	2,466.82

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

3422C4A20G	ASSY-SNSR,COND.,K1.0, 3/4	\$	2,466.82
3422C6A	SNSR-COND,1.0,SS-T 20'	\$	2,466.82
3422D3A	Conductivity Sensor,K5.0	\$	1,965.74
3422D3A50N	ASSY-SNSR,COND.,K5.0,50'	\$	2,466.82
3422D3B	ASSY-SNSR,COND., K5.0, 3/	\$	2,158.46
3422D3C	ASSY-SNSR, COND., K5.0, 3	\$	2,235.55
3422D4A	ASSY-SNSR,COND.,K5.0,	\$	2,466.82
3422D6A	SNSR-COND,5.0,SS-T,20'	\$	2,466.82
3422E3A	ASSY-SNSR,COND.,K10.0	\$	1,965.74
3422E3A20G	ASSY-SNSR,COND.,K10.0,GEN	\$	1,965.74
3422E3B	ASSY-SNSR,COND.,K10.0,	\$	2,158.46
3422E4A	ASSY-SNSR,COND.,K10.0,	\$	2,466.82
3422E4A20E	ASSY-SNSR,COND.,K10.0,	\$	2,563.18
3422E6A	ASSY-SNSR,COND.,K10.0,	\$	2,466.82
34250000	ULR FLUORESCENCE SULFITE REAGENT KIT	\$	312.57
34250001	ULR FLUORESCENCE SULFITE REAGENT REFILL KIT	\$	299.92
34251000	ULR FLUORESCENCE FREE CHLORINE REAGENT KIT	\$	507.70
34251001	ULR FLUORESCENCE FREE CHLORINE REAGENT REFILL KIT	\$	496.86
34252000	ULR FLUORESCENCE TOTAL CHLORINE REAGENT KIT	\$	552.87
34252001	ULR FLUORESCENCE TOTAL CHLORINE REAGENT REFILL KIT	\$	540.22
34332	POTASSIUM IODIDE SOLN, 100ML MDB	\$	52.20
34332-LM	POTASSIUM IODIDE SOLN, 100ML MDB	\$	52.20
3433B8A	ASSY-SNSR,COND,K0.5, 3/4"	\$	1,618.85
3433B8A30N	SNSR,COND,K0.5, 3/4" NPT	\$	2,081.38
3433B8A50N	SNSR,COND,K0.5, 3/4" NPT	\$	2,139.19
3433B8C	ASY-SNSR,COND,K0.5,AL	\$	1,871.31
3433E8A	ASY-SNSR,COND.,K10.0	\$	1,618.85
3433E8A30N	ASSY-SNSR,COND,K10.0,30'C	\$	2,736.62
3433E8B	ASY-SNSR,COND.,K10.0,PP	\$	2,312.64
3433E8D	ASY-SNSR,COND.,K10.0,SS	\$	2,139.19
34349	POTASSIUM IODIDE SLN,30G/L 500ML	\$	90.70
3444B8A	ASSY-SNSR,CON.,K0.5	\$	3,141.34
3444B8A50N	0.5 Cell Conductivity Sensor, 50 Ft.	\$	3,623.14
3444B8C	ASSY-SNSR,COND,K0.5,AL HD	\$	3,411.14
3444B8D	ASSY-SNSR,COND,K0.5,SS HD	\$	3,623.14
3444D8A	ASSY-SNSR,CON.,K5.0	\$	3,141.34
3444D8B	ASSY-,SNSR,COND.,K5.0,PP	\$	3,353.33
3444D8C	ASSY-SNSR,COND,K5.0,AL,HD	\$	3,411.14
3444D8D	ASSY-SNSR,COND,K5.0,SS HD	\$	3,623.14
3455A6A	ASSY-SNSR,OVERALL,0.05K,	\$	2,890.80
3455A6B	ASSY-SNSR,OVERALL,0.05K,	\$	3,025.70
3455A6D	ASSY-SNSR,OVERALL,0.05K,	\$	3,391.87
3455C6A	ASSY-SNSR,OVERALL,1K,	\$	2,890.80
3455C7A	ASSY-SNSR,OVERALL,1K,	\$	2,967.89
3455E6A	ASSY-SNSR OVRL,K10.0,	\$	2,890.80
3455E7A	ASSY-SNSR,OVERALL,K10.0,	\$	2,967.89
3455E7B	ASSY-SNSR,OVERALL,K10.0,	\$	3,160.61

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

34630000	DR1300 FL FLUORESCENCE STANDARDS KIT	\$	345.09
34932	STARCH INDICATOR SOLN, 100ML MDB	\$	40.83
34953	STARCH INDICATOR SOLN, 1000ML	\$	62.32
35126	SODIUM THIOSULFATE 50ML SCDB	\$	45.89
35253	SODIUM THIOSULFATE, .025N 1000ML	\$	68.10
356099,05090	O-RING FOR VENTURI VALVE, SODI	\$	38.35
356099,05220	O-ring Viton 22x1 mm (around mirror)	\$	70.54
356099,15040	O RING, VITON, FOR TEMPERATURE	\$	31.41
3563500	HOLDER, VIAL, 16mm ROUND VIAL, 4ml FILL	\$	25.63
35732	COPPER SULFATE-SULFAMIC ACID IN-	\$	42.09
35732-LM	COPPER SULFATE-SULFAMIC ACID IN-	\$	42.09
35749	COPPER SULFATE-SULFAMIC ACID	\$	66.13
359000,40780	PT100, COMPLETE	\$	356.53
359016,10000	ADAPTER, AK9 TO S7/S8	\$	118.33
359016,10040	AS7 type cable connector	\$	206.21
359016,10050	AS9 type cable connector	\$	206.21
359016,10101	oo Electrode cable, 0.5m with AS7 connector (for reference electrode only) - Use on 9073, 9186	\$	674.52
359016,10102	Electrode cable, AS7 connector, 8810 (V1, V2) reference electrode	\$	310.28
359016,10103	Electrode cable, AS9 connector for 8810 (V1, V2)	\$	402.78
359016,10110	ELECTRODE CABLE, 3 M W/ PLUG A	\$	522.27
359016,10120	ELECTRODE CABLE, 10 M W/ PLUG	\$	792.08
359016,10122	CABLE FOR PH W/AS9 PLG, 20M (8362)	\$	840.26
359016,10140	pH electrode cable, 3 m with AS7 connector (for reference electrode)	\$	375.80
359016,10141	pH electrode cable, 10 m with AS7 connector (for reference electrode)	\$	381.59
359016,10142	pH electrode cable, 20 m with AS7 connector (for reference electrode)	\$	321.84
359026,00000	ELECTROLYTE RESERVOIR W/ TUBIN	\$	1,221.84
359035,00142	Drive shaft Inox 1500 mm	\$	443.26
359035,00700	Electrode container	\$	1,337.48
359072,00120	NUT, BLACK, DO PROBE/BODY	\$	25.63
359090,00024	NOZZLE, VENTURI	\$	1,588.01
359090,20000	FILTER CARTRIDGE 12 PC SET	\$	1,416.49
359090,70015	TUBE SET TYGON (12 pcs) GREY	\$	1,331.70
359090,70020	TUBING SET, TYGON, 12PCS (8810)	\$	892.29
359090,70021	TUBE SET TYGON 12 PC PURPLE-	\$	714.99
359090,80000	TUBE CASSETTE FOR PERI PUMP	\$	1,416.49
359091,00050	Stainless steel counter electrode	\$	686.08
359092,00080	SEAL, FOR MEASURING CELL, UPPE	\$	47.99
359092,00201	oo FARADAY CELL	\$	1,349.04
359103,10077	Elbow connector (90 deg) to link 3/8 inch gas thread to 10x12mm tubing in PP-used on 8810	\$	343.04
359103,10094	Male/male fittings G3/8	\$	67.84
359103,10132	1/8" NPT nipple for 6/8 mm tube (8544)	\$	44.90



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

359110,01100	Mounting plate for analyser 8810 made of PP	\$	3,758.04
359110,01300	VESSEL	\$	10,464.70
359110,01400	SIPHON	\$	909.64
359110,01500	TUBE, OVERFLOW	\$	269.81
359110,01650	STIRRER COMPLETE	\$	3,102.79
359110,01810	SPRINKLER COMPLETE	\$	2,235.55
359110,02000	VALVE, 2/2 SOLENOID 610 NW1	\$	5,203.44
359110,02100	VALVE, 3/2-WAY SOLENOID M311 N	\$	1,169.81
359110,02200	VALVE, 2/2 SOLENOID (DRAIN)	\$	4,837.27
359110,02300	Screw in socket G3/8PP	\$	240.90
359110,02400	DETECTOR, SAMPLE LEVEL COMPLET	\$	2,331.91
359110,05000	Complete tubing level detector for 8810	\$	128.54
359110,40000	Power supply board for 8810 + 8811 - 115V	\$	5,916.50
359110,40001	Power supply board for 8810 + 8811 - 220V	\$	8,229.14
359110,40002	Power supply board for 8810 + 8811 - 240V	\$	8,036.42
359110,52012	TBNG, SILICN 12X17MM TO DRAIN (8810)	\$	194.65
35949	FLUORIDE STD SOLN, 10MG/L 500ML	\$	108.93
359500,11106	RACCORD CANNELE POUR TUBE DE 1,6mm	\$	70.54
3630000	FILM FORMING AMINE KIT	\$	858.21
363130,00500	BUFFER SOLN, PH4.00 500 ML	\$	68.99
363131,00500	BUFFER SOLN, PH6.88 500 ML	\$	76.70
363132,00500	Buffer solution pH 9.22, 500 ml	\$	125.03
363140,00500	CHEMICAL, KCL SOLUTION, 500ML	\$	200.55
363143,00500	ELECTROLYTE, REFERENCE KCL 0.	\$	166.22
363144,00500	ELECTROLYTE, REFERENCE KBR, 5	\$	411.94
363371,00050	Flow-through assembly NW50/ND10	\$	8,325.50
363633,10000	oo Electrolyte reservoir for wall mounting	\$	1,061.89
363700,74106	DIAPHRAGM ASSEMBLY	\$	1,025.27
363700,83401	HOOD COMPLETE: 2 NUTS, O-RING	\$	988.65
363700,83402	Set: 4 O-rings, 2 washers, 2 plugs, 2 sleeves	\$	298.72
363700,83403	Glass bowl: vessel, gasket	\$	1,131.27
363700,83441	Set: 5 O-rings, 2 washers, 2 blind plugs, 2 sleeves	\$	506.85
363700,83462	CLEANING BLOCK SET: BLOCK, HUB	\$	1,104.29
363700,83601	Set: 4 O-rings, washer	\$	321.84
363700,83602	Glass bowl: sample cup, gasket	\$	703.43
363700,88770	Set of beads, Chloromat + Ozonmat + Hydrastat + Nitramat	\$	595.50
363877,06000	FILTER, IN-LINE	\$	907.71
368332,00000	PROBE, INDUCTIVE , K=10, PT100	\$	8,248.42
368366,00000	PROBE, PH POLYPACK	\$	3,603.86
368402,00000	ELECTRODE, GLASS	\$	1,358.68
368404,00000	ELECTRODE, STANDARD GLASS	\$	1,370.24
368411,10000	Combination electrode (8411 B1)	\$	1,431.91
368417,00000	Combination-Elect.(8417B)	\$	2,081.38
368418,00000	ee ELECTRODE, COMBINATION W/ GEL	\$	1,393.37
368422,00000	ELECTRODE, REFERENCE	\$	1,017.56
368429,00000	ELECTRODE, REFERENCE FOR ELEC	\$	1,431.91
368430,00000	Platinum electrode (8430B)	\$	1,225.70



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

368444,00000	Standard reference electrode (8444B)	\$	1,268.10
368447,00000	ee ELECTRODE, REFERENCE	\$	1,541.76
368470,00000	FLUORIDE ION-SELECTIVE ELECTRO	\$	16,651.01
368471,00000	Chloride ion-SELECTive electrode with DIN plug	\$	10,233.43
368474,00000	ELECTRODE, CYANIDE DIN PLUG (8810)	\$	10,503.24
368481,10000	DIAPHRAGM TUBE	\$	657.18
368490,00000	ELECTRODE, SULFIDE DIN PLUG (8810)	\$	8,344.78
368495,00000	SENSOR, TEMPERATURE PT100, 3-M	\$	2,100.65
368495,10000	Temp. sensor Pt100, fixed cable for 8810 (V1, V2)	\$	238.97
368810,21240	oo 8810 Analyser ORP, 1 peristaltic pump, Silver elect. 240V/50Hz	\$	86,839.63
368810,31240	oo 8810 Analyser ISE, 1 peristaltic pump, Chloride elec. 240V/50Hz	\$	89,499.17
368810,52220	oo 8810 Analyser ORP, 1 piston pump, Copper elect. 220V/50Hz	\$	92,332.15
368810,56000	CHEMICAL CLEANING	\$	13,355.50
368810,56500	Chemical cleaning pump (w/o tubing, canister or level detector)	\$	8,537.50
368810,60000	AUTOMATIC SYSTEM CALIBRATION	\$	13,413.31
368810,70050	Additional reagent pump (3 max.) 24V, 50 Hz complete	\$	18,404.76
368810,72000	Automatic calibration pulse pump for 8810 ISE complete	\$	21,025.75
368810,76000	8810 AUTOMATIC HEATING SYSTEM ISE PH	\$	1,676.66
368810,91010	KIT 2 ANS-8810 pH SANS P.P.	\$	9,173.47
368810,92000	Spare part kit_2 years, 8810 ORP analyser w/ platinum electrode	\$	15,822.31
368810,92004	Spare parts kit_2 years, 8810 ORP analyser w/ sulfide electrode	\$	26,036.47
368810,92011	Spare part kit_2 years, 8810 ORP w/ silver stick elect., piston pump	\$	10,503.24
368810,92012	Spare part kit_2 years, 8810 ORP w/ copper elect., piston pumps	\$	17,653.15
368810,93001	Spare parts kit_2 years, 8810 ISE analyser w/ high range chloride electrode	\$	13,259.14
368810,93008	2 years kit,8810 analyser LL Chloride (0-5mg/Cl-)	\$	14,319.10
368810,95000	TUBING AND CONN.S KIT	\$	1,031.05
370=304=002	Shielded 4 conductor cable (per meter)	\$	123.73
370=506=025	6 COND. CABLE, SHIELDED, LI, PERM.	\$	67.64
3700E-SANX	Small-Bore, Sanitary Components	\$	62.63
3705E2T	ELECTRODELESS, SANITARY, POLY	\$	3,372.60
3705E2T100N	ELECT.COND.SNSR,SAN. FLNG	\$	5,164.90
3705E2T20P	FINISHED SENSOR,SANITARY,	\$	3,391.87
3705E2T50N	ELECT.COND.SNSR,SANTRY FL	\$	4,625.28
3706E2T	Electrodeless, Sanitary, Kynar	\$	3,796.58
3706E2T-1008	SMALL BORE SENSOR, ANDERSON, KYNAR	\$	3,931.49
3708E2T	Electrodeless, Sanitary, Teflon	\$	4,394.02
3708E2T20G	Generic Sanitary Sensor	\$	4,394.02
3708E2T50N	ELECT.COND SNSR,TEFLON,50	\$	5,685.24
3725-9020	SMALL BORE, MCGEAN ROHCO	\$	2,929.34
3725E2T	Electrodeless, Convertible, Poly	\$	2,948.62
3725E2T20G	Electrodeless, Convertible, Poly, Generic	\$	2,948.62
3725E2T20P	FINISHED SNSR,CONV.,PP	\$	3,160.61
3725E2T33N	ELECT.COND.SNSR (SML BOR)	\$	4,143.48



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

3725E2T50	SMALL BORE ELECT.COND CL,	\$	4,355.47
3726E2T	Electrodeless, Convertible, Kynar	\$	3,507.50
3726E2T30N	ELECTRODLS COND.SNSR,CONV	\$	4,586.74
3726E2T50N	GLI-SMALL BORE ELECT. SNS	\$	4,818.00
3727E2T	Electrodeless, Convertible, PEEK	\$	3,738.77
3727E2T20P	SENSOR-ELTRDLS,CONV,PEEK	\$	3,738.77
3727E2T50N	PEEK Elec.Cond.Sensor,50'	\$	4,952.90
3728E2T	Electrodeless, Convertible, Teflon	\$	3,989.30
3728E2T100N	ELCT.CND SNSR,CNVRT,100',	\$	5,858.69
3728E2T20G	Generic Convertible	\$	3,989.30
3728E2T20S	ELECT.COND.SNSR,CONV,PFA,	\$	4,066.39
3728E2T30N	ELECT.COND.SNSR,CONV,PFA,	\$	5,068.54
37333	PH PAPER, 3.0-5.5 PK/5	\$	118.69
37933	PH PAPER, 6.0-8.0 PK/5	\$	118.69
384126	CONN, .156 CTR PC BD EDGE 12 CKT	\$	33.34
38533	PH PAPER, 9.0-12.0 PK/5	\$	115.63
39133	PH PAPER, 1.0-11.0 PK/5	\$	103.88
39532	SOAP SOLUTION(50ML SMPL) 100ML MDB	\$	56.90
39553	SOAP SOLUTION(50ML SAMPLE) 1000ML	\$	135.51
39732	SILVER NITRATE SOLN, 100ML MDB	\$	83.47
39732-LM	SILVER NITRATE SOLN, 100ML MDB	\$	83.47
3A1031	EXPEND, BUFFER, PH 4,	\$	1,013.59
3A1032	EXPEND, BUFFER, PH 7,	\$	1,013.59
3A1033	EXPEND, BUFFER, PH 10, 20	\$	1,013.59
3A1066	BATT,NICAD,3.6V,110MAH,PC.	\$	312.21
3A2015	CHART,40-GRID/ROLL	\$	166.70
3A2017	CHART,60-GRID/ROLL	\$	166.70
3A2018	CHART,70-GRID/ROLL	\$	166.70
3A2019	CHART,50-GRID/0-100%/ROLL	\$	166.70
3A2020	CHART,50-GRID/24HR/ROLL	\$	166.70
3A2022	CHART,0-50 GRID/24 HOUR	\$	291.01
3A2023	CHART,0-60 GRID/24 HOUR	\$	291.01
3A2024	CHART,0-70 GRID/24 HOUR	\$	291.01
3A2025	CHART,0-80 GRID/24 HOUR	\$	291.01
3E0918	ELECTRODE, PLATINUM, COMB	\$	1,298.93
3E1006	ELECTRODE, PLATINUM, ORP	\$	1,073.45
3E1069	ELCTRD,SILVER,ORP,W/EPOXY	\$	855.68
3H1003	ENCL. UNILET, PVC, 1+1/4	\$	98.09
3H1020	CORD GRIP=S1121 CORD GRIP	\$	65.14
3H1091	CORD GRIP, .157 TO .354,	\$	70.34
3H1228	CORD GRIP,NYLON,.187-.312	\$	65.14
3H1230	LOCKNUT,1/2" CONDUIT,	\$	33.34
3H1252	THERMO-HEAD,ALUM,1/2X3/4,	\$	229.34
3H1263	Fitting, Single Barbed	\$	42.98
3K1104	BANANA PLUG, 6-32 X 3/8"L	\$	70.34
3P2101-105	PIPE, CPVC, 1" SCH 80	\$	6.17
3P2101-107	PIPE/FT, CPVC, 80,1+1/2	\$	6.17

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

3P2102-107	UNION, CPVC, 80,THD,1+1	\$	240.90
3P2104-107	TEE, CPVC, 80,THD,1+1	\$	127.77
3P2104-108	TEE,CPVC SCH.80, THD. 2"	\$	194.65
3P2120-103	COUPLING,CPVC, 80,THD,1/2	\$	129.32
3P2120-105	COUPLING,CPVC, 80,THD, 1	\$	102.91
3P2120-107	COUPLING,CPVC, 80,THD,1+1	\$	120.06
3P2120-122	Reducer Coupling, CPVC, 3/4" x 1/2"	\$	135.87
3P2120-125	RED.COUP,CPVC, 80,THD,1 X 1/2	\$	135.87
3P2120-132	RED.COUP,CPVC, 80,THD,1+1/2" x 1"	\$	237.05
3Q1075	MOSFET,ISOL TO220,2.2 OHM	\$	27.37
3T1001	RELAY,SS,SPST,115/230V 3A	\$	491.44
3T1003	RELAY,SS,3-32VDC IN, 25A	\$	362.31
3V0049	BATTERY AND HOLDER FOR	\$	65.14
3W1100	PLUG 2 POS.W/5.08mm SPACE	\$	40.28
3W1106	3 POS.PLUG W/5.088mm SPAC	\$	40.28
400004	Cable, RS-485	\$	37.77
4001LOA0ANN1400	GP LVL SNSR,316SS,TFE INS	\$	1,510.92
4003LOA0ANN2400	GP LVL SNSR,316SS,POLYETH	\$	1,566.81
4009800	TIE, CABLE	\$	7.52
400A1000-001	FLOW CHAMBER,KYNAR,ORBISP	\$	647.54
4012500	DRIVEBAR, PISTON	\$	443.26
4012800	PLATE, CAM	\$	635.98
4012900	ACTUATOR, PINCH VALVE/DRIVEBAR	\$	321.84
40149	POTASSIUM IODIDE-IODATE, 500ML	\$	94.48
4017600	PLATE, SLIDER	\$	29.87
401814	SEAL QUAD RING .9375 NOM ID	\$	16.57
4018300	SPACER, BEARING	\$	51.26
4018400	SHIM, REACTION PLATE	\$	110.24
402-400-400	Stainless steel DN40 nut for 8362	\$	180.00
4021100	O-RING FLCAR .103W X .299ID	\$	17.34
4021500	O-RING FLCAR .103W X .424ID	\$	21.58
4021800	FTG, TEE FITS .062 ID TUBE POLYP	\$	33.15
4022800	BEARING, LIGHT	\$	574.31
40264BB	GLI pH/DO Sensor Extension Cable, 20ft, 1W1100, 6 Wire - USABB Only	\$	321.84
4026700	BEARING, BALL	\$	69.76
4028200	FLOW METER, 50-700 ML/MIN WITH	\$	757.39
4028400	FLOW METER W/VALVE, 100-1750 ML/MIN	\$	1,032.98
4028800	FLOW METER, KIT, AIR 10-100 SCFH	\$	788.22
4029400	si COVER, VENT SS6	\$	1,032.98
4029900	THREADED DISK	\$	116.79
4030100	si FITTING, 3/4" MODIFIED ELBOW	\$	78.05
4030200	ORIFICE, FLOW MULTIPLIER	\$	381.59
4034900	SCREW 4-40X .125 +.00 -.03	\$	2.51
4035500	BULKHEAD ADAPTER, 1", PP EPDM	\$	323.77
4037200	FTG, BARB 1.00TBF X 1.0NPT TFPP	\$	39.89
4041500	CONN, AC RECEPTACLE UL/CSA 3 CKT	\$	31.22

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

4041800	IC REGULATOR VOLT +15V 340T	\$	19.46
4043900	FTG, BARB .75TBF X .75NPT TFPP	\$	39.89
4049600	FTG, ELBOW 3/4 ID TBF X 3/4NPT	\$	40.86
40502	FLUORIDE STD SOLN 0.2 MG/L 500ML	\$	73.16
40505	FLUORIDE STD SOLN 0.5 MG/L 500ML	\$	73.16
40508	FLUORIDE STD SOLN 0.8 MG/L 500ML	\$	73.16
40512	FLUORIDE STD SOLN 1.2 MG/L 500ML	\$	73.16
40515	FLUORIDE STD SOLN 1.5 MG/L 500ML	\$	73.16
40520	FLUORIDE STD SOLN 2.0 MG/L 500ML	\$	73.16
4052600	TRAP ASSY, SS6-HST	\$	556.96
4057000	SCREW 4-40X .312FLAT SST PHIL	\$	11.37
4058800	IC AMP,OP LOW POWER QUAD 324	\$	38.93
40653	OXALIC ACID RGT, 75G/L 1000ML	\$	158.80
4065800	UNION, BULKHEAD .25	\$	402.78
4066103	FILTER ASSY, INTERFERENCE 650 NM	\$	558.89
4075000	FTG, BARB .062-.093 ID TUBE POLP	\$	27.94
4078200	WASHER, BEARING SPACER	\$	8.48
4078800	O-RING FLCAR .070W X .864ID	\$	39.31
4079000	SCREW 2-56X .250FLAT SST PHIL	\$	8.67
4081100	si WASHER, BOD BOTTLE CAP	\$	28.14
4082500	TUBING HEAT SHRINK .063 MIN ID	\$	31.22
4088300	FILTER, INTERFERENCE 650 NM	\$	541.54
4088600	BEARING, BALL	\$	98.09
4093800	INSERT, SLIDER BLOCK	\$	254.39
4099700	O-RING FLCAR .070W X .239ID	\$	13.10
4099800	COVER, 1 CM / 1 IN PHOTOCCELL	\$	198.50
4100200	O-RING FLCAR .103W X 1.237ID	\$	43.94
4100300	CASE,MINI-LINE WHITE TING-301	\$	128.93
4101700	WRENCH, LENS RETAINER	\$	171.91
410316	MAGNET, BAR 2.00LX.500HX.250W	\$	31.80
4104400	si PLUG, LOWER CYLINDER	\$	179.81
4107000	si oo HOLDER,SAMPLE CELL 1CM	\$	21.39
4107700	HEATSINK, LOW PROFILE	\$	10.41
4108300	FASTNER,KNOB .30IDX.19HTX.014THK	\$	17.34
4108700	SCREW 2- X .312 PAN SST TPG	\$	5.59
4111300	JACK,SUB-MIN 2-COND SHUNTED	\$	16.57
4114388	INSTRUCTIONS, KIT MOTOR	\$	29.10
4117000	FLOC JAR, WAGNER	\$	458.67
4117008	INSTRUCTION SHT-FLOC JAR, WAGNER	\$	29.29
4117600	FILTER, INTERFERENCE 480 NM	\$	1,364.46
4125800	INSERT, FOAM UPPER DR 100	\$	29.10
4126000	oo HOLDER ASSY, SAMPLE CELL, 1CM	\$	51.46
4131600	O-RING, -114 NITRILE FDA APVD	\$	35.85
4134008	INSTRUCTIONS, END POINT ADAPTOR	\$	12.33
4134400	TBG TYGON 1/4" OD X 1/16" W	\$	31.99
4138600	TUBING, GLASS STD WALL	\$	68.03
4147000	CABLE, 3 COND 18 AWG STRD SJT	\$	7.71



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

4157800	DELIVERY TUBE,90DEG&180HOOK PK/5	\$	24.48
4158700	IC CONV POS TO NEG VOLT 7660	\$	29.10
4159900	LED, RED, T 1-3/4 SIZE SUPER BRT	\$	10.99
4165802	SAMPLE CELL, 1-CM/10-ML PK/2	\$	78.05
4168100	FAN, VENTILATION 115V 50/60HZ	\$	398.93
4168200	FILTER, FAN	\$	130.09
416831	THERMOMETER, -10 TO 100 DEG C	\$	283.30
4170200	FILTER, INTERFERENCE 810 NM	\$	181.73
4190900	CONN, .100 CTR STR POST 6 CKT	\$	19.46
4191300	CONN, .156 CTR 22AWG 3 CKT HSG	\$	8.48
4192900	CLOSED END W/PART GOLD CONTACTS	\$	44.90
420119	LENS HOLDER SPRING	\$	74.97
4207600	TBG TPE .125 OD .062 ID	\$	27.37
4208600	TBG LATEX 3/8" ID X 1/2" OD X	\$	28.33
4215200	TBG POLY .250OD .040W	\$	25.63
42200	CLAMP, 3 PRONG	\$	151.09
4224100	RELAY, DIP REED SPDT	\$	198.50
42426	BUFFER SOLN, HDNS 1 50ML SCDB	\$	39.56
42426-LM	BUFFER SOLN, HDNS 1 50ML SCDB	\$	39.56
42432	BUFFER SOLN, HDNS 1 100ML MDB	\$	56.01
42432-LM	BUFFER SOLN, HDNS 1 100ML MDB	\$	56.01
42449	BUFFER SOLN, HARDNESS 1 500ML	\$	91.23
42449-LM	BUFFER SOLN, HARDNESS 1 500ML	\$	91.23
42453	BUFFER SOLN, HARDNESS 1 1000ML	\$	159.72
42453-LM	BUFFER SOLN, HARDNESS 1 1000ML	\$	159.72
425=290=210	AMP SENSOR NUT FOR 9184sc	\$	92.31
42532	HARDNESS 2 TEST SOLN 100ML MDB	\$	58.71
42532-LM	HARDNESS 2 TEST SOLN 100ML MDB	\$	58.71
42549	HARDNESS 2 TEST SOLN 500ML	\$	109.30
42549-LM	HARDNESS 2 TEST SOLN 500ML	\$	109.30
42632	HARDNESS 3 TEST SOLN 100ML MDB	\$	33.05
42649	HARDNESS 3 TEST SOLN 500ML	\$	97.20
42653	HARDNESS 3 TEST SOLN 1000ML	\$	137.67
427=010=060	KIT FOR 8810 V2 CABLE ENDING	\$	217.77
427=010=070	Kit of caps & counter-nuts for 8810 (V1,V2) (30x)	\$	389.29
4270400	MOTOR, SYNCH,24 V 60 HZ 1.5 WATT	\$	373.88
4270500	MAGNET, U-SHAPE ROUND .50 DIA	\$	69.76
4271700	TBG TPE .193OD .063ID +/- .002	\$	39.31
4274100	si BLOCK, REAGENT	\$	178.07
4274200	si BLOCK, SAMPLE	\$	196.57
4274300	si PLATE, PRESSURE	\$	250.54
428=112=025	EPDM DN40 GASKET (8362)	\$	60.90
42849	CALCIUM CHLORIDE SOLN, 500ML	\$	79.50
429=042=019	O-RING, SMALL, 9180 SERIES PROBE	\$	20.43
429=120=030	O-RING 12X3MM, EPDM	\$	31.03
429=500=510	JOINT CLAMP EPDM 2""(DIAMETRE:51)	\$	46.64
429=600=500	EPDM gasket for 8398,3	\$	107.54

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

42953	FERRIC CHLORIDE SOLN, 1000ML	\$	60.35
43049	MAGNESIUM SULFATE SOLN, 500ML	\$	56.01
431=100=210	Colorimetric cell nut PG21 - Brass coated nikel	\$	58.97
431=210=006	CAP IN PE FOR 6MM HOLE	\$	3.85
431=210=020	CAP IN PE FOR 20MM HOLE	\$	3.85
43149	PHOSPHATE BUFFER SOLN, 500ML	\$	63.60
43149-LM	PHOSPHATE BUFFER SOLN, 500ML	\$	63.60
4319600	ADJUST, 0.1 ML VOLUME	\$	158.42
4321300	SYRINGE, 3CC LUER LOCK TIP	\$	9.64
4323800	VARIATOR ASSY, DILUTER/PIPETTER	\$	116.98
4329300	TBG TPE .250OD .125ID	\$	26.21
433519	LENS MOUNT	\$	166.70
4346400	SYRINGE, 1 CC MODIFIED	\$	210.06
4350200	oo FTG, TEE .062ID TBG 10-32TAPER	\$	28.14
4350800	LED, ORANGE, T-1 3/4 SIZE W/LEAD	\$	24.48
4351600	SCREW 10-24X .625RND SSTHUMB	\$	36.04
4351900	MOTOR ASSY, HARDNESS MONITOR	\$	826.77
4352900	si CAM, OPTICAL	\$	516.49
4353500	si BASE, STRAINER	\$	283.30
4361500	COUPLING, 1/4"-28 UNF-2B PVC	\$	130.09
4365500	FEET, PUSH RIVET BLACK VINYL	\$	14.07
437044	THERMISTOR 30K AT 25 DEG C	\$	187.90
4379400	BUSHING, STRAIN RELIEF	\$	44.90
43800	TUBE, MEASURING PLASTIC 5.83ML	\$	9.83
4381300	HOUSING, DISPLAY H-1	\$	45.67
4382200	IC AMP,OP LOW POWER 7611 CMOS	\$	20.43
4382500	COVER, DISPLAY H-1	\$	51.65
4383500	BUMPER, BUTTON, .375 DIA X.125TH	\$	12.53
4387501	CASE, CUSTOM FOR HACH ONE PH MTR	\$	406.64
43900	BOTTLE, SQ GLASS 29ML *SELL-06	\$	6.75
43906	BOTTLE, SQUARE GLASS 29ML PK/6	\$	36.81
4397500	TEST KIT, SAMPLE DEGASSING	\$	119.29
4397510	TEST KIT, FILTRATION & DEGASSING	\$	262.10
440426	SOCKET, LAMP MOUNTING	\$	36.04
4404400	HINGE ASSY, SSX	\$	161.88
4409200	IC SWITCH QUAD BILAT 74HC4066	\$	26.02
4410700	OPTICAL HEAD, 1720 SERIES	\$	321.84
4411400	LENS, ASPHERIC CONDENSER 30X25MM	\$	206.21
4411600	CLOSURE, SNAP-TITE	\$	61.48
4411700	si APERTURE, 1720C	\$	80.94
4412700	TRANSFORMER,POWER,DUAL PRIMARIES	\$	229.34
4412800	CONNECTOR, STRAIN RELIEF .400	\$	76.51
4415300	si CYLINDER ASSY, 1720C CALIBRATION	\$	310.28
4415600	KTO: CALIB CYLINDER 1720C/1720D/E	\$	1,460.82
4415600-LM	KTO: CALIB CYLINDER 1720C/1720D/E	\$	1,460.82
4416200	SPACER 1/2 OD X 9/32 ID X 1/2LG	\$	2.12
4418011	MANUAL,SERV. CHLORINE ANALYZER	\$	215.85

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

4418600	BATTERY, LITHIUM COIN TYPE 3.0V	\$	65.33
4419900	REGULATOR MODULE, VOLTAGE, 5VDC	\$	529.98
4420700	CAM, PUSHER	\$	373.88
4422200	SCREW 10-32X1.75 PAN SST PHIL	\$	2.51
4422300	STDF,M-F 10-32X3.50 3/8 HEX ALU	\$	51.65
4424700	KIT, WALL MOUNTING	\$	149.94
4425200	FTG, ELBOW .062IDTBGX10-32BLPOLY	\$	21.58
4425300	TBG TPE .093 OD .031 ID	\$	21.78
4425400	KIT, MAINT CHLORINE ANALYZER	\$	601.29
4425700	MAGNET, STIR BAR	\$	77.67
4427900	si IC EPROM CHLORINE CL17 REV P1.8	\$	227.41
4430094	MANUAL, HACH ONE PH ELCTD FRENCH	\$	31.41
4431200	FUNNEL, POWDER POLYPROPYLENE	\$	74.97
4431700	INSULATOR BLOCK	\$	724.63
4438900	PLATE, INSULATOR, 1 IN THK	\$	844.11
44417	SPADNS REAGENT SOLN, 3.78L	\$	278.24
44442	SPADNS REAGENT SOLUTION, 100ML	\$	36.86
44449	SPADNS REAGENT SOLN, 500ML	\$	79.50
44449-LM	SPADNS REAGENT SOLN, 500ML	\$	79.50
44453	SPADNS REAGENT SOLN, 1000ML	\$	114.36
44453-LM	SPADNS REAGENT SOLN, 1000ML	\$	114.36
4447169	AMMONIA IONIC STR ADJ PP PK/100	\$	103.35
4447169-LM	AMMONIA IONIC STR ADJ PP PK/100	\$	103.35
4447226	AMMONIA ELECTRODE FILL SOLN 50ML	\$	41.56
4450126	ELECTRODE FILLING SOLN, FLUORIDE	\$	50.40
4450126-LM	ELECTRODE FILLING SOLN, FLUORIDE	\$	50.40
4451569	SODIUM & POTASSIUM IONIC STR	\$	108.04
4451569-LM	SODIUM & POTASSIUM IONIC STR	\$	108.04
4451659	SODIUM ELECTRODE STOR SOLN, 50ML	\$	41.56
4453000	CKT BD ASSY, TOP BRD, 2100P BD	\$	1,580.30
4456369	NITRATE IONIC STRG ADJ PP PK/100	\$	189.71
4456369-LM	NITRATE IONIC STRG ADJ PP PK/100	\$	189.71
4457649	CALCIUM STD SOLN, 100MG/L 500ML	\$	105.88
4458200	CONN, .200 CTR PLUG 3CKT MARKED	\$	39.89
4458300	CONN, .200 CTR PLUG 2 CKT MARKED	\$	37.19
4459000	si FUSE, 1A 250V SLOW 5X20MM	\$	47.79
4459200	FUSE, 0.5A 250V SLOW 5X20MM IEC	\$	48.37
4468940	INSERT,MOLDED DIGITAL TITRATOR	\$	45.10
4472900	IC EPROM LAB PH/ISE METER	\$	258.24
4481500	PHOTOCELL, SILICON, UV-IR	\$	356.53
4490201	si INSTRUMENT CAP W/ O-RING	\$	82.87
4490300	COVER, PHOTOCELL CHAMBER	\$	82.87
4490400	MANIFOLD, PRESSURE SERIES 5000	\$	499.14
4490700	SAMPLE CELL SERIES 5000	\$	416.28
4491600	HINGE, ADJUSTABLE SWING PLASTIC	\$	118.91
4492600	si HOUSING, LAMP SERIES 5000	\$	568.52
4492800	si HOUSING, DETECTOR SERIES 5000	\$	1,576.45



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

4493600	STIR BAR, STAR HEAD,17MMDX13MMH	\$	246.68
4494200	VALVE, PRESSURE RELIEF 1/8 NPT	\$	339.19
4494500	FTG, TUBE 1/8 OD MALE KYNAR	\$	78.63
4494600	FTG, TUBE 1/4 OD MALE CONN KYNAR	\$	74.39
4494800	VALVE, SOLENOID, DIAPHRAGM, 2WAY	\$	1,183.30
4495300	FILTER, FOOT STRAINER 149 MICRON	\$	33.15
4495400	COLORIMETER ASSY SILICA	\$	6,591.02
4495401	COLORIMETER ASSY ORTHOPHOSPHATE	\$	5,319.07
4495403	COLORIMETER ASSY ALUMINUM	\$	5,473.25
4495430	Colorimeter Assembly, S5000 Silica	\$	6,591.02
4495431	COLORIMETER ASSY S5000 PHOSPHATE	\$	5,723.78
4495433	COLORIMETER ASSY S5000 ALUMINUM	\$	5,569.61
4495435	COLORIMETER ASSY S5000 LR PHOS	\$	6,957.19
4495500	MOTOR ASSY, STIRRER SERIES 5000	\$	1,102.36
4495600	THERMISTOR ASSY, SERIES 5000	\$	333.41
4495700	CAPILLARY MODULE W/3 RESISTORS	\$	1,865.53
4495703	CAPILLARY MODULE W/4 RESISTORS	\$	2,852.26
4495704	CAPILLARY MODULE W/2 RESISTORS	\$	1,015.63
4495730	CAPILLARY MODULE W/3 RESISTORS	\$	1,665.10
4495733	CAPILLARY MODULE W/4 RESISTORS	\$	2,004.29
4495734	CAPILLARY MODULE W/2 RESISTORS	\$	1,152.47
4495735	CAPILLARY MODULE W/3 RESISTORS	\$	1,665.10
4497200	CABLE ASSY, SIGNALS SERIES 5000	\$	1,505.14
4497300	CABLE ASSY, DRIVERS SERIES 5000	\$	1,204.50
4497400	CABLE ASSY, POWER & ALARMS	\$	1,125.48
4498200	CABLE ASSY, SERIAL, SERIES 5000	\$	1,574.52
4498900	TBG TEFLON .020 ID X .0625 OD	\$	17.54
4499300	REPLACEMENT LATCH, SSX ENCLOSURE	\$	211.99
4499600	CKT BD ASSY, PHOTOCCELL S5000	\$	1,279.66
4500200	TURB BODY ASSY, SS6	\$	1,462.74
4500400	LIGHT SOURCE ASSY, SS6	\$	3,006.43
4500800	BACK PANEL, SS6	\$	955.89
4501100	FTG, BULKHEAD 3/4FPTX3/4FPT PVC	\$	217.77
4501200	FTG, BULKHEAD 1FPT X1FPT PVC	\$	213.92
4501500	CLAMP, HOSE RANGE .984-1.14 PLAS	\$	29.10
4501600	CLAMP, HOSE RANGE 1.22-1.41 PLAS	\$	25.63
4502100	si CALIBRATION CUP	\$	362.31
4503200	COVER, REAR	\$	82.48
4503400	si LAMP ASSY, SS7	\$	294.86
4504100	RING, RTNG CIRCULAR .875 HSDIA CS	\$	30.45
4507300	VALVE, BALL, PVC 3/4 NPT	\$	381.59
4507600	TEMPLATE, LIGHT SOURCE ALIGNMENT	\$	140.11
4521700	FTG, ELBOW .125IDTBGX10-32BLPOLY	\$	31.22
45233	BUFFER SOLN, SULFATE TYPE 25ML DB	\$	20.60
45233-LM	BUFFER SOLN, SULFATE TYPE 25ML DB	\$	20.60
4528200	RISER CELL 10ML 1 IN SQ DR2000/2010	\$	56.27
4529900	SHIELD ASSY, SS6	\$	260.17

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

4530000	STAND, ELECTRODE HOLDER	\$	506.85
4530600	si ROD, ELECTRODE HOLDER	\$	127.39
4531500	STIR BAR, POLYGON 7/8X3/16 IN	\$	25.82
4539700	FOOT, RUBBER BLACK	\$	10.41
4542750	COVER ASSY, SAMPLE CELL REPLACE	\$	813.28
4542752	COVER ASSY, SAMPLE CELL REPLACE	\$	992.51
4542753	COVER ASSY, SAMPLE CELL REPLACE	\$	865.31
4542754	COVER ASSY, SAMPLE CELL REPLACE	\$	805.57
4542800	NUT, COMPRESSION, FTG, 1/8 TUBE	\$	43.75
4546100	O-RING FLCAR .139W X .609ID	\$	46.06
4546200	TBG TEFLON .170 ID X .250 OD	\$	55.70
4547300	KIT, TUBING REPLACE SERIES 5000	\$	776.66
4547303	KIT, TUBING REPLACE SERIES 5000	\$	714.99
4547400	FTG, TUBE 1/4 OD UNION TEE PLYPR	\$	73.81
4548600	FTG, BARB .062TBG X1/8-27NPT PLP	\$	34.11
4551401	REAGENT CAP W/O-RING	\$	16.57
4551500	si O-RING, CAPS, 3L., S-5000	\$	11.37
4552400	TBG POLYETHYL .125 OD X .062 ID	\$	12.53
4552400A	TBG POLYETHYL .125 OD X .062 ID, NOT FDA CONTROLLED	\$	7.52
4552800	O-RING FLCAR .060W X .056ID	\$	11.37
4554800	GAGE, FITTING SPACER	\$	92.51
4555100	si FUNNEL, GRAB SAMPLE SERIES 5000	\$	132.01
4555300	COVER, GRAB SAMPLE FUNNEL	\$	87.49
4557400	NUT, PANEL MOUNTING, REGULATOR	\$	59.94
4557500	REGULATOR, NON-RELIEVING	\$	838.33
4562700	KTO: SILICA REAGENT SET, S5000	\$	870.85
4562700-LM	KTO: SILICA REAGENT SET, S5000	\$	870.85
4563300	KTO: PHOSPHATE LR REAGENT SET, S5000	\$	897.95
4563900	KTO: PHOSPHATE HR REAGENT SET, S5000	\$	847.37
4565500	THERMOMETER, 0 TO 200 DEGREES C DIAL	\$	191.18
4566200	IC EPROM KIT S5000 OXYGEN SCAV	\$	705.36
4566700	FTG, TUBE 1/4 OD MALE POLYPROP	\$	43.55
4568500	oo BOX, LAB HEATER, MOLDED	\$	410.49
4568800	THERMOMETER, PLATINUM RESISTANCE	\$	231.26
4570000	DIPSTICK, POTASSIUM MEASURE	\$	95.78
4570500	DETECTOR ASSY, TRANSMITTED 2100P	\$	653.32
4573900	FLAT MIRROR HOUSING ASSEMBLY	\$	738.12
4580600	oo POWER SUPPLY 110VAC 7VDC 700MA	\$	807.50
4584500	BOX HINGED 4-9/16 SQ, 1-1/4 DEEP	\$	23.70
4584600	INSERT, FOAM- ACCUVAC	\$	47.02
4590008	MANUAL, BACTERIAL INCUBATOR	\$	196.57
4590011	MANUAL,SERV,BACTERIAL INCUBATOR	\$	133.75
45901	CAP, DISPENSER FOR 253335	\$	66.49
4591300	FUSE, 1A 250V SLO-BLO	\$	21.97
4591703	PANEL, REMOVABLE MOLDED	\$	164.78
4592200	VENT TRAP BODY, MOLDED	\$	140.11
4592300	ADAPTER, CAPILLARY FUNNEL-MOLDED	\$	116.02

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

4592400	VENT TRAP CAP, MOLDED	\$	68.03
4592500	BAFFLE, COLUMN, MOLDED	\$	82.48
4592600	RECEPTACLE, FRACTIONATING HEAD	\$	177.88
4597500	PUMP BELLOWS ASSY, REPLACEMENT	\$	1,273.88
4597800	TRANSMITTAL, HONEYWELL RECORDER	\$	85.37
4598000	REGULATOR, WATER, NON-RELIEVING	\$	618.63
4598100	VALVE, PRESSURE RELIEF, WATER	\$	512.64
4598300	KIT,SAMPLE PRES. COND'G SST	\$	11,929.37
4598400	REGULATOR ASSY, 316 SS	\$	6,860.83
4598500	FTG, CPRS 1/4 OD TBGX 1/4 NPT SS	\$	182.12
4598600	FILTER ASSY, TEE TYPE, BYPASS	\$	2,563.18
4598700	FLOW METER, 20-300 CC WATER/MIN	\$	479.87
4598800	FILTER ELEMENT KIT, REPLACEMENT	\$	749.68
4599400	si SHIELD, LIGHT SERIES 5000	\$	196.57
4599500	FTG, THRU .062 TUBING KYNAR	\$	33.34
4599600	FTG, CPRS 1/4 ODTUBE X1/8NPT SS	\$	165.55
46001	SODIUM THIOSULFATE, ACS 454G	\$	108.04
4600400	LID, SAMPLE CELL, 2100P BD	\$	101.56
4600500	COVER, BATTERY, 2100P BD	\$	54.15
4601422	DOCUMENTATION PACKAGE, DR/700	\$	429.77
4606600	SPACER, .375ODX.250IDX.250LG NYL	\$	7.71
4607500	BATTERY HOLDER, 4 CELL, SIZE AA	\$	44.90
4609093	CHEMICAL ANALYSIS INFO SPANISH	\$	49.34
4609095	CHEMICAL ANALYSIS INFO GERMAN	\$	49.34
460944	THERMISTOR 3K AT 25 DEG C	\$	211.99
460-S0467.70	BROMINE WATER 30 G/L 29ML	\$	114.36
460-S0468.70	PHENOL, 30 G/L 29ML	\$	74.97
460-S0890A.87LM	DPD FREE CHLORINE PP 25ML PK/100	\$	76.24
4611200	FILTER ELEMENT KIT, REPLACEMENT	\$	693.79
4611700	VALVE, NEEDLE, 1/8 MALE NPT TO 1/4 PUSH CONN	\$	691.86
4611800	FTG, CPRS 1/4 OD TBGX 1/8 NPT BR	\$	62.06
4612000	FTG, CPRS 1/4 OD TBGX 1/8 NPT BR	\$	23.13
4612100	VALVE, NEEDLE, 1/4 NPT, 316 SS	\$	953.96
4612200	NIPPLE, HEX PIPE, 1/4 NPT, 316SS	\$	107.54
4612300	FTG, TEE 1/4 NPT FEM ALL PORT SS	\$	427.84
4612500	REGULATOR, PRESSURE, RELIEVING	\$	1,167.88
4612600	KIT, REPLACE 3WAY W/ 2WAY VALVE	\$	2,871.53
4612700	FTG, BARB .170TBG X1/8 NPT KYNAR	\$	30.26
4613290	MANUAL,AGRICULTURAL CONB NITRATE	\$	27.37
4614000	TIP, FLEXIBLE, POLYPROPYLENE	\$	45.10
4615100	FTG, PLUG 1/4 NPT MALE STAINLESS	\$	105.23
4615500	VALVE, SOLENOID, PINCH, 2WAY NC	\$	967.45
4615900	si ASSEMBLY, SOIL SIEVE	\$	146.47
4616400	TBG SILICONE .031 ID .065 OD	\$	95.20
4620233	BINDER, CEL700, 3/4 INCH	\$	66.49
46212	SOAP SOLUTION(25ML SMPL) 50ML DB	\$	54.20
4624400	STRAP, VELCRO FASTENER, 8 IN.	\$	29.68

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

46249	SOAP SOLUTION(25ML SAMPLE) 500ML	\$	136.59
4630700	KIT, STRAIN RELIEF, 2/ .35 CABLE	\$	213.92
4637100	FUSE, .100A 250V 5MMX20MM GLASS	\$	46.06
4639100	IC EPROM KIT EC1000/PH	\$	809.42
4643600	KIT, FLOWMETER & 1/4 OD TBG FTGS	\$	849.90
4644200	STDF, M-F 4-40X.187 3/16 HEX SS	\$	10.41
4645200	GASKET, COVER ASSY, SAMPLE CELL	\$	86.34
4646800	TRANSMITTAL, SAMPLE PUMP	\$	85.37
4650200	HOUSING, BOTTOM 2100P	\$	549.25
4650300	HOUSING, TOP 2100P	\$	383.51
4650400	DISPLAY,LCD CUSTOM FOR 2100P	\$	383.51
4650500	KEYBOARD, 3X3 X-Y MATRIX, 2100P	\$	431.69
4651000	EOL FILTER ,LIGHT ,2100P	\$	148.20
4651000L1	FILTER, LIGHT, 2100Q	\$	132.59
4652700	COVER, OPTICS 2100P	\$	304.50
4653000	SHIELD, ESD, 1900C, 2100P BD	\$	83.83
4653100	LID ASSY, 2100P	\$	260.17
4653400	CKT BD ASSY, BOTTOM BRD, 2100P	\$	634.05
4653700	DETECTOR ASSY, 90 DEGREE 2100P	\$	591.65
4653900	LAMP ASSY, 2100 PORTABLE TURBS	\$	356.53
4654500	JACK, POWER, 2 CNDCT, .10 DIA	\$	68.61
4654600	SCREW 4-40X .25 PAN SST TPG	\$	8.09
4656100	BATTERY, LITHIUM 1/2 AA 3.0V	\$	116.79
4660000	TUBE,PLASTIC *SELL-04	\$	15.61
4660004	si TUBE,PLASTIC-PKG/4	\$	60.13
4660100	CASE ASSY, HC-1	\$	41.82
4660200	CASE ASSY, HC-2	\$	64.37
4660400	CASE ASSY, HC-4	\$	119.29
4660500	CASE ASSY, HC-5	\$	152.63
4660600	si CASE ASSY, HC-6	\$	153.02
4660700	si CASE ASSY, HC-7	\$	206.21
4660800	si CASE ASSY, HC-8	\$	535.76
4661500	FTG, BARB .375TBG X1/4 NPT WT PP	\$	13.30
4661600	STRAINER, SEDIMENT PVC 30 MESH	\$	992.51
4661700	si STRAINER,MODIFIED,BYPASS,SAMPLE	\$	1,025.27
4661800	Y-STRAINER KIT, 30 MESH, BYPASS	\$	1,422.27
4661900	FTG, ELBOW, .062TBGX.125TBG POLY	\$	23.90
4662400	FILTER ELEMENT	\$	120.84
4663900	CONN, AC RCPT W/ON-OFF SWITCH	\$	184.82
4666700	CLAMP, HOSE WORM DRIVE 316 SST	\$	36.23
4666800	CLAMP, HOSE WORM DRIVE 316 SST	\$	36.23
4667600	FTG, PLUG THREAD 3/4-14 NPT PVC	\$	28.14
4667800	si CAP, BUBBLE TRAP TOP	\$	343.04
4668000	si BUBBLE TRAP/ HEAD REGULATOR	\$	1,401.07
4668500	HANGER, PIPE, 3" IRON PIPE SIZE	\$	103.11
4668800	VALVE, NEEDLE 1/2" NPT PVC	\$	410.49
4669100	KIT, SSX TUBING REPLACEMENT	\$	422.06

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

4669222	KIT, SSX AUTO FLUSH 220 VAC	\$ 11,524.66
4671300	ff FILTER, INTERFERENCE 528 NM	\$ 196.57
4674300	BATTERY, ALKALINE AAA 1.5V PK/4	\$ 54.35
4679500	KIT, RETROFIT S-5000 PINCH VALVE	\$ 1,424.20
4681300	DIPSTICK, CYANURIC ACID MEASURE	\$ 92.89
4681400	DIPSTICK, SULFATE MEASURE	\$ 113.13
4683600	CORD,POWER 3X0.75MM 6A 250V 6FT	\$ 65.33
4685900	si oo CAP-SEALING	\$ 94.82
4686300	VALVE,CHECK CARTRIDGE .250 OD	\$ 125.27
468645	THERMAL COMPOUND, TYPE 120	\$ 154.95
4689400	CABLE, 24 AWG 4 COND W/SHIELDING	\$ 14.07
4692400	FTG, TUBE 1/4 OD MALE POLYPROP	\$ 31.03
4693400	FTG, Y FITS .062 ID TUBE POLYPRO	\$ 24.48
4693900	si DRAIN/INLET BLOCK	\$ 493.36
4694200	si CAP,INSERT 2-PORT FOR .062TBG	\$ 58.01
4694300	si CAP, MODIFIED	\$ 16.96
4696400	CORD, POWER 18/3 10A-125V 6FT	\$ 120.64
4697400	SENSOR ASSY, SAMPLE PRESSURE	\$ 1,742.19
4697430	SENSOR ASSY, SAMPLE PRESSURE	\$ 1,670.88
4697500	SENSOR ASSY, REAGENT PRESSURE	\$ 1,740.26
4697530	SENSOR ASSY, REAGENT PRESSURE	\$ 1,605.36
4697600	FTG, TUBE 1/4 OD TUBE/TUBE BULKH	\$ 96.75
4697900	si LAMP ASSY SERIES 5000	\$ 187.90
4698100	KIT, MAINTENANCE SERIES 5000	\$ 1,886.73
4698133	KIT, MAINTENANCE SERIES 5000	\$ 1,476.24
4698200	KIT, TUBING REPLACE SERIES 5000	\$ 693.79
4698233	KIT,TUBING REPLACE SERIES 5000	\$ 894.22
4698300	FUSE, 2.50 A 250V 5MM X 20MM	\$ 21.97
4699100	KIT,SAMPLE PRES.COND'G. S-5000	\$ 2,794.44
4699200	TBG,SILICONE 1/16 ID 1/8 OD	\$ 45.10
4700800	CABLE ASSY, 24 COND W/ 24 CKT 2	\$ 210.06
4700900	CABLE ASSY, 16 COND W/ 16 CKT 2	\$ 300.64
4701000	DETECTOR ASSY, BACK SCATTER	\$ 549.25
4701200	DETECTOR ASSY, 90 DEG, LAB TURB	\$ 888.44
4701900	COVER, OPTICS 2100N/2100AN	\$ 585.87
4702000	HOUSING, TOP 2100N	\$ 703.43
4702100	KEYBOARD OVERLAY, 2100N LAB TURB	\$ 269.81
4702300	HOUSING, BOTTOM 2100N/2100AN	\$ 491.44
4702400	VENT, POWER SUPPLY 2100N/2100AN	\$ 148.97
4702600	SHIELD, HEAT 2100N/2100AN	\$ 154.08
4702900	HOUSING, TOP 2100AN	\$ 703.43
4703000	DUST COVER, 2100N & 2100AN	\$ 275.59
4703200	PANEL, LAMP ACCESS 2100N/2100AN	\$ 149.74
4705100	LENS, PLASTIC SMALL 2100N/2100AN	\$ 335.33
4705200	LENS, PLASTIC LARGE 2100N/2100AN	\$ 410.49
4705400	PANEL, PRINTER ACCESS 2100AN	\$ 210.06
4705500	PANEL, PRINTER GLUED, 2100N	\$ 121.03

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

4706000	RETAINER, PRINTER 2100AN	\$	132.21
4706800	COVER, FLOW CELL LIGHT	\$	213.92
4707100	KEYBOARD OVERLAY, 2100AN	\$	354.60
4707300	SCREW 6-32X .5 HEX SST TPG	\$	5.59
4707600	CLOTH, OIL APPLICATION	\$	55.50
4707700	FOOT, 2100N/2100AN LAB TURB	\$	20.62
4708900	LAMP REPLACEMENT ASSY, LAB TURB	\$	618.63
4709000	PRINT PAPER, 1 PKG OF 5 ROLLS	\$	113.51
4709500	KIT, FLOW CELL GLASS	\$	686.08
4709900	ASSY, FLOW CELL STAND	\$	1,023.34
4711200	FUSE, 250MA, 250V	\$	18.31
4711600	TRANSMITTAL, CONDUCTIVITY METER	\$	82.87
4711700	FTG, TEE 3/4 NPT FEMALE MODIFIED	\$	154.18
4713288	INSTRUCTION SHT, PH SENSOR	\$	45.10
4721400	FUSE, 0.80A 250V SLOW 5X20MM IEC	\$	31.22
4721500	FUSE, 1.5A SLOW 5X20MM	\$	31.22
4721600	KIT, MAINTENANCE CL17-PRE APR '00	\$	892.29
4723200	CAP PLUG, RUBBER, 16MM	\$	12.53
4726000	KIT, PROPELLANT VALVE REPLACEMNT	\$	1,318.20
4726500	LATCH, NONKEYED	\$	420.13
4726600	LATCH, KEYED	\$	695.72
4728000	BRACKET, MOUNTING SAMPLE COND.	\$	76.32
4729000	SENSOR, PRESSURE, 0-15 PSIG 30 MAX	\$	1,406.86
4729200	KIT, SENSOR REPLACEMENT, S-5000	\$	2,062.10
4729900	FTG, TUBE 1/4 ODX 1/8 NPT FEM PP	\$	43.75
4730200	VALVE, SOLENOID PINCH, 2WAY NC	\$	1,214.14
4730400	FTG, TUBE 1/4 OD TUBE-MNPT ELBOW	\$	61.28
4730500	FTG, TUBE 1/4 OD TUBE ELBOW	\$	113.13
4730600	FTG, TUBE 1/4 OD TUBE ELBOW	\$	49.91
4730700	FTG, TUBE 1/8 OD TUBE & MALE NPT	\$	44.33
4730800	VALVE, 3 WAY PUSH BUTTON 303 SS	\$	637.90
4733600	TOOL, STIR BAR REMOVER, CL17	\$	33.15
4734200	FTG, TEE .125ID TBG, .062ID TBG	\$	21.58
4734600	HANGER, PIPE .750 PIPE SIZE	\$	45.10
4734700	FTG, TUBE 1/4 OD 1/8NPT STRAIGHT	\$	57.43
4734900	CLAMP, CONDUIT HANGER, 1 INCH	\$	36.04
4735100	REGULATOR, NON-RELIEVING	\$	609.00
4735800	REGULATOR, NON-RELIEVING, PRESET	\$	562.74
4737400	INSERT, PUNCHED TOP FF-2	\$	71.50
4739100	RING, NUMBER, BLACK	\$	19.08
4743500	TRANSMITTAL, SERIES 5000	\$	97.32
4743800	TBG, POLY .250 OD .040W	\$	17.54
4744800	KIT, FLOW CELL END CAPS	\$	832.55
4744900	KIT, MANUAL FLOW CELL	\$	5,087.81
4745300	FTG, ELBOW, SWIVEL .125 TBG NYL	\$	73.81
4746700	CABLE ASSY, CPU BD TO REMOTE	\$	129.89
4746800	CABLE ASSY, 2 COND 24 AWG STRD	\$	520.34



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

4748000	FTG, BULKHEAD 1/8 ID TBG NYLON	\$	48.76
4748200	JACK, PHONE .177 3 CONDUCTOR	\$	33.73
4748300	FTG, BARB 1/8 TBG X 10-32 BLKNYL	\$	17.73
475024	FRO SWITCH SPDT 2FB63-73 2456 EA	\$	177.30
4753300	PLUNGER, SPRING LOADED	\$	57.05
4753400	vw STRAP, VELCRO FASTENER, 30 IN	\$	27.37
4754900	si ff ASSY, STAND, HQ SERIES METER	\$	170.17
4755000	PHOTODIODE, SILICON, HIGH SPEED	\$	171.33
4760200	si CKT BD ASSY, POWER SPLY BD 5000	\$	2,910.07
4760205	si CKT BD ASSY, PWR SPLY LR ORPHOS	\$	3,199.15
4760410	aa CKT BD ASSY, MPU 5000 SILICA	\$	7,593.17
4760411	CKT BD ASSY, MPU 5000 ORPHOS	\$	7,593.17
4760412	CKT BD ASSY, MPU 5000 HYDRAZINE	\$	7,593.17
4760413	CKT BD ASSY, MPU 5000 ALUMINUM	\$	6,918.65
4760414	CKT BD ASSY, MPU 5000 OXY SCAVEN	\$	5,357.62
4760415	CKT BD ASSY, MPU 5000 LR ORPHOS	\$	6,012.86
4761000	SWITCH, TOGGLE DPDT POWER	\$	237.05
4761100	RUBBER BOOT, TOGGLE SWITCH	\$	63.40
4761400	CABLE ASSY, MOTOR/LAMP/SOLENOID	\$	427.84
4765000	KIT, S5000 CONTROL BOX SILICA	\$	47,544.02
4765001	KIT, S5000 CONTROL BOX PHOSPHATE	\$	44,383.42
4765002	KIT, S5000 CONTROL BOX HYDRAZINE	\$	46,368.43
4765003	KIT, S5000 CONTROL BOX ALUMINUM	\$	32,299.87
4765300	KIT, INSTL SERIES 5000 SILICA	\$	5,241.98
4765400	KIT, INSTL SERIES 5000 PHOSPHATE	\$	5,723.78
4765600	KIT, INSTL SERIES 5000 ALUMINUM	\$	3,623.14
4765800	KIT, INSTL SERIES 5000 LR PHOS	\$	5,338.34
4765900	KEYBOARD/COVER ASSY, SERIES 5000	\$	2,351.18
4766100	oo INSERT, FOAM, SMALL CUSTOM	\$	75.74
4766300	INSERT, FOAM, LARGE CUSTOM	\$	103.68
4766600	CASE ASSY, LARGE CUSTOM KIT	\$	612.85
4766703	oo KIT, SERIES 5K DSPLY/MPU REPL ALUM	\$	3,892.94
4769000	IC EPROM KIT 2100N CPU BD	\$	634.05
4770200	PRINTER MODULE ASSY, 2100AN	\$	1,507.07
4780800	si CASE ASSY, MEL PA	\$	1,624.63
4780900	CASE ASSY, MEL MPN	\$	1,682.45
4781000	CASE ASSY, MEL MF	\$	2,274.10
4788000	CELL HOLDER ASSEMBLY, TTHM	\$	452.89
4791200	KIT,PROM SET/MANUAL S5000 SILICA	\$	1,318.20
4791600	KIT,PROM SET/MANUAL S5000 ORPHOS	\$	1,086.94
4792400	KIT,PROM SET/MANUAL S5000 ALUMIN	\$	950.11
4793200	KIT,PROM SET/MANUAL S5K LR ORPHO	\$	1,304.71
47949	HARDWATER, 20GPG 500ML	\$	193.32
4809200	MOTOR, STEPPER UNIPOLAR 12VDC	\$	400.86
4813301	DOOR, LAMP ACCESS PAINTED	\$	277.52
4813600	SWITCH, INTERLOCK SAFETY	\$	210.06
4822800	SAMPLE CELL, 10MM QUARTZ PAIR	\$	1,167.88



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

4824900	KIT, TUBING REPLAC INLET & WASTE	\$	152.06
4825400	SENSOR TIP, PH FLAT REPLACEMENT	\$	869.17
4826300	CABLE ASSY, MALE, QUICK	\$	522.27
4827600	EOL ff ELECTRODE, INDUSTRIAL PROCESS PH	\$	3,411.14
4827700	CABLE, SHORTING	\$	69.76
4830700	KEYBOARD/COVER, SIGNET CONDUCT.	\$	1,171.74
4830800	KEYBOARD/COVER, SIGNET PH	\$	1,171.74
4831288	MANUAL, CL/PH POCKET COLORIMETER	\$	499.14
4831604	MANUAL, DR4000, SERVICE V2.0	\$	100.21
4840300	RISER, CELL 10ML 1 IN SQ. DR3000	\$	60.90
4846018	MANUAL, INSTR, DR/870 CLRMTR	\$	183.47
4846022	MANUAL, PROC, DR/870 COLORIMETER	\$	223.56
4846400	ADAPTER ASSY, COD DR/800 CLRMTR	\$	98.09
4853600	SENSOR, THERMOCOUPLE	\$	479.87
4853700	THERMOSTAT, BI-METALIC	\$	470.24
4853800	CONTROLLER, PID 12-24 VAC	\$	2,409.00
4853900	RELAY, SOLID STATE, 32VDC, 25AMP	\$	356.53
4854200	CIRCUIT BREAKER, 2 POLE, 10 AMPS	\$	834.48
4854300	CIRCUIT BREAKER, 1 POLE, 20 AMPS	\$	202.36
4855200	TRANSMITTAL, SAMPLE HEATER	\$	85.37
4855500	RETROFIT KIT, SAMPLE HEATER	\$	520.34
4864300	si SAMPLE CELL, 1-CM/10-ML, DR/800	\$	36.23
4864302	SAMPLE CELL W/CAP, 1CM/10ML PK/2	\$	75.35
4866088	ELECTRODE MANUAL, DOUBLE JCT REF	\$	82.87
4867088	MANUAL, SODIUM ELECTRODE	\$	40.66
4879900	RACK, 1CM MICRO CUVETTE	\$	184.05
4879910	RACK, 12-mm MICRO CUVETTE	\$	47.02
4883400	FUNNEL, FLOW THRU ASSY NO LOGO	\$	346.90
490=001=005	0.5L PE BOTTLE	\$	186.36
490=001=040	FLASK POLYETHYLENE 1L, SODIUM	\$	106.57
490=010=011	DIPA BOTTLE	\$	310.28
4901700	RES FHS 1/8W 5% 10 MOHM 0805	\$	6.55
491=001=001	RED CAP FOR 0.5L BOTTLE	\$	4.05
4911000	si BAND, DIFFUSER DR800	\$	46.06
49200	SPOON, MEASURING 0.05G	\$	10.60
4934800	CAP, ARSENIC TEST KIT	\$	15.03
4942500	CASE ASSY, DR/800 & DR900 CLRMTR.	\$	352.68
4943000	CEL PORTABLE CARRYING CASE	\$	903.86
4951700	GLOVES, INSPECTION, COTTON	\$	18.31
4952600	si FUSE, 2.50 A 250V 5MM X 20MM	\$	31.61
4961800	STRAINER, SEDIMENT 120 MESH	\$	753.54
4961900	FILTER ELEMENT, 120 MESH	\$	429.77
4962000	COUPLING, 3/4 NPT X 3/4 NPT	\$	46.64
4965900	KIT, ACCESSORY, ORIFLO	\$	51.46
4966000	si ASSY, ORIFLO	\$	572.38
4972388	MANUAL, SILICA, LR POCKET CLRMTR	\$	67.07
4972488	MANUAL, CYANIDE POCKET CLRMTR	\$	67.07



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

4972688	MANUAL, SILVER POCKET CLRMTR	\$	67.07
4972988	MANUAL, SULFATE POCKET CLRMTR	\$	67.07
4973077	MANUAL, IMMUNOASSAY PKT CLRMTR	\$	70.15
4973088	PROC., IMMUNOASAY PKT CLRMTR	\$	80.94
4973988	MANUAL, NITRATE POCKET CLRMTR	\$	67.07
4994	BRIDGE RECTIFIER,DF04S 400	\$	14.65
4D1003	TRANSDUCER,STAINLESS	\$	9,713.09
4F0393	PLUG, HOLE .87 DIA	\$	5.40
4F1017	CAPLUG, BLACK	\$	20.24
4F1059	CAP PLUG, BLACK VINYL	\$	6.17
4H0891	FITTING,TUBE,BRASS,1/4NPT	\$	48.76
4H1135	CONN,MALE,THERMOCOUPLE,3/	\$	734.26
4H1135-102	NUT FOR 99X4H1135,3/4" S.	\$	98.09
4H1234	CONN.MALE,POLY,1/2NPTX1/2	\$	27.37
4H1246	VALVE,3-WAY,1/4NPT,BRASS	\$	204.28
4H1252	FTG, STRAIGHT, 1/4NPT X 1/4" TUBE	\$	34.30
4H1259	FITTING,1/4"NPTx1/4" HOSE	\$	45.67
4H1278	CONN,MALE,THERMOCOUPLE,1	\$	356.53
4H1285	CONN.MALE,THERMOCOUPLE,1/	\$	493.36
4H1312	NOZZLE,SPRAY,1/8"NPT,SS.	\$	341.11
4H1313	CPLG.QUICK,FEMALE,1/4" I.D	\$	108.31
4H1314	CPLG,QUICK,MALE,1/4"NPT	\$	70.34
4H1315	CPLG,QUICK,MALE,1/4" BARB	\$	65.14
4H1316	1/2 "NOM." PKG OF 10 MSC NYLON DBL BND CLAMPS	\$	24.48
4H1331	COIL,120V(FOR USE ON 4H13	\$	258.24
4H1332	CONN,MALE,SS,1/2 BT X 3/4	\$	664.88
4H1408	NOZZLE,SPRAY,KYNAR,1/4 X	\$	183.66
4H1410	FITT,BARB,1/4"X 1/4"KYNAR	\$	65.14
4H5011	ADAPTER,1/4"NPTx1/16"Tube	\$	11.37
4S1021	SWITCH,SLIDE,DPST,QUAD	\$	65.14
5000TC-00-1-000	Cond. Sensor, POLYPRO	\$	2,948.62
5000TC-00-3-000	Cond. Sensor, PEEK	\$	3,758.04
5000TC-00-4-000	Cond. Sensor, PFA Teflon	\$	4,027.85
5000TC-00-5-000	Cond. Sensor, Sanitary	\$	3,391.87
5000TC-00-7-000	Cond.Sensor,Sanitary PFA	\$	4,432.56
5000TC-11-1-000	Cond. Sensor & Hardware	\$	3,565.32
5000TC-11-3-000	Cond. Sensor & Hardware	\$	4,394.02
5000TC-11-4-000	Cond. Sensor & Hardware	\$	4,625.28
5000TC-12-4-000	3728E2T Snsr & MH462G HDW	\$	4,586.74
5000TC-21-3-000	Cond. Sensor & Hardware	\$	6,436.85
5000TC-21-4-000	Cond. Sensor & Hardware	\$	6,899.38
5000TC-22-4-000	Cond. Sensor & Adaptor	\$	6,147.77
5000TC-24-1-000	Cond. Sensor & Hardware	\$	3,526.78
5000TC-24-4-000	Cond. Sensor & Hardware	\$	4,509.65
5000TC-25-2-000	Cond. Sensor & Hardware	\$	5,376.89
5000TC-25-4-000	Cond. Sensor & Hardware	\$	5,993.59
5000TC-41-4-000	Cond. Sensor & Hardware	\$	10,503.24

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

5000TC-42-0-000	Cond. Hardware for Mtg.	\$	2,736.62
5000TC-42-1-000	Cond. Sensor & Hardware for Mtg.	\$	3,122.06
5003000	si SAFETY SHIELD, LAB BENCH 15 X 16	\$	1,746.04
5004000	si SAFETY SHIELD, LAB BENCH 30 X 16	\$	2,274.10
50041H	BEAKER, LOW FORM 50ML	\$	23.51
50042H	BEAKER, LOW FORM 100ML	\$	24.67
50044H	BEAKER, LOW FORM 150ML	\$	25.63
50046H	BEAKER, LOW FORM 250ML	\$	26.60
50048	BEAKER, LOW FORM 400ML	\$	28.14
50052	BEAKER, LOW FORM 600ML	\$	29.68
50053	BEAKER, LOW FORM 1000ML	\$	50.69
50054	BEAKER, LOW FORM 2000ML	\$	88.65
50071H	BEAKER, LOW FORM 50ML PK/12	\$	196.57
50072H	BEAKER, LOW FORM 100ML PK/12	\$	196.57
50074	BEAKER, LOW FORM 150ML PK/12	\$	196.57
50076H	BEAKER, LOW FORM 250ML PK/12	\$	196.57
50078	BEAKER, LOW FORM 400ML PK/12	\$	225.48
50082	BEAKER, LOW FORM 600ML PK/6	\$	136.83
5008240	Delivery tubing with anti-diffusion	\$	466.38
5008260	STOPPER, TITR BTL W/DESICCANT&TUB.	\$	210.06
50083	BEAKER, LOW FORM 1000ML PK/6	\$	240.90
50084	BEAKER, LOW FORM 2000ML PK/4	\$	300.64
5008500	STIR BAR, F/AUTOCAT 9000	\$	35.85
500-P2553.88	LAMP KIT UV, 115V 60 HZ	\$	2,062.10
500-P2554.88	rr UV LAMP KIT, 230V	\$	2,351.18
504-3004	CAP,CER.,.1 MF,X7R,0.2LS,RADL,2	\$	7.52
504-3033	CAP,CER.,220PF,50V, 5%,NPO,RDL	\$	7.52
50444	BURET TEE	\$	65.72
50445	BULB & TUBE	\$	87.30
50449	STOPPER, RBR 2-HOLE F/22614 25ML	\$	15.61
50450	ff STOPPER, RBR 2-HOLE	\$	10.60
50540	ff FLASK, ERLNMEYER 25ML	\$	24.67
50541	FLASK, ERLNMEYER 50 ML	\$	26.60
50543	FLASK, ERLNMEYER 125ML	\$	25.63
50546	FLASK, ERLNMEYER 250ML	\$	26.60
50547	FLASK, ERLNMEYER 300ML	\$	28.14
5057	CLIP,TUBE HOLDER	\$	63.98
50571	FLASK, ERLNMEYER 50ML PK/12	\$	215.85
50573	FLASK, ERLNMEYER 125ML PK/12	\$	215.85
50576	FLASK, ERLNMEYER 250ML PK/12	\$	235.12
50655	FILTER PAPER, 9CM DIA PK/100	\$	55.12
50656	FILTER PAPER, 11CM DIA PK/100	\$	56.27
50657	FILTER PAPER, 12.5CM DIA PK/100	\$	63.79
50658	FILTER PAPER, 15CM DIA PK/100	\$	85.37
50659	FILTER PAPER, 18.5CM DIA PK/100	\$	123.15
50837	CYLINDER, GRADUATED 5ML TD WHITE	\$	87.30
50838	CYLINDER, GRADUATED 10ML TD WHITE	\$	76.32

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

50840	CYLINDER, GRADUATED 25ML TC WHITE	\$	79.02
50841	CYLINDER, GRADUATED 50ML TC WHITE	\$	79.02
50842	CYLINDER, GRADUATED 100ML TC WHITE	\$	94.63
50846	CYLINDER, GRADUATED 250ML TC WHITE	\$	136.83
50849	ff CYLINDER, GRADUATED 500ML TC WHITE	\$	192.14
50853	CYLINDER, GRADUATED 1000ML TC WHITE	\$	269.81
51000	SPOON, MEASURING 1.0G	\$	12.33
5101244	CARD, QUICK REF NITRATE/AMMONIA	\$	26.98
5102100	MIXING MODULE, 1ML APA 6000	\$	2,794.44
5102101	MIXING MODULE, 3ML APA6000	\$	2,620.99
5102102	MIXING MODULE 2ML, APA6000	\$	3,199.15
5102200	FACEPLATE ASSY, APA6000	\$	1,406.86
5102400	VALVE ASSEMBLY, 2 WAY APA 6000	\$	1,505.14
5103400	FILTER, SUPPORT ELEMENT	\$	292.93
5103600	GROMMET, MANIFOLD APA6000	\$	18.69
5103800	O-RING, VITON, FILTER SEAL	\$	252.46
5104000	INSTALLATION KIT, APA6000	\$	1,283.52
5104200	FILTER, 0.2 MICRON MEMBRANE	\$	223.56
5104603	TEMP CONTROL BLK ASSY, HR HARDNESS	\$	1,505.14
5104800	si REAGENT HOUSING ASSEMBLY APA	\$	479.87
510514	CLIP, PIPE HANGER, 1/2""	\$	33.73
5105200	INSTALL/MAINT KIT MICRO FILTER	\$	1,002.14
5105400	FILTER, MODEL 9150 230V 50-60HZ	\$	19,214.18
5105600	O-RING FLCAR .070W X .070 ID	\$	7.52
5107300	si CHAMBER, MIXING, 3ML VOL	\$	254.39
5107801	si CKT BD ASSY, LED 605 NM , APA 6000	\$	358.46
5107803	si CKT BD ASSY, LED 650 NM, APA 6000	\$	377.73
5109601	si CKT BD ASSY, CONTROLLER #2 BD, APA 6000	\$	4,509.65
51100	SPOON, MEASURING 0.1G	\$	10.60
5110400	BULKHEAD UNION .25 OD TBG TO TBG	\$	40.66
5110600	REAR COVER, MODULE BASE, APA6000	\$	84.99
5111900	FILTER, INTERFERENCE 600NM	\$	547.32
5112900	FTG, HEX NUT .062 OD TBG	\$	75.74
5113000	FTG, FERRULE .062 OD TBG (VALCO)	\$	90.96
5113100	FTG, FLANGELESS NUT .062 OD TBG	\$	59.17
5113200	FTG, FERRULE .062 OD TBG	\$	55.89
5114000	PLUG, FOR 1/4" TUBING	\$	18.69
5114300	si INSTRUMENT PANEL	\$	2,331.91
5115000	VALVE, ROTARY 10 PORT W/ACTUATOR	\$	7,015.01
5115500	SWITCH, ROCKER DPST PANEL MNT	\$	36.23
5115900	TBG POLYETHEYLENE .500 OD .062W	\$	11.95
5116000	PACKAGING ASSY, APA6000	\$	7.32
5116100	si CABLE ASSY, BACKPLANE/CONTROLLER	\$	470.24
5116400	si CABLE ASSY, AQUATREND APA	\$	628.27
5117800	TBG TFZL .030 IDX.062 OD 2000PSI	\$	25.05
5119200	FTG, INSERT, 1/4 OD TBG	\$	23.90
51200	FILTERING THIMBLE	\$	126.42

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

5121500	ASSY, VENT CAP	\$	175.38
5122200	SEAL, PH, APA 6000	\$	50.11
5122400	si PLUG, MANIFOLD APA6000	\$	26.98
5123900	si PIPE, PRECUT DRAIN, 1IN DIA PVC	\$	74.58
512432	LENS DOUBLE CONVEX	\$	200.43
5124400	IC PROG KIT, APA 6000 FLASH MEM	\$	238.97
5124600	FTG, TUBE, CONNECTOR,MALE 1/4TBG	\$	19.66
5125018	oo MANUAL, SIGNAL OUTPUT MODULE	\$	281.37
5125400	BRACKET, REAGENT BOX,LEFT REAR, APA 6000	\$	121.03
5126100	FTG, TUBE, CONNECTOR,MALE 1/4TBG	\$	25.25
5126200	FTG, TUBE, CONNECTOR,MALE 1/2TBG	\$	58.59
5126700	BRACKET, REAGENT BOX, RIGHT REAR, APA 6000	\$	118.72
5129100	KIT, TOOL FOR ALL APA6000	\$	2,177.74
5129300	PLUG, 1/4-28 FLAT BOTTOM TEFZEL	\$	57.43
5134500	ALLEN WRENCH, 7/64 INCH, LONG	\$	5.01
5134700	STATOR REPLACEMENT FOR 51150-00	\$	4,316.93
5134800	ROTOR REPLACEMENT FOR 51150-00	\$	892.29
5135600	TUBING CUTTER FOR 1/16" TEFZEL	\$	183.08
5135700	WRENCH, 1/4 INCH FOR VALCO VALVE	\$	131.05
5135900	FTG, UNION 1/4-28 TEFZEL	\$	56.27
5139300	FTG, HEX NUT VALVE PLUG 10-32	\$	99.44
5139500	VALVE, BALL, PVC 1/4 NPT	\$	72.08
5139600	STRAINER, HOUSING,3/4"W/O FILTER	\$	277.52
51400	TEST TUBE, SOAP	\$	21.78
5141000	U-BOLT 1/4-20 X 2.843	\$	35.85
5141700	CABLE, STOP, WITH EYELETS	\$	37.77
5143000	CAP, ALKALINITY REAGENT	\$	11.37
5143100	ASSEMBLY, GAS SCRUBBER	\$	516.49
5143200	ASSY, SCRUBBER VENT CAP	\$	653.32
5145018	MANUAL, SIGNAL INPUT MODULE	\$	139.14
51451255-001	si Interconnect Cable	\$	9.25
5150700	VALVE BODY REPLACEMENT 10 PORT	\$	5,511.79
5160500	SEAL, SENSOR CAP, LDO PROBE	\$	27.75
5167000	COLLAR, SHROUD, LAB LDO PROBE	\$	27.56
5167300	STIRRER FRAME, LAB LDO PROBE	\$	127.20
5168200	BOD ACCESSORY KIT, FUNNEL/STIRRER	\$	398.93
517937	TBG TYGON 3/16" OD X 1/32" W	\$	20.24
518037	TBG TYGON 5/16" OD X 1/16" W	\$	29.68
518137	TBG TYGON 3/8" OD X 1/16 W	\$	36.81
518237	TBG TYGON 7/16" OD X 1/16" W	\$	52.81
518637	TBG TYGON .500 OD .062 W	\$	74.97
5187500	si FILTER CAP, LAB LDO PROBE	\$	7.52
5188906	DESICCANT CARTRIDGE, PK/6	\$	58.71
518937	TBG TYGON 1/8" OD X 1/32" W	\$	22.93
5189900	CHAMBER, L.I.S. ASSEMBLY	\$	520.34
5190300	si CASE ASSY, HARD-SIDED, PORTABLE	\$	655.25
5191700	CASE, HARD SIDED, SENSION156	\$	423.98

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

5192588	MANUAL, PT SERIES SODIUM	\$	131.44
5192711	AMMONIA MEMBRANE REPLACEMENT KIT pk5	\$	528.05
5192900	BOTTLE, SOAKER	\$	19.85
5193088	MANUAL, ELECTRODE, PT SERIES PH	\$	89.23
5193588	MANUAL, ELECTRODE, GEL PH W/TEMP	\$	131.44
5194100	SOLUTION, ELECTRODE FILLING 65ML	\$	44.45
5195000	TEST POINT, SMT 105X40 TP-108	\$	8.09
5195588	MANUAL, COMB. FLUORIDE ELECTRODE	\$	75.74
5195900	SEAL, PERIMETER, LDO	\$	24.67
5196000	BOD STIRRER ACCESSORY	\$	616.70
5197088	INSTRUCTION CARD, DISSOLVED OXYG	\$	65.14
5197100	KTO: KIT, BOD ACCESSORY	\$	928.91
5197200	si ADAPTER, BOD FUNNEL	\$	397.00
5197300	MEMBRANES, DO FOR 51970-00 2/PK	\$	379.66
5197588	INSTRUCTION CARD, CONDUCTIVITY	\$	66.87
5198088	MANUAL, TEMPERATURE PROBE	\$	120.26
5198100	oo CAP, DO, GEL PH, TEMP PROBES	\$	20.62
520001	DESICCANT,SILICA GEL	\$	48.18
5200500	nn PHOTOCELL ASSY 1720D SELL 5218000	\$	1,023.34
5200600	1720D TURB BODY&BUBBLE TRAP ASSY	\$	909.64
5201200	si LID, BUBBLE TRAP, 1720D	\$	54.15
5201300	si PLATE, BOTTOM, 1720D	\$	37.39
5201700	si CAP BOTTOM, 1720D	\$	154.37
5201800	si CAP TOP, 1720D	\$	213.92
5202500	BUSHING, STRAIN RELIEF	\$	46.64
5202800	si GASKET, 1720D CAP	\$	50.88
5203300	SPRING, SAMPLE CELL 2100N	\$	42.40
5203800	IC PROG KIT, 1720D (CUSTOMER)	\$	1,083.09
5204700	OPTICAL HEAD, FT660sc	\$	202.36
52049	FUNNEL, SEPARATORY 500ML	\$	533.83
5207418	MANUAL, SERIAL I/O MODULE	\$	139.14
5208300	FUSE, 1.6 A 250V 5MM X 20MM	\$	19.85
5210000	CAP, CONNECTOR RECEPTACLE	\$	21.58
5215600	KIT, AT CABLE TERMINATION	\$	217.77
5217900	HEAD ASSY, 1720C/D TURB, REPL.	\$	697.65
5218000	PHOTOCELL ASSY, 1720D REPL. KIT	\$	1,509.00
5221500	ICE-PIC 1NTU CVD	\$	2,447.54
5222500	ICE-PIC .5NTU CVD	\$	2,447.54
5225000	ICE-PIC 20NTU CVD	\$	2,447.54
5227900	CONN, RCPT 5 CKT CIRCULAR PNL	\$	358.46
5230089	oo INST SHT, FT 660 S/W UPGRA	\$	156.30
5230689	INSTRUCTION SHT, LASER MODULE	\$	96.94
5233500	SCREW, M2(MM)X5(MM) CAP SKT	\$	0.00
5233900	REPLACEMENT BOTTOM PLATE - TURB BODY	\$	147.62
5236400	KTO: FILTERTRAK 660 CALIBRATION KIT	\$	2,081.38
5236489	INSTRUCTION SHT, CAL KIT, FT660	\$	111.30
5237100	FUNNEL, LONG STEM 60 DEG	\$	102.14

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

5237300	DUST COVER, CAL KIT FILTERTRAK 660	\$	196.57
5238600	si LIGHT TRAP, FT660	\$	33.73
5238800	MAGNET, BAR .250L X .250DIA	\$	22.93
5240018	MANUAL, DIGITAL DISPLAY MODULE	\$	166.70
52561	DISH, EVAPORATING PORC 120ML	\$	65.52
52565	DISH, EVAPORATING PORC 250ML	\$	91.35
52590	DISH, EVAPORATING PORC 70ML	\$	48.95
5260093	MANUAL, POCKET TURB, SPANISH	\$	64.18
5260094	MANUAL, POCKET TURB, FRENCH	\$	73.23
5262500	si SAMPLE CELL, POCKET TURB	\$	45.67
5262600	CAP, SAMPLE CELL	\$	22.93
5263100	ASSY, CELL W/CAP POCKET TURB PK/6	\$	406.64
5263200	si CAP, INSTRUMENT- POCKET TURB	\$	140.11
52652	FILTER PAPER, 5.5CM DIA PK/100	\$	84.03
52655	FILTER PAPER, 9CM DIA PK/100	\$	116.98
52657	FILTER PAPER, ASHLESS 12.5CM PK/100	\$	157.64
52658	FILTER PAPER, 15CM DIA PK/100	\$	196.57
5267500	1720D DRAIN VALVE ASSEMBLY	\$	753.54
527224	FUSE, 10 A 250V SLOW-BLOW	\$	13.10
527612	WASHER #6 FLAT SST	\$	2.70
5287000	LINE SURGE PROTECTOR	\$	470.24
5291100	SCREW, ENCLOSURE CONNECTING	\$	70.54
5309	FUSE, RESET, TS250, SMT	\$	16.57
53236	ff PIPET, SEROLOGICAL 2ML	\$	36.04
53237	PIPET, SEROLOGICAL 5ML	\$	36.04
53238	PIPET, SEROLOGICAL 10ML	\$	39.31
5330000	si CHAMBER, MIXING 2.0ML VOL, APA6000	\$	643.68
5410000	SP510, HARDNESS MONITOR	\$	9,578.18
5410001	KTO: SP510 HDNS MON W/RGTS 1MG/L	\$	9,578.18
5410001-LM	KTO: SP510 HDNS MON W/RGTS 1MG/L	\$	9,578.18
5410002	KTO: SP510 HDNS MON W/RGTS 2MG/L	\$	9,578.18
5410002-LM	KTO: SP510 HDNS MON W/RGTS 2MG/L	\$	9,578.18
5410003	KTO: SP510 HDNS MON W/RGTS 0.3MG/L	\$	9,578.18
5410005	KTO: SP510 HDNS MON W/RGTS 5MG/L	\$	9,578.18
5410005-LM	KTO: SP510 HDNS MON W/RGTS 5MG/L	\$	9,578.18
5410010	KTO: SP510 HDNS MON W/RGTS 10MG/L	\$	9,578.18
5410010-LM	KTO: SP510 HDNS MON W/RGTS 10MG/L	\$	9,578.18
5410020	KTO: SP510 HDNS MON W/RGTS 20MG/L	\$	9,578.18
5410050	KTO: SP510 HDNS MON W/RGTS 50MG/L	\$	9,578.18
5410050-LM	KTO: SP510 HDNS MON W/RGTS 50MG/L	\$	9,578.18
5410099	KTO: SP510 HDNS MON W/RGTS 100MG/L	\$	9,578.18
5410099-LM	KTO: SP510 HDNS MON W/RGTS 100MG/L	\$	9,578.18
5410100	SCREW 6-32X .750KNRL SST THMB	\$	40.86
5410500	CAM ASSY, SP510 PUMP	\$	1,088.87
5410800	TBG PVC 11/16"OD X 1/2"ID	\$	60.90
5410800A	TBG PVC 11/16"OD X 1/2"ID, NOT FDA CONTROLLED	\$	7.52
5411100	si COIL, PREWOUND CKT BD MT	\$	314.13

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

5411300	IC CONV LIGHT TO FREQUENCY 235	\$	95.01
5411800	si PRESSURE PLATE	\$	93.66
5412100	TBG PERISTALTIC .173 OD .045 ID	\$	74.97
5412700	si HOUSING, LOWER LED	\$	125.27
5412800	si HOUSING, UPPER LED	\$	125.27
5412900	STIR BAR, TEFLON W PIVOT RING 8X3MM	\$	60.13
5414700	KIT, 7016 PUMP, BARB REDUCED	\$	1,744.12
5415700	SWITCH, SPST PUSHBUTTON	\$	12.53
5416000	CLIP, RETAINER, DOOR HINGE	\$	20.43
5416300	si FITTING, INSTRUMENT DRAIN	\$	300.64
5417100	TUBING HARNESS, SAMPLE LINE	\$	260.17
5417200	TUBING HARNESS, REAGENT LINE	\$	356.53
5417300	TUBING HARNESS, BLANK	\$	38.16
5417400	si PIPE, DRAIN	\$	112.16
5417500	COUPLING, 1" SCH 40 PVC PIPE	\$	21.78
5417600	COUPLING, 1/2 FPT X 1/2 FPT PVC	\$	19.46
5417700	VALVE, BALL, PVC 1/2 NPT	\$	277.52
5417800	FTG, TUBE STEM ADPT 1/2 OD 1/2 NPT	\$	55.89
5417900	FTG, TUBE 1/2 OD UNION ELBOW	\$	55.50
5418000	FTG, TUBE STEM ADPT 1/2OD 1/4NPT	\$	58.97
5418100	FTG, TUBE 1/2 OD UNION STRAIGHT	\$	51.46
5418300	STRAINER, Y-BODY, 40 MESH SCREEN	\$	383.51
5418400	SCREEN, 40 MESH REPLACEMENT SST	\$	250.54
5430018	MANUAL, 1720D/L 870NM LOW RANGE	\$	233.19
5432300	FILTER, INTERFERENCE 600NM	\$	693.79
5440400	si CKT BD ASSY, MAIN BD, CL17	\$	2,582.45
5444300	KIT, MAINT, FOR CL17s SOLD AFTR 3/00	\$	634.05
5444301	KIT, CL17 PRE-ASSY MAINT, POST 3/00	\$	780.52
5445200	CAM ASSEMBLY, CL17 PUMP	\$	684.16
5445300	LABEL,CUSTOMER ACCESS COVER,CL17	\$	45.10
5448000	KIT, TUBING, NEW CL17	\$	283.30
5448800	KIT, POWER CORD W/STRAIN RELIEF	\$	181.93
5448900	KIT, PWR CORD W/STRAIN RELIEF EUROPE	\$	196.57
5449000	KTO: CL-17 CAL/VERIFICATION KIT	\$	759.32
5449089	INSTRUCTION SHT, CL-17 CAL/VER KIT	\$	36.23
5449800	KIT, VALVE, CL-17 CAL/VERIFICATION	\$	143.77
5449900	KIT, BOTTLE CAP/PRESSURIZATION CL-17	\$	281.37
54649	FLASK, FILTERING 500ML	\$	115.05
54653	FLASK, FILTERING 1000ML	\$	190.21
5467100	PROBE HOLDER, SENSION 156	\$	352.68
54679	FLASK, FILTERING 500ML PK/6	\$	703.43
54683	FLASK, FILTERING 1000ML PK/6	\$	1,032.98
5483800	WEIGHT, REAGENT, SST, 1/8 NPT	\$	136.23
5500.AMC.1.EU	5500sc AMMONIA MONOCHLORAMINE ANALYZER, AC, 1 CH, EU	\$	75,160.80
5500.AMC.1.KTO	KTO: 5500sc Ammonia Monochloramine, 1 Channel	\$	74,813.90
5500.AMC.2.KTO	KTO: 5500sc Ammonia Monochloramine, 2 Channel	\$	79,130.83



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

5500.AMC.3.EU	5500sc AMMONIA MONOCHLORAMINE ANALYZER, AC, 1 CH, EU, W/EXT FILT	\$	70,188.62
5500.AMC.3.KTO	KTO 5500sc Ammonia Monochloramine, 1 Chan, Ex Filt	\$	75,526.97
5500.AMC.4.EU	5500sc AMMONIA MONOCHLORAMINE ANALYZR, AC, 2 CH, EU, W/EXT FILT	\$	77,839.61
5500.AMC.4.KTO	KTO 5500sc Ammonia Monochloramine, 2 Chan, Ex Filt	\$	80,537.69
5500.AMC.AC.1.US	5500sc AMMONIA MONOCHLORAMINE ANALYZER, AC, 1 CH, US (9560000 FOR DWG)	\$	74,813.90
5500.AMC.AC.2.US	5500sc AMMONIA MONOCHLORAMINE ANALYZER, AC, 2 CH, US (9560000 FOR DWG)	\$	79,130.83
5500.AMC.AC.3.US	5500sc AMMONIA MONOCHLORAMINE ANALYZR, AC, 1 CH, US,EXT FILT (9560000 FOR DWG)	\$	75,526.97
5500.AMC.AC.4.US	5500sc AMMONIA MONOCHLORAMINE ANALYZR, AC, 2 CH, US,EXT FILT (9560000 FOR DWG)	\$	80,537.69
5500.KS0.A1U-LM	KTO: 5500sc SILICA, 1CH/AC/US, CHEM	\$	56,775.31
5500.KS0.A2U-LM	KTO: 5500sc SILICA, 2CH/AC/US, CHEM	\$	61,535.50
5500.KS0.A6U-LM	KTO: 5500sc SILICA, 6CH/AC/US, CHEM	\$	80,576.23
5500.KTO.PH.A1E	KTO: 5500sc PHOS HR, 1CH/AC/EU, CHEM	\$	56,775.31
5500.KTO.PH.A1U	KTO: 5500sc PHOS HR, 1CH/AC/US, CHEM	\$	56,775.31
5500.KTO.PH.A2E	KTO: 5500sc PHOS HR, 2CH/AC/EU, CHEM	\$	61,535.50
5500.KTO.PH.A2U	KTO: 5500sc PHOS HR, 2CH/AC/US, CHEM	\$	61,535.50
5500.KTO.PH.A4E	KTO: 5500sc PHOS HR, 4CH/AC/EU, CHEM	\$	71,075.14
5500.KTO.PH.A4U	KTO: 5500sc PHOS HR, 4CH/AC/US, CHEM	\$	71,075.14
5500.KTO.PH.D1U	KTO: 5500sc PHOS HR, 1CH/DC/US, CHEM	\$	56,775.31
5500.KTO.PL.A1E	KTO: 5500sc PHOS LR, 1CH/AC/EU, CHEM	\$	56,775.31
5500.KTO.PL.A1U	KTO: 5500sc PHOS LR, 1CH/AC/US, CHEM	\$	56,775.31
5500.KTO.PL.A2E	KTO: 5500sc PHOS LR, 2CH/AC/EU, CHEM	\$	61,535.50
5500.KTO.PL.A2U	KTO: 5500sc PHOS LR, 2CH/AC/US, CHEM	\$	61,535.50
5500.KTO.PL.A4U	KTO: 5500sc PHOS LR, 4CH/AC/US, CHEM	\$	71,075.14
5500.KTO.S0.A1E	KTO: 5500sc SILICA, 1CH/AC/EU, CHEM	\$	56,775.31
5500.KTO.S0.A1U	KTO: 5500sc SILICA, 1CH/AC/US, CHEM	\$	56,775.31
5500.KTO.S0.A2E	KTO: 5500sc SILICA, 2CH/AC/EU, CHEM	\$	61,535.50
5500.KTO.S0.A2U	KTO: 5500sc SILICA, 2CH/AC/US, CHEM	\$	61,535.50
5500.KTO.S0.A4E	KTO: 5500sc SILICA, 4CH/AC/EU, CHEM	\$	71,075.14
5500.KTO.S0.A4U	KTO: 5500sc SILICA, 4CH/AC/US, CHEM	\$	71,075.14
5500.KTO.S0.A6E	KTO: 5500sc SILICA, 6CH/AC/EU, CHEM	\$	80,576.23
5500.KTO.S0.A6U	KTO: 5500sc SILICA, 6CH/AC/US, CHEM	\$	80,576.23
5500.KTO.S0.D1U	KTO: 5500sc SILICA, 1CH/DC/US, CHEM	\$	56,775.31
5500.KTO.S0.D2U	KTO: 5500sc SILICA, 2CH/DC/US, CHEM	\$	61,535.50
5500.KTO.S0.D4U	KTO: 5500sc SILICA, 4CH/DC/US, CHEM	\$	71,075.14
5500.KTO.S0.D6U	KTO: 5500sc SILICA, 6CH/DC/US, CHEM	\$	80,576.23
5500.PH.A1E	5500sc PHOSPHATE HR ANALYZER, AC, 1CH, EU	\$	58,239.98
5500.PH.A1U	5500sc PHOSPHATE HR ANALYZER, AC, 1CH, US	\$	58,239.98
5500.PH.A2E	5500sc PHOSPHATE HR ANALYZER, AC, 2CH, EU	\$	63,000.17
5500.PH.A2U	5500sc PHOSPHATE HR ANALYZER, AC, 2CH, US	\$	63,000.17
5500.PH.A4E	5500sc PHOSPHATE HR ANALYZER, AC, 4CH, EU	\$	72,616.90
5500.PH.A4U	ee 5500sc PHOSPHATE HR ANALYZER, AC, 4CH, US	\$	72,616.90



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

5500.PH.D1U	5500sc PHOSPHATE HR ANALYZER, DC, 1CH, US	\$	58,239.98
5500.PH.D2U	5500sc PHOSPHATE HR ANALYZER, DC, 2CH, US	\$	63,000.17
5500.PH.D4U	5500sc PHOSPHATE HR ANALYZER, DC, 4CH, US	\$	72,616.90
5500.PL.A1E	5500sc PHOSPHATE LR ANALYZER, AC, 1CH, EU	\$	58,239.98
5500.PL.A1U	5500sc PHOSPHATE LR ANALYZER, AC, 1C, US	\$	58,239.98
5500.PL.A2E	5500sc PHOSPHATE LR ANALYZER, AC, 2CH, EU	\$	63,000.17
5500.PL.A2U	5500sc PHOSPHATE LR ANALYZER, AC, 2CH, US	\$	63,000.17
5500.PL.A4E	5500sc PHOSPHATE LR ANALYZER, AC, 4CH, EU	\$	72,616.90
5500.PL.A4U	5500sc PHOSPHATE LR ANALYZER, AC, 4CH, US	\$	72,616.90
5500.PL.D1U	5500sc PHOSPHATE LR ANALYZER, DC, 1CH, US	\$	58,239.98
5500.PL.D2E	5500sc PHOSPHATE LR ANALYZER, DC, 2CH, EU	\$	58,201.44
5500.PL.D2U	5500sc PHOSPHATE LR ANALYZER, DC, 2CH, US	\$	63,000.17
5500.PL.D4U	5500sc PHOSPHATE LR ANALYZER, DC, 4CH, US	\$	72,616.90
5500.SO.A1E	5500sc SILICA ANALYZER, AC, 1CH, EU	\$	58,239.98
5500.SO.A1U	5500sc SILICA ANALYZER, AC, 1 CH, US	\$	58,239.98
5500.SO.A2E	5500sc SILICA ANALYZER, AC, 2CH, EU	\$	63,000.17
5500.SO.A2U	5500sc SILICA ANALYZER, AC, 2CH, US	\$	63,000.17
5500.SO.A4E	5500sc SILICA ANALYZER, AC, 4CH, EU	\$	72,616.90
5500.SO.A4U	ee 5500sc SILICA ANALYZER, AC, 4CH, US	\$	72,616.90
5500.SO.A6E	5500sc SILICA ANALYZER, AC, 6CH, EU	\$	82,175.81
5500.SO.A6U	5500sc SILICA ANALYZER, AC, 6CH, US	\$	82,175.81
5500.SO.D1E	5500sc Silica Analyzer, DC/1 ch/EU	\$	53,768.88
5500.SO.D1U	5500sc SILICA ANALYZER, DC, 1CH, US	\$	58,239.98
5500.SO.D2E	5500sc Silica Analyzer, DC/2 ch/EU	\$	58,124.35
5500.SO.D2U	5500sc SILICA ANALYZER, DC, 2CH, US	\$	63,000.17
5500.SO.D4E	5500sc Silica Analyzer, DC/4 ch/EU	\$	66,758.21
5500.SO.D4U	5500sc SILICA ANALYZER, DC, 4CH, US	\$	72,616.90
5500.SO.D6E	5500sc SILICA ANALYZER, DC, 6CH, EU	\$	79,323.55
5500.SO.D6U	5500sc SILICA ANALYZER, DC, 6CH, US	\$	82,175.81
5500660	oo APA6000,AMM/MONOCHLORAMINE W/AT INST	\$	79,882.44
5500F1008-003	FAB-D.O. GUARD,SERVICE	\$	178.84
5500F1018-003	D.O. CAP,SERVICE	\$	847.97
5500KTOS0A4U-LM	KTO: 5500sc SILICA, 4CH/AC/US, CHEM	\$	71,075.14
55085	ff FUNNEL, BUCHNER PORC 56MM PLATE	\$	148.78
55087	FUNNEL, BUCHNER PORC 100MM PLATE	\$	246.68
55095	ff FUNNEL, BUCHNER PORC 126MM PLATE	\$	528.05
5516000	si CKT BD ASSY, SP510	\$	2,331.91
5516200	COLORIMETER FLOWCELL ASSY, SP510	\$	1,462.74
5516300	PLUG, FLUSH SNAP-IN BLACK	\$	4.43
5516400	KIT, INSTALLATION SP510	\$	1,719.06
5516402	KIT, INSTALLATION CL17	\$	1,927.20
5516500	KIT, MAINTENANCE SP510	\$	892.29
5516501	KIT, MAINTENANCE SP510 SPECIAL	\$	763.17
5516700	FUSE, 1.25 A 250V 5MM X 20MM	\$	25.63
5518500	WINDOW, CL17, UPPER DOOR	\$	73.81
5520000	GEARMOTOR, 1/2 RPM 115/230VAC	\$	757.39
5522300	DIAPHRAM, GRAB/BP REG. APA 6K	\$	58.01

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

5524600	IC FLASH MEM PROG, APA-CONT #2	\$	462.53
5524700	IC PROG KIT, APA 6000 FLASH MEM	\$	240.90
5526200	si CKT BD ASSY, MAIN BOARD CA610	\$	2,582.45
5526400	si CKT BD ASSY, POWER SUPPLY BD, CA610	\$	1,428.06
5527700	BODY, ISE ELECTRODE ASSY	\$	799.79
5528100	PH ELECTRODE ASSY	\$	819.06
5528400	GASKET, BNC CONN, APA	\$	4.43
55356	FILTER PAPER, 11CM DIA PK/100	\$	149.17
55358	FILTER PAPER, 15CM DIA PK/100	\$	213.92
55359	FILTER PAPER, 18.5CM DIA PK/100	\$	314.13
55480	FILTER, 480 nm 1/2"	\$	799.79
55810	FILTER, 810 NM 1/2"	\$	618.63
558111	SCREW 4-40X .25 PAN SST PHIL	\$	10.60
558311	SCREW 6-32X .25 PAN SST PHIL	\$	7.71
558411	SCREW 6-32X .312PAN SST PHIL	\$	11.95
558511	SCREW 6-32X .375PAN SST PHIL	\$	8.48
558611	SCREW 6-32X .438PAN SST PHIL	\$	8.67
558711	SCREW 6-32X .5 PAN SST PHIL	\$	11.95
55880	FILTER, 880 nm 1/2"	\$	585.87
558811	SCREW 6-32X .625PAN SST PHIL	\$	8.67
558911	SCREW 6-32X .75 PAN SST PHIL	\$	5.78
559111	SCREW 6-32X1. PAN SST PHIL	\$	3.47
559211	SCREW 6-32X1.25 PAN SST PHIL	\$	2.51
559611	SCREW 8-32X .25 PAN SST PHIL	\$	8.67
559711	SCREW 8-32X .375PAN SST PHIL	\$	8.48
559811	SCREW 8-32X .5 PAN SST PHIL	\$	8.48
560150,21957	SERINGUE 5ML	\$	77.09
56018	TUBING, RBR LATEX 6.4X2.4MM 3.6M	\$	71.50
56019	TUBING, RBR LATEX 7.9X2.4MM 3.6M	\$	72.85
560211	SCREW 8-32X1.25 PAN SST PHIL	\$	8.67
561011	SCREW 10-32X .375PAN SST PHIL	\$	3.47
561211	SCREW 10-32X .5 PAN SST PHIL	\$	1.93
561311	SCREW 10-32X .562PAN SST PHIL	\$	9.06
56162	SPATULA, 7.6CM STAINLESS STEEL	\$	37.77
561711	SCREW 10-32X1. PAN SST PHIL	\$	2.51
562600	SCREW 6-32X .25 FLAT SST PHIL	\$	1.93
56275	STOPCOCK GREASE, 150G TUBE	\$	151.09
56300	SUPPORT, RING STAND 5X8 BASE	\$	70.34
563011	SCREW 6-32X .5 FLAT SST PHIL	\$	12.14
563111	SCREW 6-32X .625FLAT SST PHIL	\$	11.37
5646300	INSERTION MTG, DIGITAL, 1 1/2" CPVC	\$	2,312.64
5646350	INSERTION MTG, DIGITAL, 1 1/2" SST	\$	8,691.67
5646400	INSERTION MTG, ANALOG, 1 1/2" CPVC	\$	2,312.64
5646450	INSERTION MTG, ANALOG, 1 1/2" SST	\$	8,691.67
5646700	Probe Adapter Assembly, Insertion, 1 1/2" CPVC	\$	1,717.14
56510	TUBE, TEST 10ML 13X100MM PK/10	\$	37.77
56528	ff TUBE, TEST 28ML 18X150MM PK/10	\$	78.05

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

5657200	si FITTING, TEFLON, MANIFOLD, TOC	\$	271.74
56800	TONGS, BEAKER	\$	35.65
56900	TONGS, CRUCIBLE	\$	26.60
5694800	si CKT BD ASSY, SP510, SPECIAL	\$	2,620.99
5700700	KTO: FLOW CONTROLLER, WEIR	\$	1,219.92
5700800	KTO: SCREEN/STRAINER, SM VOLUME, 40-MESH	\$	1,085.01
5701525	VISTA SOFTWARE VERSION 1.3.1	\$	9,732.36
5702500	CD KIT, SOFTWARE, 2200 PCX OPC	\$	5,550.34
57100	IMHOFF CONE, GLASS	\$	799.79
5710600	ADAPTER, RS-232 TO RS-485 PORT PWRED	\$	946.26
571211	SCREW 8-32X .625PAN SST PHIL	\$	8.48
5715000	PUMP MODULE ASSY, APA 6000	\$	11,852.28
5715100	si RETAINER, DOUBLE SEAL, APA PUMP	\$	225.48
5715200	si HEAD, APA PUMP MODULE	\$	637.90
5716100	ALIGNMENT TOOL, APA PUMP MODULE	\$	136.45
5716300	SEAL, FLANGED .722 O.D. X .104 LG	\$	97.13
5716400	PISTON, CERAMIC	\$	1,364.46
5716401	si PISTON, MODIFIED, APA 6000	\$	2,023.56
57200	IMHOFF CONE SUPPORT, HOLDS 3	\$	208.14
5733000	KTO: CALIBRATION KIT, OPTIQUANT SST	\$	800.39
5733000-LM	KTO: CALIBRATION KIT, OPTIQUANT SST	\$	800.39
5735300	BRACKET, OPTIQUANT CAL. IN-LINE	\$	472.16
5737600	si CYLINDER ASSY, GENERIC CALIBRATION	\$	555.03
5738400	KTO: INSERTION MOUNTING KITS FOR INLINE &	\$	10,368.34
5738900	KTO: FILTRAX SAMPLE FILTRATION SYS, 2M UNHEATED HOSE 115V	\$	26,980.80
5738901	KTO: FILTRAX SAMPLE FILTRATION SYS, 2M HEATED HOSE 230V	\$	26,980.80
5739000	KTO: FILTRAX SAMPLE FILTRATION SYS, 10M HEATED HOSE 115V	\$	28,445.47
5739001	KTO: FILTRAX SAMPLE FILTRATION SYS, 10M HEATED HOSE 230V	\$	28,445.47
5739100	KTO: FILTRAX SAMPLE FILTRATION SYS, 20M HEATED HOSE 115V	\$	30,661.75
5739101	KTO: FILTRAX SAMPLE FILTRATION SYS, 20M HEATED HOSE 230V	\$	30,661.75
5739200	KTO: FILTRAX SAMPLE FILTRATION SYS, 2M HEATED HOSE 115V	\$	28,869.46
5740000	CA610 FLUORIDE ANALYZR W/KITS	\$	25,670.30
5740001	KTO: CA610 FLUORIDE ANALYZER W/REAGENTS	\$	27,770.95
5740011	KTO: CA610 FLUORIDE W/REAGENTS EU	\$	26,075.02
5740060X	CA610 FLUORIDE ANALYZER INST ONLY	\$	24,648.89
5740300	FLOW CELL, CA610 FLUORIDE ANALYZER	\$	1,561.03
5740900	KEYBOARD, CA610 FLUORIDE ANALYZER	\$	111.01
5741500	CAM ASSEMBLY, CA610 PUMP	\$	697.65
5741800	si BLOCK, REAGENT,CA610 FLUORIDE ANALYZER	\$	122.96
5742000	KIT, INSTALLATION CA610 FLUORIDE	\$	319.92
5742100	KIT, MAINTENANCE CA610 FLUORIDE	\$	639.83
5742600	TUBING PERISTALTIC .158OD .028ID	\$	47.79
5742700	INSTRUMENT TUBING SET, CA610	\$	356.53
5743100	si BUSHING, MOTOR ALIGNMENT	\$	29.10
5743200	STAND ASSY, VERSA STAND	\$	4,663.82
5743500	CLAMP, VERSA STAND	\$	51.07



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

5743700	VALVE, SHUT-OFF, 2-WAY F/ 1/4" TBG	\$	116.98
5743700A	VALVE, SHUT-OFF, 2-WAY F/ 1/4" TBG, NOT FDA CONTROLLED	\$	65.33
5744400	PUMP TUBING SET, CA610	\$	189.83
5744800	KIT, FLUORIDE ELECTRODE, CA610	\$	3,680.95
5745100	FLUORIDE TIP ASSY. CA610	\$	2,967.89
5773000	KTO: SONATAX ^{sc} PROBE & PIVOT MOUNTING	\$	12,218.45
5773100	ff KTO: SONATAX ^{sc} 2 PROBES & 2 PIVOT	\$	24,417.62
5773400	ff KTO: SONATAX ^{sc} 4 PROBES & 4 PIVOTS	\$	46,773.14
578=401=501	STOPPER (RED) FOR QUICK FITTING, 6MM	\$	7.13
578=507=602	TOOL SET, 9240 ANALYZER QUICK FITTINGS	\$	181.35
578=507=603	Complete drain coupler with tubing kit for 924x panel (Z_0924x=A=x0xx)	\$	400.86
578=507=604	Complete drain coupler with tubing kit for 924x cabinet (Z_0924x=A=x1xx)	\$	686.08
578=601=703	FITTING, 6MM EXT. DIA TUBING	\$	17.73
578=602=703	QUICK FITTING FOR 8MM O.D. TUBING	\$	12.33
57870	WATCH GLASS, PYREX 100MM	\$	25.63
5791700	si INNER FLOAT PIPE, DO PROBE MTG	\$	213.92
5792000	FLOAT, MARINE, RED, INFLATABLE	\$	529.98
5792001	D.O. Float with Pipe Installed	\$	703.43
5792002	FLOAT WITH PIPE 6136600 INSTALLED	\$	1,522.49
5794300	KIT, MTG HRDWARE,DO PROBE BALL FLOAT	\$	2,389.73
5794389	INSTRUCTION SHT, MTG HDWR, DO PROBE	\$	31.22
5794400	KIT, MTG HRDWARE, DO PROBE RAIL MT	\$	1,746.04
5795450	KIT, HEAVY DUTY SWIVEL MOUNT, LDO	\$	1,667.03
5796000	CABLE ASSY, 7.7M PROBE EXTENSION	\$	639.83
5796001	aa CABLE ASSY, 7.7M PROBE EXTENSION	\$	1,753.75
5796002-225	CABLE assy, 225 Ft / 68.5 MTR PROBE EXTENSION	\$	2,524.63
5796100	CABLE ASSY, 15M PROBE EXTENSION	\$	969.38
5796101	CABLE ASSY, 15M PROBE EXTENSION	\$	3,102.79
5796200	CABLE ASSY, 31M PROBE EXTENSION	\$	1,298.93
5796201	CABLE ASSY, 31M PROBE EXTENSION	\$	3,218.42
5796500	HEAD ASSY, AIR BLAST CLEANING SYSTEM	\$	495.29
5796600	BAG, METALIZED MYLAR W/ BLK 6X16	\$	20.24
5796605	BAG, CALIBRATION, LDO, SET OF 5	\$	63.98
5797400	HOSE ASSY, AIR BLAST CLEANING SYSTEM	\$	738.12
5797700	GASKET, DO PROBE BODY	\$	40.66
5799200	FLOW CELL ASSY FOR GLI PHD PROBE	\$	585.87
58000	RING, SUPPORT 3 INCH	\$	43.94
58001	RING, SUPPORT 4 INCH	\$	50.69
580186	mc Fitting,Q-D,FERRULESS,THRU	\$	156.68
580222	mc Fitting, Q-DISC F Panel M	\$	136.83
580223	mc Fitting Q-DISC M FERRULES	\$	131.24
5806000	si CKT BD ASSY, CONDUCTIVITY	\$	194.65
5806400	si CKT BD ASSY, CONNECTOR, MP	\$	776.66
580800	HACH ENGINEERED SYSTEMS	\$	-
5810700	si SHROUD, HQ SERIES	\$	45.29

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

5811200	SENSOR CAP REPLACEMNT,HQd SERIES LDO	\$	450.96
5811200DKQ	SENSOR CAP REPLACEMNT,HQd SERIES LDO	\$	474.09
5811200PC	LDO CALIBRATION CAP, HQd SERIES LDO	\$	35.65
5811600	SHROUD, LOCKRING	\$	45.29
5811900	BOD FUNNEL, US	\$	16.96
5812500	si VALVE, BODY, POWDER DISPENSER	\$	16.38
5812711	AMMONIA INTELICAL MEMBRANE MODULES, 3-PK	\$	431.69
5813400	aa ADAPTER, USB/DC POWER JACK	\$	846.04
5816900	DISPLAY 240 x 160 LCD	\$	659.10
5817400	si CKT BD ASSY, CONNECTOR, LP	\$	637.90
5818101	si CABLE ASSY, CDC & PH STD PROBE, 1M	\$	63.98
5818200	LED, BLUE, 430NM, 85NM BW	\$	47.22
5818400	KIT, PROBE CLIP	\$	55.31
5818800	si SENSOR BARREL, CONDUCTIVITY PROBE	\$	16.96
5818900	si HOUSING, TOP, HQ PORTABLE METER	\$	26.02
5819000	si HOUSING, BOTTOM, HQ PORTABLE METER	\$	33.15
5819200	si COVER, BATTERY, HQ PORTABLE METER	\$	23.90
5819300	SEAL, PERIMETER, PORTABLE METER	\$	70.54
5819400	LENS, WINDOW, HQ40d	\$	87.49
5819800	CONTACT, BATTERY POST	\$	11.37
5824600	CONN, RCPT 5 CKT CIRCULAR PCB	\$	54.93
5824900	CABLE, FLAT FLEX 1.50 x 10.5mm	\$	13.49
5825700	BOX 11-3/8"x 9-1/16"x 2", HQ SERIES	\$	13.10
5825800	ff FIELD CASE, STANDARD, HQd	\$	556.96
5825900	KTO: SHROUD KIT,RGD PROBE, HQD SERIES	\$	92.89
5826300	USB/AC POWER ADAPTER KIT, 115V P-CORD	\$	422.06
5826600	KTO: BOD FUNNEL, US	\$	449.04
5826700	KTO: BOD FUNNEL, EU	\$	329.55
5827000	POWER SUPPLY EXT, 100-240VAC IN 7.5V OUT	\$	185.20
5827600	si BODY, SENSOR, PH, GEL	\$	9.25
5827700	CAP, REFILL	\$	5.78
5828600	MARKER, CABLE, RUGGED PROBE	\$	37.58
5828610	KTO: MARKER, CABLE, pk/5, RUG PROBE	\$	150.51
5828700	ff ASSEMBLY, GLOVE KIT COMPLETE	\$	184.05
5829400	PROBE HOLDER, HQ SERIES	\$	111.58
5829800	LID, PH CONDITIONING CHAMBER	\$	17.34
5830200	LENS, WINDOW, HQ11d	\$	75.55
5830400	LENS, WINDOW, HQ30d	\$	77.67
5830500	EMBLEM, HQ SERIES	\$	8.09
5831100	POWER ADPTR KIT, 230V EU P-CORD	\$	269.81
5832400	si SHROUD, HQ SERIES PROBE	\$	213.92
5832500	KTO:SHROUD KIT,STD PROBE, HQD SERIES	\$	217.77
5833700	ADAPTER KIT, USB/DC, 115V, US	\$	576.23
5834100	ADAPTER KIT, USB/DC, 230V, GERMAN	\$	576.23
5834400	ADAPTER KIT, USB/DC, 230V, SPANISH	\$	576.23
5835600	ff ASSY, GLOVE, WRIST STRAP, WRIST CORD	\$	79.79
5836000	PAPER, THERMAL PD24 4-3/8"X105'PK/5	\$	174.99

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

5838000	SENSOR CAP REPLACEMENT, LBOD	\$	450.96
5838300	si CLAMP BAR	\$	84.03
5838500	si CASE ASSY, HC-8, DREL	\$	1,591.87
5838700	PUNCH PLATE, DREL	\$	168.63
5839100	si CASE ASSY, DREL COMP. APP	\$	1,341.33
5839300	PLATE, PUNCHED, COMPL. & INDUSTRIAL	\$	138.57
5839600	si SPACER, CASE ASSY, HC-7	\$	55.12
5839700	si Carrying Case (HC-7), DR2800	\$	1,034.91
5839800	INSERT, FOAM	\$	470.24
5839900	PLATE, PUNCHED	\$	163.62
58451-0	Sodium Persulphate 0.2 M	\$	1,306.28
5845100	SODIUM PERSULFATE SOLN, 0.2M 5 GAL	\$	775.10
58451-00	SODIUM PERSULFATE SOLN, 0.2M 5	\$	630.56
5845100-LM	SODIUM PERSULFATE SOLN, 0.2M 5 GAL	\$	775.10
5845200	SODIUM PERSULFATE SOLN, 0.4M 5 GAL	\$	776.90
58452-00	Sodium Persulfate Soln. 0.4M,	\$	637.78
5845200-LM	SODIUM PERSULFATE SOLN, 0.4M 5 GAL	\$	776.90
5845300	SODIUM PERSULFATE SOLN, 0.6M 5 GAL	\$	775.10
58453-00	Sodium Persulfate Soln. 0.6M,	\$	637.78
5845300-LM	SODIUM PERSULFATE SOLN, 0.6M 5 GAL	\$	775.10
5845400	SODIUM PERSULFATE SOLN, 0.8M 5 GAL	\$	872.66
58454-00	Sodium Persulfate Soln. 0.8M,	\$	701.02
5845500	SODIUM PERSULFATE SOLN, 1.0M 5 GAL	\$	1,033.46
5845500-LM	SODIUM PERSULFATE SOLN, 1.0M 5 GAL	\$	1,033.46
5845600	SODIUM PERSULFATE SOLN, 1.2M 5 GAL	\$	1,262.92
58456-00	SODIUM PERSULFATE SOLN, 1.2M 5	\$	1,006.36
5845700	SODIUM PERSULFATE SOLN, 1.5M 5 GAL	\$	1,353.26
58457-00	SODIUM PERSULFATE SOLN, 1.5M 5	\$	1,080.44
5845800	PHOSPHORIC ACID SOLN, 0.1M 5 GALLON	\$	619.72
5845800-LM	PHOSPHORIC ACID SOLN, 0.1M 5 GALLON	\$	619.72
5845900	PHOSPHORIC ACID SOLN, 0.3M 5 GALLON	\$	619.72
5846000	PHOSPHORIC ACID SOLN, 0.6M 5 GALLON	\$	715.47
5846000-LM	PHOSPHORIC ACID SOLN, 0.6M 5 GALLON	\$	715.47
5846100	PHOSPHORIC ACID SOLN, 1.0M 5 GALLON	\$	715.47
5850800	ASSEMBLY, BEARING, LBOD	\$	705.36
5850800L	ASSEMBLY, BEARING, LBOD long Version, with short DOC	\$	651.39
5852700	si STIR PADDLE, LBOD	\$	4.05
5852800	si CLIP, HOLD DOWN, BEARING, LBOD	\$	17.15
5854000	KIT, HQD SERIES FIRMWARE UPGRADE	\$	1,171.74
5857100	SPRING, RETAINING, DR300	\$	9.64
5857305	Plastic Cap for Ph and Cond. Electrodes PK 5	\$	56.27
5861400	POWER SUPPLY, 90-264VAC, 47/63 HZ	\$	1,443.47
5866600	CD KIT, SOFTWARE, OPC PROCESS LINK	\$	319.92
5867000	TERMINATION BOX ASSY	\$	711.14
5867300	NUT M6X1.0 HEX MACHINE SCREW SST	\$	8.09
5868700	PLUG, HOLE,TPE,DOME,1/2 IN.COND	\$	36.23
5869800	SCREW DRIVER T10H TORX DRIVE	\$	151.09

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

587=004=002	Adaptor 1/4" G / DN8	\$	124.50
588800,29050	CABLE, SHIELDED 4-CONDUCTOR C	\$	78.05
588801,75008	ONE CLEANING BEAD (9186)	\$	18.50
58901	GLASS WOOL, 8UM DIA FIBERS 454G	\$	389.29
590=011=130	Set of 12 tubes grey-grey black viton	\$	2,158.46
590=020=165	Blue/blue ismaprene tube, one piece	\$	173.45
590=030=076	TUBE ELEMENT MARPRENE 0,76 x 2,36	\$	131.63
590=060=026	TUBING, PTFE, 2MM ID X 6MM OD - 0.3M Length	\$	19.27
590=121=810	Tubing, Polyethylene 8 x 10 mm (per meter)	\$	23.13
590=160=080	TUBING, SILICONE 6/8 MM	\$	146.66
590900	BOTTLE, BR POLY NARROW MOUTH 6OZ	\$	18.50
5909900	STATIC MIXER	\$	39.89
59100	PIPET BOTTLE COMPLETE	\$	206.21
59101	PIPET & BULB ONLY	\$	131.44
59102	BOTTLE ONLY	\$	63.79
591200	SPRINKLER WITH SCREEN	\$	13.30
5913800	PLATE, WINDOW, ODYSSEY	\$	38.74
5917900	LENS, CUSTOM, FIBER PICKUP, ODYSSEY	\$	44.90
5920400	KIT, CONVERTER, RS485/422 TO USB	\$	3,546.05
5920500	CONVERTER, RS485/422 TO USB	\$	5,858.69
5921000	CKT BD ASSY, MODBUS RS485	\$	2,254.82
5924000	CABLE, USB,TYPE A MALE-TYPE B MALE	\$	106.77
5932800	INSERT, BOTTOM, HQ PROBE SERIES	\$	58.97
5940400	KIT, POUR-THRU CELL	\$	1,208.35
5940500	CELL, MULTI-PATHLENGTH, 25ML	\$	39.31
5940506	CELL, MULTI-PATHLENGTH, 25ML 6/PK	\$	161.11
5941000	ASSY, POUR-THRU CELL, DR/2400 PORT.	\$	227.41
5941100	ADAPTER, POUR-THRU CELL	\$	117.17
5941300	si CAP, MULTI-PATHLENGTH CELL	\$	46.06
5941400	CELL, MULTI-PATHLENGTH, 25ML	\$	40.66
5941401	si CELL, MODIFIED FOR POUR THRU ASSY	\$	118.33
5942100	EMBLEM,HACH	\$	59.55
5946600	ADAPTER, MULTI-PATHLENGTH CELL	\$	117.17
5948400	T-BUSHING, GLIDE	\$	33.73
595=000=002	IN-LINE FILTER FOR TUBING 4MM ID	\$	112.93
595=006=000	Filter(60µ) -water - High perf. (8bars - 40°C) and transparent cartridge holder (made microcrystallin PP) - Inlet&Outlet quick-connectors (6mmOD) - H(85)Diam.(56) - Dismountable - See filtration cartridge packages at 595=006=9xx	\$	813.28
595=006=005	Filter(60µ) -water, imperial - High performance (8bars - 40°C) and transparent cartridge holder (made microcrystallin PP) - Inlet&Outlet quick-connectors (1/4" OD) - H(85)Diam.(56) - Dismountable	\$	813.28
595=006=906	Filtration cartridge(60µ) pack (6 per pack)-filter 595=006=0xx	\$	813.28
595=010=000	FILTER(100u) FOR WATER - HIGH PERFORMANCE	\$	612.85
595=010=005	FILTER, (100µ) FOR WATER-HIGH PERFORMANCE	\$	612.85
595=010=906	FILTRATION CARTRIDGE(100u) PACK (6 PER PACK)	\$	612.85



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

5951000	FILTER, INTERFERENCE 420NM	\$	545.40
5951100	FILTER, INTERFERENCE 450NM	\$	250.54
5951200	FILTER, INTERFERENCE 476 NM	\$	909.64
5951300	FILTER, INTERFERENCE 500 NM	\$	487.58
5951400	FILTER, INTERFERENCE 528 NM	\$	148.97
5951500	FILTER, INTERFERENCE 550 NM	\$	262.10
5951600	FILTER, INTERFERENCE 580 NM	\$	304.50
5951700	FILTER, INTERFERENCE 600 NM	\$	170.75
5951800	FILTER, INTERFERENCE 655 NM	\$	323.77
5952100	LENS, WINDOW	\$	52.61
5952200	BATTERY CONTACT, FIXED, DUAL AAA	\$	10.41
5952300	BATTERY CONTACT, SPRING, DUAL AAA	\$	11.18
5952400	KEYPAD, PKT II	\$	38.35
5953700	INSERT, SWIFTEST MOLDED, POCKET II	\$	54.15
5954000	CUP, OPTICAL, DENSICHEK, PKT II	\$	11.76
5954100	RETAINER, DETECTOR, DENSICHEK, PKT II	\$	23.32
5954200	LED, VIOLET T-1 SIZE W/LEADS	\$	27.56
5954300	HOUSING, LED, WHITE, PKT II	\$	9.64
5954400	HOUSING, LED, VIOLET, DENSICHEK	\$	38.93
5954500	RETAINER, LED, DENSICHEK, PKT II	\$	19.46
5954600	si ff ADAPTER, 1 CM	\$	37.39
5954610	ff ADAPTER, 12-mm	\$	140.49
5954800	COVER, CELL CHAMBER, PKT II	\$	39.70
5954900	COVER, BATTERY, PKT II	\$	23.13
5955000	HOUSING, TOP, PKT II	\$	23.32
5955100	HOUSING, BOTTOM, PKT II	\$	41.63
5955400	SEAL/GRIP, PERIMETER, DENSICHEK, PKT II	\$	38.35
5955500	LIGHTPIPE, ZEBRA FRAME, DENSICHEK, PKT II	\$	8.48
5955800	CLIP, DETECTOR, DENSICHEK, PKT II	\$	16.57
5955900	LANYARD	\$	23.51
5958000	CONTACT, BATTERY POST, DENSICHEK, PKT II	\$	3.85
5961400	si MANIFOLD, MODIFIED, ACRYLIC, TOC	\$	2,042.83
5977100	LOCK COLLAR NUT, PVC, INSERT MT, GLI	\$	80.94
5977400	ANALOG CLAMP HALF, SS, INSERT MT GLI	\$	323.77
5977600	LOCKING COLLAR NUT, SS, INSERT MT GL	\$	1,146.68
5989600	SWIVEL ELBOW, PP, 1/4" TUBE X 1/4"NPT	\$	33.34
5994800	TBG, PVC, BRAIDFLEX, 1"OD X 3/4"ID	\$	9.06
599990,04300	VALVE, SOLENOID 4/2-WAY (AIR	\$	2,910.07
5H0675	O-RING, 016, VITON, 75	\$	14.98
5H0694	O-RING, 112, VITON, 75	\$	16.00
5H0946	O-RING, 210, NITRILE, 70	\$	6.17
5H1016	O-RING, 113, VITON 70	\$	59.74
5H1017	O-RING 0.987IDx 0.103W VITRON PTFE	\$	41.63
5H1018	O-RING, 0.864IDx 0.070W VINTON PTFE	\$	6.55
5H1019	O-RING, 026, NITRILE, 70	\$	6.17
5H1090	O-RING, 114, VITON, 75	\$	33.34
5H1099-007	O-RING, 007, KALRZ, 4079	\$	362.31

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

5H1099-224	O-RING, 224, KALREZ 4079	\$	1,285.44
5H1152	O-RING, 224, VITON 75, BRN	\$	40.66
5H1177	O-RING, 117, VITON, 75	\$	33.34
5H1212	O-RING, 224, E.P.R., 515-	\$	21.39
5H1219	O-RING, 215, VITON 75	\$	42.98
5H1221	O-RING, 116, VITON, 75	\$	42.98
5H1223	O-RING, 109, VITON, 75	\$	13.10
5H1226	O-RING, 216, N674, NITRILE	\$	10.65
5H1233	O-Ring Viton, 2.25 ID X 2.5 OD X .125 C/S	\$	65.14
5H1236	O-RING, LATHE CUT	\$	59.74
5H1242	O-RING, 5/8ID 7/8OD 1/8VIT	\$	33.34
5H1253	O-RING, SILICONE, 11 X 2.5	\$	33.34
5H1254	O-RING, SILICONE, 19 X 3.5	\$	33.34
5H1258	O-RING, 127, VITON, V884-	\$	146.27
5H1261	O-Ring, Viton 70, Black	\$	48.76
5H1264	O-RING, 154, BUNA N/NITRILE	\$	13.49
5H1272	O-RING, -111, VITON, 75	\$	5.01
5H1273	O-RING, 336, VITON 70	\$	129.32
5H1276	O-RING, 127, E.P.R., 540-	\$	135.87
5H1281T	O-RING, 015 VITON, V747-7, TEFLON COATED	\$	3.85
5H1285	O-RING, 334 VITON, V884-7	\$	46.64
5H1290	O-RING, E.P.R., -111 E54	\$	33.34
5H1302	O-RING, 013, VITON 75	\$	12.14
5H1303	O-RING, 218, VITON 75	\$	70.34
5H1304	O-RING, 019, VITON 75	\$	46.64
5H1315	O-RING, -111, SILICONE, DURO	\$	11.37
5H1316	O-RING, -015, SILICONE, DURO	\$	11.37
5H5002	O-RING, 223, NITRILE 70	\$	33.34
6000000	oo KTO: ANALYZER, S-5000 SILICA W/CHEM	\$	30,912.29
6000018	MANUAL, S-5000 SILICA	\$	443.26
6000094	MANUAL, S-5000 SILICA FRENCH	\$	252.46
6000118	MANUAL, S-5000 ORTHOPHOSPHATE	\$	543.47
6000218	MANUAL, S-5000 HYDRAZINE	\$	474.09
6000318	oo MANUAL, S-5000 ALUMINUM	\$	720.77
6000518	MANUAL, S-5000 LR ORTHOPHOSPHATE	\$	443.26
6000544	CARD, QUICK REF S-5000 LR PHOS	\$	97.13
6001000	KTO: KIT, RGT, APA 6000 ALKALINITY	\$	364.96
6001000-LM	KTO: KIT, RGT, APA 6000 ALKALINITY	\$	364.96
6001100	KTO: KIT, STD, APA 6000 ALKALINITY	\$	211.39
6001400	KTO: KIT, RGT, APA 6000 LR AM/MONO	\$	576.35
6001500	KTO: KIT, STD, APA 6000 LR AM/MONO	\$	187.90
6001700	KTO: KIT, STD, APA 6000 LR COPPER	\$	444.46
6001800	KTO: KIT, STD, APA 6000 HR COPPER	\$	444.46
6001900	KTO: KIT, RGT, APA 6000 LR HARDNESS	\$	410.13
6002000	KTO: KIT, STD, APA 6000 LR HARDNESS	\$	224.04
6002100	KTO: KIT, RGT, APA 6000 HR HARDNESS	\$	594.42
6002100-LM	KTO: KIT, RGT, APA 6000 HR HARDNESS	\$	594.42

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

6002200	KTO: KIT, STD, APA 6000 HR HARDNESS	\$	224.04
6010018	MANUAL, 1720E TURBIDMETER	\$	174.99
6011600	CABLE ASSY, 4 COND SHLD W/CIRCULAR	\$	277.52
6014500	ASSY, FIBER OPTIC DETECTOR, FT660sc	\$	865.31
6015089	MANUAL INSERT, 1720E S 2 ENGLISH	\$	58.01
6015900	FT600 sc BODY/BUBBLE TRAP W/PKG	\$	896.15
6016700	CABLE ASSY, DETECTOR , SS7	\$	472.16
6022P0	SNSR-pH,RYTON,5W,10'	\$	3,565.32
6022P010000200N	PH-SNSR,ENCPS,CONV,RYTON,	\$	4,702.37
6022P020	SNSR-pH,RYTON,5W,20'	\$	4,663.82
6022P050	PH-SNSR,RYTON,50'CABLE,	\$	4,066.39
6022P1	SNSR-pH,RYTON,5W,ANTIMONY	\$	4,374.74
6022P210	PH-SNSR,ENCPS,CONV,RYTON,	\$	3,680.95
6022P510	SNSR-pH,RYTON,5W, FLUSH	\$	3,931.49
6022P610	PH-ENCPS,CONV.RYTON,LOW	\$	3,680.95
6028P0	Convertible pH LCP 5 Wire	\$	3,353.33
6028P001	SNSR-PH,5-Wire,LCP,1'	\$	4,374.74
6028P00500S000N	SNSR-pH,LCP,5W,5',SP LUGS	\$	4,066.39
6028P010000010N	SENSOR-pH,LCP,5-WIRE,10'	\$	3,565.32
6028P010000200N	Sensor, 5 Wire pH, Perflouro O-Rings	\$	4,471.10
6028P01000A000N	SNSR-PH,LCP,5W,10'CABLE,	\$	3,603.86
6028P010F00000N	PH-SNSR,CONV,LCP,GEL FLLD	\$	3,565.32
6028P010W00000N	SENSOR-pH,LCP,5-WIRE,10	\$	4,432.56
6028P012	PH SNSR,LCP,CONVERTIBLE,	\$	4,374.74
6028P014	PH SNSR,LCP,CONVERTIBLE,	\$	4,394.02
6028P015	SNSR-pH,LCP,5W,15'	\$	4,394.02
6028P01500A000N	PH-SNSR,ENCPS,CONV,LCP,	\$	4,702.37
6028P020	SNSR-pH,LCP,5W,20'	\$	4,432.56
6028P020000010N	SENSOR-pH,LCP,5-WIRE,20'	\$	4,663.82
6028P025000010N	SENSOR-pH,LCP,5-WIRE,25'	\$	3,777.31
6028P025W00000N	pH Sensor, LCP, 25 Foot Teflon Cable	\$	4,875.82
6028P025W30210N	pH Sensor, LCP, 25 ft. TFE Cable	\$	6,436.85
6028P033	SNSR-pH, LCP,5W,33'	\$	4,509.65
6028P050	SNSR-pH, LCP,5W,50'	\$	3,777.31
6028P090	SNSR-pH,LCP,5W,90'	\$	5,145.62
6028P1	SNSR-pH,LCP,5W,ANTIMONY,1	\$	4,066.39
6028P140	SNSR,PH,LCP,5W,ANTIMONY	\$	5,357.62
6028P210000010N	SNSR-PH,LCP,5W,10',"S"GLS	\$	3,680.95
6028P4	SENSOR-pH,LCP,5-WIRE,10'	\$	3,565.32
6028P420	SNSR-pH,LCP,5W,20'CBL,LO	\$	4,663.82
6028P510	SNSR-pH,LCP,FLUSH ELTRD 5	\$	3,680.95
6028P6	SENSOR-pH,LCP,5-WIRE,10'	\$	3,468.96
6028P610	SNSR,PH,LCP,5W,10',STAND,	\$	3,468.96
6028P633	SNSR-pH,LCP,5W,33'LOW PH	\$	4,663.82
6028P8	SENSOR-pH,LCP,5-WIRE,10',	\$	3,468.96
6052P0	SNSR-pH,F-THRU,RYTON,10'	\$	3,989.30
6058P0	KTO: SNSR-pH,F-THRU,LCP,10'	\$	3,738.77

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

6058P010000000N	PH SNSR,ENCPS,FLW THRU UN	\$	3,738.77
6058P01000A000N	SNSR-pH,F-THRU,LCP,10'	\$	4,066.39
6058P025	PH SNSR-UNION MOUNT,LCP,	\$	3,989.30
6058P050	SNSR-PH,FLOW-THRU UNION M	\$	4,297.66
6058P1	FINISH-SNSR,pH,LCP,UNION	\$	4,490.38
6058P4	SNSR-PH,LCP,5W,HIGH PURIT	\$	3,989.30
6058P6	SNSR-pH,F-THRU,LCP,10'	\$	3,835.13
6058P610F00010N	SNSR-pH,F-T,LCP,10',"H"GL	\$	4,394.02
6058P620	SNSR-pH,F-THRU,LCP,20',H	\$	4,933.63
6058U0	SNSR-pH,F-THRU,LCP,10',UL	\$	3,892.94
606932	Potassium Iodate, ACS, 100 g	\$	162.97
607726	TRIS-HYDROXYMETHYL-AMINO-ACS	\$	154.47
608000	DROPPER, PLASTIC UNMARKED	\$	9.64
608424	DPD OXALATE, ANALYTICAL 25G	\$	285.47
608424-LM	DPD OXALATE, ANALYTICAL 25G	\$	285.47
608493	DPD OXALATE	\$	4,191.66
60-9764-000-001	KIT-SALT BRDG,RYTON,PKG O	\$	254.39
60-9764-010-001	KIT-SALT BRDG,RYTON,CR/CR	\$	224.04
60-9764-200-001	KIT-SALT BRIDGE,PERFLUORO	\$	612.85
60-9765-000-001	PKG-LCP SALT BRIDGE,QTY	\$	281.37
60-9765-009-001	KIT -SALT BRDG,HI-PURITY	\$	249.33
60-9765-010-001	PKG-LCP SLTBRDG,CER/CER,Q	\$	249.33
60A2034-001	1" X 1.5"CPVC ADAPTER	\$	408.57
60A2035-001	ADAPTER FOR MH474 PVC	\$	449.04
60A2035-002	ADAPTER FOR MH434 CPVC	\$	260.17
60A2040-001	ASY-ADAPTER, 2" CPVC IN-	\$	489.51
60A2040-002	ASY-ADAPTER,2" 316 SS IN-	\$	813.28
60A2050-001	ASY-ADAPTER,1-1/2" INSRTN	\$	680.30
60A2053	ASSY,JUNCTION BOX,SURFACE	\$	178.84
60A2100-001	ADAPTER KIT FOR EXIST.	\$	934.69
60A2A9901	ASSY-ELTRD,pH,ACTIVE,SCHO	\$	637.90
60A2F1278	si pH/ORP Sensor Protector	\$	213.92
60A2F1278-300	si PROTECTOR-MOLDED RYTON	\$	244.75
60A2F1281	si FAB.ADAPTER,UNION	\$	489.51
60A2F1299	CAP PH ELECT. EPDM BLACK	\$	29.44
60A2F1317	FAB-SALTBRIDGE BOOT EPDM	\$	29.44
60A2F1327	FAB-RETAINING RING	\$	100.21
60A2F1334	SPONGE,LCP SALT BRIDGE	\$	29.44
60A3F1061-101	316 SS. LOCK RING FOR LCP	\$	655.25
60A4908-400	ASSY-SLT BRDG,RYT/KZ,CER/	\$	276.43
60A9753-001	ASY,UNION ADAPTER,CPVC,	\$	356.53
60A9919	ASSY,SALT BRIDGE CAP,KYNA	\$	224.04
60A9937-001	ASY-VALVE+PIPE,INSERT,	\$	1,543.69
60A9944	ASSY,JUNCTION BOX,SUBMERS	\$	321.84
60B3F1020-102	FAB-MH,LOCK RING,LARGE,	\$	973.24
60B3F1021-102	FAB-MH,LOCK RING,SMALL,31	\$	846.04
60B3F1022-102	FAB-MH,SEALING HUB,316 SS	\$	616.70

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

60F1019-005	si FAB-MH,ADAPTER,CPVC	\$	300.64
60F1025-102	si FAB-MH,BOSS,CPVC	\$	151.86
60F1059-101	FAB-UNION,TEE,PVC	\$	593.58
60F2024-001	FAB-ADAPTER,UPPER,CPVC	\$	190.02
60F2025-001	si FAB-WIPER,1-1/2"INSERTION	\$	87.11
60F2026-001	si FAB-NIPPLE,SHORT, 1-1/2"	\$	375.80
60F2027-001	si FAB-BACK TUBE,1-1/2" PVC	\$	431.69
60F2031-001	si FAB-ADAPTER,LOWER,CPVC	\$	477.95
60F2032-001	si FAB-LOCK RING,LARGE,MOD-	\$	292.93
60G2052	J-BOX,PIPE MNT,1" PVC BUS	\$	321.84
60G9753-101	BOM-UNION ADAPTER,LCP,VIT	\$	555.03
60G9753-104	BOM-UNION ADAPTER,LCP,EPR	\$	555.03
60G9753-301	BOM-UNION ADAPTER,RYTON,V	\$	585.87
611213H	KTO: SOFCHEK TOTAL HDNS (50) *SELL 12	\$	227.41
611217H	KTO: HACH SOFCHEK FOIL INDVDUAL pk/1000	\$	1,210.28
611434	MAGNESIUM CHLORIDE ACS 500G	\$	236.68
6120218	MANUAL, DIFF. PH SMART PROBE W/sc100	\$	155.53
6120500	Differential pH/ORP Digital Gateway	\$	514.56
6120500.99	Differential pH/ORP Digital Gateway	\$	514.56
6120600	Combination pH/ORP Digital Gateway	\$	423.98
6120600.99	Combination pH/ORP Digital Gateway	\$	423.98
6120700	Contacting Conductivity Gateway	\$	414.35
6120700.99	Contacting Conductivity Gateway	\$	304.50
6120800	ELECTRODELESS CONDUCTIVITY GATEWAY	\$	541.54
6120800.99	ELECTRODELESS CONDUCTIVITY GATEWAY	\$	514.56
6121300	CKT BD E&D, DIG COMBO PH/ORP/D.O.	\$	121.41
6121400	si CKT BD ASSY, DIG COMBO PH/ORP	\$	1,113.92
6121800	PIPE, INSERTION MOUNT, CPVC	\$	82.68
6122100	si CORD GRIP, DIGITAL PH PROBE, RYTON	\$	58.01
6122400	CABLE ASSY, 1.0M PROBE EXTENSION	\$	591.65
6122401	CABLE ASSY, 1.0M PROBE EXTENSION	\$	1,083.09
6122600	si FRONT, DIGITAL MUNICIPAL PH PROBE	\$	122.38
6122900	CABLE ASSY, 4 COND SHLD W/CIRCULAR, 10M	\$	346.90
6122900-030	CABLE ASSY, 4 COND SHLD W/CIRCULAR, 30M	\$	813.28
6126700	nn SENSOR, GALVANIC, MEMBRANE D.O.	\$	708.25
6126800	si PIPE, UNION MOUNT, 1 1/2" CPVC	\$	271.74
6128701	si STRAIN RELIEF, MODIFIED	\$	86.53
6130088	MANUAL, NEMA 4 PC 1 CM CHELANT	\$	87.69
6131000	KIT, DIG D.O. & PH FLOATATION MT	\$	2,543.90
6131100	ADAPTER ASSY, UNION MT, 1 1/2", PVC	\$	732.34
6131300	KIT, UNION MTG, CPVC	\$	1,219.92
6131400	KIT, UNION MOUNTING, 316SS	\$	2,177.74
6131800	CONN, RCPT 5 CKT CIRCULAR PNL, MALE	\$	198.50
6134300	SCREW DRIVER 2.0 MM WIDE BLADE	\$	72.85
6136300	KIT, 5740 DIG D.O. FLOW-THRU UNION	\$	1,129.34
6136400	KIT, DIG PH IMMERSION MTG, CPVC	\$	738.12
6136500	KIT, DIG PH IMMERSION MTG, SST	\$	1,219.92

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

6136600	INNER FLOAT PIPE, DIG. D.O. MOUNTING	\$	240.90
6138200	SEALING HUB, CPVC	\$	202.36
613826	L-ASCORBIC ACID, 100G	\$	125.56
6138500	ADAPTER, UPPER, SST	\$	323.77
613855	L-ASCORBIC ACID, 290.0G	\$	163.69
6138600	COUPLER, 3/4 NPT, SCHED 80 PVC	\$	46.83
6138700	ADAPTER, LOWER, CPVC	\$	202.36
6139900	CABLE LOCK, FOR CLASS 1 DIV 2 CABLES	\$	71.69
6146400A	FTG, TEE UNION 1/4OD TBGX1/4OD TBG, NOT FDA CONTROLLED	\$	33.73
6149600	BARRIER, RELAY, sc100	\$	7.71
6157400	KTO: AMTAX sc ANLZR 115-230V 0.05-20	\$	60,359.90
6157500	KTO: AMTAX sc ANLZR 115-230V 0.05-20	\$	62,980.90
6158000	KTO: AMTAX sc ANLZR 115-230V 1-100MG	\$	60,359.90
6158100	KTO: AMTAX sc ANLZR 115-230V 1-100MG	\$	62,980.90
6158600	KTO: AMTAX sc ANLZR 115-230V 10-1000	\$	60,359.90
6158700	KTO: AMTAX sc ANLZR 115-230V 10-1000	\$	62,980.90
6159600	KTO: PHOSPHAX sc ANLZR 115-230V .05-	\$	59,010.86
6159700	KTO:PHOSPHAX sc ANLZR 115-230V .05-	\$	61,689.67
6159800	KTO: PHOSPHAX sc ANLZR 115-230V 1-50	\$	59,010.86
6159900	KTO:PHOSPHAX sc ANLZR 115-230V 1-50	\$	61,689.67
6166300	VALVE, SPEED FIT TO MALE NPTF ADAPT.	\$	107.54
6168100	BRACKET, PANEL MOUNT, sc1000	\$	784.37
6169900	KIT, PANEL MOUNT, SC1000	\$	770.88
6170800	HOSE ASSY, LDO PROBE AIR BLAST SYS	\$	558.89
6172000	FUSE, 8A 250V SLOW-BLOW 5 X 20MM	\$	5.40
6172100	si ADAPTER, AIR BLASTER	\$	360.39
6172200	FTG, 90 DEG STREET ELL 1/4 NPT	\$	57.05
6173900	FTG, AIR, 1/4 NPT QUICK DISCONNECT	\$	53.58
6175600	POWER SUPPLY, 90-264VAC, 47-63 HZ	\$	1,537.91
6175800	CABLE ASSY, POWER SUPPLY, APA6000	\$	277.52
6175900	KIT, RETROFIT, APA6000 POWER SUPPLY	\$	2,351.18
6177800	si PIPE, MOUNTING, 2 METER, PVC	\$	323.77
6178000	8362sc HIGH PURITY PH PANEL, HACH	\$	6,147.77
6178000S	8362sc High Purity pH, Steel Tagging	\$	6,186.31
6178000S-5024	vv 6178000S with Cfl. Coated Board	\$	7,882.25
6178001	8362sc HIGH PURITY ORP PANEL, HACH	\$	6,147.77
6178002	8362sc HIGH PURITY PH PANEL, LANGE	\$	6,147.77
6178003	8362sc HIGH PURITY ORP PANEL, LANGE	\$	6,147.77
6178060	8362sc Gateway Junction Box, Hach	\$	1,709.43
6178061	8362sc Gateway Junction Box, Lange	\$	1,611.14
6178700	TUBE, BACK HIGH PURITY PH/ORP PANEL	\$	58.01
6178900	si TUBE,FRONT,HIGH PURITY PH/ORP PANEL	\$	58.01
6181100	KIT, MAINTENANCE, MOLYBDENUM ANALYZ	\$	724.63
6182500	KTO: AMTAX MEMBRANE CAPS AND ELECTR.	\$	794.97
6184200	KIT, REPLACEMENT WIPERS, ACCUCLOR	\$	237.05
6184750	LDO SENSOR GUARD	\$	254.39

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

6184900	KIT, MTG HARDWARE, ISE SENSORS, RAIL MOUNT PVC	\$	1,659.32
6188400	SENSOR CARTRIDGE, AMMONIUM PROBE	\$	3,603.86
6188401	SENSOR CARTRIDGE NO3Dsc PROBE	\$	3,411.14
6189000	aa oo SAMPLE CUP ASSY, ACCUCLOR	\$	1,965.74
6190200	HEAD, AIR BLAST CLEANING SYSTEM	\$	379.66
6190250	LDO AIR BLAST HEAD KIT	\$	520.34
6200600	COLORIMETER/HTR BLOCK ASSY, HRH	\$	3,102.79
6200900	KIT, SAMPLE SEQUENCING, APA6000	\$	2,100.65
6200989	INSTRUCTION SHT, SAMPLE SEQUENCING	\$	29.68
62011	BOTTLE, WASH 500ML	\$	26.60
6201100	FILTER, CERAMIC 0.5 MICRON	\$	341.11
6201200	SEAL, QUAD-RING .296 NOM ID	\$	61.09
62014	BOTTLE, WASH 125ML	\$	19.46
6201400	O-RING FLCAR .070 W X .676 ID	\$	61.28
62016	BOTTLE, WASH 1000ML	\$	39.31
6201900	si HOLDER, FILTER, 1/4" TBG.	\$	637.90
6202800	KIT, MAINTENANCE, MANUAL, APA6000	\$	2,698.08
62031	BOTTLE, WASH 250ML	\$	23.51
6203500	FTG, FEMALE PIPE ADPT 1/4" TO 1/8"	\$	20.43
6203600	FTG, FEMALE LUER LOCK TO 1/8"NPT	\$	258.24
6203789	INSTRUCTION SHT, GRAB SAMPLE KIT	\$	29.68
6203802	KIT, SAMPLE COND APA6000 LR HARDNESS	\$	2,717.35
6203889	INSTRUCTION SHT, SAMPLE COND KIT	\$	29.68
6204400	GREASE, HI TEMP, NON-FLAMMABLE	\$	1,104.29
6204700	DRAIN BLOCK ASSY, REG 20PSI, APA6000	\$	2,062.10
6205000	VALVE MODULE ASSY, APA6000	\$	12,912.24
6205400	TEMP CONTROL BLK ASSY, ALK	\$	1,836.62
6205800	COLORIMETER ASSY, 520NM	\$	3,796.58
62075	BOTTLE, WASH 750ML	\$	33.92
6207500	GREASE PENCIL, BLACK	\$	25.63
621=188=210	MANUAL, 8810 INSTRUCTION (ENGLISH)	\$	229.34
621=188=310	MANUAL, 8810 INSTRUCTION (ENGLISH)	\$	211.99
621=188=410	NOTICE 8810 TYPE METHANOL EN ANGLAIS	\$	225.48
621=188=510	Instruction manual 8810 COLORIMETRIC Orthophosphates	\$	225.48
621=997=101	MANUAL, 8810 INSTALL, NH4 0-100mg/l	\$	85.37
621=997=102	MANUAL, 8810 INSTALL, NH4 0-1000ug/l	\$	85.37
621=997=108	MANUAL FREE/TOTAL ALKA 0-500 MG EN	\$	97.32
621=997=116	MANUAL CHLORIDE 0.5-500 MG EN	\$	85.37
621=997=117	MANUAL SULPHIDE 1-100 MG EN	\$	91.54
621=997=119	Installation procedure 8810 in English Zinc in phosphatizing bath	\$	50.30
621=997=135	INSTAL. PROCEDURE FOR 8810 CHLORIDE 0-20PPM	\$	75.93
621=997=137	PI 8810 ISE CHLORURE 0-5PPM A4	\$	84.80
621=997=138	INSTALLATION PROCEDURE 8810 NITRITES	\$	79.02
62100	BOTTLE W/STOPPER, BOD 300ML	\$	57.24
62106	BOTTLE W/STPR, BOD 300ML PK/6	\$	323.77
62124	BOTTLE W/STPR, BOD 300ML PK/24	\$	1,092.72
621410	nn HACH AC 5-WAY WQ (50) pk/12	\$	375.80

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

621410H	KTO: AC 5-WAY WQ (50) pk/12	\$	375.80
6250700	O-RING CORD, 7/16 DIA., NEOPRENE	\$	13.49
62762	CASSEROLE, PORCELAIN 140ML	\$	114.48
62763	CASSEROLE, PORCELAIN 210ML	\$	147.05
63301	CLAMP, SCREW COMPRESSION LARGE	\$	16.00
63317	CLAMP, SCREW COMPRESSION SMALL	\$	12.72
63400	CLAMP, TEST TUBE	\$	16.00
63800	SPOON, MEASURING 0.2G	\$	16.19
641426	nn HACH AC NITRATE/NITRITE (25) PK/6	\$	260.17
641426H	KTO: AQUACHECK NITRATE/NITRITE (25) PK/6	\$	281.37
6422P0	SNSR-pH,RYTON,VITON,10'	\$	4,143.48
6422P020	SNSR-pH,RYTON,VITON,20'	\$	5,280.53
6422P040	Sensor,pH-Ryton,Viton,40', 2 Wire	\$	5,376.89
6428-9001	SNSR-pH,LCP,2W,VITON,5 FT	\$	4,818.00
6428P0	SNSR-pH,LCP,2W,VITON,10'	\$	3,835.13
6428P010010000N	SNSR-pH,LCP,2W,VITON,10'	\$	4,143.48
6428P015	PH SNSR,ENCPS,CONV,LCP,2W	\$	4,914.36
6428P020	SNSR-pH,LCP,2W,VITON,20'	\$	4,933.63
6428P025	SNSR-pH,LCP,2W,VITON,25'	\$	4,143.48
6428P030	SNSR-pH,LCP,2W,VITON,30'	\$	5,068.54
6428P040W00000N	PH-SNSR,2W TRANS,CONV,LCP	\$	6,571.75
6428P050	SNSR-pH,LCP,2W,VITON,50'	\$	4,394.02
6428P1	SNSR-pH,LCP,2W,ANTIMONY,1	\$	4,625.28
6458P0	SNSR-pH,2-WIRE,F-THRU,LCP	\$	4,355.47
651453	nn HACH AQUACHEK CHLORINE (50) pk/12	\$	294.86
651453H	KTO: AQUACHEK CHLORINE (50) pk/12	\$	294.86
657311	SCREW 8-32X .312PAN SST PHIL	\$	2.51
661454	nn HACH AQUACHEK COPPER (25) pk/6	\$	339.19
661455	nn HACH AQUACHEK IRON KIT (25) pk/6	\$	435.55
663711	SCREW 1/4 -20X1. HEX STL CAP	\$	3.47
67132	SODIUM HYDROXIDE, .01N 100ML MDB	\$	57.82
67132-LM	SODIUM HYDROXIDE, .01N 100ML MDB	\$	57.82
67153	SODIUM HYDROXIDE SOLN, 1000ML	\$	98.46
671G3120-101	ASSY-MDL.DSP SCL BD.MOD.	\$	3,122.06
672G5001-101	PCB,POWER SUPPLY BD.,1 FU	\$	1,289.30
6735500	VERIFICATION QUICK CHECK, FT660 sc	\$	1,884.80
6735600	KIT, APERTURE REPLACEMENT, FT660 sc	\$	33.15
673644	DIODE 1N914 ON TAPE, 52MM	\$	16.57
6737100	KIT, HARSH ENVIRONMENT LASER ASSY, FT660	\$	3,661.68
6745000	INSERT, HC-5 MOLDED MP METER	\$	69.57
67532	SODIUM CARBONATE SOLN, 100ML	\$	54.20
67532-LM	SODIUM CARBONATE SOLN, 100ML	\$	54.20
67549	SODIUM CARBONATE SOLN, 500ML	\$	119.05
67549-LM	SODIUM CARBONATE SOLN, 500ML	\$	119.05
6760100	COVER, REAGENT TUBING COMPARTMENT, PROCESS ANALYZER	\$	124.30
6760200	SHIPPING CARTON, PROCESS ANALYZER	\$	1,004.07
6760301	ee ASSY, REAGENT BOTTLE CAP, BLUE	\$	116.98

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

6760302	ee ASSY, REAGENT BOTTLE CAP, GREEN	\$	116.98
6760303	ee ASSY, REAGENT BOTTLE CAP, RED	\$	116.98
6760304	ASSY, REAGENT BOTTLE CAP, PURPLE	\$	117.17
6760305	ee ASSY, REAGENT BOTTLE CAP, YELLOW	\$	116.98
6760306	ee ASSY, REAGENT BOTTLE CAP, BLACK	\$	117.17
6760307	ASSY, REAGENT BOTTLE CAP, ORANGE	\$	127.00
6760800	CLAMP, ANALYZER PCB, PROCESS ANALYZER	\$	23.13
6760900	COVER, FLUIDICS, BOTTOM PANEL, PROCESS ANALYZER	\$	219.70
6761600	si CKT BD ASSY, INTERCONNECT, DC/DC, PROCESS ANALYZER	\$	8,190.60
6761800	NEEDLE VALVE, MINIATURE, BARBED	\$	499.14
6763000	ee POWER SUPPLY, 120-240V, 12VDC, 65W	\$	849.90
676312	NUT 6-32 HEX KEPS EXT LCK ZN	\$	11.37
6763600	ee si CKT BD ASSY, INTERCONNECT, AC/DC, PROCESS ANALYZER	\$	7,227.00
6763800	CKT BD ASSY, COLORIMETER, SILICA, PROCESS ANALYZER	\$	1,582.23
6764500	ASSY, ENCLOSURE, PROCESS ANALYZER	\$	8,498.95
6764701	DOOR, UPPER, ENCLOSURE, PROCESS ANALYZER	\$	373.88
6764702	DOOR, LOWER, ENCLOSURE, PROCESS ANALYZER	\$	371.95
6764800	SEAL, BULKHEAD, ENCLOSURE, PROCESS ANALYZER	\$	10.60
6764900	FOOT, ENCLOSURE, PROCESS ANALYZER	\$	69.57
6765400	CATCH PLATE, LATCH, ENCLOSURE, PROCESS ANALYZER	\$	184.24
6765600	ee KEYPAD, MEMBRANE, S5500	\$	262.10
6765800	ANALYTICS PANEL, PROCESS ANALYZER	\$	356.53
6766000	HINGE ASSY, ENCLOSURE, PROCESS ANALYZER	\$	106.38
6766103	INPUT PANEL, ELECTRICAL, PROCESS ANALYZER	\$	208.14
6766122	ASSY, INPUT PANEL, FLUIDICS, 10-PORT	\$	323.77
6766124	ASSY, INPUT PANEL, FLUIDICS, 13-PORT	\$	341.11
6766200	SCREW M5 X 8MM TORX 316 SST PANHEAD	\$	5.78
6766300	LATCH PLATE, ANALYTICS PANEL, PROCESS ANALYZER	\$	160.92
6766400	PLATE, CONDUIT, PROCESS ANALYZER	\$	92.31
6766600	COVER, INTERCONNECT PCB, LOW-VOLTAGE	\$	219.70
6766800	COVER, ANALYTICS PCB, PROCESS ANALYZER	\$	78.05
6766900	COVER, COLORIMETER, PROCESS ANALYZER	\$	137.79
6767100	FUNNEL, GRAB SAMPLE, PROCESS ANALYZER	\$	86.92
6767300	BOTTLE, REAGENT, 2-LITER, PETE	\$	4.24
6767400	SHIPPING CAP, REAGENT, QUICK DISCONNECT	\$	3.78
6767500	ASSY, SEAL, REAGENT BOTTLE, 2L	\$	32.88
6767510	SEAL ASSY, REAGENT BOTTLE, WITH TUBE	\$	55.70
6767520	SEAL ASSY, REAGENT BOTTLE, WITH TUBE, 1L	\$	46.42
6767600	BRACKET, WALL MOUNT, PROCESS ANALYZER	\$	165.35
6767800	si CELL CAP, PROCESS ANALYZER	\$	65.52
6767900	ADAPTOR, MOTOR MAGNET	\$	13.10
6768000	si CELL, COLORIMETER, PROCESS ANALYZER	\$	346.90
6768200	PANEL, BULKHEAD, ENCLOSURE, PROCESS ANALYZER	\$	277.52
6768300	BRACKET, LOCK, PROCESS ANALYZER	\$	92.31
6768500	TRAY STOP, ENCLOSURE, PROCESS ANALYZER	\$	115.05
6769200	AIR COMPRESSOR, 12VDC	\$	520.34
6769300	MAGNET, CYLINDER, 3/16" X 1/2", N42	\$	26.21



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

6769400	FTG, BARB, 1/4-28 THREAD, 3/32" ID TUBING, POLYPROPYLENE	\$	23.13
6770500	CONN, STRAIN RELIEF, SKINTOP, PG-09	\$	13.30
6770800	DIFFUSER, STATUS LED, PROCESS ANALYZER	\$	41.43
6771200	BRACKET, PRESSURE DUMP VALVE	\$	18.50
6771300	CLAMP, POLYPROPYLENE LOOP FOR 5/8" OD	\$	5.01
6772100	CORE, SAMPLE HEATER, PROCESS ANALYZER	\$	458.67
6772400	COIL, SAMPLE HEATER	\$	770.88
6772500	TUBING, POLYETHYLENE, 6MM OD X 1MM W	\$	8.48
6772600	si STIR BAR, STAR HEAD, 17.2 MM DIA X 10.2 HIGH	\$	73.62
6772800	FTG, 1/2" ELBOW, QUICK CONNECT, POLYPROPYLENE	\$	179.81
6773000	MANIFOLD, AIR, PROCESS ANALYZER	\$	277.52
6773100	SHROUD, COLORIMETER, PROCESS ANALYZER	\$	69.57
6773200	KIT, CABLE ASSY, EXTERNAL CONTROLLER, PROCESS ANALYZER	\$	585.87
6773300	CLIP, SPRING MOUNTING, SST FOR 1" VALVES	\$	63.98
6773400	GUARD, PRESSURE TRANSDUCER, PROCESS ANALYZER	\$	13.88
6773500	COVER, FUNNEL, PROCESS ANALYZER	\$	137.79
6773700	LCD MOUNT, PROCESS ANALYZER	\$	46.25
6773800	ee CABLE ASSY, MOTOR, PROCESS ANALYZER	\$	1,204.50
6774000	CABLE ASSY, OPTIC LED	\$	352.68
6774802	5500sc REAGENT 1 SILICA 2-LITER	\$	691.99
6774802-LM	5500sc REAGENT 1 SILICA 2-LITER	\$	691.99
6774902	5500sc REAGENT 2 SILICA 2-LITER	\$	290.89
6774902-LM	5500sc REAGENT 2 SILICA 2-LITER	\$	290.89
6775002	5500sc STANDARD 1 SILICA 2-LITER	\$	180.48
6775102	KTO: 5500sc REAGENT 3 SILICA PACKAGE	\$	377.61
6775102-LM	KTO: 5500sc REAGENT 3 SILICA PACKAGE	\$	377.61
6775202	5500sc REAGENT 3 SILICA 1.8-LITER	\$	209.58
6775202-LM	5500sc REAGENT 3 SILICA 1.8-LITER	\$	209.58
6775355	5500sc REAGENT 3 SILICA POWDER 205G	\$	144.35
6775355-LM	5500sc REAGENT 3 SILICA POWDER 205G	\$	144.35
6775402	5500sc REAGENT 1 LR PHOSPHATE 2,1-LITER INT	\$	350.51
6775402-LM	5500sc REAGENT 1 LR PHOSPHATE 2,1-LITER INT	\$	350.51
6775502	5500sc REAGENT 2 LR PHOSPHATE 2-LITER	\$	290.89
6775502-LM	5500sc REAGENT 2 LR PHOSPHATE 2-LITER	\$	290.89
6775602	5500sc STANDARD 2 LR PHOSPHATE	\$	623.33
6775702	KTO: 5500sc REAGENT 3 LR PHOSPHATE PACKAGE	\$	522.15
6775855	5500sc REAGENT 3 LR PHOSPHATE POWDER	\$	296.31
6775902	5500sc REAGENT 3 LR PHOSPHATE 1.8-LITER	\$	388.45
6776002	5500sc STANDARD 1 LR PHOSPHATE 2-LITER	\$	292.69
6776100	KTO: 5500SC PHOSPHATE HR RGT SET	\$	1,468.89
6776100-LM	KTO: 5500SC PHOSPHATE HR RGT SET	\$	1,468.89
6776102	5500sc REAGENT 1, HR PHOSPHATE, 2L	\$	706.44
6776102-LM	5500sc REAGENT 1, HR PHOSPHATE, 2L	\$	706.44
6776202	5500sc REAGENT 2 HR PHOSPHATE 2-LITER	\$	419.17
6776202-LM	5500sc REAGENT 2 HR PHOSPHATE 2-LITER	\$	419.17
6776302	5500sc REAGENT 3, HR PHOSPHATE, 2, 1-LITER INT	\$	392.06
6776302-LM	5500sc REAGENT 3, HR PHOSPHATE, 2, 1-LITER INT	\$	392.06

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

6776402	5500sc STANDARD 1 HR PHOSPHATE 2-LITER	\$	771.48
6776500	CABLE ASSY, LOWER DOOR, PROCESS ANALYZER	\$	346.90
6777401	HOUSING, PUSHROD, PRESSURE DUMP VALVE	\$	23.13
6777402	PINCH BOSS, PRESSURE DUMP VALVE	\$	23.13
6778600	si CKT BD ASSY, OPEN DOOR DETECTOR	\$	63.98
6779100	si CABLE ASSY, THERMISTOR, REAGENT	\$	450.96
6779200	si CABLE ASSY, LEAK DETECTOR	\$	946.26
6779900	CHECK VALVE, MINIATURE BARBED 1/8" TBG	\$	83.06
677G1050-201	RPLCMNT,DOOR,0-35 GPM w/T	\$	2,119.92
6780600	DISPLAY, LCD 5.7" COLOR	\$	2,659.54
6782100	CABLE ASSY, AC POWER, PROCESS ANALYZER	\$	277.52
6782200	CABLE ASSY, 12V POWER, PROCESS ANALYZER	\$	676.45
6782300	si CABLE ASSY, ANALYZER POWER	\$	1,160.17
6782400	si CABLE ASSY, ANALYZER DATA, PROCESS ANALYZER	\$	414.35
6782500	si CABLE ASSY, BACKLIGHT, PROCESS ANALYZER	\$	759.32
6782600	CABLE ASSY, THERMISTOR, SAMPLE HEATER	\$	574.31
6782700	ASSY, HEATER CARTRIDGE, PROCESS ANALYZER	\$	1,842.40
6782900	PRESSURE REGULATOR, PRESET 4PSI	\$	932.76
6783200	si CABLE ASSY, CHASSIS GROUND, PROCESS ANALYZER	\$	87.30
6783300	si ASSY, HEATER CARTRIDGE, DC, PROCESS ANALYZER	\$	1,813.50
6783500	KIT, INSTALLATION, PROCESS ANALYZER	\$	1,031.05
6783600	KTO: 5500SC SILICA REAGENT SET	\$	1,468.89
6783600-LM	KTO: 5500SC SILICA REAGENT SET	\$	1,468.89
6783700	si CABLE ASSY, REAGENT VALVE, PROCESS ANALYZER	\$	1,946.47
6784000	FTG, REDUCING BARB, 1/8" TO 3/32" ID TBG, POLYPRO	\$	10.60
6784100	SEAL, 15.5MM HOLE	\$	10.60
6784200	SEAL LOCK, 15.5MM HOLE	\$	10.60
6784300	SEAL, 21.4 MM HOLE	\$	10.60
6784400	SEAL LOCK, 21.4 MM HOLE	\$	10.60
6784500	ee CABLE ASSY, AIR COMPRESSOR, PROCESS ANALYZER	\$	3,353.33
6784700	CABLE, FLEX, FFC CARD FOR LCD	\$	78.05
6784800	ee ASSY, Y-STRAINER, PROCESS ANALYZER	\$	1,318.20
6785102	KIT, SEQUENCER LINE INSTALLATION 2 CHANNEL	\$	971.31
6785104	KIT, SEQUENCER LINE INSTALLATION 4 CHANNEL	\$	2,409.00
6785106	KIT, SEQUENCER LINE INSTALLATION 6 CHANNEL	\$	3,854.40
678517	CLAMP, CABLE-HOLDS 1/4 DIA NYL	\$	4.05
6785600	TUBING, PVC 11/16 OD x 1/8 W	\$	12.72
6786300	si CABLE ASSY, PINCH VALVE, 2 WAY, NC, PROCESS ANALYZER	\$	936.62
6786400	si CABLE ASSY, PINCH VALVE, 2 WAY, PROCESS ANALYZER	\$	2,717.35
6786500	si CABLE ASSY, PINCH VALVE, 3 WAY, PROCESS ANALYZER	\$	3,353.33
6786600	KIT, STAINLESS STEEL SAMPLE CONDITIONING, PROCESS ANALYZER	\$	10,638.14
6786800	ASSY, COLORIMETER, SILICA	\$	2,678.81
6786801	ASSY, COLORIMETER, LR PHOSPHATE, W/ PKG	\$	2,678.81
6786802	ASSY, COLORIMETER, HR PHOSPHATE, W/ PKG	\$	2,678.81
6786900	ASSY, CAPILLARY TBG, PROCESS ANALYZER	\$	314.13
6786901	ASSY, CAPILLARY TBG, 4 LINES, PROCESS ANALYZER	\$	528.05

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

6786902	ASSY, CAPILLARY, 3 LINES, PROCESS ANALYZER	\$	495.29
6787000	KIT, ADAPTER PLATE, S5000 TO 5500sc	\$	1,081.16
6787100	KIT, ADAPTER PLATE, 921x TO 5500sc OR 961Xsc	\$	1,081.16
6787200	PLUG, FOR 6MM TUBING	\$	37.19
6788200	si CKT BD ASSY, SD CARD BOARD PROCESS ANALYZER	\$	383.51
6788301	KIT, MAINTENANCE, 5500sc, SILICA, SINGLE CHANNEL	\$	1,780.73
6788302	KIT, MAINTENANCE, 5500sc, SILICA, TWO/FOUR CHANNEL	\$	1,817.35
6788303	KIT, MAINTENANCE, 5500sc, SILICA, SIX CHANNEL	\$	1,859.75
6788304	KIT, MAINTENANCE, 9610sc, SILICA, SINGLE CHANNEL	\$	1,780.73
6788305	KIT, MAINTENANCE, 9610sc, SILICA, TWO/FOUR CHANNEL	\$	1,817.35
6788306	KIT, MAINTENANCE, 9610sc, SILICA, SIX CHANNEL	\$	1,859.75
6788307	KIT, MAINTENANCE, 5500sc/9611sc, LR PHOSPHATE, SINGLE CHANNEL	\$	1,780.73
6788308	KIT, MAINTENANCE, 5500sc/9611sc, LR PHOSPHATE, TWO/FOUR CHANNEL	\$	1,817.35
6788309	KIT, MAINTENANCE, 5500sc/9611sc, HR PHOSPHATE, SINGLE CHANNEL	\$	1,780.73
6788310	KIT, MAINTENANCE, 5500sc/9611sc HR PHOSPHATE, TWO/FOUR CHANNEL	\$	1,715.21
67883SER1	KIT, TUBING REPLACEMENT #1	\$	105.03
6788400	FTG, REDUCING TEE, 3/8 X 3/8 X 5/32 POLYPROPYLENE	\$	13.88
6788500	ADAPTER PLATE, THERMAL FUSE, PROCESS ANALYZER	\$	178.07
6788800	TORX DRIVER, T20	\$	106.38
6788900	TORX DRIVER, T25	\$	106.38
6789100	KIT, FAN FILTER, PROCESS ANALYZER, HACH	\$	184.24
6789300	KIT, FAN PLUG, PROCESS ANALYZER, HACH	\$	160.92
6789800	si ASSY, BLOWER, PROCESS ANALYZER	\$	1,287.37
6790100	ASSY, SAMPLE VALVE, PROCESS ANALYZER	\$	188.67
6790300	BRACKET, BYPASS VALVE, PROCESS ANALYZER	\$	108.50
679045	TIE, CABLE, 5.6" LONG X .10" W	\$	12.53
6790600	ee ASSY, FAN FILTER, PROCESS ANALYZER	\$	50.30
6790700	TBG, SILICONE 1/4" OD x 1/16"W	\$	56.85
6791500	FTG, 1/4-28 UNF TO 1/8 BARB POLYPRO	\$	23.13
6791700	LABEL SET, REAGENT VALVES	\$	31.22
6792200	KIT, 4-20mA OUTPUT MODULE, HACH	\$	1,753.75
6792201	KIT, 4-20mA OUTPUT MODULE, HACH LANGE	\$	1,753.75
6792301	KIT, 4-20mA INPUT MODULE, HACH LANGE	\$	2,312.64
6792400	KIT, RELAY MODULE, HACH	\$	1,801.93
6792401	KIT, RELAY MODULE, HACH LANGE	\$	1,801.93
6792500	KIT, MODBUS MODULE, HACH	\$	819.06
6792501	KIT, MODBUS MODULE, HACH LANGE	\$	819.06
6792600	KIT, PROFIBUS MODULE, HACH	\$	2,274.10
6792601	KIT, PROFIBUS MODULE, HACH LANGE	\$	2,274.10
6793000	FTG, TEE 3/32" BARB POLYPROPYLENE	\$	18.50
6793200	FTG, ELBOW 3/32" BARB POLYPROPYLENE	\$	18.50
6793702	ASSY, SEQUENCER, 2 CHANNEL, PROCESS ANALYZER	\$	4,798.73
6793704	ASSY, SEQUENCER, 4 CHANNEL, PROCESS ANALYZER	\$	14,627.45

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

6793706	ASSY, SEQUENCER, 6 CHANNEL, PROCESS ANALYZER	\$	24,379.08
6793900	FTG, 1/8 BARB TO 1/4-18 FNPT	\$	13.88
6794000	PIN, NEEDLE VALVE, SEQUENCER, PROCESS ANALYZER	\$	61.86
6794100	si ASSY, SEQUENCER VALVE	\$	2,987.16
6794300	VALVE, TOGGLE, 1/8 FNPT TO 1/8 FNPT	\$	393.15
680G1020	BOM-ASSY,DOOR REPLACEMENT	\$	601.29
68206	CLAMP, PINCHCOCK DAY PK/6	\$	44.90
6836000	oo UPGRADE KIT,ACCUCHLOR SAMPLE CELL	\$	3,064.25
6836300	O-RING NIT .139 W X 2.359 ID	\$	13.49
6840000	si PLUG, PROBE ADAPTOR, WDMpsc	\$	200.43
6842700	ASSEMBLY, DRAIN MANIFOLD, WDMpsc	\$	520.34
6846600	PRESSURE REGULATOR ACCESSORY, WDMpsc	\$	763.17
68500	BRUSH, BURET 25ML	\$	12.72
68600	BRUSH, BURET 50ML	\$	16.00
6860000	CLEAN SYS,HI OUTPUT AIR BL HACH 115V	\$	7,149.91
6860100	HIGH OUTPUT AIR BLAST CLEANING SYSTEMS (230V)	\$	7,149.91
6861000	ENCLOSURE, MODIFIED, HI OUTPUT AIR BLAST	\$	1,428.06
6862400	COMPRESSOR, MODIFIED 230V	\$	2,197.01
6865600	STIR BAR, MICRO 4.93MM SPHERICAL	\$	34.11
6866100	si SEAL, SALT BRIDGE, OUTER, PHD SANTOPRENE	\$	16.57
6866700	si SEAL, SALTBRIDGE, INNER, PHD, SANTOPRENE	\$	16.57
6867000	COLORIMETER ASSY, 510NM, CL17	\$	2,775.17
6868500	si PLUG, CL17 COLORIMETER	\$	29.29
6868800	INSTALLATION KIT, MO42 ANALYZER	\$	1,462.74
6868900	SCREW 6-32 X 1.25 SKT HD SHOULDER SST	\$	46.25
68700	BRUSH, CYLINDER SIZE 2	\$	16.00
68800	BRUSH, IMHOFF CONE	\$	76.32
689=132=008	SOLENOID VALVE, 3/2 WAY, NPS 0.8	\$	1,406.86
689=132=024	SOLENOID VALVE ON OVERFLOW VESSEL	\$	1,393.37
68900	BRUSH, PIPET	\$	17.54
6892500	1900C PORTABLE TURBIDIMETER, EPA	\$	3,526.78
6892900	nn KEYBOARD, 3X3 X-Y MATRIX, 1900C	\$	82.48
69000	BRUSH, TEST TUBE	\$	16.00
69100	BRUSH, SEMI MICRO TUBE	\$	16.00
69257	FILTER PAPER, FOLD 12.5CM PK/100	\$	185.97
69258	FILTER PAPER, FOLDED 15CM PK/100	\$	202.36
69259	FILTER PAPER, FOLD 18.5CM PK/100	\$	231.26
6930100	Sanitary pH Sensor, Analog,	\$	4,066.39
694=000=001	FLOWMETER, 0-16L/H, W/ADJ VALVE,9186	\$	2,081.38
6940500	FILTER BACKWASH MOUNTING KIT	\$	3,044.98
6940600	ff KTO: FILTER BACKWASH PACKAGE, KTO	\$	20,852.30
6941700	EXT POLE for FILTER BACKWASH MOUNT	\$	855.68
6943000	si BRACKET CAP, PROBE MANIFOLD, WDMpsc	\$	379.66
694712	NUT 10-32 HEX KEPS EXT LCK ZN	\$	9.06
694912	NUT 10-32 HEX MACH SCREW SS	\$	5.59
695=004=004	PUMP 24VDC, MEMBRANE, 0.1L/MN	\$	957.82
695=004=006	Gas-pump 80V, resistant to DIPA	\$	3,738.77

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

695=100=001	Head of Piston Pump (spare part of 08810=A=7150, 8810_piston-pump)	\$	5,916.50
695=114=000	POMPE MICRODEBIT 24VCC A IMPULSIONS	\$	16,747.37
695=114=005	MICRO-PULSES PUMP, REGENERATION	\$	1,580.30
695=215=024	Pump for acid. unit LZY051 build after June 2014	\$	6,379.03
695612	WASHER #6 FLAT -REGULAR- SST	\$	8.67
695912	WASHER 1/4 FLAT -LIGHT- SST	\$	8.67
696=146=001	VANNE A AIGUILLE DN4/6 EN INOX	\$	3,334.06
696112	WASHER #10 FLAT -LIGHT- SST	\$	7.71
6961700	mc SAMPLER MAINTENANCE KIT, HST	\$	1,713.28
6H1184	VALVE,3-WAY,PVC 1/2" THD/	\$	1,499.36
6H1187	REPAIR KIT,D.I. INSRT FLO	\$	774.73
6M0920	GREASE, SILICONE, (2 GRAM	\$	31.41
6M1057	GREASE, BEARINGE, DOW#33	\$	33.34
6P95-1100-00	WDM TOC PREFPACK, CRS, 115V, WINDOW	\$	2,601.72
6P95-2100-00	WDM TOC PREFPACK,CRS,230V	\$	2,601.72
6Q1061	CMOS,ANALG.PROC.A/D,(LTB)	\$	162.46
6Q1081	LCD,14 DIG.,MULTIPLEXED,	\$	618.63
6Q1084	CMOS,8-BIT MCU (LTB)	\$	194.65
6Q1104	CMOS,QUAD L-H VOLT.SHIFTR	\$	16.00
6Q1156	LCD,5x7,8-CHAR.DOT MATRIX	\$	229.34
6Q1157	LCD,4-DIG.W/COLON,7-SEG.(\$	294.86
6Q1175	68HC68T1, SPI REAL-TIME C	\$	65.14
6V0010	ZERO POWDER	\$	98.09
700534	EDTA 500G	\$	184.29
700599	EDTA REAGENT PWD PLWS PK/100	\$	105.88
700599-LM	EDTA REAGENT PWD PLWS PK/100	\$	105.88
700726	CDTA, ACS GRADE 100G	\$	316.18
70100	LAB APRON, VINYL 91X69CM	\$	95.01
701137	SENTRY THERMAL SHUTOFF VALVE	\$	4,191.66
70532	SULFITE 3 REAGENT 100ML MDB	\$	37.75
70553	SULFITE 3 REAGENT SOLN, 1000ML	\$	106.04
7060824	TEFLON THREAD TAPE, .25 WIDE	\$	42.98
7061685	IC EPROM UV-ERASE 32K(4KX8)2732A	\$	99.06
7061791	NUT, 6-32 HEX KEPS EXT LCK SST	\$	9.25
706834	FERROUS SULFATE ACS 500G	\$	216.81
709901	SODIUM SULFATE, ANHYDROUS 454G	\$	113.63
709929	SODIUM SULFATE, ANHYDROUS 250G	\$	59.98
70F1037-003	si MH,SPEC.CAP,2" SANITARY,3	\$	294.86
711195	nn HACH LR QUANTAB (40) pk/6	\$	578.16
711195H	KTO: HACH LR QUANTAB (40) PK/6	\$	578.16
711196	nn HACH HR QUANTAB (40) pk/6	\$	1,131.27
711196H	KTO: HACH HR QUANTAB (40) PK/6	\$	830.62
7121649	STABLCAL 400 NTU, 500 ML	\$	220.42
7121649-LM	STABLCAL 400 NTU, 500 ML	\$	220.42
7122100	ASSEMBLY, DETECTOR, SS7	\$	3,950.76
71301	BORIC ACID, ACS 454G	\$	211.39



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

7132000	PLUG, 3.5MM, 2 CKT	\$	12.53
713400	TUBING, RUBBER LATEX F/1864	\$	14.07
713765	BOTTLE, M ROUND HDPE 1000ML PK/3	\$	29.29
713849	BOTTLE, MR HDPE 28-410 500ML	\$	1.16
7139800	FTG,MALE LUER SLIP BARB, 1/16" ID	\$	8.48
713B4A1001	ASSY,PCB.PREAMP PH,LOW PW	\$	1,023.34
714421	BOTTLE, BR AMB GLASS 473ML PK/6	\$	88.46
714424	BOTTLE, BR AMB GLASS 118ML PK/6	\$	48.76
714441	BOTTLE, BR AMB GLASS 237ML PK/6	\$	62.06
714463	BOTTLE, AMBER GLASS 1000ML PK/6	\$	115.05
714C4A1502-104N	Model 714 Preamplifier	\$	830.62
714G1502-104	ASSY-PREAMP/MTG PLATE,25'	\$	882.66
7213200	GEMS LS-3 SERIES LEVEL SWITCH	\$	535.76
7300100	si TOP, LDO FLOW CELL	\$	356.53
7300400	O-RING, 2-234 EPDM	\$	29.10
7300700	ASSY. LDO FLOW CELL	\$	2,081.38
7300800	ASSY, PHD & NH4D FLOW CELL	\$	1,491.65
732514	O-RING FLCAR .103W X .549ID	\$	21.58
74053	CDTA STD SOLN, 0.020N 1000ML	\$	106.04
74056	CDTA STD SOLN, 0.0200N 4L	\$	220.42
74226	HYDRAZINE SULFATE 100G	\$	129.36
742307	CABLE, 9 COND 2-18&7-22AWG STRD	\$	9.83
742707	CABLE, 5 COND 22 AWG STRD	\$	8.09
750907	CABLE, 2 COND, 22AWG, 7/30 STRD, NO SHIELD	\$	5.59
751111	SCREW 6- X .375PAN SST TPG	\$	11.37
75814	PHENOL, ACS 113G	\$	139.48
75814-LM	PHENOL, ACS 113G	\$	139.48
76215	PHOSPHORIC ACID, 85% ACS 2.5L	\$	278.24
76401	POTASSIUM CHLORIDE, 454G	\$	108.04
76401-LM	POTASSIUM CHLORIDE, 454G	\$	108.04
76714	POTASSIUM CYANIDE, 100G	\$	128.64
767512	WASHER 1/4 EXT LOCK TYPE A STL	\$	9.25
76A4000-001	ASSY,JUNCTION BOX,NEMA 4	\$	657.18
76F1006	FAB-MH,THREADED ADAPTER	\$	194.65
76G9902	BOM-ASSY,GUARD,DO,UNIT-SS	\$	1,337.48
76G9912	BOM-BALL FLOAT	\$	535.76
76G9913-001	BOM-SENSOR EXTENSION MOUN	\$	294.86
770013	Power Supply,+5,+/-15	\$	4,818.00
77601	SODIUM BICARBONATE 454G	\$	107.32
786712	#8 INT LOCK TYPE A SST	\$	2.51
787612	WASHER #6 INT LOCK TYPE A SST	\$	2.12
789500	SCREW 6-32X .375SET SST SLOT	\$	11.95
789700	SCREW 6-32X .5 SET SST SLOT	\$	5.78
7900SRB	SENSOR CABLE CLEANING	\$	3,777.31
79114	TANNIC ACID 113G	\$	213.20
79114-LM	TANNIC ACID 113G	\$	213.20
79501	ZINC 20 MESH ACS, 454G	\$	290.89

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

79501-LM	ZINC 20 MESH ACS, 454G	\$	290.89
798712	WASHER 3/8 INT LOCK TYPE A ZN	\$	1.93
7D0003	250 ML SHUNT FOR CHESSEL	\$	323.77
7E1202	KTO: MEASURING CELL TDO SEN THK MEM	\$	982.87
7E1222	END CAP FOR 7E1136, BLACK	\$	92.51
7E1230	WRAP NUT, GAS THREAD	\$	137.12
7E1231	PRONG PROT. RING FOR 7E12	\$	398.93
7F1053	FUSE 1 AMP 250V 5X20MM	\$	29.10
7K0735	BUSHING, CONN., .312 I.D.	\$	70.34
7W0200	JUMPER,TERM STRIP,.200"SP	\$	78.05
8000169	CAP FOR BOTTLE DIA (SODIUM)	\$	120.45
8001316	Panel Mounting Bracket	\$	37.58
80428	KTO: AMINO ACID RGT PP 25ML PK/1000	\$	706.44
80428-LM	KTO: AMINO ACID RGT PP 25ML PK/1000	\$	706.44
80499	AMINO ACID RGT PP 25ML PK/100	\$	79.50
80499-LM	AMINO ACID RGT PP 25ML PK/100	\$	79.50
81042H	CHROMIUM STD SOLN 50MG/L 100ML	\$	58.71
81842	HYDROXY HYDROCH, 100G/L 100ML	\$	57.98
81842-LM	HYDROXY HYDROCH, 100G/L 100ML	\$	57.98
819A1010	ASSY,CASE W/BUFFER+30MM V	\$	387.37
820=000=006	ELBOW TORX SCREW DRIVER T6	\$	90.39
8220-1090-002	SPARE CHAMBER W/PLUMB.HI	\$	7,862.98
8220G1320-202	SOURCE,I.R. PAIR, 10' CABLE	\$	1,031.05
82522	CALMAGITE, 5G	\$	155.55
8256600HS	VERIFICATION DYE, HIGH STANDARD, SL1000	\$	389.29
8256700MS	VERIFICATION DYE, MID STANDARD, SL1000	\$	551.18
8256800LS	VERIFICATION DYE, LOW STANDARD, SL1000	\$	416.28
8320A1190-001	VALVE ASSEMBLY, HI PRES.	\$	1,915.64
8326400	PLUG,SMALL HOLE	\$	1.35
8326500	PLUG,BIG HOLE	\$	1.35
8327400	WRENCH,PROBE	\$	44.71
8328000	BRACKET,ACM	\$	165.35
8329100	MAINTAINCE SET,CL-SINGEL	\$	30.26
8330900	PLUG, OD3/8"	\$	6.75
8332000	SPIRAL WRAP, 3/8" BLACK	\$	3.85
8332300	CABLE CLAMP, ADJUSTABLE,7.9/10.3DIA	\$	1.35
833311	SCREW 2-56X .125PAN SST PHIL	\$	8.67
8334000	MAINTAINCE SET, CH 2, CL	\$	47.41
8335500	CABLE CLAMP, 12.7MM/5.0MM	\$	1.35
833811	SCREW 8-32X .188SET SST SKT	\$	11.37
834-588	WASHER	\$	22.55
834-591	WASHER	\$	22.55
834-592	GASKET/B301	\$	92.70
834-593	GASKET/B305	\$	92.70
834-598	SEAL RUB BLK W/HOLE	\$	29.10
83514	DIPHENYLCARBAZONE IND 113G	\$	430.01
8355700	BASE, AIR-PUMP HOUSING, NA5600sc	\$	35.65

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

8355800	COVER, AIR-PUMP HOUSING, NA5600sc	\$	29.68
8355900	ASSY, AIR PUMP, NA5600sc	\$	2,331.91
8356000	ASSY, CAL TUBING, NA5600sc	\$	16.57
8356200	SEALING BUBBLE GASKET, D=2MM, L=411MM, NA5600sc	\$	10.41
8356300	CABLE GLAND, PG9 4 MM-8 MM	\$	28.14
8356500	SUPPORT,LEVEL SENSOR, NA5600sc	\$	104.65
8356600	HOUSING,CHECK VALVE, NA5600sc	\$	76.70
8357100	CABLE ASSY, LCD BACKLIGHT, NA5600sc	\$	18.12
8357200	CABLE, FLAT, SIGNAL, 33 PIN, NA5600sc	\$	3.47
8357300	CABLE ASSY, POWER SWITCH, 120-240V, NA5600sc	\$	132.59
8357400	ASSY, DRAIN PUMP W\PLUG, NA5600sc	\$	1,611.14
8357500	CABLE ASSY, POWER BD TO MAIN BD, NA5600sc	\$	46.64
8357600	PLUG 2.5MM SPRING CLAMP 6 POST	\$	14.65
8357800	ASSY, TUBING, 0.5L BOTTLE INNER, NA5600sc	\$	46.25
8357900	PLUG 2.5MM SPRING CLAMP 3 POST	\$	9.06
8358800	CABLE ASSY, LEVEL SENSOR, NA5600sc	\$	160.34
8358900	CABLE GLAND PG7 3MM-6.5MM	\$	28.33
8359300	CONDUCTIVITY ELECTRODE, NA5600sc	\$	225.48
8359900	CABLE , REF ELECTRODE, BLUE, NA5600sc	\$	329.55
8360000	CABLE , SODIUM ELECTRODE, BLACK, NA5600sc	\$	321.84
8360200	TBG, LLDPE, 6 ID x 8 OD MM	\$	54.93
8360600	TUBING, SILICON, 6 ID X 9 OD MM	\$	58.78
8361300	BASE 4 CHANNEL, SAMPLE SEQ, NA5600sc	\$	597.43
8361400	BASE 1 CHANNEL, SAMPLE SEQ, NA5600sc	\$	398.93
8361500	PIN NEEDLE VALVE, SAMPLE SEQ, NA5600sc	\$	141.84
8361600	COVER, SAMPLE SEQ, NA5600sc	\$	83.45
8362500	BRACKET LOCK TONGUE, FRAME VERSION, NA5600sc	\$	69.76
8362900	PIVOT HOUSING, FRAME VERSION, NA5600sc	\$	202.36
8363000	BODY, DRAIN CELL, NA5600sc	\$	23.13
8363100	COVER, DRAIN CELL, NA5600sc	\$	21.58
8363200	ASSY, REAC TUBING, 0.5L BOTTLE,NA5600sc	\$	16.57
8363300	PROTECT DOOR,SOLUTION, NA5600sc	\$	73.23
8363400	ASSY, ENCLOSURE W/FEET, NA5600sc	\$	5,261.26
8363500	STOP, ENCLOSURE VERSION, NA5600sc	\$	44.90
8363600	COVER, MAIN BOARD, HIGH-VOLTAGE, NA5600sc	\$	19.46
8363700	COVER, POWER BOARD, NA5600sc	\$	23.13
8363800	COVER, BOARDS, NA5600sc	\$	160.54
8364100	MOUNTING BRACKET, NA5600sc	\$	213.92
8364200	PUMP, 7604 FF, PK, 12V DC1.5W	\$	2,486.09
8364300	PUMP, 50UL, 12VDC, FKM	\$	4,278.38
8364400	ASSY, CHECK VALVE, NA5600sc	\$	202.36
8364500	VALVE, 3-WAY, 6606 FFKM PK 12VDC 3.4W	\$	1,333.62
8364600	BODY, KCL TANK, NA5600sc	\$	37.19
8364700	CAP, KCL TANK, NA5600sc	\$	22.36
8365000	KIT, UPGRADE 2 CHANNEL TO 4 CHANNEL, NA5600sc	\$	2,813.71
8365300	FRONT COVER, 3 PUMPS, NA5600sc	\$	26.98
8365400	HOUSING, 3 PUMPS, NA5600sc	\$	84.60



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

8365500	REAR COVER, 3 PUMPS, NA5600sc	\$	28.33
8365600	SUPPORTER, NA5600sc	\$	14.07
8365700	LOWER HINGE BRACKET, NA5600sc	\$	177.88
8365800	UPPER HINGE BRACKET, NA5600sc	\$	59.36
8367100	MOUNTING BRACKER, PUMP BOX, NA5600sc	\$	21.39
8367200	BLOCK PIECES, ENCLOSURE VERSION, NA5600sc	\$	28.33
8367300	MOUNTING BRACKET, EXCHANGER HEATER	\$	138.95
8367600	ASSY, 3-WAY VALVE, PH, NA5600sc	\$	1,347.11
8367800	ASSY, LEVEL SENSOR, NA5600sc	\$	302.57
8368400	CABLE GLAND, NPT1/2" 4 MM-9 MM	\$	33.92
8368600	CKT BD ASSY, CONVERSION BD, NA5600sc	\$	797.86
8368700	VERTICAL BRACKET, NA5600sc	\$	375.80
8368900	KIT, HEATER EXCHANGE, NA5600sc	\$	2,312.64
8369100	FOOT, ENCLOSURE, NA5600sc	\$	472.16
8369600	STANDOFF, M4 x 6 x 21 MALE-FEMALE,	\$	4.24
8369800	CABLE ASSY, CONVERSION BD TO CONTROL BD, NA5600sc	\$	200.43
83699	DIPHENYLCARBAZONE RGT PP PK/100	\$	156.27
8369900	CABLE, 7604 PUMP, 12V DC1.5W	\$	44.90
83699-LM	DIPHENYLCARBAZONE RGT PP PK/100	\$	156.27
8370000	CABLE ASSY, CAL PUMP TO CONVERSION BD, NA5600sc	\$	12.14
8370200	ASSY, SEQUENCE VALVE, NA5600sc	\$	647.54
8370300	VALVE, 3-WAY, 6126 EPDM P5 12VDC 3.4W	\$	634.05
8370400	FOOT LANCE, NA5600sc	\$	177.88
8370600	CABLE ASSY, VALVE TO CONTROL BD 2-CH, NA5600sc	\$	75.93
8370800	CKT BD ASSY, AIR-PUMP DRIVING BD, NA5600sc	\$	549.25
8371000	CABLE ASSY, VALVE TO CONTROL BD 4-CH, NA5600sc	\$	142.81
8371200	KIT, K-PUMP NA5600sc	\$	3,218.42
8371300	CABLE ASSY, AIR-PUMP DRIVING BD TO CONTROL BD, NA5600sc	\$	10.41
8371500	KIT, INSTALLATION, 4-CHANNEL, ENCLOSURE, NA5600sc	\$	9,674.54
8371700	SEALING PLUG, SD CARD ACCESS, NA5600sc	\$	37.58
8372000	STOP, SD CARD BD, NA5600sc	\$	35.65
8372600	CKT BD ASSY, OPEN DOOR DETECT BD, NA5600sc	\$	79.21
8372900	PIVOT, ENCLOSURE VERSION, NA5600sc	\$	78.05
8373300	CABLE ASSY, POWER BD TO POWER MODULE, NA5600sc	\$	9.06
8373700	LEVELING FEET BRACKET, M10X1.5, NA5600sc	\$	135.67
8373800	LEVELING FEET, M10X1.5, NA5600sc	\$	31.03
8374000	COVER, MIX MOTOR, NA5600sc	\$	194.65
8374400	CKT BD ASSY, SD CARD BD, NA5600sc	\$	136.25
8375200	ASSY, PUMP BOX REAR COVER, NA5600sc	\$	43.94
8375300	KIT, INSTALLATION, 2-CHANNEL, ENCLOSURE, NA5600sc	\$	8,633.86
8375400	KIT, INSTALLATION, 1-CHANNEL, ENCLOSURE, NA5600sc	\$	8,113.51
8375500	KIT, INSTALLATION, 4-CHANNEL, FRAME, NA5600sc	\$	8,441.14
8375600	KIT, INSTALLATION, 2-CHANNEL, FRAME, NA5600sc	\$	7,381.18
8375700	KIT, INSTALLATION, 1-CHANNEL, FRAME, NA5600sc	\$	6,918.65
8376000	TUBING, PVC 5/16" ID X 7/16" OD	\$	30.45
8376300	ASSY, REAC TUBING, MEA CELL, NA5600sc	\$	452.89
8377500	ADAPTER BRACKET, 924X PANEL VERSION, NA5600sc	\$	759.32



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

8380100	UPPER BRACKET, NA5600sc	\$	164.00
8380600	CABLE ASSY, VALVE TO CONTROL BD 1-CH, NA5600sc	\$	61.28
8380700	CABLE ASSY, STEPPER MOTOR TO CONTROL BD, NA5600sc	\$	51.07
8381100	CONNECTOR, DRAIN TUBE, NA5600sc	\$	26.60
8381300	LATCH BRACKET,NA5600sc	\$	28.33
8381400	LATCH HOLDER, NA5600sc	\$	83.45
8382100	ASSY, KCL TANK ,NA5600sc	\$	49.14
8382600	CATASTROPHIC HOLE, NA5600sc	\$	90.77
8383600	ASSY, SOLUTION DOOR,NA5600sc	\$	56.27
8384200	ANALYTICS PANEL, NA5600sc	\$	1,676.66
8387500	DRAIN CELL, NA5600sc	\$	87.11
8387600	LOCKER, MAIN PANEL, NA5600sc	\$	11.95
8387700	PLATE, REINFORCE, NA5600sc	\$	90.77
8390000	TL23 LCD SERVICE PACK, INTERNAL ONLY	\$	1,358.68
8393700	HOLDER, DISPLAY, 3.5", HQ SERIES	\$	6.55
8393800	GASKET, M12 CONNECTOR, HQ SERIES	\$	6.55
8394800	CKT BD ASSY, LITHIUM BATTERY, HQ SERIES	\$	38.35
8396600	CKT BD ASSY, INTELICAL PWR, 3 CH, HQ SERIES	\$	555.03
8398800	SEAL, USB MICRO, HQ SERIES	\$	6.55
8399800	ASSY, BATTERY DOOR, HQ SERIES	\$	124.11
840-200	Tubing for MNV1, 1 METER	\$	116.21
8405800	SCREW, M3.5 x 12, TORX PAN, PT, SS	\$	6.55
8419800	TBG, LLDPE, 6 ID x 8 OD MM	\$	54.93
841A1001	ASSY.PCB.-PWR SUPPLY,VAR.	\$	2,409.00
841A1001-102	ASSY.PCB.PWR SPLY.VAR.VOL	\$	2,832.98
841S0A1A0N	ASSY,POWER SUPPLY,115VAC	\$	2,698.08
841S0A2A0N	ASSY,POWER SUPPLY,230VAC	\$	2,890.80
842-156	STOPPER CONICAL 10MM	\$	83.06
842-157	STOPPER CONICAL 14MM	\$	83.06
842-158	STOPPER CONICAL 16MM	\$	83.06
8425200	PKG/10, HOLDER, CABLE TIE MOUNTING	\$	26.60
8425700	KIT, NETWORK MODULE, HART, NA5600/9600sc, ROW	\$	2,023.56
8425800	KIT, NETWORK MODULE, HART, NA5600/9600sc, EU	\$	1,921.42
8427600	TUBING, 1/16 ID X 1/8 OD	\$	16.19
8428000	PROGNOSYS NA5600sc/9600sc LICENSE KIT	\$	4,066.39
844G2001	ASSY-PCB,PWR SPLY,TWT,BAS	\$	905.78
844G2001-102	ASSY-PCB,PWR SPLY,TWT,230	\$	1,042.62
844X4X2001	ASY-PCB,PWR SPLY,TWT,BAS-	\$	973.24
845G1001-101	ASSY-PCB,PWR SPLY,TWT,115	\$	1,285.44
845G1001-102	ASSY-PCB,PWR SPLY,TWT,230	\$	1,445.40
847-737	STIRRER ROD 127MM/SAC90	\$	1,229.55
8483400	HOLDER, DISPLAY, 3.2IN, HQ SERIES	\$	6.55
8483500	BEZEL, DISPLAY, 3.2IN, HQ SERIES	\$	6.55
8483800	MICRO SD CARD, 16GB, CLASS 10	\$	39.31
8497000	ASSY, POWER SUPPLY WITH FERRITE	\$	196.57
8499000	KTO SL1000 PPA INSTRUMENT, CHEMKEYS AND PROBES	\$	15,167.06
8499100	LCD, 3.5 QVGA, SL1000	\$	780.52



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

8499300	KTO CHEMKEY, 300 PIECE, FREE CHLORINE	\$	813.04
8499400	KTO CHEMKEY, 300 PIECE, TOTAL CHLORINE	\$	813.04
8505350	KTO: HQ440D, PHC805 PH KIT	\$	8,383.32
8505350-LM	KTO: HQ440D, PHC805 PH KIT	\$	8,383.32
8505400RGT	ff RGT SET, RUGGED FIELD KIT, HQ40D w/5M PROBES (8505400)	\$	433.62
8505450	KTO: HQ440D, PHC735 RED-ROD PH KIT	\$	9,096.38
8505510	nn RUGGED FIELD CASE, MEDIUM	\$	1,129.34
8505511	nn RUGGED FIELD CASE, LARGE	\$	1,389.51
8505600	RUGGED FIELD CASE FOR 30M PROBES	\$	950.11
8505850	KTO: HQ440D, PHC725 RED-ROD PH KIT	\$	8,884.39
8506050	KTO: HQ440D, PHC705 RED-ROD PH KIT	\$	8,884.39
8506250	KTO: HQ440D, PHC108 PUNCTURE PH KIT	\$	8,730.22
8506450	KTO: HQ440D, PHC705A RED-ROD ALKALINE PH KIT	\$	8,884.39
8507000	KTO: SODIUM ISE KIT HQ44D ISENA38101	\$	9,289.10
8507000-LM	KTO: SODIUM ISE KIT HQ44D ISENA38101	\$	9,289.10
8507450	KTO: HQ440D, PHC729 RED-ROD SURFACE PH KIT	\$	9,462.55
8507560	KTO: BASIC LAB pH KIT, HQ411D, PHC201	\$	5,357.62
8507570	KTO: BASIC LAB pH/COND HQ440 PHC2 CDC4	\$	9,828.72
8507600	KTO: AMMONIA ISE KIT, HQD, ISENH3181	\$	10,811.59
8507800	KTO: AMMONIUM ISE KIT, HQD, ISENH4181	\$	11,042.86
8507900	KTO: CHLORIDE ISE KIT, HQD, ISECL181	\$	11,293.39
8508000	KTO: FLUORIDE ISE KIT, HQD, ISEF121	\$	11,293.39
8508100	KTO: NITRATE ISE IT, HQD, ISENO3181	\$	11,293.39
8508200	KTO: FLUORIDE/pH ULTRA KIT, HQ44D	\$	13,490.40
8508400	KTO: AMMONIA/NITRATE ISE, HQ44D	\$	15,051.43
8508500	KTO: BOD LAB SYSTEM, HQ440D, LBOD101	\$	11,833.01
8508600	KTO: pH ULTRA KIT, HQ440D	\$	9,154.20
8508700	KTO: WATER QUAL LAB, HQ440D	\$	13,972.20
8508850	UNIVERSAL PROBE STAND, W/ELEC HOLDER	\$	699.57
8509500	KTO: BASIC LAB pH KIT, HQ411D, PHC101	\$	5,550.34
8509600	KTO: LAB CONDUCTIVITY KIT, HQ430D	\$	7,535.35
8510100	KTO: ADVANCE PH/BOD, HQ440D	\$	13,972.20
8510200	KTO: ADVANCE pH/COND, HQ440D	\$	11,081.40
85199	MANVER 2 PWD PLWS 50ML PK/100	\$	60.35
85199-LM	MANVER 2 PWD PLWS 50ML PK/100	\$	60.35
8522400	POWER SUPPLY EXT, 100-240VAC IN 9V OUT	\$	271.74
8522400BD	POWER SUPPLY EXT, 100-240VAC IN 9V OUT	\$	302.57
8522700	ee PLUG, 2.5MM, 8 CKT	\$	61.86
8522900	PLUG, 3.5MM 9 POST	\$	35.65
85299	CALVER 2 PWD PLWS 50ML PK/100	\$	68.10
85299-LM	CALVER 2 PWD PLWS 50ML PK/100	\$	68.10
85428	FERROVER PWD PLWS 25ML PK/1000	\$	814.84
85428-LM	FERROVER PWD PLWS 25ML PK/1000	\$	814.84
8544/N02/2	Cleaning set (220 V) with air or water for loops based on 9135 + 8350 or 8351 when installed in immersion	\$	1,965.74
8544/N03/2	Chemical cleaning set (220 V) for loops based on 9135 + 8350 or 8351 when installed in immersion	\$	7,381.18



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

8544/N03/3	Chemical cleaning set (24 V) for loops based on 9135 + 8350 or 8351 when installed in immersion	\$	9,327.65
85499	FERROVER PWD PLWS 25ML PK/100	\$	99.90
85499-LM	FERROVER PWD PLWS 25ML PK/100	\$	99.90
8554100	ee PUMP ASSY, CL17sc	\$	637.90
8560200	ee ASSY, COLORIMETER, CL17sc	\$	1,289.30
8560400	TUBING KIT, CL17sc	\$	329.55
8560500	INSTALLATION KIT, CL17sc	\$	655.25
8561500	CABLE ASSEMBLY, POWER, CL17sc	\$	160.92
8562400	ASSY, CHIP, REAGENT DELIVERY, CL17sc	\$	7.52
8562700	BULKHEAD UNION BLACK, PP, 1/4 OD TBG	\$	46.06
8563500	ee FLOW METER, CL17sc	\$	659.10
8564600	PIPE, PVC, 1/2" NPT X 6"	\$	12.33
8565700	INSTALLATION KIT, W/PRESSURE REG, CL17sc	\$	1,364.46
8567900	ee ASSY, FLOW METER, CL17sc	\$	946.26
8568200	CALIBRATION VERIFICATION KIT CL17sc	\$	275.59
8570000	CAMSHAFT, PUMP, FOR CHLORINE ANALYZER CL17sc	\$	74.78
8570300	PUMP BLOCK, REAGENTS, CL17sc	\$	59.55
8570400	PUMP BLOCK, SAMPLE, CL17sc	\$	65.33
8572100	MAIN BOARD, PROGRAMMED, CL17sc	\$	990.58
8572200	MAIN BOARD, PROGRAMMED, CL17Dsc	\$	990.58
8572300	ee TOP LEVEL ASSY, CL17sc, STANDPIPE INSTALL KIT	\$	7,882.25
8572400	ee TOP LEVEL ASSY, CL17sc, PRESS.REG INSTALL KIT	\$	7,882.25
8572500	ee TOP LEVEL ASSY, CL17Dsc, STANDPIPE INSTALL KIT	\$	7,882.25
8572600	ee TOP LEVEL ASSY, CL17Dsc, PRESS.REG INSTALL KIT	\$	7,882.25
8572700	KTO: CL17sc, FREE CHLORINE WITH STANDPIPE	\$	8,055.70
8572700-LM	KTO: CL17sc, FREE CHLORINE WITH STANDPIPE	\$	8,055.70
8572800	KTO: CL17sc, TOTAL CHLORINE WITH STANDPIPE	\$	8,055.70
8572900	KTO: CL17sc, FREE CHLORINE WITH PRESS. REG.	\$	8,055.70
8572900-LM	KTO: CL17sc, FREE CHLORINE WITH PRESS. REG.	\$	8,055.70
8573000	KTO: CL17sc, TOTAL CHLORINE WITH PRESS. REG.	\$	8,055.70
8573100	CELL CLEANING KIT, CL17sc	\$	54.93
8573200	KTO: CL17sc, CALIBRATION/VERIF REFILL KIT	\$	138.95
8574400	ee TOP LEVEL ASSY, CL17sc, STANDPIPE INSTALL KIT, EU	\$	7,882.25
8574500	ee TOP LEVEL ASSY, CL17sc, PRESS.REG INSTALL KIT, EU	\$	7,882.25
8574600	ee TOP LEVEL ASSY, CL17Dsc, STANDPIPE INSTALL KIT, EU	\$	7,882.25
8574700	ee TOP LEVEL ASSY, CL17Dsc, PRESS.REG INSTALL KIT, EU	\$	7,882.25
8575000	TUBING KIT, EU, CL17sc	\$	316.06
8575100	INSTALLATION KIT, EU, CL17sc	\$	628.27
8575200	INSTALLATION KIT, EU, W/PRESSURE REG, CL17sc	\$	1,308.57
8576000	KIT, MODIFIED STANDPIPE, CL17sc	\$	122.38
8576300	AIR BREAK, DRAIN, CL17sc	\$	21.20
8576600	DRAIN KIT, CL17sc	\$	158.61
8606300	KTO: Titration Bundle, Free Chlorine	\$	15,032.16
8606400	KTO: Titration Bundle, Total Chlorine (Forward)	\$	15,089.98
8606500	KTO: Titration Bundle, Total Chlorine (Backward)	\$	15,205.61
8606600	KTO: Titration Bundle, Chlorine Dioxide	\$	15,051.43

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

8606700	KTO: Titration Bundle, Alkalinity, P&T (H2SO4)	\$	11,698.10
8606800	KTO: Titration Bundle, Alkalinity, P&T (HCl)	\$	11,698.10
8606900	KTO: Titration Bundle, Total Hardness, Ca & Mg (ISE)	\$	15,706.68
8607000	KTO: Titration Bundle, Total Hardness (Color)	\$	20,640.31
8610400	90 ML FLASK, SL1000	\$	4.82
8610500	dd RECYCLING PROGRAM, CHEMKEY, QTY 950	\$	254.75
8610501	dd RECYCLING PROGRAM, CHEMKEY, 30 gallon drum	\$	2,222.30
8626101	si ASSY, REAGENT BOTTLE CAP, BLUE	\$	148.01
8626102	si ASSY, REAGENT BOTTLE CAP, GREEN	\$	148.01
8626103	si ASSY, REAGENT BOTTLE CAP, RED	\$	148.01
8626104	si ASSY, REAGENT BOTTLE CAP, PURPLE	\$	148.01
8626105	si ASSY, REAGENT BOTTLE CAP, YELLOW	\$	148.01
8626106	si ASSY, REAGENT BOTTLE CAP, BLACK	\$	148.01
8626107	si ASSY, REAGENT BOTTLE CAP, ORANGE	\$	148.01
86261A	SERVICE ASSY SET, BOTTLE/CAP, 5500sc SILICA, 5500/9611sc HR PHOS	\$	595.50
86261B	SERVICE ASSY SET, BOTTLE/CAP, 5500/9611sc LR PHOS	\$	736.19
86261C	SERVICE ASSY SET, BOTTLE/CAP, 9610sc SILICA	\$	736.19
8626200	SENSOR, FREE CHLORINE CLF10 sc STAINLESS TIP	\$	5,550.34
8626200-LM	SENSOR, FREE CHLORINE CLF10 sc STAINLESS TIP	\$	5,550.34
8627000	ff KTO: TU5300sc, Sys Chk RFID Flow ACM w/sc200 2 CH	\$	14,685.26
8627800	HINGE BRACKET, ANALYTICS PANEL, PROCESS ANALYZER	\$	233.19
8628400	O-RING, -015, EPDM, 70 DURO	\$	7.90
8628700	SERVICE KIT, AIR MANIFOLD, PROCESS ANALYZER	\$	657.18
8628900	SENSOR, TOTAL CHLORINE CLT10 sc Stainless Tip	\$	5,280.53
8633100	KIT, MEMBRANE REPLACEMENT, CLF10sc, CLT10sc, STAINLESS TIP	\$	300.64
8635200	ASSY., PACK, CHEMKEY, 25 PIECE, PERACETIC ACID	\$	64.32
8636000	ASSY., PACK, CHEMKEY, 25 PIECE, DISSOLVED IRON	\$	71.19
8636100	ASSY, PACK, CHEMKEY, 25 PIECE, TOTAL ALKALINITY, HIGH RANGE	\$	61.96
8636200	ASSY, PACK, CHEMKEY, 25 PIECE, TOTAL ALKALINITY, LOW RANGE	\$	61.96
8636300	ASSY., PACK, CHEMKEY, 25 PIECE, HARDNESS, HIGH RANGE	\$	76.05
8636400	ASSY., PACK, CHEMKEY, 25 PIECE, HARDNESS, LOW RANGE	\$	61.96
8636500	ASSY., PACK, CHEMKEY. 25 PIECE, ORTHOPHOSPHATE, HIGH RANGE	\$	117.61
8636600	ASSY., PACK, CHEMKEY. 25 PIECE, ORTHOPHOSPHATE, LOW RANGE	\$	117.61
8699500	CLIP, EU POWER SUPPLY, R SERIES	\$	6.17
8701900	ASSY, ENCLOSURE, TOP, 4XXX, HQ SERIES	\$	187.90
8701901	ASSY, ENCLOSURE, TOP, 1110, HQ SERIES	\$	183.85
8701902	ASSY, ENCLOSURE, TOP, 1140, HQ SERIES	\$	185.78
8701903	ASSY, ENCLOSURE, TOP, 1130, HQ SERIES	\$	188.48
8701904	ASSY, ENCLOSURE, TOP, 2X00, HQ SERIES	\$	174.22
8702000	FASTENER AND SEAL KIT, HQ SERIES	\$	45.67
8702401	ASSY, ENCLOSURE, BOTTOM, 1 PORT, 21XX/11XX, HQ SERIES	\$	160.34



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

8702402	ASSY, ENCLOSURE, BOTTOM, 2 PORT, 22XX, HQ SERIES	\$	160.34
8702403	ASSY, ENCLOSURE, BOTTOM, 1 PORT, 4100, HQ SERIES	\$	164.00
8702404	ASSY, ENCLOSURE, BOTTOM, 2 PORT, 4200, HQ SERIES	\$	164.00
8702405	ASSY, ENCLOSURE, BOTTOM, 3 PORT, 4300, HQ SERIES	\$	164.00
8702900	SCREWDRIVER, REVERSIBLE, HACH LOGO	\$	13.49
87299	PHENOL REAGENT PWD PLWS PK/100	\$	252.95
87299-LM	PHENOL REAGENT PWD PLWS PK/100	\$	252.95
875011	SCREW 8-32X .438PAN SST PHIL	\$	2.51
8772800	CKT BD ASSY, DISPLAY BOARD, sc200	\$	1,985.02
8773200	POWER SUPPLY, AC/DC, W/RELAY BD, sc200	\$	846.04
8782200	PLUG 5.08MM SPRING CLAMP 3 POST	\$	42.98
8782600	aa CKT BD ASSY,4-20mA INPUT SENSOR MODULE, sc200	\$	423.98
8785600	PLUG 5.08MM SPRING CLAMP 4 POST	\$	47.79
8785800	PLUG 5.08MM SPRING CLAMP 8 POST	\$	97.71
8786900	KEYPAD, sc200 CONTROLLER	\$	74.78
8787100	COVER, SD CARD, sc200	\$	21.58
8787300	HINGE PIN, FRONT DOOR, sc200	\$	13.10
8787400	ENCLOSURE, REAR, sc200	\$	310.28
8787900	BLANKING PLATE, UPPER, sc200	\$	43.75
8788000	BLANKING PLATE, LWR, SENSOR MODULE, sc200	\$	31.99
8788100	BLANKING PLATE, LWR, NETWORK, sc200	\$	53.58
8789200	CABLE, FLAT RIBBON, 50 CKT .5MM PITCH, 6"	\$	46.25
8789900	PLUG, 3.5MM 12CKT	\$	34.11
8791500	ASSY., PACK, CHEMKEY, 25 PIECE, F&T AMMONIA	\$	496.86
8791900	KTO CHEMKEY, 25 pc MONO-CHLOR. AND 25 pc F&T AMMONIA	\$	540.22
8796100	2FF SIM Card for SC1500 - MSM	\$	17.73
8796200	3FF SIM Card for SC4200 - MSM	\$	22.36
8804300	CABLE ASSY, 4 COND W/5 CKT RECEPTACLE ONE END	\$	49.14
8804900	PLUG, 5.08MM SPRING CLAMP 2 POST	\$	10.21
8806200	KIT, MOUNTING HARDWARE, sc200, COMPLETE	\$	639.83
8806260	MOUNTING HARDWARE, SCREWS& WASHERS ONLY, sc200	\$	72.46
8809200	COVER ASSY, UV, sc200	\$	585.87
8809600	PLUG, 3.5MM, 6 CKT	\$	20.62
8812600	CKT BD ASSY, SD BOARD, sc200	\$	124.69
8814200K	KTO: CT 10000 PPM CARBON DIOXIDE CL 103L	\$	690.18
8814500	SCREW CAPTIVE 8-32X .375 PAN SST PHIL	\$	2.12
8814900	GASKET SET, sc200	\$	46.06
8828000	WQMP, AUTO CLEAN, SC1000	\$	49,933.75
88449	HYDROCHLORIC ACID, 6.0N 500ML	\$	89.78
88449-LM	HYDROCHLORIC ACID, 6.0N 500ML	\$	89.78
886-213	WINDOW RED/ABU9X	\$	110.62
8866800	WDMP , SIX SENSORS, TU5300sc W/O ACM & CL17sc	\$	56,948.76
8867000	WDMP, SIX SENSORS, TU5300sc WITH ACM & CLT10sc	\$	62,537.64
8867200	WDMP, SIX SENSORS, TU5300SC WITH ACM & CL17sc	\$	56,948.76
8867300	WDMP, SIX SENSORS, TU5300sc WITH ACM & CLF10sc	\$	62,537.64
8867400	WDMP, SIX SENSORS, TU5300sc W/O ACM & CLF10sc	\$	62,537.64
8867600	WDMP, SIX SENSORS, TU5300sc W/O ACM & CLT10sc	\$	62,537.64

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

8888001	HACH COLITAG, 16-HR PA, COMPARITOR	\$	49.14
8888006	COLITAG STERILE SAMPLE BOTTLE 100/PK	\$	217.77
89149	SULFATE STD SOLN, 100MG/L 500ML	\$	80.22
891901	TAG, METAL, SPECIAL FOR CUSTOMER	\$	124.30
89568	BUFFER, PH4.50 PP 50ML PK/25	\$	66.13
89568-LM	BUFFER, PH4.50 PP 50ML PK/25	\$	66.13
8977600	CKT BD ASSY, INTELICAL PWR, 1 CH, HQ SERIES	\$	191.18
8977900	DUST PLUG, USB MINI, HQ SERIES	\$	6.55
8978200	CKT BD ASSY, INTELICAL PWR, 2 CH, HQ SERIES	\$	333.41
8978300	CLIP, UK POWER SUPPLY, R SERIES	\$	8.29
8978400	CLIP, ARGENTINA POWER SUPPLY, R SERIES	\$	11.18
8978500	CLIP, BRAZIL POWER SUPPLY, R SERIES	\$	7.52
8978600	CLIP, CHINA POWER SUPPLY, R SERIES	\$	6.17
8978700	CLIP, INDIA POWER SUPPLY, R SERIES	\$	10.41
8978800	CLIP, KOREA POWER SUPPLY, R SERIES	\$	7.52
8978900	CLIP, AUSTRALIA POWER SUPPLY, R SERIES	\$	9.83
8980000	ee POWER SUPPLY, 5V, 2A, 100-240V, R SERIES	\$	75.93
8980100	CLIP, US POWER SUPPLY, R SERIES	\$	3.85
89868	BUFFER, PH8.30 PP 50ML PK/25	\$	67.93
89868-LM	BUFFER, PH8.30 PP 50ML PK/25	\$	67.93
8H1112	NUT,STAINLESS STEEL,#4-40	\$	5.40
8H1114	NUT,STAINLESS STEEL,#8-32	\$	11.37
8H1116	NUT,STAINLESS STEEL,1/4-2	\$	1.93
8H1214-038	SCREW, S.S., #8-32 BY 3/8	\$	5.40
8H1214-075	SCREW, S.S., #8-32 BY 3/4	\$	11.37
8H1215-038	SCREW, S.S., #10-32 BY 3/	\$	11.37
8H1293-031	SCREW, SOCKET HD CAP, 6-32	\$	2.12
8H1312	WASHER, LOCK(IN), #4, STAI	\$	5.40
8H1314	WASHER, LOCK(IN), #8, STAI	\$	11.37
8H1326	WASHER, LOCK(EX), 1/4 S.S	\$	2.12
8H1335	WASHER, SPLIT, #10, STAI	\$	11.37
8H1336	WASHER, SPLIT, 1/4", 316SS	\$	6.75
8H1337	WASHER, SPLIT, 5/16", 316SS	\$	2.51
8H1346	WASHER, FLAT, 1/4", 316SS	\$	5.40
8H1347	WASHER, FLAT, 5/16", 316SS	\$	2.51
8H1351	SCREW, SOCKET HEAD CAP	\$	8.48
8P2117-100	PLUG, PIPE, 1/8" PFA	\$	129.32
9005300	VALVE, 3/8 NPTF, SPEED FIT TO MALE NPTF	\$	100.41
9009100	KTO: AMTAX SC , 115-230V .02-5 1 CHANNEL CONT	\$	60,359.90
9009200	KTO: AMTAX SC 115-230V, .02-5 2 CHANNEL CONT	\$	62,980.90
9012700	SENSOR MODULE, FLOW, sc200, HACH	\$	2,524.63
9012700X	SENSOR MODULE, FLOW, sc200, HACH	\$	2,524.63
9012705	vv sc200 Durchfluss Sensor Modul	\$	2,524.63
9012800X	SENSOR MODULE, 4-20 mA INPUT, sc200, HACH	\$	1,570.67
9012805	vv sc200 4-20mA Input Modul	\$	1,480.09
9012900	SENSOR MODULE, pH/ORP/DO, sc200, HACH	\$	1,480.09



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

	KIT, SENSOR MODULE, pH/ORP/DO, sc200, HACH, CONFORMAL		
9012900-CFC	COATED	\$	2,158.46
9012900X	SENSOR MODULE, pH/ORP/DO, sc200, HACH	\$	1,570.67
9012905	vv sc200 pH/ORP/DO Modul	\$	1,480.09
	KIT, SENSOR MODULE, CONDUCTIVITY, sc200, HACH, CONFORMAL		
9013000-CFC	COATED	\$	2,158.46
9013000X	SENSOR MODULE, CONDUCTIVITY, sc200, HACH	\$	1,570.67
9013005	sc200 Leitfähigkeit Modul	\$	1,480.09
9013200	NETWORK MODULE, MODBUS, sc200, HACH	\$	2,062.10
9013205	vv sc200 Modbus RS232/RS485 Modul	\$	2,062.10
9015100	PACKAGING ASSEMBLY, sc200	\$	40.09
9020000	ASSY, PROBE, LDO MODEL 2, HACH	\$	7,766.62
9020000-030	ASSY, LDO MODEL 2, 30M CABLE, HACH	\$	10,426.15
9020000-060	ASSY, LDO MODEL 2, 60M CABLE, HACH	\$	12,700.25
9020000-AQ	ASSY, LDO MODEL 2 AQ, HACH	\$	3,738.77
9020000-AQS	ASSY, LDO MODEL 2 AQ, SEAWATER, HACH	\$	2,659.54
9020000-C1D2	CLASS 1 DIV 2 VERSION OF LDO	\$	8,286.96
9020000DK	ASSY, PROBE, LDO MODEL 2, HACH	\$	6,918.65
9020000-SW	LDO MODEL 2 PROBE, 10M, SEA WATER	\$	11,042.86
9020000-UPGRADE	KTO: LDO MODEL 2, 10M HACH & 9253400 ADAPTER	\$	7,882.25
9021100	KIT, SENSOR CAP REPLACEMENT, LDO MODEL 2, HACH	\$	813.28
9021100DK	KIT, SENSOR CAP REPLACEMENT, LDO MODEL 2, HACH	\$	747.75
9021150	KIT: SENSOR CAP REPLACEMENT, LDO MODEL 2, LANGE	\$	601.29
904-592	CHAMBER ABSORPTION COMPLETE	\$	316.06
904-613	ADAPTER/940540	\$	711.14
904-614	ADAPTER/940496	\$	1,167.88
904-677	TITRATION VESSEL 50PC/SAM90	\$	406.64
905-115	CLASPS 3X25PCS/SAC90	\$	127.39
90700	SPOON, MEASURING 0.5G	\$	16.19
908424	FUSE 1.5 AMP SLO-BLO	\$	24.48
9092700	KNOB, FLUTED 2-3/8 DIA, 3/8-16 X 1 STUD	\$	96.94
9092800	KNOB, FLUTED, 2-3/8 DIA, 1/2-13 X 2 STUD	\$	121.41
9094900	BATTERY, LI POLY 7.4V 5.0Ahr W/PCBA	\$	462.53
9097400	CABLE, USB A 4 PIN MALE TO MINI B 5 PIN MALE USB 2.0.	\$	113.13
9098400	si CKT BD ASSY, MAIN BD, SL1000	\$	2,235.55
9098800	si CKT BD ASSY, CONNECTOR BD, SL1000	\$	597.43
90A3030	ASSY,KIT,ENCL.MTG	\$	254.39
	Automatic flow control device for liquids (like on 8810 External-		
910=015=043	Dilution)	\$	3,199.15
9108600	si CLAMP HALF, DIGITAL, SST, INSERT.MT	\$	437.47
9108700	CLAMP HALF, ANALOG, SST, INSERT.MT	\$	466.38
9108800	CLAMP HALF, DIGITAL, PVC, INSERT.MT	\$	162.46
9108900	CLAMP HALF, ANALOG, PVC, INSERT.MT	\$	162.46
9109100	SCREW M6 X 1.0 25MM CAP SST SKT	\$	9.06
9109200	NUT, M6 HEX NYLOCK SST	\$	5.59
9110800	si BASE, PIVOT, POLE MOUNT	\$	190.21
9110900	si CLAMP, PIVOT, POLE MOUNT	\$	283.30

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

912-014	Stirring Magnet	\$	188.67
9140600	CONN, 3.5MM CTR PLUG 8CKT RT ANGLE	\$	33.73
9141500	PLUG, 3.5MM, 2 CKT W/SCREW FLANGE	\$	6.55
9150200	si CKT BOARD ASSY, OI CHLORINE	\$	1,852.04
9158900	FLOW CELL, MACHINED, AMPEROMETRIC CHLORINE	\$	1,805.79
9159000	NEEDLE VALVE, FLOW CELL, AMPEROMETRIC CHLORINE	\$	240.90
9159100	FLOAT, FLOW CELL, AMPEROMETRIC CHLORINE	\$	112.93
9159200	FITTING, SPEED FIT, 1/4 BSPT X 1/4 OD	\$	20.62
9159300	FITTING, SPEED FIT, 3/8 BSPT X 3/8 OD	\$	36.04
9159400	PLUG, 1/8 NPTF	\$	10.41
9159500	PLUG, 1/2 NPTF	\$	10.41
9159600	O-RING #2-108 , VITON	\$	6.17
9159700	BODY, FLOW CELL, MOLDED, pH/ORP SENSOR	\$	61.86
9159900	KIT, CONDITIONING, FLOW SAMPLE, PRESSURE	\$	1,453.11
9160000	ENCLOSURE, MODIFIED, DIGITAL GATEWAY	\$	643.68
9160100	GLAND, DOUBLE HOLE STRAIN RELIEF	\$	10.41
9160200	KIT, MEMBRANE REPLACEMENT, CLF10sc PLASTIC TIP	\$	358.46
9160300	si ASSY, CUSTOM CABLE, M-12, 5 PIN, 1 METER LENGTH	\$	260.17
9160600	ELECTROLYTE, CLF10sc (FREE)	\$	296.31
9160600-LM	ELECTROLYTE, CLF10sc (FREE)	\$	296.31
9160800	INDUCTIVE SENSOR, M18 , 3 WIRE DC PNP, 8mm SENSING RANGE	\$	566.60
9160900	CABLE ASSY, 4 PIN, SHIELDED	\$	138.95
9173900	KIT, NETWORK MODULE, PROFIBUS, sc200	\$	2,100.65
9177800	KIT, MTG HARDWARE SCREWS ONLY, SC200	\$	190.60
9177900	KIT, MOUNTING BRACKET INSERTS, SC200 MOUNTING	\$	202.36
9178000	KIT, STRAIN RELIEF	\$	56.85
9178400	Profibus M12 T-splitter	\$	1,113.92
9178500	Profibus M12 connector kit	\$	1,653.54
9179700	si POWER CORD, MODIFIED, NORTH AMERICAN	\$	152.06
9179800	POWER CORD, MODIFIED, E.U.	\$	143.38
9180000	ASSY, FLOW CELL, AMPEROMETRIC CHLORINE, ENGLISH	\$	3,314.78
9180100	ASSY, FLOW CELL, 1" NPT SENSOR, ENGLISH	\$	1,522.49
9180200	ASSY, FLOW CELL, 3/4" NPT SENSOR, ENGLISH	\$	1,522.49
9180300	ASSY, FLOW CELL, 1" NPT SENSOR, METRIC	\$	1,520.56
9180400	ASSY, FLOW CELL, 3/4" NPT SENSOR, METRIC	\$	1,522.49
9180500	DIGITAL GATEWAY, AMPEROMETRIC CHLORINE	\$	3,970.03
9180700	PANEL ASSY, AMPEROMETRIC CHLORINE, METRIC	\$	10,406.88
9180705	PANEL ASSY, AMPEROMETRIC CHLORINE, METRIC, EU	\$	10,406.88
9180900	KIT, MEMBRANE REPLACEMENT, CLT10sc	\$	358.46
9181400	TOTAL ELECTROLYTE FOR CLT10sc	\$	296.31
9181500	PHD, RYTON, pH-WIDE RANGE, 6 FOOT CABLE w/ADAPTER	\$	3,122.06
9181505	PHD, RYTON, pH-WIDE RANGE, 6 FOOT CABLE w/ADAPTER HACH/LANGE	\$	3,122.06
9181600	si pH COMBINATION SENSOR	\$	1,085.01
9181605	si pH COMBINATION SENSOR HACH/LANGE	\$	1,085.01
9181700	SPLIT RING, AMPEROMETRIC CHLORINE SENSOR	\$	26.21



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

9181800	O-RING, 25MM ID X 2.5MM THK, SILICONE	\$	44.90
9182700	ASSY, INLET BLOCK W/ CONDUCTIVITY	\$	2,409.00
919002	ff PIPET, SEROLOGICAL 1.00X0.01ML	\$	36.81
9196100	SEALING HUB, AMPERIOMETRIC CHLORINE	\$	102.91
9196200	SENSOR HOLDER, AMPEROMETRIC CHLORINE	\$	108.31
9196300	O-RING, #2-123, VITON	\$	6.17
9196500	FITTING, SPEED FIT, 3/8 BSPT X 10 MM OD TUBE	\$	46.83
9196900	FITTING, ELBOW, SPEED FIT, 1/2 BSPT X 1/2 OD	\$	82.48
9197000	FITTING, ELBOW, SPEED FIT, 1/2 BSPT x 12 MM OD	\$	102.91
9198100	TBG, LLDPE, 1/4" ID X 3/8" OD X 1/16" WALL	\$	6.17
9198200	TBG, LLDPE, 3/8" ID X 1/2" OD X 1/16" WALL	\$	6.17
9198300	ASSY, FLOW CELL, AMPEROMETRIC CHLORINE, METRIC	\$	3,314.78
9198400	KIT, CHLORINE FLOW CELL ASSY, ENGLISH	\$	3,025.70
9198500	KIT, CHLORINE FLOW CELL ASSY, METRIC	\$	2,563.18
9198600	KIT, NEEDLE VALVE, AMPEROMETRIC CHLORINE SENSOR	\$	283.30
9200000	KIT, POWER SUPPLY AC/DC, sc200 REPLACEMENT	\$	1,880.95
9200100	KIT, POWER SUPPLY, DC/DC, sc200 REPLACEMENT	\$	1,900.22
9200200	KIT, BARRIER, HIGH VOLTAGE, REAR, sc200 REPLACEMENT	\$	129.70
9200300	KIT, BARRIER, HIGH VOLTAGE, FRONT, sc200 REPLACEMENT	\$	182.12
9200400	KIT, BARRIER, HIGH VOLTAGE, FRONT, PROFIBUS, sc200 REPLACEMENT	\$	198.50
9200500	KIT, COVER, DISPLAY BD, sc200 REPLACEMENT	\$	134.71
9200600	KIT, DISPLAY BOARD, sc200 REPLACEMENT	\$	2,274.10
9200700	KIT, LCD DISPLAY, sc200 REPLACEMENT	\$	474.09
9200800	KIT, DOOR ASSY, BASIC FRONT, sc200 REPLACEMENT	\$	356.53
9200900	KIT, COVER ASSY, SD CARD, sc200 REPLACEMENT	\$	178.27
9200G	Portable D.O. Meter	\$	7,342.63
9201000	KIT, DIGITAL SENSOR, CONNECTOR,	\$	468.31
9201100	KIT, SCREWS, sc200 REPLACEMENT	\$	133.36
9201200	KIT, BLANKING PLATES, sc200 REPLACEMENT	\$	139.14
9201300	KIT, DOOR ASSY, COMPLETE, FRONT, sc200, HACH REPLACEMENT	\$	3,372.60
9201400	KIT, DOOR ASSY, COMPLETE, FRONT, sc200, HACH-LANGE REPLACEMENT	\$	3,430.42
9201900	SCREW, M4X0.7mmX6MM LG, PAN, #1 PHIL, W/INT. TOOTH LOCK WASHER	\$	-
9202200	KIT, HOLE PLUGS, sc200 REPLACEMENT	\$	210.06
9202400	BEAD, FERRITE SHIELD 14.27mmO.D.X 6.35mmI.D.X10mm LG	\$	25.25
9202900	KIT, POWER CORD, sc200 NORTH AMERICA ACCESSORIES	\$	227.41
9203000	KIT, POWER CORD, sc200 EUROPEAN ACCESSORIES	\$	227.41
92101	rr MERCURIC NITRATE 2.256N, CRTG	\$	80.03
92101-LM	rr MERCURIC NITRATE 2.256N, CRTG	\$	80.03
921145	GREASE, HIGH VACUUM	\$	186.55
921252	SCREWDRIVER PLASTIC GC8282 E	\$	21.78
9218100	MEMORY CARD, SECURE DIGITAL HIGH CAPACITY 4GB	\$	210.06
9218200	READER, SD CARD HIGH CAPACITY	\$	158.03
9219300	KIT, FUSE REPLACEMENT 4A, SC200 DC/DC	\$	65.52



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

9219400	KIT, FUSE REPLACEMENT 1.6A, SC200 AC/DC	\$	60.90
9220600	COVER ASSY, UV &SUN SHIELD, sc200	\$	714.99
9221200	ECLOX KIT INSERT	\$	98.87
9221700	DOOR, UV COVER, REPLACEMENT, sc200	\$	456.75
9222400	KIT, CONDUIT WIRING, sc200	\$	38.16
9222900	KIT, DOOR ASSY, BASIC, sc200 REPLACEMENT, HACH-LANGE	\$	449.04
9245300	HOUSING, TOP, HQd BENCHTOP	\$	46.06
9245500	COVER, BATTERY, HQd BENCHTOP	\$	20.62
9245600	DISPLAY WINDOW, MULTI, HQ440D	\$	53.77
9245700	DISPLAY, 440 X 160 DOTS, HQd BENCHTOP	\$	316.06
9246400	si FLEX CIRCUIT ASSY, HQ440d BENCHTOP	\$	805.57
9246500	si FLEX CIRCUIT ASSY, HQ411d/ HQ430d BENCHTOP	\$	747.75
9247000	DISPLAY WINDOW, FLEXI, HQ430d	\$	49.91
9247100	DISPLAY WINDOW, pH/Mv, HQ411d	\$	49.91
9248100	KIT, BATTERY HOLDER ASSY, HQd BENCHTOP	\$	47.79
9252800	ASSEMBLY, AIRBLAST WASHER HEAD	\$	423.98
9253000	KTO: POLE MOUNT, 1" NPT SENSOR	\$	1,892.51
9253100	KTO: FLOAT MOUNT, 1" NPT SENSOR	\$	2,601.72
9253200	si ADAPTER, 1 1/2" FNPT TO 1" FNPT	\$	132.78
9253300	KTO: POLE, 1" NPT SENSOR	\$	437.47
9253400	MOUNTING CONVERSION ADAPTER, LDO MODEL 1 TO LDO MODEL 2	\$	132.78
9253500	KIT: AIR BLAST, LDO MODEL 2	\$	1,801.93
9253800	ASSY, HEAD, LDO MODEL 2	\$	389.29
9253900	ASSY, SENSOR GUARD, LDO MODEL 2	\$	661.03
9257000	ASSEMBLY, UNION MOUNT, 1" NPT LDO	\$	1,480.09
9257500	KIT, FLOAT ASSEMBLY HARDWARE	\$	828.70
9257700	KIT, AIR BLAST HARDWARE COMPONENTS	\$	807.50
9257800	KIT, AIR BLAST HOSE ASSEMBLY	\$	620.56
9260300	aa BOTTLE, SOAKER, 0.5 oz W/GROMMET & CAP, REMINDER LABEL	\$	25.63
9261000	INDOOR LIGHT, COLOR DISC 2PC, IRON-FERROVER HI-RANGE	\$	208.14
9261100	COLOR DISC 2PC, PHENOL RED	\$	335.33
9261200	COLOR DISC 2PC, BROMTHYMOL BLUE	\$	335.33
9261300	COLOR DISC 2PC, AMMONIA NITRO-SALICYLATE	\$	210.06
9261400	COLOR DISC 2PC, NITRATE NITRO-NITRAVER V	\$	208.14
9261500	COLOR DISC 2PC, CHLORINE-DPD NALCO	\$	165.93
9261600	COLOR DISC 2PC, CHLORINE-DPD HIGH RANGE	\$	208.14
9261900	COLOR DISC 2PC, OZONE-DPD	\$	204.28
9262000	COLOR DISC 2PC, SILICA HIGH RANGE	\$	208.14
9262100	COLOR DISC 2PC, PHOSPHATE-PHOSVER III	\$	208.14
9262200	COLOR DISC 2PC, MLYBDTE MLYBDNM	\$	198.50
9262300	COLOR DISC 2PC, NITRITE NITROGEN LW-RG	\$	215.85
9262400	INDOOR LIGHT, COLOR DISC 2PC, IRON-FERROVER MD-RNG	\$	208.14
9262500	COLOR DISC 2PC, APHA-PT-CO STANDARD	\$	208.14
9262600	COLOR DISC 2PC, AMMONIA NITRO-NESSLER	\$	210.06
9262700	COLOR DISC 2PC, MANGANESE PERIODATE	\$	208.14



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

9262800	COLOR DISC 2PC, SULFIDE-METHYLNE BLUE	\$	103.30
9262900	COLOR DISC 2PC, NITRITE HIGH RANGE	\$	215.85
9263100	COLOR DISC 2PC, PHENOL HIGH RANGE	\$	210.06
9263200	COLOR DISC 2PC, TANNIN-LIGNIN	\$	210.06
9263300	COLOR DISC 2PC, COPPER-CUVER I	\$	210.06
9263400	COLOR DISC 2PC, MOLYBDATE MOLYBDENUM HR	\$	208.14
9263500	COLOR DISC 2PC, THYMOL BLUE pH	\$	210.06
9263600	COLOR DISC 2PC, IRON-FERROVER LR IL	\$	103.30
9263700	COLOR DISC 2PC, IRN-FERROVER HI-RNG OL	\$	103.30
9263800	COLOR DISC 2PC, IRN-FERROVER MR OL	\$	103.30
9263900	COLOR DISC 2PC, IRON-FERROVR LR OL	\$	103.30
9264300	COLOR DISC 2PC, TPTZ IRON HR	\$	103.30
9264400	COLOR DISC 2PC, IRON FERROZINE LR IL	\$	198.50
9264500	COLOR DISC 2PC, TOLYLTRIAZOLE	\$	202.36
9264600	COLOR DISC 2PC, NITRATE NITRO-NITRAV VI IL	\$	198.50
9264700	COLOR DISC 2PC, SLFIDE-METH BLUE LR	\$	103.30
9265300	COLOR DISC 2PC, CHROMATE-CHROMIUM	\$	103.30
9265400	COLOR DISC 2PC, TPTZ IRON LR	\$	103.30
9265500	COLOR DISC 2PC, ACCUVAC OZONE HIGH RNG	\$	198.50
9265700	COLOR DISC 2PC, DETERGENTS	\$	192.33
9265900	COLOR DISC 2PC, OXYGEN SCAVENGER	\$	215.85
9266100	COLOR DISC 2PC, ACCUVAC DISSLVD OXY HR	\$	206.21
9266200	COLOR DISC 2PC, MANGANESE-PAN 0-0.7 MG/L	\$	225.48
9266300	COLOR DISC 2PC, IRON FERROZINE HR IL	\$	198.50
9266600	COLOR DISC 2PC, PHOSPHATE-STANNAVER	\$	192.33
9266800	COLOR DISC 2PC, SILICA LOW RANGE	\$	215.85
9266900	COLOR DISC 2PC, BENZOTRIAZOLE	\$	202.36
9267000	COLOR DISC 2PC, GLUTARALDEHYDE	\$	202.36
9267100	COLOR DISC 2PC, CYANIDE-CYANIVER	\$	196.57
9267300	COLOR DISC 2PC, HYDRAVER	\$	206.21
9267400	COLOR DISC 2PC, CHROMATE-CHROMIUM OL	\$	103.30
9267700	COLOR DISC 2PC, CHLORINE DIOXIDE-DPD	\$	192.33
9267900	COLOR DISC 2PC, NITRATE NITRO-NITRAV VI OL	\$	198.50
9270000	COLOR DISC, FERROVER IRON 0-10 OL	\$	100.41
92728	KTO: FERROVER PWD PLWS 5ML PK/1000	\$	569.13
92728-LM	KTO: FERROVER PWD PLWS 5ML PK/1000	\$	569.13
92746	FERROVER PWD PLWS 5ML PK/50	\$	36.86
9275500	Profibus M12 socket	\$	784.37
9278300	KIT, sc100/sc200 BARRIERS, HOAB	\$	57.24
92799	FERROVER PWD PLWS 5ML PK/100	\$	68.30
92799-LM	FERROVER PWD PLWS 5ML PK/100	\$	68.30
92899	MANVER 2 PWD PILLOWS PK/100	\$	70.63
92899-LM	MANVER 2 PWD PILLOWS PK/100	\$	70.63
93153	SULFURIC ACID STD 10N, 1000ML	\$	111.83
93153-LM	SULFURIC ACID STD 10N, 1000ML	\$	111.83
9321000	KIT, SMART PROBE ADAPTOR,5500sc/961x	\$	1,123.56
9325900	FUSE .125A 250V SLOW BLOW 5X20MM GLASS	\$	5.40

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

9328100	KIT, NETWORK MODULE, HART, sc200, HACH	\$ 2,081.38
9328100-CFC	KIT, NETWORK MODULE, HART, sc200, HACH, CONFORMAL COATED	\$ 2,389.73
9328100X	NETWORK MODULE, HART, sc200, HACH	\$ 2,254.82
9334600	KIT, NETWORK MODULE, POWERED 4-20mA OUTPUT, sc200, HACH	\$ 1,298.93
9334600X	NETWORK MODULE, POWERED 4-20mA OUTPUT, sc200, HACH	\$ 1,289.30
934207	CABLE, 3 COND 22 AWG STRD	\$ 4.43
9345200	si ff DUST CAP, GLOVE ASSY, HQd PORTABLE	\$ 61.86
9359100	oo KIT, MAINTENANCE, 1735 TRO ANALYZER	\$ 799.79
93600	CLIPPERS, FOR SMALL PWD PLWS	\$ 12.72
9374200	CABLE ASSY., POWER ADAPTOR, AUTOMOTIVE 12 VOLT	\$ 150.90
9377700	ASSY, TRAY, SL1000	\$ 356.53
9385100	DR 900 COLORIMETER, HACH	\$ 5,357.62
9385100NB	DR 900 COLORIMETER, HACH - NO BATTERY	\$ 5,068.54
9385160	DR 900 COLORIMETER, HACH INST ONLY	\$ 3,565.32
9385200	DR 900 Robust Portable Datalogging	\$ 4,625.28
9386300	CAP, INSTRUMENT DR 900	\$ 89.61
9386400	COVER, BATTERY, DR 900	\$ 88.46
9389200	CKT BD ASSY, MiniUSB CONNECTOR BD, DR 900	\$ 262.10
9389500	ASSY, BATTERY HOLDER DR 900	\$ 49.34
9390500	CAP, REPLACEMENT, SERVICE DR 900 HACH	\$ 116.98
9390800	BATTERY COVER,REPLACEMENT, SERVICE	\$ 142.03
9390900	TOP HOUSING ASSY, REPLACEMENT, SERVICE DR 900	\$ 551.18
9391700	ASSY, SAMPLE PREHEATER, AC 120V	\$ 3,237.70
9391800	ASSY, SAMPLE PREHEATER, DC 24V	\$ 3,334.06
9392900	si CASE ASSY,CEL ADV DRINKING	\$ 992.51
9393400	si CASE ASSY, MEL POTABLE WATER KIT	\$ 988.65
9394900	REPLACEMENT TRAY SET OF FOUR SL1000	\$ 1,081.16
9395100	FITTING, 1/4-28 MALE FLAT BOTTOM, TO DOUBLE BARB	\$ 4.05
9397800	1735 TRO Analyzer Reagent Set	\$ 236.68
9400500	DPD COMPOUND, BW 1735 TRO 10.88 (+/-0.5) g	\$ 77.50
9401200	TRO SULFURIC ACID STD 5.25N, 15ML	\$ 13.73
9401300	TRO ANALYZER CLEANING SET	\$ 78.82
940-411	ADAPTER/7.5MM ELECTRODES	\$ 250.54
940-412	ADAPTER/9.5MM ELECTRODES	\$ 252.46
940-413	ADAPTER/12MM ELECTRODES	\$ 252.46
9406700	oo CKT BD ASSY, ENCAPSULATED, MAIN PCB, 1735 TRO	\$ 4,143.48
9408100	KIT, SUPPORT, CIRCUIT BOARD, CL17	\$ 568.52
9412000	CABLE, MODULE, SL1000	\$ 31.22
9418100	CUP, SAMPLE, SL1000	\$ 71.31
9422200	RETAINER, PROBE CABLE, SL1000	\$ 75.74
9425200	ASSY., PACK, CHEMKEY, 25 PIECE, TOTAL AMMONIA	\$ 240.30
9427900	ASSY., PACKAGE, INDIVIDUAL, CHEMKEY, VERIFICATION, CHLORINE	\$ 78.58
9429000	ASSY., PACK, CHEMKEY, 25 PIECE, FREE CHLORINE	\$ 68.83
9429100	ASSY., PACK, CHEMKEY, 25 PIECE, TOTAL CHLORINE	\$ 68.83



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

9429200	ASSY., PACK, CHEMKEY, 25 PIECE, COPPER	\$	151.41
9429300	ASSY., PACK, CHEMKEY, 25 PIECE, NITRITE	\$	132.42
9429400	ASSY., PACK, CHEMKEY, 25 PIECE, MONO- CHLOROMINE	\$	172.00
9429500	ASSY., PACK, CHEMKEY, 25 PIECE FREE AMMONIA	\$	168.56
9429600	KTO CHEMKEY, 25 PIECE MONO-CHLOR. AND 25 PIECE FREE AMMONIA	\$	377.61
94299	PHENOLPHTHALEIN PP PK/100	\$	60.35
94299-LM	PHENOLPHTHALEIN PP PK/100	\$	60.35
9430000	ee ASSY, INSTRUMENT, SL1000	\$	12,334.08
9430800	GASKET, CONNECTOR, PROBE, SL1000	\$	39.89
9430900	FRAME, LCD, SL1000	\$	52.81
9431000	ASSY, MODULE, PUMP, SL1000	\$	2,293.37
9431700	ASSY, NOSECONE, SL1000	\$	121.99
94323	BROMCRESOL GR-METH RED PWD, 15G	\$	170.56
9432900	LIGHTSHIELD, NOSECONE, SL1000	\$	39.89
9433100	ASSY, MODULE, OPTICS, RIGHT, SL1000	\$	1,551.40
9433200	ASSY, MODULE, OPTICS, LEFT, SL1000	\$	1,551.40
9434100	SCREW 6-32 X .375 CAPTIVE PAN PHIL SST	\$	6.55
9434600	si CKT BD ASSY, USB CIRCUIT BD, SL1000	\$	512.64
9435000	si CKT BD ASSY., SAMPLE DETECT, SL1000	\$	177.69
9435900	si CABLE ASSY, SAMPLE DETECT, SL1000	\$	194.65
9436400	ASSY, HOUSING, UPPER, SL1000	\$	493.36
9436500	ASSY, HOUSING, LOWER, SL1000	\$	466.38
9436600	ASSY, DOOR, BATTERY, SL1000	\$	129.89
9436700	ASSY, HANDSTRAP, SL1000	\$	57.82
9436800	ASSY., PACKAGE, INDIVIDUAL, CHEMKEY, VERIFICATION, SYSTEM	\$	118.33
943926	CONNECTOR HEX MIN 126 217 E	\$	141.65
94399	BROMCRESOL GR-METH RED PP PK/100	\$	60.35
94399-LM	BROMCRESOL GR-METH RED PP PK/100	\$	60.35
94799	CALVER 2 POWDER PILLOWS PK/100	\$	67.21
94799-LM	CALVER 2 POWDER PILLOWS PK/100	\$	67.21
9500.99.00004	POLYMETRON 9500 CONTROLLER, AC-DC	\$	4,586.74
9500.99.00604	POLYMETRON 9500 CONTROLLER, AC-DC, POLYCOND	\$	4,875.82
9500.99.00664	POLYMETRON 9500 CONTROLLER, AC-DC, 2 POLYCOND	\$	5,974.32
9500.99.00704	POLYMETRON 9500, AC-DC, 1POLY PH	\$	4,933.63
9500.99.00764	POLYMETRON 9500 CONTR, AC-DC, POLYCOND, POLY PH	\$	5,974.32
9500.99.00774	POLYMETRON 9500 CONTROLLER, AC-DC, 2POLY PH	\$	5,974.32
9500.99.00804	POLYMETRON 9500 AMPERO, AC-DC	\$	10,850.14
9500.99.01604	POLYMETRON 9500 CONTROLLER AC-DC, POLYCOND, MODBUS	\$	6,205.58
9500.99.01664	POLYMETRON 9500 CONTROLLER AC-DC 2POLYCOND MODBUS	\$	7,284.82
9500.99.01704	POLYMETRON 9500 CONTROLLER, AC-DC, POLY PH, MODBUS	\$	6,205.58
9500.99.01764	POLYMETRON 9500 CONTR AC-DC POLYCOND POLYPH MODBUS	\$	7,284.82
9500.99.01774	POLYMETRON 9500 CONTROLLER, AC-DC, 2POLY PH, MODBU	\$	7,284.82
9500.99.01804	POLYMETRON 9500 AMPERO, AC-DC, MODBUS	\$	12,237.72
9500.99.03604	POLYMETRON 9500 CONTROLLER AC-DC POLYCOND PROFIBUS	\$	6,205.58



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

9500.99.03664	POLYMETRON 9500 CONTR AC-DC 2POLYCOND PROFIBUS	\$	7,284.82
9500.99.03704	POLYMETRON 9500 CONTR, AC-DC, POLY PH, PROFIBUS	\$	6,205.58
9500.99.03764	POLYMETRON 9500 CONTR AC-DC POLYPH POLYCOND PROFIBUS	\$	7,284.82
9500.99.03774	POLYMETRON 9500 CONTR, AC-DC, 2POLY PH, PROFIBUS	\$	7,284.82
9500.99.03804	POLYMETRON 9500 AMPERO, AC-DC, PROFIBUS	\$	12,237.72
9500.99.05604	POLYMETRON 9500 CONTROLLER, AC-DC, POLYCOND, HART	\$	6,205.58
9500.99.05664	POLYMETRON 9500 CONTROLLER, AC-DC, 2POLYCOND, HART	\$	7,284.82
9500.99.05704	POLYMETRON 9500 CONTROLLER, AC-DC, POLY PH, HART	\$	6,205.58
9500.99.05764	POLYMETRON 9500 CONTR, AC-DC POLY PH POLYCOND HART	\$	7,284.82
9500.99.05774	POLYMETRON 9500 CONTROLLER, AC-DC, 2 POLY PH, HART	\$	7,284.82
9500.99.09604	POLYMETRON 9500 CONTR AC-DC POLYCOND 4-20MA OUT	\$	5,781.60
9500.99.09664	POLYMETRON 9500 CONTR AC-DC 2POLYCOND 4-20 MA OUT	\$	6,899.38
9500.99.09704	POLYMETRON 9500 CONTR, AC-DC, POLY PH, 4-20MA OUT	\$	5,993.59
9500.99.09764	POLY 9500 CONTR AC-DC POLYCOND POLYPH 4-20MA OUT	\$	6,899.38
9500.99.09774	POLYMETRON 9500 CONTR, AC-DC, 2POLY PH, 4-20MA OUT	\$	6,899.38
9500.99.09804	POLYMETRON 9500 AMPERO, 4-20MA OUTPUT	\$	11,736.65
9500.99.70004	POLYMETRON 9500 CONTROLLER, 24VDC	\$	6,860.83
9500.99.70604	POLYMETRON 9500 CONTROLLER, 24VDC, POLYCOND	\$	7,130.64
9500.99.70664	POLYMETRON 9500 CONTROLLER, 24VDC, 2 POLYCOND	\$	7,708.80
9500.99.70704	POLYMETRON 9500 CONTROLLER, 24VDC, POLY PH	\$	5,685.24
9500.99.70764	POLYMETRON 9500 CONTROLLER, 24VDC POLYCOND POLY PH	\$	8,132.78
9500.99.70774	POLYMETRON 9500 CONTROLLER, 24VDC, 2POLY PH	\$	6,803.02
9500.99.70804	POLYMETRON 9500 AMPERO, 24VDC	\$	11,929.37
9500.99.71604	POLYMETRON 9500 CONTROLLER, 24VDC POLYCOND MODBUS	\$	7,053.55
9500.99.71664	POLYMETRON 9500 CONTROLLER, 24VDC 2POLYCOND MODBUS	\$	8,132.78
9500.99.71704	POLYMETRON 9500 CONTR, 24VDC, 1POLY PH, MODBUS	\$	7,053.55
9500.99.71764	POLYMETRON9500 CONTR 24VDC POLY PH POLYCOND MODBUS	\$	8,132.78
9500.99.71774	POLYMETRON 9500 CONTR, 24VDC, 2POLY PH, MODBUS	\$	8,132.78
9500.99.71804	POLYMETRON 9500 AMPERO, 24VDC, MODBUS	\$	13,066.42
9500.99.73604	POLYMETRON 9500 CONTROLLER 24VDC POLYCOND PROFIBUS	\$	7,053.55
9500.99.73664	POLYMETRON 9500 CONTR 24VDC 2POLYCOND PROFIBUS	\$	8,132.78
9500.99.73704	POLYMETRON 9500 CONTR, 24VDC, 1POLY PH, PROFIBUS	\$	7,053.55
9500.99.73764	POLY 9500 CONTR 24VDC POLYPH POLYCOND PROFIBUS	\$	8,132.78
9500.99.73774	POLYMETRON 9500 CONTR, 24VDC, 2POLY PH, PROFIBUS	\$	8,132.78
9500.99.73804	POLYMETRON 9500 AMPERO, 24VDC, PROFIBUS	\$	13,066.42
9500.99.75604	POLYMETRON 9500 CONTROLLER, 24VDC, POLYCOND, HART	\$	7,053.55
9500.99.75664	POLYMETRON 9500 CONTROLLER, 24VDC, 2POLYCOND, HART	\$	8,132.78
9500.99.75704	POLYMETRON 9500 CONTROLLER, 24VDC, 1 POLY PH, HART	\$	7,053.55
9500.99.75764	POLYMETRON 9500 CONTR 24VDC POLYPH POLYCOND HART	\$	8,132.78
9500.99.75774	POLYMETRON 9500 CONTROLLER, 24VDC, 2 POLY PH, HART	\$	8,132.78
9500.99.75804	POLYMETRON 9500 AMPERO, 24VDC, HART	\$	13,066.42
9500.99.79604	POLYMETRON 9500 CONTR 24VDC POLYCOND 4-20 MA OUT	\$	6,610.30
9500.99.79664	POLYMETRON 9500 CONTR 24VDC 2POLYCOND 4-20 MA OUT	\$	7,708.80
9500.99.79704	POLYMETRON 9500 CONTR, 24VDC, 1POLY PH, 4-20MA OUT	\$	9,462.55
9500.99.79764	POLY 9500 CONTR, 24VDC POLY PH POLYCOND 4-20MA OUT	\$	7,708.80



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

9500.99.79774	POLYMETRON 9500 CONTR, 24VDC, 2POLY PH, 4-20MA OUT	\$	7,708.80
9500.99.79804	POLYMETRON 9500 AMPERO, 24VDC, 4-20MA OUTPUT	\$	12,661.70
950P.99.000D4	oo KIT, POLY 9500 CONTROLLER WITH 8350.5 PH SENSOR	\$	10,734.50
9523.99.00662	sc200 CATIONIC COND, AC-DC, HACH	\$	7,535.35
9523.99.00664	sc200 CATIONIC COND, AC-DC, MODBUS, HACH	\$	7,535.35
9523.99.01664	POLYMETRON 9523 CONTR ONLY, AC-DC	\$	8,903.66
9523.99.01P4	POLYMETRON 9523 CATIONIC COND, AC-DC, MODBUS	\$	23,126.40
9523.99.03664	POLYMETRON 9523 CONTR ONLY, AC-DC, PROFIBUS	\$	8,903.66
9523.99.03P4	POLYMETRON 9523 CATIONIC COND, AC-DC, PROFIBUS	\$	23,126.40
9523.99.05664	POLYMETRON 9523 CONTR ONLY, AC-DC, HART	\$	8,903.66
9523.99.05P4	POLYMETRON 9523 CATIONIC COND, AC-DC, HART	\$	23,126.40
9523.99.09664	POLYMETRON 9523 CONTR ONLY, AC-DC, 4-20MA OUT	\$	8,402.59
9523.99.09P2-LM	9523 CATIONIC COND, AC-DC, 4-20MA OUTPUT, HACH	\$	22,760.23
9523.99.09P4	POLYMETRON 9523 CATIONIC COND, AC-DC, 4-20MA OUT	\$	22,760.23
9523.99.70662	sc200 CATIONIC COND, 24VDC, HACH	\$	8,286.96
9523.99.70664	POLYMETRON 9523 CONTR ONLY, 24VDC	\$	8,286.96
9523.99.71664	POLYMETRON 9523 CONTR ONLY, 24VDC, MODBUS	\$	9,713.09
9523.99.71P4	POLYMETRON 9523 CATIONIC COND, 24VDC, MODBUS	\$	24,012.91
9523.99.73664	POLYMETRON 9523 CONTR ONLY, 24VDC, PROFIBUS	\$	9,713.09
9523.99.73P4	POLYMETRON 9523 CATIONIC COND, 24VDC, PROFIBUS	\$	24,012.91
9523.99.75664	POLYMETRON 9523 CONTR ONLY, 24VDC, HART	\$	9,713.09
9523.99.79664	POLYMETRON 9523 CONTR ONLY, 24VDC, 4-20MA OUT	\$	9,289.10
9523.99.79P4	POLYMETRON 9523 CATIONIC COND 24VDC, 4-20MA OUTPUT	\$	23,569.66
9524300	FILTER, INTERFERENCE 420NM	\$	215.85
9524400	FILTER, INTERFERENCE 520NM	\$	215.85
9524500	FILTER, INTERFERENCE 610NM	\$	215.85
9524600	FILTER, INTERFERENCE 560NM	\$	215.85
9525.99.7030	Degas System Heater, 120V	\$	3,141.34
9525.99.7031	Degas System Heater, 240V	\$	3,892.94
9525.99.7050	Degas Temperature Controller	\$	2,543.90
9525.99.7060	Degas System Sample Cooler	\$	3,256.97
9525.99.7070	Degas Fuse, 15.00A 600V	\$	163.81
9525.99.7071	Degas Fuse, 0.50A 600V	\$	163.81
9525.99.7072	Degas Fuse, 8.00A 600V	\$	171.71
9525700	SENSOR MODULE, PH/ORP, sc200 POLYMETRON, HACH	\$	1,464.67
9525800	SENSOR MODULE, CONDUCTIVITY, SC200 POLYMETRON, HACH	\$	1,412.64
9525900	SENSOR MODULE, AMPERO, SC200 POLYMETRON, HACH	\$	6,436.85
9526.99.00P4	POLYMETRON 9526 CONDUCTIVITY CAL SYSTEM AC-DC	\$	17,017.18
9526.A.1000	ANNUAL RECALIBRATION OF THE 9526 PURECAL	\$	3,218.42
9531000	ee POCKET PRO, PH, TESTER	\$	354.60
9531100	POCKET PRO, ORP, TESTER	\$	531.91
9531200	POCKET PRO, TDS LOW RANGE, TESTER	\$	341.11
9531300	POCKET PRO, TDS HIGH RANGE, TESTER	\$	339.19
9531400	POCKET PRO, COND LOW RANGE, TESTER	\$	339.19
9531500	POCKET PRO, COND HIGH RANGE, TESTER	\$	339.19
9531600	POCKET PRO, SALT, TESTER	\$	387.37
9531700	POCKET PRO, TEMP, TESTER	\$	316.06

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

9531701	POCKET PRO, TEMP, REPL. SENSOR	\$	110.81
9531900	Pocket Pro Tester sample cup/cap	\$	13.30
9532000	ee POCKET PRO+, PH, TESTER	\$	572.38
9532001	POCKET PRO+, PH, REPLACEMENT SENSOR	\$	296.79
9532100	POCKET PRO+, ORP, TESTER	\$	984.80
9532101	POCKET PRO+, ORP, REPLACEMENT SENSOR	\$	529.98
9532700	POCKET PRO+, MULTI-1, COND/TDS/SAL, TESTER	\$	572.38
9532701	POCKET PRO+, MULTI-1, COND/TDS/SAL, REPL. SENSOR	\$	296.79
9532800	ee POCKET PRO+, MULTI-2, PH/COND/TDS/SAL, TESTER	\$	770.88
9532801	POCKET PRO+, MULTI-2, PH/COND/TDS/SAL, REPL. SENSOR	\$	350.75
9532900	Pocket Pro Plus Tester sample cup/cap	\$	13.30
9534500	GAUGE,PRESSURE, 0-30 PSIG,1/8NPT	\$	566.60
9558600	CELL CAP, COLOR, 5500sc AMC	\$	128.93
9558900	O-RING, 0.472"IDX0.59"ODX0.059W, VITON, SHORE A75	\$	5.78
9559000	ASSY, CELL BODY W/LIGHT PIPE, 5500sc AMC	\$	362.31
9559600	FAN HOUSE, 5500sc AMC	\$	46.25
9559800	MOUNTING, CUT-OUT, COLOR, 5500sc AMC	\$	208.14
9560100	si CABLE ASSY, THERMO SENSOR, 5500sc AMC	\$	31.41
9560400	ASSY, Y-STRAINER W/FTG	\$	229.34
9560501	KIT, MAINTENANCE, 5500sc AMC, SINGLE CHANNEL	\$	3,430.42
9560502	KIT, MAINTENANCE, 5500sc AMC, TWO CHANNEL	\$	3,777.31
956-134	TITRATION VESSEL POLY 0.5-5ML	\$	537.69
9562000	si CABLE ASSY, COLORIMETER, 5500sc AMC	\$	1,568.74
9562600	SCREW CAPTIVE M3X0.5X13MM PAN SST TORX	\$	11.37
9562800	si CKT BD ASSY, LEAK DETECTOR, 5500sc AMC	\$	377.73
9563100	ANALYZE PLATE, W/MAGNET BAR, 5500sc AMC	\$	383.51
9563200	HINGE PLATE, 5500sc AMC	\$	456.75
9563300	HINGE POST, 5500sc AMC	\$	25.05
9563400	HINGE COVER, 5500sc AMC	\$	41.43
9563500	MESH, FAN GUARD, 40 X 40MM, SST	\$	12.33
9563700	si CABLE ASSY, ANALYZER DATA, 5500sc AMC	\$	402.78
9563900	STIR BAR, TEFLON 8X3MM	\$	10.21
9564800	BRACKET, PUMP, 5500sc AMC	\$	95.59
9565200	BODY, SAMPLE HOLDER, 5500sc AMC	\$	40.09
9565300	CAP, SAMPLE HOLDER, 5500sc AMC	\$	16.00
9565400	O-RING, 43.69IDX49.69ODX3W, EPDM, SHORE A70	\$	37.77
9565500	si CABLE ASSY, AC POWER, 5500sc AMC	\$	283.30
9565501	KIT, INSTALLATION, 5500SC AMC, AC, 1 CH, US	\$	1,050.32
9565502	KIT, INSTALLATION, 5500SC AMC, AC, 2 CH, US	\$	1,403.00
9565700	MANIFOLD, 7-WAY, 5500sc AMC	\$	684.16
9565800	FTG, BARB ELBOW W/1/4"-28UNF, 1/8"-1/8" TBG, HDPE	\$	8.67
9566000	FILTER, 40UM, 1/4", PP, 100PSI	\$	183.85
9566200	DRAIN CELL, 5500sc AMC	\$	17.34
9566300	DRAIN COVER, 5500sc AMC	\$	15.61
9566400	COVER, COLORIMETER PCBA, 5500sc AMC	\$	126.23
9566900	PANEL, SAMPLE PREP, 5500sc AMC	\$	853.75
9567000	COVER, ANALYTICS PCB, 5500sc AMC	\$	80.36



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

9567100	si CABLE ASSY, ANALYZER POWER, 5500sc AMC	\$	393.15
9567500	si CABLE ASSY, 0127 VALVE W/CABLE, 5500sc AMC	\$	1,004.07
9567600	FAN, 12V, 40 X 40 X 10MM	\$	66.87
9567700	CABLE ASSY, FLOAT LEVEL SENSOR, 5500sc AMC	\$	86.92
9567800	ASSY, PINCH VALVE, 2-WAY, 5500sc AMC	\$	601.29
9567900	ASSY, PINCH VALVE, 3-WAY, 5500sc AMC	\$	693.79
9568000	si CABLE ASSY, FAN, 5500sc AMC	\$	210.06
9568100	inn si ASSY, LIQUID PUMP, EXTENDED, W/BACKET, 5500sc AMC	\$	1,293.15
9568400	si CABLE ASSY, THERMO CUT-OUT, 5500sc AMC	\$	20.24
9568500	ASSY, COLORIMETER, 5500sc AMC	\$	3,546.05
9568601	ASSY, CAPILARY TBG, COLORIMETER, RED	\$	164.97
9568602	ASSY, CAPILARY TBG, COLORIMETER, BLK	\$	161.50
9576100	DISPLAY ASSY, LCD 240 X 160 DOTS	\$	54.15
9582.99.00P4	POLYMETRON 9582 DISSOLVED O2, AC-DC	\$	28,753.82
9582.99.01P4	POLYMETRON 9582 DISSOLVED O2, AC-DC, MODBUS	\$	30,488.30
9582.99.03P4	POLYMETRON 9582 DISSOLVED O2, AC-DC, PROFIBUS	\$	30,488.30
9582.99.05P4	POLYMETRON 9582 DISSOLVED O2, AC-DC, HART	\$	30,488.30
9582.99.09P4	POLYMETRON 9582 DISSOLVED O2, AC-DC, 4-20MA OUT	\$	29,910.14
9582.99.1000	O2 SENSOR, MODULE, ACCESS. FOR 9582; EXCL. CONTR.	\$	22,432.61
9582.99.70P4	POLYMETRON 9582 DISSOLVED O2, 24VDC	\$	29,813.78
9582.99.71P4	POLYMETRON 9582 DISSOLVED O2, 24VDC, MODBUS	\$	31,548.26
9582.99.73P4	POLYMETRON 9582 DISSOLVED O2, 24VDC, PROFIBUS	\$	31,548.26
9582.99.75P4	POLYMETRON 9582 DISSOLVED O2, 24VDC, HART	\$	31,548.26
9582.99.79P4	POLYMETRON 9582 DISSOLVED O2, 24VDC, 4-20MA OUT	\$	29,312.71
9586.99.00P4	POLYMETRON 9586 O2 SCAVENGER, AC-DC	\$	64,522.66
9586.99.01P4	POLYMETRON 9586 O2 SCAVENGER, AC-DC, MODBUS	\$	66,218.59
9586.99.03P4	POLYMETRON 9586 O2 SCAVENGER, AC-DC, PROFIBUS	\$	66,218.59
9586.99.05P4	POLYMETRON 9586 O2 SCAVENGER, AC-DC, HART	\$	62,576.18
9586.99.09P4	POLYMETRON 9586 O2 SCAVENGER, AC-DC, 4-20MA OUTPUT	\$	65,678.98
9586.99.70P4	POLYMETRON 9586 O2 SCAVENGER, 24VDC	\$	65,640.43
9586.99.71P4	POLYMETRON 9586 O2 SCAVENGER, 24VDC, MODBUS	\$	67,317.10
9586.99.73P4	POLYMETRON 9586 O2 SCAVENGER, 24VDC, PROFIBUS	\$	67,317.10
9586.99.75P4	POLYMETRON 9586 O2 SCAVENGER, 24VDC, HART	\$	67,317.10
9586.99.79P4	POLYMETRON 9586 O2 SCAVENGER, 24VDC, 4-20MA OUTPUT	\$	66,738.94
9595000	GRATING ASSY	\$	659.10
9595500	PROTECTIVE COVER,DR1900	\$	147.05
9596000	1/2" SPHERICAL MIRROR ASSY	\$	202.36
959614	O-RING NIT .070W X .864ID	\$	17.34
9597000	1" SPHERICAL MIRROR ASSY	\$	202.36
959714	O-RING NIT .070W X .145ID	\$	13.10
9598000	BATTERY COVER ASSY	\$	86.92
9598300	BATTERY HOLDER ASSY, DR1900	\$	58.59
9598500	PRISM,RT ANGLE, 5X5MM,587.6NM	\$	143.77
9598600	DOUBLE SLIT, PRISM RT ANGLE	\$	58.59
9598700	BEAMSPLITTER, 50-50,POLKA DOT,45	\$	143.77
9599000	TOP HOUSING ASSY	\$	319.92
9599300	HOUSING ASSM, BOTTOM, DR1900	\$	292.93

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

9599500	SPRING,PIVOT,GRATING HOLDER	\$	58.59
9599800	COMPARTMENT,SAMPLE CELL,DR1900	\$	55.12
9600000	PHOTODIODE FLEXIBLE PCB ASSY	\$	749.68
960024	Tubing, TYGON 1/8ID X 1/4OD	\$	39.89
960032	Tubing,HYTREL,.170 X .250	\$	29.49
9601000	FLASH BOARD ASSY	\$	601.29
9602000	MAIN BOARD ASSY	\$	1,734.48
9603000	MOTOR ASSY	\$	491.44
9603600	VALVE, 2-WAY SOLENOID, 0-6.5BAR	\$	620.56
9604000	CABLE SET	\$	81.91
9605000	FLAP ASSY	\$	86.92
9606000	SCREW SET	\$	58.59
9608400	SCREW, M4X120, CROSS PAN HEAD, SST	\$	2.51
9609300	DUST COVER	\$	158.03
9609800	KIT, ADAPTER SET, 4 PC., DR1900	\$	206.21
9609900	KIT, ADAPTER SET, 3 PC., DR1900	\$	244.75
9610.KTO.S0.A1E	KTO: 9610sc SILICA, 1CH/AC/EU, CHEM	\$	74,833.18
9610.KTO.S0.A1U	KTO: 9610sc SILICA, 1CH/AC/US, CHEM	\$	74,833.18
9610.KTO.S0.A2E	KTO: 9610sc SILICA, 2CH/AC/EU, CHEM	\$	79,400.64
9610.KTO.S0.A2U	KTO: 9610sc SILICA, 2CH/AC/US, CHEM	\$	79,400.64
9610.KTO.S0.A4E	KTO: 9610sc SILICA, 4CH/AC/EU, CHEM	\$	88,689.74
9610.KTO.S0.A4U	KTO: 9610sc SILICA, 4CH/AC/US, CHEM	\$	88,689.74
9610.KTO.S0.A6E	KTO: 9610sc SILICA, 6CH/AC/EU, CHEM	\$	97,978.85
9610.S0.A1E	9610sc SILICA ANALYZER, AC, 1CH, EU	\$	76,182.22
9610.S0.A1U	9610sc SILICA ANALYZER, AC, 1CH, US	\$	76,182.22
9610.S0.A2E	9610sc SILICA ANALYZER, AC, 2CH, EU	\$	80,711.14
9610.S0.A2U	9610sc SILICA ANALYZER, AC, 2CH, US	\$	80,711.14
9610.S0.A4E	9610sc SILICA ANALYZER, AC, 4CH, EU	\$	90,000.24
9610.S0.A4U	9610sc SILICA ANALYZER, AC, 4CH, US	\$	90,000.24
9610.S0.A6E	9610sc SILICA ANALYZER, AC, 6CH, EU	\$	99,212.26
9610.S0.A6U	9610sc SILICA ANALYZER, AC, 6CH, US	\$	99,212.26
9610.S0.D1U	9610sc SILICA ANALYZER, DC, 1CH, US	\$	76,182.22
9610.S0.D2U	9610sc SILICA ANALYZER, DC, 2 CH, US	\$	80,711.14
9610.S0.D4U	9610sc SILICA ANALYZER, DC, 4CH, US	\$	90,000.24
9610.S0.D6U	9610sc SILICA ANALYZER, DC, 6CH, US	\$	99,212.26
9610000	MAIN OPTICAL TUBE ASSY	\$	86.92
9611.KTO.PH.A1E	KTO: 9611sc PHOS HR, 1CH/AC/EU, CHEM	\$	74,833.18
9611.KTO.PH.A2E	KTO: 9611sc PHOS HR, 2CH/AC/EU, CHEM	\$	79,400.64
9611.KTO.PH.A4E	KTO: 9611sc PHOS HR, 4CH/AC/EU, CHEM	\$	88,689.74
9611.KTO.PH.A4U	KTO: 9611sc PHOS HR, 4CH/AC/US, CHEM	\$	88,689.74
9611.KTO.PL.A1E	KTO: 9611sc PHOS LR, 1CH/AC/EU, CHEM	\$	74,833.18
9611.KTO.PL.A2E	KTO: 9611sc PHOS LR, 2CH/AC/EU, CHEM	\$	79,400.64
9611.KTO.PL.A4E	KTO: 9611sc PHOS LR, 4CH/AC/EU, CHEM	\$	88,689.74
9611.KTO.PL.D4U	KTO: 9611sc PHOS LR, 4CH/DC/US, CHEM	\$	88,689.74
9611.PH.A1E	9611sc PHOSPHATE HR ANALYZER, AC, 1CH, EU	\$	76,182.22
9611.PH.A1U	9611sc PHOSPHATE HR ANALYZER, AC, 1CH, US	\$	76,182.22
9611.PH.A2E	9611sc PHOSPHATE HR ANALYZER, AC, 2CH, EU	\$	80,711.14



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

9611.PH.A2U	9611sc PHOSPHATE HR ANALYZER, AC, 2CH, US	\$	80,711.14
9611.PH.A4E	9611sc PHOSPHATE HR ANALYZER, AC, 4CH, EU	\$	90,000.24
9611.PH.A4U	9611sc PHOSPHATE HR ANALYZER, AC, 4CH, US	\$	90,000.24
9611.PH.D1U	9611sc PHOSPHATE HR ANALYZER, DC, 1CH, US	\$	76,182.22
9611.PH.D2U	9611sc PHOSPHATE HR ANALYZER DC, 2CH, US	\$	80,711.14
9611.PH.D4U	9611sc PHOSPHATE HR ANALYZER, DC, 4CH, US	\$	90,000.24
9611.PL.A1E	9611sc PHOSPHATE LR ANALYZER, AC, 1CH, EU	\$	76,182.22
9611.PL.A1U	9611sc PHOSPHATE LR ANALYZER, AC, 1CH, US	\$	76,182.22
9611.PL.A2E	9611sc PHOSPHATE LR ANALYZER, AC, 2CH, EU	\$	80,711.14
9611.PL.A2U	9611sc PHOSPHATE LR ANALYZER, AC, 2CH, US	\$	80,711.14
9611.PL.A4E	9611sc PHOSPHATE LR ANALYZER, AC, 4CH, EU	\$	90,000.24
9611.PL.A4U	9611sc PHOSPHATE LR ANALYZER, AC, 4CH, US	\$	90,000.24
9611.PL.D1U	9611sc PHOSPHATE LR ANALYZER, DC, 1CH, US	\$	76,182.22
9611.PL.D2U	9611sc PHOSPHATE LR ANALYZER, DC, 2CH, US	\$	80,711.14
9611.PL.D4U	9611sc PHOSPHATE LR ANALYZER, DC, 4CH, US	\$	90,000.24
9612000	Optical Bench	\$	3,507.50
9613400	PROTECTIVE COVER, ECHEM BOARD ,MS6100	\$	47.41
9616000	SENSOR, TEMP, MS6100	\$	103.30
9616700	HOLE PLUG, 3/4",ECHEM	\$	31.03
9617100	PROTECTIVE COVER, MAINBOARD ,MS6100	\$	68.03
9621000	CABLE CLAMP, 6.4MM/4.6MM	\$	0.77
9622500	si ASSY, LIQUID PUMP W/ BRACKET 5500sc AMC	\$	915.42
9623400	FILTER, DISC, 100UM, 1"NPT	\$	894.22
9624500	REPLACEMENT, FILTER DISC ELEMENT, 100 UM, W/O-RING X 2PCS	\$	431.69
9624600	KIT, CARTON SET, DR1900	\$	132.59
9624700	DR1900 HACH	\$	14,222.74
9624700NB	DR1900 HACH - NO BATTERY	\$	13,278.41
9625000	PIPE ,OD 3/8 ID 1/4, PVC ,CLEAR	\$	5.97
9625300	CABLE, 3 PAIR, 2 COND 24AWG STRD TWISTED PR	\$	16.38
96299	UNIVER 3 POWDER PILLOWS PK/100	\$	66.13
96299-LM	UNIVER 3 POWDER PILLOWS PK/100	\$	66.13
9634800	BATTERY CAP, POCKET CLIP, POCKET PRO/PRO+ TESTERS	\$	37.77
9634900	POCKET PRO/PRO+ BATTERY FLAP, inner cover with battery contacts	\$	25.63
9640400	REAGENT TRAY, PROCESS ANALYZER	\$	211.99
9646200	CKT BD ASSY, MAIN, TL2300 TURB	\$	2,177.74
9647200	ENCLOSURE, COVER, TL2300 TURB	\$	323.77
9647300	FRAME, COVER, TL2300 TURB	\$	86.53
9647400	ENCLOSURE, BOTTOM, TL2300 TURB	\$	323.77
9647500	HOLDER, SAMPLE CELL, TL2300 TURB	\$	43.17
9647600	SLIDE DOOR, TURBIDIMETER, TL2300 TURB	\$	161.11
9647700	ACCESS DOOR, TL2300 TURB	\$	166.70
9647800	FOAM, 23X7X5 MM	\$	11.95
9647900	HOLDER, FLAT FERRIT BEAD, TL2300 TURB	\$	27.56
9648000	TAPE, AL FOIL, 1/2"W, W/ADHESIVE	\$	11.95
9648300	ASSY, TOP COVER W/O MAGNET, TL2300 TURB	\$	381.59

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

9648400	KEYPAD, TL2300 TURB	\$	260.17
9648700	SUPPORTER, TL2300 TURB	\$	31.41
9648800	KEYPAD, TL2350 TURB	\$	260.17
9649000	CKT BD ASSY, MAIN, TL2350 TURB	\$	2,505.36
9649100	DUST COVER, TL2300 TURB	\$	170.17
9649300	SEAL GASKET, LCD	\$	35.85
9649400	ASSY, LCD FRAME W/BACK FOAM, TL2300 TURB	\$	94.43
9649600	BUMPER FOAM, TL2300 TURB	\$	11.95
9650600	ACCESS DOOR, TL2310 TURB	\$	133.75
9650700	ASSY, FPC W/FERRIT, TL2300 TURB	\$	55.31
9650800	CABLE ASSY, INFRA LED, 200MM LG	\$	1,146.68
9650900	CABLE ASSY, DETECTOR LED, 200MM LG	\$	356.53
9651300	FILTER ASSY, DTR, TL2310 TURB	\$	244.75
9651800	CABLE ASSY, DETECTOR FORWARD SCATTER	\$	258.24
9651900	CABLE CLAMP, NYLON, DIA 5/16"	\$	11.95
9652000	CABLE ASSY, DETECTOR, 90 DEGREE	\$	940.47
9652100	CABLE ASSY, DETECTOR BACK SCATTER	\$	400.86
9652500	CABLE ASSY, DETECTOR TRANSMISION	\$	1,786.51
9652900	KEYPAD, TL2310 TURB	\$	260.17
9653000	KEYPAD, TL2360 TURB	\$	260.17
9653500	FILTER ASSY, EPA TL2300 TURB	\$	894.22
9654900	ASSY, TOP COVER W/MAGNET, TL2310 TURB	\$	397.00
9655200	LED LIGHT ASSY, TL2310 TURB	\$	2,466.82
9655300	KIT, SCREW TL2310 OPTICAL MODULE	\$	11.95
9655400	KIT, SCREW TL2300 OPTICAL MODULE	\$	11.95
9655600	KIT, SCREW TL2310 & TL2360	\$	39.31
9655700	KIT, SCREW TL2300 & TL2350	\$	39.31
9657100	KIT, MANUAL FLOW CELL, TL23X0 TURB	\$	3,815.86
9657253	STABLCAL STD, 2000 NTU 1000ML	\$	489.63
9657800	ACCESSORY KIT, TL2300/TL2310 LAB TURB, EPA, CN	\$	2,447.54
9657900	ACCESSORY KIT, TL2350/TL2360 LAB TURB, EPA, CN	\$	2,890.80
9658300	KIT, AUTOMATED FLOW CELL, TL2300 TURB	\$	4,856.54
9658500	ASSY, BOTTOM ENCLOSURE W/FOOT, TL23 LAB TURB	\$	310.28
9658600	ASSY, TOP COVER W/O MAGNET, TL2300/50	\$	375.80
9658700	ASSY, TOP COVER W/MAGNET, TL2310/60	\$	387.37
9658900	LAMP, HALOGEN WITH LEADS, WA	\$	759.32
9659000	LAMP REPLACEMENT ASSY, WA, LAB TURB	\$	957.82
9659100	KIT, CANISTER, NA5600sc	\$	302.57
9660000	UPGRADE KIT, DRAIN MAINTENANCE, 5500sc AMC	\$	595.50
9660500	KIT, ONE YEAR SPARE PARTS, NA5600sc	\$	5,087.81
9660600	KIT, TECHNICAL, NA5600sc	\$	11,254.85
9669400	CKT BD ASSY, ADAPTER, TL2300 TURB	\$	113.70
9672000	FERRITE BEAD, 128OHMS 100MHz, 26X10MM, FOR FLAT CABLE	\$	23.32
9672200	CABLE, SENSOR BD TO MAIN BD, 130MM LG	\$	27.56
9672600	CABLE, GROUND, 150MM LG, OPTICAL BLOCK	\$	23.32
9672800	FERRITE BEAD, 220OHM 100MHz, 34X25MM, FOR FLAT CABLE	\$	23.32



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

9672900	FERRITE BEAD, 115OHMS 100MHz, 28X14.6MM, FOR FLAT CABLE	\$	23.32
9673700	AC ADAPTER, 100-240VAC, 12VDC, 3.4A 2.5MM ESTAR VI	\$	292.93
9675700	KTO: STABLCAL KIT, TL23X0, VERI & CALI, INTERNAL ONLY	\$	360.39
96799	DIPHENYLCARBAZONE PP PK/100	\$	101.90
96799-LM	DIPHENYLCARBAZONE PP PK/100	\$	101.90
96800	nn ff CLIPPER, FOR MEDIUM PWD PLWS	\$	18.31
9680000	TL2300 LAB TURB, INST ONLY	\$	8,325.50
9680100	TL2350 LAB TURB, INST ONLY	\$	13,721.66
9680800	CKT BD ASSY, MAIN, TL2310 TURB	\$	2,177.74
9681400	OPTIC BLOCK, TL2300 TURB	\$	3,102.79
9681500	OPTIC BLOCK, TL2350 TURB	\$	3,353.33
9681600	OPTIC BLOCK, TL2310 TURB	\$	3,892.94
9681700	OPTIC BLOCK, TL2360 TURB	\$	5,145.62
9682200	CKT BD ASSY, SENSOR, TL2310 TURB	\$	90.77
9684200	CKT BD ASSY, MAIN, TL2360 TURB	\$	2,505.36
9713300	ASSY, PKG OF SAMPLE PUMP REPLACEMENT, 5500sc AMC	\$	635.98
9713400	KIT, REAGENT TRAY REPLACEMENT, 5500/9610	\$	229.34
9725500A	STIR BAR MAGNET, ROD, 1/8" X 3/8", NOT FDA CONTROLLED	\$	7.52
9737500	ASSY, CASE, WORKSTATION, SL1000	\$	252.46
9737600	PLATE, ADAPTER, SL1000	\$	110.43
9738000	CLAMP, TUBING, 1/4" x 0.04" x 0.13"	\$	20.43
9744300A	FTG, 1/4" STEM X 1/8" PTC, NOT FDA CONTROLLED	\$	16.00
9759000	ASSY., PACK, CHEMKEY, 25 PIECE, pH	\$	58.90
9770700	CLEANING SWAB, 3/8" FOAM TIP, 6-1/8" L, 50/PK	\$	124.11
9775100	CLEANING SWAB, 3/8" FOAM TIP, 6-1/8" L, 10/PK	\$	55.31
9780400	TNT UV REAGENT SET	\$	439.04
9790200	ee ULTRA LOW RANGE CL17sc, STANDPIPE INSTALL KIT	\$	12,835.15
9790300	ee ULTRA LOW RANGE CL17sc, PRESS.REG INSTALL KIT	\$	12,835.15
9790400	ee ULTRA LOW RANGE CL17sc, STANDPIPE INSTALL KIT, EU	\$	10,927.22
9790500	ee ULTRA LOW RANGE CL17sc, PRESS.REG INSTALL KIT, EU	\$	10,927.22
9790600	MAIN BOARD, PROGRAMMED, ULTRA LOW RANGE CL17sc	\$	1,485.87
9790700	KIT, DOOR REPLACEMENT, ULTRA LOW RANGE CL17sc	\$	256.32
97909	SULFURIC ACID CONC, ACS 2.5L	\$	229.46
97909-LM	SULFURIC ACID CONC, ACS 2.5L	\$	229.46
9791100	aa KTO: REAGENT SET, ULTRA LOW RANGE CHLORINE TOTAL CL17	\$	180.31
9791100-LM	aa KTO: REAGENT SET, ULTRA LOW RANGE CHLORINE TOTAL CL17	\$	180.31
9791200	PERFORMANCE CHECK KIT, LOW-LEVEL CHLORINE	\$	333.41
9791700	KTO: ULTRA LOW RANGE CL17sc, STANDPIPE, w/sc200	\$	20,331.96
9791800	KTO: ULTRA LOW RANGE CL17sc, PRES.REG. w/sc200	\$	20,331.96
9791900	KTO: ULTRA LOW RANGE CL17sc, STANDPIPE,TOTAL,sc200	\$	20,524.68
9792000	KTO: ULTRA LOW RANGE CL17sc, PRES.REG,TOTAL,sc200	\$	20,524.68
9792500	KTO: ULTRA LOW RANGE CL17sc, STANDPIPE, w/sc4500	\$	20,351.23
9792600	KTO: ULTRA LOW RANGE CL17sc, PRES.REG. w/sc4500	\$	20,351.23
9792700	KTO: ULTRA LOW RANGE CL17sc,STANDPIPE,TOTAL,sc4500	\$	20,563.22

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

9792800	KTO: ULTRA LOW RANGE CL17sc, PRES.REG,TOTAL,sc4500	\$	20,563.22
9793300	KTO: ULTRA LOW RANGE CL17sc,STANDPIPE,TOTAL	\$	11,081.40
9793400	KTO: ULTRA LOW RANGE CL17sc, PRESSURE REG,TOTAL	\$	11,081.40
97949	SULFURIC ACID CONC, ACS 500ML	\$	121.77
97949-LM	SULFURIC ACID CONC, ACS 500ML	\$	121.77
98199	DISSOLVED OXYGEN 1 PP PK/100	\$	67.02
98199-LM	DISSOLVED OXYGEN 1 PP PK/100	\$	67.02
98299	DISSOLVED OXYGEN 2 PP PK/100	\$	67.02
98299-LM	DISSOLVED OXYGEN 2 PP PK/100	\$	67.02
983-030	ZERO POWDER	\$	98.09
98399	BUFFER, CITRATE PWD PLWS PK/100	\$	150.13
98399-LM	BUFFER, CITRATE PWD PLWS PK/100	\$	150.13
98499	SODIUM PERIODATE PP 25ML PK/100	\$	101.35
98499-LM	SODIUM PERIODATE PP 25ML PK/100	\$	101.35
98768	DISSOLVED OXYGEN 3 PK/25	\$	26.19
98768-LM	DISSOLVED OXYGEN 3 PK/25	\$	26.19
98799	DISSOLVED OXYGEN 3 PWD PW PK/100	\$	99.73
98799-LM	DISSOLVED OXYGEN 3 PWD PW PK/100	\$	99.73
990100	COLOR DISC 2-PIECE, WIDE RANGE pH	\$	198.50
990200	COLOR DISC 2 PC, CHLORINE-DPD	\$	198.50
991-545	GSKT FILTER DIA 63-79	\$	95.40
992004	PSL SPHERES 2.02UM	\$	3,391.87
992009	PSL SPHERES 10UM	\$	3,391.87
992010	PSL SPHERES 15.0UM	\$	3,391.87
992019	PSL SPHERES 5.0UM	\$	3,391.87
99949	WASH WATER BUFFER SOLN, 500ML	\$	74.08
99949-LM	WASH WATER BUFFER SOLN, 500ML	\$	74.08
9H0626	NUT, 3/8 NPT	\$	59.55
9H1000	U-BOLT, 5/16-18 UNC 2A THDS, 316SS	\$	21.39
9H1010	CLAMP,HINGED TOGGLE	\$	501.07
9H1097	RELEASE HINGE,PULL-PIN,LE	\$	129.32
9H1098	RELEASE HINGE,PULL-PIN,RG	\$	129.32
9H1108	SET SCREW,SOC,CUPPNT,6-32	\$	5.40
9H1132	CLAMP,HINGED,304 S.S.,2"	\$	695.72
9H1285	CLEVIS PIN,.25"DIAx3.0",3	\$	113.13
9H1285-040	CLEVIS PIN,.25"DIA.X4.0 3	\$	33.34
9H1310	TEE,SANITARY CLAMP,SS316L	\$	1,420.35
9H1327	GASKET, EPDM, 2"	\$	112.93
9H1394	LANYARD,S.S. CABLE,10"	\$	21.39
9H1442	GASKET,SILICONE,1 1/2"	\$	65.14
9P2104-108	TEE,RED KYNAR SCH.80 THD	\$	713.06
9P2109-122-1	BUSH.NAT.KNR.80,THD.3/4 X	\$	316.06
A23749	Mirror, Lamp	\$	1,420.35
A23750	Mirror, Grating	\$	1,293.15
A23753	Mirror, Sample Compartment	\$	4,394.02
A23754	MIRROR, FLAP	\$	944.33
A23755	Beam Splitter Ref. Element	\$	946.26



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

A23756	Holder for Carousel Adapter	\$	244.75
A23757	CUVETTE, DRIVE RING, ROUND	\$	501.07
A23759	Enclosure Foot	\$	31.22
A23761	Biconvex Lens,	\$	321.84
A23762	Alignment Tool for Optics	\$	108.89
A23763	KNURLED SCREW FOR CUVETTE ADAPTERS	\$	343.04
A23766	DUST FILTER PAD	\$	15.03
A23772	FUSE, T2A H-250V, 5MM X 20MM	\$	18.12
A23773	Power Supply Socket	\$	846.04
A23777	TRANSPONDER COIL	\$	104.26
A23778	LAMP, TUNGSTEN, DR 6000	\$	851.82
A23791	Temperature Sensor Lamp Cover	\$	106.38
A23792	LAMP, DEUTERIUM	\$	5,376.89
A23794	Carousel Motor	\$	695.72
A23797	MODULE, SIPPER PUMP	\$	6,513.94
A23802	CKT BD ASSY, CONNECTION, SIPPER MOD	\$	961.67
A23803	TOP ENCLOSURE F/PELTIER TEMP MODULE	\$	568.52
A23804	CELL HOLDER 1CM, INCLUDES THERMISTOR	\$	2,543.90
A23805	COVER F/PELTIER TEMPERATURE MODULE	\$	809.42
A23806	PLATE SPRING, CELL HOLDER PELTIER	\$	185.59
A23809	ENCLOSURE, BOTTOM, PELTIER TEMP MOD	\$	373.88
A23816	PELTIER TEMPERATURE ELEMENT	\$	1,287.37
A23817	Plate Spring Turbidity and	\$	66.10
A23818	KNURLED SCREW F/TURB MICRO CELL HLDR	\$	774.73
A24209	160 L Flow-Through Cell	\$	3,044.98
A25B400	KIT MAINTENANCE ST-02B/(A25B270)	\$	4,509.65
A25B410	KIT MAINTENANCE ST-05B/(A25B280)	\$	4,509.65
A25B710	BOTTLE, GLASS, BROWN, 1 LITER	\$	720.77
A46384	DISPLAY DU 720 GENERAL PURPOSE, EMC	\$	5,376.89
A46389	DISPLAY DU730 PROTEIN ANALYZER, EMC	\$	4,702.37
A46390	FIXED SURFACE OF SLIDING LID, DU700	\$	74.97
A46391	ENCLOSURE TOP DU700, EMC	\$	3,430.42
A46392	DISPLAY FRAME DU720	\$	383.51
A46394	DISPLAY FRAME DU730, PA, EMC	\$	93.47
A46395	SWITCHED MODE POWER SUPPLY, EMC	\$	2,582.45
A46396	MOTHERBOARD, EMC	\$	10,137.07
A46397	PERIPHERAL BOARD, EMC	\$	2,967.89
A46399	INTERFACE BOARD, EMC	\$	260.17
A46400	REFERENCE ELEMENT BOARD, EMC	\$	2,274.10
A46401	MEASUREMENT ELEMENT BOARD, EMC	\$	2,678.81
A46402	LAMP CONNECTOR BOARD, EMC	\$	2,235.55
A46403	DISPLAY CONNECTOR BOARD, EMC	\$	383.51
A50T520	BUBBLER FOR CNC CELL	\$	697.65
A94L136	CL136 CABLE S7/COAX	\$	531.91
AHA033NPT	si PROBE ADAPTER STRAIGHT 1 1/2FNPT	\$	516.49
AHA045	Adapter probe 90° elbow	\$	377.73
AHA045NPT	PROBE ADAPTER 90° ELBOW, 1 1/2 FNPT	\$	526.13



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

ALB-100C	ACCULIGHT BASIC TEST KIT, 100 TESTS	\$	2,439.11
ALBSURFACE-100C	AccuLight Basic Surface Complete, 100 tests	\$	1,246.66
AP0001.AT1102	pH/Alkalinity in water application pack	\$	2,293.37
AP0002.AT1102	pH/Alkalinity/Conductivity in water application pack	\$	3,584.59
AP0003.AT1112	Ca & Mg Hardness (ISE) in water application pack	\$	5,068.54
AP0004.AT1112	Total Hardness (Color) in water application pack	\$	10,406.88
AP0005.AT1222	pH/Alkalinity/Hardness (ISE) in water application pack	\$	7,072.82
AP0006.AT1102	FOS/TAK (BIOGAS) application pack	\$	6,610.30
AP0007.AT1122	Free & Total Chlorine, Chlorine Dioxide, Sulfite (AutoCAT) application pack	\$	2,736.62
AP0008.AT1102	pH, Total Acidity in Food & Beverage application pack	\$	3,044.98
AP0009.AT1112	Chlorides in water application pack	\$	2,351.18
AP0010.AT1112	Salt in food products application pack	\$	2,351.18
AP0011.AT1222	H, Total Acidity & Chlorides in Food & Beverage application pack	\$	5,261.26
AP0012.AT1122	Free & Total SO2 in wine application pack	\$	2,736.62
AP0013.AT1222	pH, Total Acidity, Free & Total SO2 in wine application pack	\$	5,588.88
AP0014.KF1121	Free Moisture Content (Karl Fischer) application pack	\$	3,006.43
AP0015.AT1102	TAN (Total Acid Number) in Petrochem application pack	\$	4,586.74
AP0016.AT1102	TBN (Total Base Number) in Petrochem application pack	\$	4,355.47
AP0017.AT1112	R-SH (Thiol) in Petrochem application pack	\$	2,351.18
AP0018.AT1102	Br2/I2 index in Petrochem application pack	\$	2,659.54
AP0019.AT1112	Generic RedOx Titration i=0 application pack	\$	3,623.14
AP0020.AT1112	Generic RedOx Titration i>0 application pack	\$	3,584.59
AP0021.AT1102	Generic Acid Base Titration application pack	\$	3,044.98
AP0025.AT1102	pH, Total Acidity in Soft Drink, Beer & Kombucha Application Pack, Single 25mL Syringe	\$	5,164.90
AP0026.AT1222	pH, Total Acidity in Soft Drink, Beer & Kombucha Application Pack, Double 25mL Syringe	\$	5,685.24
APPS10011887	Inc, Pump, flexi hose and flanges for Scrubber Monitor cabinets	\$	9,250.56
APPS10011888	Debubbler including tubing for Scrubber Monitor cabinets	\$	8,595.31
AQ14004/09	OPTIMA FLOW MONITOR (LOGG	\$	12,526.80
AS1000.82.12150	Sampler TitraLab AT1000 series, 12 beakers, 50 & 150 mL	\$	23,434.75
AS1000.82.20090	Sampler TitraLab AT1000 series, 20 beakers, 90 mL	\$	23,434.75
AS1000.82.30050	Sampler TitraLab AT1000 series, 30 beakers, 50 mL	\$	23,434.75
AS1000.97.12150	Sampler TitraLab AT1000 series, 12 beakers, 50 & 150 mL	\$	23,434.75
AS1000.97.20090	Sampler TitraLab AT1000 series, 20 beakers, 90 mL	\$	23,434.75
AS1000.97.30050	Sampler TitraLab AT1000 series, 30 beakers, 50 mL	\$	23,434.75
AS1000.98.12150	Sampler TitraLab AT1000 series, 12 beakers, 50 & 150 mL	\$	23,434.75
AS1000.98.20090	Sampler TitraLab AT1000 series, 20 beakers, 90 mL	\$	23,434.75
AS1000.98.30050	Sampler TitraLab AT1000 series, 30 beakers, 50 mL	\$	23,434.75
AT1102.97	TitraLab AT1000, 1 Syringe	\$	9,481.82
AT1112.97	TitraLab AT1000, 1 Syringe, 1 Pump	\$	10,734.50
AT1122.97	TitraLab AT1000, 1 Syringe, 2 Pumps	\$	11,948.64
AT1222.97	TitraLab AT1000, 2 Syringe, 2 Pumps	\$	13,509.67
ATS100	SULFURIC ACID TO STANDARDIZE 0.100 N BASE PK/12	\$	216.81
B90X500	FINISHING ABRASIVE BANKS	\$	281.37
BBP065	PIPET, 1.0 - 5.0 ML	\$	1,040.69



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

BBP068	PIPET TIPS, 75/PK (FOR BBP065)	\$	109.46
BBP078	Pipet, 0.1 - 1.0 mL	\$	1,040.69
BBP079	PIPET TIPS, 100/PK (FOR BBP078)	\$	98.29
BCZ822	REAGENT A (additional) FOR LCW823/840/893	\$	475.18
BRO075	oo MOUNTING PIPE WITH LARGER FLANGE	\$	470.24
BRO078	PIPE, JOINT, 3.35m	\$	716.92
BYF154	PACKAGING FOR DR3900	\$	238.97
BYF177	Packaging - TU5200	\$	114.67
BYF187	Verpackung TitraLab AT1000	\$	88.27
BYF197	SC4x00 Packaging	\$	162.46
CARTON-AC	AccuChlor2 Shipping Pack	\$	343.04
CARTON-ACEVAL	AC2000P Eval. Protection	\$	472.16
CBL2!!	CABLE LENGTH	\$	206.21
CDC40101	ee CONDUCTIVITY PROBE, STD, w/1m CABLE	\$	1,559.10
CDC40103	ee CONDUCTIVITY PROBE, STD, w/3m CABLE	\$	1,730.63
CDC40105	ee CONDUCTIVITY PROBE, RUG, w/5m CABLE	\$	3,122.06
CDC40110	ee CONDUCTIVITY PROBE, RUG, w/10m CABLE	\$	3,507.50
CDC40115	ee CONDUCTIVITY PROBE, RUG, w/15m CABLE	\$	3,603.86
CDC40130	ee CONDUCTIVITY PROBE, RUG, w/30m CABLE	\$	3,931.49
CG2ENT	CITYGUARD BASE TO ENTERPRISE, ADVANCE SYSTEM	\$	97,150.15
CGMAINTADDEM	Support of EA Add'l Event Mntr	\$	957.82
CGMAINTBASE	Support w/Updates to Base	\$	3,892.94
CM_INT-STRUP	Interface Setup and Installation	\$	9,693.82
CM_JOB-SERVICES	JOBCal Services	\$	649.47
CRI5342T/10.99	PH ELECTR 5342, 10M CABLE, PT1000	\$	1,210.28
D3422A1	Digital Conductivity Sensor, 0.05 K	\$	2,177.74
D3422A2	Digital Conductivity Sensor, 0.05 K	\$	2,563.18
D3422B3	Digital Conductivity Sensor, 0.50 K	\$	2,177.74
D3422B4	Digital Conductivity Sensor, 0.50 K	\$	2,698.08
D3422C3	Digital Conductivity Sensor, 1.00 K	\$	2,177.74
D3422D3	Digital Conductivity Sensor, 5.00 K	\$	2,177.74
D3422D4	Digital Conductivity Sensor, 5.00 K	\$	2,698.08
D3422E3	Digital Conductivity Sensor, 10.0 K	\$	2,177.74
D3433B8	KTO:DIG CONT COND SENS GEN PUR RYTON	\$	1,877.09
D3433E8	KTO: DIG CONT COND SENS GEN PUR RYTON	\$	1,877.09
D3444B8	KTO: DIG CONT COND SENS BOILER SST	\$	3,314.78
D3444D8	KTO: DIG CONT COND SENS BOILER SST	\$	3,314.78
D3455A6	KTO:DIG CONT COND SENS SANITARY SST	\$	3,006.43
D3455C7	KTO:DIG CONT COND SENS SANITARY SST	\$	3,141.34
D3455E7	KTO:DIG CONT COND SENS SANITARY SST	\$	3,141.34
D3705E2T	KTO:DIG ELCLESS COND SENS,SAN,POLY	\$	3,892.94
D3705E2T.99	Digital Gateway, Sanitary, Poly	\$	3,892.94
D3706E2T	KTO:DIG ELCLESS COND SENS,SAN, KYNAR	\$	4,432.56
D3706E2T.99	Digital Gateway, Sanitary, Kynar	\$	4,432.56
D3708E2T	KTO:DIG ELCLESS COND SENS,SAN,TEFLN	\$	4,914.36
D3708E2T.99	Digital Gateway, Sanitary, Teflon	\$	4,914.36
D3725E2T	KTO:DIG ELCLESS COND SENS,CONV,POLY	\$	3,546.05



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

D3725E2T.99	Digital Gateway, Convertible, Poly	\$	3,546.05
D3726E2T	KTO: DIG ELECLESS COND SENS, CONV, KYNA	\$	4,027.85
D3726E2T.99	Digital Gateway, Convertible, Kynar	\$	4,027.85
D3726E2T20S	Digital Inductive Conductivity Probe	\$	4,143.48
D3727E2T	KTO: DIG ELECLESS COND SENS, CONV, PEEK	\$	4,355.47
D3727E2T.99	Digital Gateway, Convertible, PEEK	\$	4,355.47
D3728E2T	KTO: DIG ELECLESS COND SENS, CONV, TEFL	\$	4,509.65
D3728E2T.99	Digital Gateway, Convertible, Teflon	\$	4,509.65
D71T006	EQ BOOK, TTL860/865	\$	601.29
D71T018	TITRALAB965 QUALIFICATION BOOKLET	\$	456.75
D71T019	EQ, TTL980, Doc-E Booklet	\$	491.44
DIS-CT1255-50	12X55MM TEST TUBES	\$	48.23
DIS-CT15-50R	dd 15mL Centrifuge Tubes, PP, 50/Rack	\$	110.81
DIS-CTTC-96	dd Centrifuge Tube, 2mL, Tether Cap, 96/Bag	\$	79.21
DISPT1000FIL96R	dd Pipet Tip, Universal ZAP Tips, 100-1000 uL Filtered, Sterile, 96/Rack	\$	47.41
DISPT200FIL96R	dd Pipet Tip, Universal ZAP Tips, 5-200 uL Filtered, Sterile, 96/Rack	\$	46.25
DISPT20FIL96R	dd Pipet Tip, Universal ZAP Tips, 2-20 uL Filtered, Sterile, 96/Rack	\$	47.41
DIS-SFQG-25	QGA SYRINGE FILTERS, BAG OF 25	\$	252.46
DIS-ST8-12	dd PCR Strip Tubes, 12/Bag	\$	79.21
DIW	CABLE, 3 METER D.O.GLI/	\$	1,534.05
DM_INT-STRTUP	Interface Setup + Installation	\$	12,642.43
DM_WIMS-HRLY-DEV	IIM Programming Services	\$	1,214.14
DM_WIMSHRLY-EXP	IIM Expert Services	\$	1,144.76
DM_WIMS-HRLY-ON	IIM Online Setup Services	\$	620.56
DM_WIMS-ON FRST	Online 1st Named User WW	\$	10,214.16
DM_WIMS-ON-50MB	WIMS Online Additional 50MB Disc Space	\$	1,165.96
DM_WIMS-ON-ADDL	WIMS Online Additional Named User	\$	2,929.34
DM_WIMS-ON-DISC	Online Addl 200MB Disc Sp	\$	3,758.04
DM_WIMS-ON-FACIL	WIMS Online Additional Facility	\$	580.09
DM_WIMS-ON-LCAL	WIMS Online LAB Cal per facility	\$	7,265.54
DM_WIMS-ON-SCDA	WIMS Online Standard SCADA Interface	\$	4,394.02
DM_WIMS-OSTRN-1DAY	Onsite Services One Day û Project Manager	\$	17,537.52
DM_WIMS-OSTRN-2DAY	Onsite Services Two Day (concurrent) û Project Manager	\$	23,974.37
DM_WIMS-OSTRN-3DAY	Onsite Services Three Dayá (concurrent) û Project Manager	\$	30,449.76
DM_WIMS-OSTRN-4DAY	Onsite Services Four Dayá (concurrent) û Project Manager	\$	36,867.34
DM_WIMS-OSTRN-5DAY	Onsite Services Five Day (concurrent) û Project Manager	\$	43,362.00
DM_WIMS-SERVICES	WIMS Services	\$	649.47
DMWIMS-SPT-BOD	WIMS BOD MGR Support & Maint	\$	479.87
DMWIMS-SPTDBSQL	Basic Support, SQL Server support	\$	4,143.48
DMWIMS-SPTFACBA	Basic Support, MU addl facility	\$	847.97
DMWIMS-SPTINT	Basic Support, Interfaces	\$	2,042.83
DMWIMS-SPTLCBA	Basic Support, MU LAB Cal	\$	3,468.96
DMWIMS-SPT-MU	Base SPT MU, 1USER, 1FACIL	\$	3,314.78
DMWIMS-SPTUSBAS	Basic Support, MU addl user	\$	657.18
DMWIMS-ON-1STLT	Online 1st Named User WW 50MB	\$	7,265.54

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

DMWMS-VIEW	WIMS Online Viewer û monthly	\$	229.34
DOC012.97.80101	HACH, MP4 QUICK REFERENCE	\$	39.89
DOC012.97.90328	QUICK REFERENCE GUIDE, SIP 10	\$	92.51
DOC012.99.80041	2100Q QUICK START GUIDE	\$	85.95
DOC012.99.80096	CARD, MODULE INFORMATION, sc200	\$	39.51
DOC013.99.80240	UNPACKING GUIDE, 5500sc, SiO2	\$	42.21
DOC013.99.80241	UNPACKING GUIDE, 5500sc, PO43 LR	\$	42.21
DOC013.99.80242	UNPACKING GUIDE, 5500sc, PO43 HR	\$	42.21
DOC013.99.80246	UNPACKING GUIDE, POLYMETRON 9610sc, SiO2	\$	42.21
DOC013.99.80247	UNPACKING GUIDE, POLYMETRON 9611sc, PO43, LR	\$	42.21
DOC013.99.80248	UNPACKING GUIDE, POLYMETRON, 9611sc, PO43, HR	\$	42.21
DOC022.53.00619	Bedienungsanleitung - DRB200	\$	356.53
DOC022.53.00782	DR3800 USER MANUAL, ENGLISH	\$	208.14
DOC022.53.80225	PROCEDURE MANUAL, HYDRAULIC FRACTURING	\$	0.00
DOC022.53.90072	BODTRAK II MANUAL	\$	43.75
DOC022.91.00751	MANUAL, USER, HQD SERIES, NA FRENCH	\$	28.14
DOC022.92.00751	MANUAL, USER, HQD SERIES LA SPANISH	\$	27.94
DOC022.97.80017	BASIC USER MANUAL, HQd PORTABLE METER, ROW	\$	115.25
DOC022.97.80041	BASIC USER MANUAL, 2100Q/2100Qis, ROW	\$	105.61
DOC022.97.80116	USER MANUAL, BASIC, HQd BENCHTOP, HACH, ROW	\$	115.25
DOC022.97.90324	BASIC USER MANUAL, DR3900	\$	215.85
DOC022.97.90367	BUM DR 6000 UV VIS SPEKTRALPHOTOMETER US	\$	82.87
DOC022.97.90372	MANUAL DRINKING WATER APPLICATION DR 6000	\$	82.87
DOC022.97.90373	MANUAL BREWERY APPLICATION DR 6000	\$	82.87
DOC022.97.90402	MANUAL ENZYMATIC APPLICATION DR 6000	\$	82.87
DOC022.97.90443	BUM DR 6000EDU UV VIS SPECTROPHOTOMETER	\$	85.57
DOC022.98.00619	MANUAL DRB200 EU VERSION	\$	346.90
DOC023.53.80040	aa USER MANUAL, SC200 CONTROLLER, ENGLISH, EU/ROW	\$	72.66
DOC023.53.90007	vv MANUAL, USER, SC1000, ENGLISH, ROW	\$	157.45
DOC023.53.90032	SODIUM ANALYZER USER MANUAL	\$	128.35
DOC023.53.90137	USER MANUAL FOR AN-ISE SC (HACH)	\$	118.91
DOC023.53.90154	USER MANUAL TSS sc, ENGLISH	\$	165.55
DOC023.53.90171	USER MANUAL TSS EX1 sc, ENGLISH	\$	165.55
DOC023.54.00051	ENGLISH 9180sc OPERATION MANUAL	\$	142.61
DOC023.54.03232	Bedienungsanleitung - SOLITAX sc	\$	143.19
DOC023.80.90154	USER MANUAL TSS sc, CHINESE	\$	165.55
DOC023.80.90171	USER MANUAL TSS EX1 sc, CHINESE	\$	165.55
DOC023.81.90154	USER MANUAL TSS sc, JAPANESE	\$	165.55
DOC023.81.90171	USER MANUAL TSS EX1 sc, JAPANESE	\$	165.55
DOC023.84.90154	Bedienungsanleitung TSS sc (KR)	\$	165.55
DOC023.84.90171	USER MANUAL TSS EX1 sc, KOREAN	\$	165.55
DOC023.90.90154	USER MANUAL TS sc, PORT FOR BRAZIL	\$	165.55
DOC023.90.90171	USER MANUAL TSS EX1 sc, PORT FOR BRAZIL	\$	165.55
DOC023.91.90137	USER MANUAL FOR AN-ISE SC (CA)	\$	118.91
DOC023.91.90154	USER MANUAL TSS sc, FRENCH CANADA	\$	165.55
DOC023.91.90171	USER MANUAL TSS EX1 sc FRENCH CANADA	\$	165.55
DOC023.92.90154	USER MANUAL TSS sc, SPANISH FOR S. AMERICA	\$	165.55



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

DOC023.92.90171	USER MANUAL TSS EX1 sc SPANISH FOR S. AMERICA	\$	165.55
DOC023.97.80040	BASIC USER MANUAL, SC200 CONTROLLER, ROW	\$	114.86
DOC023.97.80049	INSTRUCTIONS, USER, FLOW MODULE, sc200, HACH	\$	23.70
DOC023.97.80051	INSTRUCTIONS, USER, MODBUS MODULE, sc200, HACH	\$	55.89
DOC023.97.80076	BASIC USER MANUAL, pH ANALOG DIFF pH/ORP SENSORS, US/ROW	\$	63.02
DOC023.97.80196	USER MANUAL, SAMPLE HEATER, ROW	\$	152.25
DOC023.97.80240	INSTALLATION MANUAL, 5500sc, SILICA, ROW	\$	146.85
DOC023.97.80241	OPERATION MANUAL, 5500sc, PO43, LR, ROW	\$	146.85
DOC023.97.80242	OPERATION MANUAL, 5500sc, PO43, HR, ROW	\$	146.85
DOC023.97.80243	MAINTENANCE MANUAL, 5500sc, ROW	\$	105.03
DOC023.97.80246	INSTALLATION MANUAL, POLYMETRON 9610sc, SiO2, ROW	\$	146.85
DOC023.97.80247	OPERATION MANUAL, POLYMETRON 9611sc, PO43, LR, ROW	\$	146.85
DOC023.97.80248	OPERATION MANUAL, POLYMETRON 9611sc, PO43, HR, ROW	\$	146.85
DOC023.97.80249	MAINTENANCE MANUAL, POLYMETRON 9610sc/9611sc, ROW	\$	105.03
DOC023.97.80267	OPERATION MANUAL, 5500sc, SILICA, ROW	\$	146.85
DOC023.97.80268	INSTALLATION MANUAL, 5500sc, PO43, LR/HR, ROW	\$	146.85
DOC023.97.80269	OPERATION MANUAL, POLYMETRON 9610sc, SiO2, ROW	\$	146.85
DOC023.97.80270	INSTALLATION MANUAL, POLYMETRON, 9611sc, PO43, LR/HR, ROW	\$	146.85
DOC023.97.80282	BASIC USER MANUAL, CA610, ROW	\$	146.85
DOC023.97.80473	INSTALLATION MANUAL, 5500sc, AMMONIA MONOCHLORAMINE, ROW	\$	-
DOC023.97.80474	OPERATIONS MANUAL, 5500sc, AMMONIA MONOCHLORAMINE, ROW	\$	-
DOC023.97.80475	MAINTENANCE MANUAL, 5500sc, AMMONIA MONOCHLORAMINE, ROW	\$	-
DOC023.97.80482	USER MANUAL, HAZARD. LOC, 1735 TRO, (NK)	\$	30.06
DOC023.98.80240	INSTALLATION MANUAL, 5500sc, SILICA, EU	\$	146.85
DOC023.98.80241	OPERATION MANUAL, 5500sc, PO43, LR, EU	\$	146.85
DOC023.98.80242	OPERATION MANUAL, 5500sc, PO43, HR, EU	\$	146.85
DOC023.98.80243	MAINTENANCE MANUAL, 5500sc, EU	\$	105.03
DOC023.98.80246	INSTALLATION MANUAL, POLYMETRON 9610sc, SiO2, EU	\$	139.72
DOC023.98.80247	OPERATION MANUAL, POLYMETRON 9611sc, PO43, LR, EU	\$	146.85
DOC023.98.80248	OPERATION MANUAL, POLYMETRON 9611sc, PO43, HR, EU	\$	146.85
DOC023.98.80249	MAINTENANCE MANUAL, POLYMETRON 9610sc/9611sc, EU	\$	105.03
DOC023.98.80267	OPERATION MANUAL, 5500sc, SILICA, EU	\$	146.85
DOC023.98.80268	INSTALLATION MANUAL, 5500sc, PO43, LR/HR, EU	\$	146.85
DOC023.98.80269	OPERATION MANUAL, POLYMETRON 9610sc, SiO2, EU	\$	146.85
DOC023.98.80270	INSTALLATION MANUAL, POLYMETRON, 9611sc, PO43, LR/HR, EU	\$	146.85
DOC026.53.00117	SONATAX sc USER MANUAL	\$	174.60
DOC026.53.00741	POLYMETRON 9245 SODIUM ANAL MANUAL	\$	130.66
DOC026.80.00750	MANUAL, INSTR, DR/5000, CHINESE	\$	223.56
DOC026.81.00750	MANUAL, INSTR, DR/5000, JAPANESE	\$	223.56
DOC026.84.00750	MANUAL, INSTR, DR/5000, KOREAN	\$	223.56
DOC026.90.00750	MANUAL, INSTR, DR/5000, PORTUGUESE	\$	223.56



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

DOC026.91.00750	MANUAL, INSTR, DR/5000, FRENCH	\$	223.56
DOC026.92.00750	MANUAL, INSTR, DR/5000, SPANISH	\$	223.56
DOC082.97.80041	2100Q CD MANUAL, HACH	\$	128.35
DOC082.97.90444	CD DR 6000EDU UV VIS SPECTROPHOTOMETER	\$	200.43
DOC082.98.00670	MANUAL, PROCEDURE, DR/5000	\$	91.35
DOC082.98.80041	2100Q DOCUMENTATION CD-ROM	\$	122.76
DOC082.98.90072	CD ROM, BODTRAK II	\$	31.61
DOC086.53.01000	vv Water Analysis Handbook CD-ROM	\$	0.77
DOC272.52.80171	HQ440d Lab pH Meter Equipment Qualification	\$	443.26
DOC272.52.80172	HQ440d Lab Conductivity Meter Equipment Qualification	\$	443.26
DOC272.53.90575	Equipment Qualification - TL2350 & TL2360	\$	200.43
DOC272.53.90580	QUALIFICATION INSTRUCTIONS, IQ-OQ, TL23x0, Pharma Customers	\$	813.28
DOC273.97.80261	USER INSTRUCTIONS, 96XX EMPTY BOTTLE KIT, ROW	\$	29.68
DOC273.97.80351	USER INSTRUCTIONS, 96XX REAGENT KIT, ROW	\$	8.09
DOC273.97.80565	USER INSTRUCTIONS, 5500 REAGENT KIT, ROW	\$	19.08
DOC273.97.80618	USER INSTRUCTIONS, CALIBRATION VERIFICATION KIT, CL17sc	\$	16.19
DOC273.98.80351	USER INSTRUCTIONS, 96XX REAGENT KIT, EU	\$	39.31
DOC273.98.80565	USER INSTRUCTIONS, 5500 REAGENT KIT, EU	\$	19.08
DOC273.99.80258	USER INSTRUCTIONS, ANALYZER, WALL, PANEL, TABLE MOUNT	\$	42.21
DOC273.99.90155	INSTRUCTION SHEET TSS sc, TANK MOUNT	\$	165.55
DOC273.99.90156	Montageanleitung TSS(EX1)sc LZU301	\$	125.27
DOC316.53.01302	PROCEDURE, CHLORINE CAL VERIFICATION	\$	42.40
DOC326.98.00010	INST SHT, KIT, WTR TRT COMB CN-39WR	\$	13.88
DOC326.98.00050	INSTRUCTIONS, TEST KIT, IRON 10mg/L Fe	\$	19.27
DPC1R1A	KTO:DIG PH SENS,COMBO,CONV,RYTON,TEM	\$	1,753.75
DPC1R1N	KTO:DIG PH SENS,COMBO,CONV,RYTON	\$	1,445.40
DPC1R2A	KTO:DIG PH SENS,COMBO,CONV,RYTON,TEM	\$	1,865.53
DPC1R2N	KTO:DIG PH SENS,COMBO,CONV,RYTON	\$	1,753.75
DPC1R3A	KTO:DIG PH SENS,COMBO,CONV,RYTON,TEM	\$	2,100.65
DPC2K1A	KTO:DIG PH SENS,COMBO,INSER,KYNAR,TE	\$	1,985.02
DPC2K2A	KTO:DIG PH SENS,COMBO,INSER,KYNAR,TE	\$	2,100.65
DPC3K2A	KTO:DIG PH SENS,COMBO,SANIT,KYNAR,TE	\$	2,524.63
DPD1P1	Digital pH Sensor, PEEK, Convertible	\$	3,738.77
DPD1P1.1	DIGITAL pH SENSOR, PEEK, CONVERTIBLE Non-HEV	\$	3,796.58
DPD1P1.99	Digital pH Sensor, PEEK, Convertible	\$	3,738.77
DPD1P1.99M	DIGITAL pH SENSOR, PEEK, CONVERTIBLE, 1M	\$	3,565.32
DPD1P1Q	DpHD Sensor, Convertible, Quick Connect	\$	3,738.77
DPD1P1Q06	DpHD Sensor, Convertible, 6 Foot Cable	\$	3,738.77
DPD1P3	Digital pH Sensor, PEEK, Convertible	\$	4,490.38
DPD1R1	Digital pH Sensor,Ryton, Convertible	\$	3,526.78
DPD1R1.1	Digital pH Sensor,Ryton, Convertible Non-HEV	\$	3,565.32
DPD1R1.99	Digital pH Sensor,Ryton, Convertible	\$	3,526.78
DPD1R1-WDMP	pH Sensor, Ryton, WDMP Mounting, Dig	\$	3,353.33
DPD1R3	Digital pH Sensor,Ryton, Convertible	\$	4,297.66
DPD1R3.1	Digital pH Sensor,Ryton, Convertible Non-HEV	\$	3,815.86
DPD1R3.99	Digital pH Sensor,Ryton, Convertible	\$	4,297.66

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

DPD2P1	Digital pH Sensor, PEEK, Insertion	\$	3,835.13
DPD2P1.1	Digital pH Sensor, PEEK, Insertion Non-HEV	\$	3,931.49
DPD2P1.99	Digital pH Sensor, PEEK, Insertion	\$	3,738.77
DPD3P1	Digital pH Sensor, PEEK, Sanitary	\$	3,989.30
DPD3P1.99	Digital pH Sensor, PEEK, Sanitary	\$	3,989.30
DPS1	Digital pH Sensor, Stainless	\$	3,989.30
DPS1.1	Digital pH Sensor, Stainless Non-HEV	\$	3,603.86
DPS1.99	Digital pH Sensor, Stainless	\$	3,989.30
DQGOM-25C	Differential QGOM Complete Test Kit, 25 Tests	\$	1,669.44
DR1900-01H	DR1900 SPECTRO PKG, HACH	\$	13,124.23
DR1900-01HNB	DR1900 SPECTRO PKG, HACH - NO BATTERY	\$	12,314.81
DR1900-02L	DR1900 SPECTRO PKG, LANGE	\$	12,430.44
DR2800-45	PROCEDURE MANUAL, DR 2800, CD-ROM	\$	172.10
DR2827-45	PROCEDURE MANL, DR2800/2700, CD-ROM	\$	172.10
DR3900TKN01	aa KTO: DR3900 TKN PKG SINGLE BLOCK	\$	22,586.78
DR3900TKN02	aa KTO: DR3900 TKN PKG DUAL BLOCK	\$	24,012.91
DRB200-01	REACTOR 115V 9X13 + 2x20 (1 BLOCK)	\$	3,391.87
DRB200-02	REACTOR 115V 21X13/4X20 (2 BLOCKS)	\$	4,914.36
DRB200-03	REACTOR 115V 15X13/15X13 (2 BLOCKS)	\$	4,914.36
DRB200-04	REACTOR 115V 12X13 + 8X20 (2 BLOCKS)	\$	4,914.36
DRB200-05	KTO:REACTOR 230V 9X13+2X20 (1 BLOCK)	\$	3,391.87
DRB200-06	KTO:REACTOR 230V 21X13/4X20 2 BLOCKS	\$	4,914.36
DRB200-07	KTO: REACTOR 230V 15X13/15X13 2 BLKS	\$	4,914.36
DRB200-08	KTO:REACTOR 230V 12X13+8X20 2 BLOCKS	\$	4,914.36
DRC1R5N	KTO: DIG ORP SENS,COMBO,CONV,RYTON	\$	1,528.27
DRC2K5N	KTO:DIG ORP SENS,COMBO,INSER,KYNAR	\$	1,746.04
DRD1P5	Digital ORP Sensor,PEEK, Convertible	\$	3,989.30
DRD1P5.1	Digital ORP Sensor,PEEK, Convertible Non-HEV	\$	3,835.13
DRD1P5.99	Digital ORP Sensor,PEEK, Convertible	\$	3,989.30
DRD1P6	Digital ORP Sensor,PEEK, Convertible	\$	4,297.66
DRD1P6.99	Digital ORP Sensor,PEEK, Convertible	\$	4,297.66
DRD1R5	Digital ORP Sensor,Ryton,Convertible	\$	3,738.77
DRD1R5.1	Digital ORP Sensor,Ryton,Convertible Non-HEV	\$	3,295.51
DRD1R5.99	Digital ORP Sensor,Ryton,Convertible	\$	3,738.77
DRD1R6	Digital ORP Sensor,Ryton,Convertible	\$	3,989.30
DRD2P5	Digital ORP Sensor,PEEK,Insertion	\$	4,182.02
DRD2P5.99	Digital ORP Sensor,PEEK, Insertion	\$	4,182.02
DRD2P6	Digital ORP Sensor,PEEK,Insertion	\$	4,394.02
DRD3P5	Digital ORP Sensor,PEEK,Sanitary	\$	4,471.10
DRS5	Digital ORP Sensor,Stainless	\$	3,989.30
DRS5.1	pHD sc, Differential ORP Digital Non-HEV	\$	3,488.23
DRS5.99	pHD sc, Differential ORP Digital	\$	3,989.30
DR-USB	USB SOFTWARE UPDATE, PHOTOMETER	\$	103.68
DSA-100C	DEPOSIT AND SURFACE KIT, 100 TEST	\$	4,372.34
DSA-25C	DEPOSIT AND SURFACE ANALYSIS COMPLETE, 25 TEST	\$	1,365.90
E3C1N-1008	aa vv Anderson Electrodeless Module	\$	2,293.37
E71T001	PHOTOCOLORIMETRIC PROBE	\$	5,261.26



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

EAK040	CAP, USB PORT	\$	12.72
EQP-4WAY	dd Test Tube Rack, 4-Way Rotating	\$	79.21
EQP-HEAT-GCEZ	dd Heater Block, GeneCount Q-16	\$	2,601.72
EQP-MAGRACK	dd Magnetic Rack	\$	142.61
EQP-P1000	dd Adjustable Pipettor, 1000uL	\$	950.11
EQP-P200	dd Adjustable Pipettor, 200uL	\$	950.11
EQP-PBM	PHOTONMASTER BLUETOOTH MODULE	\$	1,780.73
EQP-PBM-PAC	PHOTONMASTER + PBM EQUIPMENT SET	\$	16,111.39
EQP-PBM-UPGB	PHOTONMASTER + PBM UPGRADE SET (BASIC)	\$	3,044.98
EQP-PBM-UPGF	PHOTONMASTER + PBM UPGRADE SET (FULL)	\$	3,738.77
EYT033	BUSINESS DOLLY	\$	998.29
EZD117	O-Ring 2,84 x 2,62 - NBR 70	\$	27.56
EZD169	O-RING 70.00 X 5.00	\$	52.03
EZH051	SILICONE GREASE (AMTX,PHOSX,FILTX); KORASILON-/Baysilone-Paste	\$	97.71
EZH079	GREASE, O-RING 60G	\$	237.05
EZZ052	ALLEN WRENCH, 5MM, UV NITRATE	\$	33.34
EZZ062	SPONGE, SVI OPTIQUANT	\$	78.44
F_1002002.PN-B	Plastic collar for C1100	\$	426.39
F1-1001-101	REPAIR KIT-PP PIN-LOCK(F1	\$	529.98
F1-1001-101G	Generic Pin-Lock Repair Kit	\$	495.29
F1-1001-102	REPAIR KIT-PVDF PIN-LK(F1	\$	637.90
F1-1001-103	REPAIR KIT-INSERTION/FLOW	\$	637.90
F1-1003-B2	BOM-PLUG PVDF & VITN, FOR	\$	535.76
F1A11A1G	aa vv FINISHED SENSOR, PIN-LOCK,	\$	2,004.29
F1A11A1T	SENSOR, PADDLE FLOW, T-MNT, POLY	\$	1,907.93
F1A11A1T-1001	vv K = 100FT CABLE	\$	3,738.77
F1A11A1T20P	SENSOR, PADDLE FLOW, T-MNT, POLY, 20 FT	\$	1,925.27
F1A11A1T20S	SENSOR, PADDLE FLOW, T-MNT, POLY, 20 FT	\$	2,004.29
F1A11A1T25N	SENSOR, PADDLE FLOW, T-MNT, POLY, 25 FT	\$	2,967.89
F1A11A1T30	SENSOR, PADDLE FLOW, T-MNT, POLY, 30 FT	\$	3,025.70
F1A11A1T50	SENSOR, PADDLE FLOW, T-MNT, POLY, 50 FT	\$	3,237.70
F1A11B2T	aa SENSOR, PADDLE FLOW, T-MNT, PVDF	\$	2,409.00
F1A11B2T20S	aa SENSOR, PADDLE FLOW, T-MNT, PVDF, 20 FT	\$	2,524.63
F1A11B2T25N	SENSOR, PADDLE FLOW, T-MNT, PVDF, 25 FT	\$	3,526.78
F1A11B2T-SP	vv F1A11B2T Special Promotion Upgrade	\$	1,965.74
F1A12C1T	aa SENSOR, PADDLE FLOW, PIPE THD, SS SLEEVE	\$	4,740.91
F1A12C1T20P	SENSOR, PADDLE FLOW, PIPE THD, SS SLEEVE, 20 FT	\$	4,740.91
F1A12C1T40	SENSOR, PADDLE FLOW, PIPE THD, SS SLEEVE, 40 FT	\$	5,993.59
F1A12D1T	aa SENSOR, PADDLE FLOW, PIPE THD, BRASS SLEEVE	\$	4,008.58
F1A12D1T20P	SENSOR, PADDLE FLOW, PIPE THD, BRASS SLEEVE, 20 FT	\$	3,738.77
F1A12D1T30N	SENSOR, PADDLE FLOW, PIPE THD, BRASS SLEEVE, 30 FT	\$	4,875.82
F1A13C1T	aa SENSOR, PADDLE FLOW, HT MNT, SS SLEEVE	\$	9,693.82
F1A13D1T	aa SENSOR, PADDLE FLOW, HT MNT, BRASS SLEEVE	\$	5,068.54
F1A13D1T30N	SENSOR, PADDLE FLOW, HT MNT, BRASS SLEEVE, 30 FT	\$	6,167.04
F1C2A1002	ASSY. INSERTION TOOL	\$	3,989.30
FG2005602	ANET CABLE, 165 FEET	\$	2,331.91



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

FG7001019	INSTRUMENT SCREW DRIVER, PAT700	\$	98.09
FLB3205_00	nn 4G modem upgrade for Flexy eWON unit	\$	907.71
G53F1021-001	MODEL 53 MEMBRANE SWITCH	\$	292.93
GCRNA-SARPOSRGT	dd GeneCount COVID-19 Positive Control	\$	37.58
GCRNA-SARWW-96C	dd KIT.GeneCount SARS-CoV-2 Wastewater RT-qPCR KIT - 96 samples	\$	11,870.35
GCRNASARWW-96CD	dd DEMO UNIT KIT.GeneCount SARS-CoV-2 Wastewater RT-qPCR KIT - 96 samples	\$	2,963.07
GL1T	ASSY-CABLE, COAX TINNED	\$	398.93
HDF186	MAINTENANCE SCHEDULE, FILTRAX	\$	20.43
HDH969	LABEL, HACH BADGE	\$	18.50
HHA155	Adapter (Drm.20/Drm.16 - 50 lang)	\$	71.88
H-LINK	H-LINK INTERFACE MODULE	\$	4,567.46
HLS339	SAMPLE DELIVERY HOSE	\$	154.75
HMP4	MP-4 PORTABLE METER	\$	3,931.49
HMP4CASE	CASE ASSY, TOP & BOTTOM, MP-4	\$	1,485.87
HMP4EP	KTO: MP-4 ESSENTIALS PACKAGE	\$	4,239.84
HMP4MO	RPL BOARD, MP-4	\$	3,334.06
HMP6	MP-6 PORTABLE METER	\$	5,319.07
HMP6CASE	CASE ASSY, TOP & BOTTOM, MP-6	\$	1,549.47
HMP6EP	KTO: MP-6 ESSENTIALS PACKAGE	\$	5,627.42
HMP6MO	RPL BOARD, MP-6	\$	3,546.05
HMP6P	MP-6p PORTABLE METER	\$	5,319.07
HMP6PCASE	CASE ASSY, TOP & BOTTOM, MP-6p	\$	1,549.47
HMP6PEP	KTO: MP-6p ESSENTIALS PACKAGE	\$	5,627.42
HMPCASE	CASE, HARDSIDED MP SERIES	\$	192.72
HMPSENS	REPLACEMENT PH/ORP SENSOR	\$	950.11
HMPSENS-LM	REPLACEMENT PH/ORP SENSOR	\$	950.11
HPD1P1	Gateway, Analog pH Sensor, & Cable	\$	3,758.04
HPL061	PROBE PIPE STAND, OPTIQUANT	\$	528.05
HPT210	aa FREE CHLORINE LIQUID DPD, 100 TESTS	\$	96.84
HPT251	IRON ULR, 0.01-1.0 MG/L, PK/20	\$	166.22
HPT291	Manganese LR PAN (50)	\$	392.06
HPT310	aa FREE & TOTAL CHLORINE LIQUID DPD	\$	189.71
HPT430	RGT SET,SULFITE CLM	\$	178.68
HPT487	PHOSHORUS REACTIVE ULR (0.01-0.50 MG/L PO4-P), PK/20	\$	166.22
HQ411D	HQ411d BENCHTOP METER, pH/mV	\$	3,758.04
HQ411DNB	HQ411d BENCHTOP METER, pH/mV - NO BATTERY	\$	2,967.89
HQ430D	HQ430d, BENCHTOP METER, FLEXI	\$	5,357.62
HQ430DNB	HQ430d, BENCHTOP METER, FLEXI - NO BATTERY	\$	4,355.47
HQ440D	HQ440d BENCHTOP METER, MULTI	\$	6,725.93
HQ440DNB	HQ440d BENCHTOP METER, MULTI - NO BATTERY	\$	5,858.69
HRO304	CONTROLLER INSTALLATION PIPE, OPTIQ	\$	705.36
HRPHC6	PH CAP W/FOAM MP-6/MP-6p	\$	116.79
HYH020	DUST COVER, DR3800	\$	182.12
HYT059	DEMO CASE FOR DR 6000 PHOTOMETER	\$	2,736.62
HYT061	Demo Case TU5200	\$	1,046.47

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

ISECL18101	aa CHLORIDE PROBE, STD, w/1m CABLE	\$	3,603.86
ISECL18103	aa CHLORIDE PROBE, STD, w/3m CABLE	\$	3,738.77
ISECL181AP	KTO: INTELICAL CHLORIDE ANALYSIS PKG	\$	3,680.95
ISECL181AP-LM	KTO: INTELICAL CHLORIDE ANALYSIS PKG	\$	3,680.95
ISEF12101	aa FLUORIDE PROBE, STD, w/1m CABLE	\$	3,603.86
ISEF12103	aa FLUORIDE PROBE, STD, w/3m CABLE	\$	3,738.77
ISEF121AP	KTO: INTELICAL FLUORIDE ANALYSIS PKG	\$	3,777.31
ISEF121AP-LM	KTO: INTELICAL FLUORIDE ANALYSIS PKG	\$	3,777.31
ISENA38101	aa SODIUM REFILLABLE PROBE, STD, w/1m CABLE	\$	2,543.90
ISENA38103	aa SODIUM REFILLABLE PROBE, STD, w/3m CABLE	\$	2,659.54
ISENA381AP	KTO: INTELICAL SODIUM ANALYSIS PKG	\$	2,755.90
ISENA381AP-LM	KTO: INTELICAL SODIUM ANALYSIS PKG	\$	2,755.90
ISENH318101	aa AMMONIA PROBE, STD, w/1m CABLE	\$	3,218.42
ISENH318103	aa AMMONIA PROBE, STD, w/3m CABLE	\$	3,372.60
ISENH3181AP	KTO: INTELICAL AMMONIA ANALYSIS PKG	\$	3,796.58
ISENH3181AP-LM	KTO: INTELICAL AMMONIA ANALYSIS PKG	\$	3,796.58
ISENH418101	aa AMMONIUM PROBE, STD, w/1m CABLE	\$	3,526.78
ISENH418103	aa AMMONIUM PROBE, STD, w/3m CABLE	\$	3,584.59
ISENH4181AP	KTO: INTELICAL AMMONIUM ANALYSIS PKG	\$	3,584.59
ISENH4181AP-LM	KTO: INTELICAL AMMONIUM ANALYSIS PKG	\$	3,584.59
ISENO318101	aa NITRATE PROBE, STD, w/1m CABLE	\$	3,603.86
ISENO318103	aa NITRATE PROBE, STD, w/3m CABLE	\$	3,738.77
ISENO3181AP	KTO: INTELICAL NITRATE ANALYSIS PKG	\$	3,777.31
JOBALPLUS-SU	KTO: JOBPLUS SINGLE USER LICENSE, DVD	\$	12,835.15
JOBPLUS-ADL-USR	JOBAL PLUS Addtl Concurrent User	\$	3,044.98
JOB-SCADA-INT	JOBAL STD SCADA IF-list on web	\$	12,218.45
K46500	PLUG, DROPPER 15MM FLAT TIP	\$	1.93
KF1121.97	TitraLab KF1000,1 Syringe, 2 Pumps	\$	15,089.98
L2137	UV PORTABLE FIELD PACK POWER SUPPLY	\$	1,135.12
L2389CA	KTO: DPD TOTAL PP 5ML SAMPLE PK/100	\$	79.50
L2416	POWER SUPPLY 100-240V, 9VDC, 18W, GE WATER	\$	124.50
L2636	THERMOMETER, POCKET -10 TO 110C	\$	33.15
L2668	GOGGLES, SPECTACLES UV CONTROL	\$	51.65
LAB-3YR	Connected lab device, 3 year agreement	\$	414.35
LAB-YRLY	Connected lab device, annual subscription	\$	153.02
LBCW-80ML	dd Lysis Buffer Concentrate W, 80mL	\$	222.23
LBOD10101	LBOD PROBE W/1m CABLE	\$	3,623.14
LC-250ML	LUMICLEAN 250ML	\$	713.06
LCA503	Ammonium Standard LCA503	\$	85.81
LCA505	Ammonium Standard LCA505	\$	72.17
LCA720	rr Mixed-Parameter Standard, NIST Wastewater High Range	\$	136.95
LCA721	rr Mixed-Parameter Std, NIST Wastewater Low Range	\$	136.95
LCE009	Xpress SARS-CoV-2 Plus (RUO), 10 Cartridges	\$	1,266.53
LCK307	aa db LCK BORON RGTS, PK/24	\$	408.33
LCW802	REAGENT SET, SULFURIC ACID, AMTAX INTER 2	\$	870.85
LCW803	REAGENT SET, 5MG/L NH4-N, AMTAX INTER 2	\$	381.22
LCW804	ZERO SOLUTION FOR AMTAX INTER/2 5.2L	\$	225.84



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LCW819	CLEANING SOLUTION (2*250ml)	\$	254.75
LCW820	VANADATE/MOLYBDATE RGT FOR PHOSPHAX INTER/2, 10L	\$	2,204.24
LCW821	CLEANING SOLUTION FOR PHOSPHAX INTER/2, 10L	\$	1,250.27
LCW823	REAGENT SET, PHOSPHAX SIGMA	\$	2,493.32
LCW824	CALIBRATION STD 2 MG/L PHOSPHAX SIGMA, 500ML	\$	225.84
LCW825	Nitrate Std 50,0 mg/L NO ₃ , 11,3 mg/L NO ₃ -N	\$	153.21
LCW826	Nitrate Std 100 mg/L NO ₃ , 22,6 mg/L NO ₃ -N	\$	202.36
LCW827	Nitrate Std 200 mg/L NO ₃ , 45,2 mg/L NO ₃ -N	\$	202.36
LCW828	Nitrate Std 25,0 mg/L NO ₃ , 5,65 mg/L NO ₃ -N	\$	202.36
LCW862	LCW 862 AMTAX INTER II, STANDARD SOLUTION	\$	205.97
LCW863	Nitrate Std 400 mg/L NO ₃ , 90,4 mg/L NO ₃ -N	\$	202.36
LCW868	SET OF ELECTROLYTE, AMTAX	\$	334.25
LCW883	TEST SOLUTION #1, NH ₄ Dsc PROBE	\$	377.61
LCW884	TEST SOLUTION #2, NH ₄ Dsc PROBE	\$	377.61
LCW892	NH ₄ -N COMBINATION STD	\$	229.46
LCW895	TEST SOLUTION 1 FOR AN-ISE SC / A-ISE SC / N-ISE S	\$	374.00
LCW896	TEST SOLUTION 2 FOR AN-ISE SC / A-ISE SC / N-ISE S	\$	374.00
LCW942	Nitrate Std 15,0 mg/L NO ₃ , 3,39 mg/L NO ₃ -N	\$	161.16
LCW943	Nitrate Std 40,0 mg/L NO ₃ , 9,04 mg/L NO ₃ -N	\$	161.16
LCW944	Nitrate Std 75,0 mg/L NO ₃ , 16,9 mg/L NO ₃ -N	\$	161.16
LCW945	Nitrate Std 150 mg/L NO ₃ , 33,9 mg/L NO ₃ -N	\$	161.16
LCW946	Nitrate Std 300 mg/L NO ₃ , 67,8 mg/L NO ₃ -N	\$	161.16
LCW955	KTO: PHOSPHAX sc LR Reagent Set	\$	701.02
LCW956	PHOSPHAX sc LR Reagent A	\$	364.96
LCW957	PHOSPHAX sc LR Reagent B	\$	364.96
LCW958	PHOSPHAX sc LR Standard Solution	\$	256.56
LCW959	PHOSPHAX sc LR Cleaning Solutionon	\$	303.53
LCW960	Nitrite Std 5,00 mg/L NO ₂ , 1,52 mg/L NO ₂ -N	\$	161.16
LCW961	Nitrite Std 15,0 mg/L NO ₂ , 4,57 mg/L NO ₂ -N	\$	161.16
LCW962	Nitrite Std 25,0 mg/L NO ₂ , 7,61 mg/L NO ₂ -N	\$	161.16
LCW963	Nitrite Std 40,0 mg/L NO ₂ , 12,2 mg/L NO ₂ -N	\$	161.16
LCW964	Nitrite Std 75,0 mg/L NO ₂ , 22,8 mg/L NO ₂ -N	\$	161.16
LCW965	Nitrite Std 100 mg/L NO ₂ , 30,4 mg/L NO ₂ -N	\$	161.16
LCW966	Nitrite Std 150 mg/L NO ₂ , 45,7 mg/L NO ₂ -N	\$	161.16
LCW967	Nitrite Std 200 mg/L NO ₂ , 60,9 mg/L NO ₂ -N	\$	161.16
LCW968	Nitrite Std 300 mg/L NO ₂ , 91,4 mg/L NO ₂ -N	\$	161.16
LDO10101	LDO PROBE, STD, w/1m CABLE	\$	3,122.06
LDO10101DK	LDO PROBE, STD, w/1m CABLE - DK	\$	3,122.06
LDO10103	LDO PROBE, STD, w/3m CABLE	\$	3,314.78
LDO10103DK	LDO PROBE, STD, w/3m CABLE	\$	3,237.70
LDO10105	LDO PROBE, RUGGED, w/5m CABLE	\$	3,931.49
LDO10110	LDO PROBE, RUGGED, w/10m CABLE	\$	4,355.47
LDO10115	LDO PROBE, RUGGED, w/15m CABLE	\$	4,471.10
LDO10130	LDO PROBE, RUGGED, w/30m CABLE	\$	4,818.00
LEV015.53.11101	KTO:HQ1110/PHC20101 FIELD	\$	3,546.05
LEV015.53.11102	KTO: HQ1110/PHC10101 FIELD	\$	3,623.14
LEV015.53.11104	KTO: HQ1110/MTC10105 FIELD	\$	4,952.90

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LEV015.53.1110A	ee HQ1110 PORTABLE pH/ORP METER	\$	1,882.87
LEV015.53.11301	KTO: HQ1130/LDO10101 FIELD	\$	5,319.07
LEV015.53.11302	KTO:HQ1130/LDO10105 FIELD	\$	5,781.60
LEV015.53.1130A	ee HQ1130 Portable DO Meter	\$	2,023.56
LEV015.53.11401	KTO:HQ1140/CDC40101 FIELD	\$	4,297.66
LEV015.53.11402	KTO: HQ1140/CDC40105 FIELD	\$	5,511.79
LEV015.53.1140A	ee HQ1140 Portable EC/TDS Meter	\$	2,023.56
LEV015.53.21001	KTO:HQ2100/PHC20101 FIELD	\$	4,355.47
LEV015.53.21002	KTO:HQ2100/PHC10101 FIELD	\$	4,432.56
LEV015.53.21003	KTO: HQ2100/PHC10105 FIELD	\$	5,762.33
LEV015.53.21004	KTO: HQ2100/LDO10101 FIELD	\$	5,993.59
LEV015.53.21005	KTO: HQ2100/LDO10105 FIELD	\$	6,436.85
LEV015.53.21006	KTO:HQ2100/CDC40101 FIELD	\$	4,933.63
LEV015.53.21007	KTO:HQ2100/CDC40105 FIELD	\$	6,147.77
LEV015.53.2100A	ee HQ2100 Portable pH/EC/TDS/DO Meter	\$	2,582.45
LEV015.53.22001	KTO: HQ2200/PHC10101/CDC40101 FIELD	\$	6,841.56
LEV015.53.22002	KTO: HQ2200/PHC10105/CDC40105 FIELD	\$	7,920.79
LEV015.53.22003	KTO:HQ2200/PHC20101/CDC40101 FIELD	\$	6,687.38
LEV015.53.22004	KTO:HQ2200/PHC20101/LDO10101 FIELD	\$	7,361.90
LEV015.53.22005	KTO:HQ2200/PHC10101/LDO10101 FIELD	\$	7,535.35
LEV015.53.22006	KTO:HQ2200/PHC10105/LDO10105 FIELD	\$	8,865.12
LEV015.53.22008	KTO:HQ2200/LDO10105 FIELD	\$	7,053.55
LEV015.53.22009	KTO:HQ2200/PHC10105 FIELD	\$	6,494.66
LEV015.53.2200A	ee HQ2200 Portable pH/EC/TDS/DO Meter	\$	3,584.59
LEV015.53.41001	KTO:HQ4100/PHC20101 FIELD	\$	5,415.43
LEV015.53.41002	KTO:HQ4100/PHC10101 FIELD	\$	5,550.34
LEV015.53.41003	KTO:HQ4100/PHC10105 FIELD	\$	6,841.56
LEV015.53.41004	KTO:HQ4100/CDC40101 FIELD	\$	5,762.33
LEV015.53.41005	KTO:HQ4100/CDC40105 FIELD	\$	6,899.38
LEV015.53.4100A	ee HQ4100 Portable Multi Meter	\$	3,777.31
LEV015.53.42001	KTO:HQ4200/PHC20101/CDC40101 FIELD	\$	7,496.81
LEV015.53.42002	KTO:HQ4200/PHC10101/CDC40101 FIELD	\$	7,631.71
LEV015.53.42003	KTO:HQ4200/PHC10105/CDC40105 FIELD	\$	9,886.54
LEV015.53.42005	KTO:HQ4200/PHC10101/LDO10101 FIELD	\$	8,672.40
LEV015.53.42006	KTO:HQ4200/PHC10105/LDO10105 FIELD	\$	9,944.35
LEV015.53.4200A	ee HQ4200 Portable Multi Meter	\$	4,702.37
LEV015.53.43001	KTO:HQ4300/PHC201/CDC/LDO/01 FIELD	\$	10,541.78
LEV015.53.43002	KTO:HQ4300/PHC101/CDC/LDO/05 FIELD	\$	13,432.58
LEV015.53.4300A	ee HQ4300 Portable Multi Meter	\$	5,319.07
LEV015.80.11101	KTO:HQ1110/PHC20101 FIELD	\$	3,546.05
LEV015.80.11102	KTO: HQ1110/PHC10101 FIELD	\$	3,623.14
LEV015.80.11103	KTO: HQ1110/PHC10105 FIELD	\$	4,952.90
LEV015.80.11104	KTO: HQ1110/MTC10105 FIELD	\$	4,952.90
LEV015.80.1110A	ee HQ1110 PORTABLE pH/ORP METER	\$	1,882.87
LEV015.80.11301	KTO: HQ1130/LDO10101 FIELD	\$	5,319.07
LEV015.80.11302	KTO:HQ1130/LDO10105 FIELD	\$	5,781.60
LEV015.80.1130A	ee HQ1130 Portable DO Meter	\$	2,023.56



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LEV015.80.11401	KTO:HQ1140/CDC40101 FIELD	\$	4,297.66
LEV015.80.11402	KTO: HQ1140/CDC40105 FIELD	\$	5,511.79
LEV015.80.1140A	ee HQ1140 Portable EC/TDS Meter	\$	2,023.56
LEV015.80.21001	KTO:HQ2100/PHC20101 FIELD	\$	4,355.47
LEV015.80.21002	KTO:HQ2100/PHC10101 FIELD	\$	4,432.56
LEV015.80.21003	KTO: HQ2100/PHC10105 FIELD	\$	5,762.33
LEV015.80.21004	KTO: HQ2100/LDO10101 FIELD	\$	5,993.59
LEV015.80.21005	KTO: HQ2100/LDO10105 FIELD	\$	6,436.85
LEV015.80.21007	KTO:HQ2100/CDC40105 FIELD	\$	6,147.77
LEV015.80.2100A	ee HQ2100 Portable pH/EC/TDS/DO Meter	\$	2,698.08
LEV015.80.22001	KTO: HQ2200/PHC10101/CDC40101 FIELD	\$	6,841.56
LEV015.80.22002	KTO: HQ2200/PHC10105/CDC40105 FIELD	\$	7,920.79
LEV015.80.22003	KTO:HQ2200/PHC20101/CDC40101 FIELD	\$	6,687.38
LEV015.80.22004	KTO:HQ2200/PHC20101/LDO10101 FIELD	\$	7,361.90
LEV015.80.22005	KTO:HQ2200/PHC10101/LDO10101 FIELD	\$	7,535.35
LEV015.80.22006	KTO:HQ2200/PHC10105/LDO10105 FIELD	\$	8,865.12
LEV015.80.22008	KTO:HQ2200/LDO10105 FIELD	\$	7,053.55
LEV015.80.22009	KTO:HQ2200/PHC10105 FIELD	\$	6,494.66
LEV015.80.2200A	ee HQ2200 Portable pH/EC/TDS/DO Meter	\$	3,584.59
LEV015.80.41001	KTO:HQ4100/PHC20101 FIELD	\$	5,415.43
LEV015.80.41002	KTO:HQ4100/PHC10101 FIELD	\$	5,550.34
LEV015.80.41003	KTO:HQ4100/PHC10105 FIELD	\$	6,841.56
LEV015.80.41004	KTO:HQ4100/CDC40101 FIELD	\$	5,762.33
LEV015.80.41005	KTO:HQ4100/CDC40105 FIELD	\$	6,899.38
LEV015.80.4100A	ee HQ4100 Portable Multi Meter	\$	3,777.31
LEV015.80.42001	KTO:HQ4200/PHC20101/CDC40101 FIELD	\$	7,496.81
LEV015.80.42002	KTO:HQ4200/PHC10101/CDC40101 FIELD	\$	7,631.71
LEV015.80.42003	KTO:HQ4200/PHC10105/CDC40105 FIELD	\$	9,886.54
LEV015.80.42005	KTO:HQ4200/PHC10101/LDO10101 FIELD	\$	8,672.40
LEV015.80.42006	KTO:HQ4200/PHC10105/LDO10105 FIELD	\$	9,944.35
LEV015.80.4200A	ee HQ4200 Portable Multi Meter	\$	4,702.37
LEV015.80.43001	KTO:HQ4300/PHC201/CDC/LDO/01 FIELD	\$	10,541.78
LEV015.80.43002	KTO:HQ4300/PHC101/CDC/LDO/05 FIELD	\$	13,432.58
LEV015.80.4300A	ee HQ4300 Portable Multi Meter	\$	5,319.07
LEV015.97.11101	KTO:HQ1110/PHC20101 FIELD	\$	3,546.05
LEV015.97.11102	KTO: HQ1110/PHC10101 FIELD	\$	3,623.14
LEV015.97.11103	KTO: HQ1110/PHC10105 FIELD	\$	4,952.90
LEV015.97.11104	KTO: HQ1110/MTC10105 FIELD	\$	4,952.90
LEV015.97.1110A	ee HQ1110 PORTABLE pH/ORP METER	\$	1,882.87
LEV015.97.11301	KTO: HQ1130/LDO10101 FIELD	\$	5,319.07
LEV015.97.11302	KTO:HQ1130/LDO10105 FIELD	\$	5,781.60
LEV015.97.1130A	ee HQ1130 Portable DO Meter	\$	2,023.56
LEV015.97.11401	KTO:HQ1140/CDC40101 FIELD	\$	4,297.66
LEV015.97.11402	KTO: HQ1140/CDC40105 FIELD	\$	5,511.79
LEV015.97.1140A	ee HQ1140 Portable EC/TDS Meter	\$	2,023.56
LEV015.97.21001	KTO:HQ2100/PHC20101 FIELD	\$	4,355.47
LEV015.97.21002	KTO:HQ2100/PHC10101 FIELD	\$	4,432.56



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LEV015.97.21003	KTO: HQ2100/PHC10105 FIELD	\$	5,762.33
LEV015.97.21004	KTO: HQ2100/LDO10101 FIELD	\$	5,993.59
LEV015.97.21005	KTO: HQ2100/LDO10105 FIELD	\$	6,436.85
LEV015.97.21006	KTO:HQ2100/CDC40101 FIELD	\$	4,933.63
LEV015.97.21007	KTO:HQ2100/CDC40105 FIELD	\$	6,147.77
LEV015.97.2100A	ee HQ2100 Portable pH/EC/TDS/DO Meter	\$	2,698.08
LEV015.97.22001	KTO: HQ2200/PHC10101/CDC40101 FIELD	\$	6,841.56
LEV015.97.22002	KTO: HQ2200/PHC10105/CDC40105 FIELD	\$	7,920.79
LEV015.97.22003	KTO:HQ2200/PHC20101/CDC40101 FIELD	\$	6,687.38
LEV015.97.22004	KTO:HQ2200/PHC20101/LDO10101 FIELD	\$	7,361.90
LEV015.97.22005	KTO:HQ2200/PHC10101/LDO10101 FIELD	\$	7,535.35
LEV015.97.22006	KTO:HQ2200/PHC10105/LDO10105 FIELD	\$	8,865.12
LEV015.97.22008	KTO:HQ2200/LDO10105 FIELD	\$	7,053.55
LEV015.97.22009	KTO:HQ2200/PHC10105 FIELD	\$	6,494.66
LEV015.97.2200A	ee HQ2200 Portable pH/EC/TDS/DO Meter	\$	3,584.59
LEV015.97.41001	KTO:HQ4100/PHC20101 FIELD	\$	5,415.43
LEV015.97.41002	KTO:HQ4100/PHC10101 FIELD	\$	5,550.34
LEV015.97.41003	KTO:HQ4100/PHC10105 FIELD	\$	6,841.56
LEV015.97.41004	KTO:HQ4100/CDC40101 FIELD	\$	5,762.33
LEV015.97.41005	KTO:HQ4100/CDC40105 FIELD	\$	6,899.38
LEV015.97.4100A	ee HQ4100 Portable Multi Meter	\$	3,777.31
LEV015.97.42001	KTO:HQ4200/PHC20101/CDC40101 FIELD	\$	7,496.81
LEV015.97.42002	KTO:HQ4200/PHC10101/CDC40101 FIELD	\$	7,631.71
LEV015.97.42003	KTO:HQ4200/PHC10105/CDC40105 FIELD	\$	9,886.54
LEV015.97.42005	KTO:HQ4200/PHC10101/LDO10101 FIELD	\$	8,672.40
LEV015.97.42006	KTO:HQ4200/PHC10105/LDO10105 FIELD	\$	9,944.35
LEV015.97.4200A	ee HQ4200 Portable Multi Meter	\$	4,702.37
LEV015.97.43001	KTO:HQ4300/PHC201/CDC/LDO/01 FIELD	\$	10,541.78
LEV015.97.43002	KTO:HQ4300/PHC101/CDC/LDO/05 FIELD	\$	13,432.58
LEV015.97.4300A	ee HQ4300 Portable Multi Meter	\$	5,319.07
LEV015.98.1110A	ee HQ1110 PORTABLE pH/ORP METER	\$	1,882.87
LEV015.98.1130A	ee HQ1130 Portable DO Meter	\$	2,023.56
LEV015.98.1140A	ee HQ1140 Portable EC/TDS Meter	\$	2,023.56
LEV015.98.2100A	ee HQ2100 Portable pH/EC/TDS/DO Meter	\$	2,698.08
LEV015.98.2200A	ee HQ2200 Portable pH/EC/TDS/DO Meter	\$	3,584.59
LEV015.98.4100A	ee HQ4100 Portable Multi Meter	\$	3,777.31
LEV015.98.4200A	ee HQ4200 Portable Multi Meter	\$	4,702.37
LEV015.98.4300A	ee HQ4300 Portable Multi Meter	\$	5,319.07
LEV015.99.1110A	ee HQ1110 PORTABLE pH/ORP METER	\$	1,882.87
LEV015.99.11302	KTO:HQ1130/LDO10105 FIELD	\$	5,781.60
LEV015.99.1130A	ee HQ1130 Portable DO Meter	\$	2,023.56
LEV015.99.1140A	ee HQ1140 Portable EC/TDS Meter	\$	2,023.56
LEV015.99.2100A	ee HQ2100 Portable pH/EC/TDS/DO Meter	\$	2,698.08
LEV015.99.22002	KTO: HQ2200/PHC10105/CDC40105 FIELD	\$	7,920.79
LEV015.99.2200A	ee HQ2200 Portable pH/EC/TDS/DO Meter	\$	3,584.59
LEV015.99.4100A	ee HQ4100 Portable Multi Meter	\$	3,777.31
LEV015.99.42002	KTO:HQ4200/PHC10101/CDC40101 FIELD	\$	7,631.71



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LEV015.99.4200A	ee HQ4200 Portable Multi Meter	\$	4,702.37
LEV015.99.43001	KTO:HQ4300/PHC201/CDC/LDO/01 FIELD	\$	10,541.78
LEV015.99.4300A	ee HQ4300 Portable Multi Meter	\$	5,319.07
LEV016.97.00001	GeneXpert IV, 2 Module Configuration w/Laptop	\$	122,878.27
LEV016.97.10001	GeneXpert IV, 2 Module Configuration w/Laptop	\$	104,473.51
LEZ015.53.A001A	ACCESSORY BOX, NORTH AMERICA, HQ SERIES	\$	335.33
LEZ015.80.A001A	ACCESSORY BOX, CHINA, HQ SERIES	\$	335.33
LEZ015.97.A001A	ACCESSORY BOX, ROW, HQ SERIES	\$	346.90
LEZ015.98.A001A	ACCESSORY BOX, EU AND UK, HQ SERIES	\$	352.68
LEZ015.99.00001	BATTERY, LITHIUM ION, 18650, 3200mAh, ENGLISH, PKGD	\$	179.42
LEZ015.99.00002	CABLE, USB 2.0 A TO MICRO B, HQ SERIES	\$	83.83
LEZ015.99.00005	ee POWER SUPPLY PLUS CHINA CLIP	\$	91.54
LEZ015.99.00006	ee POWER SUPPLY PLUS US CLIP	\$	92.51
LEZ015.99.00007	ee POWER SUPPLY PLUS ROW CLIPS	\$	130.86
LEZ015.99.A001A	HQ STANDARD PROBE FIELD CASE	\$	514.56
LEZ015.99.A002A	HQ RUGGED PROBE FIELD CASE	\$	618.63
LEZ015.99.A003A	HQ KICKSTAND AND HANDSTRAP	\$	99.25
LEZ015.99.A004A	HQ PROTECTIVE GLOVE	\$	217.77
LEZ015.99.A005A	WRIST STRAP, DUST PLUGS, CAPS, HQ SERIES	\$	44.71
LEZ015.99.A006A	ACCESSORY BOX, NO BATT, NO PWR SUPPLY, HQ SERIES	\$	136.45
LMK-5	LUMINOMETER MAINTENANCE KIT	\$	1,381.80
LPG439.01.00002	2100Q, INSTRUMENT ONLY	\$	3,738.77
LPG439.01.00012	2100Q IS, INSTRUMENT ONLY	\$	2,929.34
LPV2000.97.0002	KTO: sensION+ pH3 BENCHTOP	\$	2,478.86
LPV2000.97.8882	SENSION+ PH3 BENCHTOP W/TEST CERT	\$	2,074.15
LPV2000KIT	SENSION PH3, PH31, MM340	\$	312.21
LPV2000KIT-LM	SENSION PH3, PH31, MM340	\$	312.21
LPV2010T.97.002	KTO: sensION+ pH3 BENCHTOP KIT, 5010T	\$	3,220.83
LPV2010T.LM.002	KTO: sensION+ pH3 BENCHTOP KIT, 5010T	\$	3,220.83
LPV2011T.97.002	KTO: sensION+ pH3 BENCHTOP KIT, 5011T	\$	3,693.00
LPV2014T.97.002	KTO: sensION+ pH3 BENCHTOP KIT, 5014T	\$	3,524.37
LPV2014T.LM.002	KTO: sensION+ pH3 BENCHTOP KIT, 5014T	\$	3,524.37
LPV2021T.97.002	KTO: sensION+ pH3 BENCHTOP KIT, 5021T	\$	3,912.22
LPV2100.97.0002	KTO: sensION+ pH31 BENCHTOP	\$	3,085.93
LPV2100.97.8882	SENSION+ PH31 BENCHTOP W/TEST CERT	\$	2,832.98
LPV2110T.97.002	KTO: sensION+ pH31 BENCHTOP KIT, 5010T	\$	3,878.49
LPV2111T.97.002	KTO: sensION+ pH31 BENCHTOP KIT, 5011T	\$	4,283.20
LPV2114T.97.002	KTO: sensION+ pH31 BENCHTOP KIT, 5014T	\$	4,114.57
LPV2121T.97.002	KTO: sensION+ pH31 BENCHTOP KIT, 5021T	\$	4,519.28
LPV2200.97.0002	KTO: sensION+ MM 340 BENCHTOP	\$	4,266.34
LPV2200.97.8882	SENSION+ MM340 BENCHTOP W/TEST CERT	\$	4,047.12
LPV2200.LM.0002	KTO: sensION+ MM 340 BENCHTOP	\$	3,355.74
LPV2210T.97.002	KTO: sensION+ MM 340 BENCHTOP KIT, 5010T	\$	4,974.59
LPV2210T.LM.002	KTO: sensION+ MM 340 BENCHTOP KIT, 5010T	\$	4,974.59
LPV2211T.97.002	KTO: sensION+ MM 340 BENCHTOP KIT, 5011T	\$	5,396.16
LPV2214T.97.002	KTO: sensION+ MM 340 BENCHTOP KIT, 5014T	\$	5,311.85
LPV2221T.97.002	KTO: sensION+ MM 340 BENCHTOP KIT, 5021T	\$	5,665.97

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LPV3000KIT	SENSION EC7, EC71	\$	236.08
LPV3070.97.0002	KTO: sensION+ EC7 BENCHTOP KIT, 5070	\$	3,929.08
LPV4000KIT	SENSION MM374	\$	136.08
LPV4110.97.0002	KTO: sensION+ MM 374 BENCHTOP	\$	5,413.02
LPV4140.97.0002	KTO: sensION+ MM 374 BENCHTOP KIT 5014, 5070	\$	7,352.27
LPV415.99.12002	ULTRATURB SEAWATER sc W/ 0.35m CABLE	\$	12,334.08
LPV415.99.22002	db ULTRATURB SEAWATER sc w/5M CABLE	\$	11,293.39
LPV415.99.32002	ULTRATURB SEAWATER sc W/10m CABLE	\$	11,409.02
LPV415.99.82002	ULTRATURB SEAWATER sc W/1m CABLE	\$	11,062.13
LPV431.32.00002	SS7 HACH SENSOR ASSY	\$	15,937.94
LPV431.35.00002	SS7 HACH SENSOR ASSY	\$	15,803.04
LPV431.36.00002	SS7 HACH SENSOR ASSY	\$	15,803.04
LPV431.90.00002	SS7 HACH SENSOR ASSY	\$	15,668.14
LPV431.93.00002	SS7 HACH SENSOR ASSY	\$	15,937.94
LPV431.99.00002	SS7 HACH SENSOR ASSY	\$	16,863.00
LPV431.LM.00002	SS7 HACH SENSOR ASSY	\$	16,863.00
LPV440.99.00002	aa DR3900 SPECTROPHOTOMETER W/O RFID	\$	18,346.94
LPV440.99.00012	db aa DR3900 SPECTROPHOTOMETER WITH RFID	\$	18,346.94
LPV441.99.00002	DR 6000 UV VIS SPECTROPHOTOMETER W/O RFID TECH	\$	35,845.92
LPV441.99.00012	db DR 6000 UV VIS SPECTROPHOTOMETER W RFID TECHNOLOGY	\$	35,845.92
LPV441.99.10002	DR 6000EDU UV VIS SPECTROPHOTOMETER	\$	21,372.65
LPV442.32.01022	KTO: TU5200 Lab Turb, ISO	\$	12,334.08
LPV442.35.01012	KTO: TU5200 Lab Turb, EPA	\$	12,334.08
LPV442.35.01022	KTO: TU5200 Lab Turb, ISO	\$	12,334.08
LPV442.99.01012	KTO: TU5200 Lab Turb, EPA	\$	13,085.69
LPV442.99.01014	TU5200 (w/o StablCal Set) EPA w/o RFID	\$	11,659.56
LPV442.99.01022	KTO: TU5200 Lab Turb, ISO	\$	13,085.69
LPV442.99.03012	KTO: TU5200, Lab Turb with RFID, EPA	\$	13,856.57
LPV442.99.03022	KTO: TU5200, Lab Turb with RFID, ISO	\$	14,665.99
LPV442.LM.01012	KTO: TU5200 Lab Turb, EPA	\$	13,085.69
LPV442.LM.01022	KTO: TU5200 Lab Turb, ISO	\$	13,085.69
LPV442.LM.03012	KTO: TU5200, Lab Turb with RFID, EPA	\$	13,856.57
LPV442.LM.03022	KTO: TU5200, Lab Turb with RFID, ISO	\$	14,665.99
LPV4425301012	nn KIT, TU5200 Lab Turb, EPA á	\$	11,813.74
LPV4425301022	nn KIT, TU5200 LAB TURB, ISO	\$	12,199.18
LPV4425303012	nn KIT, TU5200, Lab Turb with RFID, EPA á	\$	12,064.27
LPV4425303022	nn KIT, TU5200, Lab Turb with RFID, ISO á	\$	12,912.24
LPV4429301012	nn KIT, TU5200 Lab Turb, EPA á	\$	12,334.08
LPV442LM03012	nn KIT, TU5200, Lab Turb with RFID, EPA á	\$	12,064.27
LPV443.99.10002	SL1000 - Portable Parallel Analyzer Packed in carry case with power cable and manual	\$	11,081.40
LPV443.99.20002	Full Bundle SL1000 - Portable Parallel Analyzer	\$	13,683.12
LPV444.32.00120	ee KTO: TL2310 LAB TURBIDIMETER, ISO, EU, 0-1000 NTU	\$	8,402.59
LPV444.32.00310	ee KTO: TL2350 LAB TURBIDIMETER, EPA, EU, 0-10000 NTU	\$	13,259.14
LPV444.32.00320	ee KTO: TL2360 LAB TURBIDIMETER, ISO, EU, 0-10000 NTU	\$	13,085.69
LPV444.80.00120	TL2310 LAB TURBIDIMETER, ISO, CN PKG	\$	10,618.87



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LPV444.80.00310	TL2350 LAB TURBIDIMETER, EPA, CN PKG	\$	19,214.18
LPV444.80.00320	TL2360 LAB TURBIDIMETER, ISO, CN PKG	\$	19,214.18
LPV444.93.00210	ee KTO: TL2300 LAB TURBIDIMETER, EPA, EU, 0-4000 NTU	\$	7,959.34
LPV444.99.00120	ee KTO: TL2310 LAB TURBIDIMETER, ISO, EU, 0-1000 NTU	\$	8,903.66
LPV444.99.00210	ee KTO: TL2300 LAB TURBIDIMETER, EPA, EU, 0-4000 NTU	\$	8,845.85
LPV444.99.00310	ee KTO: TL2350 LAB TURBIDIMETER, EPA, EU, 0-10000 NTU	\$	14,010.74
LPV444.99.00320	ee KTO: TL2360 LAB TURBIDIMETER, ISO, EU, 0-10000 NTU	\$	13,856.57
LPV444.LM.00120	ee KTO: TL2310 LAB TURBIDIMETER, ISO, EU, 0-1000 NTU	\$	8,903.66
LPV444.LM.00210	ee KTO: TL2300 LAB TURBIDIMETER, EPA, EU, 0-4000 NTU	\$	8,845.85
LPV444.LM.00310	ee KTO: TL2350 LAB TURBIDIMETER, EPA, EU, 0-10000 NTU	\$	14,010.74
LPV4443500120	ee KIT, TL2310 LAB TURBIDIMETER, ISO, NA	\$	8,325.50
LPV4449900120	ee KIT, TL2310 LAB TURBIDIMETER, ISO, NA	\$	8,980.75
LPV4449900210	ee KIT, TL2300 LAB TURBIDIMETER, EPA, NA	\$	8,884.39
LPV4449900310	ee KIT, TL2350 LAB TURBIDIMETER, EPA, NA	\$	14,126.38
LPV4449900320	ee KIT, TL2360 LAB TURBIDIMETER, ISO, NA	\$	13,933.66
LPV444LM00120	ee KIT, TL2310 LAB TURBIDIMETER, ISO, NA	\$	8,980.75
LPV444LM00210	ee KIT, TL2300 LAB TURBIDIMETER, EPA, NA	\$	8,884.39
LPV444LM00310	ee KIT, TL2350 LAB TURBIDIMETER, EPA, NA	\$	14,126.38
LPV445.97.00000	ee DR300 - Chlorine, Free + Total, NB, Pocket Colorimeter	\$	1,485.87
LPV445.97.00010	ee DR300 - Chlorine, Free + Total, Pocket Colorimeter	\$	1,350.97
LPV445.97.00100	ee DR300 - Chlorine, Free + Total, w.Box,NB, Pocket Colorimeter	\$	1,485.87
LPV445.97.00110	ee DR300 - Chlorine, Free + Total, w.Box, Pocket Colorimeter	\$	1,568.74
LPV445.97.01100	ee DR300 - Bromine, w. Box, NB, Pocket Colorimeter	\$	1,401.07
LPV445.97.01110	ee DR300 - Bromine, w. Box, Pocket Colorimeter	\$	1,568.74
LPV445.97.02110	ee DR300 - Nitrate, w. Box, Pocket Colorimeter	\$	1,568.74
LPV445.97.03110	ee DR300 - Dissolved Oxygen, w. Box, Pocket Colorimeter	\$	1,568.74
LPV445.97.04000	ee DR300 - Ozone, NB, Pocket Colorimeter	\$	1,485.87
LPV445.97.04100	ee DR300 - Ozone, w. Box, NB, Pocket Colorimeter	\$	1,410.71
LPV445.97.04110	ee DR300 - Ozone, w. Box, Pocket Colorimeter	\$	1,568.74
LPV445.97.06100	ee DR300 - Phosphate, w. Box, NB, Pocket Colorimeter	\$	1,410.71
LPV445.97.06110	ee DR300 - Phosphate, w. Box, Pocket Colorimeter	\$	1,568.74
LPV445.97.09110	ee DR300 - Zinc, w. Box, Pocket Colorimeter	\$	1,568.74
LPV445.97.10000	ee DR300 - Molybdenum, LR/HR, NB, Pocket Colorimeter	\$	1,485.87
LPV445.97.10110	ee DR300 - Molybdenum, LR/HR, w.Box, Pocket Colorimeter	\$	1,568.74
LPV445.97.12000	ee DR300 - Chlorine & pH, NB, Pocket Colorimeter	\$	1,485.87
LPV445.97.12110	ee DR300 - Chlorine & pH, w. Box, Pocket Colorimeter	\$	1,568.74
LPV445.97.15110	ee DR300 - Manganese, w. Box, HR, Pocket Colorimeter	\$	1,568.74
LPV445.97.16110	ee DR300 - Iron, TPTZ, w. Box, Pocket Colorimeter	\$	1,568.74
LPV445.97.22110	ee DR300 - Iron, FerroVer, w. Box, Pocket Colorimeter	\$	1,568.74
LPV445.97.25110	ee DR300 - Aluminum, w. Box, Pocket Colorimeter	\$	1,568.74
LPV445.97.26110	ee DR300 - Monochlor/Free NH3, w. Box, Pocket Colorimeter	\$	1,568.74
LPV445.97.40100	ee DR300 - Ammonia Nitrogen, w.Box, NB, Pocket Colorimeter	\$	1,410.71
LPV445.97.40110	ee DR300 - Ammonia Nitrogen, w. Box, Pocket Colorimeter	\$	1,568.74
LPV445.97.50100	ee DR300 - 500 nm, w. Box, NB, Pocket Colorimeter	\$	1,410.71
LPV445.97.50110	ee DR300 - 500 nm, w. Box, Pocket Colorimeter	\$	1,568.74
LPV445.97.51100	ee DR300 - Chlorine Dioxide, w.Box, NB, Pocket Colorimeter	\$	1,453.11
LPV445.97.51110	ee DR300 - Chlorine Dioxide, w. Box, Pocket Colorimeter	\$	1,568.74

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LPV445.97.52110	ee DR300 - 528 nm, w. Box, Pocket Colorimeter	\$	1,568.74
LPV445.97.60100	ee DR300 - 600 nm, w. Box, NB, Pocket Colorimeter	\$	1,404.93
LPV445.97.60110	ee DR300 - 600 nm, w. Box, Pocket Colorimeter	\$	1,568.74
LPV445.97.62110	ee DR300 - Cl, Free+Total, MR, w.Box, Pocket Colorimeter	\$	1,568.74
LPV445.97.65110	ee DR300 - 655 nm, w. Box, Pocket Colorimeter	\$	1,568.74
LPV446.99.00012	ee HCD (Hach Communication Dongle)	\$	281.37
LPV449.97.01002	DR1300 FL Portable Fluorometer (Bluetooth)	\$	5,453.98
LPV449.98.01002	DR1300 FL Portable Fluorometer (Bluetooth) (EU)	\$	5,453.98
LPV449.99.00002	DR1300 FL Portable Fluorometer (No Bluetooth)	\$	5,453.98
LPZ440.99.00001	Truecal - Lot Specific Calibration Software for DR3900/DR6000	\$	713.06
	Housing top DR 3900		
LPZ440.99.00002	version 2 (2018)	\$	112.93
	Housing bottom DR 3900		
LPZ440.99.00003	version 2 (2018)	\$	112.93
LPZ440.99.00004	Front cover DR 3900 version 2 (2018)	\$	71.88
	Housing set DR 3900		
LPZ440.99.00005	version 2 (2018)	\$	258.24
LPZ440.99.00007	HALOGEN LAMP - DR 3900 and LICO 6xx	\$	312.21
LPZ444.99.00001	Cell adapter TL23xx 11mm	\$	846.04
LPZ445.99.00006	Cover DR300	\$	40.28
LPZ445.99.00045	LED, WHITE	\$	8.48
LPZ445.99.00048	Spring, cell compartment	\$	8.48
LPZ449.99.00001	Sample Holder Adapter & Cover, DR1300 FL	\$	61.48
LPZ449.99.00002	USB Bluetooth Converter, DR1300 FL	\$	61.48
LQG148.99.00002	TOC-X5 TOC SHAKER FOR PURGING METHOD	\$	2,929.34
LQV148.99.00002	TOC-X5 TOC SHAKER FOR PURGING METHOD	\$	2,929.34
LQV155.99.00002	POWER BOX FOR sc ANALYZER W/O POWER CORD	\$	1,821.20
LQV155.99.00012	POWER BOX FOR sc ANALYZER W/POWER CORD	\$	1,965.74
LQV156.99.10012	LOC 100, Locator (RFID-Set)	\$	2,197.01
LQV156.99.20012	LOC 100, LOCATOR ONLY	\$	1,686.30
LQV157.99.00002	SIP 10, SIPPER DEVICE ONLY	\$	1,618.85
LQV157.99.10002	SIP 10, SIPPER KIT W/DUAL PATH CELL	\$	2,254.82
LQV157.99.20002	SIP 10 SIPPER SET - DR 6000 - 1 INCH	\$	3,122.06
LQV157.99.30001	SIP 10 SIPPER SET - DR 6000/LICO - 1 CM QUARTZ	\$	6,070.68
LQV157.99.30002	SIP 10 SIPPER SET - DR 6000 - 1 CM QUARTZ	\$	6,436.85
LQV157.99.40002	SIP 10 - Sipper for TU5200	\$	1,865.53
	SIP 10 - Sipper, Probenpumpe		
LQV157.99.50002	Set	\$	2,062.10
LQV157.99.60002	SIP 10, SIPPER PUMP, TL23	\$	2,177.74
LQV159.97.00002	db ee Automatic Cleaning Module TU5 (ROW)	\$	2,582.45
LQV160.99.00002	db ee TU5xxx sc Flow Sensor, Hach	\$	555.03
LRZ450.53.00010	Photocolorimetric module with sensor and power supply (US)	\$	5,319.07
LS-10ML-50R	LUMISOLVE 10ML, STABILIZER TUBE, 50/RACK	\$	1,214.14
LS-9ML-50R	LUMISOLVE 9ML, STABILIZER TUBE, 50/PACK	\$	990.58
LTV082.52.30001	KTO:DRB200 REACTR-230V 9X16MM 2X20MM	\$	3,391.87
LTV082.52.40001	KTO:DRB200 REACTOR-230V 15X16MM	\$	3,391.87
LTV082.52.42001	KTO:DRB200 REACTR-230V 21X16MM4X20MM	\$	4,914.36



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LTV082.52.44001	KTO: DRB200 REACTOR 230V 30 X 16MM	\$	4,740.91
LTV082.53.30001	DRB200 REACTOR, 110V 9X16MM 2X20MM	\$	3,391.87
LTV082.53.40001	DRB200 REACTOR, 110V 15X16MM	\$	3,391.87
LTV082.53.42001	DRB200 REACTOR, 110V 21X16MM, 4X20MM	\$	4,914.36
LTV082.53.44001	DRB200 REACTOR, 110V 30 X 16MM	\$	4,740.91
LTV082.80.44001	Reactor DRB 200 (CN) Hach 2 Blocks (9x16/6x16 + 9x16/6x16)	\$	5,222.71
LU-3ML-FD	3ML LUMINASE, FREEZE DRIED PLUS REHY BUFF	\$	523.96
LUL-3ML-FD	3ML LUMINASE LITE, FREEZE DRIED PLUS REHY BUFF	\$	325.22
LUW-3ML-FD	3ML LUMINASE W, FREEZE DRIED PLUS REHY BUFF	\$	301.73
LUXL-3ML-FD	3ML LUMINASE XL, FREEZE DRIED PLUS REHY BUFF	\$	981.07
LXU421.99.13002	AMTAX sc AMMONIA ANLZR 115-230V,1 CH	\$	61,535.50
LXV294.54.00000	db FILTRAX CONTROL MODULE, 115VAC	\$	21,430.46
LXV294.64.00000	db FILTRAX CONTROL MODULE, 230VAC	\$	21,430.46
LXV320.99.00002	TSS PORTABLE HAND-HELD METER	\$	2,543.90
LXV321.99.00002	TSS PORTABLE PROBE, CABLE 10M, PLUG	\$	8,132.78
LXV322.99.00002	db TSS PORTABLE COMPLETE IN CASE	\$	11,293.39
LXV323.99.10002	TSS sc, INSITU PROBE, W/O WIPER	\$	23,184.22
LXV323.99.20002	TSS sc TRICLAMP INLINE PROBE W/O WIPER	\$	25,246.32
LXV323.99.30002	TSS SC INLINE SENSOR	\$	25,708.85
LXV324.99.10002	db TSS, w sc, INSITU PROBE, WITH WIPER	\$	23,974.37
LXV324.99.20002	db TSS, w sc TRICLAMP INLINE WITH WIPER	\$	25,708.85
LXV324.99.30002	TSS W SC INLINE SENSOR	\$	26,152.10
LXV325.99.10002	TSS HT sc, INSITU PROBE, W/O WIPER	\$	23,184.22
LXV325.99.20002	db TSS HT sc, TRICLAMP INLINE, W/O WIPER	\$	25,246.32
LXV325.99.30002	TSS HT SC INLINE SENSOR	\$	25,651.03
LXV326.99.10002	TSS VARI sc, INLINE PROBE, W/O WIPER	\$	23,184.22
LXV327.99.10002	TSS XL sc, INLINE PROBE, W/O WIPER	\$	23,184.22
LXV328.99.10002	TSS EX1 SC IMMERSION SENSOR, W/O WIPER, CI D2	\$	23,974.37
LXV328.99.20002	TSS EX1 SC TRICLAMP INLINE SENSOR, W/O WIPER, CI D2	\$	25,708.85
LXV328.99.30002	TSS EX1 SC INLINE SENSOR, W/O WIPER, CI D2	\$	26,152.10
LXV329.99.10002	TSS TITANIUM2 sc INSITU, W/O WIPER	\$	30,873.74
LXV329.99.20002	TSS TITANIUM2 sc TRICLAMP, INLINE	\$	32,878.03
LXV330.99.10002	TSS TITANIUM7 sc, INSITU W/O WIPER	\$	40,085.76
LXV330.99.20002	TSS TITANIUM7 sc TRICLAMP, INLINE	\$	42,128.59
LXV400.99.00032	ee SC1000 PM 6Sens	\$	3,025.70
LXV400.99.00121	SC1000 PM 4Sens 4mA OUT	\$	4,625.28
LXV400.99.00122	ee SC1000 PM 4Sens 4mA OUT	\$	3,989.30
LXV400.99.00332	ee SC1000 PM 6Sens 8mA OUT	\$	5,319.07
LXV400.99.0A022	db ee SC1000 PM 4 SENS 4-20mA IN MODBUS	\$	3,989.30
LXV400.99.0A032	ee SC1000 PM 6Sens	\$	4,374.74
LXV400.99.0R032	ee SC1000 PM 6SENS REL 110-230V HACH	\$	4,066.39
LXV400.99.0R322	ee SC1000 PM 4Sens 8mA OUT	\$	6,205.58
LXV400.99.0R332	ee SC1000 PM 6Sens 8mA OUT	\$	6,475.39
LXV400.99.10082	ee MODULE, SC1000 PROBE 6 SENS 110-230V	\$	3,796.58
LXV400.99.10092	ee MODULE, SC1000 PROBE, 8 SENSOR	\$	4,027.85
LXV400.99.10172	db ee MODULE, PROBE, (4 SNSR) SC1000 4mA	\$	4,952.90
LXV400.99.10372	ee SC1000 PM 4 Sens Cond 8mA OUT	\$	6,398.30

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LXV400.99.10382	ee MODULE, PROBE (6 SNSR) SC1000, 8mA	\$	6,610.30
LXV400.99.10572	ee MODULE, PROBE, (4 SNSR) SC1000 4mA	\$	6,918.65
LXV400.99.1A072	ee MODULE, PROBE, (4 SNSR), SC1000	\$	4,952.90
LXV400.99.1A082	ee MODULE, PROBE (6 SNSR), SC1000	\$	5,280.53
LXV400.99.1B082	ee MODULE, PROBE, 6 SNSR, REL, 485, sc1000	\$	7,130.64
LXV400.99.1B282	SC1000 PM 6Sens Cond 4mA/digit IN	\$	9,327.65
LXV400.99.1B482	ee SC1000 PRB MOD, REL, 485, 2I, 6SNS	\$	10,021.44
LXV400.99.1B572	db ee MODULE, PROBE (4-SNSR), SC1000, 4mA	\$	9,616.73
LXV400.99.1B582	ee SC1000 PRB MOD, REL, 485, IO, 6	\$	10,734.50
LXV400.99.1E072	ee MODULE, PROBE, (4 SNSR), SC1000	\$	5,299.80
LXV400.99.1E082	MODULE, PROBE (6-SNSR), SC1000	\$	5,550.34
LXV400.99.1F572	db ee MODULE, PROBE (4 SENSOR), SC1000	\$	10,021.44
LXV400.99.1G082	ee SC1000 PM 6Sens Cond	\$	8,267.69
LXV400.99.1G092	ee SC1000 PM 8Sens Cond	\$	8,537.50
LXV400.99.1G172	ee SC1000 PM 4 SENS COND 4MA OUT PROGNOYS HACH	\$	9,616.73
LXV400.99.1G282	ee SC1000 PM 6 SENS COND PROGNOYS	\$	10,426.15
LXV400.99.1G382	ee SC1000 PM 6 SENS COND 2X 4MA OUT PROGNOYS HACH	\$	11,351.21
LXV400.99.1H072	ee SC1000 PM 4Sens Cond	\$	9,616.73
LXV400.99.1H082	ee SC1000 PM 6Sens Cond	\$	9,828.72
LXV400.99.1N182	ee SC1000 PM 6 SENS 4-20mA OUT RTC PROGNOYS REL	\$	11,351.21
LXV400.99.1N582	ee SC1000 PM 6 SENS 4-20mA IN/OUT RTC PROGNOYS REL	\$	14,704.54
LXV400.99.1R082	ee MODULE, PROBE, 6 SNSR, REL, sc1000	\$	5,511.79
LXV400.99.1R172	ee PROBE MODULE, SC1000	\$	6,398.30
LXV400.99.1R282	ee MODULE, PROBE, 6 SNSR, REL, INPUT MOD, sc1000	\$	7,535.35
LXV400.99.1R382	ee MODULE, PROBE (6-SNSR), SC1000 8mA	\$	8,055.70
LXV400.99.1R572	db ee PROBE MODULE, SC1000	\$	8,267.69
LXV400.99.1R582	db ee MODULE, PROBE (6 SENSOR), SC1000	\$	8,479.68
LXV400.99.1R682	ee SC1000 PRB MOD, COND, REL, 30, 6SNS	\$	9,520.37
LXV400.99.1S082	ee SC1000 PM 6 SENS PROFIBUS PROGNSYS HACH	\$	10,156.34
LXV400.99.1U172	ee SC1000 PM 4Sens Cond 4mA OUT	\$	11,062.13
LXV400.99.1U382	ee SC1000 PM 6 SENS 2X4-20MA OUT PROGNSYS RELAY HACH	\$	12,854.42
LXV400.99.1U572	ee SC1000 PM 4 SENS 4-20MA OUT/IN PROGNSYS RELAY HACH	\$	13,047.14
LXV400.99.1U582	ee SC1000 PM 6 SENS 4-20MA OUT/IN PROGNSYS RELAY HACH	\$	13,336.22
LXV400.99.5A282	ee SC1000 PM 6Sens Cond 4mA/digit IN MODRS485 110-230V US HACH	\$	8,903.66
LXV400.99.5A322	ee SC1000 PRB MOD, US, MOD485, 2O, 4SNS	\$	7,284.82
LXV400.99.5B522	ee PROBE MODULE, SC1000	\$	9,828.72
LXV400.99.5B532	db ee MODULE, PROBE (6-SNSR), SC1000, 4mA	\$	10,156.34
LXV400.99.5G072	ee SC1000 PM 4 SENS WITH CONDUITS PROGNOYS US POWER CORD HACH	\$	9,077.11
LXV400.99.5R622	PROBE MODULE, SC1000	\$	9,828.72
LXV400.99.5R682	ee SC1000 PM 6Sens 8mAOUT 4mA/digitIN US Cord	\$	8,055.70
LXV400.99.5U582	ee SC1000 PM 6 SENS 4-20MA OUT/IN PROGNSYS REL HACH, US POWER CORD	\$	11,524.66
LXV400.99.Z0032	ee SC1000 PRB MOD, 24V, 6SNS	\$	4,182.02
LXV400.99.Z0042	ee SC1000 PRB MOD, 24V, 8SNS	\$	4,374.74
LXV400.99.Z0172	ee SC1000 PRB MOD, 24V, 10, 4X420, 4SNS	\$	5,280.53



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LXV400.99.Z0322	ee SC1000 PRB MOD 24V, 20, 4SNS	\$	6,263.40
LXV400.99.Z0332	ee SC1000 PRB MOD, 24V, 20, 6SNS	\$	6,899.38
LXV400.99.Z0382	ee SC1000 PRB MOD, 24V, 10, 4X420, 6SNS	\$	6,494.66
LXV400.99.Z0522	SC1000 Probe Module, 24V, I/O, 4 sensors	\$	6,648.84
LXV400.99.ZA022	ee SC1000 PRB MOD, 24V, 485, 4SNS	\$	5,685.24
LXV400.99.ZA032	ee SC1000 PRB MOD, 24V, 485, 6SNS	\$	5,955.05
LXV400.99.ZB522	ee SC1000 PM 4Sens 4mAOUT 4mA/digitIN	\$	10,156.34
LXV400.99.ZB532	ee SC1000 PRB MOD, 24V, REL, 485, IO, 6	\$	9,828.72
LXV400.99.ZE022	SC1000 PRB MOD, 24V, PROF DP, 4SNS	\$	5,685.24
LXV400.99.ZE032	ee SC1000 PRB MOD, 24V, PROF DP, 6SNS	\$	5,955.05
LXV400.99.ZF522	ee SC1000 PRB MOD, 24V, REL PROF DP, 4SNS	\$	10,561.06
LXV400.99.ZR122	ee SC1000 PRB MOD, 24V, REL, 10, 4SNS	\$	7,535.35
LXV400.99.ZR332	ee SC1000 PRB MOD, 24V, REL, I/O, 6SNS	\$	8,267.69
LXV400.99.ZR522	SC1000 PRB MOD, 24V, REL, I/O, 4SNS	\$	8,460.41
LXV400.99.ZR532	SC1000 PRB MOD, 24V, REL, 10 I, 6SNS	\$	8,980.75
LXV400.99.ZR572	ee SC1000 PRB MOD, 24V, REL, I/O, 4SNS	\$	8,460.41
LXV402.99.00002	db ee MODULE, DISPLAY W/O GSM, SC1000	\$	10,946.50
LXV402.99.01001	ee MODULE,sc1000 DISPLAY W/GSM EU&AUSTR	\$	12,892.97
LXV402.99.01002	db ee MODULE, sc1000 DISPLAY W/GSM. NON-US	\$	13,297.68
LXV402.99.02002	aa db ee MODULE, sc1000 DISPLAY W/GSM,US ONLY	\$	12,892.97
LXV402.99.10002	aa db ee sc1000 DISPLAY MODULE, TCP, NO GSM	\$	15,263.42
LXV402.99.11002	aa ee sc1000 DISPLAY MODULE, TCP, GSM (EU)	\$	17,209.90
LXV402.99.12002	aa db ee sc1000 DISPLAY MODULE, TCP, GSM (US)	\$	17,209.90
LXV404.99.00001	sc200 CONTROLLER, AC-DC, HACH-LANGE	\$	4,875.82
LXV404.99.00002	sc200 CONTROLLER, AC-DC, HACH	\$	4,875.82
LXV404.99.00101	ff sc200 CONTROLLER, AC-DC, pH, HACH-LANGE	\$	5,068.54
LXV404.99.00102	ff sc200 CONTROLLER, AC-DC, 1pH, HACH	\$	5,068.54
LXV404.99.00112	ff sc200 CONTROLLER, AC-DC, 2pH, HACH	\$	6,147.77
LXV404.99.00201	ff sc200 CONTROLLER,AC-DC, COND,HACH-LANGE	\$	5,068.54
LXV404.99.00202	ff sc200 CONTROLLER, AC-DC, COND, HACH	\$	5,068.54
LXV404.99.00211	ff sc200 Controller, 1Leitf., 1pH&O2	\$	6,147.77
LXV404.99.00212	ff sc200 CONTROLLER, AC-DC, COND, pH, HACH	\$	6,147.77
LXV404.99.00221	ff sc200 CONTROLLER,AC-DC, 2 COND, HACH-LANGE	\$	6,147.77
LXV404.99.00222	ff sc200 CONTROLLER, AC-DC,2 COND, HACH	\$	6,147.77
LXV404.99.00302	ff sc200 CONTROLLER, AC-DC, FLOW, HACH	\$	5,858.69
LXV404.99.00501	nn ff sc200 CONTROLLER, AC-DC, DIG, HACH-LANGE	\$	6,552.48
LXV404.99.00502	inn sc200 CONTROLLER, AC-DC, DIG,HACH	\$	6,552.48
LXV404.99.00512	ff sc200 CONTROLLER, AC-DC, DIG, pH, HACH	\$	7,631.71
LXV404.99.00522	ff sc200 CONTROLLER, AC-DC,DIG, COND, HACH	\$	7,631.71
LXV404.99.00532	ff sc200 CONTROLLER, AC-DC, DIG, FLOW, HACH	\$	8,460.41
LXV404.99.00551	aa nn ff sc200 CONTROLLER, AC-DC, 2 DIG,HACH-LANGE	\$	7,535.35
LXV404.99.00552	inn ff sc200 CONTROLLER, AC-DC, 2 DIG, HACH	\$	7,535.35
LXV404.99.01102	sc200 CONTROLLER, AC-DC, pH,MODBUS, HACH	\$	6,552.48
LXV404.99.01112	sc200 CONTROLLER, AC-DC, 2pH, MODBUS, HACH	\$	7,593.17
LXV404.99.01202	sc200 CONTROLLER, AC-DC, COND, MOD, HACH	\$	6,552.48
LXV404.99.01212	sc200 CONTROLLER, AC-DC, COND, pH, MODBUS,HACH	\$	7,593.17
LXV404.99.01222	sc 200 CONTROLLER, AC-DC, 2 COND, MODBUS, HACH	\$	7,593.17

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LXV404.99.01502	sc200 CONTROLLER,AC-DC, DIG,MODBUS,HACH	\$	7,997.88
LXV404.99.01512	sc200 CONTROLLER, AC-DC, DIG, pH,MOD, HACH	\$	9,134.93
LXV404.99.01522	sc200 CONTROLLER, AC-DC, DIG,COND,MOD, HACH	\$	9,134.93
LXV404.99.01542	sc200 CONTR, AC-DC, DIG, 4-20mA IN, MOD,HACH	\$	9,134.93
LXV404.99.01552	sc200 CONTROLLER, AC-DC,2DIG, MOD, HACH	\$	8,980.75
LXV404.99.03102	ff sc200 CONTR, 100-240V, PROFIBUS, pH, HACH	\$	6,552.48
LXV404.99.03202	ff sc200 CONTR, 100-240V, PROFIBUS, COND, HACH	\$	6,552.48
LXV404.99.03512	ff sc200 CONTR, 100-230V, DIGITAL CONN. ,pH, PROFIBUS, HACH	\$	9,134.93
LXV404.99.05102	sc200 CONTR, 100-240V, pH, HART, HACH	\$	6,552.48
LXV404.99.05112	sc200 CONTR, 100-240V, 2 pH, HART, HACH	\$	7,593.17
LXV404.99.05202	sc200 CONTR, 100-240V, COND, HART, HACH	\$	6,552.48
LXV404.99.05502	sc200 CONTROLLER, AC-DC, DIG, HART, HACH	\$	7,997.88
LXV404.99.05512	sc200 CONTR, 100-240V, 1 DIG CONN, pH, HART, HACH	\$	9,134.93
LXV404.99.05552	sc200 CONTR, 100-240V, 2 DIG CONN, HART, HACH	\$	8,980.75
LXV404.99.09112	ff sc200 CONTROLLER, AC-DC, 4-20mA OUT, 2pH, HACH	\$	7,053.55
LXV404.99.09222	ff sc200 CONTROLLER, AC-DC, 4-20mA OUT, 2 CONDUCTIVITY, HACH	\$	7,149.91
LXV404.99.09512	ff sc200 CONTR, 100-240V, 4-20mA OUTPUT, 1 DIG CONN, pH, HACH	\$	8,826.58
LXV404.99.11202	sc200 CONTR, 100-240V, 2CORD GRIP, MODBUS, COND, HACH	\$	6,860.83
LXV404.99.11212	sc200 CONTR, 100-240V, 2 CORD GRIP, MODBUS, COND,pH, HACH	\$	7,901.52
LXV404.99.11222	sc200 CONTR, 100-240V, 2 CORD GRIPS,MODBUS, 2 COND, HACH	\$	7,901.52
LXV404.99.11402	sc200 CONTR, 100-240V, 2 CORD GRIP,MODBUS, 4-20mA INPUT, HACH	\$	6,860.83
LXV404.99.11412	sc200 CONTR, 100-240V, 2 CORD GRIPS, MODBUS, 4-20mA INPUT, pH, HACH	\$	7,901.52
LXV404.99.11422	sc200 CONTR, 100-240V, 2 CORD GRIPS,MODBUS, 4-20mA INPUT, COND, HACH	\$	7,901.52
LXV404.99.11502	sc200 CONTR, 100-240V, 2 CORD GRIPS, MODBUS,1DIG CONN, HACH	\$	8,267.69
LXV404.99.11552	sc200 CONTR, 100-240V, 2 CORD GRIPS,MODBUS, 2 DIG CONN, HACH	\$	9,289.10
LXV404.99.19522	ff sc200 CONTR, 100-240V, 2 CORD GRIPS, 4-20mA OUTPUT, 1 DIG CONN, CON	\$	8,922.94
LXV404.99.20101	vv ff sc200 Controller, 1pH&O2	\$	5,299.80
LXV404.99.20201	ff sc200 Controller, 1Leitf.	\$	5,299.80
LXV404.99.20511	ff sc200 Controller, 1Digital, 1pH&O2	\$	7,882.25
LXV404.99.23201	ff sc200 Controller, 1Leitf.	\$	6,783.74
LXV404.99.23521	ff sc200 CONTR, 100-240V, 1 DIG CONN, COND,EU PWR CORD,PROFIBUS,HACH-L	\$	9,308.38
LXV404.99.29511	ff sc200 CONTR, 100-240V,4-20 mA OUT,1 DIG CONN, pH,EU PWR CORD, HACH-	\$	8,865.12
LXV404.99.50102	ff sc200 CONTR, 100-230V, pH, N AMER PWR CORD,HACH	\$	5,299.80
LXV404.99.50112	ff sc200 CONTR, 100-230V, 2pH, N AMER PWR CORD, HACH	\$	6,263.40



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LXV404.99.50222	ff sc200 CONTR, 100-230V, 2 COND, N AMER PWR CORD,HACH	\$	6,263.40
LXV404.99.51502	sc200 CONTROLLER,AC-DC, DIG,MOD,US, HACH	\$	8,132.78
LXV404.99.51552	sc200 CONTR, 100-240V, MODBUS, 2 DIG CONN, N AMER PWR CORD, HACH	\$	9,154.20
LXV404.99.70001	sc200 CONTROLLER, DC-DC, HACH-LANGE	\$	5,781.60
LXV404.99.70002	sc200 CONTROLLER, DC-DC, HACH	\$	5,781.60
LXV404.99.70102	ff sc200 CONTROLLER, pH, 24VDC, HACH	\$	5,974.32
LXV404.99.70112	ff sc200 CONTROLLER, 2pH, 24VDC, HACH	\$	7,053.55
LXV404.99.70202	ff sc200 CONTROLLER, COND, 24VDC, HACH	\$	5,974.32
LXV404.99.70212	ff sc200 CONTR, 24VDC, COND, pH, HACH	\$	7,053.55
LXV404.99.70222	ff sc200 CONTROLLER, COND, 24VDC, HACH	\$	7,053.55
LXV404.99.70512	ff sc200 CONTROLLER, DIG, pH, 24VDC, HACH	\$	8,691.67
LXV404.99.71212	sc200 CONTR, 24VDC, MODBUS, COND, pH, HACH	\$	8,672.40
LXV404.99.71442	sc200 CONTR, 24VDC, MODBUS, 2 4-20mA INPUTS, HACH	\$	8,672.40
LXV404.99.71501	sc200 CONTR, DIG ,24VDC,MODBUS, HACH-LANGE	\$	8,903.66
LXV404.99.71502	sc200 CONTROLLER, DIG, 24VDC,MOD, HACH	\$	8,903.66
LXV404.99.71522	sc200 CONTR, 24VDC, MODBUS, 1 DIG CONN, COND, HACH	\$	10,175.62
LXV404.99.71542	sc200 CONTR, 24VDC, MODBUS, 1 DIG CONN, 4-20mA INPUT, HACH	\$	10,175.62
LXV404.99.71551	vv sc200 Controller, 2Digital	\$	9,828.72
LXV404.99.71552	sc200 CONTROLLER, 2DIG, 24VDC,MODBUS, HACH	\$	9,828.72
LXV404.99.75102	sc200 CONTR, 24VDC, HART, pH, HACH	\$	7,381.18
LXV404.99.75202	sc200 CONTR, 24VDC, HART, COND, HACH	\$	7,381.18
LXV404.99.75222	sc200 CONTR, 24VDC, HART, 2 COND, HACH	\$	8,556.77
LXV404.99.75502	sc200 CONTR, 24VDC, HART, 1 DIG CONN, HACH	\$	8,903.66
LXV404.99.75552	sc200 CONTR, 24VDC, HART, 2 DIG CONN, HACH	\$	9,886.54
LXV404.99.79212	ff sc200 CONTR, 24VDC, 4-20mA OUTPUT, COND,pH, HACH	\$	8,132.78
LXV404.99.79222	ff sc200 CONTR, 24VDC, 4-20mA OUTPUT, 2 COND, HACH	\$	8,132.78
LXV404.99.DUBEC	sc200 ELECTRODELESS ANALYZER, AC/DC, DUBOIS	\$	5,993.59
LXV404.99.DUBPH	sc200 PH ANALYZER, AC/DC, DUBOIS	\$	5,955.05
LXV404.CC.00002	sc200 CONTROLLER, 100-240V, HACH, CONFORMAL COATED	\$	6,147.77
LXV404.CC.00502	sc200 CONTROLLER, 100-240V, 1 DIG CONN, HACH, CONFORMAL COATED	\$	7,959.34
LXV404.CC.00552	sc200 CONTROLLER, 100-240V, 2 DIG CONN, HACH, CONFORMAL COATED	\$	8,922.94
LXV404.UL.00002	sc200 CONTROLLER, 100-240V, HACH, UL	\$	5,511.79
LXV404.UL.00502	sc200 CONTROLLER, 100-240V, 1 DIG CONN, HACH, UL	\$	7,246.27
LXV404.UL.00552	sc200 CONTROLLER, 100-240V, 2 DIG CONN, UL, HACH	\$	8,248.42
LXV404.US.00112	sc200 CONTROLLER, AC-DC, 2pH, HACH, US	\$	6,205.58
LXV404.US.00222	sc200 CONTROLLER, AC-DC,2 COND, HACH, US	\$	6,436.85
LXV404.US.00502	sc200 CONTROLLER, 1 DIG, HACH, US	\$	6,552.48
LXV404.US.00522	sc200 CONTROLLER, DIG, COND, HACH, US	\$	7,631.71
LXV404.US.00552	sc200 CONTROLLER, 2 DIG, HACH, US	\$	8,113.51
LXV404.US.01552	sc200 CONTROLLER, AC-DC, 2DIG, MOD, HACH, US	\$	11,254.85
LXV404.US.05502	sc200 CONTROLLER, AC-DC, DIG, HART, HACH, US	\$	8,537.50
LXV404.US.05552	sc200 CONTR, 100-240V, 2 DIG CONN, HART, HACH, US	\$	6,841.56
LXV404.US.70552	sc200 CONTROLLER, 2 DIG, 24VDC, HACH, US	\$	8,479.68



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LXV404.US.75502	sc200 CONTR, 24VDC, HART, 1 DIG CONN, HACH, US	\$	9,269.83
LXV4049900552Q	SC200 CONTROLLER, AC-DC, 2 DIG, HACH	\$	7,284.82
LXV416.99.20001	LDO MODEL 2 PROBE, 10M, LANGE	\$	6,070.68
LXV416.99.20S01	LDO MODEL 2 PROBE, 10M, SEA WATER, LANGE	\$	11,023.58
LXV416.99.2UP01	KTO: LDO MODEL 2, 10M LANGE & 9253400 ADAPTER	\$	7,882.25
LXV417.99.10002	db ff NITRATAX PLUS SC 1MM	\$	66,719.66
LXV417.99.20002	db ff NITRATAX PLUS SC 2MM	\$	66,719.66
LXV417.99.50002	db ff NITRATAX PLUS SC 5 MM	\$	66,719.66
LXV418.99.10002	db UVAS sc PROBE, 1mm	\$	70,400.62
LXV418.99.20002	db UVAS sc PROBE, 2mm	\$	70,400.62
LXV418.99.50002	db UVAS sc PROBE, 5mm	\$	70,400.62
LXV418.99.90002	db UVAS sc PROBE, 50mm	\$	70,400.62
LXV420.99.50002	db NITRATAX CLEAR SC 5MM	\$	54,597.58
LXV421.53.10012	AMTAX IDSC;0.05-20;230V;PREP.F.FP	\$	58,413.43
LXV421.53.13012	AMTAX indoor sc Ammonium Analyzer, 0,05-20 mg/l NH4-N, 110-230V 1 Channel	\$	51,706.78
LXV421.53.13112	AMTAX indoor sc Ammonium Analyzer, 0,05-20 mg/l NH4-N, 110-230V 1 Channel	\$	51,706.78
LXV421.53.14012	AMTAX indoor sc Ammonium Analyzer, 0,05-20 mg/l NH4-N, 110-230V 2 Channel	\$	54,135.05
LXV421.53.14112	AMTAX IDSC;0.05-20;115V;2CH	\$	54,135.05
LXV421.53.23012	AMTAX indoor sc Ammonium Analyzer, 1-100 mg/l NH4-N, 110-230V 1 Channel	\$	50,087.93
LXV421.53.23112	AMTAX IDSC;1-100;115V;1CH	\$	50,087.93
LXV421.53.24012	AMTAX indoor sc Ammonium Analyzer, 1-100 mg/l NH4-N, 110-230V 2 Channel	\$	52,496.93
LXV421.53.24112	AMTAX indoor sc Ammonium Analyzer, 1-100 mg/l NH4-N, 110-230V 2 Channel	\$	52,496.93
LXV421.53.33012	AMTAX IDSC;10-1000;230V;1CH	\$	51,706.78
LXV421.53.33112	AMTAX IDSC;10-1000;115V;1CH	\$	51,706.78
LXV421.53.34012	AMTAX indoor sc Ammonium Analyzer, 10-1000 mg/l NH4-N, 110-230V 2 Channel	\$	54,135.05
LXV421.53.34112	AMTAX indoor sc Ammonium Analyzer, 10-1000 mg/l NH4-N, 110-230V 2 Channel	\$	54,135.05
LXV421.53.43012	AMTAX indoor sc Ammonium Analyzer, 0,02-5mg/l NH4-N, 110-230V 1 Channel	\$	51,706.78
LXV421.53.43112	AMTAX indoor sc Ammonium Analyzer, 0,02-5 mg/l NH4-N, 110-230V 1 Channel	\$	51,706.78
LXV421.53.44012	AMTAX indoor sc Ammonium Analyzer, 0,02-5 mg/l NH4-N, 110-230V 2 Channel	\$	54,135.05
LXV421.53.44112	AMTAX IDSC;0.02-5;115V;2CH	\$	54,135.05
LXV421.99.10112	AMTAX INDOOR sc MR1 (0,05-20MG/L)	\$	77,531.26
LXV421.99.11012	AMTAX W. FILTRATION PROBE 5M, 230V	\$	73,831.03
LXV421.99.11112	AMTAX W. FILTRATION PROBE 5 M (115V)	\$	74,197.20
LXV421.99.12112	AMTAX W FILTRATION PROBE 10M, 115V	\$	82,445.62
LXV421.99.13002	db AMTAX sc AMMONIA ANLZR 115-230V,1 CH	\$	58,413.43
LXV421.99.13012	AMTAX 1 CHAN CONT SAMPLE 115-230V	\$	61,323.50



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LXV421.99.14002	AMTAX sc AMMONIA ANLZR 115-230V,2 CH	\$	61,111.51
LXV421.99.14012	AMTAX 2 CHANNEL CONT. SAMPLE 230V	\$	64,079.40
LXV421.99.14112	AMTAX INDOOR 2 CHAN, 115V	\$	64,079.40
LXV421.99.21012	AMTAX INDOOR sc MR (1-100 MG/L)	\$	79,073.02
LXV421.99.21112	AMTAX INDOOR sc MB2 (1-100 MG/L) FILTRATION 5M , 115V	\$	79,073.02
LXV421.99.43002	db AMTAX SC (0.02-5 MG/L NH4-N), 1 CHAN	\$	58,297.80
LXV422.53.33001	PHOSPHAX sc LR outdoor, 110-230V, 1 CH, continuous sample prep. (EN)	\$	59,492.66
LXV422.99.10011	PHOSPHAX INDOOR sc MB1 (0,05-15 MG/L) FOR FILTRATION PROBE	\$	59,415.58
LXV422.99.10012	PHOSPHAX INDOOR sc MB 1 (0,05-15 MG/L)	\$	59,415.58
LXV422.99.10112	PHOSPHAX INDOOR sc MR1 (0,05-15MG/L), 115V	\$	57,180.02
LXV422.99.11012	PHOSPHAX INDOOR sc MR1(0,05-15MG/L), 230V	\$	74,698.27
LXV422.99.11112	PHOSPHAX INDOOR sc , MR1 (0,05-13MG/L), 5M, 115V	\$	74,698.27
LXV422.99.12012	AMTAX INDOOR sc MR1 (0,05-20MGL) W/FILTR. PROBE 10M 230V	\$	76,047.31
LXV422.99.13002	db PHOSPHAX sc PHOSP ANLZ 115-230V,1 CH	\$	58,413.43
LXV422.99.14002	PHOSPHAX sc PHOSP ANLZ 115-230V,2 CH	\$	61,111.51
LXV422.99.23002	PHOSPHAX sc PHOSP ANLZ 115-230V,1 CH	\$	58,413.43
LXV422.99.24002	PHOSPHAX sc PHOSP ANLZ 115-230V,2 CH	\$	61,111.51
LXV422.99.33001	KTO PHOSPHAX sc LR outdoor, 110-230V, 1 CH, continuous sample preparation	\$	66,257.14
LXV422.99.33011	KTO PHOSPHAX sc LR indoor, 110-230V, 1 CH, continuous sample preparation	\$	62,075.11
LXV423.99.00100	db TS-LINE sc/IMMERSION 50g/l WIPER SS	\$	17,807.33
LXV423.99.00200	db HS-LINE sc/IMMERSION 500g/l WIPER SS	\$	19,676.71
LXV423.99.02100	TS-LINE sc/IMMER 50g/l w/o WIPER SS	\$	17,460.43
LXV423.99.02200	HS-LINE sc/IMMER 500g/l w/o WIPER SS	\$	19,349.09
LXV423.99.10000	db T-LINE sc/IMMERS 4000NTU WIPER PVC	\$	14,434.73
LXV423.99.10100	db TS-LINE sc/IMMERSION 50g/l WIPER PVC	\$	17,807.33
LXV423.99.10200	db HS-LINE sc/IMMERS 500g/l WIPER PVC	\$	19,676.71
LXV423.99.12000	T-LINE sc/IMRS 4000NTU w/o WIPER PVC	\$	14,107.10
LXV423.99.12100	TS-LINE sc/IMRS 50g/l w/o WIPER PVC	\$	16,670.28
LXV423.99.12200	HS-LINE sc/IMRS 500g/l w/o WIPER PVC	\$	18,404.76
LXV424.99.00100	db INLINE sc/INSERTION 50g/l WIPER SS	\$	20,197.06
LXV424.99.00200	db HIGHLINE sc/INSERT 500g/l WIPER SS	\$	24,552.53
LXV424.99.02100	INLINE sc/INSERT 50g/l w/o WIPER SS	\$	19,927.25
LXV424.99.02200	HIGHLINE sc/INSERT 500g/l w/oWIPER SS	\$	24,321.26
LXV426.99.10001	1200-S pH sc-Sonde	\$	3,565.32
LXV426.99.20001	1200-S sc ORP probe	\$	3,758.04
LXV428.99.00001	3798-S sc Digital inductive conductivity sensor, 10m cable	\$	3,237.70
LXV429.99.01000	ff FILTERPROBE sc,5M HEATED HOSE,230VAC	\$	17,036.45
LXV429.99.01100	db ff FILTERPROBE sc,5M HEATED HOSE,115V	\$	17,383.34
LXV429.99.02000	ff FILTERPROBE sc,10M HEATED HOSE,230V	\$	20,717.40
LXV429.99.02100	db ff FILTERPROBE sc,10M HEATED HOSE,115V	\$	20,717.40
LXV430.99.00001	aa 9184sc HOCL FREE CHLORINE SENSOR	\$	13,278.41
LXV430.LM.00001	aa 9184sc HOCL FREE CHLORINE SENSOR	\$	13,278.41



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LXV431.99.00002	db SONATAX sc SLUDGE LEVEL PROBE	\$	10,676.69
LXV432.99.00001	9184sc TFC FREE CHLORINE ANALYZER	\$	18,925.10
LXV433.99.00001	9185sc O3 OZONE SENSOR	\$	13,952.93
LXV434.99.00001	9187sc CIO2 CHLORINE DIOXIDE SENSOR	\$	13,509.67
LXV440.53.10002	AISE SC W RFID (USA)	\$	24,860.88
LXV440.53.10012	AISE SC W/O RFID (USA)	\$	24,860.88
LXV440.53.20002	NISE SC W RFID (USA)	\$	24,090.00
LXV440.53.20012	NISE SC W/O RFID (USA)	\$	24,090.00
LXV440.99.00002	aa AN-ISE sc, HACH SENSOR w/RFID	\$	40,182.12
LXV440.99.00012	db AN-ISEsc, HACH, NO RFID	\$	40,721.74
LXV440.99.10002	KTO: AISE SC W RFID W CARTRIDGE	\$	29,293.44
LXV440.99.10012	KTO: AISE SC W/O RFID W CARTRIDGE	\$	29,293.44
LXV440.99.20002	KTO: NISE SC W RFID W CARTRIDGE	\$	28,484.02
LXV440.99.20012	KTO: NISE SC W/O RFID W CARTRIDGE	\$	28,484.02
LXV441.99.11102	db FP360sc, 500ppb, SS, 10M, W/O CLEAN	\$	72,828.89
LXV441.99.11202	FP360sc, 500ppb, SS 10M, W/CLEANING	\$	74,332.10
LXV441.99.11302	FP360sc, 500ppb, SS, 1.5M, W/O CLEAN	\$	72,828.89
LXV441.99.12102	FP360sc, 500ppb, Ti, 10M, W/O CLEAN	\$	82,503.43
LXV441.99.12202	FP360sc, 500ppb, Ti, 10M, W/CLEANING	\$	83,948.83
LXV441.99.12302	FP360sc, 500ppb, Ti, 1.5M, W/O CLEAN	\$	82,503.43
LXV441.99.21102	FP360sc, 5000ppb, SS, 10M, W/O CLEAN	\$	72,828.89
LXV441.99.21202	FP360sc, 5000ppb, SS, 10M, W/CLEANING	\$	74,332.10
LXV441.99.21302	FP360sc, 5000ppb, SS, 1.5M W/O CLEAN	\$	72,828.89
LXV441.99.22102	FP360sc, 5000ppb, Ti, 10M, W/O CLEAN	\$	82,503.43
LXV441.99.22202	FP360sc, 5000ppb, Ti, 10M, W/CLEANING	\$	83,948.83
LXV441.99.22302	FP360sc, 5000ppb, Ti, 1.5M, W/O CLEAN	\$	82,503.43
LXV445.99.10112	KTO:US TU5300sc TURB,EPA	\$	7,034.28
LXV445.99.10122	db ee TU5300sc TURB,ISO	\$	7,034.28
LXV445.99.10212	db ee TU5400sc TURB,EPA	\$	20,158.51
LXV445.99.10222	db ee TU5400sc TURB,ISO	\$	20,158.51
LXV445.99.11112	db ee TU5300sc TURB,SYSCHK,EPA	\$	8,113.51
LXV445.99.11122	db ee TU5300sc TURB,SYSCHK,ISO	\$	8,113.51
LXV445.99.11212	db ee TU5400sc TURB,SYSCHK,EPA	\$	21,276.29
LXV445.99.11222	db ee TU5400sc TURB,SYSCHK,ISO	\$	21,276.29
LXV445.99.12112	db ee TU5300sc TURB,RFID,EPA	\$	8,113.51
LXV445.99.12122	db ee TU5300sc TURB,RFID,ISO	\$	8,113.51
LXV445.99.12212	db ee TU5400sc TURB,RFID,EPA	\$	21,276.29
LXV445.99.12222	db ee TU5400sc TURB,RFID,ISO	\$	21,276.29
LXV445.99.13112	db ee TU5300sc TURB,SYSCHK,RFID,EPA	\$	9,269.83
LXV445.99.13122	db ee TU5300sc TURB,SYSCHK,RFID,ISO	\$	9,269.83
LXV445.99.13212	db ee TU5400sc TURB,SYSCHK,RFID,EPA	\$	22,374.79
LXV445.99.13222	db ee TU5400sc TURB,SYSCHK,RFID,ISO	\$	22,374.79
LXV445.99.20112	ee KTO:US TU5300sc TURB,FLOW,EPA	\$	7,477.54
LXV445.99.20122	ee KTO:US TU5300sc TURB,FLOW,ISO	\$	7,477.54
LXV445.99.20212	ee KTO:US TU5400sc TURB,FLOW,EPA	\$	20,717.40
LXV445.99.20222	db ee TU5400sc TURB,FLOW,ISO	\$	20,717.40
LXV445.99.21112	ee KTO:US TU5300sc TURB,FLOW,SYSCHK,EPA	\$	8,556.77

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LXV445.99.21122	ee KTO:US TU5300sc TURB,FLOW,SYSCCHK,ISO	\$ 8,556.77
LXV445.99.21212	ee KTO:US TU5400sc TURB,FLOW,SYSCCHK,EPA	\$ 21,873.72
LXV445.99.21222	ee KTO:US TU5400sc TURB,FLOW,SYSCCHK,ISO	\$ 21,873.72
LXV445.99.22112	ee KTO:US TU5300sc TURB,FLOW,RFID,EPA	\$ 8,556.77
LXV445.99.22122	ee KTO:US TU5300sc TURB,FLOW,RFID,ISO	\$ 8,556.77
LXV445.99.22212	ee KTO:US TU5400sc TURB,FLOW,RFID,EPA	\$ 21,873.72
LXV445.99.22222	db ee TU5400sc TURB,FLOW,RFID,ISO	\$ 21,873.72
LXV445.99.23112	ee KTO: US TU5300sc TURB,FLOW,SYSCCHK,RFID,EPA	\$ 9,713.09
LXV445.99.23122	ee KTO: US TU5300sc TURB,FLOW,SYSCCHK,RFID,ISO	\$ 9,713.09
LXV445.99.23212	ee KTO:US TU5400sc TURB,FLOW,SYSCCHK,RFID,EPA	\$ 22,972.22
LXV445.99.23222	ee KTO:US TU5400sc TURB,FLOW,SYSCCHK,RFID,ISO	\$ 22,972.22
LXV445.99.30112	ee KTO: US TU5300sc TURB,CLEAN,EPA	\$ 9,308.38
LXV445.99.30122	ee KTO:US TU5300sc TURB,CLEAN,ISO	\$ 9,308.38
LXV445.99.30212	ee KTO:US TU5400sc TURB,CLEAN,EPA	\$ 22,625.33
LXV445.99.30222	db ee TU5400sc TURB,CLEAN,ISO	\$ 22,625.33
LXV445.99.31112	ee KTO:US TU5300sc TURB,CLEAN,SYSCCHK,EPA	\$ 10,464.70
LXV445.99.31122	ee KTO:US TU5300sc TURB,CLEAN,SYSCCHK,ISO	\$ 10,464.70
LXV445.99.31212	db ee TU5400sc TURB,CLEAN,SYSCCHK,EPA	\$ 23,704.56
LXV445.99.31222	ee KTO:US TU5400sc TURB,CLEAN,SYSCCHK,ISO	\$ 23,704.56
LXV445.99.32112	ee KTO:US TU5300sc TURB,CLEAN,RFID,EPA	\$ 10,464.70
LXV445.99.32122	ee KTO:US TU5300sc TURB,CLEAN,RFID,ISO	\$ 10,464.70
LXV445.99.32212	ee KTO:US TU5400sc TURB,CLEAN,RFID,EPA	\$ 23,704.56
LXV445.99.32222	db ee TU5400sc TURB,CLEAN,RFID,ISO	\$ 23,704.56
LXV445.99.33112	ee KTO:US TU5300sc TURB,CLEAN,SYSCCHK,RFID,EPA	\$ 11,524.66
LXV445.99.33122	ee KTO:US TU5300sc TURB,CLEAN,SYSCCHK,RFID,ISO	\$ 11,524.66
LXV445.99.33212	ee KTO:US TU5400sc TURB,CLEAN,SYSCCHK,RFID,EPA	\$ 24,803.06
LXV445.99.33222	db ee TU5400sc TURB,CLEAN,SYSCCHK,RFID,ISO	\$ 24,803.06
LXV445.99.50112	ee KTO: US TU5300sc TURB,FLOW,CLEAN,EPA	\$ 9,905.81
LXV445.99.50122	ee KTO:US TU5300sc TURB,FLOW,CLEAN,ISO	\$ 9,905.81
LXV445.99.50212	db ee TU5400sc TURB,FLOW,CLEAN,EPA	\$ 23,107.13
LXV445.99.50222	db ee TU5400sc TURB,FLOW,CLEAN,ISO	\$ 23,107.13
LXV445.99.51112	ee KTO:US TU5300sc TURB,FLOW,CLEAN,SYSCCHK,EPA	\$ 11,023.58
LXV445.99.51122	ee KTO:US TU5300sc TURB,FLOW,CLEAN,SYSCCHK,ISO	\$ 11,023.58
LXV445.99.51212	ee KTO:US TU5400sc TURB,FLOW,CLEAN,SYSCCHK,EPA	\$ 24,321.26
LXV445.99.51222	db ee TU5400sc TURB,FLOW,CLEAN,SYSCCHK,ISO	\$ 24,321.26
LXV445.99.52112	ee KTO:US TU5300sc TURB,FLOW,CLEAN,RFID,EPA	\$ 11,023.58
LXV445.99.52122	ee KTO:US TU5300sc TURB,FLOW,CLEAN,RFID,ISO	\$ 11,023.58
LXV445.99.52212	db ee TU5400sc TURB,FLOW,CLEAN,RFID,EPA	\$ 24,321.26
LXV445.99.52222	db ee TU5400sc TURB,FLOW,CLEAN,RFID,ISO	\$ 24,321.26
LXV445.99.53112	ee KTO: US TU5300sc TURB,FLOW,CLEAN,SYSCCHK,RFID,EPA	\$ 9,905.81
LXV445.99.53122	ee KTO: US TU5300sc TURB,FLOW,CLEAN,SYSCCHK,RFID,ISO	\$ 9,905.81
LXV445.99.53212	ee KTO: US TU5400sc TURB,FLOW,CLEAN,SYSCCHK,RFID,EPA	\$ 23,107.13
LXV445.99.53222	ee KTO:US TU5400sc TURB,FLOW,CLEAN,SYSCCHK,RFID,ISO	\$ 23,107.13
LXV446.99.101L1	SC1500; 2 SENS 4mA OUT 110V/COND EXT MOD	\$ 8,479.68
LXV446.99.101M1	SC1500; 4 SENS 4mA OUT 110V/COND EXT MOD	\$ 8,730.22
LXV446.99.103L1	SC1500; 2 SENS 8mA OUT 110V/COND EXT MOD	\$ 10,002.17
LXV446.99.103M1	SC1500; 4 SENS 8mA OUT 110V/COND EXT MOD	\$ 10,156.34



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LXV446.99.103N1	SC1500; 6 SENS 8mA OUT 110V/COND EXT MOD	\$	10,368.34
LXV446.99.1R1L1	SC1500; 2 SENS 4mA OUT 110V/COND 4 REL EXT MOD	\$	10,002.17
LXV446.99.1R1M1	SC1500; 4 SENS 4mA OUT 110V/COND 4 REL EXT MOD	\$	10,156.34
LXV446.99.1R1Q1	SC1500; 2 SENS 4mA OUT 110V/COND 4 REL/C EXT MOD	\$	10,002.17
LXV446.99.1R1R1	SC1500; 4 SENS 4mA OUT 110V/COND 4 REL/C EXT MOD	\$	10,156.34
LXV446.99.1R3L1	SC1500; 2 SENS 8mA OUT 110V/COND 4 REL EXT MOD	\$	11,409.02
LXV446.99.1R3M1	SC1500; 4 SENS 8mA OUT 110V/COND 4 REL EXT MOD	\$	11,524.66
LXV446.99.1R3N1	SC1500; 6 SENS 8mA OUT 110V/COND 4 REL EXT MOD	\$	11,736.65
LXV446.99.1R3Q1	SC1500; 2 SENS 8mA OUT 110V/COND 4 REL/C EXT MOD	\$	11,409.02
LXV446.99.1R3R1	SC1500; 4 SENS 8mA OUT 110V/COND 4 REL/C EXT MOD	\$	11,524.66
LXV446.99.1R3S1	SC1500; 6 SENS 8mA OUT 110V/COND 4 REL/C EXT MOD	\$	11,736.65
LXV446.99.201L1	SC1500; 2 SENS 4mA OUT 110-230V EU EXT MOD	\$	8,479.68
LXV446.99.201M1	SC1500; 4 SENS 4mA OUT 110-230V EU EXT MOD	\$	8,730.22
LXV446.99.203L1	SC1500; 2 SENS 8mA OUT 110-230V EU EXT MOD	\$	10,002.17
LXV446.99.203M1	SC1500; 4 SENS 8mA OUT 110-230V EU EXT MOD	\$	10,156.34
LXV446.99.203N1	SC1500; 6 SENS 8mA OUT 110-230V EU EXT MOD	\$	10,368.34
LXV446.99.2R1L1	SC1500; 2 SENS 4mA OUT 110-230V EU 4 REL EXT MOD	\$	10,002.17
LXV446.99.2R1M1	SC1500; 4 SENS 4mA OUT 110-230V EU 4 REL EXT MOD	\$	10,156.34
LXV446.99.2R3L1	SC1500; 2 SENS 8mA OUT 110-230V EU 4 REL EXT MOD	\$	11,409.02
LXV446.99.2R3M1	SC1500; 4 SENS 8mA OUT 110-230V EU 4 REL EXT MOD	\$	11,524.66
LXV446.99.2R3N1	SC1500; 6 SENS 8mA OUT 110-230V EU 4 REL EXT MOD	\$	11,736.65
LXV448.99.11001	db ee NT3100sc - 1mm (NO3) Nitrate Sensor	\$	36,269.90
LXV448.99.12001	db ee NT3200sc - 1mm (NO3+NO2) Nitrate/Nitrite Sensor	\$	44,171.42
LXV448.99.21001	db ee NT3100sc - 2mm (NO3) Nitrate Sensor	\$	36,269.90
LXV448.99.22001	db ee NT3200sc - 2mm (NO3+NO2) Nitrate/Nitrite Sensor	\$	44,171.42
LXV448.99.51001	db ee NT3100sc - 5mm (NO3) Nitrate Sensor	\$	36,269.90
LXV448.99.52001	db ee NT3200sc - 5mm (NO3+NO2) Nitrate/Nitrite Sensor	\$	44,171.42
LXV449.97.01010	CAx440EX FIELD TRANSMITTER H2S EX, WIRELESS 4G (US & CA)	\$	7,728.07
LXV449.99.02000	CDx440sc DIGITAL SC GATEWAY	\$	3,122.06
LXV449.99.10000	GS1440 SENSOR H2S	\$	17,113.54
LXV449.99.20000	GS2440EX SENSOR H2S EX	\$	20,563.22
LXV45A.99.11022	aa KTO: CLF10sc, GRAB SAMPLE, PANEL ONLY	\$	16,940.09
LXV45A.99.12022	KTO: CLF10sc, COMBO pH, PANEL ONLY	\$	17,653.15
LXV45A.99.13022	aa KTO: CLF10sc, pH, PANEL ONLY	\$	20,158.51
LXV45A.99.21022	aa KTO: CLF10sc, GRB SMPL, METRIC, PANEL ONLY	\$	14,261.28
LXV45A.99.23022	aa KTO: CLF10sc, pH, METRIC, PANEL ONLY	\$	20,158.51
LXV45A.LM.11022	aa KTO: CLF10sc, GRAB SAMPLE, PANEL ONLY	\$	16,940.09
LXV45A.LM.13022	aa KTO: CLF10sc, pH, PANEL ONLY	\$	20,158.51
LXV45B.99.11022	aa KTO: CLT10sc, GRAB SAMPLE, PANEL ONLY	\$	16,940.09
LXV45B.99.12022	KTO: CLT10sc, COMBO pH, PANEL ONLY	\$	17,653.15
LXV45B.99.13022	aa KTO: CLT10sc, pH, PANEL ONLY	\$	20,158.51
LXV510.99.10000	AF7000 SCM, 110V 60Hz, US CORD	\$	47,274.22
LXV510.99.11000	AF7000 SCM, AUTO FLSH, 110V 60 Hz, US CORD	\$	51,571.87
LXV510.99.11010	AF7000 SCM, AUTOFLSH, GRT FLTR, 110V 60 Hz, US CORD	\$	53,402.71
LXV510.99.20000	AF7000 SCM, 230V 50Hz, EU CORD	\$	47,274.22
LXV510.99.22000	AF7000 SCM, AUTOFLSH, 230V 50 Hz, EU CORD	\$	51,571.87
LXV510.99.22010	AF7000 SCM, AUTOFLSH, GRT FLTR, 230V 50Hz, EU CORD	\$	53,402.71

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LXV510.99.30000	AF7000 SCM, 110V 50HZ, UK CORD	\$	43,516.18
LXV515.53.0003B	RTC , 15" touch screen (Beckhoff)	\$	20,139.24
LXV515.99.0001B	RTC , DIN Rail IPC without UI and Basic SW (Beckhoff)	\$	12,218.45
LXV515.99.0002B	RTC , 12" touch screen (Beckhoff)	\$	19,522.54
LXV515.99.0003B	RTC , 15" touch screen (Beckhoff)	\$	20,370.50
LXV515.99.0003C	RTC , 15" touch screen (Siemens)	\$	19,850.16
LXV515.99.0004B	RTC , 19" touch screen (Beckhoff)	\$	21,141.38
LXV515.99.0004C	RTC , 19" touch screen (Siemens)	\$	20,582.50
LXV515.99.0005B	RTC , DIN Rail IPC with UI and Basic SW	\$	14,569.63
LXV515.99.0005D	RTC, 15.6" Fanless Panel PC (ADVANTECH)	\$	14,781.62
LXV524.99.000A0	SC4200c w/o plug CC	\$	2,235.55
LXV524.99.010A0	SC4200c w/o plug 3G/4G US CC	\$	2,235.55
LXV524.99.030A0	SC4200c w/o plug 3G/4G EU CC	\$	2,235.55
LXV524.99.040A0	SC4200c w/o plug USB CC	\$	2,235.55
LXV524.99.100A0	SC4200c w/ plug CC	\$	2,235.55
LXV524.99.110A0	SC4200c w/ plug 3G/4G US ATT CC	\$	2,235.55
LXV524.99.140A0	SC4200c w/ plug USB CC	\$	2,235.55
LXV525.32KTO025	KTO: CLF10sc, sc4500 DUAL INPUT, GRB SMPLE for Vietnam	\$	14,762.35
LXV525.97KTO001	KTO:áCLF10sc,ásc4500áSINGLÉáINPUT,ápHD	\$	19,715.26
LXV525.97KTO002	KTO:ááCLF10ásc,ásc4500áSGLááINPUT,áCOMBOápH	\$	18,212.04
LXV525.97KTO003	KTO:áSS7áW/sc4500,áCHá1	\$	21,179.93
LXV525.97KTO004	KTO:áCLF10sc,ásc4500áSINGLÉáINPUT,áGRBáSMPLE	\$	17,209.90
LXV525.97KTO005	KTO:áTU5300scáw/sc4500á1áCH	\$	12,006.46
LXV525.97KTO006	KTO:áTU5300sc,áSysáChk,áRFID,áFlowáw/sc4500á1áCH	\$	14,473.27
LXV525.97KTO007	KTO: TU5300sc, Sys Chk RFID Flow ACM w/sc4500 1 CH	\$	14,646.72
LXV525.97KTO008	KTO:áTU5300scá(ISO)áw/sc4500á1áCH	\$	11,370.48
LXV525.97KTO009	KTO:áTU5300scá(ISO),áSysáChk,áRFID,áFlowáw/sc4500á1áCH	\$	14,588.90
LXV525.97KTO010	KTO:áTU5300scá(ISO),áSysáChkáRFIDáFlowáACMáw/sc4500á1áCH	\$	14,781.62
LXV525.97KTO011	KTO:áTU5300scáACMáw/sc4500á1áCH	\$	14,126.38
LXV525.97KTO012	KTO:áTU5300scáFlow,áACMáw/sc4500á1áCH	\$	14,646.72
LXV525.97KTO013	KTO:áá1720Eáw/sc4500áCHá1	\$	12,449.71
LXV525.97KTO014	KTO:áá1720Eáw/sc4500áCHá1	\$	12,449.71
LXV525.97KTO015	KTO:áSS7áW/sc4500,áCHá1	\$	21,179.93
LXV525.97KTO016	KTO:áá1720Eáw/sc4500áCHá1	\$	12,449.71
LXV525.97KTO017	KTO:ááUVASáscáPROBE,á50MM,áW/sc4500	\$	66,314.95
LXV525.97KTO018	KTO: SOLITAX sc SS WIPER W/sc4500	\$	20,852.30
LXV525.97KTO019	KTO:á1720Eáw/sc4500,á2CH	\$	12,642.43
LXV525.97KTO020	KTO:áCLF10sc,ásc4500áDUALááINPUT,ápHD	\$	19,927.25
LXV525.97KTO021	KTO:ááUVASáscáPROBE,á5MM,áW/sc4500	\$	66,276.41
LXV525.97KTO022	KTO:áSS7áW/sc4500áSTANDARD	\$	21,372.65
LXV525.97KTO023	KTO:áCLT10sc,ásc4500áDUALááINPUTáCOMBOápH	\$	18,250.58
LXV525.97KTO024	KTO:áCLF10sc,ásc4500áDUALááINPUT,áCOMBOápH	\$	18,424.03
LXV525.97KTO025	KTO:áCLF10sc,ásc4500áDUALááINPUT,áGRBáSMPLE	\$	17,421.89
LXV525.97KTO026	KTO:áSOLITAXááINLINEáscáSSáWIPER/sc4500	\$	32,820.22
LXV525.97KTO027	KTO:ááUVASáscáPROBE,á1MM,áW/sc4500	\$	66,314.95
LXV525.97KTO029	KTO:áCLT10sc,ásc4500áDUALááINPUT,áGRBáSMPLE	\$	17,286.98

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LXV525.97KTO030	KTO:ááUVASáscáPROBE,á2MM,áW/sc4500	\$	66,314.95
LXV525.97KTO032	KTO:áHIGHLINEáscáSSáWIPERáW/sc4500	\$	36,385.54
LXV525.97KTO033	KTO:áLDOáMODELá2,á10MáHACHá&á9253400áADAPTERá&ásc4500á2áCHáDIGITALáCONTROLLER	\$	12,488.26
LXV525.97KTO034	KTO:ááSONATAXáscáW/sc4500áSTANDARD	\$	17,441.16
LXV525.97KTO035	KTO:áSOLITAXáHSáscáSSáWIPERáW/sc4500	\$	22,471.15
LXV525.97KTO036	KTO:áTU5300scáw/sc4500á2áCH	\$	12,218.45
LXV525.97KTO037	KTO:áTU5300sc,áSysáChkáRFIDáFlowáACMáw/sc4500á2áCH	\$	14,877.98
LXV525.97KTO038	KTO:áTU5400sc,áSysáChk,áRFID,áFlowáw/sc4500á2áCH	\$	26,287.01
LXV525.97KTO039	KTO:áTU5400sc,áSysáChkáRFIDáFlowáACMáw/sc4500á2áCH	\$	26,460.46
LXV525.97KTO040	KTO:áTU5300sc,áSysáChk,áRFID,áFlowáw/sc4500á2áCH	\$	14,665.99
LXV525.97KTO041	KTO:áTU5300scá(ISO)áw/sc4500á2áCH	\$	12,314.81
LXV525.97KTO042	KTO:áTU5300scá(ISO),áSysáChkáRFIDáFlowáACMáw/sc4500á2áCH	\$	14,993.62
LXV525.97KTO043	KTO:áTU5300scá(ISO),áSysáChk,áRFID,áFlowáw/sc4500á2áCH	\$	14,800.90
LXV525.97KTO044	KTO:áá1720EáW/sc4500á24VDCáCHá1	\$	12,584.62
LXV525.97KTO045	KTO:áTU5300scáw/sc4500á2áCHá24V	\$	12,372.62
LXV525.97KTO046	KTO:áTU5300scáSysáChkáRFIDáFlowáACMásc4500á2áCHá24V	\$	15,032.16
LXV525.97KTO047	KTO:áTU5300scá(ISO)SysáChkáRFIDáFlowáw/sc4500á2áCHá24V	\$	14,993.62
LXV525.99A01501	SC4500 w/o plug mAOutput 1DigSensor	\$	7,284.82
LXV525.99A01551	SC4500 w/oPlug mAOutput 2 digital Sensors	\$	7,496.81
LXV525.99A01601	SC4500 w/o plug mAOutput 1UPW COND	\$	7,284.82
LXV525.99A01661	SC4500 w/o plug mAOutput 2UPW COND	\$	7,959.34
LXV525.99A01701	SC4500 w/o plug mAOutput 1UPW pH/ORP	\$	7,169.18
LXV525.99A01771	SC4500 w/o plug mAOutput 2UPW pH/ORP	\$	7,747.34
LXV525.99A02501	SC4500 w/o plug ProfibusDP 1DigSensor	\$	9,115.66
LXV525.99A02551	SC4500 w/o plug ProfibusDP 2DigSensors	\$	9,327.65
LXV525.99A02601	SC4500 w/o plug ProfibusDP 1UPW COND	\$	9,000.02
LXV525.99A02661	SC4500 w/o plug ProfibusDP 2UPW COND	\$	9,674.54
LXV525.99A02701	SC4500 w/o plug ProfibusDP 1UPW pH/ORP	\$	8,884.39
LXV525.99A02771	SC4500 w/o plug ProfibusDP 2UPW pH/ORP	\$	9,462.55
LXV525.99A05501	SC4500 w/o plug Modbus TCP+LAN 1DigSensor	\$	9,115.66
LXV525.99A05551	SC4500 w/o plug Modbus TCP+LAN 2DigSensors	\$	9,327.65
LXV525.99A0B551	SC4500 w/o plug Network mA Output 2 digital Sensors	\$	9,327.65
LXV525.99A0C551	SC4500 w/oPlug LAN+Profibus 2 digital Sensors	\$	9,327.65
LXV525.99A0F541	SC4500 w/o plug LAN+Profinet 1DigSensor 1mAInput	\$	10,175.62
LXV525.99A0F551	SC4500 w/o plug LAN+Profinet 2DigSensors	\$	10,252.70
LXV525.99A0G541	SC4500 w/o plug LAN+Ethernet IP 1DigSensor 1mAInput	\$	10,175.62
LXV525.99A0G551	SC4500 w/o plug LAN+Ethernet IP 2DigSensors	\$	10,252.70
LXV525.99A11101	SC4500 w/o plug Prognosys mAOutput 1pH/ORP	\$	6,629.57
LXV525.99A11111	SC4500 w/o plug Prognosys mAOutput 2pH/ORP	\$	6,957.19
LXV525.99A11211	SC4500 w/o plug Prognosys mAOutput 1COND 1pH/ORP	\$	6,957.19
LXV525.99A11221	SC4500 w/o plug Prognosys mAOutput 2COND	\$	6,957.19
LXV525.99A11401	SC4500 w/o plug Prognosys mAOutput 1mAInput	\$	7,207.73
LXV525.99A11441	SC4500 w/o plug Prognosys mAOutput 2mAInputs	\$	7,323.36
LXV525.99A11501	SC4500 w/o plug Prognosys mAOutput 1DigSensor	\$	7,284.82
LXV525.99A11521	SC4500 w/o plug Prognosys mAOutput 1DigSensor 1COND	\$	6,668.11

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LXV525.99A11541	SC4500 w/o plug Prognosys mAOutput 1DigSensor 1mAInput	\$	7,381.18
LXV525.99A11551	SC4500 w/o plug Prognosys mAOutput 2DigSensors	\$	7,496.81
LXV525.99A11601	SC4500 w/o plug Prognosys mAOutput 1UPW COND	\$	7,284.82
LXV525.99A11661	SC4500 w/o plug Prognosys mAOutput 2UPW COND	\$	7,959.34
LXV525.99A11701	SC4500 w/o plug Prognosys mAOutput 1UPW pH/ORP	\$	7,169.18
LXV525.99A11761	SC4500 w/o plug Prognosys mAOutput 1UPW pH/ORP 1UPW COND	\$	7,862.98
LXV525.99A11771	SC4500 w/o plug Prognosys mAOutput 2UPW pH/ORP	\$	7,747.34
LXV525.99A12211	SC4500 w/o plug Prognosys ProfibusDP 1COND 1pH/ORP	\$	9,000.02
LXV525.99A12401	SC4500 w/o plug Prognosys ProfibusDP 1mAInput	\$	9,057.84
LXV525.99A12441	SC4500 w/o plug Prognosys ProfibusDP 2mAInputs	\$	9,134.93
LXV525.99A12501	SC4500 w/o plug Prognosys ProfibusDP 1DigSensor	\$	9,115.66
LXV525.99A12541	SC4500 w/o plug Prognosys ProfibusDP 1DigSensor 1mAInput	\$	9,269.83
LXV525.99A12551	SC4500 w/o plug Prognosys ProfibusDP 2DigSensors	\$	9,327.65
LXV525.99A12601	SC4500 w/o plug Prognosys ProfibusDP 1UPW COND	\$	9,000.02
LXV525.99A12661	SC4500 w/o plug Prognosys ProfibusDP 2UPW COND	\$	9,674.54
LXV525.99A12701	SC4500 w/o plug Prognosys ProfibusDP 1UPW pH/ORP	\$	8,884.39
LXV525.99A12761	SC4500 w/o plug Prognosys ProfibusDP 1UPW pH/ORP 1UPW COND	\$	9,578.18
LXV525.99A12771	SC4500 w/o plug Prognosys ProfibusDP 2UPW pH/ORP	\$	9,462.55
LXV525.99A15211	SC4500 w/o plug Prognosys Modbus TCP+LAN 1COND 1pH/ORP	\$	8,113.51
LXV525.99A15501	SC4500 w/o plug Prognosys Modbus TCP+LAN 1DigSensor	\$	9,115.66
LXV525.99A15541	SC4500 w/o plug Prognosys Modbus TCP+LAN 1DigSensor 1mAInput	\$	7,458.26
LXV525.99A15551	SC4500 w/o plug Prognosys Modbus TCP+LAN 2DigSensors	\$	9,327.65
LXV525.99A15601	SC4500 w/o plug Prognosys Modbus TCP+LAN 1UPW COND	\$	8,113.51
LXV525.99A15661	SC4500 w/o plug Prognosys Modbus TCP+LAN 2UPW COND	\$	8,788.03
LXV525.99A15701	SC4500 w/o plug Prognosys Modbus TCP+LAN 1UPW pH/ORP	\$	7,997.88
LXV525.99A15761	SC4500 w/o plug Prognosys Modbus TCP+LAN 1UPW pH/ORP 1UPW COND	\$	8,691.67
LXV525.99A15771	SC4500 w/o plug Prognosys Modbus TCP+LAN 2UPW pH/ORP	\$	8,576.04
LXV525.99A16211	SC4500 w/o plug Prognosys Profinet 1COND 1pH/ORP	\$	8,884.39
LXV525.99A16221	SC4500 w/o plug Prognosys Profinet 2Conductivity	\$	8,884.39
LXV525.99A16441	SC4500 w/o plug Prognosys Profinet 2mAInputs	\$	9,770.90
LXV525.99A16501	SC4500 w/o plug Prognosys Profinet 1DigSensor	\$	9,809.45
LXV525.99A16541	SC4500 w/o plug Prognosys Profinet 1DigSensor 1mAInput	\$	9,713.09
LXV525.99A16551	SC4500 w/o plug Prognosys Profinet 2DigSensors	\$	10,252.70
LXV525.99A17441	SC4500 w/o plug Prognosys Ethernet IP 2mAInputs	\$	9,770.90
LXV525.99A17541	SC4500 w/o plug Prognosys Ethernet IP 1DigSensor 1mAInput	\$	9,713.09
LXV525.99A17551	SC4500 w/o plug Prognosys Ethernet IP 2DigSensors	\$	9,616.73
LXV525.99A1A501	SC4500 w/o plug Prognosys LAN 1DigSensor	\$	6,803.02
LXV525.99A1B441	SC4500 w/o plug Prognosys LAN+mAOutput 2mAInputs	\$	7,959.34
LXV525.99A1B501	SC4500 w/o plug Prognosys LAN+mAOutput 1DigSensor	\$	6,803.02
LXV525.99A1B541	SC4500 w/o plug Prognosys LAN+mAOutput 1 dig Sensor 1mA Input	\$	9,269.83
LXV525.99A1B551	SC4500 w/o plug Prognosys LAN+mAOutput 2 digital Sensors	\$	9,327.65
LXV525.99A1B661	SC4500 w/o plug Prognosys LAN+mAOutput 2UPW COND	\$	8,171.33



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LXV525.99A1B771	SC4500 w/o plug Prognosys LAN+mAOutput 2UPW pH/ORP	\$	7,959.34
LXV525.99A1C501	SC4500 w/o plug Prognosys LAN+ProfibusDP 1DigSensor	\$	8,036.42
LXV525.99A1C551	SC4500 w/o plug Prognosys LAN+Profibus 2 digital Sensors	\$	9,327.65
LXV525.99A1F441	SC4500 w/o plug Prognosys LAN+Profinet 2mAInputs	\$	9,770.90
LXV525.99A1F541	SC4500 w/o plug Prognosys LAN+Profinet 1DigSensor 1mAInput	\$	10,175.62
LXV525.99A1F551	SC4500 w/o plug Prognosys LAN+Profinet 2DigSensors	\$	10,252.70
LXV525.99A1G441	SC4500 w/o plug Prognosys LAN+Ethernet IP 2mAInputs	\$	9,770.90
LXV525.99A1G541	SC4500 w/o plug Prognosys LAN+Ethernet IP 1DigSensor 1mAInput	\$	10,175.62
LXV525.99A1G551	SC4500 w/o plug Prognosys LAN+Ethernet IP 2DigSensors	\$	10,252.70
LXV525.99AA0441	SC4500 w/o plug Claros No output 2mAInputs	\$	6,841.56
LXV525.99AA0541	SC4500 w/o plug Claros No output 1DigSensor 1mAInput	\$	6,668.11
LXV525.99AA0551	SC4500 w/o plug Claros No output 2DigSensors	\$	6,475.39
LXV525.99AA1001	SC4500 w/o plug Claros mAOutput	\$	6,591.02
LXV525.99AA1101	SC4500 w/o plug Claros mAOutput 1pH/ORP	\$	6,629.57
LXV525.99AA1111	SC4500 w/o plug Claros mAOutput 2pH/ORP	\$	6,957.19
LXV525.99AA1401	SC4500 w/o plug Claros mAOutput 1mAInput	\$	6,591.02
LXV525.99AA1441	SC4500 w/o plug Claros mAOutput 2mAInputs	\$	6,706.66
LXV525.99AA1501	SC4500 w/o plug Claros mAOutput 1DigSensor	\$	6,687.38
LXV525.99AA1511	SC4500 w/o plug Claros mAOutput 1DigSensor 1pH/ORP	\$	6,706.66
LXV525.99AA1521	SC4500 w/o plug Claros mAOutput 1DigSensor 1COND	\$	6,706.66
LXV525.99AA1541	SC4500 w/o plug Claros mAOutput 1DigSensor 1mAInput	\$	6,783.74
LXV525.99AA1551	SC4500 w/o plug Claros mAOutput 2DigSensors	\$	6,899.38
LXV525.99AA1601	SC4500 w/o plug Claros mAOutput 1UPW COND	\$	7,284.82
LXV525.99AA1741	SC4500 w/o plug Claros mAOutput 1UPW pH/ORP 1mAInput	\$	7,458.26
LXV525.99AA2111	SC4500 w/o plug Claros ProfibusDP 2pH/ORP	\$	9,000.02
LXV525.99AA2221	SC4500 w/o plug Claros ProfibusDP 2Conductivity	\$	9,000.02
LXV525.99AA2501	SC4500 w/o plug Claros ProfibusDP 1DigSensor	\$	8,498.95
LXV525.99AA2521	SC4500 w/o plug Claros ProfibusDP 1DigSensor 1COND	\$	4,374.74
LXV525.99AA2551	SC4500 w/o plug Claros ProfibusDP 2DigSensors	\$	8,710.94
LXV525.99AA5111	SC4500 w/o plug Claros Modbus TCP+LAN 2pH/ORP	\$	8,113.51
LXV525.99AA5221	SC4500 w/o plug Claros Modbus TCP+LAN 2Conductivity	\$	8,113.51
LXV525.99AA5441	SC4500 w/o plug Claros Modbus TCP+LAN 2mAInputs	\$	7,959.34
LXV525.99AA5501	SC4500 w/o plug Claros Modbus TCP+LAN 1DigSensor	\$	8,498.95
LXV525.99AA5541	SC4500 w/o plug Claros Modbus TCP+LAN 1DigSensor 1mAInput	\$	7,862.98
LXV525.99AA5551	SC4500 w/o plug Claros Modbus TCP+LAN 2DigSensors	\$	8,710.94
LXV525.99AA6501	SC4500 w/o plug Claros Profinet 1DigSensor	\$	9,366.19
LXV525.99AA6541	SC4500 w/o plug Claros Profinet 1DigSensor 1mAInput	\$	9,539.64
LXV525.99AA7501	SC4500 w/o plug Claros Ethernet IP 1DigSensor	\$	8,383.32
LXV525.99AA7541	SC4500 w/o plug Claros Ethernet IP 1DigSensor 1mAInput	\$	9,539.64
LXV525.99AAA551	SC4500 w/o plug Claros LAN 2DigSensors	\$	8,633.86
LXV525.99AAB441	SC4500 w/o plug Claros LAN+mAOutput 2mAInputs	\$	7,072.82
LXV525.99AAB551	SC4500 w/o plug Claros LAN+mAOutput 2DigSensors	\$	8,710.94
LXV525.99AAB701	SC4500 w/o plug Claros LAN+mAOutput 1UPW pH/ORP	\$	7,381.18
LXV525.99AAB761	SC4500 w/o plug Claros LAN+mAOutput 1UPW pH/ORP 1UPW COND	\$	8,074.97
LXV525.99AAB771	SC4500 w/o plug Claros LAN+mAOutput 2UPW pH/ORP	\$	7,959.34



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LXV525.99AAC551	SC4500 w/o plug Claros LAN+ProfibusDP 2DigSensors	\$	8,633.86
LXV525.99AAC701	SC4500 w/o plug Claros LAN+ProfibusDP 1UPW pH/ORP	\$	8,884.39
LXV525.99AAC771	SC4500 w/o plug Claros LAN+ProfibusDP 2UPW pH/ORP	\$	9,462.55
LXV525.99AAF501	SC4500 w/o plug Claros LAN+Profinet 1DigSensor	\$	8,383.32
LXV525.99AAF541	SC4500 w/o plug Claros LAN+Profinet 1DigSensor 1mAInput	\$	8,653.13
LXV525.99AAF551	SC4500 w/o plug Claros LAN+Profinet 2DigSensors	\$	9,481.82
LXV525.99AAG551	SC4500 w/o plug Claros LAN+Ethernet IP 2DigSensors	\$	9,674.54
LXV525.99AC1501	SC4500 w/o plug RTC-N/DN(ORP) mAOutput 1DigSensor	\$	29,197.08
LXV525.99AC1541	SC4500 w/o plug RTC-N/DN(ORP) mAOutput 1DigSensor 1mAInput	\$	29,466.89
LXV525.99AC5551	SC4500 w/o plug RTC-N/DN(ORP) Modbus TCP+LAN 2DigSensors	\$	15,051.43
LXV525.99ACB501	SC4500 w/o plug RTC-N/DN(ORP) LAN+mAOutput 1DigSensor	\$	29,409.07
LXV525.99AD1551	SC4500 w/o plug RTC-N/DN(ORP)(DO) mAOutput 2DigSensors	\$	36,501.17
LXV525.99AD2551	SC4500 w/o plug RTC-N/DN(ORP)(DO) ProfibusDP 2DigSensors	\$	38,216.38
LXV525.99AD5541	SC4500 w/o plug RTC-N/DN(ORP)(DO) Modbus TCP+LAN 1DigSensor 1mAInput	\$	37,522.58
LXV525.99AD5551	SC4500 w/o plug RTC-N/DN(ORP)(DO) Modbus TCP+LAN 2DigSensors	\$	37,329.86
LXV525.99ADB551	SC4500 w/o plug RTC-N/DN(ORP)(DO) LAN+mAOutput 2DigSensors	\$	36,713.16
LXV525.99ADC551	SC4500 w/o plug RTC-N/DN(ORP)(DO) LAN+ProfibusDP 2DigSensors	\$	38,216.38
LXV525.99AE1541	SC4500 w/o plug RTC-N/DN mAOutput 1DigSensor 1mAInput	\$	29,466.89
LXV525.99AE1551	SC4500 w/o plug RTC-N/DN mAOutput 2DigSensors	\$	29,274.17
LXV525.99AE5551	SC4500 w/o plug RTC-N/DN Modbus TCP+LAN 2DigSensors	\$	30,102.86
LXV525.99AF1551	SC4500 w/o plug RTC-N/DN(DO) mAOutput 2DigSensors	\$	36,501.17
LXV525.99AF5551	SC4500 w/o plug RTC-N/DN(DO) Modbus TCP+LAN 2DigSensors	\$	18,655.30
LXV525.99AFC541	SC4500 w/o plug RTC-N/DN(DO) LAN+ProfibusDP 1DigSensor 1mAInput	\$	38,409.10
LXV525.99AG1551	SC4500 w/o plug RTC-P mAOutput 2DigSensors	\$	24,937.97
LXV525.99AG5551	SC4500 w/o plug RTC-P Modbus TCP+LAN 2DigSensors	\$	25,785.94
LXV525.99AG7551	SC4500 w/o plug RTC-P Ethernet IP 2DigSensors	\$	26,537.54
LXV525.99AH5551	SC4500 w/o plug RTC-SP Modbus TCP+LAN 2DigSensors	\$	30,102.86
LXV525.99AHB541	SC4500 w/o plug RTC-SP LAN+mAOutput 1DigSensor 1mAInput	\$	29,698.15
LXV525.99AHB551	SC4500 w/o plug RTC-SP LAN+mAOutput 2DigSensors	\$	29,505.43
LXV525.99C00201	SC4500 w/ EU plug No output 1COND	\$	7,072.82
LXV525.99C01111	SC4500 w/ EU plug mAOutput 2pH/ORP	\$	7,419.72
LXV525.99C01501	SC4500 w/ EU plug mAOutput 1DigSensor	\$	7,458.26
LXV525.99C01541	SC4500 EUplug mAOutput 1 dig Sensor 1mA Input	\$	7,573.90
LXV525.99C01551	SC4500 w/ EU plug mAOutput 2DigSensors	\$	7,670.26
LXV525.99C0B441	SC4500 w/ EU plug LAN+mAOutput 2mAInputs	\$	9,308.38
LXV525.99C0B541	SC4500 EUplug LAN+mAOutput 1 dig Sensor 1mA Input	\$	9,404.74
LXV525.99C0B551	SC4500 EU plug Network mA Output 2 digital Sensors	\$	9,520.37
LXV525.99C0C551	SC4500 w/ EU plug LAN+ProfibusDP 2DigSensors	\$	9,520.37
LXV525.99C0F551	SC4500 w/ EU plug LAN+Profinet 2DigSensors	\$	10,426.15
LXV525.99C0G551	SC4500 w/ EU plug LAN+Ethernet IP 2DigSensors	\$	10,426.15
LXV525.99C10101	SC4500 w/ EU plug Prognosys No output 1pH/ORP	\$	7,072.82



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LXV525.99C10201	SC4500 w/ EU plug Prognosys No output 1COND	\$	7,072.82
LXV525.99C10551	SC4500 w/ EU plug Prognosys No output 2DigSensors	\$	6,610.30
LXV525.99C11101	SC4500 w/ EU plug Prognosys mAOutput 1pH/ORP	\$	6,764.47
LXV525.99C11111	SC4500 w/ EU plug Prognosys mAOutput 2pH/ORP	\$	7,419.72
LXV525.99C11401	SC4500 w/ EU plug Prognosys mAOutput 1mAInput	\$	7,342.63
LXV525.99C11441	SC4500 w/ EU plug Prognosys mAOutput 2mAInputs	\$	7,477.54
LXV525.99C11501	SC4500 w/ EU plug Prognosys mAOutput 1DigSensor	\$	7,458.26
LXV525.99C11511	SC4500 w/ EU plug Prognosys mAOutput 1DigSensor 1pH/ORP	\$	7,284.82
LXV525.99C11521	SC4500 w/ EU plug Prognosys mAOutput 1DigSensor 1COND	\$	6,841.56
LXV525.99C11541	SC4500 w/ EU plug Prognosys mAOutput 1DigSensor 1mAInput	\$	7,573.90
LXV525.99C11551	SC4500 EUplug Prognosys mAOutput 2 digital Sensors	\$	7,670.26
LXV525.99C12401	SC4500 w/ EU plug Prognosys ProfibusDP 1mAInput	\$	9,231.29
LXV525.99C12441	SC4500 w/ EU plug Prognosys ProfibusDP 2mAInputs	\$	9,308.38
LXV525.99C12501	SC4500 w/ EU plug Prognosys ProfibusDP 1DigSensor	\$	9,289.10
LXV525.99C12541	SC4500 EUplug Prognosys Profibus 1 dig Sensor 1mA Input	\$	9,404.74
LXV525.99C12551	SC4500 EUplug Prognosys Profibus 2digitalSensors	\$	9,520.37
LXV525.99C15211	SC4500 w/ EU plug Prognosys Modbus TCP+LAN 1COND 1pH/ORP	\$	8,248.42
LXV525.99C15541	SC4500 EUplug Prognosys Modbus 1 dig Sensor 1mA Input	\$	9,404.74
LXV525.99C15551	SC4500 EUplug Prognosys Modbus 2 digital Sensors	\$	9,520.37
LXV525.99C16441	SC4500 EUplug Prognosys Profinet 2mA Input	\$	10,233.43
LXV525.99C16501	SC4500 w/ EU plug Prognosys Profinet 1DigSensor	\$	8,518.22
LXV525.99C16541	SC4500 EUplug Prognosys Profinet 1dig Sensor 1mA Input	\$	10,329.79
LXV525.99C16551	SC4500 EUplug Prognosys Profinet 2digitalSensors	\$	10,426.15
LXV525.99C17441	SC4500 EUplug Prognosys Ethernetá IP 2mAInput	\$	10,233.43
LXV525.99C17541	SC4500 EUplug Prognosys Ethernet IP 1DigSensor 1mAInput	\$	10,329.79
LXV525.99C17551	SC4500 EUplug Prognosys Ethernet IP 2DigSensors	\$	10,426.15
LXV525.99C1A551	SC4500 EUplug Prognosys LAN 2digitalSensors	\$	9,520.37
LXV525.99C1B441	SC4500 w/ EU plug Prognosys LAN+mAOutput 2mAInput	\$	7,188.46
LXV525.99C1B501	SC4500 w/ EU plug Prognosys LAN+mAOutput 1DigSensor	\$	6,687.38
LXV525.99C1B541	SC4500 EU plug Prognosys LAN+mAOutput 1 dig Sensor 1mA Input	\$	9,404.74
LXV525.99C1B551	SC4500 EUplug Prognosys LAN+mAOutput 2digitalSensors	\$	9,520.37
LXV525.99C1C211	SC4500 w/ EU plug Prognosys LAN+ProfibusDP 1COND 1pH/ORP	\$	8,730.22
LXV525.99C1C541	SC4500 w/ EU plug Prognosys LAN+ProfibusDP 1DigSensor 1mAInput	\$	9,404.74
LXV525.99C1C551	SC4500 EU plug Prognosys LAN+Profibus 2 digital Sensors	\$	9,520.37
LXV525.99C1C761	SC4500 w/ EU plug Prognosys LAN+ProfibusDP 1UPW pH/ORP 1UPW COND	\$	9,713.09
LXV525.99C1F441	SC4500 EUplug Prognosys LAN+Profinet 2mA Input	\$	10,233.43
LXV525.99C1F541	SC4500 EUplug Prognosys LAN+Profinet 1dig Sensor 1mA Input	\$	10,329.79
LXV525.99C1F551	SC4500 w/ EU plug Prognosys LAN+Profinet 2DigSensors	\$	10,426.15
LXV525.99C1G441	SC4500 EUplug Prognosys LAN+Ethernet IP 2mAInput	\$	10,233.43
LXV525.99C1G541	SC4500 EUplug Prognosys LAN+Ethernet IP 1DigSensor 1mAInput	\$	10,329.79
LXV525.99C1G551	SC4500 EUplug Prognosys LAN+Ethernet IP 2DigSensors	\$	10,426.15
LXV525.99CA1001	SC4500 w/ EU plug Claros mAOutput	\$	6,725.93
LXV525.99CA1101	SC4500 w/ EU plug Claros mAOutput 1pH/ORP	\$	6,764.47



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LXV525.99CA1201	SC4500 w/ EU plug Claros mAOutput 1COND	\$	7,072.82
LXV525.99CA1401	SC4500 w/ EU plug Claros mAOutput 1mAInput	\$	6,706.66
LXV525.99CA1441	SC4500 w/ EU plug Claros mAOutput 2mAInputs	\$	6,976.46
LXV525.99CA1501	SC4500 w/ EU plug Claros mAOutput 1DigSensor	\$	6,803.02
LXV525.99CA1521	SC4500 w/ EU plug Claros mAOutput 1DigSensor 1COND	\$	6,841.56
LXV525.99CA1541	SC4500 w/ EU plug Claros mAOutput 1DigSensor 1mAInput	\$	6,957.19
LXV525.99CA1551	SC4500 w/ EU plug Claros mAOutput 2DigSensors	\$	7,053.55
LXV525.99CA1601	SC4500 w/ EU plug Claros mAOutput 1UPW COND	\$	7,419.72
LXV525.99CA2521	SC4500 w/ EU plug Claros ProfibusDP 1DigSensor 1COND	\$	8,884.39
LXV525.99CA2541	SC4500 w/ EU plug Claros ProfibusDP 1DigSensor 1mAInput	\$	8,942.21
LXV525.99CA2551	SC4500 w/ EU plug Claros ProfibusDP 2DigSensors	\$	8,884.39
LXV525.99CA5541	SC4500 w/ EU plug Claros Modbus TCP+LAN 1DigSensor 1mAInput	\$	8,826.58
LXV525.99CA5551	SC4500 w/ EU plug Claros Modbus TCP+LAN 2DigSensors	\$	8,884.39
LXV525.99CA6501	SC4500 w/ EU plug Claros Profinet 1DigSensor	\$	8,518.22
LXV525.99CA6541	SC4500 w/ EU plug Claros Profinet 1DigSensor 1mAInput	\$	8,768.76
LXV525.99CA6551	SC4500 w/ EU plug Claros Profinet 2DigSensors	\$	9,578.18
LXV525.99CA7551	SC4500 w/ EU plug Claros Ethernet IP 2DigSensors	\$	10,175.62
LXV525.99CAA551	SC4500 w/ EU plug Claros LAN 2DigSensors	\$	7,920.79
LXV525.99CAB541	SC4500 w/ EU plug Claros LAN+mAOutput 1DigSensor 1mAInput	\$	8,826.58
LXV525.99CAB551	SC4500 EUplug Claros LAN+mAOutput 2digitalSensors	\$	8,884.39
LXV525.99CAB561	SC4500 w/ EU plug Claros LAN+mAOutput 1DigSensor 1UPW COND	\$	7,728.07
LXV525.99CAB571	SC4500 w/ EU plug Claros LAN+mAOutput 1DigSensor 1UPW pH/ORP	\$	7,612.44
LXV525.99CAC211	SC4500 w/ EU plug Claros LAN+ProfibusDP 1COND 1pH/ORP	\$	8,730.22
LXV525.99CAC541	SC4500 w/ EU plug Claros LAN+ProfibusDP 1DigSensor 1mAInput	\$	8,826.58
LXV525.99CAC551	SC4500 w/ EU plug Claros LAN+ProfibusDP 2DigSensors	\$	8,884.39
LXV525.99CAG551	SC4500 w/ EU plug Claros LAN+Ethernet IP 2DigSensors	\$	9,809.45
LXV525.99CC2551	SC4500 w/ EU plug RTC-N/DN(ORP) ProfibusDP 2DigSensors	\$	31,143.55
LXV525.99CCB551	SC4500 w/ EU plug RTC-N/DN(ORP) LAN+mAOutput 2DigSensors	\$	29,640.34
LXV525.99CD5551	SC4500 w/ EU plug RTC-N/DN(ORP)(DO) Modbus TCP+LAN 2DigSensors	\$	37,484.04
LXV525.99CDB551	SC4500 w/ EU plug RTC-N/DN(ORP)(DO) LAN+mAOutput 2DigSensors	\$	36,867.34
LXV525.99CDC541	SC4500 w/ EUplug RTC-N/DN(ORP)(DO) LAN+ProfibusDP 1DigSensor 1mAInput	\$	38,563.27
LXV525.99CEB551	SC4500 w/ EU plug RTC-N/DN LAN+mAOutput 2DigSensors	\$	29,640.34
LXV525.99CEC541	SC4500 w/ EU plug RTC-N/DN LAN+ProfibusDP 1DigSensor 1mAInput	\$	31,336.27
LXV525.99CF5551	SC4500 w/ EU plug RTC-N/DN(DO) Modbus TCP+LAN 2DigSensors	\$	37,484.04
LXV525.99CFB551	SC4500 w/ EU plug RTC-N/DN(DO) LAN+mAOutput 2DigSensors	\$	36,867.34
LXV525.99CGB541	SC4500 w/ EU plug RTC-P LAN+mAOutput 1DigSensor 1mAInput	\$	25,477.58
LXV525.99CGC541	SC4500 w/ EU plug RTC-P LAN+ProfibusDP 1DigSensor 1mAInput	\$	26,980.80
LXV525.99CH5551	SC4500 w/ EU plug RTC-SP Modbus TCP+LAN 2DigSensors	\$	30,257.04
LXV525.99CHB541	SC4500 w/ EU plug RTC-SP LAN+mAOutput 1DigSensor 1mAInput	\$	29,833.06



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LXV525.99CHB551	SC4500 w/ EU plug RTC-SP LAN+mAOutput 2DigSensors	\$	29,640.34
LXV525.99D0B551	SC4500 UKplug Network mAOutput 2digitalSensors	\$	9,520.37
LXV525.99D0B601	SC4500 w/ UK plug LAN+mAOutput 1UPW COND	\$	7,631.71
LXV525.99D0B701	SC4500 w/ UK plug LAN+mAOutput 1UPW pH/ORP	\$	7,516.08
LXV525.99D11501	SC4500 w/ UK plug Prognosys mAOutput 1DigSensor	\$	7,458.26
LXV525.99D11541	SC4500 w/ UK plug Prognosys mAOutput 1DigSensor 1mAInput	\$	7,573.90
LXV525.99D11551	SC4500 w/ UK plug Prognosys mAOutput 2DigSensors	\$	7,670.26
LXV525.99D12541	SC4500 w/ UK plug Prognosys ProfibusDP 1DigSensor 1mAInput	\$	9,404.74
LXV525.99D12551	SC4500 w/ UK plug Prognosys ProfibusDP 2DigSensors	\$	9,520.37
LXV525.99D17551	SC4500 w/ UK plug Prognosys Ethernet IP 2DigSensors	\$	8,132.78
LXV525.99D1B541	SC4500 UK plug Prognosys LAN+mAOutput 1 dig Sensor 1mA Input	\$	9,404.74
LXV525.99D1B551	SC4500 UK plug Prognosys LAN+mAOutput 2 digital Sensors	\$	9,520.37
LXV525.99D1C551	SC4500 UK plug Prognosys LAN+Profibus 2 digital Sensors	\$	9,520.37
LXV525.99DA1551	SC4500 w/ UK plug Claros mAOutput 2DigSensors	\$	6,610.30
LXV525.99DA7551	SC4500 w/ UK plug Claros Ethernet IP 2DigSensors	\$	8,132.78
LXV525.99DAF541	SC4500 w/ UK plug Claros LAN+Profinet 1DigSensor 1mAInput	\$	9,713.09
LXV525.99DAG551	SC4500 w/ UK plug Claros LAN+Ethernet IP 2DigSensors	\$	9,809.45
LXV525.99DCB541	SC4500 w/ UK plug RTC-N/DN(ORP) LAN+mAOutput 1DigSensor 1mAInput	\$	29,833.06
LXV525.99DCB551	SC4500 w/ UK plug RTC-N/DN(ORP) LAN+mAOutput 2DigSensors	\$	29,640.34
LXV525.99DDB541	SC4500 w/ UK plug RTC-N/DN(ORP)(DO) LAN+mAOutput 1DigSensor 1mAInput	\$	37,060.06
LXV525.99DE1501	SC4500 w/ UK plug RTC-N/DN mAOutput 1DigSensor	\$	29,331.98
LXV525.99DEC541	SC4500 w/ UK plug RTC-N/DN LAN+ProfibusDP 1DigSensor 1mAInput	\$	31,336.27
LXV525.99DF1541	SC4500 w/ UK plug RTC-N/DN(DO) mAOutput 1DigSensor 1mAInput	\$	36,848.06
LXV525.99DFB551	SC4500 w/ UK plug RTC-N/DN(DO) LAN+mAOutput 2DigSensors	\$	36,867.34
LXV525.99DG1551	SC4500 w/ UK plug RTC-P mAOutput 2DigSensors	\$	25,072.87
LXV525.99DG5551	SC4500 w/ UK plug RTC-P Modbus TCP+LAN 2DigSensors	\$	25,901.57
LXV525.99DGB541	SC4500 w/ UK plug RTC-P LAN+mAOutput 1DigSensor 1mAInput	\$	25,477.58
LXV525.99E01551	SC4500 w/ US plug mAOutput 2DigSensors	\$	7,670.26
LXV525.99E0B551	SC4500 USplug Network mAOutput 2digitalSensors	\$	9,520.37
LXV525.99E0G551	SC4500 w/ US plug LAN+Ethernet IP 2DigSensors	\$	10,426.15
LXV525.99E11401	SC4500 w/ US plug Prognosys mAOutput 1mAInput	\$	7,342.63
LXV525.99E11441	SC4500 w/ US plug Prognosys mAOutput 2mAInputs	\$	7,477.54
LXV525.99E11501	SC4500 w/ US plug Prognosys mAOutput 1DigSensor	\$	7,458.26
LXV525.99E11541	SC4500 w/ US plug Prognosys mAOutput 1DigSensor 1mAInput	\$	7,573.90
LXV525.99E11551	SC4500 USplug Prognosys mAOutput 2 digital Sensors	\$	7,670.26
LXV525.99E12401	SC4500 w/ US plug Prognosys ProfibusDP 1mAInput	\$	9,231.29
LXV525.99E12441	SC4500 w/ US plug Prognosys ProfibusDP 2mAInputs	\$	9,308.38
LXV525.99E12501	SC4500 w/ US plug Prognosys ProfibusDP 1DigSensor	\$	9,289.10
LXV525.99E12541	SC4500 w/ US plug Prognosys ProfibusDP 1DigSensor 1mAInput	\$	9,404.74
LXV525.99E12551	SC4500 w/ US plug Prognosys ProfibusDP 2DigSensors	\$	9,520.37
LXV525.99E15541	SC4500 USplug Prognosys Modbus 1 dig Sensor 1mA Input	\$	9,404.74
LXV525.99E15551	SC4500 USplug Prognosys Modbus 2 digital Sensors	\$	9,520.37

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LXV525.99E16441	SC4500 w/ US plug Prognosis Profinet 2mAInputs	\$	10,233.43
LXV525.99E16541	SC4500 w/ US plug Prognosis Profinet 1DigSensor 1mAInput	\$	10,329.79
LXV525.99E16551	SC4500 w/ US plug Prognosis Profinet 2DigSensors	\$	10,426.15
LXV525.99E17441	SC4500 w/ US plug Prognosis Ethernet IP 2mAInputs	\$	10,233.43
LXV525.99E17541	SC4500 w/ US plug Prognosis Ethernet IP 1DigSensor 1mAInput	\$	10,329.79
LXV525.99E17551	SC4500 USplug Prognosis Ethernet IP 2 digital Sensors	\$	10,426.15
LXV525.99E1A551	SC4500 USplug Prognosis LAN 2 digital Sensors	\$	9,520.37
LXV525.99E1B541	SC4500 US plug Prognosis LAN+mAOutput 1 dig Sensor 1mA Input	\$	9,404.74
LXV525.99E1B551	SC4500 US plug Prognosis LAN+mAOutput 2 digital Sensors	\$	9,520.37
LXV525.99E1B561	SC4500 w/ US plug Prognosis LAN+mAOutput 1DigSensor 1UPW COND	\$	7,728.07
LXV525.99E1B571	SC4500 w/ US plug Prognosis LAN+mAOutput 1DigSensor 1UPW pH/ORP	\$	7,612.44
LXV525.99E1C551	SC4500 US plug Prognosis LAN+Profibus 2 digital Sensors	\$	9,520.37
LXV525.99E1F441	SC4500 w/ US plug Prognosis LAN+Profinet 2mAInputs	\$	10,233.43
LXV525.99E1F541	SC4500 w/ US plug Prognosis LAN+Profinet 1DigSensor 1mAInput	\$	10,329.79
LXV525.99E1F551	SC4500 w/ US plug Prognosis LAN+Profinet 2DigSensors	\$	10,426.15
LXV525.99E1G441	SC4500 w/ US plug Prognosis LAN+Ethernet IP 2mAInputs	\$	10,233.43
LXV525.99E1G541	SC4500 w/ US plug Prognosis LAN+Ethernet IP 1DigSensor 1mAInput	\$	10,329.79
LXV525.99E1G551	SC4500 w/ US plug Prognosis LAN+Ethernet IP 2DigSensors	\$	10,426.15
LXV525.99EA1401	SC4500 w/ US plug Claros mAOutput 1mAInput	\$	6,725.93
LXV525.99EA1441	SC4500 w/ US plug Claros mAOutput 2mAInputs	\$	6,841.56
LXV525.99EA1541	SC4500 w/ US plug Claros mAOutput 1DigSensor 1mAInput	\$	6,957.19
LXV525.99EA1551	SC4500 w/ US plug Claros mAOutput 2DigSensors	\$	7,053.55
LXV525.99EA6441	SC4500 w/ US plug Claros Profinet 2mAInputs	\$	9,346.92
LXV525.99EA6541	SC4500 w/ US plug Claros Profinet 1DigSensor 1mAInput	\$	9,289.10
LXV525.99EA6551	SC4500 w/ US plug Claros Profinet 2DigSensors	\$	9,578.18
LXV525.99EA7441	SC4500 w/ US plug Claros Ethernet IP 2mAInputs	\$	9,346.92
LXV525.99EA7541	SC4500 w/ US plug Claros Ethernet IP 1DigSensor 1mAInput	\$	9,289.10
LXV525.99EA7551	SC4500 w/ US plug Claros Ethernet IP 2DigSensors	\$	9,578.18
LXV525.99EAB541	SC4500 USplug Claros LAN+mAOutput 1dig Sensor 1mA Input	\$	8,826.58
LXV525.99EAB551	SC4500 USplug Claros LAN+mAOutput 2digital Sensors	\$	8,884.39
LXV525.99EAF441	SC4500 w/ US plug Claros LAN+Profinet 2mAInputs	\$	9,346.92
LXV525.99EAF541	SC4500 w/ US plug Claros LAN+Profinet 1DigSensor 1mAInput	\$	9,289.10
LXV525.99EAF551	SC4500 w/ US plug Claros LAN+Profinet 2DigSensors	\$	9,578.18
LXV525.99EAG441	SC4500 w/ US plug Claros LAN+Ethernet IP 2mAInputs	\$	9,346.92
LXV525.99EAG541	SC4500 w/ US plug Claros LAN+Ethernet IP 1DigSensor 1mAInput	\$	9,289.10
LXV525.99EAG551	SC4500 w/ US plug Claros LAN+Ethernet IP 2DigSensors	\$	9,809.45
LXV525.99ECB551	SC4500 w/ US plug RTC-N/DN(ORP) LAN+mAOutput 2DigSensors	\$	29,640.34
LXV525.99EDB541	SC4500 w/ US plug RTC-N/DN(ORP)(DO) LAN+mAOutput 1DigSensor 1mAInput	\$	37,060.06
LXV525.99EDB551	SC4500 w/ US plug RTC-N/DN(ORP)(DO) LAN+mAOutput 2DigSensors	\$	36,867.34
LXV525.99EE1501	SC4500 w/ US plug RTC-N/DN mAOutput 1DigSensor	\$	29,331.98



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LXV525.99EEB551	SC4500 w/ US plug RTC-N/DN LAN+mAOutput 2DigSensors	\$	29,640.34
LXV525.99EF1541	SC4500 w/ US plug RTC-N/DN(DO) mAOutput 1DigSensor 1mAInput	\$	36,848.06
LXV525.99EFB541	SC4500 w/ US plug RTC-N/DN(DO) LAN+mAOutput 1DigSensor 1mAInput	\$	37,060.06
LXV525.99EG1501	SC4500 w/ US plug RTC-P mAOutput 1DigSensor	\$	24,995.78
LXV525.99EGB551	SC4500 w/ US plug RTC-P LAN+mAOutput 2DigSensors	\$	25,284.86
LXV525.99EH5551	SC4500 w/ US plug RTC-SP Modbus TCP+LAN 2DigSensors	\$	30,257.04
LXV525.99EHB541	SC4500 w/ US plug RTC-SP LAN+mAOutput 1DigSensor 1mAInput	\$	29,833.06
LXV525.99EHB551	SC4500 w/ US plug RTC-SP LAN+mAOutput 2DigSensors	\$	29,640.34
LXV525.99F11501	SC4500 w/ Switzerland Prognosis mAOutput 1DigSensor	\$	7,458.26
LXV525.99F11551	SC4500 w/ Switzerland Prognosis mAOutput 2DigSensors	\$	7,670.26
LXV525.99F12501	SC4500 w/ Switzerland Prognosis ProfibusDP 1DigSensor	\$	9,289.10
LXV525.99F12551	SC4500 w/ Switzerland Prognosis ProfibusDP 2DigSensors	\$	9,520.37
LXV525.99F15541	SC4500 w/ Switzerland Prognosis Modbus TCP+LAN 1DigSensor 1mAInput	\$	9,404.74
LXV525.99F15551	SC4500 w/ Switzerland Prognosis Modbus TCP+LAN 2DigSensors	\$	9,520.37
LXV525.99F1A551	SC4500 CHplug Prognosis LAN 2digitalSensors	\$	9,520.37
LXV525.99F1B541	SC4500 CHplug Prognosis LAN+mAOutput 1 dig Sensor 1mA Input	\$	9,404.74
LXV525.99F1B551	SC4500 CHplug Prognosis LAN+mAOutput 2digitalSensors	\$	9,520.37
LXV525.99F1B601	SC4500 w/ SwissPrognosis LAN+mAOutput 1UPW COND	\$	7,631.71
LXV525.99F1B701	SC4500 w/ Switzerland Prognosis LAN+mAOutput 1UPW pH/ORP	\$	7,516.08
LXV525.99F1C541	SC4500 w/ Switzerland Prognosis LAN+ProfibusDP 1DigSensor 1mAInput	\$	9,404.74
LXV525.99F1C551	SC4500 w/ Switzerland Prognosis LAN+ProfibusDP 2DigSensors	\$	9,520.37
LXV525.99FA5541	SC4500 w/ Switzerland Claros Modbus TCP+LAN 1DigSensor 1mAInput	\$	8,826.58
LXV525.99FA5551	SC4500 w/ Switzerland Claros Modbus TCP+LAN 2DigSensors	\$	8,884.39
LXV525.99FAB541	SC4500 w/ Switzerland Claros LAN+mAOutput 1DigSensor 1mAInput	\$	8,826.58
LXV525.99FAB551	SC4500 w/ Switzerland Claros LAN+mAOutput 2DigSensors	\$	8,884.39
LXV525.99FAC541	SC4500 w/ Switzerland Claros LAN+ProfibusDP 1DigSensor 1mAInput	\$	8,826.58
LXV525.99FAC551	SC4500 w/ Switzerland Claros LAN+ProfibusDP 2DigSensors	\$	8,884.39
LXV525.99FC2551	SC4500 w/ Swiss RTC-N/DN(ORP) ProfibusDP 2DigSensors	\$	31,143.55
LXV525.99FDB541	SC4500 w/ Swiss RTC-N/DN(ORP)(DO) LAN+mAOutput 1DigSensor 1mAInput	\$	37,060.06
LXV525.99FDC541	SC4500 w/ Swiss RTC-N/DN(ORP)(DO) LAN+ProfibusDP 1DigSensor 1mAInput	\$	38,563.27
LXV525.99KTO048	KTO: SC4500 EU WiFi Bundle, mA Output	\$	7,689.53
LXV525.99KTO049	KTO: SC4500 US WiFi Bundle, mA Output	\$	8,229.14
LXV525.99KTO050	KTO: SC4500 RoW WiFi Bundle, mA Output	\$	8,229.14
LXV525.99KTO051	KTO: SC4500 EU GSM Bundle, mA Output	\$	9,250.56
LXV525.99KTO052	KTO: SC4500 US GSM Bundle, mA Output	\$	9,250.56
LXV525.99P11201	SC4500 w/o plug C1D2 Prognosis mAOutput 1COND	\$	7,361.90

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LXV525.99P11401	SC4500 w/o plug C1D2 Prognosys mAOutput 1mAInput	\$	7,670.26
LXV525.99P11441	SC4500 w/o plug C1D2 Prognosys mAOutput 2mAInputs	\$	7,785.89
LXV525.99P11501	SC4500 w/o plug C1D2 Prognosys mAOutput 1DigSensor	\$	7,766.62
LXV525.99P11541	SC4500 w/o plug C1D2 Prognosys mAOutput 1DigSensor 1mAInput	\$	7,882.25
LXV525.99P11551	SC4500 w/o plug C1D2 Prognosys mAOutput 2DigSensors	\$	7,959.34
LXV525.99P12401	SC4500 w/o plug C1D2 Prognosys ProfibusDP 1mAInput	\$	9,520.37
LXV525.99P12441	SC4500 w/o plug C1D2 Prognosys ProfibusDP 2mAInputs	\$	9,616.73
LXV525.99P12501	SC4500 w/o plug C1D2 Prognosys ProfibusDP 1DigSensor	\$	9,578.18
LXV525.99P12541	SC4500 w/o plug C1D2 Prognosys ProfibusDP 1DigSensor 1mAInput	\$	9,713.09
LXV525.99P12551	SC4500 w/o plug C1D2 Prognosys ProfibusDP 2DigSensors	\$	9,809.45
LXV525.99P1B401	SC4500 w/o plug C1D2 Prognosys LAN+mAOutput 1mAInput	\$	9,520.37
LXV525.99P1B441	SC4500 w/o plug C1D2 Prognosys LAN+mAOutput 2mAInputs	\$	9,616.73
LXV525.99P1B501	SC4500 w/o plug C1D2 Prognosys LAN+mAOutput 1DigSensor	\$	9,578.18
LXV525.99P1B541	SC4500 w/o plug C1D2 Prognosys LAN+mAOutput 1DigSensor 1mAInput	\$	9,713.09
LXV525.99P1B551	SC4500 w/o plug C1D2 Prognosys LAN+mAOutput 2DigSensors	\$	9,809.45
LXV525.99P1C761	SC4500 w/o plug C1D2 Prognosys LAN+ProfibusDP 1UPW pH/ORP 1UPW COND	\$	10,002.17
LXV525.99PAB211	SC4500 w/o plug C1D2 Claros LAN+mAOutput 1COND 1pH/ORP	\$	9,424.01
LXV525.99PAB401	SC4500 w/o plug C1D2 Claros LAN+mAOutput 1mAInput	\$	8,884.39
LXV525.99PAB441	SC4500 w/o plug C1D2 Claros LAN+mAOutput 2mAInputs	\$	9,038.57
LXV525.99PAB501	SC4500 w/o plug C1D2 Claros LAN+mAOutput 1DigSensor	\$	8,980.75
LXV525.99PAB541	SC4500 w/o plug C1D2 Claros LAN+mAOutput 1DigSensor 1mAInput	\$	9,096.38
LXV525.99PAB551	SC4500 w/o plug C1D2 Claros LAN+mAOutput 2DigSensors	\$	9,231.29
LXV525.99PAC541	SC4500 w/o plug C1D2 Claros LAN+ProfibusDP 1DigSensor 1mAInput	\$	9,096.38
LXV525.99Y11401	SC4500 w/o plug 24VDC C1D2 Prognosys mAOutput 1mAInput	\$	7,805.16
LXV525.99Y11441	SC4500 w/o plug 24VDC C1D2 Prognosys mAOutput 2mAInputs	\$	7,920.79
LXV525.99Y11501	SC4500 w/o plug 24VDC C1D2 Prognosys mAOutput 1DigSensor	\$	7,901.52
LXV525.99Y11541	SC4500 w/o plug 24VDC C1D2 Prognosys mAOutput 1DigSensor 1mAInput	\$	7,997.88
LXV525.99Y11551	SC4500 w/o plug 24VDC C1D2 Prognosys mAOutput 2DigSensors	\$	8,113.51
LXV525.99Y12401	SC4500 w/o plug 24VDC C1D2 Prognosys ProfibusDP 1mAInput	\$	9,674.54
LXV525.99Y12441	SC4500 w/o plug 24VDC C1D2 Prognosys ProfibusDP 2mAInputs	\$	9,751.63
LXV525.99Y12501	SC4500 w/o plug 24VDC C1D2 Prognosys ProfibusDP 1DigSensor	\$	9,732.36
LXV525.99Y12541	SC4500 w/o plug 24VDC C1D2 Prognosys ProfibusDP 1DigSensor 1mAInput	\$	9,886.54
LXV525.99Y12551	SC4500 w/o plug 24VDC C1D2 Prognosys ProfibusDP 2DigSensors	\$	9,944.35
LXV525.99Y1B401	SC4500 w/o plug 24VDC C1D2 Prognosys LAN+mAOutput 1mAInput	\$	9,674.54
LXV525.99Y1B441	SC4500 w/o plug 24VDC C1D2 Prognosys LAN+mAOutput 2mAInputs	\$	9,751.63



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LXV525.99Y1B501	SC4500 w/o plug 24VDC C1D2 Prognosys LAN+mAOutput 1DigSensor	\$	9,732.36
LXV525.99Y1B541	SC4500 w/o plug 24VDC C1D2 Prognosys LAN+mAOutput 1DigSensor 1mAInput	\$	9,886.54
LXV525.99Y1B551	SC4500 w/o plug 24VDC C1D2 Prognosys LAN+mAOutput 2DigSensors	\$	9,944.35
LXV525.99Y1C761	SC4500 w/o plug 24V C1D2 Prognosys LAN+ProfibusDP UPW pH/ORP UPW COND	\$	10,002.17
LXV525.99YAB401	SC4500 w/o plug 24VDC C1D2 Claros LAN+mAOutput 1mAInput	\$	9,057.84
LXV525.99YAB441	SC4500 w/o plug 24VDC C1D2 Claros LAN+mAOutput 2mAInputs	\$	9,134.93
LXV525.99YAB501	SC4500 w/o plug 24VDC C1D2 Claros LAN+mAOutput 1DigSensor	\$	9,115.66
LXV525.99YAB541	SC4500 w/o plug 24VDC C1D2 Claros LAN+mAOutput 1DigSensor 1mAInput	\$	9,269.83
LXV525.99YAB551	SC4500 w/o plug 24VDC C1D2 Claros LAN+mAOutput 2DigSensors	\$	9,327.65
LXV525.99Z01501	SC4500 w/o plug 24VDC mAOutput 1DigSensor	\$	7,766.62
LXV525.99Z01541	SC4500 w/oPlug 24VDC mAOutput 1dig Sensor 1mA Input	\$	7,882.25
LXV525.99Z01551	SC4500 w/oPlug 24VDC mAOutput 2dig Sensor	\$	7,959.34
LXV525.99Z01601	SC4500 w/o plug 24VDC mAOutput 1UPW COND	\$	7,708.80
LXV525.99Z01661	SC4500 w/o plug 24VDC mAOutput 2UPW COND	\$	8,402.59
LXV525.99Z01701	SC4500 w/o plug 24VDC mAOutput 1UPW pH/ORP	\$	7,593.17
LXV525.99Z01771	SC4500 w/o plug 24VDC mAOutput 2UPW pH/ORP	\$	8,171.33
LXV525.99Z05541	SC4500 w/o plug 24VDC Modbus TCP+LAN 1DigSensor 1mAInput	\$	9,713.09
LXV525.99Z05551	SC4500 w/o plug 24VDC Modbus TCP+LAN 2DigSensors	\$	9,809.45
LXV525.99Z07441	SC4500 w/o plug 24VDC Ethernet IP 2mAInputs	\$	10,541.78
LXV525.99Z07541	SC4500 w/o plug 24VDC Ethernet IP 1DigSensor 1mAInput	\$	10,638.14
LXV525.99Z07551	SC4500 w/o plug 24VDC Ethernet IP 2DigSensors	\$	10,734.50
LXV525.99Z0C541	SC4500 w/o plug 24VDC LAN+ProfibusDP 1DigSensor 1mAInput	\$	9,713.09
LXV525.99Z0G441	SC4500 w/o plug 24VDC LAN+Ethernet IP 2mAInputs	\$	10,541.78
LXV525.99Z0G551	SC4500 w/o plug 24VDC LAN+Ethernet IP 2DigSensors	\$	10,734.50
LXV525.99Z11101	SC4500 w/o plug 24VDC Prognosys mAOutput 1pH/ORP	\$	7,034.28
LXV525.99Z11111	SC4500 w/o plug 24VDC Prognosys mAOutput 2pH/ORP	\$	7,708.80
LXV525.99Z11201	SC4500 w/o plug 24VDC Prognosys mAOutput 1COND	\$	7,361.90
LXV525.99Z11211	SC4500 w/o plug 24VDC Prognosys mAOutput 1COND 1pH/ORP	\$	7,708.80
LXV525.99Z11221	SC4500 w/o plug 24VDC Prognosys mAOutput 2Conductivity	\$	7,708.80
LXV525.99Z11401	SC4500 w/o plug 24VDC Prognosys mAOutput 1mAInput	\$	7,670.26
LXV525.99Z11441	SC4500 w/o plug 24VDC Prognosys mAOutput 2mAInputs	\$	7,785.89
LXV525.99Z11501	SC4500 w/o plug 24VDC Prognosys mAOutput 1DigSensor	\$	7,766.62
LXV525.99Z11511	SC4500 w/o plug 24VDC Prognosys mAOutput 1DigSensor 1pH/ORP	\$	6,051.41
LXV525.99Z11521	SC4500 w/o plug 24VDC Prognosys mAOutput 1DigSensor 1COND	\$	7,130.64
LXV525.99Z11541	SC4500 w/o plug 24VDC Prognosys mAOutput 1DigSensor 1mAInput	\$	7,882.25
LXV525.99Z11551	SC4500 w/o plug 24VDC Prognosys mAOutput 2DigSensors	\$	7,959.34
LXV525.99Z11601	SC4500 w/o plug 24VDC Prognosys mAOutput 1UPW COND	\$	7,708.80
LXV525.99Z11661	SC4500 w/o plug 24VDC Prognosys mAOutput 2UPW COND	\$	8,402.59



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LXV525.99Z11701	SC4500 w/o plug 24VDC Prognosys mAOutput 1UPW pH/ORP	\$	7,593.17
LXV525.99Z11771	SC4500 w/o plug 24VDC Prognosys mAOutput 2UPW pH/ORP	\$	8,171.33
LXV525.99Z12101	SC4500 w/o plug 24VDC Prognosys ProfibusDP 1pH/ORP	\$	9,077.11
LXV525.99Z12201	SC4500 w/o plug 24VDC Prognosys ProfibusDP 1COND	\$	9,077.11
LXV525.99Z12401	SC4500 w/o plug 24VDC Prognosys ProfibusDP 1mAInput	\$	9,520.37
LXV525.99Z12441	SC4500 w/o plug 24VDC Prognosys ProfibusDP 2mAInputs	\$	9,616.73
LXV525.99Z12501	SC4500 w/o plug 24VDC Prognosys ProfibusDP 1DigSensor	\$	9,578.18
LXV525.99Z12541	SC4500 w/o plug 24VDC Prognosys ProfibusDP 1DigSensor 1mAInput	\$	9,713.09
LXV525.99Z12551	SC4500 w/o plug 24VDC Prognosys ProfibusDP 2DigSensors	\$	9,809.45
LXV525.99Z15501	SC4500 w/o plug 24VDC Prognosys Modbus TCP+LAN 1DigSensor	\$	7,593.17
LXV525.99Z15541	SC4500 w/o plug 24VDC Prognosys Modbus TCP+LAN 1DigSensor 1mAInput	\$	9,713.09
LXV525.99Z15551	SC4500 w/o plug 24VDC Prognosys Modbus TCP+LAN 2DigSensors	\$	9,809.45
LXV525.99Z15601	SC4500 w/o plug 24VDC Prognosys Modbus TCP+LAN 1UPW COND	\$	8,537.50
LXV525.99Z15661	SC4500 w/o plug 24VDC Prognosys Modbus TCP+LAN 2UPW COND	\$	9,231.29
LXV525.99Z15701	SC4500 w/o plug 24VDC Prognosys Modbus TCP+LAN 1UPW pH/ORP	\$	8,421.86
LXV525.99Z15771	SC4500 w/o plug 24VDC Prognosys Modbus TCP+LAN 2UPW pH/ORP	\$	9,000.02
LXV525.99Z16441	SC4500 w/o plug 24VDC Prognosys Profinet 2mAInputs	\$	10,541.78
LXV525.99Z16541	SC4500 w/o plug 24VDC Prognosys Profinet 1DigSensor 1mAInput	\$	10,638.14
LXV525.99Z16551	SC4500 w/o plug 24VDC Prognosys Profinet 2DigSensors	\$	10,734.50
LXV525.99Z17441	SC4500 w/o plug 24VDC Prognosys Ethernet IP 2mAInputs	\$	10,541.78
LXV525.99Z17541	SC4500 w/o plug 24VDC Prognosys Ethernet IP 1DigSensor 1mAInput	\$	10,638.14
LXV525.99Z17551	SC4500 w/o plug 24VDC Prognosys Ethernet IP 2DigSensors	\$	10,734.50
LXV525.99Z1A501	SC4500 w/o plug 24VDC Prognosys LAN 1DigSensor	\$	6,918.65
LXV525.99Z1B501	SC4500 w/o plug 24VDC Prognosys LAN+mAOutput 2DigSensors	\$	8,980.75
LXV525.99Z1B551	SC4500 w/o Plug 24VDC Prognosys LAN+mAOutput 2dig Sensor	\$	9,809.45
LXV525.99Z1C501	SC4500 w/o plug 24VDC Prognosys LAN+ProfibusDP 1DigSensor	\$	8,441.14
LXV525.99Z1C541	SC4500 w/o plug 24VDC Prognosys LAN+ProfibusDP 1DigSensor 1mAInput	\$	9,713.09
LXV525.99Z1F441	SC4500 w/o Plug 24VDC Prognosys LAN+Profinet 2mA Input	\$	10,541.78
LXV525.99Z1F541	SC4500 w/o Plug 24VDC Prognosys LAN+Profinet 1dig Sensor 1mA Input	\$	10,638.14
LXV525.99Z1F551	SC4500 w/o Plug 24VDC Prognosys LAN+Profinet 2dig Sensor	\$	10,734.50
LXV525.99Z1G441	SC4500 w/o plug 24VDC Prognosys LAN+Ethernet IP 2mAInputs	\$	10,541.78
LXV525.99Z1G541	SC4500 w/o plug 24VDC Prognosys LAN+Ethernet IP 1DigSensor 1mAInput	\$	10,638.14
LXV525.99Z1G551	SC4500 w/o plug 24VDC Prognosys LAN+Ethernet IP 2DigSensors	\$	10,734.50
LXV525.99ZA1101	SC4500 w/o plug 24VDC Claros mAOutput 1pH/ORP	\$	7,034.28

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LXV525.99ZA1111	SC4500 w/o plug 24VDC Claros mAOutput 2pH/ORP	\$	7,708.80
LXV525.99ZA1201	SC4500 w/o plug 24VDC Claros mAOutput 1COND	\$	7,361.90
LXV525.99ZA1211	SC4500 w/o plug 24VDC Claros mAOutput 1COND 1pH/ORP	\$	7,708.80
LXV525.99ZA1221	SC4500 w/o plug 24VDC Claros mAOutput 2Conductivity	\$	7,708.80
LXV525.99ZA1501	SC4500 w/o plug 24VDC Claros mAOutput 1DigSensor	\$	7,130.64
LXV525.99ZA1511	SC4500 w/o plug 24VDC Claros mAOutput 1DigSensor 1pH/ORP	\$	6,051.41
LXV525.99ZA1521	SC4500 w/o plug 24VDC Claros mAOutput 1DigSensor 1COND	\$	7,130.64
LXV525.99ZA1541	SC4500 w/o plug 24VDC Claros mAOutput 1DigSensor 1mAInput	\$	7,882.25
LXV525.99ZA1551	SC4500 w/o plug 24VDC Claros mAOutput 2DigSensors	\$	7,342.63
LXV525.99ZA2501	SC4500 w/o plug 24VDC Claros ProfibusDP 1DigSensor	\$	9,366.19
LXV525.99ZA2551	SC4500 w/o plug 24VDC Claros ProfibusDP 2DigSensors	\$	9,578.18
LXV525.99ZA5101	SC4500 w/o plug 24VDC Claros Modbus TCP+LAN 1pH/ORP	\$	8,190.60
LXV525.99ZA5541	SC4500 w/o plug 24VDC Claros Modbus TCP+LAN 1DigSensor 1mAInput	\$	9,096.38
LXV525.99ZA5551	SC4500 w/o plug 24VDC Claros Modbus TCP+LAN 2DigSensors	\$	9,231.29
LXV525.99ZA6501	SC4500 w/o plug 24VDC Claros ProfibusDP 1DigSensor	\$	10,252.70
LXV525.99ZA6551	SC4500 w/o plug 24VDC Claros Profinet 2DigSensors	\$	10,464.70
LXV525.99ZA7551	SC4500 w/o plug 24VDC Claros Ethernet IP 2DigSensors	\$	8,421.86
LXV525.99ZAA551	SC4500 w/o plug 24VDC Claros LAN 2DigSensors	\$	8,633.86
LXV525.99ZAB541	SC4500 w/oPlug 24VDC Claros LAN+mAOutput 1dig Sensor 1mA Input	\$	9,096.38
LXV525.99ZAB551	SC4500 w/oPlug 24VDC Claros LAN+mAOutput 2dig Sensor	\$	9,231.29
LXV525.99ZAB701	SC4500 w/o plug 24VDC Claros LAN+mAOutput 1UPW pH/ORP	\$	7,805.16
LXV525.99ZAB771	SC4500 w/o plug 24VDC Claros LAN+mAOutput 2UPW pH/ORP	\$	8,383.32
LXV525.99ZAC541	SC4500 w/o plug 24VDC Claros LAN+ProfibusDP 1DigSensor 1mAInput	\$	9,096.38
LXV525.99ZAC551	SC4500 w/o plug 24VDC Claros LAN+ProfibusDP 2DigSensors	\$	9,578.18
LXV525.99ZD1541	SC4500 w/o plug 24VDC RTC-N/DN(ORP)(DO) mAOutput 1DigSensor 1mAInput	\$	37,117.87
LXV525.99ZD5541	SC4500 w/o plug 24VDC RTC-N/DN(ORP)(DO) Modbus+LAN 1DigS. 1mAInput	\$	37,946.57
LXV525.99ZD5551	SC4500 w/o plug 24VDC RTC-N/DN(ORP)(DO) Modbus TCP+LAN 2DigSensors	\$	37,753.85
LXV525.99ZDB551	SC4500 w/o plug 24VDC RTC-N/DN(ORP)(DO) LAN+mAOutput 2DigSensors	\$	37,137.14
LXV525.LMKTO020	KTO: CLF10sc, sc4500 DUAL INPUT, pHD	\$	19,349.09
LXV526.97.1011A	NA5600sc SODIUM, WALL, 1 CH, NO AUTOCALIBRATION, US	\$	48,893.06
LXV526.97.1011G	NA9600sc SODIUM, WALL, 1 CH, CN	\$	47,351.30
LXV526.97.1012A	NA5600sc SODIUM, WALL, 2 CH, NO AUTOCALIBRATION, US	\$	54,481.94
LXV526.97.1012G	NA9600sc SODIUM, WALL, 2 CH, CN	\$	52,824.55
LXV526.97.1014A	NA5600sc SODIUM, WALL, 4 CH, NO AUTOCALIBRATION, US	\$	59,608.30
LXV526.97.1014G	NA9600sc SODIUM, WALL, 4 CH, CN	\$	57,835.27
LXV526.97.1111A	NA5600sc SODIUM, WALL, AUTOCAL., 1 CH, US	\$	46,021.54
LXV526.97.1111G	NA9600sc SODIUM, WALL, AUTOCAL., 1 CH, CN	\$	48,449.81
LXV526.97.1112A	NA5600sc SODIUM, WALL, AUTOCAL., 2 CH, US	\$	51,918.77
LXV526.97.1112G	NA9600sc SODIUM, WALL, AUTOCAL., 2 CH, CN	\$	54,616.85
LXV526.97.1114A	NA5600sc SODIUM, WALL, AUTOCAL., 4 CH, US	\$	63,771.05



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LXV526.97.1114G	NA9600sc SODIUM, WALL, AUTOCAL., 4 CH, CN	\$	67,124.38
LXV526.97.1211A	NA5600sc SODIUM, WALL, K KIT, 1 CH, US	\$	50,280.65
LXV526.97.1212A	NA5600sc SODIUM, WALL, K KIT, 2 CH, US	\$	56,120.06
LXV526.97.1214A	NA5600sc SODIUM, WALL, K KIT, 4 CH, US	\$	61,419.86
LXV526.97.1311A	NA5600sc SODIUM, WALL, K KIT & AUTOCAL., 1 CH, US	\$	49,779.58
LXV526.97.1312A	NA5600sc SODIUM, WALL, K KIT & AUTOCAL., 2 CH, US	\$	55,618.99
LXV526.97.1314A	NA5600sc SODIUM, WALL, K KIT & AUTOCAL., 4 CH, US	\$	67,490.54
LXV526.97.2011A	NA5600sc SODIUM, PANEL, 1 CH, NO AUTOCALIBRATION, US	\$	46,638.24
LXV526.97.2011G	NA9600sc SODIUM, PANEL, 1 CH, CN	\$	45,289.20
LXV526.97.2012A	NA5600sc SODIUM, PANEL, 2 CH, NO AUTOCALIBRATION, US	\$	52,227.12
LXV526.97.2012G	NA9600sc SODIUM, PANEL, 2 CH, CN	\$	50,685.36
LXV526.97.2014A	NA5600sc SODIUM, PANEL, 4 CH, NO AUTOCALIBRATION, US	\$	56,929.49
LXV526.97.2014G	NA9600sc SODIUM, PANEL, 4 CH, CN	\$	55,175.74
LXV526.97.2111A	NA5600sc SODIUM, PANEL, AUTOCAL., 1 CH, US	\$	45,192.84
LXV526.97.2111G	NA9600sc SODIUM, PANEL, AUTOCAL, 1 CH, CN	\$	47,563.30
LXV526.97.2112A	NA5600sc SODIUM, PANEL, AUTOCAL., 2 CH, US	\$	51,128.62
LXV526.97.2112G	NA9600sc SODIUM, PANEL, AUTOCAL, 2 CH, CN	\$	53,807.42
LXV526.97.2114A	NA5600sc SODIUM, PANEL, AUTOCAL., 4 CH, US	\$	62,961.62
LXV526.97.2114G	NA9600sc SODIUM, PANEL, AUTOCAL, 4 CH, CN	\$	66,295.68
LXV526.97.2211A	NA5600sc SODIUM, PANEL, K KIT, 1 CH, US	\$	48,835.25
LXV526.97.2212A	NA5600sc SODIUM, PANEL, K KIT, 2 CH, US	\$	54,674.66
LXV526.97.2214A	NA5600sc SODIUM, PANEL, K KIT, 4 CH, US	\$	59,589.02
LXV526.97.2311A	NA5600sc SODIUM, PANEL, K KIT & AUTOCAL., 1 CH, US	\$	49,182.14
LXV526.97.2312A	NA5600sc SODIUM, PANEL, K KIT & AUTOCAL., 2 CH, US	\$	55,117.92
LXV526.97.2314A	NA5600sc SODIUM, PANEL, K KIT & AUTOCAL., 4 CH, US	\$	67,008.74
LXV541.99.0001H	ee PAH500, DIGITAL SENSOR	\$	52,072.94
LXV541.99.1001H	KTO: PAH500, DIGITAL SENSOR, 1.0M CABLE ASSY, BRACKET	\$	53,653.25
LXV541.99.2001H	KTO: PAH500, DIGITAL SENSOR, 7.7M CABLE ASSY, BRACKET	\$	53,865.24
LXV541.99.3001H	KTO: PAH500, DIGITAL SENSOR, 15M CABLE ASSY, BRACKET	\$	54,135.05
LXV541.99.4001H	KTO: PAH500, DIGITAL SENSOR, 31M CABLE ASSY, BRACKET	\$	55,079.38
LXZ400.99.00001	CONNECTOR POWER CONNECTION CARD w/ FERRITE; SC1000	\$	69.19
LXZ446.99.00001	GATEWAY IIoT RevPi; Basic Module; DIN rail housing	\$	1,108.14
LXZ446.99.00002	GATEWAY EtherNet/IP Slave; DIN rail housing	\$	553.11
LXZ446.99.00003	JUMPER FOR PiBridge	\$	25.25
LXZ446.99.00005	NA 4G modem Kit, SC1500	\$	1,626.56
LXZ446.99.00007	GATEWAY Profinet IRT Slave; DIN rail housing	\$	553.11
LXZ448.99.00001	USB-RS485 Connection cable with male connector 5 pin	\$	112.93
LXZ448.99.00002	Set 1mm Wiper NT3x00sc	\$	481.80
LXZ448.99.00003	Set 2mm Wiper NT3x00sc	\$	302.57
LXZ448.99.00004	EMC Cushion, self adhesive NT3x00sc	\$	31.41
LXZ448.99.00005	Cable w strain relief NT3x00sc	\$	387.37
LXZ448.99.00006	Set optical Filter, NT3200sc	\$	3,642.41
LXZ448.99.00007	Wiper Arm 5mm incl Screw NT3x00sc	\$	47.22
LXZ448.99.00008	Grub screw, 5mm Wiper NT3x00sc	\$	29.68
LXZ448.99.00009	Set enclosure screws, NT3x00sc	\$	40.09
LXZ448.99.00010	Set of screws, internal, NT3x00sc	\$	51.26
LXZ448.99.00011	Sensor cap NT3x00sc	\$	57.43



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LXZ448.99.00012	Set of Orings NT3x00sc	\$	42.40
LXZ448.99.00013	UV Cover, NT3x00sc	\$	34.69
LXZ448.99.00014	Position Detection, cplt, NT3x00sc	\$	144.73
LXZ448.99.00015	Housing, Lamp Side, NT3x00sc	\$	653.32
LXZ448.99.00016	Housing, Receiver Side, NT3x00sc	\$	686.08
LXZ448.99.00017	Chassis, NT3x00sc	\$	48.37
LXZ448.99.00018	Flash Connector, Pair, NT3x00sc	\$	139.34
LXZ448.99.00019	Grommet, incl O rings, NT3x00sc	\$	84.41
LXZ448.99.00020	EMC steel, NT3x00sc	\$	51.26
LXZ448.99.00021	Head, Lamp Side, 1mm, NT3x00sc	\$	809.42
LXZ448.99.00022	Head, Lamp Side, 2mm, NT3x00sc	\$	809.42
LXZ448.99.00023	Head, Lamp Side, 5mm, NT3x00sc	\$	809.42
LXZ448.99.00024	Lamp Sleeve w Heat Board NT3x00sc	\$	181.35
LXZ448.99.00025	Rg + Spring Lamp Seat, NT3x00sc	\$	120.06
LXZ448.99.00026	Lamp Holder incl Lense, NT3x00sc	\$	628.27
LXZ448.99.00027	Shaft, cplt, NT3x00sc	\$	152.06
LXZ448.99.00028	Shaft Stub, cplt, NT3x00sc	\$	104.45
LXZ448.99.00029	Beam Path NT3200sc	\$	503.00
LXZ448.99.00030	Beam Path NT3100sc	\$	503.00
LXZ448.99.00031	Probe Head, Rec. Sd, cplt, NT3x00sc	\$	678.37
LXZ448.99.00032	Coupling, incl Screws, NT3x00sc	\$	150.32
LXZ448.99.00033	Set 5mm Wiper NT3x00sc	\$	319.92
LXZ448.99.00034	Geared motor NT3x00sc	\$	522.27
LXZ448.99.00035	Flash lamp NT3x00sc	\$	1,678.59
LXZ448.99.00036	Set optical Filter, NT3100sc	\$	1,466.60
LXZ448.99.00037	Measuring window NT3x00sc	\$	1,254.61
LXZ448.99.00038	Pressure ring NT3x00sc	\$	144.93
LXZ448.99.00039	O Ring Set for Shaft and Coupling	\$	41.82
LXZ448.99.00040	Mounting Ring, Lense NT3x00sc	\$	27.56
LXZ448.99.00041	Lense NT3x00sc	\$	204.28
LXZ449.99.00001	TEE FTG SS, 1 1/4 in BSPT, GS1440/GS2440EX	\$	182.12
LXZ449.99.00002	SUSPENSION CHAIN AND CARABINERS, 3 m/10 ft, GS1440/GS2440EX	\$	377.73
LXZ449.99.00003	BATTERIES FOR CAX440EX FIELD TRANSMITTER H2S EX, 8 COUNT	\$	1,131.27
LXZ449.99.00004	ADAPTER FOR INLINE MOUNTING, SS, 1 1/4 in BSPT, GS1440/GS2440EX	\$	753.54
LXZ449.99.00005	ADAPTER, INLINE MOUNTING WITH TEE, SS, 1 1/4 in BSPT, GS1440/GS2440EX	\$	940.47
LXZ449.99.00006	CABLE GUARD EXTENSION, 1.4 m/4.6 ft, GS1440/GS2440EX	\$	1,350.97
LXZ449.99.00007	CROSSBEAM, GS1440/GS2440EX	\$	1,350.97
LXZ449.99.00008	CROSSBEAM CLAMPS, GS1440/GS2440EX	\$	377.73
LXZ449.99.00009	CAX440EX FIELD TRANSMITTER H2S EX EXTERNAL ANTENNA	\$	503.00
LXZ449.99.00010	FLUSH FRONT SENSOR CAP FOR SUBMERGED SENSORS, GS1440/GS2440EX	\$	1,171.74
LXZ449.99.00011	AIR FLOW CELL, GS1440/GS2440EX	\$	1,350.97
LXZ449.99.00012	INSTALLATION SPARE PARTS KIT, GS1440/GS2440EX	\$	129.51



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LXZ449.99.00013	POWER/4-20 mA BARRIER, GS2440EX	\$	2,948.62
LXZ449.99.00015	ADAPTER FOR INLINE MOUNTING, SS, 1 1/4 in NPT, GS1440/GS2440EX	\$	753.54
LXZ449.99.00017	CABLE GUARD, WITH CHAIN, GS1440/GS2440EX	\$	1,350.97
LXZ449.99.00018	CAX440EX FIELD TRANSMITTER H2S EX DIN RAIL AND WALL BRACKET (US & CA)	\$	753.54
LXZ449.99.00019	ADAPTER, INLINE MOUNTING WITH TEE, SS, 1 1/4 in NPT, GS1440/GS2440EX	\$	940.47
LXZ449.99.00022	TEE FTG SS, 1 1/4 in NPT, GS1440/GS2440EX	\$	182.12
LXZ449.99.02G05	CABLE ASSY, 5 METER, GAS MEASUREMENT	\$	397.00
LXZ449.99.02G10	CABLE ASSY, 10 METER, GAS MEASUREMENT	\$	487.58
LXZ449.99.02G30	CABLE ASSY, 30 METER, GAS MEASUREMENT	\$	732.34
LXZ449.99.02L05	CABLE ASSY, W/GUARD, 5 METER, LIQUID MEASUREMENT	\$	1,751.82
LXZ449.99.02L10	CABLE ASSY, W/GUARD, 10 METER, LIQUID MEASUREMENT	\$	1,840.48
LXZ449.99.02L30	CABLE ASSY, W/GUARD, 30 METER, LIQUID MEASUREMENT	\$	2,081.38
LXZ449.99.D0001	CAX440EX FIELD TRANSMITTER DATA LICENSE, 1 YEAR	\$	2,004.29
LXZ449.99.D0002	CAX440EX FIELD TRANSMITTER DATA LICENSE, 2 YEAR	\$	3,237.70
LXZ514.53.B0101	Software-RTC DOSING 1 channel	\$	25,477.58
LXZ514.53.B0111	Software-RTC DOSING 2 channel	\$	43,362.00
LXZ514.99.B0101	RTC C/N/RTC-P (1C):Organic load based nutrient dosing to adjust C/N/P combined with effluent nutrient control (1Channel)	\$	25,284.86
LXZ514.99.B0111	RTC C/N/RTC-P (2C):Organic load based nutrient dosing to adjust C/N/P combined with effluent nutrient control (2 channel)	\$	42,938.02
LXZ515.53.A1010	Software-RTC P RTC 1 channel	\$	20,563.22
LXZ515.53.A1011	Software-RTC P RTC 2 channel	\$	34,940.14
LXZ515.53.A1110	RTC P (1COL): Open loop PO4 control, adjusting precipitant flow (1Channel)	\$	20,601.77
LXZ515.53.A1210	RTC P (1CCLCL): Closed loop PO4 control, adjusting precipitant flow considering P_tot / TSS in effluent (1Channel)	\$	30,873.74
LXZ515.53.A1310	RTC P (1COLCL): Combination open/closed loop PO4 control, adjustig precipitant flow (1Channel)	\$	30,873.74
LXZ515.53.A1311	RTC P (2CCLOL): Closed loop PO4 control adjust. precip. flow 1 combined with open loop PO4 control adjust precip. flow 2 (2Channel)	\$	35,017.22
LXZ515.53.A1411	RTC P (2COLCL):Combination open / closed loop PO4 control, adjustig precipitant flow (2Channel)	\$	37,484.04
LXZ515.99.00001	RTC SW upgrade from std. single to std. comb	\$	3,507.50
LXZ515.99.00002	RTC Std. combined Extension	\$	6,089.95
LXZ515.99.00003	RTC basic software configuration on existing hardware (Visualization and data base)	\$	10,040.71
LXZ515.99.00005	RTC Software adoption or modification or extension	\$	2,582.45
LXZ515.99.A1010	RTC P(1CCL):Closed loop PO4 control, adjusting precipitant flow (1Channel)	\$	20,370.50
LXZ515.99.A1011	RTC P(2CCLCL):Closed loop PO4 control, adjusting precipitant flow (2Channel)	\$	34,593.24
LXZ515.99.A1110	RTC P (1COL): Open loop PO4 control, adjusting precipitant flow (1Channel)	\$	20,139.24



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LXZ515.99.A1111	RTC P (2COLCL): Open loop PO4 control, adjusting precipitant flow (2Channel)	\$	33,244.20
LXZ515.99.A1210	RTC P (1CCLCL): Closed loop PO4 control, adjusting precipitant flow considering P_tot / TSS in effluent (1Channel)	\$	30,141.41
LXZ515.99.A1310	RTC P (1COLCL): Combination open/closed loop PO4 control, adjustig precipitant flow (1Channel)	\$	30,141.41
LXZ515.99.A1311	RTC P (2CCLOL): Closed loop PO4 control adjust. precip. flow 1 combined with open loop PO4 control adjust precip. flow 2 (2Channel)	\$	38,813.81
LXZ515.99.A1411	RTC P (2COLCL):Combination open / closed loop PO4 control, adjustig precipitant flow (2Channel)	\$	45,848.09
LXZ515.99.B0000	SOFTWARE- New: Opc Server	\$	12,218.45
LXZ516.53.A0101	Software-RTC SD-RTC 1 channel	\$	25,477.58
LXZ516.53.A0111	Software-RTC SD-RTC 2 channel	\$	43,362.00
LXZ516.99.A0101	RTC SD (1C): Open loop contr. of TSS in dewatered sludge by adjustment of polymer flow or inflow (1Channel)	\$	25,284.86
LXZ516.99.A0111	RTC SD(2C):Open loop contr. of TSS in dewatered sludge by adjustment of polymer flow or inflow (2 Channel)	\$	42,938.02
LXZ517.53.A0101	Software-RTC ST-RTC 1 channel	\$	25,477.58
LXZ517.53.A0111	Software-RTC ST-RTC 2 channel	\$	43,362.00
LXZ517.53.B0101	RTC DAF(1C):Control of TSS in flotated sludge and TSS in clear water by adjustment of coagulant C/N/P, polymer C/N/Ping and pH adjustment (1Channel)	\$	30,989.38
LXZ517.53.B0111	RTC DAF (2C):Control of TSS in flotated sludge and TSS in clear water by adjustment of coagulant C/N/P, polymer C/N/P and pH adjustment (2 Channel)	\$	52,766.74
LXZ517.53.B0201	RTC DAF+ (1C):Control of effluent TSS	\$	35,479.75
LXZ517.99.A0101	RTC ST (1C): Open and closed loop contr. of TSS in thicken sludge and/or filtrate by adj. of polymer flow and/or feed flow (1Channel)	\$	25,284.86
LXZ517.99.A0111	RTC ST (2C):Open and closed loop contr. of TSS in thicken sludge and/or filtrate by adj. of polymer flow and/or feed flow (2Channel)	\$	42,938.02
LXZ517.99.B0101	RTC DAF(1C):Control of TSS in flotated sludge and TSS in clear water by adjustment of coagulant C/N/P, polymer C/N/Ping and pH adjustment (1Channel)	\$	30,738.84
LXZ517.99.B0111	RTC DAF (2C):Control of TSS in flotated sludge and TSS in clear water by adjustment of coagulant C/N/P, polymer C/N/P and pH adjustment (2 Channel)	\$	52,323.48
LXZ517.99.B0201	RTC DAF+ (1C):Control of effluent TSS	\$	35,479.75
LXZ517.99.K0101	SC4500 RTC-SP module upgrade	\$	19,869.43
LXZ518.53.A0101	Software-RTC SLUDGE RETENTION TIME (SRT-RTC) 1 channel	\$	25,477.58
LXZ518.53.A0111	Software-RTC SLUDGE RETENTION TIME (SRT-RTC) 2 channel	\$	43,362.00



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LXZ518.99.A0101	RTC SRT(1C): Adjustment of sludge retention time according to temperature by controlling surplus activated sludge flow (1Channel)	\$ 25,284.86
LXZ518.99.A0111	RTC SRT (2C): Adjustment of sludge retention time according to temperature by controlling surplus activated sludge flow (2 Channel)	\$ 42,938.02
LXZ519.53.B0101	Software-RTC N RTC STANDARD 1 channel	\$ 25,477.58
LXZ519.53.B0111	Software-RTC N RTC STANDARD 2 channel	\$ 43,362.00
LXZ519.53.B1101	Software-RTC N RTC STANDARD WITH O2 CONTROL 1 channel	\$ 46,021.54
LXZ519.53.B1111	Software-RTC N RTC STANDARD WITH O2 CONTROL 2 channel	\$ 78,225.05
LXZ519.53.B2101	Software-RTC N RTC STANDARD WITH O2 AND VARIABLE FREQUENCE DRIVE CONTROL 1 channel	\$ 50,916.62
LXZ519.53.B2111	Software-RTC N RTC STANDARD WITH O2 AND VARIABLE FREQUENCE DRIVE CONTROL 2 channel	\$ 86,608.37
LXZ519.53.D0101	Software-RTC N-RTC WITH 4 DO SETPOINT PROFILE 1 channel	\$ 35,344.85
LXZ519.53.D0111	Software-RTC N-RTC WITH 4 DO SETPOINT PROFILE 2 channel	\$ 60,032.28
LXZ519.53.D1101	RTC N_STERTC-P(1C):Combining open / closed loop NH4 control on Step Feed reactors, providing 3 O2 set point (1Channel)	\$ 49,837.39
LXZ519.53.D1111	RTC N_STERTC-P(2C): Combining open / closed loop NH4 control on Step Feed reactors, providing 3 O2 set point (2 Channel)	\$ 77,839.61
LXZ519.99.B0101	RTC N (1C):Combination open/closed loop NH4 control, providing O2 setpoint (1Channel)	\$ 25,284.86
LXZ519.99.B0111	RTC N (2C):Combination open/closed loop NH4 control, providing O2 setpoint (2Channel)	\$ 42,938.02
LXZ519.99.B1101	RTC N_DO (1C):Combination open/closed loop NH4 control, with O2 control (Step) (1Channel)	\$ 45,616.82
LXZ519.99.B1111	RTC N_DO(2C):Combination open/closed loop NH4 control, with O2 control (Step) (2Channel)	\$ 77,569.80
LXZ519.99.B2101	RTC N_DOVFD (1C):Combination open/closed loop NH4 control, with O2 control (Step & VSD) (1Channel)	\$ 50,531.18
LXZ519.99.B2111	RTC N_DOVFD(2C):Combination open/closed loop NH4 control, with O2 control (Step & VSD)(2 Channel)	\$ 85,895.30
LXZ519.99.D0101	RTC N_4Z (1C): Combination open/closed loop NH4 control, providing profile of 4 O2 setpoints (1Channel)	\$ 35,036.50
LXZ519.99.D0111	RTC N_4Z (2C): Combination open/closed loop NH4 control, providing profile of 4 O2 setpoints (2 Channel)	\$ 59,589.02
LXZ519.99.D1101	RTC N_STERTC-P(1C):Combining open / closed loop NH4 control on Step Feed reactors, providing 3 O2 set point (1Channel)	\$ 49,413.41
LXZ519.99.D1111	RTC N_STERTC-P(2C): Combining open / closed loop NH4 control on Step Feed reactors, providing 3 O2 set point (2 Channel)	\$ 77,165.09
LXZ520.53.C0101	Software-RTC N/DN RTC STANDARD 1 channel	\$ 25,477.58
LXZ520.53.C0111	Software-RTC N/DN RTC STANDARD 2 channel	\$ 43,362.00
LXZ520.53.C1101	Software-RTC N/DN RTC STANDARD WITH O2 CONTROL 1 channel	\$ 46,021.54



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LXZ520.53.C1111	Software-RTC N/DN RTC STANDARD WITH O2 CONTROL 2 channel	\$	78,225.05
LXZ520.53.C2101	Software-RTC N/DN RTC STANDARD WITH O2 AND FREQUENCE DRIVE CONTROL 1 channel	\$	50,916.62
LXZ520.53.C2111	Software-RTC N/DN RTC STANDARD WITH O2 AND FREQUENCE DRIVE CONTROL 2 channel	\$	86,608.37
LXZ520.53.D0101	Software-RTC N/DN-SBR RTC STANDARD 1 channel	\$	26,672.45
LXZ520.53.D0111	Software-RTC N/DN-SBR RTC STANDARD 2 channel	\$	45,443.38
LXZ520.53.D1101	Software-RTC N/DN-SBR RTC STANDARD WITH O2 CONTROL 1 channel	\$	47,235.67
LXZ520.53.D1111	Software-RTC N/DN-SBR RTC STANDARD WITH O2 CONTROL 2 channel	\$	80,344.97
LXZ520.53.D2101	Software-RTC N/DN-SBR RTC STANDARD WITH O2 AND FREQUENCE DRIVE CONTROL 1 channel	\$	52,111.49
LXZ520.53.D2111	Software-RTC N/DN-SBR RTC STANDARD WITH O2 AND FREQUENCE DRIVE CONTROL 2 channel	\$	88,689.74
LXZ520.99.C0101	RTC N/DN (1C): Intermittent aeration control (1Channel)	\$	25,284.86
LXZ520.99.C0111	RTC N/DN (2C): Intermittent aeration control (2 Channel)	\$	42,938.02
LXZ520.99.C1101	RTC N/DN_DO (1C):Intermittent aeration control & O2 control (Step) (1 Channel)	\$	45,616.82
LXZ520.99.C1111	RTC N/DN_DO(2C):Intermittent aeration control & O2 control (Step) (2 Channel)	\$	77,569.80
LXZ520.99.C2101	RTC N/DN_DOVFD1C:Intermittent aeration control & O2 control (Step & VSD) (1Channel)	\$	50,531.18
LXZ520.99.C2111	RTC N/DN_DOVFD2C: Intermittent aeration control & O2 control (Step & VSD) (2 Channel)	\$	85,895.30
LXZ520.99.D0101	RTC N/DNSBR (1C):Intermittent aeration control (SBR) (1Channel)	\$	26,499.00
LXZ520.99.D0111	RTC N/DNSBR (2C):Intermittent aeration control (SBR) (2 Channel)	\$	45,019.39
LXZ520.99.D1101	RTC N/DNSBR_DO (1C):Intermittent aeration control (SBR) & O2 control (Step) (1Channel)	\$	46,830.96
LXZ520.99.D1111	RTC N/DNSBR_DO (2C):Intermittent aeration control (SBR) & O2 control (Step) (2Channel)	\$	79,689.72
LXZ520.99.D2101	RTC N/DNSBR_DOVFD (1C):Intermittent aeration control (SBR) & O2 control (Step & VSD) (1Channel)	\$	51,764.59
LXZ520.99.D2111	RTC N/DNSBR_DOVFD (2C):Intermittent aeration control (SBR) & O2 control (Step & VSD) (2 Channel)	\$	87,976.68
LXZ521.53.A0101	Software-RTC DN-RTC INTERNAL RECIRCULATION 1 channel	\$	20,563.22
LXZ521.53.A0111	Software-RTC DN-RTC INTERNAL RECIRCULATION 2 channel	\$	43,362.00
LXZ521.53.B0101	Software-RTC DN-RTC INTERNAL RECIRCULATION AND EXTERN C 1 channel	\$	25,477.58
LXZ521.53.B0111	Software-RTC DN-RTC INTERNAL RECIRCULATION AND EXTERN C 2 channel	\$	78,225.05
LXZ521.53.C0101	Software-RTC SWINGZONE RTC (ONLY IN COMBINATION WITH N RTC) 1 channel	\$	20,563.22



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LXZ521.53.C0111	Software-RTC SWINGZONE RTC (ONLY IN COMBINATION WITH N RTC) 2 channel	\$	86,608.37
LXZ521.53.D0101	Software-RTC DN-RTC EXTERN C 1 channel	\$	20,563.22
LXZ521.53.D0111	Software-RTC DN-RTC EXTERN C 2 channel	\$	60,032.28
LXZ521.53.E0101	SOFTWARE RTC DN Extern C sec Anox 1 Channel	\$	22,837.32
LXZ521.53.E0111	SOFTWARE RTC DN Extern C sec Anox 2 Channel	\$	38,833.08
LXZ521.53.F0101	SOFTWARE RTC DN IRC Extern C sec Anox 1 Channel	\$	25,207.78
LXZ521.53.F0111	SOFTWARE RTC DN IRC Extern C sec Anox 2 Channel	\$	45,790.27
LXZ521.99.A0101	RTC DN_IRC(1C): Closed loop control NO3, adjusting internal recirculation (1Channel)	\$	20,370.50
LXZ521.99.A0111	RTC DN_IRC (2C):Closed loop control NO3, adjusting internal recirculation (2Channel)	\$	34,593.24
LXZ521.99.B0101	RTC DN_IRC_C (1C):Closed loop control NO3, adjusting internal recirculation and external carbon (1Channel)	\$	25,284.86
LXZ521.99.B0111	RTC DN_IRC_C (2C):Closed loop control NO3, adjusting internal recirculation and external carbon (2Channel)	\$	42,938.02
LXZ521.99.C0101	RTC N_SZ1C(1C):Control of swingzone (only in combination with N_4Z) (1Channel)	\$	20,370.50
LXZ521.99.D0101	RTC DN_C(1C):Closed loop control NO3, adjusting external carbon (1Channel)	\$	20,563.22
LXZ521.99.D0111	RTC DN_C (2C):Closed loop control NO3, adjusting external carbon (2 Channel)	\$	60,032.28
LXZ521.99.E0101	SOFTWARE RTC DN Extern C sec Anox 1 Channel	\$	21,565.37
LXZ521.99.E0111	SOFTWARE RTC DN Extern C sec Anox 2 Channel	\$	36,655.34
LXZ521.99.F0101	SOFTWARE RTC DN IRC Extern C sec Anox 1 Channel	\$	24,880.15
LXZ521.99.F0111	SOFTWARE RTC DN IRC Extern C sec Anox 2 Channel	\$	45,385.56
LXZ522.53.A0101	Software-RTC OXIDATION DITCH (INTENSITY) 1 channel	\$	25,477.58
LXZ522.53.A0111	Software-RTC OXIDATION DITCH (INTENSITY) 2 channel	\$	43,362.00
LXZ522.53.B0101	Software-RTC OXIDATION DITCH (DIFFERENT ZONES) 1 channel	\$	30,854.47
LXZ522.53.B0111	Software-RTC OXIDATION DITCH (DIFFERENT ZONES) 2 channel	\$	52,439.11
LXZ522.99.A0101	RTC SND (1C): NH4 & NO3 control, adjusting aerated volume (intensity outp. converted by PLC) (1Channel)	\$	25,284.86
LXZ522.99.A0111	RTC SND (2C): NH4 & NO3 control, adjusting aerated volume (intensity outp. converted by PLC) Double channel) (2 Channel)	\$	42,938.02
LXZ522.99.B0101	RTC SND (1C6Z): NH4 & NO3 control, adjusting aerated volume by valve or blower adjustment (1Channel)	\$	30,584.66
LXZ522.99.B0111	RTC SND (2C6Z): NH4 & NO3 control, adjusting aerated volume by valve or blower adjustment (2 Channel)	\$	51,918.77
LXZ524.97.00006	ee SC4x00c mA output module, 5 outputs; RoW	\$	1,441.55
LXZ524.97.00007	ee SC4x00c Profibus DP network module, RoW	\$	2,274.10
LXZ524.97.00042	SC4x00c input module, RoW	\$	1,724.84
LXZ524.98.00042	SC4x00c input module, EU	\$	1,724.84
LXZ524.99.00004	SC4x00 UV protection screen	\$	510.71

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LXZ524.99.00005	SC4x00 UV protection screen w/ sun roof	\$	1,015.63
LXZ524.99.00006	SC4x00c mA output module, 5 outputs	\$	1,524.42
LXZ524.99.00007	SC4x00 Profibus DP network module	\$	2,274.10
LXZ524.99.00009	SC4x00 M12/RJ45 connector kit (Eth. Cable)	\$	585.87
LXZ524.99.00010	SC4x00 M12/M12 connector kit (Eth. Cable)	\$	381.59
LXZ524.99.00011	10m Ethernet cable M12 to M12	\$	713.06
LXZ524.99.00012	5m Ethernet cable M12 to RJ45	\$	610.92
LXZ524.99.00013	SC4x00 Holder for Mobiles	\$	713.06
LXZ524.99.00016	SC4200c External USB Box V2	\$	670.67
LXZ524.99.00017	SC4x00 Modules holder	\$	231.26
LXZ524.99.00018	SC4x00 Fuse 1,5A Time Lag	\$	81.91
LXZ524.99.00019	SC4200c Mainboard cover	\$	93.47
LXZ524.99.00021	SC4x00 Door Hinge pin	\$	37.77
LXZ524.99.00022	SC4200c Front door with keypad	\$	695.72
LXZ524.99.00023	SC4200c Front door assembly	\$	368.10
LXZ524.99.00024	SC4200c Keypad	\$	167.09
LXZ524.99.00025	SC4200c Hach badge	\$	34.30
LXZ524.99.00026	SC4200c Back enclosure	\$	373.88
LXZ524.99.00027	SC4200c SC connector kit	\$	389.29
LXZ524.99.00028	SC4200c Ethernet connector kit	\$	389.29
LXZ524.99.00029	SC4200c SC and Ethernet conector cap	\$	32.18
LXZ524.99.00030	SC4200c SC and Ethernet conector hole plug M16	\$	32.18
LXZ524.99.00031	SC4x00 Hole Plug Kit	\$	22.36
LXZ524.99.00032	SC4200c USB cover with screws	\$	107.15
LXZ524.99.00033	SC4x00 Sunroof Visor	\$	300.64
LXZ524.99.00034	SC4200c US-Power Cord Kit, strain relief, 125 VAC	\$	304.50
LXZ524.99.00036	SC4x00 Mounting HW Sunroof w/ visor	\$	1,420.35
LXZ524.99.00037	SC4x00 Sunroof w/ visor	\$	713.06
LXZ524.99.00039	US dongle Verizon kit	\$	2,254.82
LXZ524.99.00041	Mounting screws pack for modem	\$	3.85
LXZ524.99.00043	Remote service connector for sc1000	\$	595.50
LXZ524.99.00046	Remote service connector for sc200	\$	456.75
LXZ524.99.A0008	Claros Connector w/ LAN; 1YR	\$	3,546.05
LXZ524.99.A0009	Claros Connector Internal 4G US AT&T; 1YR	\$	5,357.62
LXZ524.99.A0010	Claros Connector Internal 4G EU; 1YR	\$	4,663.82
LXZ524.99.A0011	Claros Connector Ext USB Box Verizon; 1YR	\$	5,203.44
LXZ524.99.A0012	Claros Connector Ext USB Box Wi-Fi; 1YR	\$	3,758.04
LXZ524.99.A0013	Claros Connector w/ LAN; 3YR	\$	6,263.40
LXZ524.99.A0014	Claros Connector Internal 4G US AT&T; 3YR	\$	11,274.12
LXZ524.99.A0015	Claros Connector Internal 4G EU; 3YR	\$	8,132.78
LXZ524.99.A0016	Claros Connector Ext USB Box Verizon; 3YR	\$	6,263.40
LXZ524.99.A0017	Claros Connector Ext USB Box Wi-Fi; 3YR	\$	6,263.40
LXZ525.99.00001	SC4500 Connector incl. cable nut and cap	\$	132.21
LXZ525.99.00002	SC4500 Ethernet connector incl. cable nut and cap	\$	133.75
LXZ525.99.00003	SC4x00 Profinet+IEP connect w/cable,nut,cap,gasket	\$	139.14
LXZ525.99.00004	SC4500 USB box, nut for cable gland	\$	40.86
LXZ525.99.00005	SC4x00 Ribbon Cable	\$	33.53



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LXZ525.99.00006	SC4500 DC Fuse Kit	\$	89.61
LXZ525.99.00008	SC4500 M20 Hole Plug	\$	63.40
LXZ525.99.00009	SC4500 Module holder for Profinet and Ethernet IP	\$	189.06
LXZ525.99.00010	SC4500 Front door assembly incl. Touchscreen + screws	\$	373.88
LXZ525.99.00011	SC4500 Back enclosure	\$	158.99
LXZ525.99.00012	SC4500 USB cover (incl. screws)	\$	59.36
LXZ525.99.00013	SC4500 Main board cover	\$	92.51
LXZ525.99.00014	SC4500 Set of screws	\$	37.19
LXZ525.99.00015	SC4x00 and Ethernet connector cap	\$	60.90
LXZ525.99.00016	SC4500 Hach Badge	\$	27.56
LXZ525.99.00017	SC4x00 USB Stick	\$	92.51
LXZ525.99.00018	SC4500 Wi-Fi USB Box RoW	\$	1,225.70
LXZ525.99.00019	SC4500 USB-Box cable alone (4.5m) - WiFi	\$	516.49
LXZ525.99.00020	SC4500 external antenna for USB box	\$	343.04
LXZ525.99.00021	SC4x00 Plugs kit for openings	\$	53.38
LXZ525.99.00022	Sc-Controller M16 Cable Gland Kit	\$	70.34
LXZ525.99.00023	SC4x00 Set of front door screws	\$	57.05
LXZ525.99.00024	1/2 NPT Strain relief/Cable gland kit	\$	74.97
LXZ525.99.00025	SC4500 Mounting Plate, USB Box	\$	56.27
LXZ525.99.00026	sc4500 10m Ethernet cable M12 to M12 / C1D2	\$	553.11
LXZ525.99.00027	SC4500 KIT-Mounting Bracket, WIFI - USB Box	\$	888.44
LXZ525.99.00028	SC4500 KIT-Mounting Screws TORX, USB Box	\$	17.34
LXZ525.99.00029	SC4500 KIT-Mounting Bracket, GSM - USB Box	\$	1,759.53
LXZ525.99.00030	SC4500 Kit Bolt C1D2	\$	1,645.83
LXZ525.99.00031	SC4500 USB-Box cable alone (2m) - GSM	\$	637.90
LXZ525.99.A0001	SC4500 External USB box EU WiFi	\$	1,225.70
LXZ525.99.A0002	SC4500 External USB box NA WiFi	\$	1,225.70
LXZ525.99.A0011	SC4500 External USB box EU GSM	\$	2,505.36
LXZ525.99.A0012	SC4500 External USB box NA GSM	\$	2,505.36
LXZ525.99.B0001	SC4500 Power Cord Kit EU plug	\$	223.56
LXZ525.99.B0002	SC4500 Power Cord Kit US plug	\$	202.36
LXZ525.99.B0003	SC4500 Power Cord Kit UK plug	\$	227.41
LXZ525.99.C0001	SC4500 Profinet IO Upgrade Kit	\$	2,331.91
LXZ525.99.C0002	SC4500 Ethernet IP Upgrade Kit	\$	2,331.91
LXZ525.99.C0003	SC4500 LAN Ethernet Upgrade Kit	\$	312.21
LXZ525.99.C0004	SC4500 Prognosys Upgrade Kit	\$	1,645.83
LXZ525.99.C0005	dd SC4500 Modbus TCP/IP Upgrade Kit	\$	1,270.02
LXZ525.99.C0006	dd SC4500 Modbus TCP/IP Upgrade Kit - License Only	\$	1,021.42
LXZ525.99.D0002	SC4500 mA Output module (5 Outputs)	\$	487.58
LXZ525.99.D0003	SC4500 pH/ORP module	\$	1,312.42
LXZ525.99.D0004	SC4500 COND module	\$	1,312.42
LXZ530.53.A0111	RTC MOV (12C):Closed loop zone DO controller acting on corresponding air valve position (12 valves)	\$	42,822.38
LXZ530.53.B0101	RTC MOV (8C):Closed loop zone DO controller acting on corresponding air valve position (8 valves)	\$	33,745.27
LXZ530.53.B0111	RTC MOV (16C):Closed loop zone DO controller acting on corresponding air valve position (16 valves)	\$	51,109.34



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LXZ530.99.A0101	RTC MOV (4C): Closed loop zone DO controller acting on corresponding air valve position (4 valves)	\$	24,880.15
LXZ530.99.A0111	RTC MOV (12C):Closed loop zone DO controller acting on corresponding air valve position (12 valves)	\$	41,878.06
LXZ530.99.B0101	RTC MOV (8C):Closed loop zone DO controller acting on corresponding air valve position (8 valves)	\$	32,993.66
LXZ530.99.B0111	RTC MOV (16C):Closed loop zone DO controller acting on corresponding air valve position (16 valves)	\$	49,895.21
LXZ530.99.C0101	RTC DO (4C):Closed loop zone DO, controller (4 actors)	\$	22,818.05
LXZ530.99.D0101	RTC DO (8C):Closed loop zone DO, controller (8 actors)	\$	27,096.43
LXZ531.53.A1010	RTC-C/DC (1C): Chlor/dechlor RTC for wastewater	\$	27,751.68
LXZ531.53.A1011	RTC-C/DC (2C): Chlor/dechlor RTC for wastewater	\$	46,484.06
LXZ531.99.A1010	RTC-C/DC (1C): Chlor/dechlor RTC for wastewater	\$	28,060.03
LXZ531.99.A1011	RTC-C/DC (2C): Chlor/dechlor RTC for wastewater	\$	47,042.95
LXZ532.53.A1010	Coagulation RTC for Source Water small, Software	\$	22,162.80
LXZ532.53.A1011	Coagulation RTC for Source Water medium, Software	\$	31,394.09
LXZ532.53.A1012	Coagulation RTC for Source Water large, Software	\$	40,567.56
LXZ532.53.B1010	Pre-oxidation RTC for Source Water small, Software	\$	9,828.72
LXZ532.53.B1011	Pre-oxidation RTC for Source Water medium, Software	\$	15,976.49
LXZ532.53.B1012	Pre-oxidation RTC for Source Water large, Software	\$	22,162.80
LXZ541.99.0001H	ee KIT, MOUNTING BRACKETS, PAH500	\$	360.39
LYA0020	STANDARD FAT ORBISPHERE AND POLYMETRON ANALYZERS	\$	3,218.42
LYP106.99.01003	Sample Conditioning Panel, <300 PSIG(20.7 bar), <600F (315C)	\$	23,916.55
LYP106.99.01004	Sample Conditioning Panel, <600 PSIG(41.4 bar), <600F (315C)	\$	19,773.07
LYP106.99.01005	Sample Conditioning Panel, <3000 PSIG(206.8 bar), <1000F (537C)	\$	24,109.27
LYP514.99.00079	RTC CUSTOM PROJECT LINCOLN NE - 15.6" wide touch screen controller	\$	260,788.70
LYW185	K ⁿ vetten-Test Menthol	\$	225.84
LYW964	PIPET STAND	\$	225.48
LZE100	Syringe 2,5mL (Glass, PTFE Luer L.)	\$	450.96
LZE101	Syringe 5mL (Glass, PTFE Luer L.)	\$	481.80
LZE102	Syringe 10mL (Glass, PTFE Luer L.)	\$	481.80
LZE103	Syringe 25mL (Glass, PTFE Luer L.)	\$	1,422.27
LZE104	Ring for 2,5mL Syringe	\$	8.09
LZE105	Ring for 5mL Syringe	\$	8.09
LZE106	Ring for 10mL Syringe	\$	8.09
LZE107	Ring for 25mL Syringe	\$	8.09
LZE108	Set of Beaker 50mL (5x) for AT	\$	35.85
LZE109	Set of Beaker 150mL (5x) for AT	\$	39.31
LZE110	KF beaker, 150mL, Glass	\$	310.28
LZE112	Delivery Tubing with Anti Diffusion	\$	223.56
LZE113	Delivery Tubing with Anti Diffusion	\$	281.37
LZE114	Set of 2 tubing for Syringe,	\$	150.51
LZE115	Set of 2 tubing for Pump, Suction &	\$	150.51
LZE116	Bottle Stopper GL25	\$	87.30
LZE117	Bottle Stopper GL45	\$	87.30



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LZE118	Bottle Stopper S40	\$	114.86
LZE119	Glass Bottle, 1L	\$	51.65
LZE120	Dessicant Tube Empty with cotton	\$	8.09
LZE121	Dessicant Tube Filled with Molecular Sieve	\$	35.85
LZE122	Conical Adapter NS14.5/23	\$	31.80
LZE123	Stopper 20 mm, KF Cell	\$	154.56
LZE124	Tube Gland (2 pcs), KF Cell	\$	31.03
LZE125	Spare pump cassette, AT model	\$	472.16
LZE126	Spare pump cassette, KF model	\$	472.16
LZE127	TitraMaster 1000 Series, TitraLab	\$	2,543.90
LZE128	USB key Apps Water	\$	256.32
LZE129	USB key Apps Biogas	\$	4,586.74
LZE130	USB key Apps Amperometric	\$	358.46
LZE131	USB key Food & Beverage Apps	\$	256.32
LZE132	USB key Apps Karl Fischer	\$	356.53
LZE133	USB key Apps PetroChem	\$	256.32
LZE134	USB key Apps General Chemistry &	\$	256.32
LZE135	USB key Apps Plating	\$	256.32
LZE136	Magnetic Bar, 20x6mm (x5)	\$	55.12
LZE137	Set of Accessories, AT, 10 mL	\$	1,067.67
LZE138	Set of Accessories, AT, 5 mL	\$	950.11
LZE139	Set of Accessories, AT, 25 mL	\$	2,023.56
LZE140	Set of Accessories, KF, 5 mL	\$	1,306.64
LZE141	Tube holder, AT, 4 positions.	\$	10.41
LZE142	Sample Leveling Pump	\$	1,738.33
LZE143	Propeller stirrer, 80mm shaft	\$	1,029.12
LZE144	Power Supply for AT1000 Tisch-Netzgerat	\$	206.21
LZE145	TitraLab RS-232 Adapter cable	\$	126.62
LZE146	Top cover AT kit	\$	1,173.66
LZE147	Top cover KF kit	\$	1,173.66
LZE148	Bottom Cover Kit	\$	468.31
LZE149	Feet Serie SJ5000 (Set of 4pcs)	\$	35.65
LZE150	USB cover (2pcs)	\$	35.65
LZE151	AT Pump cassette: 100ml/min ? PharMed	\$	738.12
	KF Pump cassette: 170ml/min KF Pump cassette: 170ml/min ?		
LZE152	Novoprene	\$	898.08
LZE153	Display kit	\$	2,659.54
LZE154	Front cover + badge	\$	55.12
LZE155	Stickers set (Side USB r panel)	\$	55.12
LZE156	Z axis protection window	\$	55.12
LZE157	Z axis syringe support 1 seringue	\$	206.21
LZE158	Z axis syringes support 2 seringues	\$	316.06
LZE159	Syringe screw	\$	67.64
LZE161	Belt Synchroflex 10T5/215	\$	63.40
LZE162	Opto sensor OPB829C + Connector	\$	100.79
LZE163	Z axis Silicone Gasket	\$	83.83
LZE164	Display Gasket	\$	83.83



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LZE165	Valve+connector	\$	1,719.06
LZE166	Stirrer assembly complete w magnet, motor, holder and cable	\$	655.25
LZE167	Keyboard	\$	377.73
LZE168	Cable - Controlboard to Motherboard	\$	252.46
LZE169	Set of screws + clips for boards +	\$	80.56
LZE170	Set of Cables	\$	96.55
LZE171	Z-Axis assembly complete with motor	\$	2,717.35
LZE173	Probe holder parts AT1xxx	\$	418.20
LZE174	Probe holder parts KF1xxx	\$	1,135.12
LZE175	Pump connector + circlips	\$	83.83
LZE176	Pump door w/sticker	\$	167.86
LZE177	Z-Axis assembly complete with motor 2 syringes	\$	2,910.07
LZE179	Pump door area gasket	\$	48.95
LZE180	Motherboard updater kit	\$	316.06
LZE181	Application development kit	\$	126.62
LZE182	Conical Adapter with large cut (for diameter 7.5)	\$	31.80
LZE183	KF Conical adapter	\$	31.80
LZE184	Adapter BNC male to Banana	\$	79.02
LZE185	Storage Tube Electrode	\$	24.09
LZE186	Bottle stopper GL32	\$	87.30
LZE187	Sample Leveling Pump stand-alone	\$	1,285.44
LZE188	Set of tubes for sample leveling pump	\$	597.43
LZE189	Spare pump cassette, sample leveling pump	\$	472.16
LZE191	Holder Hach sampler, one probe	\$	1,233.41
LZE192	Holder Hach sampler, two probes	\$	1,042.62
LZE193	Set of beakers 50 mL PP, Sampler 30	\$	568.52
LZE194	Set of beakers 90 mL PP, Sampler 20	\$	377.73
LZE195	TitraLab 1000 Series, RS232 cable with USB adapter	\$	643.68
LZE196	Titralab RS232 cable, Balance connection	\$	134.71
LZE197	Zaxis assembly complete, AS1000	\$	4,278.38
LZE198	Set of housing for sampler, AS1000	\$	809.42
LZE199	Electronic board, AS1000	\$	1,009.85
LZE200	Screw to fix Tray, AS1000	\$	329.55
LZE201	Carrousel 20 positions, 90mL beakers, AS1000	\$	601.29
LZE202	Carrousel 30 positions, 50 mL beakers, AS1000	\$	1,761.46
LZE203	Carrousel 12 positions, 150ml beaker, AS1000	\$	2,023.56
LZE204	Carrousel 12 positions, 50ml beaker, AS1000	\$	3,122.06
LZE205	Rotating module 20 positions, AS1000	\$	6,012.86
LZE206	Rotating module 30 positions, AS1000	\$	6,012.86
LZE207	Rotating module 12 positions, AS1000	\$	6,012.86
LZE208	Bottom base plate with stirrer, AS1000	\$	1,310.50
LZE211	Set of cables, AS1000	\$	346.90
LZE212	Stepper motor, AS1000	\$	601.29
LZE213	Beaker ring adapter, 50 mL, AS1000 (x12pcs)	\$	1,491.65
LZE217	Magnetic bar, 12x4.5mm (x30)	\$	988.65
LZE219	Set of Accessories, AT, 1 syringe common - prepackage	\$	626.34



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LZE220	Set of Accessories, AT, 1 syringe hardness + petrochem - prepackage	\$	618.63
LZE221	Set of Accessories, AT, 2 syringes - prepackage	\$	1,262.32
LZE222	Set of Accessories, AT, AutoCAT - prepackage	\$	1,104.29
LZE223	Set of Accessories, AT, Karl Fischer - prepackage	\$	1,137.05
LZE225	Set for AT1000, Single 25mL Syringe	\$	1,780.73
LZE227	Replacement O-Ring for KF1000 Volumetric Karl Fischer Titrator	\$	16.57
LZE228	AT1000 Syringe Volume Verification KIT	\$	1,509.00
LZH165	DIN RAIL, NS 35/15, 35 CM	\$	61.86
LZH166	POWER SUPPLY, 90-240 VAC/24 V DIN RAIL MOUNT	\$	356.53
LZH167	TERMINAL BLOCK, FEED-THROUGH, W/O POWER SUPPLY	\$	44.90
LZH168	GROUNDING TERMINAL BLOCK	\$	58.59
LZH169	SUB-D CONNECTOR, FOR RTC	\$	346.90
LZH171	CPU MODULE, WITH ETHERNET, RS422/RS485	\$	11,389.75
LZH172	POWER SUPPLY MODULE, I/O INTERFACE	\$	2,678.81
LZH173	DIGITAL OUTPUT MODULE, 24 V DC, 2 OUTPUTS	\$	260.17
LZH174	DIGITAL OUTPUT MODULE, 24 V DC, 4 OUTPUTS	\$	300.64
LZH175	ANALOG OUTPUT MODULE, 1 OUTPUT	\$	1,137.05
LZH176	ANALOG OUTPUT MODULE, 2 OUTPUTS	\$	1,603.43
LZH177	ANALOG INPUT MODULE, 1 INPUT	\$	1,067.67
LZH178	BUS TERMINATION MODULE, FOR RTC	\$	129.12
LZH204	DIGITAL INPUT MODULE, 24 V DC, 2 INPUTS	\$	191.37
LZH205	DIGITAL OUTPUT MODULE, 24 V DC, 8 OUTPUTS	\$	493.36
LZH206	DIGITAL OUTPUT MODULE, 24 V DC, 16 OUTPUTS	\$	795.93
LZH216	FERRITE, RTC MODULE	\$	65.33
LZH409.99.00100	dd Scrubber Monitor SM100	\$	49,471.22
LZH409.99.00110	dd Scrubber Monitor SM100 no display	\$	41,800.97
LZH409.99.00210	dd Scrubber Monitor SM200 no display	\$	45,366.29
LZH409.99.00300	dd Scrubber Monitor SM300	\$	63,636.14
LZH409.99.00310	dd Scrubber Monitor SM300 no display	\$	55,985.16
LZH409.99.00400	dd Scrubber Monitor SM400	\$	81,308.57
LZH409.99.00410	dd Scrubber Monitor SM400 no display	\$	73,638.31
LZH409.99.00500	dd Scrubber Monitor SM500	\$	88,400.66
LZH409.99.00510	dd Scrubber Monitor SM500 no display	\$	80,711.14
LZH409.99.00600	dd Scrubber Monitor SM600	\$	56,563.32
LZH409.99.00610	dd Scrubber Monitor SM600 no display	\$	48,912.34
LZH409.99.00700	dd Scrubber Monitor SM700	\$	84,835.34
LZH409.99.00710	dd Scrubber Monitor SM700 no display	\$	77,184.36
LZH409.99.00800	dd Scrubber Monitor SM800	\$	45,944.45
LZH409.99.00810	dd Scrubber Monitor SM800 no display	\$	38,274.19
LZI-12012	aa OPTICAL HOLDER, TSS TITANIUM7 sc	\$	7,824.43
LZI-12013	aa aa OPTICAL HOLDER, TSS sc, NO WIPER	\$	5,261.26
LZI-12014	aa aa OPTICAL HOLDER, TSS W sc, WITH WIPER	\$	6,610.30
LZI-12015	aa OPTICAL HOLDER, TSS XL sc	\$	6,610.30
LZI-12016	aa OPTICALHOLDER, TSS VARI sc	\$	7,612.44
LZI-12017	aa MOTOR FOR TSS sc	\$	479.87
LZI-12018	sc Kabel 10m	\$	678.37

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LZI-12021	Kabel 10 m mit Markierung	\$	551.18
LZI-66.11478	aa Schraubenset	\$	23.13
LZI-66.11479	aa Dichtungsset	\$	9.83
LZI-66.11480	Montageteileset	\$	319.92
LZI-66.11485	aa VARIVENT SEAL, TSS sc	\$	42.01
LZI-66.11487	aa Schutzbecher	\$	9.83
LZI-66.11492	aa Sendep rint	\$	1,131.27
LZI-66.11493	aa Sendep rint	\$	1,487.80
LZI-66.11494	aa PCB TSS/COSMOS-25, Recipient print	\$	830.62
LZP045	Sample cell, 10 mm rectangular, glass, pk/3	\$	316.06
LZP065	20mm DIGESTION VIALS pk/5	\$	59.98
LZP273	Thermal Paper Roll for Thermal Printer	\$	112.93
LZP277	HALOGEN LAMP, AMTAX INTER 2, 1 piece	\$	501.07
LZP303	LOCK SETS, QUICK RELEASE, FILTRAX	\$	383.51
LZP320	KTO: PIPET SET	\$	2,139.19
LZP335	Flow-Through Cell 370 =	\$	2,543.90
LZP336	FLOW-THROUGH CELL/QUARTZ/50MM	\$	3,218.42
LZP341	CUVETTE, SEMIMICRO 50MM PLASTIC W/CAP	\$	178.07
LZP361	SAMPLE OVERFLOW VESSEL,AMTAX/PHOSPHX	\$	711.14
LZP365	SOLENOID STIRRING STICK (8x3 mm)	\$	78.44
LZP368	Optikblock	\$	2,659.54
LZP373	FLOW THROUGH CUVETTE FOR PHOSPHAX INTER 2	\$	2,370.46
LZP376	SET OF WEARING PARTS, AMTAX INTER	\$	2,890.80
LZP452	CUVETTE, NEPHLA TURBIDIMETER	\$	551.18
LZP538	ELBOW FITTING, OUTLET NITRATAX FLOW CELL	\$	174.60
LZP618	DOUBLE FITTING	\$	144.73
LZP777	PUMP CARTRIDGE, FILTRAX	\$	194.65
LZP782	SET OF SEALS, ULTRA TURB	\$	291.01
LZP847	TEFLON VALVE PISTON FOR PHOSPHAX	\$	174.60
LZP865	CONNECTOR FOR STEP MOTOR	\$	141.07
LZP882	ff TUBE SET, PSIGMA, FIRST DELIVERY	\$	759.32
LZQ066	OPERATOR RFID TAG, PK/2	\$	68.03
LZQ067	SAMPLE RFID TAG, BLACK, PK/3	\$	138.18
LZQ068	SAMPLE RFID TAG, YELLOW, PK/3	\$	138.18
LZQ069	SAMPLE RFID TAG, RED, PK/3	\$	138.18
LZQ070	SAMPLE RFID TAG, BLUE, PK/3	\$	138.18
LZQ071	SAMPLE RFID TAG, GREEN, PK/3	\$	138.18
LZQ072	SAMPLE LOCATION RFID TAG, PK/5	\$	210.06
LZQ090	FLEXIBLE-TUBE PUMP/MOTOR	\$	1,125.48
LZQ091	HOUSING SET SIPPER SIP 10	\$	310.28
LZQ092	SCREW SET FOR SIPPER SIP 10	\$	45.67
LZQ097	TUBE SET SIPPER DR 6000 ACID RESISTENT	\$	763.17
LZQ098	TUBE SET SIPPER DR 6000 ACID RESISTENT	\$	763.17
LZQ100	TUBE SET SIPPER DR 6000 DRINKING WATER	\$	763.17
LZQ102	TUBE SET SIPPER DR 6000 1 INCH POUR THRU CELL	\$	377.73
LZQ104	USB CONNECTION CABLE WITH FERRIT SIPPER SIP 10	\$	101.37
LZQ105	LID FOR SIPPER APPLICATION DR 6000	\$	289.08

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LZQ107	RACK FOR SIPPER SIP 10	\$	289.08
LZQ109	Folientastatur Locator (LOC 100)	\$	135.29
LZQ110	Gehäuse komplett (LOC 100)	\$	337.26
LZQ111	Verbinder Batteriekontakte	\$	47.99
LZQ170	Stage Assy ACM TU5	\$	1,331.70
LZQ171	Plug ACM TU5	\$	91.73
LZQ176	Fiber Wiper ACM TU5x00sc	\$	93.28
LZQ177	Cover Clip Set ACM TU5x00sc	\$	59.94
LZT032	HEATER BLOCK DRB200, 9X13MM & 2X20MM	\$	2,678.81
LZT033	HEATER BLOCK, DRB200, 15 X 13 MM	\$	2,871.53
LZT034	Heizblock 6x13 u. 4x20 / Lange	\$	2,775.17
LZT036	HEATING BLOCK 15X16	\$	3,044.98
LZT039	TOP COVER ENCLOSURE W/KEYPAD	\$	1,742.19
LZT041	COVER, DRB200, 9X13MM & 2X20MM	\$	147.82
LZT042	COVER, DRB200, 15 X 13MM	\$	397.00
LZT043	Abdeckung 6x13 und 4x20	\$	408.57
LZT047	Abdeckung ohne Öffnungen	\$	292.93
LZT048	Schutzklappe, transparent	\$	260.17
LZT115	INSULATION SET, DRB200	\$	174.60
LZT116	BOTTOM HOUSING, DRB200	\$	281.37
LZT130	Heizplatte / Heating plate	\$	400.86
LZT144	ADAPTER SLEEVE DRB 200 FOR MICRO DIST	\$	528.05
LZU1003.99	COAXIAL CABLE, AD9 CONNECTOR, 3M	\$	316.06
LZU1005.99	COAXIAL CABLE, AS9 CONNECTOR, 5M	\$	373.88
LZU1010.99	COAXIAL CABLE, AS9 CONNECTOR, 10M	\$	481.80
LZU1015.99	COAXIAL CABLE, AS9 CONNECTOR, 15M	\$	597.43
LZU1020.99	COAXIAL CABLE AS9 CONNECTOR, 20M	\$	639.83
LZU1025.99	COAXIAL CABLE, AS9 CONNECTOR, 25M	\$	757.39
LZU1099.99	AS9 connector, for coaxial cable °5	\$	210.06
LZU1103.99	MULTIPLE CABLE, VP CONNECTOR, 3M	\$	713.06
LZU1130.99	MULTIPLE CABLE, MP4, FOR AMP CI, 3M	\$	514.56
LZU1131.99	MULTIPLE CABLE, MP4, FOR AMP CI, 5M	\$	535.76
LZU1132.99	Multiple cable with MP-4 con., 10 m	\$	597.43
LZU1133.99	MULTIPLE CABLE, MP4, FOR AMP CI, 15M	\$	639.83
LZU215.99.1E120	Flow through fitting; PP-H; 3/4" inner thread; cond. sensor; DN25 coupling	\$	3,141.34
LZU300.99.00000	db BALL VALVE ARMATURE FOR TSS sc TRICLAMP	\$	8,402.59
LZU301.99.00000	BALL VALVE ARMATURE FOR TSS EX1 sc TRICLAMP	\$	8,402.59
LZU302.99.00000	WELDED BRANCH PIECE, TSS sc TRICLAMP	\$	506.85
LZU302.99.00060	ff WELDED BRANCH DN200, TSS sc TRICLAMP	\$	1,283.52
LZU302.99.10000	WELDED BRANCH FOR TSS XL sc	\$	506.85
LZU303.99.00000	Screwed Branch SS for TSS XL (Open tank)	\$	506.85
LZU304.99.00020	ff MEASURING TUBE DN80 FOR TSS VARI sc	\$	5,627.42
LZU304.99.00030	MEASURING TUBE DN100 FOR TSS VARI sc	\$	6,128.50
LZU4420.99	CONDUCTIVITY MEASUREMENT BOARD	\$	1,695.94
LZU4494.99	PROFIBUS DP OUTPUT, MM44	\$	2,042.83
LZU5300.97.0002	5300 pH ELECTRODE, GEL, PLASTIC	\$	535.76

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LZU5303.97.0002	5303 pH ELECTRODE, GEL	\$	1,029.12
LZU5330.97.0002	5330 pH ELECTRODE, GEL	\$	1,298.93
LZU5334T.97.002	5334T pH ELECTRODE, PT1000, MP4	\$	1,464.67
LZU5335.97.0002	5335 pH ELECTRODE, SOLID ELECTROLYTE	\$	1,298.93
LZU5336T.97.002	5336T, pH ELECTRODE, Pt1000, MP4	\$	1,298.93
LZU5338.97.0002	5338 pH ELECTRODE, GEL, PLASTIC	\$	639.83
LZU5350.97.0002	5350 REDOX ELECTRODE, GEL	\$	1,298.93
LZU5354.97.0002	5354 REDOX ELECTRODE, GOLD, GEL	\$	925.06
LZU5355.97.0002	5355 REDOX ELECTRODE, Pt, SOLID	\$	1,298.93
LZU5358.97.0002	ORP Pt electrode 5358, plastic, gel	\$	639.83
LZU5359.99	Galvanic free chlorine cell with	\$	2,794.44
LZU5361.97.0002	5361 REDOX ELECTRODE, GEL	\$	535.76
LZU5362.97.0002	5362 REDOX ELECTRODE, GEL	\$	1,029.12
LZU5364.97.0002	5364 pH ELECTRODE, DIFF MEDIA	\$	1,626.56
LZU5388.97.0002	5388 COND CELL, Pt1000, MP5	\$	1,965.74
LZU5390.97.0002	Cond cell 5390 with Pt1000, C=0.1,	\$	2,004.29
LZU5392.97.0002	Cond cell 5392 with Pt1000, C=0.3,	\$	2,004.29
LZU5395.97.0002	5395 COND CELL, S8	\$	1,676.66
LZU5396.97.0002	5396 COND CELL, S8	\$	1,140.90
LZU5398.97.0002	5396 COND CELL, Pt1000, MP5	\$	2,004.29
LZU5399.97.0002	Cond cell 5399 with Pt1000, C=10,	\$	2,062.10
LZU5525.99	DOUBLE TEMP SENSOR, Pt100/Pt1000	\$	703.43
LZU5526.99	TEMP SENSOR, Pt1000	\$	703.43
LZU6000.99	FLOATATON HOUSING, FI 60	\$	1,306.64
LZU6001.99	FLOW CHAMBER CR60, 2 SENSORS, W/O THREAD	\$	925.06
LZU6002.99	FLOW CHAMBER CR60, 2 SNENSORS, W/THREAD	\$	925.06
LZU6003.99	FLOW CHAMBER CR60, 3 SENSORS, W/THREAD	\$	1,081.16
LZU6006.99	FLOW CHAM CR60, W/FLOW MON, 2 SENSORS PG	\$	1,356.75
LZU6048.99	MEMBRANE 6050 AMP DO SENSOR	\$	634.05
LZU6071.99	Spare membrane for amperometric	\$	595.50
LZU6101.99	IMMERSION HOUSING IN 61, PVD, 800MM	\$	790.15
LZU6108.99	THREE-PIECE ADAPTER DN50, PVC	\$	46.64
LZU6109.99	WELD-IN-SOCKET, PVC, STRAIGHT	\$	265.95
LZU6110.99	WELD-IN-SOCKET, PVC, INCLINED	\$	331.48
LZU6111.99	THREE-PIECE ADAPTER DN75, PVC	\$	102.91
LZU6114.99	WELD-IN-SOCKET, PP, STRAIGHT	\$	316.06
LZU6115.99	WELD-IN-SOCKET, PP, INCLINED	\$	316.06
LZU6118.99	THREE-PIECE ADAPTER CN63, PVC	\$	46.64
LZU6121.99	WELD-IN-SOCKET, SS, STRAIGHT	\$	373.88
LZU6122.99	WELD-IN-SOCKET, SS, INCLINED	\$	435.55
LZU6125.99	ADAPTER, SENSORS, PG13, 3 PIECE DN20	\$	46.64
LZU6126.99	SCREW CAP, 3-PIECE ADAPTER DN20	\$	46.64
LZU6131.99	IMMERSION HOUSING IN 61, PVC, 1600MM	\$	1,085.01
LZU6191.99	CABLE GLANDS, COAXIAL CABLE 7, PK/2	\$	46.64
LZU63.99	TUBE TO TUBE JOINT	\$	46.64
LZU6301.99	IMMERSION FITTING IN 63, PVC, 800MM	\$	1,468.53
LZU6331.99	IMMERSION FITTING IN 63, PVC, 1600MM	\$	1,798.08



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LZU6899.99	PG13.5 SCREW SEALING, PVC	\$	79.59
LZU7500.97.0002	si075, pH SIMULATOR, AS7/AS9 CONNECTORS	\$	1,611.14
LZU7501.97.0002	si075 CASE, pH SIM, AS7/AS9 CONNECTOR	\$	1,907.93
LZU7504.97.0002	si075 CASE, pH SIM, AS7/AS7/AS9, CAL CERT	\$	2,717.35
LZU7600.99	Sensor fitting, PVC, PG13.5 for	\$	435.55
LZU7601.99	Insertion fitting In76, PVC, 3/4"	\$	759.32
LZU7801.99	Insertion fitting In 78, SS,	\$	2,447.54
LZU9044.99	Multiple cable with MP-5 connector,	\$	373.88
LZU9045.99	Multiple cable with MP-5 connector,	\$	429.77
LZU9046.99	MULTIPLE CABLE, MP5, 10M	\$	535.76
LZU9047.99	Multiple cable with MP-5 connector,	\$	639.83
LZU9055.99	AS7 CABLE, BNC, 1M	\$	210.06
LZU9059.99	AS7-MP5 CONNECTION CABLE, 1M	\$	316.06
LZU9188.99	BNC CONNECTOR, L-SHAPED, F COAXIAL	\$	46.64
LZU9251.99	Multiple cable with MP-5 connector,	\$	686.08
LZU9256.99	MP4 CONNECTOR, COAXIAL, 2 TEMP WIRE	\$	300.64
LZU9300.99	MULTIPLE CABLE, MP4, 3M	\$	481.80
LZU9301.99	MULTIPLE CABLE, MP4, 5M	\$	535.76
LZU9302.99	MULTIPLE CABLE, MP4, 10M	\$	639.83
LZU9303.99	Multiple cable with MP-4 connector,	\$	703.43
LZU9340.99	GLASS PEARLS, 3MM DIAM., 100/PK	\$	46.64
LZU9450.99	AS7 CABLE, W/OPEN STRIPS, 1M S175	\$	210.06
LZU962.99	962 CONNECTOR, COAXIAL	\$	354.60
LZU9810.99	ELECTROLYTE, 6070 DO, 30 ML FLASK	\$	153.98
LZV009	STIRRER MOTOR, AMTAX INTER 2	\$	664.88
LZV148	PROFIBUS CARD	\$	2,409.00
LZV149	MAINTENANCE KIT, AMTAX	\$	4,856.54
LZV150	MAINTENANCE KIT, PHOSPHAX	\$	2,505.36
LZV168	PUMP HEAD 1-CHANNEL	\$	1,848.18
LZV169	PUMP HEAD 1-CHANNEL	\$	1,362.53
LZV172	LARGE SCREW FOR CUVETTE (PHOSPHAX)	\$	231.26
LZV173	SMALL SCREW FOR CUVETTE (AMTX,PHOSX)	\$	352.68
LZV174	Kleinteilesatz	\$	200.43
LZV177	12V DIAPHRAGM PUMP FOR AMTAX COMPACT	\$	1,840.48
LZV178	CUVETTE, MEASURING (AMTAX)	\$	1,225.70
LZV179	CUVETTE, MEASURING (PHOSPHAX)	\$	2,042.83
LZV180	EXPULSION CUVETTE (AMTAX)	\$	1,129.34
LZV181	PUMP CARTRIDGE (AMTAX,PHOSPHAX)	\$	302.57
LZV182	PUMP ROLLERS, PK/5 (AMTAX,PHOSPHAX)	\$	319.92
LZV184	Jahresverschleisteilsatz, 2-Kanal	\$	589.72
LZV205	Luftventil 2/2-Wege	\$	851.82
LZV275	WIPER BLADES, 4 SETS, ULTRATURB SC	\$	289.08
LZV281	AMTAX INTER 2, SPARE PART KIT, 1 YEAR	\$	3,083.52
LZV299	WIPER HOLDER, ULTRA TURB PLUS	\$	1,225.70
LZV302	MEASURING CHAMBER, ULTRA TURB PLUS	\$	8,614.58
LZV325	FILTERS, ZERO CALIBRATION, ULTRATURB	\$	1,343.26
LZV330	FTG, O-RING, SET	\$	281.37



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LZV393	POWER SWITCH, DR3800	\$	59.17
LZV508	1 Inch Square Pour Through	\$	292.93
LZV510	QUARTZ POUR THRU CELL 450uL 1CM SQR	\$	3,006.43
LZV537	FILTER SET, DR5000 & DR2800	\$	4,740.91
LZV553	FOOT, DRB200	\$	10.41
LZV563	Stecker-Netzgerät	\$	358.46
LZV565	HALOGEN LAMP, DR 3800, DR 2800, DR 2700, LICO 500	\$	458.67
LZV566	BARCODE SCANNER - USB - DR5000	\$	2,910.07
LZV569	TUBING, POUR THRU CELL DR3900	\$	558.89
LZV579	Doppelnippel	\$	112.16
LZV582	KEYBOARD - USB,QWERTY	\$	957.82
LZV587	SET OF LENSES, DR2800	\$	381.59
LZV589	"HACH" HOUSING TOP, DR2800	\$	971.31
LZV591	FILTER WHEEL W/OPTICAL FILTERS	\$	1,653.54
LZV592	Konkavgitter	\$	3,950.76
LZV593	Tooth belt for filter changer DR3900	\$	51.65
LZV594	Tooth belt for barcode unit DR3900	\$	51.65
LZV595	Akkudeckel	\$	29.87
LZV596	HOUSING BOTTOM, DR2800	\$	414.35
LZV597	HACH NAMEPLATE	\$	31.03
LZV599	CELL COMPARTMENT, ROUND	\$	20.24
LZV600	CELLCOMPARTMENT, 50MM	\$	25.63
LZV601	TUBING, DR2800	\$	80.75
LZV603	BELT PULLEY, BARCODE UNIT	\$	40.09
LZV606	SEAL FOR DISPLAY FRAME, DR2800	\$	17.73
LZV609	DISPLAY FRAME, DR2800	\$	150.13
LZV611	TOUCHSCREEN, COMPLETE	\$	2,736.62
LZV612	MEASUREMENT SENSOR	\$	595.50
LZV613	PUSH BUTTON, DR2800	\$	314.13
LZV614	Schrittmotor-Filterwechsler	\$	319.92
LZV615	Griffmuldenabdeckung	\$	29.10
LZV616	Griffmulde "Hach"	\$	59.17
LZV618	Dichtung-Displayrahmen	\$	87.69
LZV621	BOTTOM HOUSING, DR3800	\$	1,855.89
LZV622	Displayrahmen DR 3800	\$	38.74
LZV632	USB Connection Cable	\$	169.98
LZV643	SEAL FOR DISPLAY, DR2800	\$	67.64
LZV646	Light Shield F/DR2800	\$	109.46
LZV652	CONNECTOR ASSEMBY, ULTRA TURB PLUS	\$	312.21
LZV658	Barcodeleseinheit (Redesign)	\$	2,042.83
LZV662	LEGACY ADAPTER, HQD 5PIN TO BNC	\$	699.57
LZV663	O-Ring 10,00 x 1,50 EPDM 70	\$	283.30
LZV664	Barcodeleseinheit	\$	2,216.28
LZV-702-520.110	HOSE SCREW JOINT M5 f. TUBE 6x1	\$	75.93
LZV-702-808.379	LINE FILTER, IL5XX	\$	967.45
LZV-702-823.049	DRYER MEMBRANE, PTFE, 24"	\$	6,610.30
LZV-702-886.179	Schlauchverschraubung M5	\$	39.89



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LZV-702-886.180	Schwenkverschraubung M5	\$	100.21
LZV709	FAN, DR3800	\$	487.58
LZV743	ADAPTER BOX, DR3800	\$	312.21
LZV747	LINE FILTER, DRB200	\$	258.24
LZV753	Gittermotor (Redesign)	\$	1,724.84
LZV754	Lampenstecker (Redesign)	\$	252.46
LZV755	BEAM SPLITTER MIRROR	\$	637.90
LZV756	PUSH-BUTTON	\$	312.21
LZV758	INTELLICAL PROBE STAND REPLACEMENT CLIPS	\$	198.50
LZV762	Fixed Surface of Sliding Lid	\$	82.87
LZV764	ff DISPLAY DR 5000, EMC	\$	3,044.98
LZV771	Displayrahmen DR 2700	\$	134.33
LZV774	LINK2SC SOFTWARE, DR3800	\$	1,298.93
LZV779	HOUSING BOTTOM EMC, DR5000	\$	1,462.74
LZV783	Adapter 13 mm f. Turbidimeter	\$	1,376.02
LZV784	Filter Barcodeleseinheit	\$	377.73
LZV786	FILTER KG2	\$	528.05
LZV797	Blank Module, 2100Q	\$	60.90
LZV803	aa POWER SUPPLY, 2100Q, 110-240VAC	\$	279.44
LZV804.99.00002	Module, Power/Stromversorgungsmodul	\$	472.16
LZV813.99.00002	USB+POWER MODULE, 2100Q & DR1900 Hach	\$	1,223.77
LZV818	STD USB CABLE W/MINI USB CONNECTOR	\$	112.93
LZV819	SPECIAL USB CABLE ADAPTER	\$	74.78
LZV821	RUBBER FOOT SET, 2100Q	\$	39.89
LZV822	Lamp Cover, 2100Q, ASSY	\$	60.90
LZV823	Battery Cover, 2100Q, Set	\$	60.90
LZV824	Module Cover	\$	21.97
LZV825	Connector Cover, USB+Power Module	\$	21.97
LZV826	Connector Cover, Power Module	\$	21.97
LZV827	Cap, 2100Q, ASSY	\$	40.86
LZV828	Display Set, 2100Q	\$	329.55
LZV829	Enclosure Bottom, 2100Q EPA	\$	140.11
LZV830	Enclosure Bottom, 2100Q ISO	\$	128.93
LZV831	Enclosure Top Assy, 2100Q EPA	\$	146.27
LZV832	Enclosure Top Assy, 2100Q ISO	\$	178.27
LZV833	OPTIC ASSY, 2100Q, EPA	\$	1,501.29
LZV834	Optic Assy, 2100Q ISO	\$	1,349.04
LZV835	ESD Shield, 2100Q	\$	57.05
LZV836	Service Kit, 2100Q	\$	89.23
LZV837	LED Lamp Assy, 2100Q	\$	114.28
LZV842	Wischerhalter (Seawater)	\$	2,659.54
LZV844	POWER SUPPLY, 100-240V, DR3900	\$	672.59
LZV845	DUST COVER, DR3900	\$	98.29
LZV846	ADAPTER A, DR3900	\$	175.76
LZV847	ADAPTER B, DR3900	\$	175.76
LZV848	REPLACEMENT CELL COMPARTMENT, DR3900	\$	186.17
LZV849	LIGHT SHIELD, DR3900	\$	186.17



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LZV850	STEPPING DRIVE FOR BARCODE	\$	398.93
LZV855	REAR COVER, DR3900	\$	80.75
LZV856	SAMPLE COMPARTMENT LID, DR3900	\$	45.10
LZV857	LAMP COVER, DR3900	\$	45.10
LZV859	TRAVERSE, HACH, DR3900	\$	80.75
LZV861	CELL COMPARTMENT, ROUND, DR3900	\$	120.45
LZV862	CAMERA, PREADJUSTED	\$	1,204.50
LZV863	GRATING HOLDER W/SCREWS, DR3900	\$	138.37
LZV864	DISPLAY WITH TOUCHSCREEN	\$	2,004.29
LZV865	Lichtleiter 0.Ordnung (DR 3900)	\$	80.75
LZV866	LIGHT PIPE, ROUND VIAL RECOGNITION	\$	80.75
LZV867	LIGHT PIPE, CAMERA ILLUMINATION	\$	80.75
LZV868	RECTANGULAR CELL RECOGNITION, DR3900	\$	80.75
LZV870	Kabelsatz (DR 3900)	\$	160.15
LZV871	LABEL SET, DR3900	\$	80.75
LZV873	ETHERNET CABLE, 2m	\$	90.58
LZV874	PROTECTIVE FILM FOR DR3900 DISPLAY	\$	77.09
LZV875	SIP 10 TUBING SET, LAGOPRENE + TYGON	\$	339.19
LZV876	POUR-THRU CELL, DUAL PATH, 1 in/1 cm	\$	196.57
LZV877	SIP 10 PUMP TUBING, LAGOPRENE	\$	121.41
LZV882	LUMBERG TONG	\$	167.47
LZV884	Schraubensatz (DR 3900)	\$	80.56
LZV885	DISPLAY AND TOUCHSCREEN, DR 6000 - SPARE	\$	2,370.46
LZV886	DUST COVER FOR DR 6000	\$	173.64
LZV887	ADAPTER FOR 10 CM CELL FOR DR 6000	\$	1,946.47
LZV897	Spring for glass lens DR 5000/6000	\$	63.60
LZV898	PHILLIPS COUNTERSUNK SCREW	\$	23.13
LZV899	POUR-THRU CELL KIT, DR3900	\$	967.45
LZV901	Pr ⁿ ffilter UV-Bereich (BG36)	\$	1,248.83
LZV902.99.00002	CAROUSEL ADAPTER, 1CM CELLS, FOR DR 6000	\$	1,212.21
LZV902.99.00012	CAROUSEL ADAPTER, 1 INCH CELLS, FOR DR 6000	\$	1,212.21
LZV902.99.00020	UNIVERSAL ADAPTER FOR DR 6000	\$	1,212.21
LZV904	Screw set DR 6000	\$	142.03
LZV906	Round cuvette holder	\$	237.05
LZV907	Grating with holder	\$	1,123.56
LZV908	FILTER WHEEL WITH FILTER GLASSES - SPARE	\$	643.68
LZV909	Motor, filter wheel	\$	169.59
LZV911	Detent plate of the carousel holder	\$	74.58
LZV912	Sliding Lid Cuvette Compartment	\$	215.85
LZV914	Fixed Surface of Sliding Lid	\$	101.56
LZV915	FILTER PAD FOR DR 6000	\$	38.16
LZV916	Housing upper part	\$	745.83
LZV918	Chassiseinsatz	\$	142.03
LZV919	Housing bottom DR 6000	\$	714.99
LZV920	Rear cover	\$	169.40
LZV921	Badge Logo Set	\$	194.65
LZV922	Dust cover optical assembly DR6000	\$	302.57



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LZV923	Light pipe camera illumination	\$	65.72
LZV924	Mirror Flap Motor	\$	202.36
LZV925	Temperature Sensor Lamp Cover	\$	74.58
LZV926	Fan 120mm	\$	443.26
LZV927	Fan 80mm	\$	393.15
LZV929	SAMPLE COMPARTMENT HACH - SPARE	\$	101.56
LZV930	Lamp Cover with screw	\$	101.56
LZV931	Fan cover	\$	101.56
LZV932	Front cover	\$	169.40
LZV933	CABLE GUIDE - SPARE	\$	52.61
LZV934	Cable set DR 6000	\$	635.98
LZV938	REMOTE OPERATING SOFTWARE, PHOTOMETER	\$	1,626.56
LZV940	Tubing for SIP 10	\$	42.98
LZV941	APPLICATION SOFTWARE DRINKING WATER ANALYSIS	\$	555.03
LZV942	APPLICATION SOFTWARE BREWERY ANALYSIS	\$	555.03
LZV943	APPLICATION SOFTWARE ENZYMATIC FOOD ANALYSIS	\$	555.03
LZV944	Touch Pen, Silver	\$	13.30
LZV945	Holder for Carousel Adapter	\$	148.01
LZV946	Sample Vials	\$	198.50
LZV947	Dust cover (TU5200-LAB Benchtop)	\$	102.33
LZV952	Vial Stand	\$	118.91
LZV954	Display mit Touch und Traverse	\$	1,676.66
LZV958	Display mit Touch und Traverse DR3900	\$	1,676.66
LZV967	Sipper Kit, Tubing Adapter	\$	705.36
LZV968	vv Lid Assembly for TU5200	\$	156.68
LZV969	Seal, Lid, TU5xxx	\$	19.85
LZV995	Filter Set Retrofit Kit	\$	265.95
LZW2320.99	MAGNETIC STIRRER WITH STAND, TITRALAB SO2	\$	909.64
LZW2325.99	ELECTRODE HOLDER FOR STAND, TITRALAB SO2	\$	967.45
LZW2326.99	ELECTRODE HOLDER FOR STAND, PH/TA TITRALAB	\$	967.45
LZW2328.99	MAGNETIC STIRRER FOR STAND, SO2 PH/TA &SALT TITRLB	\$	341.11
LZW2329.99	MAGNETIC STIRRER FOR STAND, KARL FISCHER TITRALAB	\$	967.45
LZW5010T.97.002	ee 5010T, pH ELECTRODE	\$	886.51
LZW5011T.97.002	5011T, pH ELECTRODE	\$	1,370.24
LZW5012.97.0002	pH electrode sensION+ 5012, glass body, liquid electrolyte	\$	551.18
LZW5012T.97.002	pH electrode sensION+ ATC 5012T , glass body, liquid electrolyte	\$	714.99
LZW5014.97.0002	5014, pH ELECTRODE	\$	967.45
LZW5014T.97.002	5014T, pH ELECTRODE	\$	1,219.92
LZW5015T.97.002	pH electrode sensION+ ATC 5015T , for chemical & pharma	\$	1,056.11
LZW5021T.97.002	5021, pH ELECTRODE	\$	1,614.99
LZW5044.97.0002	5044, REFERENCE ELECTRODE	\$	1,219.92
LZW5045.97.0002	5045, MULTI ELECTRODE, pH, ORP	\$	1,589.94
LZW5046.97.0002	5046, ORP, SILVER ELECTRODE, COMBINATION, TITRALAB	\$	1,441.55
LZW5048.97.0002	5048, MULTI ELECTRODE, pH, ORP, COND	\$	1,985.02
LZW5049.97.0002	5049, MULTI ELECTRODE, PH, COND, DO	\$	4,066.39
LZW5050T.97.002	ee 5050T, pH ELECTRODE	\$	932.76
LZW5051T.97.002	5051T, pH ELECTRODE	\$	1,489.73

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LZW5052T.97.002	5052T, pH ELECTRODE	\$	1,456.96
LZW5053T.97.002	pH field electrode sensION+ 5053T with ATC + stainless-steel sheath	\$	1,131.27
LZW5055.97.0002	5055, REDOX ELECTRODE	\$	967.45
LZW5056.97.0002	5056, REDOX ELECTRODE	\$	1,312.42
LZW5057.97.0002	5057, REDOX ELECTRODE	\$	1,347.11
LZW5059.97.0002	5059 PT, MULTI ELECTRODE, pH, COND	\$	1,499.36
LZW5060.97.0002	5060, COND ELECTRODE	\$	1,347.11
LZW5062.97.0002	5062, COND ELECTRODE	\$	1,426.13
LZW5064.97.0002	Conductivity cell sensION+ ATC 5064, field use, Ti cell, ABS body	\$	695.72
LZW5070.97.0002	5070, COND ELECTRODE	\$	1,381.80
LZW5071.97.0002	Conductivity cell sensION+ ATC 5071, for lab, Pt cell, glass body	\$	1,283.52
LZW5080.97.0002	5080 Double platinum electrode (Pt)	\$	1,805.79
LZW5123.99	PROTECTOR FLASK 5120 DO PROBE, LIGHT	\$	62.06
LZW5125.99	SPARE MEMBRANE, 5120 DO PROBE	\$	470.24
LZW5130.97.0002	5130 PORTABLE DO PROBE, MP5	\$	1,753.75
LZW5131.97.0002	5131 BENCHTOP DO PROBE, BNC	\$	1,626.56
LZW5202.97.0002	pH electrode, for general use, large membrane	\$	803.64
LZW5208.97.0002	pH electrode, for micro samples, ° 3 mm	\$	1,017.56
LZW5209.97.0002	pH electrode, for micro samples, ° 6 mm	\$	961.67
LZW5211.97.0002	pH electrode, food, 3 diaphragms, pH 0...12, 0...100 C.	\$	1,231.48
LZW5221.97.0002	pH electrode, for difficult media, pH 0...12, 0...60 C	\$	1,243.04
LZW5231.97.0002	pH electrode, puncture, stainless-steel, solid electrolyte	\$	888.44
LZW5233.97.0002	pH electrode, puncture, POM, solid electrolyte, pH 2...14, 0...60 C	\$	788.22
LZW5243.99	CERAMIC DIAPHRAGM, LZW5044 REF ELECTRODE	\$	139.72
LZW5264.97.0002	5264, DOUBLE PLATINUM ELECTRODE, TITRALAB	\$	1,676.66
LZW5292.97.0002	Conductivity cell, PT cell, glass body, for general use	\$	903.86
LZW5293.97.0002	Conductivity cell, PT cell, glass body, banana + telephonic connector	\$	977.09
LZW8204.99	THERMAL PRINTER FOR SENSION+ INSTRUMENTS	\$	1,191.01
LZW8205.99	CABLE, THERMAL PRINTER TO SENSION+	\$	102.91
LZW8206.99	Thermal Paper for KyoMouse RS	\$	179.42
LZW8705.99	INLET TUBE, BURETTE PUMP TITRALAB	\$	208.14
LZW8710.99	TUBE, EMPTIES/DISPENSES PUMP, CONIC, TITRALAB	\$	291.01
LZW8741.99	CYLINDER, POSITIONING, 5 ML SYRINGE, TITRALAB	\$	57.62
LZW8805.99	REAR PANEL HOUSING, COMPLETE, PH/TA TITRALAB	\$	184.82
LZW8806.99	REAR PANEL HOUSING, COMPLETE, SO2 TITRALAB	\$	184.82
LZW8807.99	REAR PANEL HOUSING, COMPLETE, SALT TITRALAB	\$	184.82
LZW8808.99	MOTOR, PERISTALTIC PUMP, 12V, 170ML/MIN	\$	1,113.92
LZW8810.99	FRONT ENCLOSURE, ONE BURETTE, TITRALAB	\$	749.68
LZW8811.99	Lexan/Keyboard TitraLab SO2	\$	111.39
LZW8812.99	KEYBOARD, TITRALAB PH/TA	\$	111.39
LZW8813.99	KEYBOARD, TITRALAB SALT	\$	111.39
LZW8815.99	REAR PANEL W/ CIRCUIT BDS, TITRALAB KARL FISCHER	\$	3,334.06
LZW8817.99	REAR PANEL W/ CIRCUIT BOARDS, TITRALAB PH/TA	\$	3,507.50
LZW8818.99	REAR PANEL W/ CIRCUIT BOARDS, TITRALAB SALT	\$	3,507.50



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LZW8819.99	KIT, SCREW SET, PISTON SYRINGE CONN, TITRALAB	\$	31.03
LZW8820.99	BURETTE SHAFT OPTICAL POSITN DETECTN ASSY, TITRALAB	\$	148.59
LZW8821.99	BURETTE MECHANICAL ASSEMBLY, TITRALAB	\$	6,976.46
LZW8822.99	BURETTE MOTOR ASSY W/PULLEY & BELT, TITRALAB	\$	1,113.92
LZW8823.99	RING, TUBING HOLDER FOR TITRALAB	\$	148.59
LZW8824.99	KF SMPL CELL HEAD INCL STOPPERS & O-RING, TITRALAB	\$	1,624.63
LZW8826.99	KIT, SCREWS & SML PARTS FOR FRNT ENCL, TITRALAB	\$	42.21
LZW8827.99	KIT, SCREWS & SML PARTS FOR 2 PUMP MOTORS, TITRALB	\$	74.78
LZW8828.99	REAR PANEL HSG, COMPLETE, KARL FISCHER	\$	184.82
LZW9008.99	POWER SUPPLY, SENSION+ BT, 230-115VAC, 12V	\$	283.30
LZW9055.99	AS7 CABLE, 1M, BNC, TITRALAB	\$	240.90
LZW9110.97	FLASKS FOR BENCHTOP pH CAL, 3X50 ML	\$	35.65
LZW9111.99	FLASKS FOR BENCHTOP COND CAL, 3X50	\$	35.65
LZW9118.99	Pyrex glass chamber, sensION+	\$	537.69
LZW9135.99	RS232 cable for sensION+	\$	265.95
LZW9135USB.99	RS232 cable for sensION+	\$	531.91
LZW9136.99	TUBE REDOX CAL, 1 X 10 ml	\$	17.73
LZW9137.97	TUBES FOR PORTABLE pH CAL, 3X 10 ML	\$	35.65
LZW9138.99	TUBES FOR PORTABLE COND CAL, 3X10 ML	\$	35.65
LZW9140.99	TUBE PORTABLE pH 10.01 CAL 1 x 10ML	\$	17.73
LZW9147.99	TUBE PORTABLE PH 7.00, CAL 1X10	\$	17.73
LZW9155.99	Holder and clamp for three sensors	\$	443.26
LZW9161.99	PP protector, electrode storage	\$	16.19
LZW9162.99	PP protector, for measurement	\$	25.05
LZW9169.99	Replacement protective shroud for sensION+ 5049 multi probe	\$	219.70
LZW9171.99	REPLACEMENT PH ELECTRODE, 5049 MULTISENSOR	\$	609.00
LZW9172.99	REPLACEMENT COND CELL, 5049 MULTISENSOR	\$	1,127.41
LZW9173.99	REPLACEMENT DO SENSOR, 5049 MULTISENSOR	\$	1,219.92
LZW9179.99	REPLACEMENT HANDLE, 5049 MULTISENSOR	\$	1,352.89
LZW9200.99	DESICCANT CARTRIDGE, MOLECULAR SIEVE, TITRALAB	\$	67.26
LZW9267.99	Holder and clamp for one sensor	\$	447.11
LZW9308.99	Plug for electrolyte refill hole,	\$	17.73
LZW9314.99	90 ML FLASK, SensION	\$	17.15
LZW9321.99	SENSOR HOLDER FOR 3, senSION+ BENCH	\$	84.03
LZW9325.99	ELECTRODE HOLDER, SENSION+, RA ELECT.	\$	148.01
LZW9400.99	STANDARD REDOX SOLN, 220 mV, 250 ML	\$	94.48
LZW9401.99	STANDARD REDOX SOLN, 220 mV, 1 L	\$	281.85
LZW9402.99	STANDARD REDOX SOLN, 220 mV, 125 ML	\$	65.57
LZW9460.99	pH STD SOLN, 4.01, 125 ML	\$	48.42
LZW9461.97	pH STD SOLN, 7.00, 125 ML	\$	48.42
LZW9461.97-LM	pH STD SOLN, 7.00, 125 ML	\$	48.42
LZW9462.99	pH STD SOLN, 9.21, W/CERT, 135ML	\$	35.05
LZW9463.99	pH STD SOLN, 4.01 W/CERT, 250ML	\$	59.98
LZW9464.97	pH STD SOLN, 7.00, W/CERT, 250ML	\$	59.98
LZW9464.97-LM	pH STD SOLN, 7.00, W/CERT, 250ML	\$	59.98
LZW9465.99	pH STD SOLN, 9.21 W/CERT, 250 ml	\$	52.40
LZW9466.99	pH STD SOLN, 4.01, W/CERT, 1L	\$	114.55

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LZW9467.97	pH STD SOLN, 7.00, W/ CERT, 1 L	\$	114.55
LZW9467.97-LM	pH STD SOLN, 7.00, W/ CERT, 1 L	\$	114.55
LZW9470.99	pH STD SOLN., 10.01, 125 ML	\$	48.42
LZW9471.99	pH STD SOLN., 10.00 W/CERT, 250ML	\$	59.98
LZW9472.99	pH STD SOLN., 10.00 W/CERT, 1 L	\$	114.55
LZW9500.99	KCl 3M SOLN, 250ML	\$	84.72
LZW9509.99	KCl 3M SOLN, 50ML	\$	35.58
LZW9509.99-LM	KCl 3M SOLN, 50ML	\$	35.58
LZW9510.99	KCl 3M SOLN, 125ML	\$	54.20
LZW9521.99	ELECTROLYTE SOLN, KNO3 250 ML	\$	101.71
LZW9650.97.0002	9650, SODIUM ISE	\$	1,786.51
LZW9652.97.0002	9652, CHLORIDE ISE	\$	2,351.18
LZW9652AP.97.02	KTO: SENSION+ CHLORIDE ANALYSIS PKG	\$	3,025.70
LZW9652C.97.002	9652C, CHLORIDE ISE, COMBINATION	\$	3,179.88
LZW9655.97.0002	9655, FLUORIDE ISE	\$	2,351.18
LZW9655AP.97.02	KTO: SENSION+ FLUORIDE ANALYSIS PKG	\$	3,141.34
LZW9655C.97.002	9655C, FLUORIDE ISE, COMBINATION	\$	3,179.88
LZW9660C.97.002	SENSION+ 9660C Combination Calcium	\$	2,832.98
LZW9662.97.0002	9662, NITRATE ISE	\$	2,351.18
LZW9662C.97.002	9662C, NITRATE ISE, COMBINATION	\$	3,179.88
LZW9663C.97.002	9663C, AMMONIUM ISE, COMBINATION	\$	3,044.98
LZW9700.99	COND STD SOLN, 147uS/cm, CERT, 250ML	\$	72.27
LZW9701.99	COND STD SOLN, 147uS/cm, 125ML	\$	54.20
LZW9701.99-LM	COND STD SOLN, 147uS/cm, 125ML	\$	54.20
LZW9710.99	COND STD SOLN, 1413 uS/cm, CERT, 250ML	\$	73.52
LZW9710.99-LM	Conductivity Standard Solution 1413 uS/cm, 250ml	\$	73.52
LZW9711.99	COND STD SOLN, 1413uS/cm, 125ML	\$	54.20
LZW9711.99-LM	COND STD SOLN, 1413uS/cm, 125ML	\$	54.20
LZW9720.99	COND STD SOLN, 12.88mS/cm, CERT, 250ML	\$	72.27
LZW9721.99	COND STD SOLN, 12.88 mS/cm, 125ML	\$	54.20
LZW9721.99-LM	COND STD SOLN, 12.88 mS/cm, 125ML	\$	54.20
LZW9811.99	ELECTROLYTIC SOLN, sensION+ DO, 25ML	\$	151.41
LZW9901.99	ELECTROLYTIC SOLN, sensION+ ISE, 125ML	\$	59.98
LZX012	WIPERS, REPLACEMENT, UV NIT, ORG2MM	\$	1,063.81
LZX017	FILTER ELEMENTS, PK/8 (FILTRAX)	\$	138.37
LZX018	ANNUAL MAINTENANCE KIT (FILTRAX)	\$	4,297.66
LZX019	PUMP ROLLERS, SET OF 5, FILTRAX	\$	233.19
LZX024	COMPRESSOR, 230V FILTRAX	\$	2,062.10
LZX025	COMPRESSOR, 115V, FILTRAX	\$	2,139.19
LZX026	PHOTOMETER REPLACEMENT FOR FILTRAX	\$	1,640.05
LZX027	FAN, COMPLETE (FILTRAX)	\$	896.15
LZX029	CONNECTOR SET	\$	728.48
LZX030	AIR FILTER FOR COMPRESSOR (FILTRAX)	\$	102.14
LZX050	WIPERS, REPLACEMENT FOR TURBIDITY	\$	265.95
LZX056	OPTICS HOLDER, PLASTIC	\$	5,877.96
LZX059	Wischerwelle	\$	755.46
LZX077	GEARED MOTOR (SOLITAX, SONATAX)	\$	795.93



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LZX107	SAMPLE HOSE, 5M HEATED 115 VAC, FILTRAX	\$	6,417.58
LZX117	WIPERS, REPLACEMENT, UV NIT,ORG 5MM	\$	616.70
LZX119	WIPERS, REPLACEMENT, UV,ORG 50MM	\$	506.85
LZX148	SPARE WIPERS FOR 1MM NITRATAX	\$	2,274.10
LZX182	GEARED MOTOR (UVAS, NITRATAX)	\$	1,778.81
LZX216	CLEANING VESSEL, FILTRAX FILTER MOD	\$	699.57
LZX217	KIT, FILTRAX CLEANING	\$	173.83
LZX223	WIPER COMPLETE	\$	487.58
LZX278	DISCHARGE HOSE, 2M (AMTAX,PHOSPHAX)	\$	68.61
LZX302	DESICCANT, 20X40MM PROBE, SET OF 2	\$	819.06
LZX303	DESICCANT 1, ULTRA TURB PLUS	\$	20.24
LZX304	DESICCANT 2, ULTRA TURB PLUS	\$	25.05
LZX328	WIPERS, REPLACEMENT FOR SLUDGE PROBE	\$	456.75
LZX337	BALL VALVE SAFETY ARM. FOR SOLITAX IN- & HIGHLINE	\$	11,890.82
LZX341	MODULE BODY, UVAS PLUS 50MM LXG170	\$	23,338.39
LZX354	DISPLAY	\$	1,630.41
LZX372	MODULE BODY, UV NITRATE PROBE LXV245	\$	36,173.54
LZX396	TEST VERIFICATION FILTER, UV ORGA	\$	568.52
LZX407	TUBING, CLEAR, UV NITRATE BIPASS PNL	\$	56.27
LZX414.00.60000	KIT, MTG HRDWARE,RAIL MT,FILTERPROBE	\$	3,141.34
LZX414.00.73000	RAIL MOUNTING ASSEMBLY, SONATAX	\$	2,910.07
LZX416	MISC HARDWARE F/PROBE INST KIT OPTIQ	\$	358.46
LZX417	KIT, SEAL FOR SEALED PROBE ADAPTERS	\$	260.17
LZX419	MODULE BODY, NITRATAX PLUS 1MM LXG	\$	17,094.26
LZX421	SET OF SEALS FOR SOLITAX PROBES	\$	116.02
LZX423	OPTICS HOLDERS SET F/SOLITAX-FILTERS	\$	1,044.54
LZX426	SET OF SEALS	\$	450.96
LZX428	2 O-RINGS FOR FLOW-THROUGH LZX867 LZX868 LZX869	\$	17.73
LZX450	SEDIMENTER	\$	8,113.51
LZX453	Blitzlampe	\$	3,546.05
LZX455	OPTICS ASSEMBLY, SOLITAX SC	\$	5,955.05
LZX461	Einbauarmatur	\$	1,902.15
LZX509	MODULE BODY, UVAS PLUS 1 MM LXG211	\$	25,188.50
LZX510	MODULE BODY, NITRATAX PLUS 5 MM	\$	23,338.39
LZX567	Set of seals for clear- probes	\$	194.65
LZX572	SEALS FOR NITRATAX/UVAS CLEAR FLOW	\$	108.31
LZX578	SET OF WIPER BLADES IN VITON MAT'L	\$	738.12
LZX657	INTERFERENCE FILTER SET, OPTIQUANT N	\$	3,661.68
LZX658	INTERFERENCE FILTER SET-UV PROBE	\$	3,507.50
LZX660	WELDING FLANGE IN-LINE PROBES SS 1.4571	\$	686.08
LZX661	WELDING FLANGE,IN-LINE PROBES,C-STEEL, C22.8, 1.0460	\$	350.75
LZX667	MAINTENANCE KIT, FILTRAX	\$	1,667.03
LZX671	FILTRAX HEATED HOSE, 10M, 115V	\$	2,312.64
LZX672	FILTRAX HEATED HOSE, 10M, 230V	\$	2,312.64
LZX673	FILTRAX HEATED HOSE, 20M, 115V	\$	4,644.55
LZX674	FILTRAX HEATED HOSE, 20M, 230V	\$	4,644.55
LZX675	FILTRAX DELIVERY HOSE, 2M, UNHEATED	\$	643.68

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LZX676	HOLDER, CONTROL UNIT (FILTRAX)	\$	1,528.27
LZX677	FILTER MODULE, FILTRAX	\$	1,753.75
LZX701	TUBE ADAPTER KIT (AMTX,PHOSX,FILTX)	\$	429.77
LZX736	WIPER HOLDER, UVAS 5MM	\$	1,285.44
LZX737	Wischerarm 50 mm UVASplus-Sonde	\$	1,605.36
LZX738	Wischerarm 50 mm UVASplus bypass	\$	1,927.20
LZX739	FILTRAX, PUMP HEAD, 2 CHAN, 5 ROOLS	\$	1,086.94
LZX741	DISPLAY MODULE, FILTRAX	\$	770.88
LZX763	Kleinteileset	\$	797.86
LZX768	Schaltschrankheizung	\$	1,144.76
LZX769	Adapter SS for LDO; 1 7/8-12UN	\$	256.32
LZX771	Set 2 -Schrauben	\$	129.70
LZX774	Kabelverschraubung (PVDF)	\$	146.27
LZX823	REPLACEMENT GATEWAY CIRCUIT BOARD	\$	1,778.81
LZX847	SENSOR CABLE, 0.35 METERS, ULTRA TURB	\$	379.66
LZX848	EXTENSION CABLE, 5 METERS, SC-SENSORS	\$	616.70
LZX849	EXTENSION CABLE, 10 METERS, SC-SENSORS	\$	533.83
LZX850	EXTENSION CABLE, 15 METERS, SC-SENSORS	\$	736.19
LZX851	EXTENSION CABLE, 20 METERS, SC-SENSORS	\$	925.06
LZX852	EXTENSION CABLE, 30 METERS, SC-SENSORS	\$	1,129.34
LZX853	EXTENSION CABLE, 50 METERS, SC-SENSORS	\$	1,530.20
LZX858	Messfenstersatz ohne Gewinde	\$	1,927.20
LZX866	db FLOWTHRU CELL NITRATAX CLEAR SC 5MM	\$	9,057.84
LZX867	db FLOWTHRU CELL NT3x00sc 5MM, NITRATAX PLUS SC 5MM, UVAS PLUS SC 5MM	\$	9,057.84
LZX868	db FLOWTHRU CELL UVAS PLUS SC 50MM	\$	8,537.50
LZX869	db FLOWTHRU CELL NT3x00sc 1/2MM, NITRATAX PLUS SC 2MM, UVAS PLUS SC 2MM	\$	9,057.84
LZX885	pH-Elektrode	\$	803.64
LZX887	Interfacekabel RS232 f. sc100/sc200	\$	965.53
LZX888	ORPpH combination electrodefor LXV426	\$	1,233.41
LZX889	ee pH ELECTRODE FOR 1200SSC	\$	1,092.72
LZX890	Replacement ORP electrode	\$	1,244.97
LZX893	Dig. Cable assembly, cable and tube adaptor, 1200Ssc + 3798Ssc	\$	776.66
LZX894	Kleinteilesatz f'r 1200Ssc+3798Ssc	\$	143.96
LZX895	Electronic board for 3798 S sc	\$	1,162.10
LZX896	Replacement sensor for 3798 S sc	\$	1,588.01
LZX897	Replacement pH electrode cablewith AS/ connector, 1200 S sc	\$	339.19
LZX898	Electronic board for 1200 S sc	\$	1,171.74
LZX899	Electrode protection for 1200 S sc	\$	126.42
LZX902	OPTICS CARRIER	\$	2,274.10
LZX903	PROBE HOUSING, WELDS	\$	1,343.26
LZX904	SondengehZuse, Edelstahl	\$	1,580.30
LZX905	HOUSING, SC SOLITAX	\$	1,081.16
LZX908	PLUG, SEALING	\$	14.65
LZX914.99.11110	CHAIN MOUNT KIT, FP360sc, 5M	\$	1,443.47
LZX914.99.11200	CHAIN MOUNT KIT, STAINLESS STEEL, LDO2	\$	2,100.65

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LZX914.99.11300	Kettenhalterung - Edelstahl, 5m	\$	2,312.64
LZX914.99.12400	KIT, CHAIN MOUNT PVC, ISE PROBES	\$	1,570.67
LZX915	MODULE, SC1000 BASE DIN RAIL	\$	738.12
LZX918	CONNECTOR, INTERNAL NTWRK FOR SC1000	\$	204.28
LZX919	Externes Ausgangs-Modul analog	\$	1,104.29
LZX920	MODULE, RELAY (4) DIN RAIL, SC1000	\$	1,618.85
LZX921	MODULE, DIN RAIL SC1000 mA INPUTS-2	\$	1,327.84
LZX923	SCREWS (4PC) F/FLOW THROUGH CELL	\$	53.58
LZX925	HOUSING, FRONT SC1000 DISPLAY (HACH)	\$	1,900.22
LZX927	DISPLAY MODULE, BACKHOUSING, SC1000 DISPLAY	\$	664.88
LZX928	FRAME, INTERNAL SC1000 DISPLAY MOD	\$	106.77
LZX929	GASKET, SC1000 DISPLAY MODULE	\$	35.85
LZX930	Set of Screws SC1000 Display Module	\$	59.55
LZX931	CONNECTOR, SC1000 ANTENNA	\$	237.05
LZX932	STAIN RELIEF, RELAYS, SC1000	\$	53.38
LZX933	Verbindersatz (Connector Set)	\$	119.10
LZX934	Kabel	\$	703.43
LZX935	STRAP, CARRYING, SC1000 DISP MODULE	\$	33.53
LZX936	BALL VALVE ARMATURE FOR SOLITAX INLINE & HIGHLINE	\$	5,627.42
LZX937	SPRINGS, CONTACT SC1000	\$	82.87
LZX938	COVER, SIM CARD, SC1000	\$	87.69
LZX948	KIT, SMALL PARTS, SC1000 SUNSHIELD	\$	62.63
LZX949	HOUSING, FRONT SC1000 PROBE MD HACH	\$	2,062.10
LZX953	HOUSING, BACK SC1000 PROBE MD HACH	\$	1,277.73
LZX954	GASKET, SC1000 PROBE MODULE	\$	176.53
LZX955	CABLE, SC1000 ANTENNA EXT 10M	\$	439.40
LZX956	GSM Quadband Antenne	\$	302.57
LZX957	HARDWARE, MTG. CONTROLLER W SUNROOF	\$	1,877.09
LZX958	SUNROOF, CONTROLLER	\$	612.85
LZX961	SUNSHIELD SC200	\$	639.83
LZX961.54	SC100 PLASTIC SUNSHIELD	\$	597.43
LZX962	FAN, SC1000 PROBE MODULE	\$	273.66
LZX963	CARD, SC1000 MULTIMEDIA (16MB)	\$	350.75
LZX967	KIT, CONNECTOR, SC1000 MODULE	\$	429.77
LZX968	COVER, INTERNAL RELAY, SC1000	\$	235.12
LZX969	KIT, SC SENSOR CONNECTOR (2) SC1000	\$	240.90
LZX970	CONNECTORS (2), POWER OUTLET, SC1000	\$	304.50
LZX971	PLUG SC-PROBE CABLE, UVAS	\$	206.21
LZX972	SLEEVE, SC-SENSOR CABLE	\$	310.28
LZX973	SCREWS (4) EXTERNAL, SC1000 PROBE MD	\$	41.43
LZX974	SCREWS (6), INTERNAL SC1000 PROBE MD	\$	20.24
LZX975	SCREW, SET, SC1000	\$	20.24
LZX976	Replacement fuses for SC1000 (2 large, 2 small fuses)	\$	65.14
LZX977	COVER, INTERNAL SC1000 NTWRK CONNEC	\$	87.69
LZX978	STRAIN RELIEFS (2), SC1000	\$	77.09
LZX979	Satz Blindkappen	\$	170.94
LZX980	STRAIN RELIEF, SC1000 POWER SUPPLY	\$	34.30

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LZX981	CONNECTOR, CONDUIT, SC1000	\$ 77.09
LZX982	COVER, SC1000 POWER OUTPUT	\$ 47.02
LZX983	COVER, INTERNAL SC1000 POWER SUPPLY	\$ 123.73
LZX985	Calibration tool for 3798-S sc	\$ 316.06
LZX986	Buchse sc-Sondenkabel	\$ 167.67
LZX987	ET FEHL TR CARD SC ECO/CLEAR	\$ 187.90
LZX995	Luftfilter, rund	\$ 141.84
LZX998	CABLE, ETHERNET CROSSOVER	\$ 79.59
LZY001	BRACKET, SUNSHIELD, SC1000	\$ 191.37
LZY007	PLUG SET, RUBBER, F. SC-ANALYZER	\$ 39.89
LZY008	Trigger-Fassung	\$ 3,218.42
LZY044	KIT, MOUNTING, SC-ANALYZER	\$ 360.39
LZY045	Schraubkappensatz - AMTAX sc	\$ 948.18
LZY046	Schraubkappensatz - PHOSPHAX sc	\$ 873.02
LZY048	SCREW CAP, COMPLETE FOR 1.5" PVC PIPE	\$ 294.86
LZY051	9180sc ACIDIFICATION UNIT	\$ 10,098.53
LZY053	FLOW CELL, AMPEROMETRIC SC ANALYZER	\$ 5,357.62
LZY054	FLOW CELL MAINTENANCE KIT FOR 9184sc	\$ 481.80
LZY059	PANEL, AMPEROMETRIC SC ANALYZER	\$ 1,252.68
LZY061	AMPEROMETRIC ELECTRODE PARTS KIT	\$ 136.25
LZY069	NH4-N ELECTRODE, AMTAX ANALYZER	\$ 3,256.97
LZY070	ELECTRODE, SPECIAL, AMTAX SC	\$ 3,449.69
LZY080	Stainless steel flowcell for PG 13.5 120 mm long sensors	\$ 3,758.04
LZY105	CABLE TO TRANSMITTER, DIGITAL AMP	\$ 1,059.96
LZY111	FITTING SET, 3.2MM	\$ 168.05
LZY119	SEALING FOR TOP COVER	\$ 28.14
LZY123	VALVE EXHAUST, 230V	\$ 1,462.74
LZY124	VALVE EXHAUST, 115V	\$ 1,462.74
LZY125	SAMPLE PUMP, 230V	\$ 4,606.01
LZY126	SAMPLE PUMP, 115V	\$ 4,606.01
LZY127	SLEEVE NUT, WITH SEALING	\$ 283.30
LZY130	WEAR PARTS, SET FOR SAMPLE PUMP	\$ 512.64
LZY135	ELBOW FITTING 4/6MM	\$ 106.19
LZY136	FITTING FOR EXHAUST VALVE	\$ 106.19
LZY137	HUMIDITY SENSOR W/CONNECTOR	\$ 35.65
LZY138	EXHAUST, FOR AIR CLEANING	\$ 141.07
LZY139	EXHAUST, FILTRATION PROBE, COPPER	\$ 106.19
LZY140	FILTER MODULE, FILTERPROBE sc	\$ 1,433.84
LZY141	ADAPTER, FOR FILTRATION MODULE	\$ 246.68
LZY142	SEAL, FOR FILTER MODULE	\$ 16.00
LZY143	DOOR, FOR AMTAX SC OUTDOOR ENCLOSURE	\$ 2,254.82
LZY145	ENCLOSURE, SC ANALYZER	\$ 5,955.05
LZY146	Auffangwanne sc-Analyser	\$ 522.27
LZY147	CLAMP LOCK, SC-ANALYZER	\$ 246.68
LZY149	COMPRESSOR, SWITCHABLE 115/230V	\$ 2,177.74
LZY151	AIR TUBING, FOR COMPRESSOR	\$ 114.86
LZY152	FAN, INCOMING AIR	\$ 805.57



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LZY153	FAN, AIR CIRCULATION	\$	316.06
LZY154	FILTER SET, WITH FILTER PADS	\$	28.14
LZY156	HEATNIG FOR ANALYZER ENCLOSURE	\$	882.66
LZY157	Gehäusegitter rechts	\$	38.16
LZY160	VALVEBLOCK, FOR PHOSPHAX SC MR:1-50mg/l	\$	9,867.26
LZY165	BOTTOM, OVERFLOW VESSEL	\$	316.06
LZY166	OVERFLOW VESSEL, UPPER	\$	238.97
LZY167	OVERFLOW VESSEL UPPER W/VALVE	\$	1,509.00
LZY168	Ventil 2/2 Wege einzeln	\$	1,063.81
LZY169	VALVE BLOCK, AMTAX SC, ALL RANGES	\$	9,867.26
LZY170	CONVERSION KIT, FORM 1-CHANNEL	\$	4,606.01
LZY171	VALVE, 3/2 WAY	\$	1,125.48
LZY172	Ventilblock Umschaltung 2 Kanal	\$	909.64
LZY173	VALVE BLOCK, AMTAX SC	\$	1,202.57
LZY174	VALVE BLOCK, UPPER, AMTAX ANALYZER	\$	225.48
LZY175	VALVE BLOCK, UPPER WITH VALVE,AMTAX	\$	1,297.01
LZY177	PISTON PUMP, SC ANALYZER	\$	6,764.47
LZY183	VALVE BLOCK, PHOSPHAX SC MR: 0.05-15mg/l	\$	9,867.26
LZY187	GASKET, DOOR, SC-ANALYZER	\$	168.05
LZY188	LOCK, DOOR, SC-ANALYZER	\$	56.27
LZY189	KIT, ACCESSORIES SC ANALYZER	\$	298.72
LZY198	SEALING FOR OVERFLOW VESSEL	\$	10.41
LZY200	TOOL FOR ELECTRODE, SPECIAL	\$	49.14
LZY201	CUTTER, FOR TUBING	\$	100.21
LZY203	MASTER LINK SET FOR SOLITAX sc	\$	122.18
LZY218	#Kabelverschraubung PG 9	\$	333.41
LZY220	SCREW SET, RAIL MOUNTING	\$	143.77
LZY221	Signalleuchte LED	\$	100.41
LZY222	Signalleuchte LED	\$	100.41
LZY223	SCREW SET, SC ANALYZER	\$	15.03
LZY224	COVER, INSULATION, FOR GSE AMTAX SC	\$	265.95
LZY225	COVER, INSULATION FOR PHOTMETER	\$	100.21
LZY226	COVER, INSULATION FOR PHOTMETER	\$	106.38
LZY228	RING W/O THREADS	\$	591.65
LZY229	RING WITH THREADS FOR FLOW THROUGH	\$	616.70
LZY230	SEALING PLUG, NITRATAX FLOW CELL	\$	17.34
LZY232	Extender Chain with master link,SS (5m)	\$	331.48
LZY241	CONVERSION KIT, FORM FILTER PROBE	\$	1,761.46
LZY242	CONVERSION KIT, FORM 1-CHANNEL	\$	2,871.53
LZY261	SCREW SET, PROBE ADAPTER, SOLITAX sc MOUNTING	\$	56.27
LZY267	VALVE BLOCK, 2 CHANNEL SWITCH	\$	2,042.83
LZY268	UPPER PART OF OVERFLOW VESSEL, A/Psc	\$	254.39
LZY269	OVERFLOW VESSEL FOR 2-CHANNEL INSTR.	\$	2,563.18
LZY270	COVER PLATE FOR SC ANALYZER	\$	163.04
LZY271	VALVE BLOCK FOR PHOSPHAX	\$	1,202.57
LZY272	VALVE BLOCK FOR PHOSPHAX SC	\$	1,202.57
LZY285	KIT, MTG HRDWARE,RAIL MT,AMTAX sc W/	\$	3,102.79



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LZY286	KIT, MTG HRDWARE,STAND MT,AMTAX sc	\$	3,218.42
LZY287	KIT, MTG HRDWARE,STAND MT,AMTAX sc	\$	2,948.62
LZY290	Schraubkappe mit 1/2"se für Armatur 1"	\$	262.10
LZY300	ACIDIFICATION TUBING REPLACEMENT KIT	\$	1,412.64
LZY302	CONNECTING HOSE, HEATED DRAIN	\$	1,668.96
LZY303	Beheizter Ablauf-/Verbindungs-	\$	3,141.34
LZY307	SET OF 3 MEMBRANE CAPS USED W/ AMTAX	\$	237.05
LZY316	KIT, MTG HRDWARE,RAIL MT, AMTAX sc	\$	2,948.62
LZY325	POWER SUPPLY FOR FILTRAX	\$	2,023.56
LZY330	ELECTRODE KEY, AMTAXsc	\$	79.21
LZY332	AIRFILTER FOR AMTAX sc COMPRESSOR	\$	65.91
LZY334	SPACE BOLT, L=15MM, M3	\$	11.18
LZY335	STANDOFF L 18 MM, LXG424	\$	35.85
LZY336	SPACER BOLT L=10MM, M3	\$	47.79
LZY344	Wischerarm	\$	637.90
LZY345	GUIDANCE SCREW (TO ARRET WIPER ARM)	\$	69.76
LZY346	HOUSING, BTM W/TRANSDUCER & PT100	\$	2,158.46
LZY347	HOUSING TOP, SONATAX	\$	1,655.46
LZY348	SET SPACER FOR PCBA, SONATAX	\$	56.27
LZY349	O-Ring 146 x 4,5 NBR	\$	52.03
LZY350	Gasket	\$	47.79
LZY351	SEALING KIT FOR CABLE GLAND, SONATAX	\$	52.03
LZY362	KIT, SCREW, M5 WITH ALLEN KEY	\$	31.41
LZY367	CLEANING SET, SC ANALYZER	\$	840.26
LZY376	SPECIAL VASELINE, ULTRA TURB PLUS	\$	116.60
LZY379	Druckring 1	\$	501.07
LZY403	Verbinder	\$	258.24
LZY410	LIGHT DUCT, SONATAX	\$	182.31
LZY411	SEALING INSERT, CABLE GLAND, SONATAX	\$	41.43
LZY418	WINDOW FOR MEASURING CHAMBER PHOSPHAX sc	\$	215.85
LZY421	10cm cell holder, DR 6000	\$	1,304.71
LZY429	HRDW SET, 4 SCKT-HD SCREWS W WSHRS	\$	79.59
LZY431	EXTENSION CABLE, 5M SC POWER	\$	714.99
LZY440	KEYS FOR ENCLOSURE SC ANALYZER	\$	59.74
LZY441	HAND WHEEL FOR IN LINE ARMATURE	\$	177.88
LZY449	BRACKET FOR PIPES	\$	437.47
LZY457	MIDDLE PART PROBE ENCLOSURE	\$	2,447.54
LZY463	ENCLOSURE SCREWS, NITRATAX AND UVAS	\$	15.61
LZY465	Ammonium Analyzer AMTAX sc System Kit (2nd year of Operation)	\$	3,758.04
LZY467	Phosphax sc System Kit (2nd year of Operation)	\$	886.51
LZY472-01	SC4x00 USB to Ethernet Adapter	\$	217.77
LZY474	LED, BLOCK PHOSPHAX SC 0.05 - 15 mg/L	\$	1,676.66
LZY488	sc1000 bus cable-flexible (per foot)	\$	8.29
LZY489	sc1000 bus cable-fix/outdoor (per ft)	\$	16.00
LZY513	COUPLING WITH POSITION MARK	\$	112.55
LZY520	SD CARD, SC1000 Display Module	\$	165.93



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LZY522	MOBILEMATE CARD READER, SD/MMC	\$	161.11
LZY525	HOSE CLAMP, FILTERPROBE	\$	192.72
LZY540	MOTOR COUPLING COMPLETE FOR WIPER	\$	431.69
LZY542	SET OF SCREWS M3X10	\$	16.96
LZY547	SENSOR HEAD NITRATA/UVAS SC, 2MM	\$	734.26
LZY548	FASTENING PARTS FOR PROTECTIONPLATE	\$	198.50
LZY553	SC1000 ETHERNET OUTDOOR CONNECT KIT	\$	300.64
LZY560	WIPER COUPLING FOR UVAS sc 50mm	\$	63.79
LZY568	TRANSACTION RELIEF CABLE	\$	356.53
LZY590	Fitting, Luft u. Probe	\$	647.54
LZY591	Cord, Power, Amtax and Phosphax sc	\$	310.28
LZY592	Verbinder Datenleitung sc	\$	408.57
LZY598	aa SC1000 Modbus TCP/IP Lizenz Kit	\$	3,796.58
LZY604	BATTERIES, RECHARGEABLE, NiMH PK/6	\$	162.46
LZY605	STANDARD CASE FOR TSS PORTABLE	\$	612.85
LZY606	BATTERY HOLDER	\$	53.58
LZY607	Plug charger for NiMH battery, for EU, UK, USA, AUS	\$	458.67
LZY613	90 DEGREE ELBOW FOR OVERFLOW	\$	66.10
LZY616	Flanschset Edelstahl f"r FP 360 sc	\$	1,356.75
LZY617	FLANGE SET,CS FOR FP 360 SC INSERTION MOUNTING HW	\$	585.87
LZY619	TUBING, FP360sc TO HOAB, 5M	\$	289.08
LZY620	TUBING, FP360sc TO HOAB, 10M	\$	655.25
LZY621	TUBING, FP360sc TO HOAB, 25M	\$	828.70
LZY623	EXTENSION CABLE, FP360sc, 1.5M	\$	6,012.86
LZY624	EXTENSION CABLE, FP360sc, 10M	\$	6,070.68
LZY625	O-RING, SCREW SET, FP360sc FLOW CELL	\$	450.96
LZY626	FITTING KIT FOR FP360sc FLOW CELL	\$	468.31
LZY627	Ersatz O-Ringe 68,00 x 4,00 - NBR70	\$	56.27
LZY630.00.10000	SS BALL VALVE SAFETY ARM. TSS sc INLINE, SS KIT	\$	13,374.77
LZY630.00.11000	SS BALL VALVE SAFETY ARM. TSS sc INLINE, CS KIT	\$	13,047.14
LZY630.00.12000	Sicherheitseinbauarmatur /o.Flansch	\$	12,835.15
LZY630.00.20000	SS Ball VALVE SAFE ARM. TSS EX1 sc INLINE, SS KIT	\$	13,374.77
LZY630.00.21000	SS Ball VALVE SAFE.ARM. TSS EX1 sc INLINE, CS KIT	\$	13,047.14
LZY630.00.22000	BALL VALVE SAFETY ARMATURE FOR TSS EX1 sc INLINE	\$	12,835.15
LZY630.00.30000	ff BALL VALVE MNT. FP360 W SS FL. SET	\$	15,205.61
LZY630.00.31000	ff BALL VALVE MNT. FP360 W C-S FL. SET	\$	14,685.26
LZY630.00.32000	ff BALL VALVE MOUNTING HARDWARE, FP360sc	\$	14,107.10
LZY634	WIPER SET TSS sc	\$	487.58
LZY635	MAINTENANCE WIPER SET TSS sc	\$	489.51
LZY652	PUMPHEAD, MAGNET DOSING PUMP	\$	2,775.17
LZY653	GASKET, SILICON, TSS sc, TRICLAMP MT	\$	53.58
LZY654	GASKET, PTFE, TSS sc, TRICLAMP MT	\$	135.87
LZY655	GASKET, FPM, TSS sc, TRICLAMP MT	\$	135.87
LZY656	CLAMP 2-PARTS, TS sc, SCREW TYPE	\$	809.42
LZY657	Schelle 3-teilig mit Schraubknebel	\$	1,085.01
LZY660	TUBE SET FILTRAX PHOTOMETER	\$	185.59
LZY661	PROBE HEAD W/PCB, FP360sc, 500ppB	\$	68,897.40



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LZY662	PROBE HEAD W/PCB, FP360sc, 5000ppb	\$	68,897.40
LZY664	TI REAR COVER W/PCBs, FP360sc	\$	29,004.36
LZY669	db FLOW CELLKIT, FP360sc, PANEL, HOSES	\$	9,327.65
LZY672	TUBING, FP360sc FOR FLOW CELL, 5M	\$	271.74
LZY673	TUBING, FP360sc FOR FLOW CELL, 10M	\$	381.59
LZY674	O-RING, BRACKET, FP360sc FLOW CELL	\$	2,717.35
LZY675	TFT 5,7" Color Display	\$	591.65
LZY676	sc1000 DISPLAY MODULE INNER FRAME	\$	82.48
LZY678	FILTER HOLDER SET 5M/230 VAC, FILTRAX ECO	\$	4,278.38
LZY679	SAMPLE HOSE, 2M HEATED 230 VAC, FILTRAX (ECO)	\$	1,478.16
LZY680	DELIVERY HOSE 2M,115V FOR FILTRAX (ECO)	\$	1,564.89
LZY691	WIPER PLUG (STAINLEES STEEL)	\$	102.53
LZY692	WIPER PLUG (PVC)	\$	142.42
LZY693	O-RING FOR MEASURING WINDOW, FP360sc	\$	1,088.87
LZY694	SENSOR CARTRIDGE INCL. SHIPPING BOOT, POLISHING STRIPS, AN-ISE sc / A-ISE sc / N-ISE sc	\$	4,818.00
LZY697	Kabelbinder Satz NH4D sc	\$	14.65
LZY698	SET OF CABLE TIES NISE sc & NO3D sc, YELLOW	\$	12.33
LZY703	Safety ring w/cable for armature LZU301.99.00000	\$	678.37
LZY706	CLEANING UNIT AN-ISEsc,A-ISEsc,N-ISEsc	\$	1,204.50
LZY707	Safety ring w/ cable for armature LZU300.99.00000	\$	576.23
LZY709	CLEANING CAP WITH NOZZLES	\$	620.56
LZY710	SET OF TUBE ADAPTER	\$	190.21
LZY711	CHECK VALVE	\$	75.74
LZY712	ADAPTER FOR CLEANING DEVICE	\$	5,723.78
LZY713	O-RING FOR SENSOR CARTRIDGE LZY694	\$	40.28
LZY714.99.00020	2 M SS INSTALLATION PIPE WITH FLANGE	\$	855.68
LZY714.99.00030	EXTENSION POLE 1M, STAINLESS STEEL, FOR LZY714.99.	\$	1,092.72
LZY714.99.00040	EXTENSION POLE 1.8M, STAINLESS STEEL, FOR LZY714	\$	1,306.64
LZY714.99.00050	IMMERSION POLE WITH SIDE HOLE, 2M, STAINLESS STEEL	\$	913.49
LZY714.99.00060	INSTALLATION PIPE 2M W WELDED HINGE	\$	988.65
LZY714.99.00100	END CAP BLACK FOR 40MM MOUNTING TUBES	\$	55.12
LZY714.99.02000	BRACKET 24CM, INCL. PLATE AND 2 SETS POLE CLAMPS	\$	1,397.22
LZY714.99.03000	BRACKET 10CM, INCL. PLATE AND 2 SETS POLE CLAMPS	\$	1,252.68
LZY714.99.32220	Beckenrandbefestigung - NISE	\$	2,351.18
LZY714.99.32320	POLE MOUNTING HARDWARE AN-ISE, 24CM BRACKET, SS PO	\$	3,006.43
LZY714.99.32620	POLE MOUNTING HARDWARE AISE, 24CM BRACKET, SS POLE	\$	3,006.43
LZY714.99.42020	POLE MOUNTING HARDWARE FILTERPROBE/FILTRAX, 24CM B	\$	2,042.83
LZY714.99.42120	POLE MOUNTING HARDWARE SONATAx, 24CM BRACKET, SS P	\$	2,100.65
LZY714.99.43050	POLE MOUNTING HARDWARE FILTRAX, 10CM BRACKET, SS P	\$	1,884.80
LZY714.99.50000	Probe adapter 90 Degree SS	\$	383.51
LZY714.99.52120	POLE MOUNTING HARDWARE SOLIDS, 24CM BRACKET, SS PO	\$	2,119.92
LZY714.99.52220	POLE MOUNTING HARDWARE NITRAT, 24CM BRACKET, SS PO	\$	2,081.38
LZY714.99.52520	POLE MOUNTING HARDWARE ORGANIC, 24CM BRACKET, SS P	\$	2,119.92
LZY714.99.53120	POLE MOUNTING HARDWARE SOLIDS, 10CM BRACKET, SS PO	\$	1,838.55
LZY714.99.53220	POLE MOUNTING HARDWARE NITRATAx, 10CM BRACKET, SS	\$	1,838.55
LZY714.99.53520	POLE MOUNTING HARDWARE UVAS, 10CM BRACKET, SS POLE	\$	1,838.55

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LZY714.99.60000	PIVOT EXTENSION, 35CM POLE & PROTECTION SHIELD	\$	799.79
LZY714.99.62160	POLE MOUNTING HW - SONATAX PIVOT, SS POLE 2M+0.35M	\$	2,524.63
LZY714.99.70000	PIVOT EXTENSION, 100CM POLE & PROTECTION SHIELD	\$	882.66
LZY714.99.72160	POLE MOUNTING HW - SONATAX PIVOT, SS POLE 2M + 1M	\$	2,524.63
LZY714.99.90000	Tube adapter for swivel mount LZY714	\$	142.42
LZY715	SCREWS FOR SENSOR CARTRIDGE (4 PCS)	\$	86.92
LZY716	COUPLING NUT AN/A/N-ISEsc	\$	398.93
LZY720	TEST CARTRIDGE FOR AN/A/N-ISEsc	\$	1,264.24
LZY721	HOUSING, GASKET SET, AN/A/N-ISEsc	\$	3,102.79
LZY722	PROBE ADAPTER FOR ANISE	\$	3,353.33
LZY723	Dichtungssatz	\$	210.06
LZY724	Goldkontakte, 3 Stk.	\$	379.66
LZY729	PROBE HEAD, FP360sc, 500, RECALIBRATED	\$	11,678.83
LZY736	SET OF SCREWS FOR PIPE AND ADAPTOR	\$	44.33
LZY739	SOLID STANDARD SET, 2 STANDARDS, FP360	\$	20,910.12
LZY743	SOLID STANDARD, MID-RANGE, FP360, 500PPB	\$	9,250.56
LZY744	SOLID STANDARD, HIGH-RANGE, FP360, 5000PPB	\$	9,250.56
LZY748-00	CF CARD WITH PROGRAMMING, RTC	\$	2,948.62
LZY749	SENSOR CABLE 10 M WITH FERRITE ESP. SOLITAX SC	\$	635.98
LZY755	KIT, CAM ROD, PISTON LINKAGE, AF7000	\$	618.63
LZY756	STREAMING CURRENT SENSOR, AF7000	\$	10,233.43
LZY757	PISTON, AF7000	\$	1,360.60
LZY758	KIT, SENSOR BEARING/ WASHER, AF7000	\$	1,121.63
LZY759	KIT, PART REPLACEMENT TOOL, AF7000	\$	1,231.48
LZY763	RELAY OUTPUT, AF7000	\$	493.36
LZY764	4-20mA PID OUTPUT CARD, AF7000	\$	2,235.55
LZY765	KIT, WTR CONNw/AUTO FLUSH, AF7000	\$	6,783.74
LZY766	KIT, WTR CONN w/o AUTO FLUSH AF7000	\$	3,160.61
LZY770	Befestigungsglasche mit Gewinde	\$	570.45
LZY771	SEDIMENTER FOR AISE, NISE AND AN-ISE sc	\$	8,286.96
LZY772	Flange Seal for UltraTurb	\$	37.19
LZY789	ASSEMBLING AID FOR O-RING, sc CABLE	\$	33.15
LZY791	AIR STRAINER FOR HOAB	\$	279.44
LZY804	MOUNTING PLATE HPL061	\$	556.96
LZY813	Nachr"stsatz Druckluftreinigung	\$	969.38
LZY814	SPARE PART SET AIR CLEANING FP 360 SC	\$	277.52
LZY815	Sondenkabel mit Ferrit 10m	\$	1,017.56
LZY820	ff ISOLATION KIT 5M SAMPLING TUBE FILTRAX/FILTERPROBE	\$	2,119.92
LZY821	ff ISOLATION KIT 10MSAMPLING TUBE FILTRAX/FILTERPROBE	\$	3,507.50
LZY822	Screw and sealing set	\$	55.89
LZY823	SCREW SET	\$	105.03
LZY827	Sockel 10cm	\$	398.93
LZY828.99.00002	TU5xxx sc Bubble Trap	\$	585.87
LZY830	VALVE PUMP CONVERSION KIT PSC & ASC	\$	2,717.35
LZY831	VALVE PUMP CONVERSION KIT AMTAX SC, 1 PUMP	\$	2,832.98
LZY832	VALVE PUMP CONVERSION KIT AMTAX SC,2 PUMPS	\$	5,588.88
LZY834	ee Process Vial	\$	128.93



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LZY835	StablCal Calibration Set w. RFID	\$	1,295.44
LZY835-LM	StablCal Calibration Set w. RFID	\$	1,295.44
LZY837	StablCal Cal. Vial 20 NTU w. RFID	\$	1,109.34
LZY837-LM	StablCal Cal. Vial 20 NTU w. RFID	\$	1,109.34
LZY838	StablCal Cal. Vial 600 NTU w. RFID	\$	1,109.34
LZY838-LM	StablCal Cal. Vial 600 NTU w. RFID	\$	1,109.34
LZY846	SCREW SET FOR LZY832	\$	28.33
LZY868	GLAND NUT FOR TSS sc / COSMOS	\$	55.12
LZY870	Set for Wall Bracket	\$	35.46
LZY871	Set for Flow Sensor Wall Bracket	\$	102.91
LZY873	Kit, Service Bracket, TU5xxx sc	\$	59.74
LZY874	NON RETURN VALVE A&P SC FOR CONT. MODE OPERATION	\$	61.86
LZY876	Desiccant cartridge TU5300/TU5400sc	\$	131.05
LZY877	StablCal Verf. Vial 10 NTU w. RFID	\$	1,024.43
LZY878	StablCal Verf. Vial 10 NTU w/o RFID	\$	852.79
LZY885.99.00002	PROGNOSYS COMMUNICATION CARD, W/MANUAL, FOR SC1000	\$	3,931.49
LZY886	Ersatz O-Ringe 59,92 x 3,53 - NBR70	\$	49.14
LZY887	Replacement O-Rings 40,00 x 3,50 -EPDM70	\$	49.14
LZY888	Replacem. o-rings TSS EX1 Inline armatures (NBR70)	\$	49.14
LZY894	KIT, MAINTENANCE TUBING FILTERSONDE sc	\$	86.92
LZY895	SET, O-RINGS, FILTERSONDE/SOLITAX sc	\$	13.49
LZY898	StablCal Calibration Set w/o RFID	\$	1,109.34
LZY898-LM	StablCal Calibration Set w/o RFID	\$	1,109.34
LZY899	StablCal Cal. Vial 20 NTU w/o RFID	\$	921.44
LZY899-LM	StablCal Cal. Vial 20 NTU w/o RFID	\$	921.44
LZY900	StablCal Cal. Vial 600 NTU w/o RFID	\$	921.44
LZY901	Glass Rod Sec. Turb. Stand. <0,1NTU	\$	3,197.95
LZY903	Vial Wiper TU5xxx	\$	59.74
LZY904.97.00002	Calibration Lid TU5xxx (RoW)	\$	74.78
LZY904.98.00002	Calibration Lid TU5xxx (EU)	\$	74.78
LZY906	Tool, Vial Removal, TU5xxx	\$	38.93
LZY907.35.00002	Maintenance Case TU5xxx; RoW	\$	3,892.94
LZY907.36.00002	Maintenance Case TU5xxx; RoW	\$	3,989.30
LZY907.91.00002	Maintenance Case TU5xxx; RoW	\$	4,143.48
LZY907.97.00002	Maintenance Case TU5xxx; RoW	\$	4,182.02
LZY907.98.00002	Maintenance Case TU5xxx; EU	\$	4,182.02
LZY907.LM.00002	Maintenance Case TU5xxx; RoW	\$	4,182.02
LZY910	Vial Compartment Wiper TU5xxx	\$	75.74
LZY912	Tubing Connect.set ULTRATURB to TU5	\$	79.02
LZY914	Seal for Automatic Cleaning Module TU5xxx	\$	80.56
LZY915	Silicone Wiper, Cleaning TU5xxx Silikonwischer, Reinigung TU5xxx	\$	138.95
LZY918	Seal for online vial TU5xxx	\$	59.74
LZY945	Micro Fiber cloth	\$	59.74
LZY946	Process Head Holder, TU5xxx sc	\$	269.81
LZY953	Calibration Syringe TU5xxx	\$	59.74



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

LZY954	Reducer Set 1/4" (6mm)	\$	109.27
LZY963	KIT, Flow Regulator TU5x00sc	\$	233.19
LZY971	ff Modem Kit	\$	1,946.47
LZY982	Socket for ethernet	\$	237.05
LZY990	Network Cable, M12 to RJ45	\$	235.12
LZY996	WiFi USB Adapter 2.4/5 GHz (US)	\$	609.00
LZY998	Sealing insert probe cable, slotted	\$	31.22
LZZ003	StablCal 1 NTU KIT W/ RFID TAG	\$	910.60
LZZ004	StablCal 1 NTU KIT W/O RFID TAG	\$	767.87
LZZ005	StablCal LR CALIBRATION SET w. RFID	\$	1,503.22
LZZ006	StablCal LR CALIBRATION SET w/o RFID	\$	1,418.30
M_2401_064_FF	TBG, 0.8ID X 1.6OD MM, PTFE	\$	52.61
M_2401_316_FF	PTFE TUBING 2 ID x 6 OD MM	\$	38.93
M00X002	SWITCH-ON-OFF/TIM900,CDM210,CDM230	\$	69.57
M00Y001	RS232 FLAT CABLE	\$	200.43
M02P003	Electrovalve,MDE150	\$	1,634.27
M02P510	Keyboard, MDE150 FROM R7	\$	2,235.55
M02T003	VENT FAN/ABU900	\$	458.67
M02T301	CASING, ABU901, WHITE	\$	529.98
M02T511	RACKSHAFT MOTOR ASSEMBLY/ABU901	\$	13,586.76
M02T512	STOPCOCK OPTICAL CELLS/ABU901	\$	470.24
M02T513	RACKSHAFT OPT. CELLS/ ABU901	\$	477.95
M035J35	CLIP CAP ELECTRODE 12mm	\$	8.29
M03M501	MAIN BOARD/PHM210	\$	1,848.18
M04M501	MAIN BOARD/PHM220	\$	3,989.30
M04M503	MAIN BOARD PHM240	\$	5,588.88
M04V501	MAIN BOARD CTV101	\$	2,775.17
M05M503	MAIN BOARD/CDM230	\$	9,385.46
M05V101	ON/OFF SWITCH, PGZ SERIES (S/NS >314593)	\$	229.34
M09T317	Window & Lexan ID, BIOGAS	\$	720.77
M11T301	CASING, AUTOCAT 9000	\$	2,081.38
M11T302	WINDOW, AUTOCAT 9000	\$	1,393.37
M11T503	KEYPAD REPLACEMENT, AUTOCAT	\$	1,002.14
M14T308	BURETTE LOCKING ASSEMBLY	\$	3,141.34
M14T501	RFID BOARD, PO66	\$	1,268.10
M16T207	RFID MODULE, SACX50	\$	8,036.42
M16T210	ULTRASONIC TRANSDUCER, SACx50	\$	2,254.82
M16T211	Z ARM STEPPER MOTOR W/CONNECTORS, SACX50	\$	2,023.56
M16T212	MOTOR & OPTOCOUPLER, SACX50	\$	2,081.38
M16T302	NEEDLE ROLLER BEARING, SACX50	\$	695.72
M16T309	DRIVING BELTS, Z ARM, SACX50	\$	1,050.32
M16T310	KEYBOARD LEXAN, SACX50	\$	520.34
M16T314	MAINT KIT, SACX50 SAMPLE TRAY REPAIR	\$	653.32
M16T315	UPDATE KIT, SACx50, Z ARM MECHANIC	\$	6,224.86
M16T317	PTFE BEARING, 4PCS, SACX50	\$	684.16
M16T320	SHOWERHEAD UPGD KIT, SACX50	\$	238.97
M16T321	Ballbearings, shower head, SACx50,	\$	770.88



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

M16T322	Z ARM, SACX50, 2 LATERAL SIDES	\$	9,693.82
M16T501	MAIN PCB BOARD (P069), SACX50	\$	13,509.67
M16T502	POWER SUPPLY, SACX50	\$	1,431.91
M16T504	KEYBOARD PCB, SACX50	\$	1,212.21
MAN-3444	MANUAL-COND/RESIST.SNSRS,	\$	122.18
MAN-3455	MANUAL-COND/RESIST.SNSRS,	\$	122.18
MAN-E53	Electrodeless 1/2 DIN	\$	192.72
MBW-5ML	dd Binding Beads W, 5 mL	\$	2,945.00
MH018S8SZ	Sanitary Mounting Kit, 2"	\$	2,100.65
MH1012S	ASSY-ADAPTER, INSERTION,	\$	894.22
MH113M2C	CONDUCTIVITY, INSERTION, 0.5" NPT	\$	2,274.10
MH114M3C	CONDUCTIVITY, INSERTION, 0.75" NPT	\$	2,409.00
MH116M3MZ	316SS Insertion Hardware	\$	4,509.65
MH118M9NZ	INSERTION,316SS,2" PROC.C	\$	12,623.16
MH138	INSERTION,CPVC,2" PROC.CO	\$	2,312.64
MH138M9NZ	INSRTN,CPVC,2"NPT,NO J-BX	\$	3,179.88
MH236B00Z	MH,1" pH,HANDRAIL,CPVC	\$	1,736.41
MH236B00Z-1001	MH, 1"pH, Handrail, 10 ft. Length	\$	2,004.29
MH276H00F	D.O. Mounting, Handrail	\$	1,788.44
MH276J00F	D.O. Mounting, Floatation	\$	2,524.63
MH276R00FP	MTG.HDWR-D.O.PVC,PAPR TAG	\$	661.03
MH313N3NZ	MH,TEE,316 SS,3/4", FLOW	\$	618.63
MH314N4MZ	Flow-Thru Hardware, 316SS	\$	813.28
MH333N3NZ	MH-TEE,CPVC,3/4", FLOW TH	\$	393.15
MH334N4NZ	FINISHED MTG HARDWARE	\$	489.51
MH336	MH,FLO-THRU,CPVC,1+1/2"	\$	555.03
MH373N3NZ	MH,TEE,PVC,3/4",FLOW-THRU	\$	292.93
MH376	pH/ORP Mounting, ,PVC	\$	489.51
MH378W4PZ	BOM-ASSY,DO,MH,FLOW-THRU	\$	1,061.89
MH414A00B	Submersion Hardware,316SS	\$	1,285.44
MH432G	CPVC PIPE MOUNT W/ JUNCTION BOX	\$	753.54
MH434A00B	Submersion Hardware, CPVC	\$	813.28
MH434B	SUBM.W/PIPE J-BOX,CPVC,1"	\$	813.28
MH462G	FINISHED MOUNTING HARDWAR	\$	1,798.08
MH474B	pH/ORP Immersion Mounting	\$	691.86
MH516N9Z	Union Mounting Kit, 316SS	\$	1,634.27
MH518M3NZ	MTG HARDWARE,316SS,W/O TE	\$	2,447.54
MH518N3NZ	Cond. Union Mounting	\$	3,141.34
MH518N9A	FINISHED MOUNTING HARDWAR	\$	2,794.44
MH538M3NZ	MTG HARDWARE,CPVC,W/O TEE	\$	614.78
MH538N3N	FINISHED MTG HARDWARE	\$	728.48
MH538N3NZ	Cond. Union Mounting	\$	828.70
MH538N9A	FINISHED MOUNTING HARDWAR	\$	788.22
MH568N3NZ	MTG HARDWARE,PVDF,W/O J-B	\$	2,235.55
MH576N3MZ	ASY-MTG HARDWARE	\$	616.70
MH578N9A	pH/ORP Mounting, PVC	\$	691.86
MH578N9AZ	FINISHED Mtg Hdwr,No JBox	\$	423.98



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

MH736M4MZ-1000	vv 10.5" Insert Depth & Mod. Adapter	\$	2,929.34
MHF15B2	Mounting Tee, 1" PVC	\$	855.68
MHF15D2	Mounting Tee, 1.5" PVC	\$	824.84
MHF15E2	Mounting Tee, 2" PVC	\$	753.54
MHF15G2	Mounting Tee, 3" PVC	\$	948.18
MHF15H2	Mounting Tee, 4" PVC	\$	1,135.12
MHF15J2	Mounting Tee, 0.5" PVC	\$	855.68
MHF15K2	Mounting Tee, 0.75" PVC	\$	855.68
MHF31B2	MH-TEE,CPVC,1",FLW-THRU	\$	1,264.24
MHF31D2	MH-TEE,CPVC,1.5",FLW-THRU	\$	1,050.32
MHF31E2	MH-TEE,CPVC,2",FLW-THRU	\$	1,135.12
MHF31G2	MH-TEE,CPVC,3",FLW-THRU	\$	1,509.00
MHF31H2	MH-TEE,CPVC,4",FLW-THRU	\$	1,668.96
MP000147-01	Adapter, 6 mm Ferrule: 1/4" compression fitting; W	\$	208.14
MP000344-01	Overlay: 4-Key Membrane Switch; 3 LEDs; ARTI; WPC-	\$	3,642.41
MP4K	nn aa MP-4 CASE W/STANDARDS	\$	118.91
MP6K	nn aa MP-6 & MP-6p CASE W/STANDARDS	\$	439.40
MSM_PHD	oo offering no longer available	\$	805.57
MTC10101	aa ORP GEL-FILLED PROBE, STD w/1m CABLE	\$	1,516.71
MTC10103	aa ORP GEL-FILLED PROBE, STD w/3m CABLE	\$	1,626.56
MTC10105	aa ORP GEL-FILLED PROBE, RUG w/5m CABLE	\$	2,698.08
MTC10110	aa ORP GEL-FILLED PROBE, RUG w/10m CABLE	\$	3,025.70
MTC10115	aa ORP GEL-FILLED PROBE, RUG w/15m CABLE	\$	3,218.42
MTC10130	aa ORP GEL-FILLED PROBE, RUG w/30m CABLE	\$	3,391.87
MTC30101	aa ORP REFILLABLE PROBE STD, w/1m CABLE	\$	1,694.01
MTC30103	aa ORP REFILLABLE PROBE STD, w/3m CABLE	\$	1,801.93
MTC30601	ee MTC306 METAL PROBE, STD	\$	1,108.14
MTC69501	ee MTC695 METAL PROBE, STD	\$	1,194.86
NFW-55ML	dd Nuclease-Free Water, 55mL	\$	29.63
PC1R1A	Convertible pH Sensor	\$	1,067.67
PC1R1A-33T	PC1R1A Electrode with 33 Foot Cable	\$	1,252.68
PC1R1A-STC	vv pH Sensor for Cold Temp.Operation	\$	1,965.74
PC1R1A-V12	Convertible pH Sensor with pt100	\$	1,067.67
PC1R1N	Convertible pH Sensor	\$	952.04
PC1R2A	Convertible pH Sensor	\$	1,189.08
PC1R2N	Convertible pH Sensor	\$	1,067.67
PC1R3A	Convertible pH Sensor	\$	1,410.71
PC2K1A	Insertion pH Sensor	\$	1,298.93
PC2K2A	Insertion pH Sensor, Flat Glass	\$	1,410.71
PC3K2A	Sanitary pH Sensor	\$	2,023.56
PD1P1	PHD, PEEK, pH-Wide Range	\$	3,044.98
PD1P1.99	PHD, PEEK, pH-Wide Range, Generic	\$	3,044.98
PD1P1A15P1S1N	PHD, PEEK, pH-Wide Range	\$	4,374.74
PD1P1A15P2S1N	PHD, PEEK, pH-Wide Range, Perflouro	\$	4,586.74
PD1P1A15V2S1N	PHD, PEEK, pH-Wide Range	\$	3,353.33
PD1P1A25	ff PHD, PEEK, pH-Wide Range	\$	4,297.66
PD1P1A33	PHD, PEEK, pH-Wide Range	\$	4,297.66

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

PD1P1A50	PHD, PEEK, pH-Wide Range	\$	4,490.38
PD1P3	PHD, PEEK, pH-HF Glass	\$	3,777.31
PD1P3A30	PHD,HF Resist.Elect,PEEK	\$	4,933.63
PD1P3A50	PHD,HF Resistant Electrode, PEEK	\$	4,374.74
PD1R1	PHD,RYTON, pH-Wide Range	\$	2,794.44
PD1R1.99	PHD, Ryton, pH-Wide Range, Generic	\$	2,794.44
PD1R3	PHD, Ryton, pH-HF Glass	\$	3,565.32
PD1R3.99	PHD, Ryton, HF Resistant	\$	3,565.32
PD1R3A25	ELECTRODE, PHD, HF RES RYTON, 25FT	\$	3,777.31
PD2P1	PHD, PEEK, pH-Wide Range	\$	3,160.61
PD2P1.99	PHD, PEEK, pH-Wide Range, Generic	\$	3,160.61
PD2P1A15P2S1N	PHD, PEEK, pH-Wide Range, Insertion	\$	4,702.37
PD2P1A30	PHD, PEEK, pH-Wide Range	\$	4,374.74
PD2P1A50	PHD, PEEK, pH-Wide Range	\$	3,680.95
PD2P3	PHD, PEEK, pH-HF Glass	\$	3,892.94
PD3P1	PHD, PEEK, pH-Wide Range	\$	3,353.33
PD3P1A50	PHD, PEEK, pH-Wide Range	\$	3,777.31
PEH30N	M.H.-HORIZ.MNT,ALUM,3 HOL	\$	982.87
PHC10101	pH GEL PROBE, STD, w/1m CABLE	\$	928.91
PHC10103	pH GEL PROBE, STD, w/3m CABLE	\$	1,104.29
PHC10105	pH GEL PROBE, RUGGED, w/5m CABLE	\$	2,698.08
PHC10110	pH GEL PROBE, RUGGED, w/10m CABLE	\$	3,025.70
PHC10115	pH GEL PROBE, RUGGED, w/15m CABLE	\$	3,218.42
PHC10130	pH GEL PROBE, RUGGED, w/30m CABLE	\$	3,546.05
PHC10801	ee PHC108 pH PROBE, PUNCTURE, STD	\$	1,555.25
PHC20101	pH GEL-FILLED PROBE, GENERAL PURPOSE, STD w/1m CABLE	\$	1,013.71
PHC20103	pH GEL-FILLED PROBE, GENERAL PURPOSE, STD w/3m CABLE	\$	1,165.96
PHC28101	aa ULTRA pH REFILLABLE PROBE, w/1m CABLE	\$	1,985.02
PHC28103	aa ULTRA pH REFILLABLE PROBE, w/3m CABLE	\$	2,081.38
PHC30101	pH LIQUID PROBE, STD, w/1m CABLE	\$	1,104.29
PHC30101-LM	pH LIQUID PROBE, STD, w/1m CABLE	\$	1,104.29
PHC30103	pH LIQUID PROBE, STD, w/3m CABLE	\$	1,266.17
PHC70501	ee PHC705 pH PROBE, RED-ROD, GLASS, GEN PURP, STD	\$	1,659.32
PHC70501AP	KTO: PHC705 RedRod Lab pH Electrode with Cal. Pack	\$	1,530.20
PHC705A01	ee PHC705A pH PROBE, RED-ROD, GLASS, ALKALINE, STD	\$	1,659.32
PHC705A01AP	KTO: PHC705A RedRod Lab pH Electrode with Cal. Pack	\$	1,530.20
PHC72501	ee PHC725 pH PROBE, RED-ROD, GLASS, GEN PURP, STD	\$	1,659.32
PHC72501AP	KTO: PHC725 RedRod Lab pH Electrode with Cal. Pack	\$	1,530.20
PHC72901	ee PHC729 pH PROBE, RED-ROD, GLASS, FLAT, STD	\$	2,177.74
PHC72901AP	KTO: PHC729 RedRod Lab pH Electrode with Cal. Pack	\$	2,004.29
PHC73501	ee PHC735 pH PROBE, RED-ROD, GLASS, HIGH FLOW, STD	\$	1,927.20
PHC73501AP	KTO: PHC735 RedRod Lab pH Electrode with Cal. Pack	\$	1,757.61
PHC74501	ee PHC745 pH PROBE, RED-ROD, GLASS, SLEEVE, STD	\$	2,119.92
PHC74501AP	KTO: PHC745 RedRod Lab pH Electrode with Cal. Pack	\$	1,946.47
PHC80501	ee PHC805 pH PROBE, GLASS, GEN PURP, STD	\$	1,131.27
PHC80501AP	KTO: PHC805 Lab pH Electrode with Cal. Pack	\$	1,073.45
PRO-C3A1N	PRO Series Transmitter	\$	2,466.82



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

PRO-E3A1N	PRO Series Transmitter	\$	2,466.82
PRO-P3A1N	PRO Series Transmitter	\$	2,466.82
PS4ABASECHLOR	WDM PIPESONDE IN-PIPE PROBE,CHLORINE	\$	-
QG21IT-100C	QUENCHGONE21 INDUSTRIAL COMPLETE, tATP ONLY, 100 TESTS	\$	5,203.44
QG21ST-100C	QUENCH-GONE SPECIALTY ATP, 100 TESTS	\$	6,901.79
QG21W-50C	QUENCHGONE21 WW COMPLETE, 50 SAMPLES	\$	5,203.44
QG21WA-25C	QUENCHGONE21 WW ADVANCED COMPLETE, 25 SAMPLES	\$	3,541.23
QGA-100	QUENCH-GONE AQUEOUS REAGENT-ONLY, 100 TESTS	\$	2,547.52
QGA-100C	QUENCH-GONE AQUEOUS KIT, 100 TESTS	\$	4,751.75
QGA-25C	QUENCH-GONE AQUEOUS KIT, 25 TESTS	\$	1,374.94
QGOM-100C	QUENCH-GONE ORG MOD COMPLTE, 100TESTS	\$	5,329.91
QGOM-5	QUENCH-GONE ORG MOD DEMO, 5 TESTS	\$	231.26
R7540	OXYGEN CHEMETS REFILL	\$	118.69
RC1R5N	Convertible ORP Sensor	\$	1,040.69
RC1R5N-HF	ORP COMB SENSOR W/ ULTEM	\$	1,526.34
RD1P5	PHD, PEEK, ORP-Platinum	\$	3,353.33
RD1P5.99	ORP, PEEK, Pt Electrode, Generic	\$	3,353.33
RD1P5A15V2S1N	PHD, PEEK, ORP-CONVERTBLE	\$	3,565.32
RD1P5A50	PHD, PEEK, ORP-Platinum, Convertible	\$	3,777.31
RD1P5A80	PHD, PEEK, ORP-Platinum	\$	4,933.63
RD1P6	PHD, PEEK, ORP-Gold	\$	3,565.32
RD1R5	PHD, RYTON,ORP-Platinum	\$	3,044.98
RD1R5.99	ORP, Ryton, Pt Electrode, Generic	\$	3,044.98
RD1R5A-1000	K = 100 Foot Cable Length	\$	4,702.37
RD1R5A25	PHD, RYTON,ORP-Platinum	\$	3,353.33
RD1R5A50	PHD, RYTON,ORP-Platinum	\$	3,565.32
RD1R6	PHD, Ryton, ORP-Gold	\$	3,353.33
RD2P5	PHD, PEEK, ORP-Platinum	\$	3,468.96
RD2P5A50	PHD, PEEK, ORP-Platinum	\$	4,413.29
RD2P6	PHD, PEEK, ORP-Gold	\$	3,680.95
RD3P5	PHD, PEEK, ORP-Platinum	\$	3,565.32
RP6300A4100-001	ASSY-CABLE W/HALL EFFECT	\$	589.72
RTCMISC	ff RTC MISCELLANEOUS	\$	0.00
S41M006	ELECTRODE FILLING SOLN, NITRATE, RA	\$	83.64
SA000120-01	mc Brush Kit: Flow Cell Cleaning	\$	69.76
SA000120-02	Brush Kit: Flow Cell Cleaning Brush; 1UM; WPC	\$	364.24
SARSMM-0.825ML	dd SARS-CoV-2 Advanced Mastermix, 0.825mL	\$	2,908.87
SB-P1SP	PEEK Salt Bridge, Kynar	\$	851.82
SB-P1SP-LM	PEEK Salt Bridge, Kynar	\$	851.82
SB-P1SV	PEEK Salt Bridge, Kynar	\$	249.33
SB-P1SV-LM	PEEK Salt Bridge, Kynar	\$	249.33
SB-P1SV-NL	PEEK Salt Bridge, Kynar	\$	234.88
SB-P2SV	PEEK Salt Bridge, Ceramic	\$	328.83
SB-P2SV-LM	PEEK Salt Bridge, Ceramic	\$	328.83
SB-R1SV	Ryton Salt Bridge, Kynar	\$	249.33
SB-R1SV-LM	Ryton Salt Bridge, Kynar	\$	249.33

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

SCP001	SCP Proximity Switch Kit, Thermal Shutoff Valve	\$	1,081.16
SE 244	SCREW,4-20,TWIN HELIX PLASTITE	\$	34.11
SM5010	pH Beakers, natural	\$	31.80
SM5011	pH Beakers, red	\$	31.80
SM5012	pH Beakers, yellow	\$	31.80
SM5013	pH Beakers, blue	\$	31.80
SM5014	pH Beakers, green	\$	31.80
SM730-1033	Ultraviolet Blocking Eyewear	\$	124.88
SMART-LINK	SMART-LINK INTERFACE MOD.	\$	15,244.15
SPMP-CL	SPMP- SC4500, CL17SC	\$	17,286.98
SPMP-TU	SPMP - SC4500, TU5300	\$	21,391.92
SRW6321	CASE SCREWS, MP, EA	\$	12.53
SWCUSTOM	SOURCE WATER PANEL, CUSTOM	\$	-
TL2300.NA	KIT, TL2300 LAB TURBIDIMETER, EPA, NA	\$	8,845.85
TL2310.NA	KIT, TL2310 LAB TURBIDIMETER, ISO, NA	\$	8,903.66
TL2350.NA	KIT, TL2350 LAB TURBIDIMETER, EPA, NA	\$	14,010.74
TL2360.NA	KIT, TL2360 LAB TURBIDIMETER, ISO, NA	\$	13,856.57
TNT810	TOC, TNT+, LR(1.5-30MG/L), PK/25	\$	1,179.81
TNT810-LM	TOC, TNT+, LR(1.5-30MG/L), PK/25	\$	1,179.81
TNT811	TOC, TNT+, HR (30 - 300 MG/L)	\$	1,179.81
TNT811-LM	TOC, TNT+, HR (30 - 300 MG/L)	\$	1,179.81
TNT815	rr COD IN SALTWATER,TNT+, LR (7-70 MG/L), PK/25	\$	176.51
TNT816	rr COD IN SALTWATER ,TNT+, HR (70-1000 MG/L), PK/25	\$	176.51
TNT817	TNTplus International Bitterness Units	\$	339.67
TNT819	TNTplus VDK (vicinal diketones) Vial Tests	\$	334.25
TNT819-LM	TNTplus VDK (vicinal diketones) Vial Tests	\$	334.25
TNT820	rr COD,TNT+,ULR(1.0-60 MG/L) PK/25	\$	195.13
TNT82006	rr KTO: COD,TNT+ ULR(150 VIALS)PK/6	\$	1,060.56
TNT820-LM	rr COD,TNT+,ULR(1.0-60 MG/L) PK/25	\$	195.13
TNT821	rr COD,TNT+, LR (3-150 MG/L), PK/25	\$	195.13
TNT82106	rr KTO: TNT+,LR COD, 6-PK(150 VIALS) (3-150MG/L)	\$	1,060.56
TNT82106-LM	rr KTO: TNT+,LR COD, 6-PK(150 VIALS) (3-150MG/L)	\$	1,060.56
TNT821-LM	rr COD,TNT+, LR (3-150 MG/L), PK/25	\$	195.13
TNT821-PH	rr COD,TNT+, LR (3-150 MG/L), PK/25	\$	174.89
TNT822	rr COD,TNT+, HR(20-1500 MG/L), PK/25	\$	195.13
TNT82206	rr KTO: COD TNT+,HR 6-PK(150 VIALS)	\$	1,060.56
TNT82206-LM	rr KTO: COD TNT+,HR 6-PK(150 VIALS)	\$	1,060.56
TNT822-LM	rr COD,TNT+, HR(20-1500 MG/L), PK/25	\$	195.13
TNT822-PH	COD,TNT+ ,HR (20-1500 MG/L), PK/25	\$	174.89
TNT823	rr COD TNT+,HR+(250-15000 MG/L)PK/25	\$	195.13
TNT82306	rr KTO: TNT+COD HR+,(150VIALS)6PK/25	\$	1,060.56
TNT82306-LM	rr KTO: TNT+COD HR+,(150VIALS)6PK/25	\$	1,060.56
TNT823-LM	rr COD TNT+,HR+(250-15000 MG/L)PK/25	\$	195.13
TNT823-PH	COD TNT+,HR+(250-15000 MG/L)PK/25	\$	174.89
TNT824	rr COD TNT+ UHR (5000-60000 MG/L) PK/25	\$	184.29
TNT824-LM	rr COD TNT+ UHR (5000-60000 MG/L) PK/25	\$	184.29
TNT825	HG-FREE COD TNT+ HR(25-1000MG/L) PK/25	\$	195.13

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

TNT82506	KTO: COD2, TNT+, HR (MERC FREE) PK/6	\$	1,060.56
TNT825-LM	HG-FREE COD TNT+ HR(25-1000MG/L) PK/25	\$	195.13
TNT826	NITROGEN TOT LR TNT+, 1-16MG/L PK/25	\$	252.95
TNT826-LM	NITROGEN TOT LR TNT+, 1-16MG/L PK/25	\$	252.95
TNT827	NITROGEN TOT HR TNT+, 5-40MG/L PK/25	\$	242.10
TNT827-LM	NITROGEN TOT HR TNT+, 5-40MG/L PK/25	\$	242.10
TNT827-PH	NITROGEN TOT HR TNT+, 5-40MG/L PK/25	\$	252.95
TNT828	NITRGN TOT UHR TNT+ 20-100MG/L PK/25	\$	252.95
TNT828-LM	NITRGN TOT UHR TNT+ 20-100MG/L PK/25	\$	252.95
TNT829	TNTplus Hach Ammonia ULR+ 0,5 - 5 mg/L NH3-N	\$	180.48
TNT830	AMMONIA TNT+, ULR, 0.015-2.0 MG/L PK/25	\$	184.29
TNT830-LM	AMMONIA TNT+, ULR, .015-2.0MG/L PK/25	\$	184.29
TNT831	AMMONIA, TNT+, LR (1-12 MG/L) PK/25	\$	184.29
TNT831-LM	AMMONIA, TNT+, LR (1-12 MG/L) PK/25	\$	184.29
TNT831-PH	AMMONIA, TNT+, LR (1-12 MG/L) PK/25	\$	182.48
TNT832	AMMONIA, TNT+, HR (2-47 MG/L) PK/25	\$	184.29
TNT832-LM	AMMONIA, TNT+, HR (2-47 MG/L) PK/25	\$	184.29
TNT832-PH	AMMONIA, TNT+, HR (2-47 MG/L) PK/25	\$	182.48
TNT833	TNTplus Ammonia Vial Test	\$	184.29
TNT833-LM	TNTplus Ammonia Vial Test	\$	184.29
TNT834	AMMONIA, TNT+, UHR (100-1800 MG/L) PK/25	\$	180.12
TNT835	NITRATE,TNT+ LR, 0.2-13.5MG/L PK/25	\$	146.35
TNT835-LM	NITRATE,TNT+ LR, 0.2-13.5MG/L PK/25	\$	146.35
TNT835-PH	NITRATE,TNT+ LR, 0.2-13.5MG/L PK/25	\$	146.71
TNT836	NITRATE, TNT+ HR (5-35MG/L -N) PK/25	\$	151.57
TNT836-LM	NITRATE, TNT+ HR (5-35MG/L -N) PK/25	\$	151.57
TNT836-PH	NITRATE, TNT+ HR (5-35MG/L -N) PK/25	\$	146.71
TNT837	TNTplus Hach Ammonia HR+ 10-100 mg/L NH3-N	\$	180.48
TNT838	NITRATE, TNT+ UHR (15-150 MG/L-N) PK/25	\$	151.57
TNT839	db NITRITE, TNT+ LR 0.015-0.6MG/L PK/25	\$	126.28
TNT839-LM	db NITRITE, TNT+ LR 0.015-0.6MG/L PK/25	\$	126.28
TNT840	db HR NITRITE TNT+ (0.6-6.0 MG/L, NO2-N), PK/25	\$	126.28
TNT840-LM	db HR NITRITE TNT+ (0.6-6.0 MG/L, NO2-N), PK/25	\$	126.28
TNT841	UHR NITRITE TNT+ (2-90 MG/L, NO2-N), PK/25	\$	125.03
TNT843	aa PHOSPHORUS TNT+ LR PK/25 (0.05-1.5MG/L PO4-P)	\$	184.29
TNT843-LM	aa PHOSPHORUS TNT+ LR PK/25 (0.05-1.5MG/L PO4-P)	\$	184.29
TNT844	aa PHOSPHORUS TNT+ HR PK/25 (0.5-5 MG/L PO4-P)	\$	184.29
TNT844-LM	aa PHOSPHORUS TNT+ HR PK/25 (0.5-5 MG/L PO4-P)	\$	184.29
TNT844-PH	aa PHOSPHORUS TNT+ HR PK/25 (0.5-5 MG/L PO4-P)	\$	187.90
TNT845	aa PHOSPHORUS TNT+ UHR PK/25 (2-20MG/L PO4-P)	\$	184.29
TNT845-LM	aa PHOSPHORUS TNT+ UHR PK/25 (2-20MG/L PO4-P)	\$	184.29
TNT846	PHOS, TNT+, ORTHO, 1.67-30MG/L PK/25	\$	126.28
TNT846-LM	PHOS, TNT+, ORTHO, 1.67-30MG/L PK/25	\$	126.28
TNT848	ALUMINUM, TNT+, 0.02-0.5MG/L, PK/25	\$	252.95
TNT848-LM	ALUMINUM, TNT+, 0.02-0.5MG/L, PK/25	\$	252.95
TNT849	MAGNESIUM TNT+ (0.5-50 MG/L) PK/25	\$	254.75
TNT850	LEAD, TNT+, 0.1 - 2.0 MG/L, PK/25	\$	252.95



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

TNT850-LM	TNTplus HACH Lead Vial Test (LM), Measuring Range 0,1-2,0 mg/L	\$	252.95
TNT852	CADMIUM, TNT+, (0.02-0.3 MG/L) PK/25	\$	252.95
TNT852-LM	CADMIUM, TNT+, (0.02-0.3 MG/L) PK/25	\$	252.95
TNT854	CHROMIUM, TNT+, (0.03-1 MG/L) PK/25	\$	252.95
TNT854-LM	CHROMIUM, TNT+, (0.03-1 MG/L) PK/25	\$	252.95
TNT856	NICKEL, TNT+, 0.1-6 MG/L, PK/25	\$	252.95
TNT856-LM	NICKEL, TNT+, 0.1-6 MG/L, PK/25	\$	252.95
TNT858	db IRON TNT+ (0.2-6.0 MG/L), PK/25	\$	184.29
TNT858-LM	db IRON TNT+ (0.2-6.0 MG/L), PK/25	\$	184.29
TNT860	db COPPER TNT+ (0.1-8.0 MG/L), PK/25	\$	252.95
TNT860-LM	db COPPER TNT+ (0.1-8.0 MG/L), PK/25	\$	252.95
TNT861	SULFIDE TNT+ (0.1-2.0 MG/L) PK/25	\$	234.88
TNT861-LM	SULFIDE TNT+ (0.1-2.0 MG/L) PK/25	\$	234.88
TNT862	TNTplus HACH Cyanide,free Vial Test	\$	229.46
TNT864	SULFATE TNT+ LR 40-150 MG/L, PK/25	\$	238.49
TNT864-LM	SULFATE TNT+ LR 40-150 MG/L, PK/25	\$	238.49
TNT865	SULFATE TNT+ HR, 150-900 MG/L, PK/25	\$	238.49
TNT865-LM	SULFATE TNT+ HR, 150-900 MG/L, PK/25	\$	238.49
TNT866	db FREE CHLORINE TNT+ , PK/24	\$	138.76
TNT866-LM	db TNTplus Free Chlorine VialTest (LM), Measuring Range 0,05-2,0 mg/L	\$	138.76
TNT867	db TOTAL & FREE CHLORINE TNT+, PK/24	\$	152.49
TNT867-LM	TOTAL + FREE CHLORINE TNT+, PK/24	\$	152.49
TNT868	TNTplus HACH Phenol Vial Test	\$	258.37
TNT868-LM	TNTplus HACH Phenol Vial Test	\$	231.26
TNT869	WATER HARDNESS TNT+ (20-350 MG/L) PK/25	\$	234.88
TNT870	TOTAL ALKALINITY TNT+ (25-400MG/L), PK/25	\$	146.15
TNT870-LM	TOTAL ALKALINITY TNT+ (25-400MG/L), PK/25	\$	146.15
TNT870-PH	TOTAL ALKALINITY TNT+ (25-400MG/L), PK/25	\$	143.46
TNT871	FORMALDEHYDE TNT+ (0.5-10 MG/L) PK/25	\$	234.88
TNT872	VOLATILE ACIDS TNT+ (50-2500 MG/L), PK/25	\$	240.30
TNT872-LM	TNTplus HACH Volat. Acids Test (LM), Measuring Range 50-2500 mg/L	\$	240.30
TNT872-PH	VOLATILE ACIDS TNT+ (50-2500 MG/L), PK/25	\$	238.49
TNT874	TNTplus HACH Anionic Surfactants 0.1-4.0 mg/L	\$	243.91
TNT874-LM	TNTplus HACH Anion.Surfactants (LM), Measuring Range 0,1-4,0 mg/L	\$	243.91
TNT875	NONIONIC SURFACTANTS TNT+ LR(0.2-6.0 MG/L) PK25	\$	242.10
TNT876	NONIONIC SURFACTANTS TNT+ HR (6-200 MG/L) PK25	\$	242.10
TNT877	aa db TNTplus Hach Boron Vial Test	\$	245.72
TNT877-LM	aa db TNTplus Hach Boron Vial Test	\$	245.72
TNT878	TNTplus Hach Fluoride Vial Test	\$	184.29
TNT878-LM	TNTplus Hach Fluoride Vial Test	\$	184.29
TNT879	rr CHLORIDE TNT+ (1-1000 MG/L) PK/24	\$	234.88
TNT879-LM	CHLORIDE TNT+ (1-1000 MG/L) PK/24	\$	234.88
TNT880	s-TKN TNT+ (0-16 MG/L N), PK/25	\$	482.40



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

TNT880-LM	s-TKN TNT+ (0-16 MG/L N), PK/25	\$	482.40
TNT882	ANAMMOX ACTIVITY, TNT+, PK/25	\$	179.40
TNT882KTO	KTO: ANAMMOX ACTIVITY, TNT+, PK/25	\$	263.79
TNT885	CAT. SURFACTANTS, TNT+, (0.2-2.0 MG/L), PK/25	\$	224.04
TNT890	TNT+ METALS PREP SET (50-DIGESTS)	\$	278.24
TNT890-LM	TNT+ METALS PREP SET (50-DIGESTS)	\$	278.24
TNT892	CALCIUM SEPARATION SET FOR TNT852	\$	202.36
TNT892-LM	CALCIUM SEPARATION SET FOR TNT852	\$	202.36
TNT919	db SAMPLE BLANK VIALS, TNTPLUS, 5/PK	\$	50.23
TNT999E	aa vv TNTPLUS DEMO KIT	\$	106.04
U53BRACKET	oo FIXED U53 SENSOR SUPPORT	\$	977.09
UC1-5ML	ULTRACHECK 1, 2 x 2.5ML DROPPER BOTTLES	\$	262.10
UCCASTB	RPL SNAP-IN COND CELL W/THERM ASSY	\$	474.09
UIIDD MOS	DIG DISP W/BD REPL W/SUPPORTS	\$	728.48
UL7-125ML	ULTRALYSE 7, 125ML BOTTLE	\$	283.30
UL7-5ML-25R	ULTRALYSE 7, 5ML EXTRACTION TUBE, 25/RACK	\$	549.25
ULU-9ML-50R	ULTRALUTE, 9ML, DILUTION TUBE, 50/RACK	\$	780.52
ULUR-8ML-50R	ULTRALUTE/RESIN 8ML DILUTION TUBE 50/RACK	\$	1,138.25
VAA056	TOOL FOR SEAL REPLACEMENT ON	\$	124.30
VAA315	TOOL, WINDOW REMOVAL	\$	381.59
VAA356	nn Mounting stand	\$	287.15
VAA406	SPECIAL LENS REMOVAL TOOL	\$	211.99
VAA585	ADJUSTMENT SCREEN 2	\$	1,543.69
VAA586	ADJUSTMENT SCREEN 3	\$	1,543.69
VAA587	ADJUSTMENT SCREEN 4	\$	1,543.69
VAA591	aa Validation Kit	\$	3,835.13
VAA598	ADJUSTMENT SCREEN, ROUND CELL HOLDER	\$	1,551.40
VAA599	MIRROR APERTURE, LIGHTING MIRROR1	\$	1,543.69
VAA600	ADJUSTMENT SCREEN, REF. CELL	\$	1,688.23
VAA608	USB STICK WITH PROOF PROGRAM FOR SPECS	\$	296.79
VAA630	TESTFILTER KIT-STRAYLIGHT	\$	3,006.43
VAA631	TEST CUVETTE KIT-NEUTRAL GLASS	\$	2,389.73
VAA632	TESTFILTER KIT-WAVELENGTH ACCURACY	\$	2,119.92
VAA633	TEST SOLUTION KIT FOR SPECTROPHOTOMETER	\$	1,185.23
VAA634	COD-CUVETTE	\$	695.72
VAA643	aa Abstandslehre	\$	526.13
VAA654	Montagehilfe f"r SONATAX sc	\$	1,913.71
VAA752	db Druckkontrolleinheit	\$	948.18
VAA761-001	BLADE FOR TORQUE SCREW DRIVER	\$	271.74
VAA761-010	TORQUE SCREW DRIVER /	\$	680.30
VAA840	aa CALIBRATION SCREEN, 50mm	\$	2,042.83
VAA880	aa aa CAMERA TEST CELL	\$	724.63
VAA896	aa SERVICE RFID TAGS	\$	70.54
VAA897	aa aa ADAPTER FOR VAA880	\$	443.26
VAA901	PROOF FILTER UV RANGE, DR 5000/6000	\$	2,119.92
VAA903	MOUNTING AID CELL COMPARTMENT DR 6000	\$	1,085.01
VAA910	Montagevorrichtung /Justagehilfe	\$	1,015.63

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

VAA959	Test Cuvette for TU5 Series / Side LED vial	\$	680.30
VAA961	aa StablCal? Cali Set - for Hach Service Only	\$	2,486.09
VAA962	Spacing gauge Diffusor	\$	331.48
VAA964	Hall Sensor Tester	\$	418.20
VAA966	Secondary Turb. Stand. < 0.1 NTU	\$	599.84
VAA967	Secondary Turbidity Standard >20 NTU	\$	599.84
VAA976	Flow Sensor Test Set	\$	680.30
VAM084	Insertion Tool for TU5 Diffuser motor	\$	431.69
VP752401	Strain Relief: Liquid Tight; .181 to .312 wire dia	\$	26.02
W_3838_007_X	Compression spring 10 x 4.0 x 0.67 71/2.	\$	37.39
W_4107_002_PN	FTG, TBG, 8MM, HALF-CARTRIDGE, NITRILE O-RING	\$	15.03
W_4107_007_PN	ADAPTER, REDUCER 10/6 MM	\$	21.20
W_4108_131_99	FTG, BULKHEAD, TBG 6MM,QUICK CONNECT	\$	15.61
W_4509_002_X	CHECK VALVE, CRYSTAL SODIUM	\$	21.01
W_4510_001_FF	STIR BAR, MAGNETIC, DIAMETER 4.5MM LENGTH 15MM	\$	33.15
W_5109_001_00	BUTTON MAGNET, MIX MOTOR, DIAMETER 12.5MM	\$	10.60
WIMS_CONF_ST	Standard Claros Summit annual conference registration	\$	2,543.90
WM-BOD_MGR	KTO: WIMS BOD MGR LICENSE, DVD	\$	3,353.33
WM-BOD-YSI	WM YSI DO Meter Int	\$	11,659.56
WM-DBORACLE	WIMS for Oracle DB	\$	31,528.99
WM-DBSQLSERVER	WIMS for MS SQL Server DB	\$	21,873.72
WM-FACILITY	WIMS Additional Facilities	\$	4,875.82
WM-LABCAL-ADDL	WM LabCal for addl facility	\$	4,875.82
WM-LBCALC	WIMS LOAD BALANCING TOOL	\$	46,696.06
WM-LIMS-INT	WIMS STD LIMS IF-list on web	\$	12,218.45
WM-MOBILE-INT	WIMS STD MOBILE IF	\$	10,869.41
WM-MU	KTO: WM Multi Computer License	\$	16,997.90
WM-MU-CONUSR	WIMS CONCURRENT USER LICENSE	\$	3,758.04
WM-MU-LABCAL	WIMS LABCAL SOFTWARE MODULE	\$	18,385.49
WM-MU-LC	KTO: SOFTWARE, WIMS STANDALONE LAB CAL	\$	23,878.01
WM-OL-VIEW	KTO: HACH WIMS ONLINE VIEWER	\$	2,466.82
WM-SCADA-INT	WIMS STD SCADA/HIST IF- list on web	\$	12,218.45
WS1CW-180ML	dd Wash Solution 1 Concentrate	\$	229.46
WS2CW-200ML	dd Wash Solution 2 Concentrate	\$	129.36
X16T020	PERISTALTIC PUMP ASSEMBLY, SACXXX	\$	5,646.70
X31T053	Suction Tubing,35 cm	\$	217.77
X31T054	Omega clips, D=30mm 10pcs	\$	112.55
X31T071	oo KF CELL BAYONET HEAD ADAP F/HOMOGNZR	\$	202.36
X31T075	TUBING ASSEMBLY, MPA250 TURB RESERVR	\$	439.40
X31T087	WATERLAB BEAKERS W/COVER 125ML 50PCS	\$	427.84
X31T118	BEAKER HOLDER FOR TITRATION VESSEL, 0.5-5ML	\$	113.51
X31T128	BOTTLE TOP BURETTE HOLDER, B6XX	\$	2,871.53
X31T161	BEAKER COVER, D=48 MM, 42 PCS	\$	1,528.27
X41T003	Stirrer Rod, 100 mm, w/propellers	\$	1,568.74
X61T002	PERISTALTIC PUMP HEAD EXTENSION	\$	3,449.69
X61T003	PERISTALTIC PUMP HEAD	\$	3,449.69
X61T005	PHOTOCOLORIMETRIC TITRATION MODULE	\$	4,355.47

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

X91T091	KIT 25-POSITIONS, SAC9xx, 250ML	\$	8,479.68
XLH057	POWER CORD WITH UK PLUG	\$	174.60
XLH926	USB PC INTERFACE CABLE	\$	481.80
XMJ425	FUSE PLUG, (5x20) UL 248-14	\$	52.03
XMJ901	FUSE PLUG, (6,3x32) UL 248-14, 4.0A	\$	60.90
XMJ902	FUSE PLUG, (6,3x32) UL 248-14, 7.0A	\$	60.90
XQB091	Display 7" TFT LCD	\$	1,562.96
XQF024	LED LAMP, RED (FILTRAX)	\$	93.28
XQF025	LED LAMP, GREEN (FILTRAX)	\$	93.28
YAA249	PHOTOMETER CARD FOR PHOSPHAX INTER/2	\$	4,490.38
YAA368	TERMINAL BRD F/OPTIQUANT CONTROLLER	\$	1,297.01
YAA373	LIGHT BARRIER CARD	\$	1,985.02
YAA416	BOARD W/MESS ELEMENT, ULTRA TURB+	\$	300.64
YAA417	BOARD WITH REF. ELEMENT, UT plus	\$	300.64
YAA554	TERMINAL/CIRCUIT BOARD CARD (FILTRX)	\$	1,428.06
YAA555	INTERFACE BOARD, FILTRAX	\$	393.15
YAA857	MODBUS CARD FOR OPITQUANT CONTROLLER	\$	2,100.65
YAA981	FLASH SUPPLY CARD	\$	2,832.98
YAA982	FLASH SUPPLY CARD 2	\$	3,237.70
YAA985	CKT BD ASSY, PROCESSOR BOARD, FILTRAX	\$	15,301.97
YAB002	Netzteilplatte, doppelt/power pack b	\$	3,025.70
YAB014	Powerplatte SOLITAX sc	\$	3,661.68
YAB017	DISPLAY	\$	1,137.05
YAB018	CARD, INTERNAL SC1000, 4 mA INPUTS	\$	3,237.70
YAB019	CARD, 4mA OUTPUTS, SC1000 INTERNAL	\$	2,447.54
YAB020	CARD, INTERNAL PROFIBUS DP, SC1000	\$	1,667.03
YAB021	CARD, INTRNAL MODBUS (RS485), SC1000	\$	2,524.63
YAB023	Netzanschlusskarte 100-240VAC	\$	732.34
YAB024	ee Klemmenkarte SC1000	\$	1,927.20
YAB025	Leuchtdiodenkarte SC1000	\$	300.64
YAB031	Processor boards NITRATAX plus sc/UVAS plus sc	\$	4,625.28
YAB032	PCB, MAIN PROCESSOR SC1000 DSPLY MD	\$	6,860.83
YAB036	ANALOG CARD	\$	5,376.89
YAB037	WIPER CARD	\$	379.66
YAB042	CARD FILTRATION PROBE	\$	2,235.55
YAB043	ff PROCESSOR CARD FOR SC-ANALYZER	\$	3,815.86
YAB044	AMPLIFIER CARD, AMTAX sc	\$	2,312.64
YAB045	AMPLIFIER CARD FOR PHOSPHAX SC	\$	2,062.10
YAB048	Ersatzteilesatz Lampenplatte	\$	1,534.05
YAB049	WIPER MOTOR SOLITAX sc w/MOTOR PCB	\$	2,563.18
YAB051	MICROCOMPUTER PCB	\$	5,492.52
YAB052	POWER SUPPLY PCB	\$	2,254.82
YAB054	MOTHER BOARD DR3800	\$	7,824.43
YAB055	Data phone GSM (EU) Quad-Band	\$	7,805.16
YAB056	Data phone GSM (US) Quad-Band	\$	8,286.96
YAB058	DISPLAY INVERTER BOARD, DR3800	\$	568.52
YAB059	HALL SENSOR PCB	\$	294.86

OFICINA: 55 7575 1661 - 55 7575 1662

55 7160 2454 - 55 7158 7031

<http://laborwasserdemexico.com/>

LABOR WASSER
DE MÉXICO S.A. DE C.V.



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

YAB076	RELAY CARD, 4 POS. TYPE C, SC1000	\$	1,339.40
YAB077	MICROCOMPUTER BOARD	\$	5,704.51
YAB078	BATTERY POWER SUPPLY BOARD	\$	2,254.82
YAB081	LAMP SOCKET BOARD	\$	194.65
YAB082	PLUG-IN POWER SUPPLY BOARD	\$	291.01
YAB083	REFERENCE BOARD	\$	1,426.13
YAB084	DISPLAY, DR2800	\$	3,102.79
YAB089	CARD W/TEMP SENSOR AMTAX sc/PHOSPHAX	\$	377.73
YAB090	MAIN CIRCUIT BOARD, UT plus	\$	2,235.55
YAB091	MAIN BOARD W/MOTOR, SONATAX	\$	2,158.46
YAB092	LIGHT BARRIER BOARD. SONATAX	\$	277.52
YAB096	MMC ADAPTER WITH SD CARD	\$	227.41
YAB097	LED BOARD, DR2800	\$	371.95
YAB098	sc1000 DISPLAY MODULE LED ADAPTOR	\$	329.55
YAB099	PROCESSOR BOARD SC ANALYZER	\$	4,297.66
YAB100	KEYPAD, DRB200	\$	106.77
YAB102	NETWORK MODULE, PROFIBUS, sc200	\$	1,923.35
YAB103	MODULE, PROFIBUS DP, SC1000	\$	1,840.48
YAB105	MODULE, PROFIBUS DP, SC1000 WITH DOC	\$	1,840.48
YAB110	MAIN BOARD, 2100Q, EPA	\$	1,586.09
YAB111	MAIN BOARD, 2100Q, ISO	\$	1,927.20
YAB113	Prozessorkarte mit Adapter AN-ISE	\$	3,603.86
YAB117	inn RTC SYSTEMS COMMUNICATION CARD	\$	3,931.49
YAB119	ee MAIN BOARD W/PROCESSOR, DR3900	\$	5,627.42
YAB120	RFID MODULE	\$	807.50
YAB122	PRINTED CIRCUIT BOARD SIP 10	\$	529.98
YAB123	nn PROGNOISYS COMMUNICATION CARD, FOR SC1000	\$	3,931.49
YAB124	Mainboard DR 6000	\$	5,588.88
YAB125	HALL SENSOR BOARD - SPARE	\$	601.29
YAB126	Reference Element Board (PCBA)	\$	1,046.47
YAB127	Measurement Element Board	\$	1,377.95
YAB128	LAMP CONNECTOR BOARD - SPARE	\$	2,832.98
YAB129	Grating Motor with PCB	\$	1,040.69
YAB145	Switched Mode Power Supply Unit	\$	1,015.63
YAB146	CIRCUIT BOARD ASSY W/ DISPLAY, LOC 100 LOCATOR	\$	809.42
YAB147	PCBA, RFID SENSOR, LOC 100 & TU5x00SC	\$	271.74
YAB152	ANALOG BOARD UVAS PLUS sc, 1, 2MM	\$	4,625.28
YAB154	ANALOG BOARD UVAS PLUS sc, 50MM	\$	4,625.28
YAB155	ANALOG BOARD NITRATAX PLUS SC, 1, 2MM	\$	4,625.28
YAB156	ANALOG BOARD NITRATAX PLUS sc, 5MM	\$	4,625.28
YAB171	Display AP Comp1, Inter, TOctax, CODMax1 & P Sigma	\$	1,169.81
YAB172	LCD board for A/P Inter, P Sigma, TOctax & CODMax 1.	\$	1,379.88
YAB173	Motherboard + USB cover TitraLab AT1000	\$	2,813.71
YAB174	Controlboard + USB cover TitraLab AT1000	\$	1,509.00
YAB175	sc board w/o RFID online TU5400 sc board ohne RFID Pozess TU5400	\$	429.77
YAB189	Probe connector board, assembly with accessory	\$	312.21



PRECIOS LISTA HACH 2024

LABOR WASSER
DE MÉXICO S.A. DE C.V.

YAB190	Carrier board incl. processor module	\$	647.54
YAB191	Power Supply Card assembled with accessory	\$	2,062.10
YAB192	Impedance Adapterboard for AMTAX and PHOSPHAX on sc1500	\$	204.28
YAB200	See SC4200c power board and carrier board assembly	\$	1,755.68
YAB203	See SC4200c RF board US AT&T; assembly	\$	3,044.98
YAB205	See SC4200c Main board and VF61 COM board and Antenna board; assembly	\$	4,374.74
YAB212	Powersupply board 24V DC SC1000	\$	780.52
YAB217	SC4500 AC with carrier board,ribbon cable,screws	\$	795.93
YAB218	SC4500 DC with carrier board,ribbon cable,screws	\$	882.66
YAB219	SC4500 Main and COM board incl. screw	\$	1,325.91
YAB220	Sensor Flexboard NT3100sc	\$	131.82
YAB221	Sensor Flexboard NT3200sc	\$	131.82
YAB223	Flash Module, NT3x00sc	\$	736.19
YAB224	Mainboardgroup NT 3200sc	\$	1,086.94
YAB225	Mainboardgroup NT 3100sc	\$	1,086.94
ZDA417	THERMAL LINK	\$	33.34